

PPAR信号通路研究-产品清单

高级验证	产品名	货号	反应种属	验证应用
	ACADL Rabbit pAb	ER1901-11	Human, Mouse	WB,IF-Cell,IHC-P,FC
KD	ACADM Rabbit pAb	ER1804-01	Human, Mouse	WB,IF-Cell,IHC-P,FC
KD	ACADM Recombinant Rabbit mAb [JE48-09]	ET7109-74	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	ACBP Rabbit pAb	ER65484	Human, Mouse, Rat	WB
	ACOX2 Rabbit pAb	ER60439	Human, Mouse, Rat	WB
	ACOX3 Rabbit pAb	ER60440	Human	WB
	ACSL1 Rabbit pAb	ER60807	Human, Rat, Mouse	WB
KD	ACSL1 Mouse mAb [A9B5]	HA601112	Human, Mouse, Rat	WB,IHC-P
	ACSL3 Recombinant Mouse mAb [A8F6-R] - BSA and Azide free	HA610091	Human	WB,IHC-P
	ACSL3 Recombinant Mouse mAb [A8F6-R]	HA601213	Human	WB,IHC-P
	ACSL3 Mouse mAb [A8F6]	HA601069	Human	WB
	ACSL3 Rabbit pAb	ER60808	Human, Rat, Mouse	WB
	ACSL3 Mouse mAb [A8F7]	HA601070	Human	WB,IHC-P
KD	ACSL4 / FACL4 Recombinant Rabbit mAb [JE56-33] - BSA and Azide free	HA750469	Human, Mouse, Rat	WB,IHC-P
KD	ACSL4 / FACL4 Recombinant Rabbit mAb [JE56-33]	ET7111-43	Human, Mouse, Rat	WB,IHC-P
	ACSL5 Recombinant Mouse mAb [A9D4-R] - BSA and Azide free	HA610050	Mouse, Rat	WB,IHC-P
	ACSL5 Rabbit pAb	ER60809	Human, Rat, Mouse	WB
	ACSL5 Recombinant Mouse mAb [A9D3-R]	HA601166	Human, Mouse, Rat	WB,IHC-P
	ACSL5 Recombinant Mouse mAb [A9D4-R]	HA601149	Mouse, Rat	WB,IHC-P
	ACSL5 Recombinant Mouse mAb [A9D3-R] - BSA and Azide free	HA610053	Human, Mouse, Rat	WB,IHC-P
	ACSL6 Rabbit pAb	ER62582	Human, Mouse, Rat	WB,IHC-P
	ACSVL4 Rabbit pAb	ER62584	Human	IF-Cell
KD	ADFP / PLIN2 Recombinant Rabbit mAb [JA31-81] - BSA and Azide free	HA750407	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,FC,IHC-P
KD	ADFP / PLIN2 Recombinant Rabbit mAb [JA31-81]	ET1704-17	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,FC,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Adiponectin Recombinant Rabbit mAb [PSH05-48]	HA722286	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Adiponectin Recombinant Rabbit mAb [PSH05-48] - BSA and Azide free	HA751007	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Adiponectin Recombinant Rabbit mAb [JE65-49]	HA721064	Mouse	WB,IF-Cell,IHC-P
	Angiopoietin-like 4 Rabbit pAb	ER1903-68	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
KD	ANGPTL4 Recombinant Rabbit mAb [PSH07-35] - BSA and Azide free	HA751156	Human, Monkey	WB
KD	ANGPTL4 Recombinant Rabbit mAb [PSH07-35]	HA722823	Human, Monkey	WB
	ApoA-V Rabbit pAb	ER62685	Human	WB
	APOC3 Rabbit pAb	ER1802-47	Human	WB,IHC-P
	Apolipoprotein A1 Rabbit pAb	RT1046	Human, Mouse, Rat	WB,IP,IF
	Apolipoprotein A1 Recombinant Rabbit mAb [JF0548]	ET1702-23	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Apolipoprotein A1 Rabbit pAb	ER1802-29	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Apolipoprotein A1 Mouse mAb	M0809-13	Human	WB,IF-Cell,IHC-P
	Apolipoprotein A II Recombinant Rabbit mAb [JB70-35]	ET7107-47	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
Relative expression	BLBP Recombinant Rabbit mAb [PSH0-88]	HA721523	Human	WB
Relative expression	BLBP Recombinant Rabbit mAb [PSH0-88] - BSA and Azide free	HA750668	Human	WB
Relative expression	BLBP Recombinant Rabbit mAb [PSH0-89]	HA721524	Human	WB,IHC-P
Relative expression	BLBP Recombinant Rabbit mAb [PSH0-89] - BSA and Azide free	HA750669	Human	WB,IHC-P
Cell treatment	CD36 Recombinant Rabbit mAb [JJ2005] - BSA and Azide free	HA750307	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	CD36 Rabbit pAb	ER1802-7	Human, Mouse, Rat	IF-Cell,IHC-P,FC
Relative expression	CD36 Recombinant Rabbit mAb [PSH18-47]	HA724000	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Cell treatment	CD36 Recombinant Rabbit mAb [JJ2005]	ET1701-24	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	CD36 Recombinant Rabbit mAb	IRS118RB	Human	mIHC
	CPT1A Recombinant Rabbit mAb [PSH10-94]	HA723261	Mouse, Rat	WB
	CPT1A Recombinant Rabbit mAb [PSH04-00]	HA722081	Human, Mouse	WB,IHC-P,IP
	CPT1A Recombinant Rabbit mAb [PSH04-00] - BSA and Azide free	HA750923	Human, Mouse	WB,IHC-P,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CPT1A Recombinant Rabbit mAb [PSH10-94] - BSA and Azide free	HA751384	Mouse, Rat	WB
KD	CPT2 Recombinant Rabbit mAb [SN06-70]	ET1611-64	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	CPT2 Recombinant Rabbit mAb [SN06-70] - BSA and Azide free	HA750265	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	CPTI-M Rabbit pAb	ER63116	Human, Mouse, Rat	WB,IHC-P
	CYP27A1 Recombinant Rabbit mAb [JG40-36]	ET7109-05	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	CYP4A11 Recombinant Rabbit mAb [JE40-81]	HA722348	Human, Mouse, Rat	WB,IHC-P
	CYP7A1 Rabbit pAb	HA500207	Human, Rat	WB,IF-Cell,FC
	CYP8B1 Rabbit pAb	ER1907-09	Human	WB
	FABP3 Recombinant Rabbit mAb [JU73-03]	ET1706-50	Mouse, Rat	WB,IHC-P
	FABP3 Recombinant Rabbit mAb [JE53-37]	HA721859	Human, Mouse, Rat	WB,IHC-P
Relative expression	FABP4 Recombinant Rabbit mAb [JM10-99]	ET1703-98	Human, Mouse, Rat	WB,IF-Cell,IHC-Fr
Relative expression	FABP4 Recombinant Rabbit mAb [JM10-99] - BSA and Azide free	HA750398	Human, Mouse, Rat	WB,IF-Cell,IHC-Fr
	FABP5 Rabbit pAb	ER62026	Human, Rat, Mouse	WB
	FABP6 Rabbit pAb	ER62027	Human, Rat, Mouse	WB
	FATP2 Mouse mAb [A9-B3]	EM1708-57	Human	WB,IHC-P
	GK2 Rabbit pAb	ER63540	Human, Mouse	WB,IHC-P,IF-Cell
KD	Glycerol Kinase Recombinant Rabbit mAb [JE55-12]	ET7110-96	Human, Mouse, Rat	WB,IHC-P,FC
	Biotin Conjugated H-FABP Recombinant Rabbit mAb [PSH16-99] - Detector	HA725178B	Human, Mouse	ELISA(Det),ELISA
	H-FABP Recombinant Rabbit mAb [PSH16-98] - BSA and Azide free (Capture)	HA725177	Human, Mouse	ELISA(Cap)
	H-FABP Recombinant Rabbit mAb [PSH16-99] - BSA and Azide free (Detector)	HA725178	Human, Mouse	ELISA(Det)
Cell treatment	HMGCS1 Recombinant Rabbit mAb [PSH08-52] - BSA and Azide free	HA751227	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IP
	HMGCS1 Recombinant Rabbit mAb [PSH08-51] - BSA and Azide free	HA751226	Human, Mouse, Rat, Monkey	WB,IP
	HMGCS1 Recombinant Rabbit mAb [PSH08-51]	HA722997	Human, Mouse, Rat, Monkey	WB,IP
Cell treatment	HMGCS1 Recombinant Rabbit mAb [PSH08-52]	HA722998	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	HMGCS2 Recombinant Rabbit mAb [JE54-76]	HA720083	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Human Adiponectin Recombinant Rabbit mAb [PSH05-31] - BSA and Azide free (Capture)	HA722260	Human	ELISA(Cap)
	Human Adiponectin Recombinant Rabbit mAb [PSH05-32] - BSA and Azide free (Detector)	HA722261	Human	ELISA(Det)
	Biotin Conjugated Human Adiponectin Recombinant Rabbit mAb [PSH05-32] - Detector	HA722595B	Human	ELISA(Det),ELISA
	Human ANGPTL4 Recombinant Rabbit mAb [PSH20-15] - BSA and Azide free (Detector)	HA725390	Human	ELISA(Det)
	Biotin Conjugated Human ANGPTL4 Recombinant Rabbit mAb [PSH20-15] - Detector	HA725390B	Human	ELISA(Det),ELISA
	Human ANGPTL4 Recombinant Rabbit mAb [PSH20-14] - BSA and Azide free (Capture)	HA725389	Human	ELISA(Cap)
	Biotin Conjugated Human Apolipoprotein A1 Recombinant Rabbit mAb [PSH18-90] - Detector	HA724027B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human Apolipoprotein A1 Recombinant Rabbit mAb [PSH18-91] - Detector	HA724028B	Human	ELISA(Det),ELISA
	Human Apolipoprotein A1 Recombinant Rabbit mAb [PSH18-91] - BSA and Azide free (Detector)	HA724028	Human	ELISA(Det)
	Human Apolipoprotein A1 Recombinant Rabbit mAb [PSH18-90] - BSA and Azide free (Detector)	HA724027	Human	ELISA(Det)
	Human Apolipoprotein A1 Recombinant Rabbit mAb [PSH18-89] - BSA and Azide free (Capture)	HA724026	Human	ELISA(Cap)
	Biotin Conjugated Human MMP-1 Recombinant Rabbit mAb [PSH07-60]	HA722876B	Human	ELISA(Det)
	Human MMP-1 Recombinant Rabbit mAb [PSH07-60] - BSA and Azide free (Detector)	HA722875	Human	ELISA(Det)
	Human MMP-1 Recombinant Rabbit mAb [PSH07-59] - BSA and Azide free (Capture)	HA722874	Human	ELISA(Cap)
	ILK Rabbit pAb	RT1330	Human, Mouse, Rat	WB,IP,IF
KD	Integrin linked ILK Recombinant Rabbit mAb [SC68-04]	ET1610-76	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC
	Integrin linked ILK Rabbit pAb	ER1803-78	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	intestinal FABP Rabbit pAb	ER1911-53	Human, Mouse	WB,IHC-P,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	K27-linkage Specific Ubiquitin Recombinant Rabbit mAb [JE75-30]	HA723066	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,Dot Blot
	K63-linkage Specific Ubiquitin Recombinant Rabbit mAb [JM09-67]	ET1703-07	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	L-FABP Mouse mAb [AH54-31]	EM170403	Human, Mouse	ELISA(Cap),WB,ICC,IHC-P
	L-FABP Mouse mAb [AH59-31]	EM170404	Human, Mouse	ELISA(Det),WB,ICC,IHC-P
	Lipoprotein lipase Recombinant Rabbit mAb [JA22-02]	ET1704-48	Human	WB,IHC-P
	liver FABP Recombinant Rabbit mAb [JE10-77]	HA722321	Human, Mouse	WB
	LOX 1 Rabbit pAb	ER1912-42	Human, Mouse, Rat, Rabbit	WB,IHC-P,IF-Tissue
	LOX-1 Rabbit pAb	ER1706-53	Human, Mouse	WB,IHC-P,IF-Cell,FC
	LXR alpha Recombinant Rabbit mAb [JA20-38]	ET1704-51	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	LXR alpha Recombinant Rabbit mAb [JA20-38] - BSA and Azide free	HA750414	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Malic enzyme 1 Rabbit pAb	ER63901	Human, Mouse, Rat	WB
	ME3 Rabbit pAb	ER63936	Human, Mouse	WB,IHC-P
Relative expression	MMP-1 Recombinant Rabbit mAb [PSH09-56] - BSA and Azide free	HA751303	Human	WB,IHC-P
Relative expression	MMP-1 Recombinant Rabbit mAb [PSH09-56]	HA723118	Human	WB,IHC-P
	MMP-1 Rabbit pAb	ER31211	Human, Mouse, Monkey	WB,IHC-P
	Mouse Adiponectin Recombinant Rabbit mAb [PSH14-84] - BSA and Azide free (Capture)	HA723682	Mouse	ELISA(Cap)
	Mouse Adiponectin Recombinant Rabbit mAb [PSH14-82] - BSA and Azide free (Capture)	HA723679	Mouse	ELISA(Cap)
	Mouse Adiponectin Recombinant Rabbit mAb [PSH14-83] - BSA and Azide free (Detector)	HA723680	Mouse	ELISA(Det)
	Biotin Conjugated Mouse Adiponectin Recombinant Rabbit mAb [PSH14-83] - Detector	HA723681B	Mouse	ELISA(Det),ELISA
	PBFE Rabbit pAb	ER64529	Human, Rat	WB,IHC-P
	PCK1 Recombinant Rabbit mAb [JU84-39]	ET7106-81	Human	IF-Cell,IHC-P,WB
KD	PCK2 Recombinant Rabbit mAb [JB52-39]	ET7107-29	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	PDPK1 Recombinant Rabbit mAb [SD204-9]	ET1612-27	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	Perilipin-5 Rabbit pAb	ER1902-53	Human, Mouse	WB,IHC-P,FC
	Perilipin A Recombinant Rabbit mAb [JM10-96]	ET1703-38	Human, Mouse	WB,IF-Cell
	PPAR Rabbit pAb	HA500569	Mouse, Rat	WB
	PPAR alpha Mouse mAb	EM1707-71	Human, Monkey	FC
	PPAR delta Rabbit pAb	ER1902-24	Human, Mouse	WB,FC
	PPAR Delta Mouse mAb [2F9]	EM60127	Human, Rat, Mouse	IF-Cell,IHC-P
	Retinoid X Receptor alpha Recombinant Rabbit mAb [JG99-38]	ET7108-99	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
	RPS27A Rabbit pAb	HA500371	Human, Rat, Mouse	WB,IHC-P,IF-Cell
	RXR gamma Rabbit pAb	ER64896	Human, Mouse	WB,IHC-P,IF-Cell
	SCD1 Rabbit pAb	HA500119	Human, Mouse, Rat	WB
	SCD1 Recombinant Rabbit mAb [PSH15-47] - BSA and Azide free	HA751562	Mouse, Rat	WB,IHC-Fr,IHC-P,FC
	SCD1 Rabbit pAb	ER1916-26	Human, Mouse, Rat	IHC-P,IF-Cell,IF-Tissue
	SCD1 Recombinant Mouse mAb [A9A11-R]	HA601180	Human, Mouse, Rat	WB
	SCD1 Recombinant Rabbit mAb [PSH15-47]	HA723753	Mouse, Rat	WB,IHC-Fr,IHC-P,FC
	SCD1 Recombinant Mouse mAb [A9A11-R] - BSA and Azide free	HA610063	Human, Mouse, Rat	WB
KD	SLC27A4 Recombinant Rabbit mAb [JE58-29]	HA720001	Human, Mouse, Rat	WB,IHC-P,FC
	SLC27A5 Mouse mAb [C10-B11]	EM1708-95	Human, Mouse	WB,IHC-P
	Sterol carrier protein 2 Recombinant Rabbit mAb [JE62-79]	HA720099	Human, Mouse	WB,IHC-P
	Ub Rabbit pAb	ER65617	Human, Mouse, Rat	IF-Cell,WB,IHC-P
	Ubiquitin Mouse mAb [6G6]	EM60128	Human, Rat, Mouse	WB,IHC-P
Cell treatment	Ubiquitin Recombinant Rabbit mAb [ST46-03]	ET1609-21	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
Cell treatment	Ubiquitin Rabbit pAb	ER31212	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Ubiquitin Recombinant Rabbit mAb [ST46-03] - BSA and Azide free	HA750165	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	UCP-1 Rabbit pAb	ER1918-37	Mouse	WB,FC
	UCP1 Rabbit pAb	HA500118	Human, Mouse, Rat	WB,IHC-P
Relative expression	UCP1 Recombinant Rabbit mAb [PSH15-73]	HA723776	Mouse, Rat	WB,IHC-Fr,IHC-P,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	UCP1 Recombinant Rabbit mAb [PSH15-73] - BSA and Azide free	HA751580	Mouse, Rat	WB,IHC-Fr,IHC-P,IP
	UCP1 Rabbit pAb	ER62425	Human, Mouse, Rat	WB