

猴种属反应性产品清单

高级验证	产品名	货号	反应种属	验证应用
KD	14-3-3 epsilon Recombinant Rabbit mAb [JJ08-40]	ET1701-1	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IF-Tissue,IHC-P
	4E-BP1 Rabbit pAb	ER62532	Mouse, Rat, Monkey	WB,IHC-P
	AAK1 Rabbit pAb	ER62539	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell
KD	Acetyl CoA synthetase Recombinant Rabbit mAb [JF0917]	ET1702-21	Human, Mouse, Rat, Green monkey	WB,IF-Cell
	Acid phosphatase/ACP1 Recombinant Rabbit mAb [PSH04-57] - BSA and Azide free	HA750964	Human, Green monkey	WB,IHC-P
	Acid phosphatase/ACP1 Recombinant Rabbit mAb [PSH04-57]	HA722156	Human, Green monkey	WB,IHC-P
Cell treatment	Acid sphingomyelinase Recombinant Rabbit mAb [PSH15-61] - BSA and Azide free	HA751572	Human, Mouse, Rat, Green monkey	WB,IHC-P,IP
Cell treatment	Acid sphingomyelinase Recombinant Rabbit mAb [PSH15-61]	HA723768	Human, Mouse, Rat, Green monkey	WB,IHC-P,IP
	Actin Recombinant Rabbit mAb [JF47-01]	ET1702-52	Human, Mouse, Rat, Green monkey, Zebrafish	WB,IHC-P,IF-Cell,FC
	Biotin Conjugated Activin E Recombinant Rabbit mAb [PSH22-84] - Detector	HA724329B	Human, Mouse, Rat, Cynomolgus monkey	ELISA(Det),ELISA
	Activin E Recombinant Rabbit mAb [PSH22-83] - BSA and Azide free (Capture)	HA724328	Human, Mouse, Rat, Cynomolgus monkey	ELISA(Cap)
	Activin E Recombinant Rabbit mAb [PSH22-84] - BSA and Azide free (Detector)	HA724329	Human, Mouse, Rat, Cynomolgus monkey	ELISA(Det)
	ADRM1 Recombinant Rabbit mAb [PSH05-41]	HA722270	Human, Mouse, Rat, Green monkey	WB,IF-Cell,FC,IP
	ADRM1 Recombinant Rabbit mAb [PSH05-41] - BSA and Azide free	HA751000	Human, Mouse, Rat, Green monkey	WB,IF-Cell,FC,IP
	AGO1 Recombinant Rabbit mAb [PSH19-77] - BSA and Azide free	HA751732	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P
	AGO1 Recombinant Rabbit mAb [PSH19-77]	HA724101	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P
	AKT1 Recombinant Rabbit mAb [PSH0-50]	HA721364	Human, Mouse, Rat, Green monkey	WB,IF-Cell,FC

	AKT1 Recombinant Rabbit mAb [PSH0-50] - BSA and Azide free	HA750627	Human, Mouse, Rat, Green monkey	WB,IF-Cell,FC
	AKT1/2/3 Recombinant Rabbit mAb [ST48-09]	ET1609-51	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
	AKT1/2/3 Recombinant Rabbit mAb [ST48-09] - BSA and Azide free	HA750179	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
KD	AKT1/2/3 Recombinant Rabbit mAb [JE75-09]	HA721870	Human, Mouse, Rat, Green monkey	WB,IF-Cell,FC,IP
KD	AKT1/2/3 Recombinant Rabbit mAb [JE75-09] - BSA and Azide free	HA750840	Human, Mouse, Rat, Green monkey	WB,IF-Cell,FC,IP
	ALDH1B1 Rabbit pAb	ER62640	Human, Monkey	WB
	alpha Actinin Rabbit pAb	ER1803-60	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC
	AMID Rabbit pAb	ER62655	Human, Mouse, Monkey	WB,IHC-P,IF-Cell
KD	ANGPTL4 Recombinant Rabbit mAb [PSH07-35] - BSA and Azide free	HA751156	Human, Green monkey	WB
KD	ANGPTL4 Recombinant Rabbit mAb [PSH07-35]	HA722823	Human, Green monkey	WB
	ARHGDI1A Mouse mAb [G8-C3]	EM1710-11	Human, Mouse, Monkey	WB,IF-Cell,FC
	ARHGEF10 Rabbit pAb	ER62720	Human, Monkey	WB,IHC-P
	ASH2L Recombinant Rabbit mAb [ST47-01]	ET1609-24	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IF-Tissue,IHC-P
	Asparagine synthetase Recombinant Rabbit mAb [SC67-07]	ET1610-50	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P
	ATG4B Recombinant Rabbit mAb [JE59-97]	HA721334	Human, Mouse, Rat, Green monkey	WB
KD	ATG4B Recombinant Rabbit mAb [PSH02-53]	HA721829	Human, Mouse, Rat, Green monkey	WB
KD	ATG4B Recombinant Rabbit mAb [PSH02-53] - BSA and Azide free	HA750816	Human, Mouse, Rat, Green monkey	WB
	ATG4B Rabbit pAb	ER62759	Human, Mouse, Rat, Green monkey	WB,IHC-P
	ATG4C Recombinant Rabbit mAb [PSH05-23] - BSA and Azide free	HA750994	Human, Mouse, Rat, Green monkey	WB,IP
	ATG4C Recombinant Rabbit mAb [PSH05-23]	HA722252	Human, Mouse, Rat, Green monkey	WB,IP
KO	ATG5 Recombinant Rabbit mAb [SN73-07] - BSA and Azide free	HA750251	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,FC(Intra)
KO	ATG5 Recombinant Rabbit mAb [SN73-07]	ET1611-38	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,FC(Intra)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	ATL3 Recombinant Rabbit mAb [PSH02-47] - BSA and Azide free	HA750812	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC,IP
	ATL3 Recombinant Rabbit mAb [PSH02-47]	HA721824	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC,IP
	ATP13A2 Recombinant Rabbit mAb [PSH06-69] - BSA and Azide free	HA751098	Human, Mouse, Rat, Green monkey, Cynomolgus monkey	WB,IHC-P,IHC-Fr
	ATP13A2 Recombinant Rabbit mAb [PSH06-69]	HA722741	Human, Mouse, Rat, Green monkey, Cynomolgus monkey	WB,IHC-P,IHC-Fr
	ATP6V1A Recombinant Rabbit mAb [PSH21-14] - BSA and Azide free	HA751786	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-Fr,IHC-P,FC,IP
	ATP6V1A Recombinant Rabbit mAb [PSH21-14]	HA724200	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-Fr,IHC-P,FC,IP
	ATPase Inhibitory Factor 1 / IF1 Recombinant Rabbit mAb [PSH16-94] - BSA and Azide free	HA751608	Human, Green monkey	WB,IF-Cell,IHC-P,FC,IP
	ATPase Inhibitory Factor 1 / IF1 Recombinant Rabbit mAb [PSH16-94]	HA723890	Human, Green monkey	WB,IF-Cell,IHC-P,FC,IP
	ATXN1 Mouse mAb [4C7B11]	EM1711-76	Human, Mouse, Rat, Monkey	WB,IHC-P,FC
Relative expression	Axl Recombinant Rabbit mAb [PSH06-68] - BSA and Azide free	HA751097	Human, Green monkey	WB,IF-Cell,FC
Relative expression	Axl Recombinant Rabbit mAb [PSH06-68]	HA722740	Human, Green monkey	WB,IF-Cell,FC
KD	BAF180 Recombinant Rabbit mAb [PSH02-34]	HA721811	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Tissue
KD	BAF180 Recombinant Rabbit mAb [PSH02-34] - BSA and Azide free	HA750799	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Tissue
	BANF1 Recombinant Rabbit mAb [JG40-43]	ET7109-08	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC
	BCKDHA Recombinant Rabbit mAb [PSH15-62]	HA723769	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC
	BCKDHA Recombinant Rabbit mAb [PSH15-62] - BSA and Azide free	HA751573	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC
	Bcl2-L-13 Recombinant Rabbit mAb [PSH17-12]	HA723899	Human, Mouse, Rat, Green monkey	WB,IHC-P,IP
	Bcl2-L-13 Recombinant Rabbit mAb [PSH17-12] - BSA and Azide free	HA751616	Human, Mouse, Rat, Green monkey	WB,IHC-P,IP
	BDH1 Recombinant Rabbit mAb [PSH14-69]	HA723667	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P
	BDH1 Recombinant Rabbit mAb [PSH14-69] - BSA and Azide free	HA751532	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	HRP Conjugated beta Actin Recombinant Mouse mAb [A2-F6-R] - Loading control	HA601585H	Human, Mouse, Rat, Green monkey, Pig	WB
	HRP Conjugated beta Actin Recombinant Rabbit mAb [PSH03- 63] - Loading control	HA724235H	Human, Mouse, Rat, Green monkey, Pig	WB
	beta COP Recombinant Rabbit mAb [PSH19-01] - BSA and Azide free	HA751691	Human, Mouse, Rat, Green monkey	WB,IF-Cell,FC
	beta COP Recombinant Rabbit mAb [PSH19-01]	HA724037	Human, Mouse, Rat, Green monkey	WB,IF-Cell,FC
	beta I Tubulin Recombinant Rabbit mAb [JE01-23]	HA722290	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC- P,FC,IF-Tissue
	beta Tubulin Recombinant Mouse mAb [A1-A4-R]	HA601187	Human, Mouse, Rat, Green monkey	WB,IF-Cell,FC,IHC-P
	beta Tubulin Mouse mAb [A1-A4]	M1305-2	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC
	beta Tubulin Recombinant Mouse mAb [A1-A4-R] - BSA and Azide free	HA610067	Human, Mouse, Rat, Green monkey	WB,IF-Cell,FC,IHC-P
	HRP Conjugated beta Tubulin Recombinant Rabbit mAb [SR25-04] - Loading control	HA724236H	Human, Mouse, Rat, Green monkey	WB
	BMPR2 Mouse mAb [F8-G6]	EM1709-92	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P
	Bmx Rabbit pAb	ER62821	Human, Mouse, Monkey	WB,IHC-P,IF-Cell
	BNIP3L Recombinant Rabbit mAb [JE59-98]	HA722056	Human, Mouse, Rat, Green monkey	WB,IHC-P
KD	BRAF Mouse mAb [F1-F4-E10]	EM1701-32	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,FC
	BRCC36 Recombinant Rabbit mAb [PSH05-10]	HA722226	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC
	BRCC36 Recombinant Rabbit mAb [PSH05-10] - BSA and Azide free	HA750983	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC
	BRD9 Recombinant Rabbit mAb [PSH15-95] - BSA and Azide free	HA751586	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,IP
	BRD9 Recombinant Rabbit mAb [PSH15-95]	HA723800	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,IP
	Bub3 Recombinant Rabbit mAb [JG38-81]	ET7108-82	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC- P,FC,IF-Tissue
	C3orf38 Recombinant Rabbit mAb [JE64-28]	HA722360	Human, Green monkey	WB,IF-Cell
	CAB39 Recombinant Rabbit mAb [PSH15-16] - BSA and Azide free	HA751545	Human, Mouse, Rat, Green monkey	WB,IHC-P,IP
	CAB39 Rabbit pAb	ER60546	Human, Mouse, Rat, Green monkey	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CAB39 Recombinant Rabbit mAb [PSH15-16]	HA723720	Human, Mouse, Rat, Green monkey	WB,IHC-P,IP
	Calpain 1 Recombinant Rabbit mAb [PSH24-37] - BSA and Azide free	HA751938	Human, Mouse, Rat, Green monkey	WB
	Calpain 1 Recombinant Rabbit mAb [PSH24-37]	HA724444	Human, Mouse, Rat, Green monkey	WB
	Calreticulin Recombinant Antibody [SU37-03] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610363	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell
	Calreticulin Recombinant Antibody [SU37-03] - Mouse IgG1 (Chimeric)	HA601547	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell
	CARM1 Mouse mAb [B10-H2]	EM1709-69	Human, Monkey, Rat, Mouse	WB,IF-Cell,IHC-P,FC
	Casein Kinase I gamma 1 Rabbit pAb	ER62902	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell
KD	Cathepsin L Recombinant Rabbit mAb [PSH03-91]	HA722063	Human, Mouse, Rat, Green monkey	WB,IHC-P,FC
KD	Cathepsin L Recombinant Rabbit mAb [PSH03-91] - BSA and Azide free	HA750916	Human, Mouse, Rat, Green monkey	WB,IHC-P,FC
KD	Cathepsin L Recombinant Rabbit mAb [PSH03-92] - BSA and Azide free	HA750917	Human, Mouse, Rat, Green monkey	WB,IHC-P
KD	Cathepsin L Recombinant Rabbit mAb [PSH03-92]	HA722064	Human, Mouse, Rat, Green monkey	WB,IHC-P
	CBL Recombinant Rabbit mAb [PSH17-61] - BSA and Azide free	HA751635	Human, Mouse, Rat, Green monkey	WB,IF-Cell,FC,IP
	CBL Recombinant Rabbit mAb [PSH17-61]	HA723936	Human, Mouse, Rat, Green monkey	WB,IF-Cell,FC,IP
	Cbx8 Recombinant Rabbit mAb [PSH06-61]	HA722734	Human, Mouse, Rat, Green monkey	WB,IF-Cell,FC,IP
	Cbx8 Recombinant Rabbit mAb [PSH06-61] - BSA and Azide free	HA751091	Human, Mouse, Rat, Green monkey	WB,IF-Cell,FC,IP
Relative expression	CCK4 Mouse mAb [A11A9]	HA601344	Human, Mouse, Rat, Green monkey	WB,IF-Cell
	CCL2 Mouse mAb [E11-D8]	EM1710-22	Human, Mouse, Monkey	WB,IF-Cell,IHC-P,FC
	CD2BP2 Rabbit pAb	ER62942	Human, Monkey	WB,IHC-P
	CD354 Mouse mAb [6B12D1]	EM1711-44	Human, Rat, Monkey	WB,FC
	Cdk6 Rabbit pAb	ER40101	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P
	CELF1 Recombinant Rabbit mAb [JE40-85]	HA722552	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,FC,IP
	CELF1 Mouse mAb [7A1]	EM1701-77	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CELF1 Recombinant Mouse mAb [7A1-R]	HA601327	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P
	CELF1 Recombinant Mouse mAb [7A1-R] - BSA and Azide free	HA610199	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P
	CEP70 Rabbit pAb	ER62986	Human, Monkey	WB,IHC-P
	Cerebellin 2 Rabbit pAb	ER62991	Human, Mouse, Rat, Monkey	WB
KD	CHD1 Mouse mAb [2-2]	M1211-5	Human, Mouse, Green monkey	WB,IHC-P
	Clathrin heavy chain Recombinant Rabbit mAb [JA53-02]	ET1704-50	Human, Mouse, Rat, Green monkey	WB,IHC-P
	CLC-6 Rabbit pAb	ER63029	Human, Mouse, Rat, Monkey	WB
	Cleaved-Cathepsin L2 (L114) Rabbit pAb	ER60017	Human, Monkey	WB
Cell treatment	Cleaved PARP1 Recombinant Rabbit mAb [JE00-15]	HA722218	Human, Mouse, Rat, Green monkey	WB
	CMTM3 Rabbit pAb	ER63052	Human, Monkey	WB
	CMTM6 Recombinant Rabbit mAb [PSH15-02] - BSA and Azide free	HA751540	Human, Green monkey	WB,IP
	CMTM6 Recombinant Rabbit mAb [PSH15-02]	HA723702	Human, Green monkey	WB,IP
	c-Myc Mouse mAb [6G6]	RT1149	Human, Mouse, Rat, Monkey	WB,IP,IF,IHC-P
	c-Myc Recombinant Rabbit mAb [JE01-20]	HA722895	Human, Mouse, Rat, Green monkey	WB,IF-Cell
	Cofilin Recombinant Rabbit mAb [PSH0-33]	HA721326	Human, Mouse, Rat, Green monkey	WB,IF-Cell,FC
	Cofilin Rabbit pAb	HA500523	Human, Mouse, Rat, Green monkey	WB,IF-Cell,FC
	Cofilin Recombinant Rabbit mAb [PSH0-38]	HA721340	Human, Mouse, Rat, Green monkey	WB,FC
	Cofilin Recombinant Rabbit mAb [PSH0-38] - BSA and Azide free	HA750616	Human, Mouse, Rat, Green monkey	WB,FC
	Cofilin Recombinant Rabbit mAb [PSH01-44]	HA721682	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,FC
	Cofilin Recombinant Rabbit mAb [PSH0-33] - BSA and Azide free	HA750612	Human, Mouse, Rat, Green monkey	WB,IF-Cell,FC
	Cofilin Recombinant Rabbit mAb [PSH01-44] - BSA and Azide free	HA750721	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,FC
	Cofilin Recombinant Rabbit mAb [PSH0-34]	HA721327	Human, Mouse, Rat, Green monkey	WB,FC
	Cofilin Recombinant Rabbit mAb [PSH0-34] - BSA and Azide free	HA750613	Human, Mouse, Rat, Green monkey	WB,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Cofilin Rabbit pAb	ER63062	Human, Mouse, Rat, Green monkey	WB,IHC-P
	COL6A2 Rabbit pAb	ER63079	Human, Mouse, Monkey	WB,IHC-P,IF-Cell
	COP zeta1 Rabbit pAb	ER63093	Human, Mouse, Monkey	WB,IHC-P
	CREB Recombinant Rabbit mAb [PSH15-48]	HA723754	Human, Mouse, Rat, Green monkey	WB,IHC-P,IHC-Fr,IF-Cell,FC
	CREB Recombinant Rabbit mAb [PSH15-48] - BSA and Azide free	HA751563	Human, Mouse, Rat, Green monkey	WB,IHC-P,IHC-Fr,IF-Cell,FC
	CREB Recombinant Rabbit mAb [SA04-04] - BSA and Azide free	HA750007	Human, Mouse, Green monkey, Zebrafish, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	CREB Recombinant Rabbit mAb [SA04-04]	ET1601-15	Human, Mouse, Green monkey, Zebrafish, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	CRTC2 Mouse mAb [D2-B10]	EM1709-78	Human, Monkey	WB,IF-Cell,IHC-P,FC
	CRTC3 Recombinant Rabbit mAb [JE63-18]	HA721887	Human, Mouse, Rat, Green monkey	WB,IHC-P,FC
	CRTC3 Mouse mAb [G8-A9]	EM1709-65	Human, Monkey	WB,IF-Cell,IHC-P
	CSNK2A1 Recombinant Rabbit mAb [PSH04-54] - BSA and Azide free	HA750961	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	CSNK2A1 Recombinant Rabbit mAb [PSH04-54]	HA722153	Human, Mouse, Rat, Cynomolgus monkey	WB,IF-Cell,IHC-P,FC
	CT A Recombinant Rabbit mAb [JE48-23]	HA722459	Human, Green monkey	WB,IF-Cell,FC
	CTU2 Recombinant Rabbit mAb [PSH04-05] - BSA and Azide free	HA750928	Human, Mouse, Green monkey	WB,IHC-P
	CTU2 Recombinant Rabbit mAb [PSH04-05]	HA722086	Human, Mouse, Green monkey	WB,IHC-P
	Cx40 / GJA5 Recombinant Rabbit mAb [PSH10-82]	HA723249	Human, Mouse, Rat, Green monkey	WB,IP
	Cx40 / GJA5 Recombinant Rabbit mAb [PSH10-82] - BSA and Azide free	HA751374	Human, Mouse, Rat, Green monkey	WB,IP
	Cyclin A Rabbit pAb	ER63153	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell
	Cyclophilin 40 Recombinant Rabbit mAb [PSH02-71] - BSA and Azide free	HA750828	Human, Green monkey	WB,IHC-P,FC,IF-Cell
	Cyclophilin 40 Recombinant Rabbit mAb [PSH02-71]	HA721849	Human, Green monkey	WB,IHC-P,FC,IF-Cell
	Cyclophilin 40 Rabbit pAb	ER65466	Human, Mouse, Rat, Green monkey	WB,IHC-P
KO	CYLD Recombinant Rabbit mAb [PSH11-21]	HA723303	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KO	CYLD Recombinant Rabbit mAb [PSH11-22] - BSA and Azide free	HA751631	Human, Mouse, Green monkey	WB,IF-Cell
KO	CYLD Recombinant Rabbit mAb [PSH11-21] - BSA and Azide free	HA751630	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell
KO	CYLD Recombinant Rabbit mAb [PSH11-22]	HA723304	Human, Mouse, Green monkey	WB,IF-Cell
	Cynomolgus IFN-gamma, Tag Free Protein	HA211458	Monkey	
	Cynomolgus PCSK9, C-His Tag Protein	HA211416	Monkey	
	CYP1A1/2 Rabbit pAb	ER63159	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell
	CYP2R1 Rabbit pAb	ER63172	Human, Mouse, Monkey	WB,IHC-P
	CysLTR1 Rabbit pAb	ER63188	Human, Monkey	WB,IHC-P,IF-Cell
	Cytokeratin 5 Mouse mAb [6D-4D]	EM1706-36	Human, Mouse, Monkey	WB,IHC-P,IF-Cell,FC
	DDX3X Mouse mAb	EM1707-67	Human, Mouse, Monkey, Rat	WB,IF-Cell,IHC-P,FC
	DHODH Recombinant Rabbit mAb [PSH09-69]	HA723134	Human, Mouse, Rat, Green monkey	WB,IHC-P
	DHODH Recombinant Rabbit mAb [PSH09-69] - BSA and Azide free	HA751311	Human, Mouse, Rat, Green monkey	WB,IHC-P
	Dishevelled 3 Recombinant Rabbit mAb [PSH15-93] - BSA and Azide free	HA751584	Human, Mouse, Rat, Green monkey	WB,IP
	Dishevelled 3 Recombinant Rabbit mAb [PSH15-93]	HA723798	Human, Mouse, Rat, Green monkey	WB,IP
	DLST Recombinant Rabbit mAb [PSH10-07] - BSA and Azide free	HA751337	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P
	DLST Recombinant Rabbit mAb [PSH10-07]	HA723176	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P
Relative expression	DNAJC5 Recombinant Rabbit mAb [PSH24-38] - BSA and Azide free	HA751939	Human, Mouse, Rat, Green monkey, Pig	WB,IF-Cell,IHC-Fr,IHC-P,FC(Intra),IP
Relative expression	DNAJC5 Rabbit pAb	ER63138	Human, Mouse, Rat, Green monkey, Pig	WB,IF-Cell,FC(Intra)
Relative expression	DNAJC5 Recombinant Rabbit mAb [PSH24-38]	HA724445	Human, Mouse, Rat, Green monkey, Pig	WB,IF-Cell,IHC-Fr,IHC-P,FC(Intra),IP
	DNAM-1 Rabbit pAb	ER63273	Human, Mouse, Monkey	WB,IHC-P,IF-Cell
	DNA Polymerase lambda Recombinant Rabbit mAb [PSH05-24] - BSA and Azide free	HA750995	Human, Mouse, Green monkey	WB,IHC-P
	DNA Polymerase lambda Recombinant Rabbit mAb [PSH05-24]	HA722253	Human, Mouse, Green monkey	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	Dnmt2 Recombinant Rabbit mAb [PSH03-43]	HA721988	Human, Green monkey	WB,IF-Cell,IHC-P,FC
KD	Dnmt2 Recombinant Rabbit mAb [PSH03-43] - BSA and Azide free	HA750877	Human, Green monkey	WB,IF-Cell,IHC-P,FC
Relative expression	Drebrin Recombinant Rabbit mAb [JE30-20]	HA722062	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC
	E2F-4/5 Rabbit pAb	ER63306	Human, Mouse, Rat, Monkey	WB,IHC-P
	EDC4 Recombinant Rabbit mAb [PSH22-01]	HA724270	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,FC
	EDC4 Recombinant Rabbit mAb [PSH22-01] - BSA and Azide free	HA751817	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,FC
KD	EEA1 Recombinant Rabbit mAb [JE59-34]	HA722147	Human, Mouse, Rat, Green monkey	WB,IF-Cell,FC
Relative expression, KD	EGFR Mouse mAb [A8E12]	HA601073	Human, Green monkey	WB,IHC-P
	EGFR Rabbit pAb	R1511-18	Human, Green monkey	WB,IF-Cell,IHC-P
	EGFR Rabbit pAb	R1608-3	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC
Relative expression, KD	EGFR Recombinant Rabbit mAb [PD00-75] - BSA and Azide free	HA750550	Human, Green monkey	WB,IHC-P,IF-Cell,FC
Relative expression, KO	EGFR Mouse mAb [7-F2-F8]	EM1901-67	Human, Green monkey	WB,IHC-P,FC
Relative expression	EGFR Rabbit pAb	ER1512-6	Human, Mouse, Green monkey	WB,IF-Cell,IHC-P,FC
Relative expression, KO	EGFR Recombinant Rabbit mAb [SP00-86] - BSA and Azide free	HA750078	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IP,FC
	EGFR Rabbit pAb	0407-21	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P
Relative expression, KD	EGFR Recombinant Rabbit mAb [PD00-75]	HA721220	Human, Green monkey	WB,IHC-P,IF-Cell,FC
Relative expression, KD	EGFR Recombinant Rabbit mAb [SZ40-19]	ET1603-37	Human, Mouse, Green monkey, Rat	WB,IF-Cell,IP,FC,mIHC
Relative expression, KO	EGFR Recombinant Rabbit mAb [SP00-86]	ET1604-44	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IP,FC
	EGFR Mouse mAb [B6-E5-D9]	EM1701-22	Human, Green monkey	WB,IHC-P,IF-Cell,FC(Intra)
	EIF2S1 Recombinant Rabbit mAb [PSH04-29]	HA722112	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	EIF2S1 Recombinant Rabbit mAb [PSH04-29] - BSA and Azide free	HA750938	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC,IP
	EIF3S1 / EIF3J Recombinant Rabbit mAb [PSH06-97] - BSA and Azide free	HA751124	Human, Rat, Green monkey	WB,IF-Cell,IHC-P,FC
	EIF3S1 / EIF3J Recombinant Rabbit mAb [PSH06-97]	HA722775	Human, Rat, Green monkey	WB,IF-Cell,IHC-P,FC
	eIF4A2 Recombinant Rabbit mAb [PSH01-89]	HA721747	Human, Mouse, Rat, Green monkey	WB,IP
	eIF4A2 Recombinant Rabbit mAb [PSH01-88]	HA721746	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,FC
	eIF4A2 Recombinant Rabbit mAb [PSH01-88] - BSA and Azide free	HA750762	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,FC
	eIF4A2 Recombinant Rabbit mAb [PSH01-89] - BSA and Azide free	HA750763	Human, Mouse, Rat, Green monkey	WB,IP
	eIF4G1 Recombinant Rabbit mAb [PSH06-28]	HA722660	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,FC
	eIF4G1 Recombinant Rabbit mAb [PSH06-28] - BSA and Azide free	HA751065	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,FC
	Eph receptor B1 Recombinant Rabbit mAb [JE39-46]	HA722345	Human, Mouse, Rat, Green monkey	WB
	Ephrin-B1 Rabbit pAb	ER63382	Human, Mouse, Rat, Monkey	WB
KD	EPS15 Recombinant Rabbit mAb [PSH23-03] - BSA and Azide free	HA751885	Human, Mouse, Rat, Green monkey	WB,IF-Cell,FC
KD	EPS15 Recombinant Rabbit mAb [PSH23-03]	HA724343	Human, Mouse, Rat, Green monkey	WB,IF-Cell,FC
	ETAR Rabbit pAb	ER63399	Human, Mouse, Monkey	WB,IHC-P,IF-Cell
	FAF1 Recombinant Rabbit mAb [PSH22-51]	HA724301	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC,IP
	FAF1 Recombinant Rabbit mAb [PSH22-51] - BSA and Azide free	HA751864	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC,IP
	FAM80A Rabbit pAb	ER63409	Human, Mouse, Monkey	WB
	FAM80B Rabbit pAb	ER63410	Human, Mouse, Rat, Monkey	WB
	FASTKD1 Rabbit pAb	ER63414	Human, Mouse, Monkey	WB
KD	Filamin A Recombinant Rabbit mAb [SA30-08]	ET1601-3	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	Flotillin-2 Recombinant Rabbit mAb [PSH03-41] - BSA and Azide free	HA750875	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,FC
KD	Flotillin-2 Recombinant Rabbit mAb [PSH03-41]	HA721986	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,FC
	FoxL1 Rabbit pAb	ER63462	Human, Monkey	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Frizzled-9 Rabbit pAb	ER63472	Human, Mouse, Monkey	WB,IF-Cell
	FUK Mouse mAb [E2-C3]	EM1710-03	Human, Mouse, Rat, Monkey	WB,FC
	Fumarate hydratase Recombinant Rabbit mAb [PSH08-80]	HA723032	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,FC
	FUS/TLS Recombinant Rabbit mAb [JJ09-31]	ET1701-86	Human, Mouse, Green monkey, Cynomolgus monkey, Pig	WB,IHC-P,FC,IF-Cell
	FUS/TLS Recombinant Rabbit mAb [JJ09-31] - BSA and Azide free	HA750326	Human, Mouse, Green monkey, Cynomolgus monkey, Pig	WB,IHC-P,FC,IF-Cell
	FXR1 Recombinant Rabbit mAb [PSH01-67] - BSA and Azide free	HA750743	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,FC,IP
	FXR1 Recombinant Rabbit mAb [PSH01-67]	HA721710	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,FC,IP
	FXR1 Recombinant Rabbit mAb [PSH01-35] - BSA and Azide free	HA750713	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,IP
	FXR1 Rabbit pAb	ER1909-14	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,FC
	FXR1 Recombinant Rabbit mAb [PSH01-35]	HA721669	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,IP
	G3BP1 Rabbit pAb	ER63481	Human, Mouse, Monkey	WB,IHC-P
	gamma Tubulin Recombinant Rabbit mAb [JF0972]	ET1702-32	Human, Mouse, Rat, Green monkey, Chicken, Hamster, Cow, Dog, Fish, Xenopus tropicalis	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Gankyrin Recombinant Rabbit mAb [PSH17-80] - BSA and Azide free	HA751639	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC,IF-Tissue
	Gankyrin Recombinant Rabbit mAb [PSH17-80]	HA723949	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC,IF-Tissue
	GAPDH Rabbit pAb	R1210-1	Human, Mouse, Rat, Green monkey, Oryza sativa	WB,IF-Cell,IHC-P,FC
	HRP Conjugated GAPDH Recombinant Rabbit mAb [JF81-04] - Loading control	ET1702-66	Human, Mouse, Rat, Green monkey, Chicken, Fish	WB
	HRP Conjugated GAPDH Recombinant Rabbit mAb [SA30-01] - Loading control	HA724237H	Human, Mouse, Rat, Green monkey	WB
KD	Gelsolin Recombinant Rabbit mAb [JB36-68]	ET7107-45	Human, Mouse, Green monkey	WB,IF-Tissue,IHC-P
	Giantin Recombinant Rabbit mAb [PSH03-86]	HA722052	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC
	Giantin Recombinant Rabbit mAb [PSH03-86] - BSA and Azide free	HA750911	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	GPI Recombinant Rabbit mAb [PSH25-26] - BSA and Azide free	HA751972	Human, Mouse, Rat, Green monkey	WB,IHC-P,IP
	GPI Recombinant Rabbit mAb [PSH25-26]	HA724519	Human, Mouse, Rat, Green monkey	WB,IHC-P,IP
	GPR101 Rabbit pAb	ER63560	Human, Mouse, Monkey	WB,IF-Cell
	GPR133 Rabbit pAb	ER63569	Human, Monkey	WB,IF-Cell
	GPR34 Rabbit pAb	ER63601	Human, Mouse, Rat, Monkey	WB,IF-Cell
	GPR40 Rabbit pAb	ER63605	Human, Monkey	WB,IF-Cell
	GPRC5C Rabbit pAb	ER63620	Human, Mouse, Rat, Monkey	WB,IF-Cell
	GRB2 Recombinant Rabbit mAb [JE49-73]	HA721923	Human, Mouse, Rat, Green monkey	WB,IF-Cell,FC
KD	GRK2 Recombinant Rabbit mAb [ST05-60]	ET1609-72	Human, Mouse, Rat, Cynomolgus monkey	WB,IHC-P
	GRP94 Recombinant Rabbit mAb [PSH25-59] - BSA and Azide free	HA751978	Human, Mouse, Rat, Green monkey, Pig	WB,IF-Cell,IHC-P,FC(Intra),IP
	GRP94 Recombinant Rabbit mAb [PSH25-59]	HA724533	Human, Mouse, Rat, Green monkey, Pig	WB,IF-Cell,IHC-P,FC(Intra),IP
	GSC Recombinant Rabbit mAb [JE51-35]	HA722460	Human, Rat, Green monkey	WB
	GSK3 alpha Recombinant Rabbit mAb [JE33-27]	HA722217	Human, Mouse, Rat, Green monkey	WB,IHC-P
	GSTM1 Mouse mAb [A4-H4]	EM1709-33	Human, Monkey	WB,IF-Cell,FC
	GWL Recombinant Rabbit mAb [PSH07-69] - BSA and Azide free	HA751181	Human, Green monkey	WB,IF-Cell,IP
	GWL Recombinant Rabbit mAb [PSH07-69]	HA722884	Human, Green monkey	WB,IF-Cell,IP
	HADHA Recombinant Rabbit mAb [PSH01-29] - BSA and Azide free	HA750707	Human, Mouse, Rat, Green monkey	WB
	HADHA Recombinant Rabbit mAb [PSH01-28] - BSA and Azide free	HA750706	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell
	HADHA Recombinant Rabbit mAb [PSH01-28]	HA721652	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell
	HADHA Recombinant Rabbit mAb [PSH01-29]	HA721653	Human, Mouse, Rat, Green monkey	WB
	Hamartin Recombinant Rabbit mAb [PSH14-67] - BSA and Azide free	HA751528	Human, Mouse, Rat, Green monkey	WB,IP
	Hamartin Recombinant Rabbit mAb [PSH14-67]	HA723665	Human, Mouse, Rat, Green monkey	WB,IP
	HCCS Rabbit pAb	ER63666	Human, Mouse, Monkey	WB,IHC-P,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	HDAC3 Recombinant Rabbit mAb [SC55-02] - BSA and Azide free	HA750221	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP
KD	HDAC3 Recombinant Rabbit mAb [SC55-02]	ET1610-5	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP
KD	HDAC6 Recombinant Rabbit mAb [JJ09-09]	ET1701-66	Human, Green monkey	WB,IF-Tissue,IHC-P,IP
	Hexokinase II Recombinant Rabbit mAb [PSH07-97] - BSA and Azide free	HA751199	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,FC,IP
	Hexokinase II Recombinant Rabbit mAb [PSH07-97]	HA722933	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,FC,IP
	H-FABP Recombinant Rabbit mAb [PSH16-99] - BSA and Azide free (Detector)	HA725178	Human, Mouse, Cynomolgus monkey	ELISA(Det)
	Biotin Conjugated H-FABP Recombinant Rabbit mAb [PSH16-99] - Detector	HA725178B	Human, Mouse, Cynomolgus monkey	ELISA(Det),ELISA
	H-FABP Recombinant Rabbit mAb [PSH16-98] - BSA and Azide free (Capture)	HA725177	Human, Mouse, Cynomolgus monkey	ELISA(Cap)
	hHR23A Recombinant Rabbit mAb [JE38-49]	HA722215	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC
	HIP2 / LIG Recombinant Rabbit mAb [PSH11-24]	HA723306	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC
	HIP2 / LIG Recombinant Rabbit mAb [PSH11-24] - BSA and Azide free	HA751403	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC
	Histone H1 (Acetyl K25) Rabbit pAb	ER60044	Human, Monkey	WB,IHC-P
	Histone H2B (Acetyl K12) Rabbit pAb	ER60047	Human, Mouse, Monkey	WB,IHC-P,IF-Cell
	HRP Conjugated Histone H3 Recombinant Rabbit mAb [JJ090-07] - Loading control	HA724276H	Human, Mouse, Rat, Green monkey	WB
	Histone H3 (di methyl K4) Recombinant Rabbit mAb [JE00-98]	HA722135	Human, Mouse, Rat, Green monkey	WB,IF-Cell,ChIP,Dot Blot,IHC-P
	Histone H3 (mono+di+tri methyl K79) Recombinant Rabbit mAb [PSH09-10] - BSA and Azide free	HA751264	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,ChIP
	Histone H3 (mono+di+tri methyl K79) Recombinant Rabbit mAb [PSH09-10]	HA723068	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,ChIP
	Histone H3 (mono methyl K14) Recombinant Rabbit mAb [JE43-29]	HA722232	Human, Mouse, Rat, Green monkey	WB,IF-Cell,ChIP,Dot Blot
KD	Histone H3 (mono methyl K79) Recombinant Rabbit mAb [PSH07-66] - BSA and Azide free	HA751175	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC,ChIP
KD	Histone H3 (mono methyl K79) Recombinant Rabbit mAb [PSH07-66]	HA722848	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC,ChIP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Histone H3 (mono methyl K9) Recombinant Rabbit mAb [JE43-27]	HA722350	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC- P,FC,ChIP,IF- Tissue,Dot Blot
	Histone H3 (tri methyl K36) Recombinant Rabbit mAb [PSH15- 01] - BSA and Azide free	HA751539	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF- Cell,IP,Dot Blot,ChIP
	Histone H3 (tri methyl K36) Recombinant Rabbit mAb [PSH15- 01]	HA723701	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF- Cell,IP,Dot Blot,ChIP
	Histone H3 (tri methyl K4) Recombinant Rabbit mAb [PSH25- 68] - BSA and Azide free	HA751987	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC- P,FC(Intra),ChIP
	Histone H3 (tri methyl K4) Recombinant Rabbit mAb [PSH25- 68]	HA724555	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC- P,FC(Intra),ChIP
	Histone H3 (tri methyl K9) Recombinant Rabbit mAb [PSH19- 58]	HA724084	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,ChIP
	Histone H3 (tri methyl K9) Recombinant Rabbit mAb [PSH19- 58] - BSA and Azide free	HA751722	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,ChIP
	HMGCS1 Recombinant Rabbit mAb [PSH08-51]	HA722997	Human, Mouse, Rat, Green monkey	WB,IP
	HMGCS1 Recombinant Rabbit mAb [PSH08-51] - BSA and Azide free	HA751226	Human, Mouse, Rat, Green monkey	WB,IP
KD	HNRNPA0 Recombinant Rabbit mAb [PSH02-77] - BSA and Azide free	HA750832	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC
KD	HNRNPA0 Recombinant Rabbit mAb [PSH02-77]	HA721855	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC
	HRI Rabbit pAb	ER63714	Human, Mouse, Monkey	WB,IHC-P
	HSC70 Interacting Protein HIP Recombinant Rabbit mAb [JE33-17]	HA722274	Human, Mouse, Green monkey	WB
	Hsp40 Recombinant Rabbit mAb [PSH12-10]	HA723426	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC
	Hsp40 Recombinant Rabbit mAb [PSH12-10] - BSA and Azide free	HA751439	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC
	Hsp90 Recombinant Rabbit mAb [JE05-06]	HA722689	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC,IP
	Hsp90 alpha Recombinant Rabbit mAb [SY14-06] - BSA and Azide free	HA750091	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IF- Tissue,IHC-P
	Hsp90 alpha Recombinant Rabbit mAb [SY14-06]	ET1605-57	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IF- Tissue,IHC-P
	HSPA14 Recombinant Rabbit mAb [SD08-48]	ET1612-93	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IF- Tissue,IHC-P
	Human Adiponectin Recombinant Rabbit mAb [PSH05-31] - BSA and Azide free (Capture)	HA722260	Human, Cynomolgus monkey	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human Adiponectin Recombinant Rabbit mAb [PSH05-32] - BSA and Azide free (Detector)	HA722261	Human, Cynomolgus monkey	ELISA(Det)
	Biotin Conjugated Human Adiponectin Recombinant Rabbit mAb [PSH05-32] - Detector	HA722595B	Human, Cynomolgus monkey	ELISA(Det),ELISA
	HRP Conjugated Human Adiponectin Recombinant Rabbit mAb [PSH05-32] - Detector	HA722595H	Human, Cynomolgus monkey	ELISA(Det),ELISA
	Human alpha Fetoprotein Recombinant Rabbit mAb [PSH14-54] - BSA and Azide free (Detector)	HA723645	Human, Cynomolgus monkey	ELISA(Det)
	Human alpha Fetoprotein Recombinant Rabbit mAb [PSH14-55] - BSA and Azide free (Detector)	HA723647	Human, Cynomolgus monkey	ELISA(Det)
	Biotin Conjugated Human alpha Fetoprotein Recombinant Rabbit mAb [PSH14-54] - Detector	HA723646B	Human, Cynomolgus monkey	ELISA(Det)
	Biotin Conjugated Human alpha Fetoprotein Recombinant Rabbit mAb [PSH14-55] - Detector	HA723648B	Human, Cynomolgus monkey	ELISA(Det)
	Human alpha Fetoprotein Recombinant Rabbit mAb [PSH14-53] - BSA and Azide free (Capture)	HA723644	Human, Cynomolgus monkey	ELISA(Cap)
	Human CCL11/Eotaxin Recombinant Rabbit mAb [PSH16-24] - BSA and Azide free (Capture)	HA725197	Human, Cynomolgus monkey	ELISA(Cap)
	Biotin Conjugated Human CCL11/Eotaxin Recombinant Rabbit mAb [PSH16-25] - Detector	HA725198B	Human, Cynomolgus monkey	ELISA(Det),ELISA
	Human CCL11/Eotaxin Recombinant Rabbit mAb [PSH16-25] - BSA and Azide free (Detector)	HA725198	Human, Cynomolgus monkey	ELISA(Det)
	Human CCL2 / MCP-1 Recombinant Rabbit mAb [PSH18-21] - BSA and Azide free (Capture)	HA725195	Human, Cynomolgus monkey	ELISA(Cap)
	Biotin Conjugated Human CCL2 / MCP-1 Recombinant Rabbit mAb [PSH18-22] - Detector	HA725196B	Human, Cynomolgus monkey	ELISA(Det)
	Human CCL2 / MCP-1 Recombinant Rabbit mAb [PSH18-22] - BSA and Azide free (Detector)	HA725196	Human, Cynomolgus monkey	ELISA(Det)
	Human CCL4 / MIP-1 beta Recombinant Rabbit mAb [PSH19-48] - BSA and Azide free (Capture)	HA725355	Human, Cynomolgus monkey	ELISA(Cap)
	Human CCL4 / MIP-1 beta Recombinant Rabbit mAb [PSH19-49] - BSA and Azide free (Detector)	HA725356	Human, Cynomolgus monkey	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human CCL4 / MIP-1 beta Recombinant Rabbit mAb [PSH19-49] - Detector	HA725356B	Human, Cynomolgus monkey	ELISA(Det),ELISA
	Human CCL5 Recombinant Rabbit mAb [PSH13-16] - BSA and Azide free (Capture)	HA725055	Human, Dog, Cynomolgus monkey	ELISA(Cap)
	HRP Conjugated Human CCL5 Recombinant Rabbit mAb [PSH13-17] - Detector	HA725057H	Human, Dog, Cynomolgus monkey	ELISA(Det),ELISA
	Biotin Conjugated Human CCL5 Recombinant Rabbit mAb [PSH13-17] - Detector	HA725057B	Human, Dog, Cynomolgus monkey	ELISA(Det),ELISA
	Human CCL5 Recombinant Rabbit mAb [PSH13-17] - BSA and Azide free (Detector)	HA725056	Human, Dog, Cynomolgus monkey	ELISA(Det)
	Biotin Conjugated Human Clusterin Recombinant Rabbit mAb [PSH15-56] - Detector	HA725093B	Human, Cynomolgus monkey	ELISA(Det),ELISA
	Human Clusterin Recombinant Rabbit mAb [PSH15-56] - BSA and Azide free (Detector)	HA725092	Human, Cynomolgus monkey	ELISA(Det)
	HRP Conjugated Human Clusterin Recombinant Rabbit mAb [PSH15-56] - Detector	HA725093H	Human, Cynomolgus monkey	ELISA(Det),ELISA
	Human Clusterin Recombinant Rabbit mAb [PSH15-55] - BSA and Azide free (Capture)	HA725091	Human, Cynomolgus monkey	ELISA(Cap)
	HRP Conjugated Human Complement C3 Recombinant Rabbit mAb [PSH15-08] - Detector	HA725142H	Human, Cynomolgus monkey	ELISA(Det),ELISA
	Human Complement C3 Recombinant Rabbit mAb [PSH15-07] - BSA and Azide free (Capture)	HA725140	Human, Cynomolgus monkey	ELISA(Cap)
	Human Complement C3 Recombinant Rabbit mAb [PSH15-08] - BSA and Azide free (Detector)	HA725141	Human, Cynomolgus monkey	ELISA(Det)
	Biotin Conjugated Human Complement C3 Recombinant Rabbit mAb [PSH15-08] - Detector	HA725142B	Human, Cynomolgus monkey	ELISA(Det),ELISA
	Human CRP Recombinant Rabbit mAb [PSH04-25] - BSA and Azide free (Capture)	HA722108	Human, Dog, Cynomolgus monkey	ELISA(Cap)
	Biotin Conjugated Human CRP Recombinant Rabbit mAb [PSH04-26] - Detector	HA722616B	Human, Dog, Cynomolgus monkey	ELISA(Det)
	Human CRP Recombinant Rabbit mAb [PSH04-26] - BSA and Azide free (Detector)	HA722109	Human, Dog, Cynomolgus monkey	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human CXCL10 / IP-10 Recombinant Rabbit mAb [PSH12-36] - BSA and Azide free (Capture)	HA725007	Human, Cynomolgus monkey	ELISA(Cap)
	Human CXCL10 / IP-10 Recombinant Rabbit mAb [PSH12-37] - BSA and Azide free (Detector)	HA725008	Human, Cynomolgus monkey	ELISA(Det)
	Biotin Conjugated Human CXCL10 / IP-10 Recombinant Rabbit mAb [PSH12-37] - Detector	HA725009B	Human, Cynomolgus monkey	ELISA(Det),ELISA
	Human EGF Recombinant Rabbit mAb [PSH02-20] - BSA and Azide free (Detector)	HA721783	Human, Cynomolgus monkey	ELISA(Det)
	Biotin Conjugated Human EGF Recombinant Rabbit mAb [PSH02-20] - Detector	HA722638B	Human, Cynomolgus monkey	ELISA(Det)
	Human EGF Recombinant Rabbit mAb [PSH02-19] - BSA and Azide free (Capture)	HA721782	Human, Cynomolgus monkey	ELISA(Cap)
	Human GM-CSF Recombinant Rabbit mAb [PS01-38] - BSA and Azide free (Detector)	HA721289	Human, Dog, Cynomolgus monkey	ELISA(Det)
	Biotin Conjugated Human GM-CSF Recombinant Rabbit mAb [PS01-38] - Detector	HA722648B	Human, Dog, Cynomolgus monkey	ELISA(Det)
	Human GM-CSF Recombinant Rabbit mAb [PS01-52] - BSA and Azide free (Capture)	HA721288	Human, Dog, Cynomolgus monkey	ELISA(Cap)
	Human IL-10 Recombinant Rabbit mAb [PS00-81] - BSA and Azide free (Detector)	HA721237	Human, Dog, Cynomolgus monkey	ELISA(Det)
	Human IL-10 Recombinant Rabbit mAb [PS00-98] - BSA and Azide free (Capture)	HA721236	Human, Dog, Cynomolgus monkey	ELISA(Cap)
	Biotin Conjugated Human IL-10 Recombinant Rabbit mAb [PS00-81] - Detector	HA722652B	Human, Dog, Cynomolgus monkey	ELISA(Det)
	Biotin Conjugated Human IL-12B p40 Recombinant Rabbit mAb [PSH02-58] - Detector	HA722633B	Human, Dog, Cynomolgus monkey	ELISA(Det)
	Human IL-12B p40 Recombinant Rabbit mAb [PSH02-57] - BSA and Azide free (Capture)	HA721834	Human, Dog, Cynomolgus monkey	ELISA(Cap)
	Human IL-12B p40 Recombinant Rabbit mAb [PSH02-58] - BSA and Azide free (Detector)	HA721835	Human, Dog, Cynomolgus monkey	ELISA(Det)
	Human IL-17A Recombinant Rabbit mAb [PS01-24] - BSA and Azide free (Detector)	HA721233	Human, Dog, Cynomolgus monkey	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human IL-17A Recombinant Rabbit mAb [PS01-20] - BSA and Azide free (Capture)	HA721232	Human, Dog, Cynomolgus monkey	ELISA(Cap)
	HRP Conjugated Human IL-17A Recombinant Rabbit mAb [PS01-24] - Detector	HA721296	Human, Dog, Cynomolgus monkey	ELISA(Det)
	Biotin Conjugated Human IL-17A Recombinant Rabbit mAb [PS01-24] - Detector	HA721631	Human, Dog, Cynomolgus monkey	ELISA(Det)
	Human IL-18 Recombinant Rabbit mAb [PSH01-70] - BSA and Azide free (Detector)	HA721723	Human, Dog, Cynomolgus monkey	ELISA(Det)
	Human IL-18 Recombinant Rabbit mAb [PSH01-71] - BSA and Azide free (Capture)	HA721722	Human, Dog, Cynomolgus monkey	ELISA(Cap)
	Biotin Conjugated Human IL-18 Recombinant Rabbit mAb [PSH01-70] - Detector	HA722644B	Human, Dog, Cynomolgus monkey	ELISA(Det)
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-85] - BSA and Azide free (Capture)	HA723881	Human, Dog, Cynomolgus monkey	ELISA(Cap)
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-86] - BSA and Azide free (Detector)	HA723882	Human, Dog, Cynomolgus monkey	ELISA(Det)
	Biotin Conjugated Human IL-1 beta Recombinant Rabbit mAb [PSH16-86] - Detector	HA723882B	Human, Dog, Cynomolgus monkey	ELISA(Det),ELISA
	Biotin Conjugated Human IL1RL1 / ST2 Recombinant Rabbit mAb [PSH19-93] - Detector	HA725376B	Human, Cynomolgus monkey	ELISA(Det),ELISA
	Human IL1RL1 / ST2 Recombinant Rabbit mAb [PSH19-93] - BSA and Azide free (Detector)	HA725376	Human, Cynomolgus monkey	ELISA(Det)
	Human IL1RL1 / ST2 Recombinant Rabbit mAb [PSH19-92] - BSA and Azide free (Capture)	HA725375	Human, Cynomolgus monkey	ELISA(Cap)
	Human IL-4 Recombinant Rabbit mAb [PS00-67] - BSA and Azide free (Capture)	HA721230	Human, Cynomolgus monkey	ELISA(Cap)
	Biotin Conjugated Human IL-4 Recombinant Rabbit mAb [PS00-66] - Detector	HA721368	Human, Cynomolgus monkey	ELISA(Det)
	HRP Conjugated Human IL-4 Recombinant Rabbit mAb [PS00-66] - Detector	HA721295	Human, Cynomolgus monkey	ELISA(Det)
	Human IL-4 Recombinant Rabbit mAb [PS00-66] - BSA and Azide free (Detector)	HA721231	Human, Cynomolgus monkey	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human IL-6 Recombinant Rabbit mAb [PS00-21] - BSA and Azide free (Capture)	HA721167	Human, Dog, Cynomolgus monkey	ELISA(Cap)
	Human IL-6 Recombinant Rabbit mAb [PS00-15] - BSA and Azide free (Detector)	HA721168	Human, Dog, Cynomolgus monkey	ELISA(Det)
	HRP Conjugated Human IL-6 Recombinant Rabbit mAb [PS00-15] - Detector	HA721293	Human, Dog, Cynomolgus monkey	ELISA(Det)
	Biotin Conjugated Human IL-6 Recombinant Rabbit mAb [PS00-15] - Detector	HA722233B	Human, Dog, Cynomolgus monkey	ELISA(Det)
	Human IL-8 Recombinant Rabbit mAb [PS00-83] - BSA and Azide free (Capture)	HA722098	Human, Dog, Cynomolgus monkey	ELISA(Cap)
	Biotin Conjugated Human IL-8 Recombinant Rabbit mAb [PS01-22] - Detector	HA722611B	Human, Dog, Cynomolgus monkey	ELISA(Det)
	Human IL-8 Recombinant Rabbit mAb [PS01-22] - BSA and Azide free (Detector)	HA722099	Human, Dog, Cynomolgus monkey	ELISA(Det)
	Human Interferon gamma Recombinant Rabbit mAb [PSH01-45] - BSA and Azide free (Capture)	HA721683	Human, Dog, Cynomolgus monkey	ELISA(Cap)
	Human Interferon gamma Recombinant Rabbit mAb [PSH01-47] - BSA and Azide free (Detector)	HA721685	Human, Dog, Cynomolgus monkey	ELISA(Det)
	Biotin Conjugated Human Interferon gamma Recombinant Rabbit mAb [PSH01-47] - Detector	HA722645B	Human, Dog, Cynomolgus monkey	ELISA(Det),ELISA
	Biotin Conjugated Human Lipocalin-2 Recombinant Rabbit mAb [PSH13-68] - Detector	HA725111B	Human, Cynomolgus monkey	ELISA(Det),ELISA
	Human Lipocalin-2 Recombinant Rabbit mAb [PSH13-67] - BSA and Azide free (Capture)	HA725109	Human, Cynomolgus monkey	ELISA(Cap)
	Human Lipocalin-2 Recombinant Rabbit mAb [PSH13-68] - BSA and Azide free (Detector)	HA725110	Human, Cynomolgus monkey	ELISA(Det)
	HRP Conjugated Human Lipocalin-2 Recombinant Rabbit mAb [PSH13-68] - Detector	HA725111H	Human, Cynomolgus monkey	ELISA(Det),ELISA
	HRP Conjugated Human Lipoprotein A Recombinant Rabbit mAb [PSH13-83] - Detector	HA723574H	Human, Cynomolgus monkey	ELISA(Det),ELISA
	Human Lipoprotein A Recombinant Rabbit mAb [PSH13-82] - BSA and Azide free (Detector)	HA723571	Human, Cynomolgus monkey	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human Lipoprotein A Recombinant Rabbit mAb [PSH13-83] - Detector	HA723574B	Human, Cynomolgus monkey	ELISA(Det),ELISA
	Biotin Conjugated Human Lipoprotein A Recombinant Rabbit mAb [PSH13-82] - Detector	HA723572B	Human, Cynomolgus monkey	ELISA(Det),ELISA
	Human Lipoprotein A Recombinant Rabbit mAb [PSH13-81] - BSA and Azide free (Capture)	HA723570	Human, Cynomolgus monkey	ELISA(Cap)
	Human Lipoprotein A Recombinant Rabbit mAb [PSH13-83] - BSA and Azide free (Detector)	HA723573	Human, Cynomolgus monkey	ELISA(Det)
	Human/Monkey ACE Recombinant Rabbit mAb [PSH22-81] - BSA and Azide free (Capture)	HA724326	Human, Cynomolgus monkey	ELISA(Cap)
	Human/Monkey ACE Recombinant Rabbit mAb [PSH22-82] - BSA and Azide free (Detector)	HA724327	Human, Cynomolgus monkey	ELISA(Det)
	Biotin Conjugated Human/Monkey ACE Recombinant Rabbit mAb [PSH22-82] - Detector	HA724327B	Human, Cynomolgus monkey	ELISA(Det),ELISA
	Human/Monkey Albumin Recombinant Rabbit mAb [PSH20-17] - BSA and Azide free (Detector)	HA724128	Human, Dog, Cynomolgus monkey	ELISA(Det)
	Biotin Conjugated Human/Monkey Albumin Recombinant Rabbit mAb [PSH20-17] - Detector	HA724128B	Human, Dog, Cynomolgus monkey	ELISA(Det),ELISA
	HRP Conjugated Human/Monkey Albumin Recombinant Rabbit mAb [PSH20-17] - Detector	HA724128H	Human, Dog, Cynomolgus monkey	ELISA(Det),ELISA
	Human/Monkey Albumin Recombinant Rabbit mAb [PSH20-16] - BSA and Azide free (Capture)	HA724127	Human, Dog, Cynomolgus monkey	ELISA(Cap)
	Human/Monkey Alpha-1-Microglobulin Recombinant Rabbit mAb [PSH25-81] - BSA and Azide free (Capture)	HA724573	Human, Cynomolgus monkey	ELISA(Cap)
	Human/Monkey Alpha-1-Microglobulin Recombinant Rabbit mAb [PSH25-82] - BSA and Azide free (Detector)	HA724574	Human, Cynomolgus monkey	ELISA(Det)
	Biotin Conjugated Human/Monkey Alpha-1-Microglobulin Recombinant Rabbit mAb [PSH25-82] - Detector	HA724574B	Human, Cynomolgus monkey	ELISA(Det),ELISA
	HRP Conjugated Human/Monkey Apolipoprotein CIII Recombinant Rabbit mAb [PSH24-27] - Detector	HA724434H	Human, Rhesus monkey	ELISA(Det),ELISA
	Human/Monkey Apolipoprotein CIII Recombinant Rabbit mAb [PSH24-26] - BSA and Azide free (Capture)	HA724433	Human, Rhesus monkey	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human/Monkey Apolipoprotein CIII Recombinant Rabbit mAb [PSH24-27] - Detector	HA724434B	Human, Rhesus monkey	ELISA(Det),ELISA
	Human/Monkey Apolipoprotein CIII Recombinant Rabbit mAb [PSH24-27] - BSA and Azide free (Detector)	HA724434	Human, Rhesus monkey	ELISA(Det)
	Human/Monkey Apo-M Recombinant Rabbit mAb [PSH24-30] - BSA and Azide free (Capture)	HA724437	Human, Rhesus monkey	ELISA(Cap)
	Human/Monkey Apo-M Recombinant Rabbit mAb [PSH24-29] - BSA and Azide free (Detector)	HA724436	Human, Rhesus monkey	ELISA(Det)
	Biotin Conjugated Human/Monkey Apo-M Recombinant Rabbit mAb [PSH24-29] - Detector	HA724436B	Human, Rhesus monkey	ELISA(Det),ELISA
	Human/Monkey Apo-M Recombinant Rabbit mAb [PSH24-28] - BSA and Azide free (Capture)	HA724435	Human, Rhesus monkey	ELISA(Cap)
	Human/Monkey Coagulation factor VIII Recombinant Rabbit mAb [PSH21-70] - BSA and Azide free (Detector)	HA724241	Human, Cynomolgus monkey	ELISA(Det)
	Biotin Conjugated Human/Monkey Coagulation factor VIII Recombinant Rabbit mAb [PSH21-70] - Detector	HA724241B	Human, Cynomolgus monkey	ELISA(Det),ELISA
	Human/Monkey Coagulation factor VIII Recombinant Rabbit mAb [PSH21-69] - BSA and Azide free (Capture)	HA724240	Human, Cynomolgus monkey	ELISA(Cap)
	Biotin Conjugated Human/Monkey ESM1 Recombinant Rabbit mAb [PSH24-32] - Detector	HA724439B	Human, Rhesus monkey	ELISA(Det),ELISA
	Human/Monkey ESM1 Recombinant Rabbit mAb [PSH24-32] - BSA and Azide free (Detector)	HA724439	Human, Rhesus monkey	ELISA(Det)
	Human/Monkey ESM1 Recombinant Rabbit mAb [PSH24-31] - BSA and Azide free (Capture)	HA724438	Human, Rhesus monkey	ELISA(Cap)
	Human/Monkey IGFBP2 Recombinant Rabbit mAb [PSH25-34] - BSA and Azide free (Detector)	HA724525	Human, Cynomolgus monkey	ELISA(Det)
	Biotin Conjugated Human/Monkey IGFBP2 Recombinant Rabbit mAb [PSH25-34] - Detector	HA724525B	Human, Cynomolgus monkey	ELISA(Det),ELISA
	Human/Monkey IGFBP2 Recombinant Rabbit mAb [PSH25-33] - BSA and Azide free (Capture)	HA724524	Human, Cynomolgus monkey	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human/Monkey Mannose Receptor(CD206) Recombinant Rabbit mAb [PSH23-47] - BSA and Azide free (Detector)	HA724378	Human, Cynomolgus monkey	ELISA(Det)
	Biotin Conjugated Human/Monkey Mannose Receptor(CD206) Recombinant Rabbit mAb [PSH23-47] - Detector	HA724378B	Human, Cynomolgus monkey	ELISA(Det),ELISA
	Human/Monkey Mannose Receptor(CD206) Recombinant Rabbit mAb [PSH23-46] - BSA and Azide free (Capture)	HA724377	Human, Cynomolgus monkey	ELISA(Cap)
	Human/Monkey Pro-COL1A1 Recombinant Rabbit mAb [PSH25-15] - BSA and Azide free (Capture)	HA724509	Human, Rhesus monkey	ELISA(Cap)
	Human/Monkey Pro-COL1A1 Recombinant Rabbit mAb [PSH25-16] - BSA and Azide free (Detector)	HA724510	Human, Rhesus monkey	ELISA(Det)
	Biotin Conjugated Human/Monkey Pro-COL1A1 Recombinant Rabbit mAb [PSH25-16] - Detector	HA724510B	Human, Rhesus monkey	ELISA(Det),ELISA
	Human/Monkey TIMP4 Recombinant Rabbit mAb [PSH21-49] - BSA and Azide free (Capture)	HA724220	Human, Rhesus monkey	ELISA(Cap)
	Biotin Conjugated Human/Monkey TIMP4 Recombinant Rabbit mAb [PSH21-50] - Detector	HA724221B	Human, Rhesus monkey	ELISA(Det),ELISA
	Human/Monkey TIMP4 Recombinant Rabbit mAb [PSH21-50] - BSA and Azide free (Detector)	HA724221	Human, Rhesus monkey	ELISA(Det)
	Biotin Conjugated Human/Monkey VEGFR3 / FLT4 Recombinant Rabbit mAb [PSH22-17] - Detector	HA724278B	Human, Cynomolgus monkey	ELISA(Det),ELISA
	Human/Monkey VEGFR3 / FLT4 Recombinant Rabbit mAb [PSH22-17] - BSA and Azide free (Detector)	HA724278	Human, Cynomolgus monkey	ELISA(Det)
	Human/Monkey VEGFR3 / FLT4 Recombinant Rabbit mAb [PSH22-16] - BSA and Azide free (Capture)	HA724277	Human, Cynomolgus monkey	ELISA(Cap)
	HRP Conjugated Human/Monkey VEGFR3 / FLT4 Recombinant Rabbit mAb [PSH22-17] - Detector	HA724278H	Human, Cynomolgus monkey	ELISA(Det),ELISA
	Biotin Conjugated Human/Mouse MASP2 Recombinant Rabbit mAb [PSH20-52] - Detector	HA725392B	Human, Mouse, Rhesus monkey	ELISA(Det),ELISA
	Human/Mouse MASP2 Recombinant Rabbit mAb [PSH20-52] - BSA and Azide free (Detector)	HA725392	Human, Mouse, Rhesus monkey	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human/Mouse MASP2 Recombinant Rabbit mAb [PSH20-51] - BSA and Azide free (Capture)	HA725391	Human, Mouse, Rhesus monkey	ELISA(Cap)
	HRP Conjugated Human/Mouse MASP2 Recombinant Rabbit mAb [PSH20-52] - Detector	HA725392H	Human, Mouse, Rhesus monkey	ELISA(Det),ELISA
	Biotin Conjugated Human/Mouse/Rat FGF1 Recombinant Rabbit mAb [PSH18-35] - Detector	HA723988B	Human, Mouse, Rat, Cynomolgus monkey	ELISA(Det),ELISA
	Human/Mouse/Rat FGF1 Recombinant Rabbit mAb [PSH18-34] - BSA and Azide free (Capture)	HA723987	Human, Mouse, Rat, Cynomolgus monkey	ELISA(Cap)
	Human/Mouse/Rat FGF1 Recombinant Rabbit mAb [PSH18-35] - BSA and Azide free (Detector)	HA723988	Human, Mouse, Rat, Cynomolgus monkey	ELISA(Det)
	Biotin Conjugated Human NCAM1 Recombinant Rabbit mAb [PSH08-09] - Detector	HA722950B	Human, Cynomolgus monkey	ELISA(Det)
	Human NCAM1 Recombinant Rabbit mAb [PSH08-07] - BSA and Azide free (Capture)	HA722946	Human, Cynomolgus monkey	ELISA(Cap)
	Biotin Conjugated Human NCAM1 Recombinant Rabbit mAb [PSH08-08] - Detector	HA722948B	Human, Cynomolgus monkey	ELISA(Det)
	Human NCAM1 Recombinant Rabbit mAb [PSH08-08] - BSA and Azide free (Capture/Detector)	HA722947	Human, Cynomolgus monkey	ELISA(Det),ELISA(Cap)
	Human NCAM1 Recombinant Rabbit mAb [PSH08-09] - BSA and Azide free (Detector)	HA722949	Human, Cynomolgus monkey	ELISA(Det)
	Human Osteopontin Recombinant Rabbit mAb [PSH08-58] - BSA and Azide free (Capture)	HA723004	Human, Cynomolgus monkey	ELISA(Cap)
	Human Osteopontin Recombinant Rabbit mAb [PSH08-59] - BSA and Azide free (Detector)	HA723005	Human, Cynomolgus monkey	ELISA(Det)
	Biotin Conjugated Human Osteopontin Recombinant Rabbit mAb [PSH08-59] - Detector	HA723007B	Human, Cynomolgus monkey	ELISA(Det),ELISA
	Human PD-L1 Recombinant Rabbit mAb [PSH15-37] - BSA and Azide free (Detector)	HA723741	Human, Cynomolgus monkey	ELISA(Det)
	Human PD-L1 Recombinant Rabbit mAb [PSH15-36] - BSA and Azide free (Capture)	HA723740	Human, Cynomolgus monkey	ELISA(Cap)
	Biotin Conjugated Human PD-L1 Recombinant Rabbit mAb [PSH15-37] - Detector	HA723742B	Human, Cynomolgus monkey	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human Prolactin Recombinant Rabbit mAb [PSH03-14] - Detector	HA722627B	Human, Cynomolgus monkey	ELISA(Det)
	Human Prolactin Recombinant Rabbit mAb [PSH03-13] - BSA and Azide free (Capture)	HA721941	Human, Cynomolgus monkey	ELISA(Cap)
	Human Prolactin Recombinant Rabbit mAb [PSH03-14] - BSA and Azide free (Detector)	HA721942	Human, Cynomolgus monkey	ELISA(Det)
	Biotin Conjugated Human Prostate Specific Antigen Recombinant Rabbit mAb [PSH20-05] - Detector	HA725258B	Human, Cynomolgus monkey	ELISA(Det),ELISA
	Human Prostate Specific Antigen Recombinant Rabbit mAb [PSH20-05] - BSA and Azide free (Detector)	HA725258	Human, Cynomolgus monkey	ELISA(Det)
	Human Prostate Specific Antigen Recombinant Rabbit mAb [PSH20-04] - BSA and Azide free (Capture)	HA725257	Human, Cynomolgus monkey	ELISA(Cap)
	Human Resistin Recombinant Rabbit mAb [PSH14-62] - BSA and Azide free (Detector)	HA723659	Human, Cynomolgus monkey	ELISA(Det)
	Human Resistin Recombinant Rabbit mAb [PSH14-63] - BSA and Azide free (Capture)	HA723661	Human, Cynomolgus monkey	ELISA(Cap)
	Biotin Conjugated Human Resistin Recombinant Rabbit mAb [PSH14-62] - Detector	HA723660B	Human, Cynomolgus monkey	ELISA(Det),ELISA
	Human Resistin Recombinant Rabbit mAb [PSH14-61] - BSA and Azide free (Capture)	HA723658	Human, Cynomolgus monkey	ELISA(Cap)
	Human TNF alpha Recombinant Rabbit mAb [PS00-40] - BSA and Azide free (Capture)	HA721265	Human, Dog, Cynomolgus monkey	ELISA(Cap)
	Human TNF alpha Recombinant Rabbit mAb [PS01-07] - BSA and Azide free (Detector)	HA721266	Human, Dog, Cynomolgus monkey	ELISA(Det)
	Biotin Conjugated Human TNF alpha Recombinant Rabbit mAb [PS01-07] - Detector	HA721559	Human, Dog, Cynomolgus monkey	ELISA(Det)
	Biotin Conjugated Human VEGF Receptor 1 Recombinant Rabbit mAb [PSH19-23] - Detector	HA725366B	Human, Cynomolgus monkey	ELISA(Det),ELISA
	Human VEGF Receptor 1 Recombinant Rabbit mAb [PSH19-22] - BSA and Azide free (Capture)	HA725365	Human, Cynomolgus monkey	ELISA(Cap)
	Human VEGF Receptor 1 Recombinant Rabbit mAb [PSH19-23] - BSA and Azide free (Detector)	HA725366	Human, Cynomolgus monkey	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human YKL-40 Recombinant Rabbit mAb [PSH13-09] - BSA and Azide free (Detector)	HA725089	Human, Cynomolgus monkey	ELISA(Det)
	Biotin Conjugated Human YKL-40 Recombinant Rabbit mAb [PSH13-09] - Detector	HA725090B	Human, Cynomolgus monkey	ELISA(Det),ELISA
	HRP Conjugated Human YKL-40 Recombinant Rabbit mAb [PSH13-09] - Detector	HA725090H	Human, Cynomolgus monkey	ELISA(Det),ELISA
	Human YKL-40 Recombinant Rabbit mAb [PSH13-08] - BSA and Azide free (Capture)	HA725088	Human, Cynomolgus monkey	ELISA(Cap)
	IARS2 Rabbit pAb	ER63724	Human, Mouse, Monkey	WB,IHC-P,IF-Cell
	IDH1 Mouse mAb [A11-B4]	EM40705	Human, Mouse, Rat, Green monkey, Zebrafish	WB,IHC-P
	IFT88 Recombinant Rabbit mAb [PSH11-61]	HA723357	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell
	IFT88 Recombinant Rabbit mAb [PSH11-61] - BSA and Azide free	HA751417	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell
Relative expression	IGFBP2 Recombinant Rabbit mAb [JU37-46]	ET7106-63	Human, Green monkey	WB,IHC-P
Relative expression	IGFBP2 Recombinant Rabbit mAb [PSH21-02] - BSA and Azide free	HA751784	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-Fr,IHC-P,FC,IP
Relative expression	IGFBP2 Recombinant Rabbit mAb [PSH21-02]	HA724197	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-Fr,IHC-P,FC,IP
	IKB beta Recombinant Rabbit mAb [SN0648]	ET1611-1	Human, Green monkey	WB,IF-Cell,IHC-P,IP,FC
KD	Integrin alpha 5 Recombinant Rabbit mAb [JJ08-94]	ET1701-58	Human, Mouse, Green monkey, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Integrin beta 5 Rabbit pAb	ER63768	Human, Mouse, Rat, Monkey	WB
	Integrin beta 5 Recombinant Rabbit mAb [PSH10-13]	HA723182	Human, Rat, Green monkey	WB,IHC-P
	Integrin beta 5 Recombinant Rabbit mAb [PSH10-13] - BSA and Azide free	HA751342	Human, Rat, Green monkey	WB,IHC-P
Relative expression	IRAK4 Recombinant Rabbit mAb [JE41-77]	HA721347	Human, Mouse, Green monkey	WB,IHC-P
	IRAK4 Mouse mAb [B9-H2]	EM1709-86	Human, Mouse, Monkey	WB,IHC-P,FC
	IRE1 Recombinant Rabbit mAb [PSH10-60] - BSA and Azide free	HA751363	Human, Mouse, Rat, Green monkey	WB
	IRE1 Recombinant Rabbit mAb [PSH10-60]	HA723225	Human, Mouse, Rat, Green monkey	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment, Relative expression	ISG15 Mouse mAb	RT1674	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC
KD	JMJD6 Recombinant Rabbit mAb [PSH15-03] - BSA and Azide free	HA751541	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC,IP,ChIP
KD	JMJD6 Recombinant Rabbit mAb [PSH15-03]	HA723703	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC,IP,ChIP
	KAT1 / HAT1 Recombinant Rabbit mAb [JE36-62]	HA722332	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,IP,IF-Tissue
	KAT3B / p300 Rabbit pAb	ER1908-18	Human, Mouse, Rat, Green monkey	WB
	KAT3B / p300 Recombinant Rabbit mAb [PSH15-31] - BSA and Azide free	HA751553	Human, Mouse, Rat, Green monkey	WB,IF-Cell
	KAT3B / p300 Recombinant Rabbit mAb [PSH15-31]	HA723737	Human, Mouse, Rat, Green monkey	WB,IF-Cell
	KDM1 / LSD1 Recombinant Rabbit mAb [JE04-60]	HA722723	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,IP
	KGF Recombinant Rabbit mAb [JE40-14]	HA722404	Human, Mouse, Rat, Green monkey	WB
	KIR3DL1 Mouse mAb [6D9F6]	EM1710-62	Human, Rat, Monkey	WB,IHC-P,FC
	KMT1E / SETDB1 Recombinant Rabbit mAb [PSH20-40]	HA724151	Human, Mouse, Rat, Green monkey	WB,IHC-P,IP
	KMT1E / SETDB1 Recombinant Rabbit mAb [PSH20-39]	HA724150	Human, Mouse, Rat, Green monkey	WB,IHC-P,IP
	KMT1E / SETDB1 Recombinant Rabbit mAb [PSH20-40] - BSA and Azide free	HA751756	Human, Mouse, Rat, Green monkey	WB,IHC-P,IP
	KMT1E / SETDB1 Recombinant Rabbit mAb [PSH20-39] - BSA and Azide free	HA751755	Human, Mouse, Rat, Green monkey	WB,IHC-P,IP
	KMT6 / EZH2 Recombinant Rabbit mAb [PSH04-14] - BSA and Azide free	HA750937	Human, Mouse, Rat, Green monkey	WB,IF-Cell
	KMT6 / EZH2 Recombinant Rabbit mAb [PSH04-14]	HA722095	Human, Mouse, Rat, Green monkey	WB,IF-Cell
	KSR1 Recombinant Rabbit mAb [JE57-75]	HA721963	Human, Mouse, Rat, Green monkey	WB,IHC-P,FC
	Ku70 Rabbit pAb	R1401-5	Human, Green monkey	WB,IHC-P,IF-Cell,FC
	Ku70 Mouse mAb [D1-A6]	EM50109	Human, Mouse, Rat, Green monkey	WB,IHC-P
	Ku80 Recombinant Rabbit mAb [SC06-14]	ET1610-40	Human, Green monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Kv3.4 Rabbit pAb	ER63826	Human, Mouse, Monkey	WB,IHC-P,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Lamin B2 Recombinant Mouse mAb [J9-C4-R] - BSA and Azide free	HA610162	Human, Green monkey	WB,IF-Cell,IHC-P,IF-Tissue
	Lamin B2 Recombinant Mouse mAb [J9-C4-R]	HA601289	Human, Green monkey	WB,IF-Cell,IHC-P,IF-Tissue
	Laminin gamma-1 Rabbit pAb	ER63838	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell
	LARP1 Recombinant Rabbit mAb [PSH12-09] - BSA and Azide free	HA751438	Human, Green monkey	WB,IF-Cell
	LARP1 Recombinant Rabbit mAb [PSH12-09]	HA723425	Human, Green monkey	WB,IF-Cell
	LARP1 Rabbit pAb	ER60797	Human, Rat, Green monkey	WB,IF-Cell,IHC-P,FC
	LATS1 / WARTS Recombinant Rabbit mAb [PSH07-20] - BSA and Azide free	HA751141	Human, Mouse, Green monkey	WB
	LATS1 / WARTS Recombinant Rabbit mAb [PSH07-20]	HA722808	Human, Mouse, Green monkey	WB
	LDHB Recombinant Rabbit mAb [PSH22-91] - BSA and Azide free	HA751880	Human, Mouse, Rat, Green monkey, Pig	WB,IF-Cell,FC,IP
	LDHB Recombinant Rabbit mAb [PSH22-91]	HA724336	Human, Mouse, Rat, Green monkey, Pig	WB,IF-Cell,FC,IP
	LIMK2 Recombinant Rabbit mAb [PSH08-82] - BSA and Azide free	HA751247	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,IP
	LIMK2 Recombinant Rabbit mAb [PSH08-82]	HA723034	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,IP
Relative expression	LIM Kinase 1 Recombinant Rabbit mAb [PSH08-68] - BSA and Azide free	HA751233	Human, Mouse, Rat, Green monkey	WB,IHC-P,FC
Relative expression	LIM Kinase 1 Recombinant Rabbit mAb [PSH08-68]	HA723020	Human, Mouse, Rat, Green monkey	WB,IHC-P,FC
KD	LIMPII Recombinant Rabbit mAb [PSH03-25] - BSA and Azide free	HA750865	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell
KD	LIMPII Recombinant Rabbit mAb [PSH03-25]	HA721967	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell
	LIS1 Recombinant Rabbit mAb [JE37-43]	HA722143	Human, Mouse, Rat, Green monkey	WB
Relative expression	LKB1 Recombinant Rabbit mAb [PSH09-29] - BSA and Azide free	HA751283	Human, Mouse, Rat, Green monkey	WB,IHC-P,IP
Relative expression	LKB1 Recombinant Rabbit mAb [PSH09-29]	HA723087	Human, Mouse, Rat, Green monkey	WB,IHC-P,IP
	LMO3 Rabbit pAb	ER63861	Human, Mouse, Rat, Monkey	WB,IHC-P
	LPCAT3 Recombinant Rabbit mAb [PSH10-02]	HA723171	Human, Green monkey	WB,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	LPCAT3 Recombinant Rabbit mAb [PSH10-02] - BSA and Azide free	HA751332	Human, Green monkey	WB,IF-Cell
	LRPPRC Recombinant Rabbit mAb [PSH24-58]	HA724454	Human, Mouse, Rat, Green monkey, Pig	WB,IF-Cell,IHC-P,FC(Intra)
	LRPPRC Recombinant Rabbit mAb [PSH24-58] - BSA and Azide free	HA751948	Human, Mouse, Rat, Green monkey, Pig	WB,IF-Cell,IHC-P,FC(Intra)
KD	Lyn Recombinant Rabbit mAb [JE44-75]	HA721991	Human, Mouse, Rat, Cynomolgus monkey	WB,FC
Relative expression	MAPRE1 / EB1 Recombinant Rabbit mAb [PSH25-89]	HA724581	Human, Mouse, Rat, Pig, Green monkey	WB,IF-Cell,IHC-P,FC(Intra),IP
Relative expression	MAPRE1 / EB1 Recombinant Rabbit mAb [PSH25-89] - BSA and Azide free	HA752001	Human, Mouse, Rat, Pig, Green monkey	WB,IF-Cell,IHC-P,FC(Intra),IP
	MCM2 Recombinant Antibody [PD00-90] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610350	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P
	MCM2 Recombinant Antibody [PD00-90] - Mouse IgG1 (Chimeric)	HA601523	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P
	MCM2 Recombinant Rabbit mAb [PD00-90]	HA721177	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
	MCM7 Recombinant Mouse mAb [15E1-R]	HA601288	Human, Mouse, Rat, Green monkey	WB
	MCM7 Mouse mAb [15E1]	EM1901-12	Human, Mouse, Rat, Green monkey	WB,IHC-P,FC
	MCM7 Recombinant Mouse mAb [15E1-R] - BSA and Azide free	HA610161	Human, Mouse, Rat, Green monkey	WB
	MCU Recombinant Rabbit mAb [PSH13-74]	HA723563	Human, Mouse, Rat, Green monkey	WB,IHC-Fr,IHC-P,IF-Cell,IP
	MCU Mouse mAb [A2E11]	EM1901-85	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell
	MCU Recombinant Rabbit mAb [PSH13-74] - BSA and Azide free	HA751487	Human, Mouse, Rat, Green monkey	WB,IHC-Fr,IHC-P,IF-Cell,IP
Cell treatment	MDM2 Recombinant Mouse mAb [A10E11-R]	HA601310	Human, Green monkey	WB,IF-Cell,IHC-P
	MEK4 / MKK4 Recombinant Rabbit mAb [PSH11-99]	HA723412	Human, Mouse, Rat, Green monkey	WB,IHC-P,FC
	MEK4 / MKK4 Recombinant Rabbit mAb [PSH11-99] - BSA and Azide free	HA751432	Human, Mouse, Rat, Green monkey	WB,IHC-P,FC
Relative expression, KD	Merlin Recombinant Rabbit mAb [PSH02-52] - BSA and Azide free	HA750815	Human, Mouse, Rat, Green monkey	WB,IF-Cell,FC
Relative expression, KD	Merlin Recombinant Rabbit mAb [PSH02-52]	HA721828	Human, Mouse, Rat, Green monkey	WB,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	METTL14 Recombinant Rabbit mAb [PSH07-85] - BSA and Azide free	HA751192	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC,IP
KD	METTL14 Recombinant Rabbit mAb [PSH07-85]	HA722910	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC,IP
	MFF Recombinant Rabbit mAb [PSH14-36] - BSA and Azide free	HA751518	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,FC,IP
	MFF Recombinant Rabbit mAb [PSH14-36]	HA723630	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,FC,IP
	Mindin Rabbit pAb	ER63958	Human, Mouse, Rat, Monkey	WB
	MMP-1 Rabbit pAb	ER31211	Human, Mouse, Green monkey	WB,IHC-P
	MOB1 Recombinant Rabbit mAb [PSH17-36] - BSA and Azide free	HA751624	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC,IP
	MOB1 Recombinant Rabbit mAb [PSH17-36]	HA723919	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC,IP
	Monkey Albumin/ALB, C-His Tag Protein	HA211229	Monkey	
	Monkey CRP/C-Reactive Protein enzyme-linked immunoassay kit	EH0281	Monkey	ELISA
	Biotin Conjugated Monkey PCSK9 Recombinant Mouse mAb [PSH23-55] - Detector	HA601603B	Cynomolgus monkey	ELISA(Det),ELISA
	HRP Conjugated Monkey PCSK9 Recombinant Mouse mAb [PSH23-55] - Detector	HA601603H	Cynomolgus monkey	ELISA(Det),ELISA
	Monkey PCSK9 Recombinant Mouse mAb [PSH23-55] - BSA and Azide free (Detector)	HA601603	Cynomolgus monkey	ELISA(Det)
	Monkey PCSK9 Recombinant Rabbit mAb [PSH23-54] - BSA and Azide free (Capture)	HA724385	Cynomolgus monkey	ELISA(Cap)
	MRP7 Rabbit pAb	ER64010	Human, Monkey	WB
	MRP-S21 Rabbit pAb	ER64048	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell
Relative expression	MSK1 Recombinant Rabbit mAb [PSH13-72] - BSA and Azide free	HA751485	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC
Relative expression	MSK1 Recombinant Rabbit mAb [PSH13-72]	HA723561	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC
	MSP / MST1 Recombinant Rabbit mAb [JE39-26]	HA722342	Human, Mouse, Rat, Green monkey	WB
	MTA1 Recombinant Rabbit mAb [PSH21-25] - BSA and Azide free	HA751791	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC,IP
	MTA1 Recombinant Rabbit mAb [PSH21-25]	HA724208	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	MVK Mouse mAb [A9G12]	HA601151	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell
	MyD88 Recombinant Rabbit mAb [PSH09-74]	HA723139	Human, Mouse, Rat, Green monkey	WB
	MyD88 Recombinant Rabbit mAb [PSH09-74] - BSA and Azide free	HA751316	Human, Mouse, Rat, Green monkey	WB
	MyD88 Recombinant Rabbit mAb [PSH05-58] - BSA and Azide free	HA751021	Human, Green monkey	WB,IF-Cell,IHC-P,FC,IP
	MyD88 Recombinant Rabbit mAb [PSH05-58]	HA722378	Human, Green monkey	WB,IF-Cell,IHC-P,FC,IP
	Myf-5 Rabbit pAb	ER64072	Human, Mouse, Monkey	WB,IHC-P,IF-Cell
	MYH14 Recombinant Rabbit mAb [PSH20-41] - BSA and Azide free	HA751757	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,FC,IP
	MYH14 Recombinant Rabbit mAb [PSH20-41]	HA724152	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,FC,IP
	MYL7 Rabbit pAb	ER64080	Human, Mouse, Rat, Monkey	WB
	Na+ CP type VII alpha Rabbit pAb	ER64101	Human, Monkey	WB,IF-Cell
	NDUFA8 Recombinant Rabbit mAb [PSH21-96]	HA724265	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,FC,IP
	NDUFA8 Recombinant Rabbit mAb [PSH21-96] - BSA and Azide free	HA751812	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,FC,IP
	NDUFA9 Recombinant Rabbit mAb [PSH14-07] - BSA and Azide free	HA751502	Human, Mouse, Rat, Green monkey	WB,IHC-P,IP
	NDUFA9 Recombinant Rabbit mAb [PSH14-07]	HA723597	Human, Mouse, Rat, Green monkey	WB,IHC-P,IP
	NDUFA9 Recombinant Rabbit mAb [PSH14-08]	HA723598	Human, Mouse, Rat, Green monkey	WB,IP
	NDUFA9 Recombinant Rabbit mAb [PSH14-08] - BSA and Azide free	HA751503	Human, Mouse, Rat, Green monkey	WB,IP
	NDUFS1 Recombinant Rabbit mAb [PSH14-45] - BSA and Azide free	HA751524	Human, Mouse, Rat, Green monkey	WB,IHC-P
	NDUFS1 Recombinant Rabbit mAb [PSH14-45]	HA723641	Human, Mouse, Rat, Green monkey	WB,IHC-P
	NDUFS1 Rabbit pAb	ER64138	Human, Mouse, Rat, Green monkey	WB
	NDUFV1 Rabbit pAb	ER60893	Human, Mouse, Rat, Green monkey	WB
	NDUFV1 Recombinant Rabbit mAb [PSH21-97] - BSA and Azide free	HA751813	Human, Mouse, Rat, Green monkey, Pig	WB,IF-Cell,IHC-P,IP
	NDUFV1 Recombinant Rabbit mAb [PSH21-97]	HA724266	Human, Mouse, Rat, Green monkey, Pig	WB,IF-Cell,IHC-P,IP
	NEDD8 Recombinant Rabbit mAb [JF69-09]	ET1702-84	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	NELF-B Recombinant Rabbit mAb [JE77-56]	HA722320	Human, Mouse, Rat, Green monkey	WB
KD	Neurofibromin Recombinant Rabbit mAb [PSH25-01]	HA724496	Human, Mouse, Rat, Green monkey	WB
KD	Neurofibromin Recombinant Rabbit mAb [PSH25-01] - BSA and Azide free	HA751964	Human, Mouse, Rat, Green monkey	WB
KO, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [PSH0-27]	HA721307	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,FC
KO, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [PSH0-27] - BSA and Azide free	HA750601	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,FC
	NFR-kappa B Rabbit pAb	ER64181	Human, Mouse, Monkey	WB,IHC-P
	Nkx-6.3 Rabbit pAb	ER64205	Human, Mouse, Monkey	WB,IF-Cell
	NMDA epsilon 4 Rabbit pAb	ER64208	Human, Mouse, Rat, Monkey	WB
	NSUN2 Recombinant Rabbit mAb [PSH05-05] - BSA and Azide free	HA750978	Human, Mouse, Green monkey	WB,IF-Cell,IHC-P,IF-Tissue,IP
	NSUN2 Recombinant Rabbit mAb [PSH05-05]	HA722219	Human, Mouse, Green monkey	WB,IF-Cell,IHC-P,IF-Tissue,IP
Relative expression	NUMB Recombinant Rabbit mAb [PSH25-86]	HA724578	Human, Rat, Green monkey	WB,FC(Intra),IP
Relative expression	NUMB Recombinant Rabbit mAb [PSH25-86] - BSA and Azide free	HA751998	Human, Rat, Green monkey	WB,FC(Intra),IP
	Olfactory receptor 10S1 Rabbit pAb	ER64284	Human, Monkey	WB,IF-Cell
	OTUB1 Recombinant Rabbit mAb [PSH18-50] - BSA and Azide free	HA751674	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,FC,IP
	OTUB1 Recombinant Rabbit mAb [PSH18-50]	HA724003	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,FC,IP
	OTULIN Rabbit pAb	HA500587	Human, Mouse, Rat, Green monkey	WB
	p18 Rabbit pAb	ER64492	Human, Mouse, Monkey	WB,IHC-P
	p38 alpha/MAPK14 + p38 beta/MAPK11 + p38 gamma/MAPK12 Recombinant Rabbit mAb [PSH03-60] - BSA and Azide free	HA750894	Human, Mouse, Rat, Green monkey	WB
	p38 alpha/MAPK14 + p38 beta/MAPK11 + p38 gamma/MAPK12 Recombinant Rabbit mAb [PSH03-60]	HA722019	Human, Mouse, Rat, Green monkey	WB
	p38 beta/MAPK11 + p38 gamma/MAPK12 + p38 delta/MAPK13 Recombinant Rabbit mAb [PSH03-59]	HA722018	Human, Mouse, Rat, Green monkey	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	p38 beta/MAPK11 + p38 gamma/MAPK12 + p38 delta/MAPK13 Recombinant Rabbit mAb [PSH03-59] - BSA and Azide free	HA750893	Human, Mouse, Rat, Green monkey	WB
	p38 MAPK Recombinant Rabbit Multiclonal Antibody [PSH03-72]	HA722036	Human, Mouse, Rat, Green monkey	WB
	p38 MAPK Recombinant Rabbit Multiclonal Antibody [PSH03-72] - BSA and Azide free	HA750906	Human, Mouse, Rat, Green monkey	WB
KD	P4HB Recombinant Rabbit mAb [JE54-98]	ET7110-92	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC
Relative expression	p53 Mouse mAb [6B2]	EM1701-91	Human, Mouse, Green monkey	WB,IF-Cell
KD	p53 Recombinant Mouse mAb [6B2-R]	HA601322	Human, Green monkey	WB,IF-Cell
Cell treatment, Relative expression, KD	p53 Recombinant Mouse mAb [6B5-R]	HA601315	Human, Mouse, Green monkey	WB,IF-Cell
Cell treatment, KD, Relative expression	p53 Recombinant Mouse mAb [6B5-R] - BSA and Azide free	HA610187	Human, Mouse, Green monkey	WB,IF-Cell
KD	p53 Recombinant Mouse mAb [6B2-R] - BSA and Azide free	HA610194	Human, Green monkey	WB,IF-Cell
	p70 S6 Kinase Recombinant Rabbit mAb [PSH13-50] - BSA and Azide free	HA751469	Human, Mouse, Rat, Green monkey	WB,IF-Cell,FC,IP
	p70 S6 Kinase Recombinant Rabbit mAb [PSH13-50]	HA723539	Human, Mouse, Rat, Green monkey	WB,IF-Cell,FC,IP
	PABP Recombinant Rabbit mAb [PSH05-51]	HA722289	Human, Mouse, Green monkey	WB,IF-Cell,IHC-P,FC,IF-Tissue
	PABP Recombinant Rabbit mAb [PSH05-51] - BSA and Azide free	HA751010	Human, Mouse, Green monkey	WB,IF-Cell,IHC-P,FC,IF-Tissue
Relative expression	PADI2 Recombinant Rabbit mAb [PSH18-07]	HA723970	Human, Green monkey	WB,IP
Relative expression	PADI2 Recombinant Rabbit mAb [PSH18-07] - BSA and Azide free	HA751653	Human, Green monkey	WB,IP
Relative expression, Cell treatment	PAI1 Recombinant Rabbit mAb [PSH24-33] - BSA and Azide free	HA751934	Human, Mouse, Green monkey	WB,IP
Relative expression, Cell treatment	PAI1 Recombinant Rabbit mAb [PSH24-33]	HA724440	Human, Mouse, Green monkey	WB,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	PAK2 Recombinant Rabbit mAb [JB41-84]	ET7107-79	Human, Mouse, Rat, Green monkey	WB,IHC-P,FC,IF-Cell,IF-Tissue
	Parkin Recombinant Rabbit mAb [PSH08-11]	HA722952	Human, Mouse, Rat, Green monkey, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IF-Cell,FC,IP,IHC-Fr
	Parkin Recombinant Rabbit mAb [PSH08-11] - BSA and Azide free	HA751206	Human, Mouse, Rat, Green monkey, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IF-Cell,FC,IP,IHC-Fr
Relative expression	PAX8 Recombinant Mouse mAb [6G5-R]	HA601270	Human, Mouse, Green monkey	WB,IHC-P,IF-Cell
Relative expression	PAX8 Recombinant Mouse mAb [6G5-R] - BSA and Azide free	HA610145	Human, Mouse, Green monkey	WB,IHC-P,IF-Cell
	PEA-15 Rabbit pAb	ER64545	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell
KD	PGM1 Recombinant Rabbit mAb [JE65-06]	HA721110	Human, Mouse, Rat, Green monkey	WB,IF-Cell,FC
KD	PHF6 Recombinant Rabbit mAb [PSH02-35] - BSA and Azide free	HA750800	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell
KD	PHF6 Recombinant Rabbit mAb [PSH02-35]	HA721812	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell
Cell treatment	Phospho-BCKDHA (S293) Recombinant Rabbit mAb [PSH14-02] - BSA and Azide free	HA751498	Human, Mouse, Rat, Green monkey	WB
Cell treatment	Phospho-BCKDHA (S293) Recombinant Rabbit mAb [PSH14-02]	HA723592	Human, Mouse, Rat, Green monkey	WB
Cell treatment	Phospho-Cofilin (S3) Recombinant Rabbit mAb [PSH07-18]	HA722806	Human, Rat, Green monkey	WB
Cell treatment	Phospho-Cofilin (S3) Recombinant Rabbit mAb [PSH07-18] - BSA and Azide free	HA751139	Human, Rat, Green monkey	WB
Cell treatment	Phospho-Histone H2A.X (S139) Recombinant Rabbit mAb [SR33-09]	ET1602-2	Human, Mouse, Rat, Cynomolgus monkey	WB,IHC-P,IF-Cell,IF-Tissue
Cell treatment	Phospho-Histone H2A.X (S139) Recombinant Rabbit mAb [SR33-09] - BSA and Azide free	HA750034	Human, Mouse, Rat, Cynomolgus monkey	WB,IHC-P,IF-Cell,IF-Tissue
Cell treatment	Phospho-Rb (S608) Recombinant Rabbit mAb [PSH11-04]	HA723271	Human, Green monkey	WB,IF-Cell
Cell treatment	Phospho-Rb (S608) Recombinant Rabbit mAb [PSH11-04] - BSA and Azide free	HA751394	Human, Green monkey	WB,IF-Cell
Cell treatment	Phospho-Src (S17) Recombinant Rabbit mAb [PSH11-03] - BSA and Azide free	HA751393	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC,IP
Cell treatment	Phospho-Src (S17) Recombinant Rabbit mAb [PSH11-03]	HA723270	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	PI 3 Kinase p85 alpha/gamma Rabbit pAb	ER64588	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Tissue
	PIASx Rabbit pAb	ER64592	Human, Mouse, Rat, Monkey	WB
	PKA C-alpha Recombinant Rabbit mAb [PSH08-57]	HA723003	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,FC,IP
	PKA C-alpha Recombinant Rabbit mAb [PSH08-57] - BSA and Azide free	HA751232	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,FC,IP
	PKC eta Recombinant Rabbit mAb [JE32-41]	HA722389	Human, Mouse, Rat, Green monkey	WB
Relative expression	PKM2 Recombinant Rabbit mAb [PSH11-74]	HA723370	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC,IP
Relative expression	PKM2 Recombinant Rabbit mAb [PSH11-74] - BSA and Azide free	HA751429	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC,IP
	PLCG1 Mouse mAb [D11-C2]	EM1707-02	Human, Rat, Monkey	WB,IF-Cell,IHC-P,FC
	PLK1 Recombinant Rabbit mAb [PSH04-35]	HA722118	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P
	PLK1 Recombinant Rabbit mAb [PSH04-35] - BSA and Azide free	HA750944	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P
	Pmp24 Rabbit pAb	ER64634	Human, Monkey	WB,IHC-P
	POLR2A Recombinant Rabbit mAb [PSH18-30]	HA723983	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC,IP,ChIP
	POLR2A Recombinant Rabbit mAb [PSH18-30] - BSA and Azide free	HA751666	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC,IP,ChIP
	POLR3A Rabbit pAb	ER64642	Human, Mouse, Monkey	WB,IHC-P
	PP4R2 Rabbit pAb	ER64658	Human, Monkey	WB,IHC-P
	PPAR alpha Mouse mAb	EM1707-71	Human, Monkey	FC
	PRA1 Recombinant Rabbit mAb [PSH01-96] - BSA and Azide free	HA750770	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell
	PRA1 Recombinant Rabbit mAb [PSH01-96]	HA721754	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell
	PRAK / MK5 Recombinant Rabbit mAb [PSH14-06]	HA723596	Human, Mouse, Rat, Green monkey	WB,IP
	PRAK / MK5 Recombinant Rabbit mAb [PSH14-06] - BSA and Azide free	HA751501	Human, Mouse, Rat, Green monkey	WB,IP
	PRKAA2 Mouse mAb [H5-E11]	EM1708-49	Human, Monkey, Mouse, Rat	WB,IHC-P,IF-Cell
	PRMT5 Recombinant Rabbit mAb [ST51-06] - BSA and Azide free	HA750174	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	PRMT5 Recombinant Rabbit mAb [ST51-06]	ET1609-43	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment, Relative expression, KO	Pro Caspase-3 Recombinant Rabbit mAb [SR03-01] - BSA and Azide free	HA750048	Human, Green monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Cell treatment, Relative expression, KO	Pro Caspase-3 Recombinant Rabbit mAb [SR03-01]	ET1602-39	Human, Green monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	Pro Caspase-3 Recombinant Rabbit mAb [SZ02-08]	ET1603-26	Human, Green monkey, Mouse	WB,IF-Tissue,IHC-P
	Prohibitin(Mitochondrial Marker) Rabbit pAb	R1510-30	Human, Mouse, Rat, Green monkey	WB,IHC-P
	PSMA2 Recombinant Rabbit mAb [JE51-51]	ET7110-19	Human, Mouse, Rat, Green monkey	WB,IHC-P
	PSMA7 / HSPC Rabbit pAb	HA500004	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,FC
	PSMA7 / HSPC Recombinant Rabbit mAb [PSH14-09]	HA723599	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC
	PSMA7 / HSPC Recombinant Rabbit mAb [PSH14-09] - BSA and Azide free	HA751504	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC
Cell treatment	PSMB5 / MB1 Rabbit pAb	HA500533	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,FC
Cell treatment	PSMB5 / MB1 Recombinant Rabbit mAb [PSH15-20]	HA723724	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC
Cell treatment	PSMB5 / MB1 Recombinant Rabbit mAb [PSH15-20] - BSA and Azide free	HA751549	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC
	PSMB7 Rabbit pAb	HA500345	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,FC
	PSMC2 Rabbit pAb	HA500404	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC
	PSMD4 Rabbit pAb	HA500390	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P
	PTG Rabbit pAb	ER64697	Human, Monkey	WB,IHC-P
KD	RAB10 Recombinant Rabbit mAb [PSH03-24]	HA721966	Human, Mouse, Rat, Green monkey	WB
KD	RAB10 Recombinant Rabbit mAb [PSH03-24] - BSA and Azide free	HA750864	Human, Mouse, Rat, Green monkey	WB
	Rab4A Recombinant Rabbit mAb [SR27-08]	ET1602-15	Human, Rat, Green monkey, Mouse	WB,IF-Cell,IF-Tissue,IP
	Rad51C Rabbit pAb	ER64745	Human, Monkey	WB,IHC-P,IF-Cell
	RALA Mouse mAb	EM1707-77	Human, Mouse, Monkey	WB,IF-Cell,IHC-P,FC
	RALB Mouse mAb [B4-F8]	EM1707-30	Human, Mouse, Monkey	WB,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	RALBP1 Recombinant Rabbit mAb [SC06-49]	ET1610-88	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,FC
	RAN Mouse mAb [A1-6D]	EM1706-18	Human, Mouse, Monkey, Rat	WB,IHC-P,IF-Cell,FC
	RBBP7 Rabbit pAb	HA500155	Human, Mouse, Rat, Green monkey	WB,IHC-P
	Rbm15 Mouse mAb [A10H4]	HA601347	Human, Mouse, Rat, Green monkey	WB,IHC-P
	RCC1 Recombinant Rabbit mAb [JE39-10]	HA722311	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,IP
	RHOG Recombinant Rabbit mAb [PSH02-70] - BSA and Azide free	HA750827	Human, Green monkey	WB,FC
	RHOG Recombinant Rabbit mAb [PSH02-70]	HA721848	Human, Green monkey	WB,FC
	Ribosomal Protein L10L Rabbit pAb	ER64821	Human, Monkey	WB,IHC-P
	RICTOR Mouse mAb [E9-B3]	EM1709-50	Human, Mouse, Monkey	WB,IF-Cell,IHC-P,FC
	RING2 Recombinant Rabbit mAb [JB38-41]	ET7107-07	Human, Mouse, Rat, Green monkey	WB,IHC-P
	RIP Recombinant Rabbit mAb [PSH19-57] - BSA and Azide free	HA751721	Human, Mouse, Rat, Green monkey	WB
	RIP Recombinant Rabbit mAb [PSH19-57]	HA724083	Human, Mouse, Rat, Green monkey	WB
	RPB3 Recombinant Rabbit mAb [JE34-81]	HA722448	Human, Mouse, Rat, Green monkey	WB
	RPL26 Rabbit pAb	HA500532	Human, Mouse, Green monkey, Rat	WB,IF-Cell,IHC-P,FC
	RPS19 Recombinant Rabbit mAb [JU33-43]	ET7106-94	Human, Mouse, Rat, Green monkey	WB,IHC-P,IP,IF-Cell,FC
	RRAS Rabbit pAb	HA500397	Human, Mouse, Rat, Green monkey	WB
	Rubicon / Baron Recombinant Rabbit mAb [PSH08-25] - BSA and Azide free	HA751216	Human, Mouse, Rat, Green monkey	WB
	Rubicon / Baron Recombinant Rabbit mAb [PSH08-25]	HA722969	Human, Mouse, Rat, Green monkey	WB
	RUVB2 Recombinant Mouse mAb [A1E12-R] - BSA and Azide free	HA610114	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P
	RUVB2 Recombinant Mouse mAb [A1E12-R]	HA601239	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P
	RUVB2 Mouse mAb [A1E12]	EM1901-58	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P
	SA2 Recombinant Rabbit mAb [JE57-34]	HA721908	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	SAFB Recombinant Rabbit mAb [JE64-94]	HA721085	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P
	SCNN1G Recombinant Rabbit mAb [PSH16-29]	HA723827	Human, Mouse, Rat, Green monkey	WB,IF-Cell
	SCNN1G Recombinant Rabbit mAb [PSH16-29] - BSA and Azide free	HA751597	Human, Mouse, Rat, Green monkey	WB,IF-Cell
	SEC24A Recombinant Rabbit mAb [PSH04-36]	HA722119	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC
	SEC24A Recombinant Rabbit mAb [PSH04-36] - BSA and Azide free	HA750945	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC
	SENP5 Rabbit pAb	ER64925	Human, Mouse, Monkey	WB,IHC-P
KD	SET7 Recombinant Rabbit mAb [PSH02-90] - BSA and Azide free	HA750845	Human, Mouse, Rat, Green monkey	WB,IF-Cell,FC
KD	SET7 Recombinant Rabbit mAb [PSH02-90]	HA721909	Human, Mouse, Rat, Green monkey	WB,IF-Cell,FC
	SETDB1 Mouse mAb [A12-H6]	EM1709-76	Human, Mouse, Green monkey	WB,IF-Cell
Relative expression	SHMT2 / SHMT Recombinant Rabbit mAb [PSH19-37]	HA724069	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,FC,IP
Relative expression	SHMT2 / SHMT Recombinant Rabbit mAb [PSH19-37] - BSA and Azide free	HA751707	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,FC,IP
	SIRT4 Recombinant Rabbit mAb [PSH15-46] - BSA and Azide free	HA751561	Human, Mouse, Rat, Green monkey	WB,IHC-P,IP
	SIRT4 Recombinant Rabbit mAb [PSH15-46]	HA723752	Human, Mouse, Rat, Green monkey	WB,IHC-P,IP
	SIRT6 Recombinant Rabbit mAb [PSH10-98]	HA723265	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P
KO	SIRT6 Recombinant Rabbit mAb [PSH10-99]	HA723266	Human, Mouse, Rat, Green monkey	WB,IHC-P
	SIRT6 Recombinant Rabbit mAb [PSH10-98] - BSA and Azide free	HA751388	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P
KO	SIRT6 Recombinant Rabbit mAb [PSH10-99] - BSA and Azide free	HA751389	Human, Mouse, Rat, Green monkey	WB,IHC-P
	SKP2 Recombinant Rabbit mAb [PSH06-22]	HA722654	Human, Green monkey	WB
	SKP2 Recombinant Rabbit mAb [PSH06-22] - BSA and Azide free	HA751060	Human, Green monkey	WB
Relative expression	SLC16A3 / MCT 4 Recombinant Rabbit mAb [PSH06-91] - BSA and Azide free	HA751115	Human, Green monkey	WB,IF-Cell,IHC-P,FC
Relative expression	SLC16A3 / MCT 4 Recombinant Rabbit mAb [PSH06-91]	HA722766	Human, Green monkey	WB,IF-Cell,IHC-P,FC
	SLC22A1 Recombinant Rabbit mAb [JE37-67]	HA722335	Human, Green monkey	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Smad2/3 Rabbit pAb	ER64981	Human, Mouse, Rat, Monkey	WB,IHC-P
	Smad7 Recombinant Mouse mAb [PSH21-17] - BSA and Azide free	HA610387	Human, Mouse, Rat, Green monkey	WB
	Smad7 Recombinant Mouse mAb [PSH21-17]	HA601573	Human, Mouse, Rat, Green monkey	WB
Relative expression	SMARCA2 Recombinant Rabbit mAb [PSH01-15] - BSA and Azide free	HA750691	Human, Green monkey	WB,IHC-P,FC
Relative expression	SMARCA2 Recombinant Rabbit mAb [PSH01-15]	HA721627	Human, Green monkey	WB,IHC-P,FC
	SMC2 Recombinant Rabbit mAb [PSH0-62] - BSA and Azide free	HA750639	Human, Green monkey	WB,IHC-P,IF-Cell
	SMC2 Recombinant Rabbit mAb [PSH0-62]	HA721450	Human, Green monkey	WB,IHC-P,IF-Cell
	SMC2 Recombinant Rabbit mAb [PSH0-63]	HA721451	Human, Green monkey	WB,IHC-P,IF-Cell
	SMC2 Recombinant Rabbit mAb [PSH0-63] - BSA and Azide free	HA750640	Human, Green monkey	WB,IHC-P,IF-Cell
	SMC5 Recombinant Rabbit mAb [PSH15-29] - BSA and Azide free	HA751551	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P
	SMC5 Recombinant Rabbit mAb [PSH15-29]	HA723735	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P
	SMCR7L / MID51 Recombinant Rabbit mAb [PSH24-77]	HA724474	Human, Mouse, Rat, Green monkey, Pig	WB,IP
	SMCR7L / MID51 Recombinant Rabbit mAb [PSH24-77] - BSA and Azide free	HA751955	Human, Mouse, Rat, Green monkey, Pig	WB,IP
	Smg1 Recombinant Rabbit mAb [JE38-38]	HA722452	Human, Mouse, Rat, Green monkey	WB
	SMUG1 Rabbit pAb	ER64988	Human, Mouse, Rat, Monkey	WB,IHC-P
Relative expression	SOCS1 Recombinant Rabbit mAb [PSH09-68]	HA723133	Human, Mouse, Rat, Green monkey	WB
Relative expression	SOCS1 Recombinant Rabbit mAb [PSH09-68] - BSA and Azide free	HA751310	Human, Mouse, Rat, Green monkey	WB
	SP1 Recombinant Rabbit mAb [JF0950] - BSA and Azide free	HA750330	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IF-Tissue,IHC-P,FC,ChIP,IP
	SP1 Recombinant Rabbit mAb [JF0950]	ET1702-02	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IF-Tissue,IHC-P,FC,ChIP,IP
KO	Src Recombinant Rabbit mAb [JF0947]	ET1702-03	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KO	Src Recombinant Rabbit mAb [JF0947] - BSA and Azide free	HA750331	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	SSBP1 Recombinant Rabbit mAb [PSH21-31] - BSA and Azide free	HA751794	Human, Mouse, Rat, Green monkey	WB,IHC-P
	SSBP1 Recombinant Rabbit mAb [PSH21-31]	HA724214	Human, Mouse, Rat, Green monkey	WB,IHC-P
	SSRP1 Recombinant Rabbit mAb [PSH09-75]	HA723140	Human, Mouse, Rat, Green monkey	WB,IHC-P
	SSRP1 Recombinant Rabbit mAb [PSH09-75] - BSA and Azide free	HA751317	Human, Mouse, Rat, Green monkey	WB,IHC-P
	STX17 Recombinant Rabbit mAb [PSH23-50]	HA724381	Human, Mouse, Rat, Green monkey	WB,IHC-P,FC(Intra)
	STX17 Recombinant Rabbit mAb [PSH23-50] - BSA and Azide free	HA751898	Human, Mouse, Rat, Green monkey	WB,IHC-P,FC(Intra)
	TAF5L Rabbit pAb	ER65092	Human, Monkey	WB,IHC-P
KD	TAK1 Recombinant Rabbit mAb [JM73-19]	ET1705-14	Human, Mouse, Rat, Green monkey	WB
KD	TAK1 Recombinant Rabbit mAb [JM73-19] - BSA and Azide free	HA750424	Human, Mouse, Rat, Green monkey	WB
	TAPP-1 Recombinant Rabbit mAb [JE64-38]	HA722361	Human, Mouse, Rat, Green monkey	WB,FC
	TATA binding protein TBP Recombinant Rabbit mAb [PSH22-13] - BSA and Azide free	HA751855	Human, Mouse, Rat, Green monkey	WB,IF-Cell,FC,IP
	TATA binding protein TBP Recombinant Rabbit mAb [PSH22-13]	HA724274	Human, Mouse, Rat, Green monkey	WB,IF-Cell,FC,IP
	TAZ / WWTR1 Recombinant Rabbit mAb [PSH15-57] - BSA and Azide free	HA751568	Human, Mouse, Rat, Green monkey	WB,IP
	TAZ / WWTR1 Recombinant Rabbit mAb [PSH15-57]	HA723764	Human, Mouse, Rat, Green monkey	WB,IP
	TCPTP Recombinant Rabbit mAb [PSH24-14] - BSA and Azide free	HA751928	Human, Mouse, Rat, Green monkey	WB,IHC-P
	TCPTP Recombinant Rabbit mAb [PSH24-14]	HA724425	Human, Mouse, Rat, Green monkey	WB,IHC-P
Relative expression	TGF beta 1 Recombinant Rabbit mAb [PD00-17]	HA721143	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Tissue
	TGF beta 1 Recombinant Rabbit mAb [PSH18-59] - BSA and Azide free	HA751683	Human, Mouse, Green monkey	WB,IF-Cell,IHC-P,FC
Relative expression	TGF beta 1 Recombinant Rabbit mAb [PD00-17] - BSA and Azide free	HA750485	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Tissue
	TGF beta 1 Recombinant Rabbit mAb [PSH18-59]	HA724012	Human, Mouse, Green monkey	WB,IF-Cell,IHC-P,FC
	TGF beta 1 Recombinant Rabbit mAb [PSH15-22] - BSA and Azide free (Detector)	HA723726	Human, Mouse, Rat, Dog, Cynomolgus monkey, Pig	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	TGF beta 1 Recombinant Rabbit mAb [PSH15-21] - BSA and Azide free (Capture)	HA723725	Human, Mouse, Rat, Dog, Cynomolgus monkey, Pig	ELISA(Cap)
	Biotin Conjugated TGF beta 1 Recombinant Rabbit mAb [PSH15-22] - Detector	HA723727B	Human, Mouse, Rat, Dog, Cynomolgus monkey, Pig	ELISA(Det),ELISA
	THAP4 Rabbit pAb	ER65158	Human, Mouse, Rat, Monkey	WB,IHC-P
	Thioredoxin Recombinant Rabbit mAb [PSH23-52]	HA724383	Human, Mouse, Rat, Green monkey	WB,IP
	Thioredoxin Recombinant Rabbit mAb [PSH23-52] - BSA and Azide free	HA751900	Human, Mouse, Rat, Green monkey	WB,IP
	TIMP1 Recombinant Rabbit mAb [PSH21-64]	HA724231	Human, Green monkey	WB,IHC-P
	TIMP1 Recombinant Rabbit mAb [PSH21-64] - BSA and Azide free	HA751803	Human, Green monkey	WB,IHC-P
Relative expression	TIMP2 Recombinant Rabbit mAb [PSH11-62]	HA723358	Human, Green monkey	WB
Relative expression	TIMP2 Recombinant Rabbit mAb [PSH11-62] - BSA and Azide free	HA751418	Human, Green monkey	WB
	TIMP4 Recombinant Rabbit mAb [PSH20-91] - BSA and Azide free	HA751773	Human, Mouse, Rat, Green monkey	WB,FC
	TIMP4 Recombinant Rabbit mAb [PSH20-91]	HA724186	Human, Mouse, Rat, Green monkey	WB,FC
	TNFRSF10D Mouse mAb [7G7H9]	EM1711-12	Human, Monkey	WB
	TNFSF11 Mouse mAb [8A7B9]	EM1710-42	Human, Monkey	WB,IHC-P,IF-Cell,FC
	TOM70 Recombinant Rabbit mAb [PSH13-00]	HA723505	Human, Mouse, Rat, Green monkey	WB,IHC-P,FC,IF-Cell,IP
	TOM70 Recombinant Rabbit mAb [PSH13-00] - BSA and Azide free	HA751455	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,IP,FC(Intra)
	TRAF6 Recombinant Rabbit mAb [JE01-42]	HA722857	Human, Mouse, Rat, Green monkey	WB,IF-Cell,FC
	Transketolase Recombinant Rabbit mAb [PSH13-70] - BSA and Azide free	HA751483	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC
	Transketolase Recombinant Rabbit mAb [PSH13-70]	HA723559	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC
	Transketolase Rabbit pAb	ER62380	Human, Mouse, Rat, Green monkey	WB
	TRH-R1 Rabbit pAb	ER65215	Human, Mouse, Rat, Monkey	WB,IF-Cell
KD	TRIM56 Recombinant Rabbit mAb [PSH02-17]	HA721776	Human, Mouse, Rat, Green monkey	WB,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	TRIM56 Recombinant Rabbit mAb [PSH02-17] - BSA and Azide free	HA750787	Human, Mouse, Rat, Green monkey	WB,FC
	Triosephosphate isomerase Recombinant Rabbit mAb [PSH14-38]	HA723632	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,FC
	Triosephosphate isomerase Recombinant Rabbit mAb [PSH14-38] - BSA and Azide free	HA751520	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,FC
	Triosephosphate isomerase Rabbit pAb	HA500283	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,FC
KD	TRPA1 Recombinant Rabbit mAb [PSH23-56]	HA724319	Human, Green monkey	WB
KD	TRPA1 Recombinant Rabbit mAb [PSH23-56] - BSA and Azide free	HA751875	Human, Green monkey	WB
	TYK2 Recombinant Rabbit mAb [PSH18-53] - BSA and Azide free	HA751677	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,FC(Intra)
	TYK2 Recombinant Rabbit mAb [PSH18-53]	HA724006	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,FC(Intra)
	UAP56 Recombinant Rabbit mAb [JE34-73]	HA722726	Human, Mouse, Rat, Green monkey	WB,IHC-P,IF-Cell,FC,IF-Tissue
	UCP2 Recombinant Rabbit mAb [PSH19-75]	HA724099	Human, Mouse, Rat, Green monkey	WB
	UCP2 Recombinant Rabbit mAb [PSH19-75] - BSA and Azide free	HA751730	Human, Mouse, Rat, Green monkey	WB
	UHRF1 Recombinant Rabbit mAb [PSH10-63]	HA723228	Human, Mouse, Green monkey	WB,IF-Cell,FC,IP,ChIP
	UHRF1 Recombinant Rabbit mAb [PSH10-63] - BSA and Azide free	HA751366	Human, Mouse, Green monkey	WB,IF-Cell,FC,IP,ChIP
	UQCRB Recombinant Rabbit mAb [PSH19-74] - BSA and Azide free	HA751729	Human, Mouse, Rat, Green monkey	WB,IHC-P,IP
	UQCRB Recombinant Rabbit mAb [PSH19-74]	HA724098	Human, Mouse, Rat, Green monkey	WB,IHC-P,IP
	USP15 Recombinant Rabbit mAb [PSH15-94]	HA723799	Human, Mouse, Rat, Green monkey	WB,IHC-P,IP
	USP15 Recombinant Rabbit mAb [PSH15-94] - BSA and Azide free	HA751585	Human, Mouse, Rat, Green monkey	WB,IHC-P,IP
	USP15 Rabbit pAb	HA500497	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC
	USP25 Recombinant Rabbit mAb [PSH25-91] - BSA and Azide free	HA752003	Human, Mouse, Rat, Green monkey	WB,IF-Cell,FC(Intra),IP
	USP25 Recombinant Rabbit mAb [PSH25-91]	HA724583	Human, Mouse, Rat, Green monkey	WB,IF-Cell,FC(Intra),IP
	VCP Recombinant Rabbit mAb [JM11-15]	ET1703-56	Human, Mouse, Rat, Zebrafish, Green monkey	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	VEGFA Recombinant Rabbit mAb [PSH19-17]	HA724051	Human, Mouse, Rat, Green monkey	WB
	VEGFA Recombinant Rabbit mAb [PSH19-17] - BSA and Azide free	HA751696	Human, Mouse, Rat, Green monkey	WB
KD	Vinculin Recombinant Rabbit mAb [JM42-43] - BSA and Azide free	HA750448	Human, Mouse, Rat, Green monkey, Pig	WB,IF-Cell,IF-Tissue,IHC-P
KD	Vinculin Recombinant Rabbit mAb [JM42-43]	ET1705-94	Human, Mouse, Rat, Green monkey, Pig	WB,IF-Cell,IF-Tissue,IHC-P
	HRP Conjugated Vinculin Recombinant Rabbit mAb [JM42-43]	HA724470H	Human, Mouse, Rat, Green monkey, Pig	WB
	VIRMA Recombinant Rabbit mAb [PSH08-28] - BSA and Azide free	HA751219	Human, Rat, Green monkey	WB,IF-Cell,IHC-P,FC,IP
	VIRMA Recombinant Rabbit mAb [PSH08-28]	HA722972	Human, Rat, Green monkey	WB,IF-Cell,IHC-P,FC,IP
	VPAC1 Recombinant Rabbit mAb [PD00-96]	HA721209	Human, Mouse, Rat, Green monkey	WB,IHC-P
KD	VPS35 Recombinant Rabbit mAb [JB33-82]	ET7107-01	Human, Mouse, Rat, Green monkey, Cynomolgus monkey, Pig	WB,IHC-P,IF-Cell
KD	VPS35 Recombinant Rabbit mAb [PSH16-54] - BSA and Azide free	HA751603	Human, Mouse, Rat, Green monkey	WB,IF-Cell,FC,IHC-P,IP
KD	VPS35 Recombinant Rabbit mAb [PSH16-54]	HA723858	Human, Mouse, Rat, Green monkey	WB,IF-Cell,FC,IHC-P,IP
	Wee1 Recombinant Rabbit mAb [PSH05-11]	HA722227	Human, Green monkey	WB,IHC-P
	Wee1 Recombinant Rabbit mAb [PSH05-11] - BSA and Azide free	HA750984	Human, Green monkey	WB,IHC-P
	WTAP Mouse mAb [B6-B6]	EM1706-68	Human, Monkey	WB,IF-Cell,IHC-P,FC
	WWOX Rabbit pAb	HA500592	Human, Mouse, Rat, Green monkey	WB,IF-Cell,FC
	WWOX Recombinant Rabbit mAb [PSH15-43] - BSA and Azide free	HA751558	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IP
	WWOX Recombinant Rabbit mAb [PSH15-43]	HA723749	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IP
	WWOX Recombinant Rabbit mAb [PSH15-44]	HA723750	Human, Mouse, Rat, Green monkey	WB,IP
	WWOX Recombinant Rabbit mAb [PSH15-44] - BSA and Azide free	HA751559	Human, Mouse, Rat, Green monkey	WB,IP
KD	XIAP Recombinant Rabbit mAb [PSH04-30] - BSA and Azide free	HA750939	Human, Mouse, Rat, Green monkey	WB,IHC-P
KD	XIAP Recombinant Rabbit mAb [PSH04-30]	HA722113	Human, Mouse, Rat, Green monkey	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	XPD Recombinant Rabbit mAb [JE55-02]	ET7110-93	Human, Mouse, Green monkey	WB,IHC-P
	YTHDF2 Recombinant Rabbit mAb [PSH06-17]	HA722514	Human, Mouse, Rat, Green monkey	WB,IHC-P,FC
	YTHDF2 Recombinant Rabbit mAb [PSH06-17] - BSA and Azide free	HA751055	Human, Mouse, Rat, Green monkey	WB,IHC-P,FC
	ZDHHC2 Rabbit pAb	ER65333	Human, Mouse, Rat, Monkey	WB
	ZFP106 Rabbit pAb	ER65335	Human, Mouse, Monkey	WB
Relative expression	ZHX2 Recombinant Rabbit mAb [JE63-11]	HA721886	Human, Mouse, Rat, Green monkey	WB,IF-Cell
	ZO1 Recombinant Rabbit mAb [PSH07-09]	HA722797	Human, Mouse, Green monkey, Cynomolgus monkey	WB,IF-Cell,mIHC
	ZO1 Recombinant Rabbit mAb [PSH21-37]	HA722851	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P
	ZO1 Recombinant Rabbit mAb - BSA and Azide free	HA751628	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P
	ZO1 Recombinant Rabbit mAb [PSH07-09] - BSA and Azide free	HA751131	Human, Mouse, Green monkey, Cynomolgus monkey	WB,IF-Cell
	ZO3 Rabbit pAb	ER61305	Human, Mouse, Monkey	WB,IF-Cell
	ZW10 Recombinant Rabbit mAb [PSH05-80]	HA722479	Human, Mouse, Rat, Green monkey	WB,IHC-P
	ZW10 Recombinant Rabbit mAb [PSH05-80] - BSA and Azide free	HA751028	Human, Mouse, Rat, Green monkey	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层