

2026年7月—新上市产品清单

高级验证	产品名	货号	反应种属	验证应用
KD	ADAM10 Recombinant Rabbit mAb [JM32-11] - BSA and Azide free	HA752014	Human, Mouse, Rat	WB,IHC-P,IP
	ADX Recombinant Rabbit mAb [JE63-56] - BSA and Azide free	HA752023	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Alpha-Synuclein Recombinant Rabbit mAb [JB42-33] - BSA and Azide free	HA752018	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Beta III Tubulin Recombinant Antibody [PSH15-32] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610444	Human, Mouse, Rat	IHC-P,IF-Tissue,IHC-Fr,WB,IF-Cell
	Beta III Tubulin Recombinant Antibody [PSH15-32] - Rat IgG1 (Chimeric)	HA601655	Human, Mouse, Rat	IHC-P,IF-Tissue,IHC-Fr,WB,IF-Cell
	Beta III Tubulin Recombinant Antibody [PSH15-32] - Mouse IgG1 (Chimeric)	HA601654	Human, Mouse, Rat	IHC-P,IF-Tissue,IHC-Fr,WB,IF-Cell
	Beta III Tubulin Recombinant Antibody [PSH15-32] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610445	Human, Mouse, Rat	IHC-P,IF-Tissue,IHC-Fr,WB,IF-Cell
Cell treatment	BNIP3 Recombinant Rabbit mAb [JA71-10] - BSA and Azide free	HA752015	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	BSA Recombinant Rabbit mAb [PSH26-82] - BSA and Azide free	HA724647	Species independent	ELISA
Cell treatment	CCL1 Recombinant Rabbit mAb [PSH26-39] - BSA and Azide free	HA752032	Human	WB
Cell treatment	CCL1 Recombinant Rabbit mAb [PSH26-39]	HA724618	Human	WB
KD	CD127 Recombinant Rabbit mAb [JE37-10] - BSA and Azide free	HA752021	Human, Mouse	WB,IF-Cell,IHC-P,FC
Relative expression	CD133 Recombinant Rabbit mAb [PSH26-71] - BSA and Azide free	HA752039	Human	WB,IHC-P,FC
Relative expression	CD133 Recombinant Rabbit mAb [PSH26-71]	HA724636	Human	WB,IHC-P,FC
	CDC20 Recombinant Rabbit mAb [PSH26-87] - BSA and Azide free	HA752049	Human, Green monkey	WB
	CDC20 Recombinant Rabbit mAb [PSH26-87]	HA724652	Human, Green monkey	WB
Relative expression	CDCP1 Recombinant Rabbit mAb [PSH26-74]	HA724639	Human, Mouse, Rat	WB,IHC-P,FC
Relative expression	CDCP1 Recombinant Rabbit mAb [PSH26-74] - BSA and Azide free	HA752042	Human, Mouse, Rat	WB,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Cell-Surface Vimentin (CSV) Recombinant Mouse mAb [PSH26-07]	HA601653	Human	IF-Cell,FC
	Cell-Surface Vimentin (CSV) Recombinant Mouse mAb [PSH26-07] - BSA and Azide free	HA610443	Human	IF-Cell,FC
Relative expression	CIP2A Recombinant Rabbit mAb [PSH16-11] - BSA and Azide free	HA752006	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Relative expression	CIP2A Recombinant Rabbit mAb [PSH16-11]	HA723816	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	COL1A1 Recombinant Antibody [PSH06-20] - Rat IgG1 (Chimeric)	HA601656	Human	mIHC
	Biotin Conjugated CREB1 Recombinant Rabbit mAb [PSH26-25] - Detector	HA724611B	Human, Mouse, Rat	ELISA(Det),ELISA
	CREB1 Recombinant Rabbit mAb [PSH26-25] - BSA and Azide free (Detector)	HA724611	Human, Mouse, Rat	ELISA(Det)
	CREB1 Recombinant Rabbit mAb [SA04-04] - BSA and Azide free (Capture)	HA724610	Human, Mouse, Rat	ELISA(Cap)
	CTHRC1 Recombinant Mouse mAb [PSH26-63] - BSA and Azide free	HA610446	Human	ELISA
	CTHRC1 Recombinant Mouse mAb [PSH26-63]	HA601657	Human	ELISA
	Cullin 5 / CUL-5 Recombinant Rabbit mAb [PSH25-83] - BSA and Azide free	HA751995	Human, Mouse, Rat	WB,IHC-P,IP
	Cullin 5 / CUL-5 Recombinant Rabbit mAb [PSH25-83]	HA724575	Human, Mouse, Rat	WB,IHC-P,IP
Relative expression, Cell treatment	CX3CL1 Recombinant Rabbit mAb [PSH26-38] - BSA and Azide free	HA752031	Mouse, Rat, Pig	WB,FC,IP
Cell treatment, Relative expression	CX3CL1 Recombinant Rabbit mAb [PSH26-38]	HA724617	Mouse, Rat, Pig	WB,FC,IP
Relative expression	CX3CL1 Recombinant Rabbit mAb [PSH26-37] - BSA and Azide free	HA752030	Mouse, Rat	WB,IHC-Fr,IHC-P,IP
Relative expression	CX3CL1 Recombinant Rabbit mAb [PSH26-37]	HA724616	Mouse, Rat	WB,IHC-Fr,IHC-P,IP
	HRP Conjugated DYKDDDDK Tag (FLAG) Recombinant Mouse mAb [A2- A4-R]	HA601652H	Species independent	WB
	Dystrophin Recombinant Rabbit mAb [PSH26-84] - BSA and Azide free (Detector)	HA724649	Human, Mouse, Rat	ELISA(Det)
	Dystrophin Recombinant Rabbit mAb [PSH26-83] - BSA and Azide free (Capture)	HA724648	Human, Mouse, Rat	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Dystrophin Recombinant Rabbit mAb [PSH26-84] - Detector	HA724649B	Human, Mouse, Rat	ELISA(Det),ELISA
Cell treatment	EGR1 Recombinant Rabbit mAb [JE57-62] - BSA and Azide free	HA752019	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P
	Estrogen Receptor Recombinant Antibody [PD00-04] - Mouse IgG1 (Chimeric)	HA601648	Human	IHC-P
	iFluor™ 647 Conjugated F4/80 Recombinant Rabbit mAb [PSH01-87]	HA724595F	Mouse	IHC-Fr,IF-Tissue
	Formyl-Histone H2B type 2E (K120) Recombinant Rabbit mAb [PSH26-76]	HA724641	Human, Mouse, Rat	WB,IHC-P
	Histone H3 (tri methyl K4) Recombinant Rabbit mAb [PSH25-68] - BSA and Azide free	HA751987	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC(Intra),ChIP
	Histone H3 (tri methyl K4) Recombinant Rabbit mAb [PSH25-68]	HA724555	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC(Intra),ChIP
Cell treatment	Histone H4 (acetyl K12) Recombinant Rabbit mAb [PSH25-67] - BSA and Azide free	HA751986	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,ChIP
Cell treatment	Histone H4 (acetyl K12) Recombinant Rabbit mAb [PSH25-67]	HA724554	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,ChIP
Relative expression	HLA-DR Recombinant Rabbit mAb [PSH26-51] - BSA and Azide free	HA752035	Human	WB,IF-Cell,IHC-P
Relative expression	HLA-DR Recombinant Rabbit mAb [PSH26-51]	HA724628	Human	WB,IF-Cell,IHC-P
Relative expression	Homer1 Recombinant Rabbit mAb [JE65-19] - BSA and Azide free	HA752020	Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr
	Human Acetylcholinesterase/ACHE enzyme-linked immunoassay kit	EH0226	Human	ELISA
	Human Amphiregulin/ARGE enzyme-linked immunoassay kit	EH0273	Human	ELISA
	Human APOC3/Apolipoprotein CIII enzyme-linked immunoassay kit(one step)	EHY0260	Human	ELISA
	Human Apolipoprotein E/Apo E enzyme-linked immunoassay kit	EH0237	Human	ELISA
	Human Apolipoprotein E/Apo E enzyme-linked immunoassay kit(one step)	EHY0237	Human	ELISA
	Human BDNF, Tag Free Protein	HA211599	Human	
	Human CCL18/PARC/MIP-4 enzyme-linked immunoassay kit	EH0271	Human	ELISA
	Human CD127/IL-7 R alpha enzyme-linked immunoassay kit	EH0259	Human	ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human CD30 Ligand/TNFSF8/CD 30-L enzyme-linked immunoassay kit	EH0274	Human	ELISA
	Human CD4, Tag Free Protein	HA211593	Human	
	PE Conjugated Human CD63 Recombinant Rabbit mAb [PSH06-64]	HA724635F	Human	FC
	Human CD70 Recombinant Rabbit mAb [PSH26-00] - BSA and Azide free (Capture)	HA724590	Human	ELISA(Cap)
	Biotin Conjugated Human CD70 Recombinant Rabbit mAb [PSH26-01] - Detector	HA724591B	Human	ELISA(Det),ELISA
	Human CD70 Recombinant Rabbit mAb [PSH26-01] - BSA and Azide free (Detector)	HA724591	Human	ELISA(Det)
	Human CD70 Recombinant Rabbit mAb [PSH26-02] - BSA and Azide free (Capture)	HA724592	Human	ELISA(Cap)
	Human Chemerin Recombinant Rabbit mAb [PSH26-16] - BSA and Azide free (Detector)	HA724601	Human	ELISA(Det)
	Biotin Conjugated Human Chemerin Recombinant Rabbit mAb [PSH26-17] - Detector	HA724602B	Human	ELISA(Det),ELISA
	Human Chemerin Recombinant Rabbit mAb [PSH26-17] - BSA and Azide free (Detector)	HA724602	Human	ELISA(Det)
	Biotin Conjugated Human Chemerin Recombinant Rabbit mAb [PSH26-16] - Detector	HA724601B	Human	ELISA(Det),ELISA
	Human Chemerin Recombinant Rabbit mAb [PSH26-15] - BSA and Azide free (Capture)	HA724600	Human	ELISA(Cap)
	Biotin Conjugated Human c-Kit Recombinant Rabbit mAb [PSH26-41] - Detector	HA724620B	Human	ELISA(Det),ELISA
	Human c-Kit Recombinant Rabbit mAb [PSH26-42] - BSA and Azide free (Detector)	HA724621	Human	ELISA(Det)
	Human c-Kit Recombinant Rabbit mAb [PSH26-41] - BSA and Azide free (Detector)	HA724620	Human	ELISA(Det)
	Biotin Conjugated Human c-Kit Recombinant Rabbit mAb [PSH26-42] - Detector	HA724621B	Human	ELISA(Det),ELISA
	Human c-Kit Recombinant Rabbit mAb [PSH26-40] - BSA and Azide free (Capture)	HA724619	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human Complement C2 Recombinant Rabbit mAb [PSH26-57] - BSA and Azide free (Capture)	HA725525	Human	ELISA(Cap)
	Biotin Conjugated Human Complement C2 Recombinant Rabbit mAb [PSH26-58] - Detector	HA725526B	Human	ELISA(Det),ELISA
	Human Complement C2 Recombinant Rabbit mAb [PSH26-58] - BSA and Azide free (Detector)	HA725526	Human	ELISA(Det)
	Human Complement factor D enzyme-linked immunoassay kit(one step)	EHY0242	Human	ELISA
	Human Complement factor H enzyme-linked immunoassay kit	EH0246	Human	ELISA
	Human CXADR/CAR enzyme-linked immunoassay kit(one step)	EHY0262	Human	ELISA
	Human DKK1 enzyme-linked immunoassay kit	EH0233	Human	ELISA
	Human DLL4, Tag Free Protein	HA211588	Human	
	Human DMP1, C-His Tag Protein	HA211515	Human	
	Biotin Conjugated Human DPT Recombinant Rabbit mAb [PSH26-69] - Detector	HA725498B	Human	ELISA(Det),ELISA
	Human DPT Recombinant Rabbit mAb [PSH26-69] - BSA and Azide free (Detector)	HA725498	Human	ELISA(Det)
	Human DPT Recombinant Rabbit mAb [PSH26-68] - BSA and Azide free (Capture)	HA725497	Human	ELISA(Cap)
	Human EGFR enzyme-linked immunoassay kit	EH0254	Human	ELISA
	Human Elafin/PI3 enzyme-linked immunoassay kit	EH0236	Human	ELISA
	Human Endostatin/COL18A1 enzyme-linked immunoassay kit(one step)	EHY0232	Human	ELISA
	Human FABP1/L-FABP enzyme-linked immunoassay kit	EH0257	Human	ELISA
	Human FABP2 enzyme-linked immunoassay kit	EH0256	Human	ELISA
	Human FAP Recombinant Rabbit mAb [PSH26-56] - BSA and Azide free (Detector)	HA725500	Human	ELISA(Det)
	Human FAP Recombinant Rabbit mAb [PSH26-55] - BSA and Azide free (Capture)	HA725499	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human FAP Recombinant Rabbit mAb [PSH26-56] - Detector	HA725500B	Human	ELISA(Det),ELISA
	Human FGF-18, C-Strep II Tag Protein	HA211584	Human	
	Human FGF-21 enzyme-linked immunoassay kit	EH0261	Human	ELISA
	Human FGL1 enzyme-linked immunoassay kit	EH0263	Human	ELISA
	Human Flt3 Ligand/FLT3LG enzyme-linked immunoassay kit	EH0255	Human	ELISA
	Biotin Conjugated Human GC / Vitamin D binding protein Recombinant Rabbit mAb [PSH26-30] - Detector	HA725374B	Human	ELISA(Det),ELISA
	Human GC / Vitamin D binding protein Recombinant Rabbit mAb [PSH26-30] - BSA and Azide free (Detector)	HA725374	Human	ELISA(Det)
	Human GC / Vitamin D binding protein Recombinant Rabbit mAb [PSH26-29] - BSA and Azide free (Capture)	HA725373	Human	ELISA(Cap)
	HRP Conjugated Human GDF15 Recombinant Rabbit mAb [PSH16-02] - Detector	HA723807H	Human	ELISA(Det),ELISA
	HRP Conjugated Human GDF15 Recombinant Rabbit mAb [PSH16-01] - Detector	HA723806H	Human	ELISA(Det),ELISA
	Human GDF-15 enzyme-linked immunoassay kit	EH0267	Human	ELISA
	Human Granzyme B Recombinant Rabbit mAb [PSH26-06] - BSA and Azide free (Detector)	HA725512	Human	ELISA(Det)
	Biotin Conjugated Human Granzyme B Recombinant Rabbit mAb [PSH26-06] - Detector	HA725512B	Human	ELISA(Det),ELISA
	Human Granzyme B Recombinant Rabbit mAb [PSH26-05] - BSA and Azide free (Capture)	HA725511	Human	ELISA(Cap)
	Human Hemoglobin subunit alpha/HBA1 enzyme-linked immunoassay kit	EH0247	Human	ELISA
	Human HIF-1 alpha Recombinant Rabbit mAb [PSH26-18] - BSA and Azide free (Capture/Detector)	HA724603	Human	ELISA(Cap),ELISA(Det)
	Human HIF-1 alpha Recombinant Rabbit mAb [PSH26-19] - BSA and Azide free (Detector)	HA724604	Human	ELISA(Det)
	Human HIF-1 alpha Recombinant Rabbit mAb [PSH26-20] - BSA and Azide free (Capture)	HA724605	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human HIF-1 alpha Recombinant Rabbit mAb [PSH26-19] - Detector	HA724604B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human HIF-1 alpha Recombinant Rabbit mAb [PSH26-18] - Detector	HA724603B	Human	ELISA(Det),ELISA
	Human IGFBP1 enzyme-linked immunoassay kit(one step)	EHY0258	Human	ELISA
	Human IL-10RB enzyme-linked immunoassay kit	EH0223	Human	ELISA
	Human IL-17F enzyme-linked immunoassay kit	EH0249	Human	ELISA
	Human IL-17RA/CD217 enzyme-linked immunoassay kit	EH0278	Human	ELISA
	Human IL-1 alpha enzyme-linked immunoassay kit	EH0243	Human	ELISA
	Human IL-1R2 Recombinant Rabbit mAb [PSH26-09] - BSA and Azide free (Detector)	HA725252	Human	ELISA(Det)
	Biotin Conjugated Human IL-1R2 Recombinant Rabbit mAb [PSH26-09] - Detector	HA725252B	Human	ELISA(Det),ELISA
	Human IL-1R2 Recombinant Rabbit mAb [PSH26-08] - BSA and Azide free (Capture)	HA725251	Human	ELISA(Cap)
	Human IL-23R, C-His Tag Protein	HA211579	Human	
	Human INHBE, N-His Tag Protein	HA211580	Human	
	Human Interferon alpha-2/IFNA2 enzyme-linked immunoassay kit	EHSS0272	Human	ELISA
	Human Interferon omega-1/IFNW1 enzyme-linked immunoassay kit	EH0265	Human	ELISA
	Human KIM-1/HAVCR1/TIM-1 enzyme-linked immunoassay kit	EH0241	Human	ELISA
	Human LAMC2 enzyme-linked immunoassay kit	EH0220	Human	ELISA
	Biotin Conjugated Human LAMP3 Recombinant Rabbit mAb [PSH26-32] - Detector	HA725312B	Human	ELISA(Det),ELISA
	Human LAMP3 Recombinant Rabbit mAb [PSH26-32] - BSA and Azide free (Detector)	HA725312	Human	ELISA(Det)
	Human LAMP3 Recombinant Rabbit mAb [PSH26-31] - BSA and Azide free (Capture)	HA725311	Human	ELISA(Cap)
	Human MAPK1/ERK2, C-His Tag Protein	HA211598	Human	

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human MCP-3/CCL7 enzyme-linked immunoassay kit	EH0269	Human	ELISA
	Human MKI67/Ki67, C-His Tag Protein	HA211578	Human	
	Human/Monkey VEGFR3/FLT4 enzyme-linked immunoassay kit(one step)	EHY0276	Human, Monkey	ELISA
	Human/Mouse IL-20 enzyme-linked immunoassay kit	EH0275	Human, Mouse	ELISA
	Human NLRP3, Tag Free Protein	HA211600	Human	
	Human PARK7/DJ-1 enzyme-linked immunoassay kit(one step)	EHY0230	Human	ELISA
	Biotin Conjugated Human PDGFR alpha Recombinant Rabbit mAb [PSH26-34] - Detector	HA725516B	Human	ELISA(Det),ELISA
	Human PDGFR alpha Recombinant Rabbit mAb [PSH26-33] - BSA and Azide free (Capture)	HA725515	Human	ELISA(Cap)
	Human PDGFR alpha Recombinant Rabbit mAb [PSH26-34] - BSA and Azide free (Detector)	HA725516	Human	ELISA(Det)
	Human PDLIM5 Recombinant Rabbit mAb [PSH26-13] - BSA and Azide free (Capture)	HA724598	Human	ELISA(Cap)
	Biotin Conjugated Human PDLIM5 Recombinant Rabbit mAb [PSH26-14] - Detector	HA724599B	Human	ELISA(Det),ELISA
	Human PDLIM5 Recombinant Rabbit mAb [PSH26-14] - BSA and Azide free (Detector)	HA724599	Human	ELISA(Det)
	Human Prealbumin/Transthyretin enzyme-linked immunoassay kit(one step)	EHY0250	Human	ELISA
	Human PSMA / FOLH1, Tag Free Protein	HA211553	Human	
	Human S100A12 enzyme-linked immunoassay kit(one step)	EHY0229	Human	ELISA
	Human SLAMF1/SLAM/CD150, C-His Tag Protein	HA211554	Human	
	Human TACI/TNFRSF13B enzyme-linked immunoassay kit	EH0264	Human	ELISA
	Human TNFSF15 Recombinant Rabbit mAb [PSH26-47] - BSA and Azide free (Capture)	HA725277	Human	ELISA(Cap)
	Human TNFSF15 Recombinant Rabbit mAb [PSH26-48] - BSA and Azide free (Detector)	HA725278	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human TNFSF15 Recombinant Rabbit mAb [PSH26-48] - Detector	HA725278B	Human	ELISA(Det),ELISA
	Human VEGFR1/VEGF Receptor/FLT1 enzyme-linked immunoassay kit	EH0270	Human	ELISA
	IP3 Receptor Recombinant Rabbit mAb [JA11-35] - BSA and Azide free	HA752016	Human, Mouse, Rat	WB,IHC-P,IHC-Fr,IF-Tissue
	K48-linkage Specific Ubiquitin Recombinant Rabbit mAb [PSH25-69]	HA724556	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC(Intra)
	K48-linkage Specific Ubiquitin Recombinant Rabbit mAb [PSH25-69] - BSA and Azide free	HA751988	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC(Intra)
Relative expression, KD	LEF1 Recombinant Rabbit mAb [JE52-28] - BSA and Azide free	HA752026	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Malachite green Recombinant Mouse mAb [E0-A4-R] - BSA and Azide free	HA601658	Species independent	ELISA
	Melanoma gp100 Recombinant Rabbit mAb [JE42-47] - BSA and Azide free	HA752024	Human, Mouse	WB,IHC-P,IF-Cell
Relative expression, KD	MGMT Recombinant Rabbit mAb [JJ089-6] - BSA and Azide free	HA752010	Human	WB,IHC-P,IF-Cell,FC(Intra)
	Monkey CRP/C-Reactive Protein enzyme-linked immunoassay kit	EH0281	Monkey	ELISA
	Mouse Angiopoietin-2/ANGPT2, C-His Tag Protein	HA211583	Mouse	
	HRP Conjugated Mouse Apolipoprotein CIII Recombinant Rabbit mAb [PSH25-28] - Detector	HA725156H	Mouse	ELISA(Det),ELISA
	Mouse Apolipoprotein E enzyme-linked immunoassay kit(one step)	EMY0065	Mouse	ELISA
	Mouse CD16/32 Recombinant Rabbit mAb [PSH06-60] - BSA and Azide free	HA724562	Mouse	FC
	HRP Conjugated Mouse CD48 Recombinant Rabbit mAb [PSH25-10] - Detector	HA724504H	Mouse	ELISA(Det),ELISA
	Mouse Chil3 / Ym-1 Recombinant Rabbit mAb [PSH26-28] - BSA and Azide free (Detector)	HA724614	Mouse	ELISA(Det)
	Biotin Conjugated Mouse Chil3 / Ym-1 Recombinant Rabbit mAb [PSH26-27] - Detector	HA724613B	Mouse	ELISA(Det),ELISA
	Mouse Chil3 / Ym-1 Recombinant Rabbit mAb [PSH26-27] - BSA and Azide free (Detector)	HA724613	Mouse	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse Chil3 / Ym-1 Recombinant Rabbit mAb [PSH26-26] - BSA and Azide free (Capture)	HA724612	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse Chil3 / Ym-1 Recombinant Rabbit mAb [PSH26-28] - Detector	HA724614B	Mouse	ELISA(Det),ELISA
	Mouse CREB1, C-His Tag Protein	HA211601	Mouse	
	Mouse DLL4, Tag Free Protein	HA211581	Human	
	Mouse GP130/IL16ST enzyme-linked immunoassay kit	EM0076	Mouse	ELISA
	Mouse IL-16 enzyme-linked immunoassay kit	EM0070	Mouse	ELISA
	Mouse IL1RL1/ST2 enzyme-linked immunoassay kit	EM0077	Mouse	ELISA
	Mouse MMP-9 enzyme-linked immunoassay kit	EM0068	Mouse	ELISA
	Mouse Osteopontin/OPN enzyme-linked immunoassay kit(one step)	EMY0073	Mouse	ELISA
	Biotin Conjugated Mouse Osteoprotegerin Recombinant Rabbit mAb [PSH26-46] - Detector	HA724625B	Mouse	ELISA(Det),ELISA
	Mouse Osteoprotegerin Recombinant Rabbit mAb [PSH26-45] - BSA and Azide free (Capture)	HA724624	Mouse	ELISA(Cap)
	Mouse Osteoprotegerin Recombinant Rabbit mAb [PSH26-46] - BSA and Azide free (Detector)	HA724625	Mouse	ELISA(Det)
	Mouse Periostin/POSTN enzyme-linked immunoassay kit	EM0075	Mouse	ELISA
	Biotin Conjugated Mouse/Rat Angiopoietin 2 Recombinant Rabbit mAb [PSH25-97] - Detector	HA724587B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat Angiopoietin 2 Recombinant Rabbit mAb [PSH25-97] - BSA and Azide free (Detector)	HA724587	Mouse, Rat	ELISA(Det)
	Mouse/Rat Angiopoietin 2 Recombinant Rabbit mAb [PSH25-96] - BSA and Azide free (Capture)	HA724586	Mouse, Rat	ELISA(Cap)
	Mouse/Rat B7-H4 Recombinant Rabbit mAb [PSH26-11] - BSA and Azide free (Capture)	HA724596	Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Mouse/Rat B7-H4 Recombinant Rabbit mAb [PSH26-12] - Detector	HA724597B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat B7-H4 Recombinant Rabbit mAb [PSH26-12] - BSA and Azide free (Detector)	HA724597	Mouse, Rat	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Mouse/Rat EGF Recombinant Rabbit mAb [PSH26-62] - Detector	HA725504B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat EGF Recombinant Rabbit mAb [PSH26-62] - BSA and Azide free (Detector)	HA725504	Mouse, Rat	ELISA(Det)
	Mouse/Rat EGF Recombinant Rabbit mAb [PSH26-61] - BSA and Azide free (Capture)	HA725503	Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Mouse/Rat FSH beta Recombinant Rabbit mAb [PSH26-86] - Detector	HA724651B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat FSH beta Recombinant Rabbit mAb [PSH26-86] - BSA and Azide free (Detector)	HA724651	Mouse, Rat	ELISA(Det)
	Mouse/Rat FSH beta Recombinant Rabbit mAb [PSH26-85] - BSA and Azide free (Capture)	HA724650	Mouse, Rat	ELISA(Cap)
	Mouse/Rat HGF enzyme-linked immunoassay kit	EM0072	Mouse, Rat	ELISA
	Mouse/Rat IGFBP1 Recombinant Rabbit mAb [PSH26-21] - BSA and Azide free (Capture)	HA724606	Mouse, Rat	ELISA(Cap)
	Mouse/Rat IGFBP1 Recombinant Rabbit mAb [PSH26-22] - BSA and Azide free (Detector)	HA724607	Mouse, Rat	ELISA(Det)
	Biotin Conjugated Mouse/Rat IGFBP1 Recombinant Rabbit mAb [PSH26-22] - Detector	HA724607B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat PDGF-BB enzyme-linked immunoassay kit	EM0062	Mouse, Rat	ELISA
	Mouse/Rat VEGFA/VEGF164 enzyme-linked immunoassay kit	EM0071	Mouse	ELISA
	Mouse/Rat VEGFR3/FLT4 enzyme-linked immunoassay kit(one step)	EMY0074	Mouse, Rat	ELISA
	MUC5B Recombinant Rabbit mAb [PSH26-66]	HA724632	Mouse, Rat	IHC-P,IF-Tissue
	MUC5B Recombinant Rabbit mAb [PSH26-66] - BSA and Azide free	HA752037	Mouse, Rat	IHC-P,IF-Tissue
	Multi-Tag Protein	HA211568		
Relative expression	MYH2 Recombinant Rabbit mAb [JE53-30] - BSA and Azide free	HA752025	Human, Mouse, Rat	WB,IHC-P
KD	NRG1 Recombinant Rabbit mAb [PSH26-78]	HA724643	Human, Mouse, Rat, Pig	WB,IHC-P,IP
KD	NRG1 Recombinant Rabbit mAb [PSH26-78] - BSA and Azide free	HA752045	Human, Mouse, Rat, Pig	WB,IHC-P,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	Occludin Recombinant Rabbit mAb [JJ091-08] - BSA and Azide free	HA752011	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IP,FC,IHC-P,IHC-Fr,IF-Tissue,mIHC
Cell treatment	O-Linked N-Acetylglucosamine (O-GlcNAc) Recombinant Rabbit mAb [PSH25-76] - BSA and Azide free	HA751994	Human, Mouse, Rat, Green monkey, Pig	WB,IF-Cell,IHC-P,FC(Intra),ELISA
Cell treatment	O-Linked N-Acetylglucosamine (O-GlcNAc) Recombinant Rabbit mAb [PSH25-75] - BSA and Azide free	HA751993	Human, Mouse, Rat, Green monkey, Pig	WB,IF-Cell,IHC-P,FC(Intra),IP,ELISA
Cell treatment	O-Linked N-Acetylglucosamine (O-GlcNAc) Recombinant Rabbit mAb [PSH25-76]	HA724568	Human, Mouse, Rat, Green monkey, Pig	WB,IF-Cell,IHC-P,FC(Intra),ELISA
Cell treatment	O-Linked N-Acetylglucosamine (O-GlcNAc) Recombinant Rabbit mAb [PSH25-75]	HA724567	Human, Mouse, Rat, Green monkey, Pig	WB,IF-Cell,IHC-P,FC(Intra),IP,ELISA
Relative expression	Osteoprotegerin Recombinant Rabbit mAb [PSH26-52]	HA724629	Human, Mouse, Rat	WB,IHC-P,FC(Intra),IHC-Fr
Relative expression	Osteoprotegerin Recombinant Rabbit mAb [PSH26-52] - BSA and Azide free	HA752036	Human, Mouse, Rat	WB,IHC-P,FC(Intra),IHC-Fr
Relative expression	OXCT1 Recombinant Rabbit mAb [PSH26-75]	HA724640	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC(Intra)
Relative expression	OXCT1 Recombinant Rabbit mAb [PSH26-75] - BSA and Azide free	HA752043	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC(Intra)
KD	PARK7/DJ1 Recombinant Rabbit mAb [SN07-21] - BSA and Azide free	HA752009	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Pax2 Recombinant Rabbit mAb [PSH26-04]	HA724593	Human, Mouse, Rat	IHC-P
	PGC1 alpha+beta Recombinant Rabbit mAb [JF09-71] - BSA and Azide free	HA752013	Human, Mouse, Rat	WB
	Phospho-alpha Synuclein (S129) Recombinant Rabbit mAb [JB22-44] - BSA and Azide free	HA752017	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC(Intra),IF-Tissue
Cell treatment	Phospho-Aurora A (T288) Recombinant Rabbit mAb [PSH26-72] - BSA and Azide free	HA752040	Human	WB
Cell treatment	Phospho-Aurora A (T288) Recombinant Rabbit mAb [PSH26-72]	HA724637	Human	WB
Cell treatment	Phospho-MEK1/2 (S218 + S222) Recombinant Rabbit mAb [ST0490] - BSA and Azide free	HA752008	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
Cell treatment	Phospho-Smad1/5/9 (S463/S465/S467) Recombinant Rabbit mAb [JE59-46] - BSA and Azide free	HA752027	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Progesterone Receptor Recombinant Antibody [JF0549] - Mouse IgG1 (Chimeric)	HA601647	Human	IF-Cell,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	Pyruvate Dehydrogenase E2 Recombinant Rabbit mAb [JE34-38] - BSA and Azide free	HA752022	Human, Mouse, Rat	WB,IHC-P
	RAB35 Recombinant Rabbit mAb [PSH26-50] - BSA and Azide free	HA752034	Human, Mouse, Rat, Green monkey, Pig	WB,IF-Cell,IHC- P,FC(Intra),IP
	RAB35 Recombinant Rabbit mAb [PSH26-50]	HA724627	Human, Mouse, Rat, Green monkey, Pig	WB,IF-Cell,IHC- P,FC(Intra),IP
Relative expression	Rabbit IgG heavy chain(Fc fragment) Recombinant Mouse mAb [PSH23-92] - BSA and Azide free	HA601651	Rabbit	WB,ELISA
	Rabbit IL-10, Tag Free Protein	HA211561	Oryctolagus cuniculus	
	Rat CCL2/MCP-1 enzyme-linked immunoassay kit	ER0016	Rat	ELISA
	Rat CXCL2/MIP-2 enzyme-linked immunoassay kit	ER0013	Rat	ELISA
	Rat ICAM-1 enzyme-linked immunoassay kit(one step)	ERY0011	Rat	ELISA
	Rat IL-3 Recombinant Rabbit mAb [PSH26-53] - BSA and Azide free (Capture)	HA725437	Rat	ELISA(Cap)
	Rat IL-3 Recombinant Rabbit mAb [PSH26-54] - BSA and Azide free (Detector)	HA725438	Rat	ELISA(Det)
	Biotin Conjugated Rat IL-3 Recombinant Rabbit mAb [PSH26-54] - Detector	HA725438B	Rat	ELISA(Det),ELISA
	Rat IL-5 enzyme-linked immunoassay kit	ER0014	Rat	ELISA
	Rat KIM-1/HAVCR1 enzyme-linked immunoassay kit	ER0012	Rat	ELISA
	Rat LIF Recombinant Rabbit mAb [PSH26-24] - BSA and Azide free (Detector)	HA724609	Rat	ELISA(Det)
	Rat LIF Recombinant Rabbit mAb [PSH26-23] - BSA and Azide free (Capture)	HA724608	Rat	ELISA(Cap)
	Biotin Conjugated Rat LIF Recombinant Rabbit mAb [PSH26-24] - Detector	HA724609B	Rat	ELISA(Det),ELISA
	Rat/Mouse IL-18 enzyme-linked immunoassay kit	ER0015	Rat, Mouse	ELISA
	Biotin Conjugated Rat Renin Recombinant Rabbit mAb [PSH26-44] - Detector	HA724623B	Rat	ELISA(Det),ELISA
	Rat Renin Recombinant Rabbit mAb [PSH26-44] - BSA and Azide free (Detector)	HA724623	Rat	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Rat Renin Recombinant Rabbit mAb [PSH26-43] - BSA and Azide free (Capture)	HA724622	Rat	ELISA(Cap)
	Rat SOD1 Recombinant Rabbit mAb [PSH26-64] - BSA and Azide free (Capture)	HA724630	Rat	ELISA(Cap)
	Biotin Conjugated Rat SOD1 Recombinant Rabbit mAb [PSH26-65] - Detector	HA724631B	Rat	ELISA(Det),ELISA
	Rat SOD1 Recombinant Rabbit mAb [PSH26-65] - BSA and Azide free (Detector)	HA724631	Rat	ELISA(Det)
	Rat TIMP1 enzyme-linked immunoassay kit	ER0017	Rat	ELISA
Relative expression	RHOB Recombinant Rabbit mAb [PSH26-73] - BSA and Azide free	HA752041	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC(Intra),IP
Relative expression	RHOB Recombinant Rabbit mAb [PSH26-73]	HA724638	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC(Intra),IP
	SDHB Recombinant Rabbit mAb [JE77-70] - BSA and Azide free	HA752028	Human, Mouse, Rat	WB,IHC-P,IP,IF-Cell
	SFRP1 Recombinant Rabbit mAb [PSH26-59] - BSA and Azide free (Capture)	HA725173	Human, Mouse, Rat	ELISA(Cap)
	SFRP1 Recombinant Rabbit mAb [PSH26-60] - BSA and Azide free (Detector)	HA725174	Human, Mouse, Rat	ELISA(Det)
	Biotin Conjugated SFRP1 Recombinant Rabbit mAb [PSH26-60] - Detector	HA725174B	Human, Mouse, Rat	ELISA(Det),ELISA
Cell treatment	SOD1 Recombinant Rabbit mAb [PSH26-10]	HA724615	Human, Mouse, Rat, Pig	WB,IF-Cell,FC(Intra)
Cell treatment	SOD1 Recombinant Rabbit mAb [PSH26-10] - BSA and Azide free	HA752029	Human, Mouse, Rat, Pig	WB,IF-Cell,FC(Intra)
	SOX9 Recombinant Antibody [SN74-09] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610441	Human, Mouse, Rat	IHC-P
	SOX9 Recombinant Antibody [SN74-09] - Rat IgG1 (Chimeric)	HA601650	Human, Mouse, Rat	IHC-P
	SOX9 Recombinant Antibody [SN74-09] - Mouse IgG1 (Chimeric)	HA601649	Human, Mouse, Rat	IHC-P
	SOX9 Recombinant Antibody [SN74-09] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610442	Human, Mouse, Rat	IHC-P
Relative expression	STAT4 Recombinant Rabbit mAb [PSH26-81]	HA724646	Human, Mouse, Rat	WB,IHC-P,IP
Relative expression	STAT4 Recombinant Rabbit mAb [PSH26-81] - BSA and Azide free	HA752048	Human, Mouse, Rat	WB,IHC-P,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Superoxide Dismutase 1 Recombinant Rabbit mAb [JF1005] - BSA and Azide free	HA752012	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IHC-P,FC
	Tenascin C Recombinant Rabbit mAb [SU36-01] - BSA and Azide free	HA752007	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue
Relative expression	TET2 Recombinant Rabbit mAb [PSH26-49]	HA724626	Human, Rat	WB,IHC-P,IP
Relative expression	TET2 Recombinant Rabbit mAb [PSH26-49] - BSA and Azide free	HA752033	Human, Rat	WB,IHC-P,IP
	UAP1 Recombinant Mouse mAb [4H1-R] - BSA and Azide free	HA601660	Human	WB
	VP40 Recombinant Mouse mAb [A2-G7-R] - BSA and Azide free	HA601659	Species independent	ELISA
Relative expression	XPC Recombinant Rabbit mAb [PSH26-70] - BSA and Azide free	HA752038	Human	WB,IF-Cell,IHC-P,IP
Relative expression	XPC Recombinant Rabbit mAb [PSH26-70]	HA724634	Human	WB,IF-Cell,IHC-P,IP
Relative expression	YTHDC2 Recombinant Rabbit mAb [PSH26-80] - BSA and Azide free	HA752047	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC(Intra),IP
Relative expression	YTHDC2 Recombinant Rabbit mAb [PSH26-80]	HA724645	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC(Intra),IP
	ZO1 Recombinant Rabbit mAb [PSH26-36] - BSA and Azide free (Detector)	HA725518	Human, Mouse, Rat	ELISA(Det)
	Biotin Conjugated ZO1 Recombinant Rabbit mAb [PSH26-36] - Detector	HA725518B	Human, Mouse, Rat	ELISA(Det),ELISA
	ZO1 Recombinant Rabbit mAb [PSH26-35] - BSA and Azide free (Capture)	HA725517	Human, Mouse, Rat	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层