

Ras信号通路研究-产品清单

高级验证	产品名	货号	反应种属	验证应用
	ABL1 Recombinant Rabbit mAb [PSH04-08]	HA722089	Human, Mouse	WB,IF-Cell,IHC-P,FC
	ABL1 Recombinant Rabbit mAb [PSH04-08] - BSA and Azide free	HA750931	Human, Mouse	WB,IF-Cell,IHC-P,FC
	Abl1/2 Rabbit pAb	ER62565	Human, Mouse	WB,IHC-P
KD	ABL2 Recombinant Rabbit mAb [JF0951]	ET1702-07	Human, Rat, Mouse	WB,FC
KD	AKT1 Recombinant Rabbit mAb [ST05-09]	ET1609-47	Human, Mouse, Rat	WB,IHC-P,IP,FC,IHC-Fr,IF-Cell
	AKT1 Mouse mAb [D9-9-C9]	EM40507	Human, Mouse, Rat	WB,IHC-P
	AKT1 Recombinant Rabbit mAb [PSH0-50] - BSA and Azide free	HA750627	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC
KD	AKT1 Recombinant Rabbit mAb [ST05-09] - BSA and Azide free	HA750177	Human, Mouse, Rat	WB,IHC-P,IP,FC,IHC-Fr,IF-Cell
	AKT1 Recombinant Rabbit mAb [PSH0-50]	HA721364	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC
KD	AKT1/2/3 Recombinant Rabbit mAb [JE75-09]	HA721870	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IP
	AKT1/2/3 Rabbit pAb	ER1901-33	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	AKT1/2/3 Recombinant Rabbit mAb [JE75-09] - BSA and Azide free	HA750840	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IP
	AKT1/2/3 Recombinant Rabbit mAb [ST48-09] - BSA and Azide free	HA750179	Human, Mouse, Rat, Monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
	AKT1/2/3 Recombinant Rabbit mAb [ST48-09]	ET1609-51	Human, Mouse, Rat, Monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
	Akt1/3 Rabbit pAb	ER62638	Human, Mouse, Rat	WB,IHC-P
	Akt2 Mouse mAb [2A1]	RT1027	Human, Mouse, Rat	WB,IP,IF,IHC-P
	AKT2 Rabbit pAb	HA500091	Human, Mouse, Rat	WB,IHC-P
KD	Angiopoietin 1 Recombinant Rabbit mAb [SN67-02]	ET1611-28	Human, Mouse, Rat	WB
	Angiopoietin 2 Recombinant Rabbit mAb [JM71-34] - BSA and Azide free	HA750436	Human, Mouse, Rat	WB,IHC-P,FC
	Angiopoietin 2 Recombinant Rabbit mAb [JM71-34]	ET1705-6	Human, Mouse, Rat	WB,IHC-P,FC
	Angiopoietin-2 Mouse mAb [4A2]	RT1042	Human, Mouse, Rat	WB,IP,IF,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	ARF6 Recombinant Rabbit mAb [JF101-2]	ET1702-91	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
	AS250 Rabbit pAb	ER62747	Human, Mouse	WB,IHC-P,IF-Cell
	Bad Mouse mAb [2H4]	RT1067	Human, Mouse, Rat	WB,IP,IF,IHC-P
	Bad Recombinant Rabbit mAb [SA32-01] - BSA and Azide free	HA750025	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	Bad Recombinant Rabbit mAb [SA32-01]	ET1601-7	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	Bad Rabbit pAb	R1511-32	Human	WB,IF-Cell,IHC-P,FC
	BCL2L1 Rabbit pAb	R1308-16	Human, Mouse, Rat	WB,IHC-P
	Bcl-x Recombinant Rabbit mAb	RE6009	Human, Mouse, Rat	WB,ICC,IHC-P,FC,IP
	Bcl-XL Recombinant Rabbit mAb [SZ3-03]	ET1603-28	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	Bcl-XL Recombinant Rabbit mAb [SZ3-03] - BSA and Azide free	HA750063	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	BDNF Recombinant Rabbit mAb [JE59-35]	HA722912	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IHC-P,IHC-Fr,FC
KD	BDNF Recombinant Rabbit mAb [SJ12-09] - BSA and Azide free	HA750104	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,FC
	BDNF Rabbit pAb	ER130915	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	BDNF Recombinant Rabbit mAb [SJ12-09]	ET1606-42	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,FC
	Biotin Conjugated Beta Nerve Growth Factor Recombinant Rabbit mAb [PSH07-81]	HA722905B	Human	ELISA(Det)
	bFGF Rabbit pAb	ER1904-22	Human, Mouse	WB,FC
	CALL3 Rabbit pAb	ER60555	Human, Mouse, Rat	WB
	CALL6 Rabbit pAb	ER60556	Human	WB
	CALML5 Rabbit pAb	ER62879	Human	IHC-P
	Calmodulin Recombinant Rabbit mAb [SJ16-09] - BSA and Azide free	HA750106	Human, Mouse, Dog, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Calmodulin Recombinant Rabbit mAb [SJ16-09]	ET1606-46	Human, Mouse, Dog, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	CD271 Rabbit pAb	1003-9	Mouse, Rat	WB,IF-Cell,IHC-P
	CDC42 Recombinant Rabbit mAb [JJ086-04]	ET1701-7	Human, Mouse, Rat	WB,FC,IF-Cell
	CDC42 Recombinant Rabbit mAb [JJ086-04] - BSA and Azide free	HA750320	Human, Mouse, Rat	WB,FC,IF-Cell
	c-kit Rabbit pAb	0108-4	Human	WB,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression, KD	c-Kit Recombinant Rabbit mAb [PD00-24]	HA721154	Human	IHC-P,WB,mIHC
	c-Kit Recombinant Rabbit mAb	IRS013	Human	mIHC
Relative expression	c-Kit Recombinant Rabbit mAb [PSH17-96]	HA723961	Mouse, Rat	WB,IHC-Fr,IHC-P
Relative expression, KD	c-Kit Recombinant Rabbit mAb [ST04-99]	ET1609-60	Human	WB,IHC-P
Relative expression	c-Kit Recombinant Rabbit mAb [PSH17-79]	HA723948	Mouse	IHC-Fr,IF-Cell,FC
	Cleaved-Bad (D71) Rabbit pAb	ER60015	Mouse, Rat	WB,IHC-P
	Cleaved-Ephrin-A2 (N188) Rabbit pAb	ER60022	Human, Mouse	WB
Relative expression	c-Met Recombinant Mouse mAb [A9A4-R]	HA601122	Human	WB,IHC-P
Relative expression	c-Met Recombinant Rabbit mAb [PSH0-19]	HA721287	Human	WB,IHC-P,FC
	c-Met Mouse mAb [A9A5]	HA601093	Human	WB,IHC-P
	c-Met Mouse mAb [A9A4]	HA601092	Human	WB,IHC-P
Relative expression	c-Met Recombinant Rabbit mAb [PSH0-19] - BSA and Azide free	HA750594	Human	WB,IHC-P,FC
Relative expression	c-Met Recombinant Mouse mAb [A9A4-R] - BSA and Azide free	HA610037	Human	WB,IHC-P
	cMet Rabbit pAb	ER40501	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	cPLA2 Rabbit pAb	ER63109	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	cPLA2-delta Rabbit pAb	ER63111	Human	IHC-P
	cPLA2-epsilon Rabbit pAb	ER63112	Human	IHC-P
	cPLA2-gamma Rabbit pAb	ER63110	Human	WB
KD	c-Rel Recombinant Rabbit mAb [JM72-93]	ET1705-44	Human	WB,IHC-P,IP
	c-Rel Mouse mAb [1E7]	EM1712-44	Human, Mouse	WB,IHC,IF-Cell
	CSF-1-R Recombinant Rabbit mAb [PD00-45]	HA721180	Human	IHC-P,FC
Relative expression	CSF-1-R Recombinant Rabbit mAb [PSH09-77] - BSA and Azide free	HA751319	Human, Mouse, Rat	WB,IHC-P
Relative expression	CSF-1-R Recombinant Rabbit mAb [PSH09-77]	HA723142	Human, Mouse, Rat	WB,IHC-P
	CSF1R Mouse mAb [B9-B9]	EM1708-56	Human	WB,IF-Cell,IHC-P
	DNA Ligase I Recombinant Rabbit mAb [JB43-39]	ET7107-14	Human, Rat	WB,IHC-P
	EGF Recombinant Rabbit mAb [JA14-15]	ET1704-33	Human, Mouse, Rat	WB,IHC-P,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression, KO	EGFR Mouse mAb [7-F2-F8]	EM1901-67	Human, Monkey	WB,IHC-P,FC
	EGFR Rabbit pAb	R1511-18	Human, Monkey	WB,IF-Cell,IHC-P
Relative expression, KO	EGFR Recombinant Rabbit mAb [SP00-86] - BSA and Azide free	HA750078	Human, Mouse, Rat, Monkey	WB,IF-Cell,IP,FC
Relative expression, KD	EGFR Mouse mAb [A8E12]	HA601073	Human, Monkey	WB,IHC-P
Relative expression, KD	EGFR Recombinant Rabbit mAb [PD00-75]	HA721220	Human, Monkey	WB,IHC-P,IF-Cell,FC
	EGFR Rabbit pAb	ER1512-6	Human, Mouse, Monkey	WB,IF-Cell,IHC-P,FC
	EGFR Rabbit pAb	0407-21	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P
	EGFR Mouse mAb [B6-E5-D9]	EM1701-22	Human, Mouse, Monkey, Rat	WB,IHC-P,IF-Cell,FC
	EGFR Rabbit pAb	R1608-3	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC
	EGFR Recombinant Rabbit mAb	IRS046	Human	mIHC
Relative expression, KD	EGFR Recombinant Rabbit mAb [PD00-75] - BSA and Azide free	HA750550	Human, Monkey	WB,IHC-P,IF-Cell,FC
Relative expression, KO	EGFR Recombinant Rabbit mAb [SP00-86]	ET1604-44	Human, Mouse, Rat, Monkey	WB,IF-Cell,IP,FC
Relative expression, KD	EGFR Recombinant Rabbit mAb [SZ40-19]	ET1603-37	Human, Mouse, Monkey, Rat	WB,IF-Cell,IP,FC,mIHC
	ELK1 Recombinant Rabbit mAb [SU32-05]	ET1608-19	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
	EphA2/3/4 Rabbit pAb	ER63372	Human, Rat	WB,IF-Cell
	EphA2 (N-term) Rabbit pAb	ER2001-31	Human, Mouse	WB,IHC-P,FC
	EphA3 Rabbit pAb	ER1908-23	Mouse, Rat	WB,IHC-P,IF-Tissue
	Eph receptor A2 Recombinant Rabbit mAb [PSH09-21] - BSA and Azide free	HA751275	Human	WB,IF-Cell,IP,FC
	Eph receptor A2 Recombinant Rabbit mAb [PSH09-21]	HA723079	Human	WB,IF-Cell,IP,FC
	Ephrin A1 Mouse mAb [A7B10]	HA601038	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Ephrin A1 Recombinant Mouse mAb [A7B10-R] - BSA and Azide free	HA610173	Human	WB
	Ephrin A1 Recombinant Mouse mAb [A7B10-R]	HA601300	Human	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Ephrin-A2 Rabbit pAb	ER63379	Human, Mouse	WB
	Ephrin-A3 Rabbit pAb	ER63380	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Ephrin A5 Rabbit pAb	ER1908-27	Human, Mouse	WB
KD	ERK1 Recombinant Rabbit mAb [SP00-74]	ET1604-32	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	ERK1 Recombinant Rabbit mAb [SP05-09]	ET1604-16	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	ERK1 Mouse mAb [7-D6-E5]	M1505-9	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Erk1/2 Mouse mAb [2G1]	RT1484	Human, Mouse, Rat	WB,IP,IF
	ERK1/2 Rabbit pAb	ER131011	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	ERK1/2 Recombinant Rabbit mAb [SA43-03] - BSA and Azide free	HA750019	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IP,FC,IHC-P,IHC-Fr
KD	ERK1/2 Recombinant Rabbit mAb [SA43-03]	ET1601-29	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IP,FC,IHC-P,IHC-Fr
	ERK2 Mouse mAb [5-D2]	M1407-3	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	ERK2 Recombinant Mouse mAb [C3B12-R] - BSA and Azide free	HA610079	Human, Mouse, Rat	WB
KD	ERK2 Recombinant Rabbit mAb [SZ25-01]	ET1603-23	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
KD	ERK2 Recombinant Mouse mAb [C3B12-R]	HA601201	Human, Mouse, Rat	WB
KD	ERK2 Mouse mAb [B10-8]	EM1901-53	Human, Mouse, Rat	WB,IHC-P,FC
KD	ERK2 Mouse mAb [A10D2]	EM1901-54	Human, Mouse, Rat	WB,IHC-P,FC
KD	ERK2 Rabbit pAb	ER131218	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	ERK2 Recombinant Mouse mAb [5-D2-R]	HA601202	Human, Mouse, Rat	WB
KD	ERK2 Recombinant Mouse mAb [5-D2-R] - BSA and Azide free	HA610080	Human, Mouse, Rat	WB
	ERK2 Mouse mAb [C3B12]	EM1901-55	Human, Mouse, Rat	WB,IHC-P,FC
	ETS1 Recombinant Rabbit mAb [JM92-32]	ET1705-23	Human, Mouse, Rat	WB,IHC-P
	ETS2 Rabbit pAb	ER62455	Human, Mouse	WB
	EXOC2 Recombinant Rabbit mAb [JE65-98]	HA721386	Human, Mouse, Rat	WB,IHC-P
	Fas Ligand Rabbit pAb	ER1908-69	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
	FGF-1 Rabbit pAb	ER63425	Human, Mouse, Rat	IHC-P
	FGF1 Recombinant Rabbit mAb [JE33-93]	HA722695	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	FGF1 Recombinant Rabbit mAb [PSH18-61]	HA724014	Human, Mouse, Rat	WB,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	FGF10 Rabbit pAb	ER60274	Human, Mouse, Rat	WB
	FGF-16 Rabbit pAb	ER65578	Human, Mouse, Rat	WB,IHC-P
	FGF-17 Rabbit pAb	ER65622	Human, Mouse, Rat	WB
	FGF-18 Rabbit pAb	ER60001	Human, Mouse, Rat	WB,IHC-P
Cell treatment	FGF18 Recombinant Rabbit mAb [PSH01-24]	HA721641	Human, Mouse	WB
Cell treatment	FGF18 Recombinant Rabbit mAb [PSH01-25]	HA721642	Mouse	WB
Cell treatment	FGF18 Recombinant Rabbit mAb [PSH01-25] - BSA and Azide free	HA750701	Mouse	WB
Cell treatment	FGF18 Recombinant Rabbit mAb [PSH01-24] - BSA and Azide free	HA750700	Human, Mouse	WB
	FGF2 Recombinant Rabbit mAb [JM28-10]	ET1703-18	Human, Mouse	WB,IP,FC,IHC-P
	FGF-20 Rabbit pAb	ER65477	Human, Mouse, Rat	WB,IHC-P
	FGF21 Recombinant Rabbit mAb [JA10-97]	ET1704-04	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	FGF21 Recombinant Rabbit mAb [JA10-97] - BSA and Azide free	HA750399	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	FGF-22 Rabbit pAb	ER63427	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	FGF23 Rabbit pAb	ER1908-82	Human, Rat	WB
Relative expression	FGF23 Recombinant Rabbit mAb [PSH19-00]	HA724036	Human, Mouse, Rat	WB
	FGF3 Rabbit pAb	ER61973	Human, Mouse	WB
	FGF5 Rabbit pAb	ER1908-83	Human, Mouse, Rat	IHC-P,IF-Tissue
	FGF-6 Rabbit pAb	ER65475	Human, Mouse, Rat	WB
	FGF-8 Rabbit pAb	ER65576	Human, Mouse, Rat	WB
	FGF-9 Rabbit pAb	ER65478	Human, Mouse, Rat	WB,IHC-P
Relative expression, KD	FGFR1 Recombinant Rabbit mAb [SD08-25]	ET1612-62	Human	WB,IF-Cell
	FGFR2 Recombinant Rabbit mAb [JM10-60]	ET1703-78	Human, Mouse, Rat	WB,IP
	FGFR2/CD332 Recombinant Mouse mAb [A5-E10-G7-R] - BSA and Azide free	HA610003	Human	WB
	FGFR2/CD332 Rabbit pAb	R1311-7	Human	WB,IHC-P
	FGFR2/CD332 Rabbit pAb	R1310-10	Human	WB,IF-Cell,IHC-P,FC
	FGFR2/CD332 Recombinant Mouse mAb [A5-E10-G7-R]	EM50103	Human	WB
KD	FGFR3 Recombinant Rabbit mAb [JM81-10]	ET1703-36	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	FGFR3 Recombinant Rabbit mAb [JM110-33] - BSA and Azide free	HA750388	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
KD	FGFR3 Recombinant Rabbit mAb [JM110-33]	ET1703-65	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
KD	FGFR3 Recombinant Rabbit mAb [JM81-10] - BSA and Azide free	HA750379	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	FGFR-4 Rabbit pAb	ER65688	Human	WB
	FGFR4 Recombinant Rabbit mAb [JE59-96]	HA721578	Human	WB,FC
	Flk-1/Flt-4 Rabbit pAb	ER63443	Human, Mouse, Rat	IHC-P,IF-Cell
	Flt3 Rabbit pAb	ER2001-25	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell
	FLT3 Mouse mAb [C3-B7]	EM1709-22	Human	WB,IHC-P
	FLT3L Rabbit pAb	ER1908-92	Human, Mouse	WB
	Flt-4 Rabbit pAb	ER65750	Human, Mouse, Rat	WB,IHC-P
	FoxO1/3/4 Rabbit pAb	ER63465	Human, Mouse, Rat	WB,IHC-P
	iFluor™ 488 Conjugated FOXO1A Recombinant Rabbit mAb [SU33-01]	HA720161F	Human, Mouse, Rat	IF-Cell,FC
KD	FOXO1A Recombinant Rabbit mAb [SU33-01]	ET1608-25	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC
KD	FOXO1A Recombinant Rabbit mAb [SU33-01] - BSA and Azide free	HA750140	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC
	FOXO3A Rabbit pAb	ER1706-79	Human, Mouse	WB,IF-Cell,IHC-P,FC
	FOXO3A Recombinant Rabbit mAb [SP0054]	ET1604-11	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue
	FOXO4 Recombinant Rabbit mAb [SD0817]	ET1612-66	Human, Rat	WB,IHC-P
	FOXO4 Recombinant Rabbit mAb [JJ09-11]	ET1701-67	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,FC,IP
	G3BP Recombinant Rabbit mAb [JE35-20]	HA721314	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	G3BP1 Rabbit pAb	ER63481	Human, Mouse, Monkey	WB,IHC-P
	GAB1 Recombinant Rabbit mAb [JU47-19]	ET7106-70	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Gab2 Rabbit pAb	HA500292	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	GAP1-InsP4 BP Rabbit pAb	ER63502	Human, Mouse, Rat	WB
	Gbeta5 Rabbit pAb	HA500513	Human, Mouse, Rat	WB,IHC-P,FC
	G gamma 5 Rabbit pAb	ER63661	Human, Mouse, Rat	IF-Cell
	GNB1 Rabbit pAb	HA500392	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	GNB2 Recombinant Rabbit mAb [JE47-16]	ET7109-53	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	GNB3 Rabbit pAb	HA500127	Human, Mouse, Rat	WB,IF-Cell,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	GNB3 Rabbit pAb	HA500153	Human, Mouse, Rat	WB,IHC-P
	GNB4 Rabbit pAb	HA500501	Human, Mouse, Rat	WB
	GNG11 Rabbit pAb	HA500256	Human, Mouse, Rat	WB,IF-Cell,FC
	GNG12 Rabbit pAb	ER65918	Human, Mouse	IHC-P
	GNG13 Rabbit pAb	ER65917	Human, Mouse	IHC-P
	GNG2 Rabbit pAb	ER65936	Human, Mouse	IHC-P
	GNG3 Rabbit pAb	HA500380	Human, Mouse, Rat	WB,IHC-P
	GNG4 Rabbit pAb	ER65937	Human, Mouse	IHC-P
	GNG4 Rabbit pAb	HA500259	Human, Mouse	WB,IF-Cell,IHC-P
	GNG7 Rabbit pAb	ER65938	Human, Mouse, Rat	IHC-P
	GNG8 Rabbit pAb	ER65919	Human, Mouse, Rat	IHC-P
	GRB2 Rabbit pAb	ER31203	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	GRB2 Recombinant Rabbit mAb [JE49-73]	HA721923	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC
	GRB2 Mouse mAb [11-B11]	EM1701-41	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Group VI iPLA2 Rabbit pAb	ER63644	Human, Mouse, Rat	WB
	GTPase HRAS Recombinant Rabbit mAb [JE42-15]	HA721355	Human, Mouse, Rat	WB,IHC-P
	HGF Recombinant Rabbit mAb [PSH20-29]	HA724140	Human	IHC-P
	HGF Recombinant Rabbit mAb [PSH03-79] - BSA and Azide free	HA751614	Human	WB
	HGF Rabbit pAb	ER65684	Human	WB,IHC-P
	HGF Recombinant Rabbit mAb [PSH03-79]	HA723896	Human	WB
	Human Beta Nerve Growth Factor Recombinant Rabbit mAb [PSH07-79] - BSA and Azide free (Capture)	HA722902	Human	ELISA(Cap)
	Human Beta Nerve Growth Factor Recombinant Rabbit mAb [PSH07-80] - BSA and Azide free (Capture)	HA722903	Human	ELISA(Cap)
	Human Beta Nerve Growth Factor Recombinant Rabbit mAb [PSH07-81] - BSA and Azide free (Detector)	HA722904	Human	ELISA(Det)
	Biotin Conjugated Human CSF-1-R Recombinant Rabbit mAb [PSH12-66] - Detector	HA723478B	Human	ELISA(Det),ELISA
	Human CSF-1-R Recombinant Rabbit mAb [PSH12-67] - BSA and Azide free (Detector)	HA723475	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human CSF-1-R Recombinant Rabbit mAb [PSH12-68] - BSA and Azide free (Capture)	HA723477	Human	ELISA(Cap)
	Biotin Conjugated Human CSF-1-R Recombinant Rabbit mAb [PSH12-67] - Detector	HA723476B	Human	ELISA(Det),ELISA
	HRP Conjugated Human CSF-1-R Recombinant Rabbit mAb [PSH12-66] - Detector	HA723478H	Human	ELISA(Det),ELISA
	Human CSF-1-R Recombinant Rabbit mAb [PSH12-66] - BSA and Azide free (Capture/Detector)	HA723474	Human	ELISA(Cap),ELISA(Det)
	Biotin Conjugated Human EGF Recombinant Rabbit mAb [PSH02-20] - Detector	HA722638B	Human	ELISA(Det)
	Human EGF Recombinant Rabbit mAb [PSH02-19] - BSA and Azide free (Capture)	HA721782	Human	ELISA(Cap)
	Human EGF Recombinant Rabbit mAb [PSH02-20] - BSA and Azide free (Detector)	HA721783	Human	ELISA(Det)
	Human EphA2 Recombinant Rabbit mAb [PSH15-13] - BSA and Azide free (Detector)	HA723715	Human	ELISA(Det)
	Biotin Conjugated Human EphA2 Recombinant Rabbit mAb [PSH15-13] - Detector	HA723716B	Human	ELISA(Det),ELISA
	Human EphA2 Recombinant Rabbit mAb [PSH15-12] - BSA and Azide free (Capture)	HA723714	Human	ELISA(Cap)
	Human Fas Ligand Recombinant Rabbit mAb [PSH18-87] - BSA and Azide free (Detector)	HA724024	Human	ELISA(Det)
	Biotin Conjugated Human Fas Ligand Recombinant Rabbit mAb [PSH18-88] - Detector	HA724025B	Human	ELISA(Det),ELISA
	Human Fas Ligand Recombinant Rabbit mAb [PSH18-88] - BSA and Azide free (Detector)	HA724025	Human	ELISA(Det)
	Biotin Conjugated Human Fas Ligand Recombinant Rabbit mAb [PSH18-87] - Detector	HA724024B	Human	ELISA(Det),ELISA
	Human Fas Ligand Recombinant Rabbit mAb [PSH18-86] - BSA and Azide free (Capture)	HA724023	Human	ELISA(Cap)
	Biotin Conjugated Human FGF-2 Recombinant Rabbit mAb [PSH02-09]	HA722639B	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human FGF-2 Recombinant Rabbit mAb [PSH02-09] - BSA and Azide free (Detector)	HA721781	Human	ELISA(Det)
	Human FGF-2 Recombinant Rabbit mAb [PSH02-08] - BSA and Azide free (Capture)	HA721780	Human	ELISA(Cap)
	Biotin Conjugated Human Flt3 Ligand / FLT3LG Recombinant Rabbit mAb [PSH20-11] - Detector	HA725326B	Human	ELISA(Det),ELISA
	Human Flt3 Ligand / FLT3LG Recombinant Rabbit mAb [PSH20-10] - BSA and Azide free (Capture)	HA725325	Human	ELISA(Cap)
	Human Flt3 Ligand / FLT3LG Recombinant Rabbit mAb [PSH20-11] - BSA and Azide free (Detector)	HA725326	Human	ELISA(Det)
	Biotin Conjugated Human HGF Recombinant Rabbit mAb [PSH02-07]	HA722640B	Human	ELISA(Det)
	Human HGF Recombinant Rabbit mAb [PSH03-78] - BSA and Azide free (Capture)	HA722042	Human	ELISA(Cap)
	Human HGF Recombinant Rabbit mAb [PSH02-07] - BSA and Azide free (Detector)	HA721779	Human	ELISA(Det)
	Human HGF Recombinant Rabbit mAb [PSH03-79] - BSA and Azide free (Detector)	HA722043	Human	ELISA(Det)
	Biotin Conjugated Human HGF Recombinant Rabbit mAb [PSH03-79]	HA722619B	Human	ELISA(Det)
	Human HGF Recombinant Rabbit mAb [PSH02-06] - BSA and Azide free (Capture)	HA721778	Human	ELISA(Cap)
	Human M-CSF Recombinant Rabbit mAb [PSH11-36] - BSA and Azide free (Capture)	HA723323	Human	ELISA(Cap)
	Biotin Conjugated Human M-CSF Recombinant Rabbit mAb [PSH11-38] - Detector	HA723326B	Human	ELISA(Det),ELISA
	Human M-CSF Recombinant Rabbit mAb [PSH11-38] - BSA and Azide free (Detector)	HA723325	Human	ELISA(Det)
	Human M-CSF Recombinant Rabbit mAb [PSH11-37] - BSA and Azide free (Capture)	HA723324	Human	ELISA(Cap)
	Biotin Conjugated Human/Mouse/Rat FGF1 Recombinant Rabbit mAb [PSH18-35] - Detector	HA723988B	Human, Mouse, Rat	ELISA(Det),ELISA
	Human/Mouse/Rat FGF1 Recombinant Rabbit mAb [PSH18-34] - BSA and Azide free (Capture)	HA723987	Human, Mouse, Rat	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human/Mouse/Rat FGF1 Recombinant Rabbit mAb [PSH18-35] - BSA and Azide free (Detector)	HA723988	Human, Mouse, Rat	ELISA(Det)
	Biotin Conjugated Human/Mouse VEGFC Recombinant Rabbit mAb [PSH10-76] - Detector	HA723243B	Human, Mouse	ELISA(Det),ELISA
	Human/Mouse VEGFC Recombinant Rabbit mAb [PSH10-75] - BSA and Azide free (Capture)	HA723241	Human, Mouse	ELISA(Cap)
	Human/Mouse VEGFC Recombinant Rabbit mAb [PSH10-76] - BSA and Azide free (Detector)	HA723242	Human, Mouse	ELISA(Det)
	Human PDGF-BB Recombinant Rabbit mAb [PSH05-98] - BSA and Azide free (Capture)	HA722497	Human	ELISA(Cap)
	Biotin Conjugated Human PDGF-BB Recombinant Rabbit mAb [PSH05-99] - Detector	HA722585B	Human	ELISA(Det)
	Human PDGF-BB Recombinant Rabbit mAb [PSH05-99] - BSA and Azide free (Detector)	HA722498	Human	ELISA(Det)
	Human SCF Recombinant Rabbit mAb [PSH07-57] - BSA and Azide free (Capture)	HA722871	Human	ELISA(Cap)
	Human SCF Recombinant Rabbit mAb [PSH07-58] - BSA and Azide free (Detector)	HA722872	Human	ELISA(Det)
	Biotin Conjugated Human SCF Recombinant Rabbit mAb [PSH07-58]	HA722873B	Human	ELISA(Det)
	Biotin Conjugated Human VEGFA Recombinant Rabbit mAb [PSH12-92] - Detector	HA723495B	Human	ELISA(Det),ELISA
	Human VEGFA Recombinant Rabbit mAb [PSH12-91] - BSA and Azide free (Capture)	HA723493	Human	ELISA(Cap)
	Human VEGFA Recombinant Rabbit mAb [PSH12-92] - BSA and Azide free (Detector)	HA723494	Human	ELISA(Det)
	Biotin Conjugated Human VEGFA Recombinant Rabbit mAb [PSH12-93] - Detector	HA723497B	Human	ELISA(Det),ELISA
	Human VEGFA Recombinant Rabbit mAb [PSH12-93] - BSA and Azide free (Detector)	HA723496	Human	ELISA(Det)
	Human VEGF Receptor 1 Recombinant Rabbit mAb [PSH19-22] - BSA and Azide free (Capture)	HA725365	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human VEGF Receptor 1 Recombinant Rabbit mAb [PSH19-23] - BSA and Azide free (Detector)	HA725366	Human	ELISA(Det)
	Biotin Conjugated Human VEGF Receptor 1 Recombinant Rabbit mAb [PSH19-23] - BSA and Azide free (Detector)	HA725366B	Human	ELISA(Det),ELISA
	Human VEGF Receptor 2 Recombinant Rabbit mAb [PSH12-83] - BSA and Azide free (Capture)	HA725121	Human	ELISA(Cap)
	Biotin Conjugated Human VEGF Receptor 2 Recombinant Rabbit mAb [PSH12-84] - Detector	HA725123B	Human	ELISA(Det),ELISA
	Human VEGF Receptor 2 Recombinant Rabbit mAb [PSH12-84] - BSA and Azide free (Detector)	HA725122	Human	ELISA(Det)
Cell treatment	IGF1 Recombinant Rabbit mAb [JE30-01]	HA722527	Human	WB
	IGF1 (E peptide) Rabbit pAb	HA500107	Human	WB,IHC-P,FC
	IGF1R Mouse mAb [D5-G2]	HA600033	Human	WB,FC
	IGF1R Rabbit pAb	R1310-6	Human	WB,IHC-P
Relative expression	IGF1 Receptor beta Recombinant Rabbit mAb [PSH11-72] - BSA and Azide free	HA751427	Human, Mouse, Rat	WB,IF-Cell
Relative expression	IGF1 Receptor beta Recombinant Rabbit mAb [PSH11-72]	HA723368	Human, Mouse, Rat	WB,IF-Cell
	IGF2 Recombinant Rabbit mAb [JJ092-3]	ET1701-88	Human	WB,IF-Cell,IF-Tissue,IHC-P
	IGF I Rabbit pAb	ER1911-10	Human, Mouse, Rat, Dog	WB,IHC-P,IF-Tissue
	IGF-IR Rabbit pAb	ER63734	Human, Mouse, Rat	IF-Cell,WB,IHC-P
	IKK alpha Rabbit pAb	ER30911	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	IKK alpha Recombinant Rabbit mAb [SN06-99]	ET1611-15	Human, Mouse, Rat	WB,IHC-P,IP,FC
	IKK alpha Recombinant Rabbit mAb [SN06-99] - BSA and Azide free	HA750246	Human, Mouse, Rat	WB,IHC-P,IP,FC
KD	IKK alpha + IKK beta Recombinant Rabbit mAb [SN63-02] - BSA and Azide free	HA750248	Human, Mouse, Rat	WB,IF-Cell,IP
KD	IKK alpha + IKK beta Recombinant Rabbit mAb [SN63-02]	ET1611-23	Human, Mouse, Rat	WB,IF-Cell,IP
	IKK beta Rabbit pAb	R1706-13	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	IKK beta Recombinant Rabbit mAb [JE32-81]	HA721715	Human	WB,FC
KD	IKK gamma Recombinant Rabbit mAb [JA11-83]	ET1704-40	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,IP,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Insulin Recombinant Rabbit mAb [PSH04-16] - Detector	HA722610B	Human, Mouse, Rat	ELISA(Det),ELISA
	Insulin Mouse mAb [5E1]	EM1701-33	Human	IHC-P,ELISA
	Insulin Mouse mAb [A6-6]	EM80714	Human, Mouse, Rat	IHC-P,ELISA,IF-Tissue,IHC-Fr
	Insulin Recombinant Rabbit mAb [PSH04-15] - BSA and Azide free (Capture)	HA722096	Human, Mouse, Rat	ELISA(Cap)
Relative expression	Insulin Recombinant Rabbit mAb [SA0410]	ET1601-12	Human, Mouse, Rat	IF-Cell,IF-Tissue,IHC-P,mIHC,IHC-Fr,WB
	Insulin Recombinant Rabbit mAb	IRS161RB	Human	mIHC
	Insulin Recombinant Rabbit mAb	IRS233RB	Mouse	mIHC
Relative expression	Insulin Recombinant Rabbit mAb [SA0410] - BSA and Azide free	HA750004	Human, Mouse, Rat	IF-Cell,IF-Tissue,IHC-P,IHC-Fr,WB
	Insulin Recombinant Rabbit mAb [PSH04-16] - BSA and Azide free (Detector)	HA722097	Human, Mouse, Rat	ELISA(Det)
	Insulin B Chain Rabbit pAb	0807-11	Human, Mouse, Rat	IHC-P,ELISA,IF-Tissue
KD, Relative expression	Insulin Receptor Beta Recombinant Rabbit mAb [SN20-13] - BSA and Azide free	HA750276	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC
KD, Relative expression	Insulin Receptor Beta Recombinant Rabbit mAb [SN20-13]	ET1611-90	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC
	JNK1 Rabbit pAb	R1309-1	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC
	JNK1 Recombinant Rabbit mAb [JE63-04]	HA721899	Human, Mouse, Rat	WB,IHC-P
	JNK1+JNK2+JNK3 Recombinant Rabbit mAb [SA43-06] - BSA and Azide free	HA750018	Human, Mouse, Rat, Cow, Monkey, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP
	JNK1+JNK2+JNK3 Recombinant Rabbit mAb [SA43-06]	ET1601-28	Human, Mouse, Rat, Cow, Monkey, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP
	JNK1+JNK3 Recombinant Rabbit mAb [JJ08-86]	ET1701-47	Human, Mouse, Rat	WB,IP,IHC-P
	JNK2 Rabbit pAb	HA500110	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	JNK2 Recombinant Rabbit mAb [SC05-85]	ET1610-11	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
	JNK3 Recombinant Rabbit mAb [SD082-09]	ET1612-68	Human, Mouse, Rat	WB,IF-Tissue,IHC-P
	KDR / VEGFR-2 Mouse mAb [1G1]	RT1347	Human, Mouse, Rat	WB,IP,IF,IHC-P
	K-Ras Rabbit pAb	ER65706	Human, Mouse, Rat	WB,IHC-P
	KRAS Recombinant Mouse mAb [A8E6]	HA601060	Human, Mouse, Rat	WB,IHC-P
	KRAS Mouse mAb [A8E4]	HA601058	Human, Mouse, Rat	WB,IF-Cell,FC
	KRAS Mouse mAb [A8E5]	HA601059	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	KRAS Recombinant Mouse mAb [A8E6] - BSA and Azide free	HA610008	Human, Mouse, Rat	WB,IHC-P
	Krs-1/2 Rabbit pAb	ER63820	Human, Mouse	WB,IHC-P
	Ksr-1 Rabbit pAb	ER63821	Human, Mouse	WB,IHC-P,IF-Cell
	KSR1 Recombinant Rabbit mAb [JE57-75]	HA721963	Human, Mouse, Rat, Monkey	WB,IHC-P,FC
	Ksr2 Rabbit pAb	ER63822	Human, Mouse	WB,IHC-P,IF-Cell
Relative expression	LAT Recombinant Rabbit mAb [PSH09-76] - BSA and Azide free	HA751318	Human	WB,IF-Cell,FC
	LAT Rabbit pAb	ER63843	Human, Mouse, Rat	WB
Relative expression	LAT Recombinant Rabbit mAb [PSH09-76]	HA723141	Human	WB,IF-Cell,FC
	M-CSF Recombinant Rabbit mAb [SU0413]	ET1609-1	Human, Mouse	WB,IHC-P,IF-Cell,IP
KD	MEK1 Recombinant Rabbit mAb [SZ22-01]	ET1603-20	Human, Mouse, Rat, Cow, Dog	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	MEK1/2 Recombinant Rabbit mAb [SR13-07]	ET1602-3	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	MEK1/2 Recombinant Rabbit mAb [SR13-07] - BSA and Azide free	HA750042	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	MEK2 Recombinant Rabbit mAb [SD207-08]	ET1612-6	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Met (c-Met) Rabbit pAb	HA500242	Human	WB,IF-Cell,IHC-P
KD	Met (C-Met) Recombinant Rabbit mAb [SJ19-05]	ET1606-45	Human	WB,IF-Cell,IF-Tissue,IHC-P
	Mouse FGF-21 Recombinant Rabbit mAb [PSH15-27] - BSA and Azide free (Detector)	HA723732	Mouse	ELISA(Det)
	Mouse FGF-21 Recombinant Rabbit mAb [PSH15-26] - BSA and Azide free (Capture)	HA723731	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse FGF-21 Recombinant Rabbit mAb [PSH15-27] - Detector	HA723733B	Mouse	ELISA(Det),ELISA
	Biotin Conjugated Mouse M-CSF/CSF1 Recombinant Rabbit mAb [PSH15-85] - Detector	HA723794B	Mouse	ELISA(Det),ELISA
	Mouse M-CSF/CSF1 Recombinant Rabbit mAb [PSH15-84] - BSA and Azide free (Capture)	HA723792	Mouse	ELISA(Cap)
	Mouse M-CSF/CSF1 Recombinant Rabbit mAb [PSH15-85] - BSA and Azide free (Detector)	HA723793	Mouse	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse/Rat PDGF-BB Recombinant Rabbit mAb [PSH16-05] - BSA and Azide free (Capture)	HA723810	Mouse, Rat	ELISA(Cap)
	Mouse/Rat PDGF-BB Recombinant Rabbit mAb [PSH16-06] - BSA and Azide free (Detector)	HA723811	Mouse, Rat	ELISA(Det)
	Mouse/Rat PDGF-BB Recombinant Rabbit mAb [PSH16-07] - BSA and Azide free (Capture)	HA723812	Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Mouse/Rat PDGF-BB Recombinant Rabbit mAb [PSH16-06] - Detector	HA723811B	Mouse, Rat	ELISA(Det),ELISA
	Mouse VEGF Recombinant Rabbit mAb [PSH02-72] - BSA and Azide free (Capture)	HA721850	Mouse	ELISA(Cap)
	Mouse VEGF Recombinant Rabbit mAb [PSH02-73] - BSA and Azide free (Detector)	HA721851	Mouse	ELISA(Det)
	Biotin Conjugated Mouse VEGF Recombinant Rabbit mAb [PSH02-73]	HA722636B	Mouse	ELISA(Det)
	MRas Recombinant Rabbit mAb [JE64-57]	HA722579	Human, Mouse, Rat	WB
	MST1 Rabbit pAb	ER1913-22	Human, Rat	IHC-P,IF-Cell,IF-Tissue,FC
KD	MST2 Recombinant Rabbit mAb [SC05-83]	ET1610-8	Human, Mouse, Rat	WB,IHC-P,IP
	NAK / TBK1 Recombinant Rabbit mAb [PSH11-69]	HA723365	Human, Mouse, Rat	WB,IF-Cell,FC,IP
KD	NAK / TBK1 Mouse mAb [A8D9]	HA601045	Human, Mouse, Rat	WB
KD	NAK / TBK1 Mouse mAb [A8D10]	HA601044	Human, Mouse, Rat	WB
	NAK / TBK1 Recombinant Rabbit mAb [PSH11-69] - BSA and Azide free	HA751425	Human, Mouse, Rat	WB,IF-Cell,FC,IP
	NAK / TBK1 Mouse mAb [2G1]	RT1413	Human, Mouse, Rat	WB,IP
	Neurofibromin Rabbit pAb	ER64168	Human, Mouse, Rat	WB,IF-Cell
	Neurotrophin 4 Rabbit pAb	ER1913-76	Human, Rat	WB,IHC-P,IF-Tissue
	NF-kB p105 / p50 Mouse mAb [A2-7]	M1505-8	Human, Mouse	WB,IF-Cell,IHC-P
	NF-kB p105 / p50 Rabbit pAb	R1309-9	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	NF-kB p105 / p50 Recombinant Rabbit mAb [SZ20-01]	ET1603-18	Human, Mouse, Rat	WB,IHC-P,ICC,FC,IF-Tissue
KD	NF-kB p105 / p50 Recombinant Rabbit mAb [SZ20-01] - BSA and Azide free	HA750060	Human, Mouse, Rat	WB,IHC-P,ICC,FC,IF-Tissue
	NF-kB p65 Rabbit pAb	HA500402	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
	NF-kB p65 Rabbit pAb	ER0815	Human, Mouse, Rat, Zebrafish	WB,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KO, KD, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [SZ10-04] - BSA and Azide free	HA750058	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KO, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [PSH0-27] - BSA and Azide free	HA750601	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
KO, KD, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [SZ10-04]	ET1603-12	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KO, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [PSH0-27]	HA721307	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
	NGF Rabbit pAb	ER1802-59	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	NGF Recombinant Rabbit mAb [SI79-01]	ET1606-29	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
	NGF Recombinant Rabbit mAb [SI79-01] - BSA and Azide free	HA750099	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
	NGF beta Mouse mAb [1G4]	RT1428	Human, Mouse, Rat	WB,IP,IF
	NGF Receptr Mouse mAb [2G2]	RT1429	Human, Mouse, Rat	WB,IP,IF
	N/H/K-Ras Rabbit pAb	ER64095	Human, Mouse, Rat	WB,IHC-P
Relative expression	NMDAR1 Recombinant Rabbit mAb [JM11-26] - BSA and Azide free	HA750391	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Tissue,IHC-P,IHC-Fr
Relative expression	NMDAR1 Recombinant Rabbit mAb [JM11-26]	ET1703-75	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Tissue,IHC-P,IHC-Fr
	NMDAR2A Recombinant Rabbit mAb [JA31-20]	ET1704-80	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue
Relative expression	NMDAR2B Recombinant Rabbit mAb [PSH05-46]	HA722284	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr
	NMDAR2B Rabbit pAb	ER1913-90	Human, Rat	WB,IHC-P,IF-Tissue
	NMDAR2B Rabbit pAb	HA500208	Mouse, Rat, Human	WB,IF-Cell,IHC-P,FC
Relative expression	NMDAR2B Recombinant Rabbit mAb [PSH05-46] - BSA and Azide free	HA751005	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr
	N-Shc Rabbit pAb	ER64229	Human, Mouse, Rat	IHC-P
	NT-3 Rabbit pAb	ER65762	Human, Mouse, Rat	IHC-P
	p75 NGF Receptor Recombinant Antibody [PSH13-26] - Rabbit IgG (Chimeric) - BSA and Azide free	HA751462	Mouse	IHC-Fr
	p75 NGF Receptor Recombinant Rat mAb [PSH13-26] - BSA and Azide free	HA610274	Mouse, Rat	IHC-Fr,WB
KD	p75 NGF Receptor Recombinant Rabbit mAb [SA39-02]	ET1601-22	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	p75 NGF Receptor Recombinant Rat mAb [PSH13-26]	HA601436	Mouse, Rat	IHC-Fr, WB
	iFluor™ 488 Conjugated p75 NGF Receptor Recombinant Rabbit mAb [SA39-02]	HA720190F	Human, Mouse, Rat	IF-Cell
KD	p75 NGF Receptor Recombinant Rabbit mAb [SA39-02] - BSA and Azide free	HA750014	Human, Mouse, Rat	WB, IF-Cell, IF-Tissue, IHC-P, IP, FC
	p75 NGF Receptor Recombinant Antibody [PSH13-26] - Rabbit IgG (Chimeric)	HA723532	Mouse	IHC-Fr
	PA21B Rabbit pAb	ER60996	Human	WB
	PAK1 Recombinant Rabbit mAb [JE57-36]	ET7111-33	Human, Mouse, Rat	WB, IF-Cell, IHC-P, FC
	PAK1 + PAK2 + PAK3 Rabbit pAb	ER64516	Human, Mouse, Rat	WB, IHC-P
KD	PAK2 Recombinant Rabbit mAb [JB41-84]	ET7107-79	Human, Mouse, Rat, Monkey	WB, IHC-P, FC, IF-Cell, IF-Tissue
	PAK3 Mouse mAb [F11-G8]	EM1708-48	Human, Monkey	WB
	PAK3 Recombinant Rabbit mAb [JE57-39]	ET7111-35	Human, Mouse, Rat	WB
	PAK5 Rabbit pAb	ER64513	Human, Mouse, Rat	WB, IHC-P
	PAK5/6 Rabbit pAb	ER64514	Human, Mouse	WB
	PAK6 Rabbit pAb	ER60159	Human, Mouse	WB
	PC-PLD1 Rabbit pAb	ER64532	Human, Mouse, Rat	IHC-P
	PDGF B Recombinant Rabbit mAb [JE75-16]	HA722003	Human, Mouse, Rat	WB, IHC-P
	PDGF-B Rabbit pAb	ER1914-92	Human, Mouse, Rat	WB, IHC-P, IF-Tissue, FC
KD	PDGFR alpha Recombinant Rabbit mAb [JF104-6]	ET1702-49	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB, IHC-Fr
	PDGFR beta Rabbit pAb	R1510-44	Human, Mouse, Rat	WB, IF-Cell, IHC-P, FC
	PDGFR beta Recombinant Rabbit mAb	IRS058	Mouse	mIHC
	PDGFR beta Recombinant Antibody - Rat IgG1 (Chimeric)	IRS283RT	Mouse	mIHC
	PDGFR beta Recombinant Rabbit mAb	IRS038	Human	mIHC
	PDGFR beta Recombinant Rabbit mAb [SY10-08] - BSA and Azide free	HA750082	Human, Mouse, Rat	WB, IHC-P, IP, FC, IF-Tissue, IF-Cell, IHC-Fr
	PDGFR beta Recombinant Rabbit mAb [SY10-08]	ET1605-20	Human, Mouse, Rat	WB, IHC-P, IP, FC, IF-Tissue, IF-Cell, mIHC, IHC-Fr
Cell treatment	Phospho-AKT1 (S124) Recombinant Rabbit mAb [JJ08-46]	ET1701-36	Human, Mouse, Rat	WB, IF-Cell, IHC-P, IP
Cell treatment	Phospho-AKT1 (S129) Recombinant Rabbit mAb [JE59-50]	HA721677	Human, Mouse, Rat	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-AKT1 (T450) Recombinant Rabbit mAb [SD08-12]	ET1612-73	Human, Mouse, Rat	WB,IHC-P,IP,IF-Cell,FC
Cell treatment	Phospho-AKT2 (S474) Recombinant Rabbit mAb [PSH09-46]	HA723108	Human, Mouse, Rat	WB,FC
Cell treatment	Phospho-AKT2 (S474) Recombinant Rabbit mAb [PSH09-46] - BSA and Azide free	HA751293	Human, Mouse, Rat	WB,FC
Cell treatment	Phospho-AKT2 (S474) Rabbit pAb	HA500116	Human, Mouse	WB,IHC-P
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [PSH04-44] - BSA and Azide free	HA750953	Human, Mouse, Rat	WB,IHC-P,FC
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [SY28-05]	ET1607-73	Human, Mouse, Rat, Dog	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [SY28-05] - BSA and Azide free	HA750132	Human, Mouse, Rat, Dog	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [PSH04-44]	HA722129	Human, Mouse, Rat	WB,IHC-P,FC
Cell treatment	Phospho-AKT (T308) Recombinant Rabbit mAb [PSH08-10] - BSA and Azide free	HA751205	Human	WB,IF-Cell,FC
Cell treatment	Phospho-AKT (T308) Recombinant Rabbit mAb [PSH08-10]	HA722951	Human	WB,IF-Cell,FC
Cell treatment	Phospho-Cytosolic Phospholipase A2 (S505) Recombinant Rabbit mAb [PSH07-56]	HA722849	Human, Mouse	WB,Dot Blot
Cell treatment	Phospho-Cytosolic Phospholipase A2 (S505) Recombinant Rabbit mAb [PSH07-16]	HA722804	Human, Mouse	WB
Cell treatment	Phospho-Cytosolic Phospholipase A2 (S505) Recombinant Rabbit mAb [PSH07-56] - BSA and Azide free	HA751176	Human, Mouse	WB,Dot Blot
Cell treatment	Phospho-Cytosolic Phospholipase A2 (S505) Recombinant Rabbit mAb [PSH07-16] - BSA and Azide free	HA751137	Human, Mouse	WB
	Phospho-EGFR (S695) Recombinant Rabbit mAb [SN07-36]	ET1611-50	Human	WB,IF-Cell,IF-Tissue
Cell treatment	Phospho-EGFR (T669) Recombinant Rabbit mAb [PSH10-03]	HA723172	Human, Mouse, Rat	WB,IP
Cell treatment	Phospho-EGFR (T669) Recombinant Rabbit mAb [PSH10-03] - BSA and Azide free	HA751333	Human, Mouse, Rat	WB,IP
Cell treatment	Phospho-EGFR (Y1045) Recombinant Rabbit mAb [PSH03-28] - BSA and Azide free	HA750868	Human	WB,IF-Cell,FC
Cell treatment	Phospho-EGFR (Y1045) Recombinant Rabbit mAb [PSH03-28]	HA721970	Human	WB,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-EGFR (Y1068) Recombinant Rabbit mAb [SD2055] - BSA and Azide free	HA750290	Human	WB,IF-Cell,IF-Tissue,IP,FC
Cell treatment	Phospho-EGFR (Y1068) Recombinant Rabbit mAb [SD2055]	ET1612-30	Human	WB,IF-Cell,IF-Tissue,IP,FC
Cell treatment	Phospho-EGFR (Y1092) Recombinant Rabbit mAb [SJ0194]	ET1606-44	Human, Mouse	WB,IF-Tissue,IHC-P
	Phospho-EGFR (Y1148) Rabbit pAb	HA500580	Human	WB,IF-Cell
Cell treatment	Phospho-EGFR (Y1173) Recombinant Rabbit mAb [SC57-04]	ET1610-4	Human, Mouse	WB,IF-Cell,IHC-P,IP
Cell treatment	Phospho-EGFR (Y845) Recombinant Rabbit mAb [JE46-77]	HA721488	Human	WB
Cell treatment	Phospho-ELK1 (S383) Recombinant Rabbit mAb [PSH14-65] - BSA and Azide free	HA751526	Human, Mouse	WB,IHC-P,IP
	Phospho-ELK1 (S383) Rabbit pAb	ER65967	Human, Mouse, Rat	WB,IHC
Cell treatment	Phospho-ELK1 (S383) Recombinant Rabbit mAb [PSH14-65]	HA723663	Human, Mouse	WB,IHC-P,IP
Cell treatment	Phospho-Erk1 (T202)+Erk2 (T185) Recombinant Rabbit mAb [SZ2-4]	ET1603-22	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
Cell treatment	Phospho-Erk1 (T202)+Erk2 (T185) Recombinant Rabbit mAb [SZ2-4] - BSA and Azide free	HA750061	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
Cell treatment	Phospho-Erk1 (T202 + Y204) + Erk2 (T185 + Y187) Recombinant Rabbit mAb [SC58-01] - BSA and Azide free	HA750199	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP
Cell treatment	Phospho-Erk1 (T202 + Y204) + Erk2 (T185 + Y187) Recombinant Rabbit mAb [SC58-01]	ET1610-13	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP
Cell treatment	Phospho-FOXO1A (T24)+ FOXO3A (T32)+ FOXO4 (T28)+ FOXO6 (T26) Recombinant Rabbit mAb [PSH06-00] - BSA and Azide free	HA751040	Human, Mouse	WB,IF-Cell,FC
Cell treatment	Phospho-FOXO1A (T24)+ FOXO3A (T32)+ FOXO4 (T28)+ FOXO6 (T26) Recombinant Rabbit mAb [PSH06-00]	HA722499	Human, Mouse	WB,IF-Cell,FC
Cell treatment	Phospho-FOXO3a (S253) Recombinant Rabbit mAb [ST49-01]	ET1609-49	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-FOXO3a (S253) Recombinant Rabbit mAb [ST49-01] - BSA and Azide free	HA750178	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-JNK1/2/3 (T183 + T183 + T221) Recombinant Rabbit mAb [ST500] - BSA and Azide free	HA750173	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
Cell treatment	Phospho-JNK1/2/3 (T183 + T183 + T221) Recombinant Rabbit mAb [ST500]	ET1609-42	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Phospholipase A2 Recombinant Rabbit mAb [JE41-02]	HA721565	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-MEK1/2 (S218 + S222) Recombinant Rabbit mAb [ST0490]	ET1609-50	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Phospho-MEK1 (S298) Recombinant Rabbit mAb [SD206-7]	ET1612-40	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	Phospho-MEK1 (T292) Recombinant Rabbit mAb [SD2088]	ET1612-42	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-Met (Y1234/Y1235) Recombinant Rabbit mAb [PSH06-90] - BSA and Azide free	HA751114	Human	WB
Cell treatment	Phospho-Met (Y1234/Y1235) Recombinant Rabbit mAb [PSH06-90]	HA722765	Human	WB
Cell treatment	Phospho-Met (Y1349) Recombinant Rabbit mAb [JE51-68]	ET7110-14	Human, Rat	WB,IHC-P
Cell treatment	Phospho-MST1/MST2 (T180 + T183) Recombinant Rabbit mAb [JE46-64]	HA721737	Human, Mouse, Rat	WB
Cell treatment	Phospho-NF-kB p65 (S529) Recombinant Rabbit mAb [SP07-00]	ET1604-27	Human	WB,IHC-P,FC,IP
Cell treatment	Phospho-NF-kB p65 (S536) Recombinant Rabbit mAb [PSH10-58] - BSA and Azide free	HA751361	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-NF-kB p65 (S536) Recombinant Rabbit mAb [PSH10-58]	HA723223	Human, Mouse, Rat	WB,IF-Cell
	Phospho-PAK1 (S144)+PAK2 (S141)+PAK3 (S139) Recombinant Rabbit mAb [SU03-70]	ET1608-58	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
Cell treatment	Phospho-PAK4+PAK5+PAK6 (S474 + S560 + S602) Recombinant Rabbit mAb [JE57-42]	HA722146	Human	WB
Cell treatment	Phospho-PDGFR beta (Y1009) Recombinant Rabbit mAb [JE46-72]	HA721736	Human, Mouse, Rat	WB
Cell treatment	Phospho-PDGFR beta (Y1021) Recombinant Rabbit mAb [JE48-54]	HA721833	Mouse, Rat, Human	WB
Cell treatment	Phospho-PDGFR beta (Y740) Recombinant Rabbit mAb [SN0646]	ET1611-29	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue
Cell treatment	Phospho-PDGFR beta (Y751) Recombinant Rabbit mAb [PSH13-61] - BSA and Azide free	HA751480	Mouse, Rat, Human	WB
Cell treatment	Phospho-PDGFR beta (Y751) Recombinant Rabbit mAb [PSH13-61]	HA723550	Mouse, Rat, Human	WB
Cell treatment	Phospho-PI3K p85 (Y467) + PI3K p55 (Y199) Recombinant Rabbit mAb [PSH01-39]	HA721673	Human, Mouse	WB,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-PI3K p85 (Y467) + PI3K p55 (Y199) Recombinant Rabbit mAb [PSH01-38]	HA721672	Human, Mouse	WB,IF-Cell
Cell treatment	Phospho-PI3K p85 (Y467) + PI3K p55 (Y199) Recombinant Rabbit mAb [PSH01-38] - BSA and Azide free	HA750716	Human, Mouse	WB,IF-Cell
Cell treatment	Phospho-PI3K p85 (Y467) + PI3K p55 (Y199) Recombinant Rabbit mAb [PSH01-39] - BSA and Azide free	HA750717	Human, Mouse	WB,IF-Cell
Cell treatment	Phospho-PKA alpha/beta/gamma (catalytic subunit) (T197) Recombinant Rabbit mAb [JE45-17]	HA721864	Human, Mouse, Rat	WB
Cell treatment	Phospho-PKC alpha/beta II (T638/T641) Recombinant Rabbit mAb [PSH15-91]	HA723796	Human, Mouse, Rat	WB,IF-Cell,Dot Blot
Cell treatment	Phospho-PKC alpha/beta II (T638/T641) Recombinant Rabbit mAb [PSH16-27] - BSA and Azide free	HA751595	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,Dot Blot
Cell treatment	Phospho-PKC alpha/beta II (T638/T641) Recombinant Rabbit mAb [PSH16-27]	HA723825	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,Dot Blot
Cell treatment	Phospho-PKC alpha/beta II (T638/T641) Recombinant Rabbit mAb [PSH15-91] - BSA and Azide free	HA751582	Human, Mouse, Rat	WB,IF-Cell,Dot Blot
Cell treatment	Phospho-PKC alpha (T497) Recombinant Rabbit mAb [JE44-81]	HA722278	Human, Mouse, Rat	WB,IF-Cell,FC
Cell treatment	Phospho-PKC alpha (T638) Recombinant Rabbit mAb [JF0964]	ET1702-17	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,Dot Blot
Cell treatment	Phospho-PKC delta (Y311) Recombinant Rabbit mAb [JE44-20]	HA722458	Human, Mouse, Rat	WB
Cell treatment	Phospho-PKC (S729) Recombinant Rabbit mAb [JE77-54]	HA722846	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-PKC zeta (T560) Recombinant Rabbit mAb [SN0753]	ET1611-78	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-PLC-gamma-1 (S1248) Recombinant Rabbit mAb [PSH13-62] - BSA and Azide free	HA751481	Human	WB,IP
Cell treatment	Phospho-PLC-gamma-1 (S1248) Recombinant Rabbit mAb [PSH13-62]	HA723551	Human	WB,IP
Cell treatment	Phospho-PLC-gamma-1 (Y783) Recombinant Rabbit mAb [PSH06-29] - BSA and Azide free	HA751066	Human, Mouse, Rat	WB
Cell treatment	Phospho-PLC-gamma-1 (Y783) Recombinant Rabbit mAb [PSH06-29]	HA722661	Human, Mouse, Rat	WB
Cell treatment	Phospho-Rac1+Cdc42 (S71) Recombinant Rabbit mAb [JJ086-9]	ET1701-40	Human	WB,IF-Cell,IHC-P
Cell treatment	Phospho-Raf1 (S259) Recombinant Rabbit mAb [SD85-07] - BSA and Azide free	HA750304	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-Raf1 (S259) Recombinant Rabbit mAb [SD85-07]	ET1612-87	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	Phospho-Raf1 (S43) Recombinant Rabbit mAb [SD087-05]	ET1612-75	Human, Mouse, Rat	WB
	Phospho-Raf1 (S621) Recombinant Rabbit mAb [JJ085-05]	ET1701-3	Human, Rat, Mouse	WB,IHC-P
Cell treatment	Phospho-SHP2 (Y542) Recombinant Rabbit mAb [SN61-01]	ET1611-22	Mouse, Human	WB,IP
Relative expression	Phospho-TrkB (Y817) Recombinant Rabbit mAb [SC0556] - BSA and Azide free	HA750213	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
Relative expression	Phospho-TrkB (Y817) Recombinant Rabbit mAb [SC0556]	ET1610-35	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	PI3 Kinase p110 beta Recombinant Rabbit mAb [JE65-38]	HA722474	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	PI 3 Kinase p110 delta Recombinant Rabbit mAb [JM83-24]	ET1705-46	Human, Mouse, Rat	WB,IHC-P
KD	PI3-kinase p110 subunit beta Recombinant Rabbit mAb [SC0580]	ET1610-36	Human	WB,IP,IHC-P,FC
	PI3 Kinase p110 α Recombinant Rabbit mAb [SJ0186]	ET1606-36	Human, Mouse	WB,IF-Cell,IF-Tissue,IP,IHC-P
	PI 3 Kinase p55 gamma Recombinant Rabbit mAb [SR44-06]	ET1602-27	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	PI 3 Kinase p85 alpha Rabbit pAb	R1312-6	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	PI 3 Kinase p85 alpha Recombinant Rabbit mAb [SU04-07] - BSA and Azide free	HA750162	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	PI 3 Kinase p85 alpha Recombinant Mouse mAb [A3-D0-R]	HA601206	Human	WB,IF-Cell
KD	PI 3 Kinase p85 alpha Mouse mAb [L3-D6]	EM1701-62	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue
KD	PI 3 Kinase p85 alpha Recombinant Mouse mAb [A3-D0-R] - BSA and Azide free	HA610084	Human	WB,IF-Cell
KD	PI 3 Kinase p85 alpha Mouse mAb [A3-D0]	M1510-2	Human	WB,IF-Cell
	PI 3 Kinase p85 alpha Recombinant Rabbit mAb [SU04-07]	ET1608-70	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	PI3 Kinase p85 alpha Recombinant Rabbit mAb [JE04-11]	HA722522	Human, Mouse, Rat	WB
	PI 3 Kinase p85 alpha/gamma Rabbit pAb	ER64588	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Tissue
	PI 3 Kinase p85 beta Recombinant Rabbit mAb [ST04-77]	ET1609-30	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	PKA alpha/beta cat Rabbit pAb	ER64617	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	PKA alpha/beta/gamma cat Rabbit pAb	ER64618	Human, Mouse, Rat	WB,IHC-P,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	PKA beta cat Rabbit pAb	ER64619	Human, Mouse, Rat	WB,IHC-P
	PKA C alpha Rabbit pAb	ER1706-65	Human, Mouse, Rat	IF-Cell,IHC-P,FC
	PKA C-alpha Recombinant Rabbit mAb [PSH08-57]	HA723003	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC,IP
	PKA C-alpha Recombinant Rabbit mAb [PSH08-57] - BSA and Azide free	HA751232	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC,IP
	PKC Rabbit pAb	ER64620	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
KD	PKC alpha Recombinant Rabbit mAb [SU31-08]	ET1608-15	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	PKC beta Rabbit pAb	ER1802-14	Human, Mouse	WB,IF-Cell,IHC-P
	PKC beta 1 Recombinant Rabbit mAb [JM43-26]	ET1705-30	Human, Mouse, Rat	WB,IHC-P,FC
	PKC beta 2 Recombinant Rabbit mAb [ST48-06]	ET1609-44	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	PKC delta Recombinant Rabbit mAb [JJ093-02]	ET1701-85	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	PKC delta Recombinant Rabbit mAb [JJ093-02] - BSA and Azide free	HA750325	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	PKC epsilon Recombinant Rabbit mAb [JA38-21]	ET1704-85	Human, Mouse, Rat	WB,IHC-P
	PKC eta Recombinant Rabbit mAb [JE32-41]	HA722389	Human, Mouse, Rat, Monkey	WB
	PKC gamma Recombinant Rabbit mAb [JG69-32]	ET7108-42	Human, Mouse, Rat	WB,IHC-P
	PKC theta Rabbit pAb	ER64621	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Relative expression	PKC zeta Recombinant Rabbit mAb [PSH06-74]	HA722773	Human	WB
Relative expression	PKC zeta Recombinant Rabbit mAb [PSH06-74] - BSA and Azide free	HA751122	Human	WB
	PLA1A Rabbit pAb	ER64627	Human	WB,IHC-P
	PLA2G12A Mouse mAb [C11-H2]	EM1706-96	Human	WB,IHC-P,FC
	PLCE1 Rabbit pAb	ER61619	Human	IHC-P
	PLCG1 Mouse mAb [D11-C2]	EM1707-02	Human, Rat, Monkey	WB,IF-Cell,IHC-P,FC
	PLCG 2 Recombinant Rabbit mAb [JE39-12]	HA721477	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	PLCG2 Rabbit pAb	R1512-4	Human	WB,IF-Cell,IHC-P,FC
	PLC-gamma-1 Recombinant Rabbit mAb [SN06-10]	ET1611-25	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
	PLD2 Mouse mAb [D9-E4]	EM1707-42	Human	IF-Cell,IHC-P,FC
	PLGF Mouse mAb [3G3]	EM1701-88	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	PLGF Mouse mAb [3G2]	EM1701-90	Human, Mouse, Rat	WB,IHC-P,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	PLGF Recombinant Rabbit mAb [JA63-15]	ET1704-99	Human, Mouse, Rat	WB,IHC-P
	PLGF Rabbit pAb	ER1706-88	Human, Mouse	IHC-P,FC
	PRKACG Mouse mAb [B3-E4]	EM1709-59	Human	WB,IHC-P,FC
KD	Rab5 Recombinant Rabbit mAb [ST47-04] - BSA and Azide free	HA750167	Human, Mouse, Rat	WB,IHC-P,IF-Cell
KD	Rab5 Recombinant Rabbit mAb [ST47-04]	ET1609-27	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Rab5a Mouse mAb [C11-H5]	EM1708-08	Human	WB,IF-Cell,IHC-P,FC
	RAB5B Rabbit pAb	ER61113	Human, Mouse	WB
KD	RAB5C Recombinant Rabbit mAb [JE63-14]	HA721012	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	RAC1 Recombinant Rabbit mAb [PSH05-09] - BSA and Azide free	HA750982	Human, Rat	WB,IF-Cell
	RAC1 Recombinant Rabbit mAb [PSH05-09]	HA722225	Human, Rat	WB,IF-Cell
	Rac1-2-3 Recombinant Rabbit mAb [JM11-29]	ET1703-80	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	RAC1/2/3 Rabbit pAb	HA500424	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	RAC2 Rabbit pAb	ER61092	Human, Mouse	WB
	RAC3 Rabbit pAb	ER61093	Human, Mouse	WB
KD	RAC3 Recombinant Rabbit mAb [JE56-48]	HA721890	Human, Mouse, Rat	WB
KD	Raf1 Recombinant Rabbit mAb [JJ086-06]	ET1701-21	Human, Mouse, Rat	WB
KD	RALA Rabbit pAb	HA500075	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	RALA Mouse mAb	EM1707-77	Human, Mouse, Monkey	WB,IF-Cell,IHC-P,FC
KD	RALA Recombinant Rabbit mAb [JE54-12]	HA721215	Human	WB,IHC-P,FC
	RALB Mouse mAb [B4-F8]	EM1707-30	Human, Mouse, Monkey	WB,IHC-P,FC
KD	RALBP1 Recombinant Rabbit mAb [SC06-49]	ET1610-88	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
	RAP1A + RAP1B Recombinant Rabbit mAb [JE35-61]	HA721335	Human, Mouse	WB,FC
	RAP1B Rabbit pAb	ER61929	Human, Rat, Mouse	WB
	Ras Recombinant Rabbit mAb [JE44-13]	HA721883	Human, Mouse, Rat	WB,IF-Cell,FC
	Ras Recombinant Rabbit mAb [JF10-11]	ET1702-94	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IP,FC
	Ras Rabbit pAb	ER40115	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Ras Recombinant Rabbit mAb [SA39-05]	ET1601-16	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,FC,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	RASA1 Recombinant Rabbit mAb [JE60-89]	HA720052	Human, Mouse, Rat	WB,IHC-P
	RASA2 Rabbit pAb	ER61669	Human, Mouse, Rat	WB
	Repac Rabbit pAb	ER64789	Human, Mouse, Rat	WB
	RhoA Recombinant Rabbit mAb [SN0612] - BSA and Azide free	HA750242	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	RhoA Recombinant Rabbit mAb [SN0612]	ET1611-10	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Rho A+B+C Recombinant Rabbit mAb [SR38-00]	ET1602-23	Human, Mouse, Rat	WB,IF-Cell,FC
	Rho B Rabbit pAb	ER64810	Human, Mouse, Rat	WB
	RIN1 Rabbit pAb	ER64867	Human, Mouse, Rat	WB
	RRAS Rabbit pAb	HA500397	Human, Mouse, Rat, Monkey	WB
	RRAS2 Rabbit pAb	HA500379	Human, Mouse, Rat	WB,IHC-P
	SCF Recombinant Rabbit mAb [SC61-06]	ET1610-77	Human, Mouse, Rat	WB,IF-Tissue,IHC-P
	SCF Mouse mAb [4G5]	RT1553	Human, Mouse, Rat	WB,IP,IF,IHC-P
	Sck Rabbit pAb	ER64915	Human, Rat	WB
	SHC Recombinant Rabbit mAb [JE57-45]	HA722564	Human, Mouse, Rat	WB,IHC-P
	SHC1 Mouse mAb [7C-7F]	EM1706-40	Human, Mouse	WB,IF-Cell,FC
	SHC4 Rabbit pAb	ER1916-54	Mouse	WB
	SHP2 Rabbit pAb	ER1706-93	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	SHP2 Recombinant Rabbit mAb [SN72-02]	ET1611-35	Human, Mouse, Rat	WB,IP,IHC-P
KD	SHP2 Recombinant Rabbit mAb [SN07-32]	ET1611-80	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP
	SOS1 Recombinant Rabbit mAb [JE30-06]	HA720077	Human, Mouse, Rat	WB,IHC-P
	Sos 2 Rabbit pAb	ER64998	Human, Mouse	WB,IHC-P
	SR-7 Rabbit pAb	ER65021	Human, Mouse, Rat	IF-Cell
	STK4 / MST-1 Rabbit pAb	HA500043	Human, Mouse, Rat	WB,IHC-P
	STK4 / MST-1 Rabbit pAb	HA500031	Human, Mouse	WB,IHC-P
	SynGAP Recombinant Rabbit mAb [JE37-00]	HA722728	Human, Mouse, Rat	WB,IHC-P
	SynGAP Rabbit pAb	ER1902-57	Human, Mouse, Rat	WB,IHC-P,FC
	SynGAP Recombinant Rabbit mAb [JE57-79]	HA721860	Mouse, Rat	WB,IHC-P,IF-Tissue
	TGF alpha Recombinant Rabbit mAb [JB58-36]	ET7107-40	Human	WB,IF-Tissue,IHC-P,IP
	TIAM1 Rabbit pAb	ER62423	Human, Mouse	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	TIE2 Recombinant Rabbit mAb [PSH06-63]	HA722732	Human	WB,FC
Relative expression	TIE2 Recombinant Rabbit mAb [PSH06-63] - BSA and Azide free	HA751089	Human	WB,FC
Relative expression	TIE2 Rabbit pAb	ER65864	Human	WB,IF-Cell,FC
Relative expression	TrkA Recombinant Rabbit mAb [SU0354]	ET1608-44	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Relative expression	TrkA Recombinant Rabbit mAb [SU0354] - BSA and Azide free	HA750148	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	TrkA Rabbit pAb	ER1706-86	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Relative expression	TrkA+B+C Recombinant Rabbit mAb [JJ084-04]	ET1701-16	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Relative expression	TrkB Recombinant Rabbit mAb [JE58-22]	ET7111-38	Mouse, Rat	WB,IHC-P
	Trk C Rabbit pAb	ER65473	Human	WB
Relative expression	TrkC Recombinant Rabbit mAb [PSH05-40] - BSA and Azide free	HA750999	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
Relative expression	TrkC Recombinant Rabbit mAb [PSH05-40]	HA722269	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
KD	VEGF Rabbit pAb	ER30607	Human	WB,IF-Cell,IHC-P,FC
	VEGF Recombinant Rabbit mAb [SP07-01] - BSA and Azide free	HA750074	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	VEGF Recombinant Rabbit mAb [SP07-01]	ET1604-28	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	VEGFA Recombinant Rabbit mAb [PSH19-17]	HA724051	Human, Mouse, Rat, Monkey	WB
	VEGFB167 Rabbit pAb	ER1918-50	Rat	WB
	VEGFC Rabbit pAb	R1307-4	Human	WB
	VEGF-D Rabbit pAb	ER65298	Human, Mouse, Rat, Monkey	WB,IHC-P
	VEGFD Recombinant Rabbit mAb [JE54-75]	ET7110-81	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	VEGF Receptor 1 Rabbit pAb	RT1651	Human, Mouse, Rat	WB,IP,IF,IHC-P
KD	VEGF Receptor 1 Rabbit pAb	ER50803	Human, Mouse, Rat	WB,IHC-P
KD	VEGF Receptor 1 Recombinant Rabbit mAb [SY09-09]	ET1605-11	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	VEGF Receptor 1 Recombinant Rabbit mAb [SY09-09] - BSA and Azide free	HA750081	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Relative expression	VEGF Receptor 2 Recombinant Rabbit mAb [SU03-42]	ET1608-33	Human, Rat	WB,IF-Cell,IHC-P,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	ZAP70 Recombinant Rabbit mAb [PSH11-52]	HA723344	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Relative expression	ZAP70 Recombinant Rabbit mAb [JU08-39]	ET1706-42	Human	WB,IF-Cell,IHC-P,FC,IP
Relative expression	ZAP70 Recombinant Mouse mAb [A1B5-R]	HA601218	Human, Mouse, Rat	WB,IHC-P
Relative expression	ZAP70 Mouse mAb [A1B7]	EM1902-28	Human, Mouse, Rat	WB,IF-Cell,FC
Relative expression	ZAP70 Recombinant Mouse mAb [A1B5-R] - BSA and Azide free	HA610096	Human, Mouse, Rat	WB,IHC-P
Relative expression	ZAP70 Mouse mAb [A1B5]	EM1901-27	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层