

Sirtuin信号通路研究-抗体清单

高级验证	产品名	货号	反应种属	验证应用
	5 Lipoxygenase Recombinant Rabbit mAb [PSH18-05]	HA723968	Human, Mouse, Rat	WB
	5 Lipoxygenase Rabbit pAb	R1512-14	Human, Mouse, Rat	WB
	ABCA1 Rabbit pAb	R1510-42	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	ABCA1 Mouse mAb [4G2]	RT1009	Human, Mouse	WB,IP,IF,IHC-P
KD	ADAM10 Recombinant Rabbit mAb [JM32-11]	ET1703-60	Human, Mouse, Rat	WB,IHC-P,IP
	ADAM10 Rabbit pAb	ER1706-97	Human	WB,IHC-P
	ADAM17 Recombinant Rabbit mAb [JM10-35]	ET1703-06	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IF-Tissue,FC,IP,IHC-Fr,IHC-P
	ADAM17 Recombinant Rabbit mAb [JM10-35] - BSA and Azide free	HA750370	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IF-Tissue,FC,IP,IHC-Fr,IHC-P
	ADAM9 Rabbit pAb	ER60448	Human	WB
	Anoctamin-9 Rabbit pAb	ER62673	Human	IHC-P,IF-Cell
KD	Arg2 Recombinant Rabbit mAb [JE54-97]	ET7110-91	Human	WB,IHC-P
	Arginase1 Recombinant Rabbit mAb	IRS081RB	Human	mIHC
	Arginase1 Recombinant Rabbit mAb	IRS082RB	Mouse	mIHC
	Arginase II Rabbit pAb	ER1903-91	Human, Mouse	WB,IHC-P,IF-Cell,IF-Tissue,FC
Relative expression	ATP2A1 / SERCA1 Rabbit pAb	ER1803-34	Human, Mouse, Rat	WB,IHC-P,FC
	ATP2A1 / SERCA1 Recombinant Rabbit mAb [JE45-75]	ET7109-43	Human	WB,IP,IHC-P
	ATP8A1 Rabbit pAb	ER1902-61	Human, Rat	WB,FC
Relative expression	Axl Rabbit pAb	R1406-3	Human, Mouse	WB,IF-Cell,IHC-P,FC
Relative expression	Axl Recombinant Rabbit mAb [PSH06-68] - BSA and Azide free	HA751097	Human, Monkey	WB,IF-Cell,FC
Relative expression	Axl Recombinant Rabbit mAb [PSH06-68]	HA722740	Human, Monkey	WB,IF-Cell,FC
	Axl Mouse mAb [C4-A8]	M1310-5	Human	WB,IF-Cell
	BAI1 Rabbit pAb	ER60093	Human, Rat, Mouse	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	BCAR1 Recombinant Rabbit mAb [JE36-56]	HA721465	Human, Mouse, Rat	WB
	C1QA Mouse mAb [G3-C7]	EM1707-50	Human	IHC-P,FC
	C1QA Recombinant Rabbit mAb [JE66-59]	HA721439	Human	WB,IHC-P
	C1q-B Rabbit pAb	ER62847	Human	WB,IHC-P,IF-Cell
	Calreticulin Recombinant Antibody [SU37-03] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610363	Human, Mouse, Rat, Monkey	WB,IHC-P
KD	Calreticulin Recombinant Rabbit mAb [SU37-03]	ET1608-60	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
	Calreticulin Recombinant Antibody [SU37-03] - Mouse IgG1 (Chimeric)	HA601547	Human, Mouse, Rat, Monkey	WB,IHC-P
KD	Calreticulin Recombinant Rabbit mAb [SU37-03] - BSA and Azide free	HA750156	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
	CAMK2D Rabbit pAb	R1305-2	Human	WB,IF-Cell,IHC-P
	CAMK2G Mouse mAb [10C-G1]	EM1706-16	Human, Rat	WB,IHC-P,FC
	CaMKII Recombinant Rabbit mAb [SU03-57]	ET1608-47	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr,FC
	CaMKII Recombinant Rabbit mAb [SU03-57] - BSA and Azide free	HA750150	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr,FC
	CaMKII alpha Recombinant Rabbit mAb [JM11-07]	ET1703-54	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr,IF-Cell,FC,IP
	CaMKII alpha/beta/delta Rabbit pAb	ER62885	Human, Mouse, Rat	IHC-P,WB
	CaMKII alpha/delta Rabbit pAb	ER62886	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	CaMKII beta/gamma Rabbit pAb	ER62887	Human, Mouse, Rat	WB,IF-Cell
	CaMKII beta/gamma/delta Rabbit pAb	ER62888	Human, Mouse, Rat	WB,IHC-P
KD	CaMKII delta Recombinant Rabbit mAb [JA30-07]	ET1704-37	Human, Mouse, Rat, Zebrafish	WB,IHC-P
	Caspase-1 Rabbit pAb	R1510-23	Human	IF-Cell,IHC-P,FC
	Caspase-1 Recombinant Rabbit mAb [SU40-07]	ET1608-69	Human, Mouse, Rat	WB,FC,IF-Cell
	Caspase-1 Recombinant Rabbit mAb [SU40-07] - BSA and Azide free	HA750161	Human, Mouse, Rat	WB,FC,IF-Cell
	Caspase-1 Rabbit pAb	R1510-22	Human	IHC-P,FC
Relative expression	Caspase-1 + p10 + p12 Recombinant Rabbit mAb [JE56-35]	HA722222	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Caspase-1 p20 Rabbit pAb	ER1905-47	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment, KD	Caspase-3 p12 Recombinant Rabbit mAb [JE07-51]	HA722235	Human, Mouse, Rat	WB,IHC-P
	Caspase-7 Recombinant Rabbit mAb [JE59-36]	HA722317	Mouse, Rat	WB,IHC-P,IF-Tissue
Cell treatment	Caspase-7 Recombinant Rabbit mAb [PSH13-69] - BSA and Azide free	HA751482	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Caspase-7 Recombinant Rabbit mAb [PSH13-69]	HA723558	Human, Mouse, Rat	WB,IHC-P
KD	CD147 Recombinant Rabbit mAb [JF1-045]	ET1702-58	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	CD147 Rabbit pAb	1003-10	Mouse, Rat	WB,IF-Cell,IHC-P,FC
	CD172b Rabbit pAb	ER62937	Human	WB
	CD172g Rabbit pAb	ER62938	Human, Monkey	WB
	CD24 Rabbit pAb	R1303-1	Human	WB,IHC-P
Cell treatment, Relative expression	CD24 Recombinant Rabbit mAb [PSH10-85] - BSA and Azide free	HA751377	Human	IF-Cell,FC,WB
	CD24 Rabbit pAb	0804-4	Human	WB,IF-Cell,FC
	CD24 Rabbit pAb	0804-3	Mouse, Human, Rat	WB,IHC-P,FC
Cell treatment, Relative expression	CD24 Recombinant Rabbit mAb [PSH10-85]	HA723252	Human	IF-Cell,FC,WB
	CD31 Recombinant Antibody [7-A1-R] - Rat IgG1 (Chimeric)	HA601522	Human	IHC-P
	CD31 Recombinant Rabbit mAb [PSH20-35]	HA724146	Human, Mouse, Rat	IHC-P
	CD31 Recombinant Rabbit mAb [SU03-59] - BSA and Azide free	HA750151	Human	WB,IHC-P,FC,IP,IF-Cell,IF-Tissue
	CD31 Recombinant Antibody - Rabbit IgG (Chimeric)	IRS116RB	Human	mIHC
Relative expression	CD31 Recombinant Antibody [7-A1-R] - Rabbit IgG (Chimeric) - BSA and Azide free	HA751615	Human	WB,IF-Cell,IHC-P,FC
	CD31 Recombinant Rabbit mAb [SU03-59]	ET1608-48	Human	WB,IHC-P,FC,IP,IF-Cell,IF-Tissue
	CD31 Mouse mAb	IRS023	Human	mIHC
Relative expression	CD31 Recombinant Mouse mAb [7-A1-R]	HA601240	Human	WB,IF-Cell,IHC-P
Relative expression	CD31 Mouse mAb [7-A1]	M1511-8	Human, Mouse	WB,IF-Cell,IHC-P,mIHC
	CD31 Recombinant Rabbit mAb	IRS117RB	Mouse	mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	iFluor™ 488 Conjugated CD31 Recombinant Rabbit mAb [SU03-59]	HA720220F	Human	IF-Tissue
Relative expression	CD31 Rabbit pAb	ER31219	Human, Mouse, Rat	WB,IF-Cell,FC,IHC-P
	CD31 Recombinant Rabbit mAb	IRS075	Mouse	mIHC
Relative expression	CD31 Recombinant Antibody [7-A1-R] - Rabbit IgG (Chimeric)	HA723897	Human	WB,IF-Cell,IHC-P,FC
Relative expression	CD31 Recombinant Mouse mAb [7-A1-R] - BSA and Azide free	HA610115	Human	WB,IF-Cell,IHC-P
	CD31 Recombinant Antibody [7-A1-R] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610349	Human	IHC-P
Relative expression	CD36 Recombinant Rabbit mAb [PSH18-47]	HA724000	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Cell treatment	CD36 Recombinant Rabbit mAb [JJ2005] - BSA and Azide free	HA750307	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	CD36 Rabbit pAb	ER1802-7	Human, Mouse, Rat	IF-Cell,IHC-P,FC
Cell treatment	CD36 Recombinant Rabbit mAb [JJ2005]	ET1701-24	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	CD36 Recombinant Rabbit mAb	IRS118RB	Human	mIHC
	CD47 Mouse mAb [1A9]	EM1701-38	Human	IF-Cell,IHC-P,FC
Relative expression	CD47 Recombinant Rabbit mAb [PSH10-31]	HA723198	Human	WB,IF-Cell,FC
	iFluor™ 488 Conjugated CD47 Mouse mAb [1B1]	HA600104F	Human	FC
Relative expression	CD47 Recombinant Rabbit mAb [PSH10-31] - BSA and Azide free	HA751348	Human	WB,IF-Cell,FC
	CD47 Recombinant Rabbit mAb [PD00-33]	HA721158	Human	IHC-P,FC
	CDC42 Recombinant Rabbit mAb [JJ086-04] - BSA and Azide free	HA750320	Human, Mouse, Rat	WB,FC,IF-Cell
	CDC42 Recombinant Rabbit mAb [JJ086-04]	ET1701-7	Human, Mouse, Rat	WB,FC,IF-Cell
	C/EBP Beta Recombinant Rabbit mAb [SC55-05] - BSA and Azide free	HA750235	Human, Rat, Mouse	WB,FC,IP
	C/EBP Beta Recombinant Rabbit mAb [SC55-05]	ET1610-9	Human, Rat, Mouse	WB,FC,IP
	CH25H Rabbit pAb	ER1919-26	Human	IHC-P,IF-Tissue
	CH25H Rabbit pAb	ER1919-29	Human, Rat	WB,IHC-P,IF-Tissue
Cell treatment	Cleaved Caspase-3 Recombinant Rabbit mAb [SR01-02] - BSA and Azide free	HA750053	Human, Pig	WB,IF-Cell,IHC-P
Cell treatment	Cleaved Caspase-3 Recombinant Rabbit mAb [SR01-02]	ET1602-47	Human, Pig	WB,IF-Cell,IHC-P
	Cleaved Caspase-3 Rabbit pAb	IRS227RB	Mouse	mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Cleaved Caspase-3 Recombinant Rabbit mAb	IRS136RB	Human	mIHC
	Cleaved Caspase-7 (Ser199) Rabbit pAb	ER60002	Human, Mouse	WB,IHC-P,IF-Cell
Cell treatment, KO	Cleaved+pro Caspase-3 Recombinant Rabbit mAb [SU38-04] - BSA and Azide free	HA750159	Human, Mouse, Rat	WB,IHC-P,IP
Cell treatment, KD	Cleaved+pro Caspase-3 Recombinant Rabbit mAb [JE75-05]	HA722367	Human, Mouse	WB
Cell treatment, KO	Cleaved+pro Caspase-3 Recombinant Rabbit mAb [SU38-04]	ET1608-64	Human, Mouse, Rat	WB,IHC-P,IP
KD	COX2 Mouse mAb [A3F7]	EM1902-12	Human, Rat	WB,IHC-P,FC
	COX2 Mouse mAb [A3F8]	EM1902-11	Human, Rat	WB,IHC-P,FC
	COX2 Recombinant Rabbit mAb [6J8]	RE6064	Human, Mouse	WB,ICC
	COX2 Mouse mAb	RT1159	Human, Mouse, Rat	WB,IP,IF,IHC-P,ELISA
KD	COX2/Cyclooxygenase 2 Recombinant Rabbit mAb [SC56-06]	ET1610-23	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
KD	COX2/Cyclooxygenase 2 Recombinant Rabbit mAb [SC56-06] - BSA and Azide free	HA750207	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	CPT1A Recombinant Rabbit mAb [PSH04-00] - BSA and Azide free	HA750923	Human, Mouse	WB,IHC-P,IP
	CPT1A Recombinant Rabbit mAb [PSH04-00]	HA722081	Human, Mouse	WB,IHC-P,IP
	CPT1A Recombinant Rabbit mAb [PSH10-94]	HA723261	Mouse, Rat	WB
	CPT1A Recombinant Rabbit mAb [PSH10-94] - BSA and Azide free	HA751384	Mouse, Rat	WB
	CPTI-M Rabbit pAb	ER63116	Human, Mouse, Rat	WB,IHC-P
	CREB Recombinant Rabbit mAb [SA04-04] - BSA and Azide free	HA750007	Human, Mouse, Zebrafish, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	CREB Recombinant Rabbit mAb [PSH15-48] - BSA and Azide free	HA751563	Human, Mouse, Rat, Monkey	WB,IHC-P,IHC-Fr,IF-Cell,FC
	CREB Recombinant Rabbit mAb [PSH15-48]	HA723754	Human, Mouse, Rat, Monkey	WB,IHC-P,IHC-Fr,IF-Cell,FC
	CREB Recombinant Rabbit mAb [SA04-04]	ET1601-15	Human, Mouse, Zebrafish, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	CREB1 Rabbit pAb	R1309-14	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Crk-L Rabbit pAb	ER63126	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell
	CrkL Recombinant Rabbit mAb [JE50-21]	HA721740	Human, Mouse, Rat	WB,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	CrkL Rabbit pAb	HA500536	Human, Mouse, Rat	WB,IF-Cell,FC,IHC-P
	Crk p38 Recombinant Rabbit mAb [JE38-60]	HA722453	Human, Rat	WB
	CX3CL1 Rabbit pAb	HA500304	Human	WB,IHC-P
	CX3CR1 Mouse mAb [A5C7]	HA600074	Human	IF-Cell
	DBC-1 Rabbit pAb	ER63209	Human	WB,IHC-P
	Dnmt3a Rabbit pAb	ER1902-62	Human, Mouse, Rat	WB,IHC-P,FC
KD	Dnmt3a Recombinant Rabbit mAb [ST04-78]	ET1609-31	Human, Mouse, Rat	WB,IF-Tissue,IHC-P
KD	Dnmt3a Recombinant Rabbit mAb [ST04-78] - BSA and Azide free	HA750169	Human, Mouse, Rat	WB,IF-Tissue,IHC-P
	DOCK 180 Rabbit pAb	ER63276	Human, Mouse	WB
	DUSP1 Recombinant Rabbit mAb [JJ0930] - BSA and Azide free	HA750324	Human, Mouse, Rat	WB,IHC-P,IP,FC,IF-Cell,IF-Tissue
	DUSP1 Recombinant Rabbit mAb [JJ0930]	ET1701-82	Human, Mouse, Rat	WB,IHC-P,IP,FC,IF-Cell,IF-Tissue
	DUSP4 Rabbit pAb	ER1907-84	Human, Mouse	WB,FC
	DUSP5 Rabbit pAb	R1109-1	Human, Mouse	WB,IF-Cell
KD	EDG1 Recombinant Rabbit mAb [JM10-66] - BSA and Azide free	HA750376	Human, Mouse, Rat	WB,IP
KD	EDG1 Recombinant Rabbit mAb [JM10-66]	ET1703-27	Human, Mouse, Rat	WB,IP
	ELMO1 Recombinant Rabbit mAb [JG34-59]	ET7108-06	Human, Mouse, Rat	WB,IHC-P,IP
	EP4 Rabbit pAb	ER1908-19	Human, Mouse	WB,IHC-P,IF-Tissue
Cell treatment	EPO Recombinant Rabbit mAb [PSH07-33]	HA722821	Human	WB
Cell treatment	EPO Recombinant Rabbit mAb [PSH07-33] - BSA and Azide free	HA751154	Human	WB
	EPO-R Rabbit pAb	HA500529	Human	WB
	EPOR Rabbit pAb	ER1908-28	Human, Rat	WB,IHC-P,IF-Tissue
	ERK1 Mouse mAb [7-D6-E5]	M1505-9	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	ERK1 Recombinant Rabbit mAb [SP00-74]	ET1604-32	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	ERK1 Recombinant Rabbit mAb [SP05-09]	ET1604-16	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Erk1/2 Mouse mAb [2G1]	RT1484	Human, Mouse, Rat	WB,IP,IF
KD	ERK1/2 Recombinant Rabbit mAb [SA43-03] - BSA and Azide free	HA750019	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IP,FC,IHC-P,IHC-Fr
	ERK1/2 Rabbit pAb	ER131011	Human, Mouse, Rat	WB,IF-Cell,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	ERK1/2 Recombinant Rabbit mAb [SA43-03]	ET1601-29	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IP,FC,IHC-P,IHC-Fr
KD	ERK2 Recombinant Mouse mAb [5-D2-R] - BSA and Azide free	HA610080	Human, Mouse, Rat	WB
KD	ERK2 Mouse mAb [B10-8]	EM1901-53	Human, Mouse, Rat	WB,IHC-P,FC
KD	ERK2 Recombinant Mouse mAb [C3B12-R] - BSA and Azide free	HA610079	Human, Mouse, Rat	WB
KD	ERK2 Recombinant Mouse mAb [C3B12-R]	HA601201	Human, Mouse, Rat	WB
KD	ERK2 Mouse mAb [A10D2]	EM1901-54	Human, Mouse, Rat	WB,IHC-P,FC
KD	ERK2 Rabbit pAb	ER131218	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	ERK2 Recombinant Rabbit mAb [SZ25-01]	ET1603-23	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
KD	ERK2 Recombinant Mouse mAb [5-D2-R]	HA601202	Human, Mouse, Rat	WB
	ERK2 Mouse mAb [C3B12]	EM1901-55	Human, Mouse, Rat	WB,IHC-P,FC
	ERK2 Mouse mAb [5-D2]	M1407-3	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	FAK Rabbit pAb	RT1215	Human, Mouse, Rat	WB,IP,ICC,IHC-P
KD	FAK Recombinant Rabbit mAb [SR46-04] - BSA and Azide free	HA750039	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
KD	FAK Recombinant Rabbit mAb [SR46-04]	ET1602-25	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	FoxO1/3/4 Rabbit pAb	ER63465	Human, Mouse, Rat	WB,IHC-P
	iFluor™ 488 Conjugated FOXO1A Recombinant Rabbit mAb [SU33-01]	HA720161F	Human, Mouse, Rat	IF-Cell,FC
KD	FOXO1A Recombinant Rabbit mAb [SU33-01] - BSA and Azide free	HA750140	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC
KD	FOXO1A Recombinant Rabbit mAb [SU33-01]	ET1608-25	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC
	FOXO3A Rabbit pAb	ER1706-79	Human, Mouse	WB,IF-Cell,IHC-P,FC
	FOXO3A Recombinant Rabbit mAb [SP0054]	ET1604-11	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue
	FOXO4 Recombinant Rabbit mAb [JJ09-11]	ET1701-67	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,FC,IP
	FOXO4 Recombinant Rabbit mAb [SD0817]	ET1612-66	Human, Rat	WB,IHC-P
	Gas6 Rabbit pAb	ER63503	Human	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Glucose Transporter GLUT1 Rabbit pAb	ER1510-11	Human, Mouse	WB,IF-Cell,IHC-P,FC
	Glucose Transporter GLUT1 Recombinant Antibody [SA0377] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610348	Human, Mouse, Rat	IHC-Fr,IHC-P,IF-Tissue,WB
	iFluor™ 647 Conjugated Glucose Transporter GLUT1 Recombinant Rabbit mAb [SA0377]	HA720170F	Human, Mouse	IF-Tissue
	Glucose Transporter GLUT1 Recombinant Rabbit mAb	IRS151RB	Human	mIHC
KD	Glucose Transporter GLUT1 Recombinant Rabbit mAb [SA0377]	ET1601-10	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Glucose Transporter GLUT1 Recombinant Antibody [SA0377] - Rat IgG1 (Chimeric)	HA601521	Human, Mouse, Rat	IHC-Fr,IHC-P,IF-Tissue,WB
KD	Glucose Transporter GLUT1 Recombinant Rabbit mAb [SA0377] - BSA and Azide free	HA750002	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	iFluor™ 488 Conjugated Glucose Transporter GLUT1 Recombinant Rabbit mAb [SA0377]	HA720154F	Human, Mouse, Rat	IF-Cell,IF-Tissue
	GPR132 Rabbit pAb	R1308-8	Human, Mouse, Rat	WB,IF-Cell
	Group VI iPLA2 Rabbit pAb	ER63644	Human, Mouse, Rat	WB
	HIF-1 alpha Mouse mAb [5H7]	RT1278	Human, Mouse, Rat	WB,IP,ICC,IHC-P
Cell treatment	HIF-1 alpha Recombinant Rabbit mAb [JE75-33]	HA721997	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	HIF-1 alpha Recombinant Rabbit mAb [PSH07-00]	HA722778	Human, Mouse	WB,IF-Cell,IHC-P
Cell treatment	HIF-1 alpha Recombinant Rabbit mAb [PSH07-00] - BSA and Azide free	HA751127	Human, Mouse	WB,IF-Cell,IHC-P
Cell treatment	HIF-1 alpha Rabbit pAb	ER1802-41	Human	WB,IHC-P
KD	HIF-1 beta Recombinant Rabbit mAb [ST05-25] - BSA and Azide free	HA750187	Human, Mouse, Rat	WB,IHC-P
KD	HIF-1 beta Recombinant Rabbit mAb [ST05-25]	ET1609-67	Human, Mouse, Rat	WB,IHC-P
	Human CD147 Recombinant Rabbit mAb [PSH13-06] - BSA and Azide free (Capture)	HA725085	Human	ELISA(Cap)
	Human CD147 Recombinant Rabbit mAb [PSH13-07] - BSA and Azide free (Detector)	HA725086	Human	ELISA(Det)
	Biotin Conjugated Human CD147 Recombinant Rabbit mAb [PSH13-07] - Detector	HA725087B	Human	ELISA(Det),ELISA
	Human CD31 Recombinant Rabbit mAb [PSH07-82] - BSA and Azide free (Capture)	HA722906	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human CD31 Recombinant Rabbit mAb [PSH07-83] - BSA and Azide free (Detector)	HA722907	Human	ELISA(Det)
	Biotin Conjugated Human CD31 Recombinant Rabbit mAb [PSH07-83] - Detector	HA722908B	Human	ELISA(Det)
	Biotin Conjugated Human CD47 Recombinant Rabbit mAb [PSH11-06] - Detector	HA723274B	Human	ELISA(Det),ELISA
	Human CD47 Recombinant Rabbit mAb [PSH11-05] - BSA and Azide free (Capture)	HA723272	Human	ELISA(Cap)
	Human CD47 Recombinant Rabbit mAb [PSH11-06] - BSA and Azide free (Detector)	HA723273	Human	ELISA(Det)
	Human CX3CL1 Recombinant Rabbit mAb [PSH16-09] - BSA and Azide free (Detector)	HA723814	Human	ELISA(Det)
	Human CX3CL1 Recombinant Rabbit mAb [PSH16-08] - BSA and Azide free (Capture)	HA723813	Human	ELISA(Cap)
	Biotin Conjugated Human CX3CL1 Recombinant Rabbit mAb [PSH16-09] - Detector	HA723814B	Human	ELISA(Det),ELISA
	Human IL-10 Recombinant Rabbit mAb [PS00-98] - BSA and Azide free (Capture)	HA721236	Human	ELISA(Cap)
	Human IL-10 Recombinant Rabbit mAb [PS00-81] - BSA and Azide free (Detector)	HA721237	Human	ELISA(Det)
	Biotin Conjugated Human IL-10 Recombinant Rabbit mAb [PS00-81]	HA722652B	Human	ELISA(Det)
	Biotin Conjugated Human KIM-1/HAVCR1 Recombinant Rabbit mAb [PSH18-81]	HA725250B	Human	ELISA(Det),ELISA
	Human KIM-1/HAVCR1 Recombinant Rabbit mAb [PSH18-80] - BSA and Azide free (Capture)	HA725249	Human	ELISA(Cap)
	Human KIM-1/HAVCR1 Recombinant Rabbit mAb [PSH18-81] - BSA and Azide free (Detector)	HA725250	Human	ELISA(Det)
	Human MERTK Recombinant Rabbit mAb [PSH12-70] - BSA and Azide free (Detector)	HA723480	Human	ELISA(Det)
	Biotin Conjugated Human MERTK Recombinant Rabbit mAb [PSH12-70] - Detector	HA723481B	Human	ELISA(Det),ELISA
	HRP Conjugated Human MERTK Recombinant Rabbit mAb [PSH12-70] - Detector	HA723481H	Human	ELISA(Det),ELISA
	Human MERTK Recombinant Rabbit mAb [PSH12-69] - BSA and Azide free (Capture)	HA723479	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human SIRP alpha Recombinant Rabbit mAb [PSH17-64] - BSA and Azide free (Detector)	HA725348	Human	ELISA(Det)
	Biotin Conjugated Human SIRP alpha Recombinant Rabbit mAb [PSH17-64] - Detector	HA725348B	Human	ELISA(Det),ELISA
	Human SIRP alpha Recombinant Rabbit mAb [PSH17-63] - BSA and Azide free (Capture)	HA725347	Human	ELISA(Cap)
	Human TIM 4 Recombinant Rabbit mAb [PSH14-01] - BSA and Azide free (Detector)	HA723590	Human	ELISA(Det)
	Biotin Conjugated Human TIM 4 Recombinant Rabbit mAb [PSH14-01] - Detector	HA723591B	Human	ELISA(Det),ELISA
	HRP Conjugated Human TIM 4 Recombinant Rabbit mAb [PSH14-01] - Detector	HA723591H	Human	ELISA(Det),ELISA
	Human TIM 4 Recombinant Rabbit mAb [PSH14-00] - BSA and Azide free (Capture)	HA723589	Human	ELISA(Cap)
	IL-10 Rabbit pAb	ER1911-19	Human, Mouse	WB,IHC-P,IF-Tissue
	IL-10 Rabbit pAb	ER65464	Human, Mouse, Rat	WB,IHC-P
Cell treatment	IL-10 Recombinant Rabbit mAb [JE59-71]	HA722032	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	Integrin alpha V Recombinant Rabbit mAb [SC56-07] - BSA and Azide free	HA750201	Human, Mouse, Rat	WB,IF-Cell,FC,IHC-P
KD	Integrin alpha V Recombinant Rabbit mAb [SC56-07]	ET1610-15	Human, Mouse, Rat	WB,IF-Cell,FC,IHC-P
	Integrin Alpha V Rabbit pAb	ER1911-52	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,IP
	Integrin beta 3 Rabbit pAb	ER1601-4	Human, Mouse, Rat	WB,IHC-P
	Integrin beta 3 Recombinant Rabbit mAb [JE22-64]	ET7110-61	Human, Mouse	WB,IHC-P
KD	Integrin beta 3 Recombinant Rabbit mAb [SJ19-09]	ET1606-49	Human, Mouse	WB,IHC-P,IF-Tissue
Relative expression	Integrin beta 3 Recombinant Rabbit mAb [PSH09-27]	HA723085	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Integrin beta 5 Recombinant Rabbit mAb [PSH11-53]	HA723345	Human	IF-Cell,FC
	Integrin beta 5 Recombinant Rabbit mAb [PSH10-13] - BSA and Azide free	HA751342	Human, Rat, Monkey	WB,IHC-P
	Integrin beta 5 Rabbit pAb	ER63768	Human, Mouse, Rat, Monkey	WB
	Integrin beta 5 Recombinant Rabbit mAb [PSH10-13]	HA723182	Human, Rat, Monkey	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Integrin beta 5 Recombinant Rabbit mAb [PSH11-53] - BSA and Azide free	HA751415	Human	IF-Cell,FC
	JAK2 Rabbit pAb	RT1343	Human, Mouse, Rat	WB,IP,IF,IHC-P
	JAK2 Rabbit pAb	ER1706-58	Human, Mouse, Rat	IF-Cell,IHC-P
KD	JAK2 Mouse mAb [6-D3]	M1501-8	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	JAK2 Recombinant Rabbit mAb [SY0245] - BSA and Azide free	HA750118	Human, Mouse, Rat	WB,IF-Cell
KD	JAK2 Recombinant Rabbit mAb [SY0245]	ET1607-35	Human, Mouse, Rat	WB,IF-Cell
Relative expression	Liver Arginase Recombinant Rabbit mAb [PSH09-78]	HA723143	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Liver Arginase Recombinant Rabbit mAb [PD00-05]	HA721147	Human	WB,IHC-P,IF-Tissue
Relative expression	Liver Arginase Recombinant Rabbit mAb [JB21-49]	ET7107-90	Human	WB,IHC-P,IP
Cell treatment	Liver Arginase Recombinant Rabbit mAb [SY09-06]	ET1605-8	Human, Mouse, Rat	WB
	LOX15 Rabbit pAb	ER65960	Human, Mouse, Rat	WB,IHC
	LRP-1 Mouse mAb [A4-F5]	M1211-4	Mouse	WB,IHC-P
KD	LRP1 Recombinant Rabbit mAb [SA0290]	ET1601-1	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,IP,IHC-Fr
	LXR alpha Recombinant Rabbit mAb [JA20-38]	ET1704-51	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	LXR alpha Recombinant Rabbit mAb [JA20-38] - BSA and Azide free	HA750414	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	LXR alpha + LXR beta Rabbit pAb	ER1912-55	Human, Mouse	WB,FC
KD	LXR beta Recombinant Rabbit mAb [PSH0-56]	HA721403	Human, Mouse, Rat	WB,IF-Cell
	LXR beta Rabbit pAb	ER1912-56	Mouse	WB
KD	LXR beta Recombinant Rabbit mAb [PSH0-56] - BSA and Azide free	HA750633	Human, Mouse, Rat	WB,IF-Cell
KD	MAPKAP Kinase 2 Recombinant Rabbit mAb [JU30-31]	ET7106-83	Human, Mouse	WB,IP
	MCT1 Recombinant Rabbit mAb [JE41-72]	HA722556	Human	WB
	MCT1 Rabbit pAb	ER1912-87	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
	MCT1 Rabbit pAb	HA500190	Human	WB,IF-Cell,IHC-P,FC
	MEGF11 Rabbit pAb	ER1912-95	Mouse	WB
KD	MEK1 Recombinant Rabbit mAb [SZ22-01]	ET1603-20	Human, Mouse, Rat, Cow, Dog	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	MEK1/2 Recombinant Rabbit mAb [SR13-07] - BSA and Azide free	HA750042	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	MEK1/2 Recombinant Rabbit mAb [SR13-07]	ET1602-3	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	MEK2 Recombinant Rabbit mAb [SD207-08]	ET1612-6	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	MERTK Recombinant Rabbit mAb [SR29-07] - BSA and Azide free	HA750036	Human	WB,IHC-P,IP,FC,IF-Tissue
KD	MERTK Recombinant Rabbit mAb [SR29-07]	ET1602-21	Human	WB,IHC-P,FC,IF-Tissue
	MFGE8 Recombinant Rabbit mAb [PSH05-75] - BSA and Azide free	HA751025	Human	WB
	MFGE8 Recombinant Rabbit mAb [PSH05-75]	HA722442	Human	WB
	MK13 Rabbit pAb	ER60146	Human, Mouse, Rat	WB
	MKP-1/2 Rabbit pAb	ER63964	Human, Mouse, Rat	WB,IHC-P
	MKP-7 Rabbit pAb	ER63967	Human, Mouse	WB,IHC-P,IF-Cell
	Biotin Conjugated Mouse IL-10 Recombinant Rabbit mAb [PS01-93] - Detector	HA721519	Mouse	ELISA(Det)
	Mouse IL-10 Recombinant Rabbit mAb [PS01-93] - BSA and Azide free (Detector)	HA721518	Mouse	ELISA(Det)
	Mouse IL-10 Recombinant Rabbit mAb [PS01-85] - BSA and Azide free (Capture)	HA721517	Mouse	ELISA(Cap)
Relative expression	NFAT1 Recombinant Rabbit mAb [JA11-08]	ET1704-14	Human	WB,IF-Tissue,IHC-P
Relative expression	NFAT1 Recombinant Rabbit mAb [JA11-08] - BSA and Azide free	HA750405	Human	WB,IF-Tissue,IHC-P
	NFAT2 Recombinant Rabbit mAb [JA81-03]	ET1704-45	Human	WB,IHC-P,FC
	NFATc3 Rabbit pAb	ER64179	Human, Mouse	WB,IHC-P
	NFATc4 Rabbit pAb	ER64180	Human, Mouse	WB,IHC-P,IF-Cell
	NF-kB p105 / p50 Rabbit pAb	R1309-9	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	NF-kB p105 / p50 Mouse mAb [A2-7]	M1505-8	Human, Mouse	WB,IF-Cell,IHC-P
KD	NF-kB p105 / p50 Recombinant Rabbit mAb [SZ20-01] - BSA and Azide free	HA750060	Human, Mouse, Rat	WB,IHC-P,ICC,FC,IF-Tissue
KD	NF-kB p105 / p50 Recombinant Rabbit mAb [SZ20-01]	ET1603-18	Human, Mouse, Rat	WB,IHC-P,ICC,FC,IF-Tissue
KO, KD, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [SZ10-04]	ET1603-12	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KO, KD, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [SZ10-04] - BSA and Azide free	HA750058	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	NF-kB p65 Rabbit pAb	ER0815	Human, Mouse, Rat, Zebrafish	WB,IHC-P,FC
KO, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [PSH0-27] - BSA and Azide free	HA750601	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
	NF-kB p65 Rabbit pAb	HA500402	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
KO, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [PSH0-27]	HA721307	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
	NGEP Rabbit pAb	ER64185	Human, Mouse	WB,IHC-P,IF-Cell
	p130 Cas Rabbit pAb	ER64487	Human, Mouse, Rat	WB,IHC-P
	P2RY2 Rabbit pAb	ER62266	Human, Rat, Mouse	WB
	P2Y12 Recombinant Rabbit mAb [PSH04-48]	HA722133	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-P,IHC-Fr,IF-Tissue
	P2Y12 Recombinant Antibody [PSH04-48] - Rat IgG1 (Chimeric)	HA601388	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P
	P2Y12 Recombinant Antibody [PSH04-48] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610231	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P
	P2Y12 Recombinant Rabbit mAb [PSH04-48] - BSA and Azide free	HA750957	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-P,IHC-Fr,IF-Tissue
	P2Y12 Recombinant Rabbit mAb [JE63-23]	HA721011	Human, Mouse	WB
	P2Y6 Recombinant Rabbit mAb [JG37-77]	ET7108-52	Human	WB,IP,IHC-P,IF-Cell
KO	p38 alpha / MAPK14 Recombinant Rabbit mAb [JF55-07]	ET1702-65	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	p38 alpha / MAPK14 Mouse mAb [9-B9]	EM40307	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Cell treatment	p38 alpha / MAPK14 Rabbit pAb	R1308-3	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD, Cell treatment	p38 alpha / MAPK14 Recombinant Rabbit mAb [SR43-04]	ET1602-26	Human, Mouse, Rat	WB,IF-Cell,FC
KD, Cell treatment	p38 alpha / MAPK14 Recombinant Rabbit mAb [SR43-04] - BSA and Azide free	HA750040	Human, Mouse, Rat	WB,IF-Cell,FC
	p38 alpha/MAPK14 + p38 beta/MAPK11 + p38 gamma/MAPK12 Recombinant Rabbit mAb [PSH03-60] - BSA and Azide free	HA750894	Human, Mouse, Rat, Monkey	WB
	p38 alpha/MAPK14 + p38 beta/MAPK11 + p38 gamma/MAPK12 Recombinant Rabbit mAb [PSH03-60]	HA722019	Human, Mouse, Rat, Monkey	WB
	p38 beta/MAPK11 + p38 gamma/MAPK12 + p38 delta/MAPK13 Recombinant Rabbit mAb [PSH03-59] - BSA and Azide free	HA750893	Human, Mouse, Rat, Monkey	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	p38 beta/MAPK11 + p38 gamma/MAPK12 + p38 delta/MAPK13 Recombinant Rabbit mAb [PSH03-59]	HA722018	Human, Mouse, Rat, Monkey	WB
	p38 delta / MAPK13 Rabbit pAb	HA500582	Human	WB
	p38 gamma/MAPK12 Recombinant Rabbit mAb [PSH03-58]	HA722017	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	p38 gamma/MAPK12 Recombinant Rabbit mAb [PSH03-58] - BSA and Azide free	HA750892	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	p38 MAPK Recombinant Rabbit Multiclonal Antibody [PSH03-72]	HA722036	Human, Mouse, Rat, Monkey	WB
	p38 MAPK Recombinant Rabbit Multiclonal Antibody [PSH03-72] - BSA and Azide free	HA750906	Human, Mouse, Rat, Monkey	WB
	p53 Recombinant Rabbit mAb [PD01-32] - BSA and Azide free	HA750545	Human	IHC-P
Cell treatment, Relative expression, KD	p53 Recombinant Mouse mAb [6B5-R] - BSA and Azide free	HA610187	Human, Mouse, Monkey	WB,IF-Cell
Relative expression, KD, Cell treatment	p53 Recombinant Rabbit mAb [SA39-07]	ET1601-13	Human	WB,IF-Cell,ChIP
	p53 Mouse mAb [C8-A11]	M1312-2	Human	WB,IF-Cell,IHC-P
Cell treatment	p53 Recombinant Rabbit mAb [PSH03-95] - BSA and Azide free	HA750920	Mouse	WB,IF-Cell,IP
Cell treatment	p53 Recombinant Rabbit mAb [PSH03-95]	HA722074	Mouse	WB,IF-Cell,IP
KD	p53 Rabbit pAb	R1308-7	Human	WB,IF-Cell
KO	p53 Mouse mAb [7-5]	EM20603	Human	WB,IF-Cell,IHC-P,FC
	p53 Recombinant Rabbit mAb [PD01-32]	HA721213	Human	IHC-P
Relative expression, KD, Cell treatment	p53 Recombinant Rabbit mAb [SA39-07] - BSA and Azide free	HA750005	Human	WB,IF-Cell,ChIP
Cell treatment, KD, Relative expression	p53 Recombinant Mouse mAb [6B5-R]	HA601315	Human, Mouse, Monkey	WB,IF-Cell
KD	p53 Recombinant Mouse mAb [6B2-R] - BSA and Azide free	HA610194	Human, Monkey	WB,IF-Cell
	p53 Mouse mAb [6B5]	EM1701-92	Human	WB
KD	p53 Recombinant Mouse mAb [6B2-R]	HA601322	Human, Monkey	WB,IF-Cell
KD	p53 Recombinant Rabbit mAb [JB42-26]	ET7107-33	Human	WB,IHC-P,FC,IP,IF-Cell,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression, KD	p53 Recombinant Rabbit mAb [SY010-6]	ET1605-16	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Relative expression	p53 Mouse mAb [6B2]	EM1701-91	Human, Mouse	WB,IF-Cell
Cell treatment	p53 (acetyl K370) Recombinant Rabbit mAb [SR40-09] - BSA and Azide free	HA750047	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,FC
Cell treatment	p53 (acetyl K370) Recombinant Rabbit mAb [SR40-09]	ET1602-38	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,FC
Cell treatment	p53 (acetyl K382) Recombinant Rabbit mAb [PSH06-88] - BSA and Azide free	HA751112	Human	WB,IF-Cell
Cell treatment	p53 (acetyl K382) Recombinant Rabbit mAb [SD0801]	ET1612-71	Human, Mouse	WB,IF-Cell,IHC-P
Cell treatment	p53 (acetyl K382) Recombinant Rabbit mAb [PSH06-87] - BSA and Azide free	HA751111	Human, Mouse	WB,IF-Cell
Cell treatment	p53 (acetyl K382) Recombinant Rabbit mAb [PSH06-88]	HA722763	Human	WB,IF-Cell
Cell treatment	p53 (acetyl K382) Recombinant Rabbit mAb [PSH06-87]	HA722762	Human, Mouse	WB,IF-Cell
	Pannexin 1 Recombinant Rabbit mAb [JE38-97]	HA721629	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	Paxillin Recombinant Rabbit mAb [SY23-02] - BSA and Azide free	HA750114	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Paxillin Rabbit pAb	RT1466	Human, Mouse, Rat	WB,IP,IF
	Paxillin Recombinant Rabbit mAb [SY23-02]	ET1607-22	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	PBX1 Recombinant Rabbit mAb [PSH0-24] - BSA and Azide free	HA750598	Human, Mouse, Rat	WB,IHC
	PBX1 Recombinant Rabbit mAb [PSH0-24]	HA721303	Human, Mouse, Rat	WB,IHC
	PBX1 Rabbit pAb	ER1914-79	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,FC
	PGC1 alpha+beta Recombinant Rabbit mAb [JF09-71]	ET1702-96	Human, Mouse, Rat	WB
	PGC1 beta Recombinant Rabbit mAb [JM33-24]	ET1705-65	Human, Rat	WB,IP,FC
Cell treatment	Phospho-Axl (Y702) Recombinant Rabbit mAb [PSH17-98]	HA723963	Human	WB
Cell treatment	Phospho-BCAR1 (Y410) Recombinant Rabbit mAb [JE47-97]	HA722124	Human	WB
Cell treatment	Phospho-CaMKII (T286) Recombinant Rabbit mAb [JE51-81]	HA721794	Human, Mouse, Rat	WB
	Phospho-CaMKII (T286) Mouse mAb [2G1]	RT1092	Human, Mouse, Rat	WB,IP,IF
Cell treatment	Phospho-Creb (S133) Recombinant Rabbit mAb [JB25-40] - BSA and Azide free	HA750465	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-Creb (S133) Recombinant Rabbit mAb [JB25-40]	ET7107-93	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC,IF-Tissue
Cell treatment	Phospho-CrkL (Y207) Recombinant Rabbit mAb [JE55-58]	HA721675	Human, Mouse, Rat	WB
Cell treatment	Phospho-Erk1 (T202)+Erk2 (T185) Recombinant Rabbit mAb [SZ2-4] - BSA and Azide free	HA750061	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
Cell treatment	Phospho-Erk1 (T202)+Erk2 (T185) Recombinant Rabbit mAb [SZ2-4]	ET1603-22	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
Cell treatment	Phospho-Erk1 (T202 + Y204) + Erk2 (T185 + Y187) Recombinant Rabbit mAb [SC58-01] - BSA and Azide free	HA750199	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP
Cell treatment	Phospho-Erk1 (T202 + Y204) + Erk2 (T185 + Y187) Recombinant Rabbit mAb [SC58-01]	ET1610-13	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP
	Phospho-FAK (S722) Mouse mAb [4G1]	RT1216	Human, Mouse, Rat	WB,IP,IF,IHC-P
Cell treatment	Phospho-FAK (Y397) Recombinant Rabbit mAb [SC54-07] - BSA and Azide free	HA750212	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	Phospho-FAK (Y397) Recombinant Rabbit mAb [SC54-07]	ET1610-34	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	Phospho-FOXO1A (T24)+ FOXO3A (T32)+ FOXO4 (T28)+ FOXO6 (T26) Recombinant Rabbit mAb [PSH06-00]	HA722499	Human, Mouse	WB,IF-Cell,FC
Cell treatment	Phospho-FOXO1A (T24)+ FOXO3A (T32)+ FOXO4 (T28)+ FOXO6 (T26) Recombinant Rabbit mAb [PSH06-00] - BSA and Azide free	HA751040	Human, Mouse	WB,IF-Cell,FC
Cell treatment	Phospho-FOXO3a (S253) Recombinant Rabbit mAb [ST49-01]	ET1609-49	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-FOXO3a (S253) Recombinant Rabbit mAb [ST49-01] - BSA and Azide free	HA750178	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-JAK2 (Y1007 + Y1008) Recombinant Rabbit mAb [SY24-03] - BSA and Azide free	HA750117	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
Cell treatment	Phospho-JAK2 (Y1007 + Y1008) Recombinant Rabbit mAb [SY24-03]	ET1607-34	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
Cell treatment	Phospho-MEK1/2 (S218 + S222) Recombinant Rabbit mAb [ST0490]	ET1609-50	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Phospho-MEK1 (S298) Recombinant Rabbit mAb [SD206-7]	ET1612-40	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	Phospho-MEK1 (T292) Recombinant Rabbit mAb [SD2088]	ET1612-42	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-NF-kB p65 (S529) Recombinant Rabbit mAb [SP07-00]	ET1604-27	Human	WB,IHC-P,FC,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-NF-kB p65 (S536) Recombinant Rabbit mAb [PSH10-58] - BSA and Azide free	HA751361	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-NF-kB p65 (S536) Recombinant Rabbit mAb [PSH10-58]	HA723223	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-p38 alpha (T180) Recombinant Rabbit mAb [PSH06-25]	HA722657	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-p38 alpha (T180) Recombinant Rabbit mAb [PSH06-25] - BSA and Azide free	HA751062	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-p38 alpha (T180 + Y182) Recombinant Rabbit mAb [PSH06-35]	HA722669	Human, Mouse, Rat	WB,IF-Cell,FC
Cell treatment	Phospho-p38 alpha (T180 + Y182) Recombinant Rabbit mAb [PSH06-35] - BSA and Azide free	HA751074	Human, Mouse, Rat	WB,IF-Cell,FC
Cell treatment	Phospho-p38 (T180) Rabbit pAb	ER2001-51	Human, Mouse, Rat	WB
Cell treatment	Phospho-p38 (T180) Recombinant Rabbit mAb [JE77-37]	HA722150	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-p38 (T180 + Y182) Rabbit pAb	ER2001-52	Human, Mouse, Rat	WB
	Phospho-p38 (T180 + Y182) Rabbit pAb	ER1903-01	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
Cell treatment	Phospho-p53 (S15) Recombinant Rabbit mAb [PSH01-98] - BSA and Azide free	HA750772	Human	WB,IF-Cell
Cell treatment	Phospho-p53 (S15) Recombinant Rabbit mAb [PSH01-98]	HA721756	Human	WB,IF-Cell
Cell treatment	Phospho-p53 (S20) Recombinant Rabbit mAb [PSH06-86]	HA722761	Human	WB
Cell treatment	Phospho-p53 (S20) Recombinant Rabbit mAb [PSH06-86] - BSA and Azide free	HA751110	Human	WB
Cell treatment	Phospho-p53 (S33) Recombinant Rabbit mAb [JE41-97]	HA722703	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-p53 (S376) Recombinant Rabbit mAb [ST0440]	ET1609-14	Human, Mouse	WB
Cell treatment	Phospho-p53 (S392) Recombinant Rabbit mAb [SI17-04]	ET1606-24	Human, Mouse, Rat	WB,IHC-P,IP
Cell treatment	Phospho-p53 (S392) Recombinant Rabbit mAb [SI17-04] - BSA and Azide free	HA750097	Human, Mouse, Rat	WB,IHC-P,IP
Cell treatment	Phospho-p53 (T55) Recombinant Rabbit mAb [ST0439]	ET1609-13	Human	WB,IHC-P
Cell treatment	Phospho-Paxillin (Y118) Recombinant Rabbit mAb [JE09-10]	HA722526	Human	WB
Cell treatment	Phospho-Rac1+Cdc42 (S71) Recombinant Rabbit mAb [JJ086-9]	ET1701-40	Human	WB,IF-Cell,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-SHP2 (Y542) Recombinant Rabbit mAb [SN61-01]	ET1611-22	Mouse, Human	WB,IP
Cell treatment	Phospho-SIRT1 (S47) Recombinant Rabbit mAb [SN06-40]	ET1611-31	Human	WB,IF-Cell,IHC-P
	Phospho-SIRT1 (T530) Recombinant Rabbit mAb [JJ206-6]	ET1701-27	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
Cell treatment, KD	Phospho-STAT3 (S727) Recombinant Rabbit mAb [SY24-09] - BSA and Azide free	HA750121	Human, Rat, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP
Cell treatment, KD	Phospho-STAT3 (S727) Recombinant Rabbit mAb [SY24-09]	ET1607-39	Human, Rat, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP
Cell treatment	Phospho-STAT3 (Y705) Recombinant Rabbit mAb [SZ43-01]	ET1603-40	Human, Mouse, Rat	WB,IHC-P,IP
Cell treatment	Phospho-STAT3 (Y705) Recombinant Rabbit mAb [SZ43-01] - BSA and Azide free	HA750067	Human, Mouse, Rat	WB,IHC-P,IP
	PPAR Rabbit pAb	HA500569	Mouse, Rat	WB
	PPAR delta Rabbit pAb	ER1902-24	Human, Mouse	WB,FC
	PPAR Delta Mouse mAb [2F9]	EM60127	Human, Rat, Mouse	IF-Cell,IHC-P
KD	Pro Caspase-3 Recombinant Rabbit mAb [SZ02-08]	ET1603-26	Human, Monkey, Mouse	WB,IF-Tissue,IHC-P
	Pro Caspase-3 Rabbit pAb	ER30804	Human, Mouse	WB,IF-Cell,IHC-P,FC,IHC-Fr
Cell treatment, Relative expression, KO	Pro Caspase-3 Recombinant Rabbit mAb [SR03-01] - BSA and Azide free	HA750048	Human, Monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Cell treatment, Relative expression, KO	Pro Caspase-3 Recombinant Rabbit mAb [SR03-01]	ET1602-39	Human, Monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	pro Caspase-7 Recombinant Rabbit mAb [SD20-57]	ET1612-28	Human	WB,IF-Cell,IF-Tissue,IHC-P
	Protein S Rabbit pAb	0812-4	Human, Mouse, Rat	WB,IHC-P
	PTGER2 Rabbit pAb	ER1915-56	Human, Mouse, Rat	WB
	PTGER2 Recombinant Rabbit mAb [JE30-72]	HA721380	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell
	RAB14 Rabbit pAb	ER62485	Human, Mouse, Rat	WB
	RAB17 Rabbit pAb	ER61096	Human, Mouse	WB
KD	Rab5 Recombinant Rabbit mAb [ST47-04] - BSA and Azide free	HA750167	Human, Mouse, Rat	WB,IHC-P,IF-Cell
KD	Rab5 Recombinant Rabbit mAb [ST47-04]	ET1609-27	Human, Mouse, Rat	WB,IHC-P,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Rab5a Mouse mAb [C11-H5]	EM1708-08	Human	WB,IF-Cell,IHC-P,FC
	RAB5B Rabbit pAb	ER61113	Human, Mouse	WB
KD	RAB5C Recombinant Rabbit mAb [JE63-14]	HA721012	Human, Mouse, Rat	WB,IHC-P,IF-Cell
KD	RAB7 Recombinant Rabbit mAb [SN202-03]	ET1611-96	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	RAB7 Rabbit pAb	ER1802-53	Human, Mouse, Rat	WB,IHC-P
KD	RAB7 Recombinant Rabbit mAb [SN202-03] - BSA and Azide free	HA750279	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	RAC1 Recombinant Rabbit mAb [PSH05-09] - BSA and Azide free	HA750982	Human, Rat	WB,IF-Cell
	RAC1 Recombinant Rabbit mAb [PSH05-09]	HA722225	Human, Rat	WB,IF-Cell
	Rac1-2-3 Recombinant Rabbit mAb [JM11-29]	ET1703-80	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	RAC1/2/3 Rabbit pAb	HA500424	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	RAC2 Rabbit pAb	ER61092	Human, Mouse	WB
	RAC3 Rabbit pAb	ER61093	Human, Mouse	WB
KD	RAC3 Recombinant Rabbit mAb [JE56-48]	HA721890	Human, Mouse, Rat	WB
Relative expression	RAGE Recombinant Rabbit mAb [JF0975]	ET1702-27	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,mIHC
Relative expression	RAGE Recombinant Rabbit mAb [JF0975] - BSA and Azide free	HA750341	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	RAGE Recombinant Rabbit mAb	IRS195RB	Human	mIHC
	RAGE Recombinant Rabbit mAb	IRS073	Mouse	mIHC
	Rat IL-10 Recombinant Rabbit mAb [PSH02-21] - BSA and Azide free (Capture)	HA721786	Rat	ELISA(Cap)
	Rat IL-10 Recombinant Rabbit mAb [PSH02-22] - BSA and Azide free (Detector)	HA721787	Rat	ELISA(Det)
	Biotin Conjugated Rat IL-10 Recombinant Rabbit mAb [PSH02-22]	HA722637B	Rat	ELISA(Det)
	Rat KIM-1 / HAVCR1 Recombinant Rabbit mAb [PSH20-49] - BSA and Azide free (Capture)	HA725381	Rat	ELISA(Cap)
	Rat KIM-1 / HAVCR1 Recombinant Rabbit mAb [PSH20-50] - BSA and Azide free (Detector)	HA725382	Rat	ELISA(Det)
	Biotin Conjugated Rat KIM-1 / HAVCR1 Recombinant Rabbit mAb [PSH20-50] - Detector	HA725382B	Rat	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Retinoid X Receptor alpha Recombinant Rabbit mAb [JG99-38]	ET7108-99	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
	Ribosomal Protein L3 Rabbit pAb	ER64833	Human, Mouse, Rat	WB,IHC-P
	Ribosomal Protein S19BP Rabbit pAb	ER64854	Human	IHC-P
	SERCA2 Recombinant Rabbit mAb [JM10-20]	ET1703-01	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	SERCA3 Rabbit pAb	ER1916-41	Human	WB
	SGK1 Recombinant Rabbit mAb [SC05-71]	ET1610-19	Human, Mouse, Rat	WB,IHC-P,IP,FC,IF-Tissue
	SHP1 Recombinant Rabbit mAb [SR41-02]	ET1602-19	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
Relative expression	SHP1 Rabbit pAb	ER1802-74	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	SHP1 Rabbit pAb	R1511-24	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	SHP1 Recombinant Rabbit mAb [SR41-02] - BSA and Azide free	HA750033	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	SHP1 Mouse mAb [B5-A4]	M1505-4	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	SHP2 Rabbit pAb	ER1706-93	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	SHP2 Recombinant Rabbit mAb [SN07-32]	ET1611-80	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP
KD	SHP2 Recombinant Rabbit mAb [SN72-02]	ET1611-35	Human, Mouse, Rat	WB,IP,IHC-P
	SIGLEC10 Rabbit pAb	ER1916-59	Human, Mouse, Rat	WB
KD	SIRP alpha Recombinant Rabbit mAb [JE36-20]	HA721316	Human, Mouse, Rat	WB,IHC-P,FC
	SIRP alpha Mouse mAb [A3E3]	EM1902-37	Human	WB,IF-Cell,IHC-P,FC
	SIRP Alpha Rabbit pAb	ER1916-61	Human, Mouse, Rat	WB,FC
KD	SIRT1 Recombinant Rabbit mAb [SZ04-01] - BSA and Azide free	HA750064	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,ChIP
KD	SIRT1 Mouse mAb [7-C5-B2]	M1506-3	Human, Mouse	WB,IF-Cell,IHC-P
KD	SIRT1 Rabbit pAb	ER130811	Human, Mouse	WB,IHC-P,FC,IF-Cell
KD	SIRT1 Recombinant Rabbit mAb [JE11-04]	ET7110-49	Human	WB,IF-Cell
KD	SIRT1 Recombinant Rabbit mAb [SZ04-01]	ET1603-3	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,ChIP
KD	SIRT2 Recombinant Rabbit mAb [JJ08-83]	ET1701-43	Human	WB
	SIRT2 Recombinant Antibody [SN70-04] - Mouse IgG1 (Chimeric)	HA601505	Mouse, Rat	IHC-Fr,IHC-P,WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	SIRT2 Recombinant Antibody [SN70-04] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610335	Mouse, Rat	IHC-Fr,IHC-P,WB
	SIRT2 Recombinant Antibody [SN70-04] - Rat IgG1 (Chimeric)	HA601506	Mouse, Rat	IHC-Fr,IHC-P,WB
	SIRT2 Recombinant Antibody [SN70-04] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610336	Mouse, Rat	IHC-Fr,IHC-P,WB
	SIRT2 Recombinant Rabbit mAb [SN70-04]	ET1611-72	Human, Cynomolgus monkey	WB,IP,IHC-P
	SIRT2 Recombinant Rabbit mAb [PSH17-48]	HA723923	Human, Mouse, Rat, Cynomolgus monkey	IHC-Fr,IHC-P,WB
KD	SIRT3 Recombinant Rabbit mAb [PSH05-22]	HA722251	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IF-Tissue,IP,ChIP
KD	SIRT3 Recombinant Rabbit mAb [PSH05-22] - BSA and Azide free	HA750993	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IF-Tissue,IP,ChIP
	SIRT3 Rabbit pAb	R1511-3	Human, Mouse, Rat	WB
	SIRT5 Recombinant Rabbit mAb [PSH15-71]	HA723774	Human, Mouse, Rat	WB,IHC-P
	SIRT5 Recombinant Rabbit mAb [JJ084-01]	ET1701-6	Human, Mouse	WB
	SIRT5 Recombinant Rabbit mAb [PSH15-71] - BSA and Azide free	HA751578	Human, Mouse, Rat	WB,IHC-P
KD	SIRT6 Recombinant Rabbit mAb [JJ20-69]	ET1701-29	Human	WB
KO	SIRT6 Recombinant Rabbit mAb [PSH10-99] - BSA and Azide free	HA751389	Human, Mouse, Rat, Monkey	WB,IHC-P
	SIRT6 Recombinant Rabbit mAb [PSH10-98]	HA723265	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P
KO	SIRT6 Recombinant Rabbit mAb [PSH10-99]	HA723266	Human, Mouse, Rat, Monkey	WB,IHC-P
	SIRT6 Rabbit pAb	R1511-1	Human, Rat	WB,IF-Cell,FC
	SIRT6 Recombinant Rabbit mAb [PSH10-98] - BSA and Azide free	HA751388	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P
	SIRT7 Recombinant Rabbit mAb [PSH09-17] - BSA and Azide free	HA751271	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,ChIP
	SIRT7 Mouse mAb [G10-E2]	EM1707-09	Human	WB,IHC-P
	SIRT7 Recombinant Rabbit mAb [PSH09-17]	HA723075	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,ChIP
	SK2 Mouse mAb [E1-C5]	EM1707-15	Human	IF-Cell,IHC-P,FC
	SOD2 Recombinant Rabbit mAb [JJ089-02]	ET1701-54	Human, Mouse, Rat	WB,IHC-P
	SOD2 Recombinant Rabbit mAb [JJ089-02] - BSA and Azide free	HA750314	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	SOD2 (acetyl K68) Recombinant Rabbit mAb [JE62-28]	HA722463	Human, Mouse, Rat	WB,IHC-P
Relative expression	SPHK1 Recombinant Rabbit mAb [JA31-14]	ET1704-76	Human, Mouse	WB
	STAB1 Rabbit pAb	ER61883	Human, Mouse	IHC-P
	STAB2 Rabbit pAb	ER61939	Human, Rat, Mouse	IHC-P
KD	STAT3 Recombinant Rabbit mAb [SY34-01] - BSA and Azide free	HA750087	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	STAT3 Recombinant Rabbit mAb [SY24-08]	ET1607-38	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC
KD	STAT3 Recombinant Rabbit mAb [SY34-01]	ET1605-45	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	STAT3 Rabbit pAb	RT1591	Human, Mouse, Rat	WB,IP,IF,IHC-P
KD	STAT3 Recombinant Rabbit mAb [SY24-08] - BSA and Azide free	HA750120	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC
	TGF beta 1 Mouse mAb [B2-F6]	EM1709-93	Human	WB,IHC-P,FC
	TGF beta 1 Rabbit pAb	HA500496	Human	WB
Relative expression	TGF beta 1 Recombinant Rabbit mAb [PD00-17] - BSA and Azide free	HA750485	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Tissue
	TGF beta 1 Recombinant Rabbit mAb [PSH15-21] - BSA and Azide free (Capture)	HA723725	Human, Mouse, Rat	ELISA(Cap)
	Biotin Conjugated TGF beta 1 Recombinant Rabbit mAb [PSH15-22] - Detector	HA723727B	Human, Mouse, Rat	ELISA(Det),ELISA
	TGF beta 1 Recombinant Rabbit mAb [PSH15-22] - BSA and Azide free (Detector)	HA723726	Human, Mouse, Rat	ELISA(Det)
	TGF beta 1 Recombinant Rabbit mAb [PSH18-59]	HA724012	Human, Mouse, Monkey	WB,IF-Cell,IHC-P,FC
	TGF beta 1 Recombinant Rabbit mAb	IRS205RB	Human	mIHC
Relative expression	TGF beta 1 Recombinant Rabbit mAb [PD00-17]	HA721143	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Tissue
	Thrombospondin 1 Recombinant Rabbit mAb [PSH02-97] - BSA and Azide free	HA750852	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC
	Thrombospondin 1 Recombinant Rabbit mAb [PSH02-97]	HA721916	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC
	TIM 1 Recombinant Rabbit mAb [JE39-41]	HA721535	Human, Mouse, Rat	WB
	TIM 3 Mouse mAb [1E6]	EM1701-18	Human	WB,IF-Cell,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	TIM 3 Recombinant Antibody [PSH13-95] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610277	Human	WB,IF-Cell,FC
Relative expression	TIM 3 Recombinant Rabbit mAb [PSH16-50] - BSA and Azide free	HA751599	Human, Mouse, Rat	WB,IHC-P,IP
	TIM 3 Mouse mAb [1E3]	EM1701-19	Human	WB,IHC-P,IF-Cell,FC
Relative expression	TIM 3 Recombinant Rabbit mAb [PSH07-75] - BSA and Azide free	HA751187	Human	WB
Relative expression	TIM 3 Recombinant Rabbit mAb [PSH16-50]	HA723854	Human, Mouse, Rat	WB,IHC-P,IP
Relative expression	TIM 3 Recombinant Rabbit mAb [PSH10-00]	HA723169	Human	WB,IF-Cell,FC,IP
	TIM 3 Recombinant Rabbit mAb	IRS077	Human	mIHC
Relative expression	TIM 3 Recombinant Rabbit mAb [PSH07-75]	HA722891	Human	WB
	TIM 3 Mouse mAb [1E5]	EM1701-20	Human	WB,IF-Cell,IHC-P,FC
Relative expression	TIM 3 Recombinant Antibody [PSH13-95] - Mouse IgG1 (Chimeric)	HA601442	Human	WB,IF-Cell,FC
	TIMD4 Rabbit pAb	ER61913	Human, Mouse	WB
	Tin2 Recombinant Rabbit mAb [JE35-82]	HA722432	Human, Mouse, Rat	WB
	TPP1 Recombinant Rabbit mAb [PSH18-09]	HA723972	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
	Transmembrane protein 30A Rabbit pAb	ER1918-01	Mouse, Rat	IHC-P,IF-Tissue
	TYRO3 Recombinant Rabbit mAb [JE38-13]	HA721539	Human, Rat	WB
	UQCRFS1 Recombinant Rabbit mAb [JG62-31]	ET7108-28	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	VPS11 Rabbit pAb	ER62498	Human, Mouse	WB
	VPS11 Recombinant Rabbit mAb [JE33-36]	HA722447	Human	WB
	VPS18 Recombinant Rabbit mAb [JE34-83]	HA722325	Human, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层