

2026年8月—新上市产品清单

产品名	货号	Purity	Source
Aequorea victoria GFP, C-His Tag Protein	HA211504	>95% as determined by SDS-PAGE.	E.coli
AIM2 Recombinant Rabbit mAb [PSH27-15]	HA724676	Protein A affinity purified.	
AIM2 Recombinant Rabbit mAb [PSH27-15] - BSA and Azide free	HA752068	Protein A affinity purified.	
ATP6V1C1 Recombinant Rabbit mAb [PSH27-11] - BSA and Azide free	HA752064	Protein A affinity purified.	
ATP6V1C1 Recombinant Rabbit mAb [PSH27-11]	HA724672	Protein A affinity purified.	
BubR1 Recombinant Rabbit mAb [PSH27-14]	HA724675	Protein A affinity purified.	
BubR1 Recombinant Rabbit mAb [PSH27-14] - BSA and Azide free	HA752067	Protein A affinity purified.	
Cathepsin S Recombinant Rabbit mAb [PSH27-55] - BSA and Azide free	HA752078	Protein A affinity purified.	
Cathepsin S Recombinant Rabbit mAb [PSH27-55]	HA724704	Protein A affinity purified.	
CCP110 Recombinant Rabbit mAb [PSH27-13]	HA724674	Protein A affinity purified.	
CCP110 Recombinant Rabbit mAb [PSH27-13] - BSA and Azide free	HA752066	Protein A affinity purified.	
CD34 Recombinant Rabbit mAb [PSH27-02] - BSA and Azide free	HA752059	Protein A affinity purified.	
CD34 Recombinant Rabbit mAb [PSH27-02]	HA724663	Protein A affinity purified.	
CD8 alpha Recombinant Rabbit mAb [PSH27-01] - BSA and Azide free	HA752058	Protein A affinity purified.	
CD8 alpha Recombinant Rabbit mAb [PSH27-01]	HA724662	Protein A affinity purified.	
Choline Acetyltransferase Recombinant Antibody [PSH15-05] - Mouse IgG1 (Chimeric)	HA601668	Protein A affinity purified.	
Choline Acetyltransferase Recombinant Antibody [PSH15-05] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610450	Protein A affinity purified.	
c-Myb Recombinant Rabbit mAb [PSH27-17]	HA724678	Protein A affinity purified.	
c-Myb Recombinant Rabbit mAb [PSH27-17] - BSA and Azide free	HA752070	Protein A affinity purified.	

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cynomolgus Apolipoprotein E, Tag Free	HA211543	>95% as determined by SDS-PAGE.	HEK293
Cystatin C Recombinant Rabbit mAb [PSH27-39] - BSA and Azide free (Detector)	HA725450	Protein A affinity purified.	
Cystatin C Recombinant Rabbit mAb [PSH27-38] - BSA and Azide free (Capture)	HA725449	Protein A affinity purified.	
Biotin Conjugated Cystatin C Recombinant Rabbit mAb [PSH27-39] - Detector	HA725450B	Protein A affinity purified.	
Desmoglein 3 Recombinant Rabbit mAb [PSH27-06]	HA724667	Protein A affinity purified.	
Desmoglein 3 Recombinant Rabbit mAb [PSH27-06] - BSA and Azide free	HA752063	Protein A affinity purified.	
Desmoglein 3 Recombinant Rabbit mAb [PSH27-07]	HA724668	Protein A affinity purified.	
ENO1 Recombinant Rabbit mAb [PSH26-93] - BSA and Azide free	HA752054	Protein A affinity purified.	
ENO1 Recombinant Rabbit mAb [PSH26-93]	HA724658	Protein A affinity purified.	
ENO1/ENO2 Recombinant Rabbit mAb [PSH26-92]	HA724657	Protein A affinity purified.	
ENO1/ENO2 Recombinant Rabbit mAb [PSH26-92] - BSA and Azide free	HA752053	Protein A affinity purified.	
EpCAM Recombinant Antibody [PS01-69] - Mouse IgG1 (Chimeric)	HA601661	Protein A affinity purified.	
EpCAM Recombinant Antibody [PS01-69] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610447	Protein A affinity purified.	
ERCC1 Recombinant Rabbit mAb [PSH27-52]	HA724701	Protein A affinity purified.	
FGFR2 Recombinant Rabbit mAb [PSH26-90]	HA724655	Protein A affinity purified.	
FGFR2 Recombinant Rabbit mAb [PSH26-90] - BSA and Azide free	HA752052	Protein A affinity purified.	
FOXF1 Recombinant Rabbit mAb [PSH27-47]	HA724698	Protein A affinity purified.	
FOXF1 Recombinant Rabbit mAb [PSH27-47] - BSA and Azide free	HA752074	Protein A affinity purified.	
GLP-1 (7-36, amide) peptide	HA211605		
Helicobacter pylori Recombinant Rabbit mAb [PSH27-22]	HA724681	Protein A affinity purified.	
Histone H2A (mono methyl K118) Recombinant Rabbit mAb [PSH27-21]	HA724680	Protein A affinity purified.	
Human ADAMTS13, Tag Free Protein	HA211544	>95% as determined by SDS-PAGE.	HEK293

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Human Albumin Recombinant Rabbit mAb [PSH27-50] - BSA and Azide free (Detector)	HA725336	Protein A affinity purified.	
Biotin Conjugated Human Albumin Recombinant Rabbit mAb [PSH27-50] - Detector	HA725336B	Protein A affinity purified.	
Human Albumin Recombinant Rabbit mAb [PSH27-49] - BSA and Azide free (Capture)	HA725335	Protein A affinity purified.	
Human ALK-7, Tag Free Protein	HA211563	>95% as determined by SDS-PAGE.	HEK293
Biotin Conjugated Human ATG7 Recombinant Rabbit mAb [PSH26-95] - Detector	HA725494B	Protein A affinity purified.	
Human ATG7 Recombinant Rabbit mAb [PSH26-94] - BSA and Azide free (Capture)	HA725493	Protein A affinity purified.	
Human ATG7 Recombinant Rabbit mAb [PSH26-95] - BSA and Azide free (Detector)	HA725494	Protein A affinity purified.	
Human ATG7, Tag Free Protein	HA211607	>95% as determined by SDS-PAGE.	E.coli
Human CA125/MUC16 enzyme-linked immunoassay kit	EH0282		
Human Caspase-3/CASP3, C-His Tag Protein	HA211516	>95% as determined by SDS-PAGE.	E.coli
Human CCL16/LEC enzyme-linked immunoassay kit (one step)	EHY0294		
Human CCL24/Eotaxin-2/MPIF-2 enzyme-linked immunoassay kit	EH0283		
HRP Conjugated Human CD163 Recombinant Rabbit mAb [PSH13-79] - Detector	HA723568H	Protein A affinity purified.	
Human CD40 enzyme-linked immunoassay kit	EH0280		
Human CD5L, Tag Free Protein	HA211560	>95% as determined by SDS-PAGE.	HEK293
Human CNP/CNPase, Tag Free Protein	HA211541	>95% as determined by SDS-PAGE.	HEK293
HRP Conjugated Human Complement C2 Recombinant Rabbit mAb [PSH26-58] - Detector	HA725526H	Protein A affinity purified.	
HRP Conjugated Human CXCL16 Recombinant Rabbit mAb [PSH21-55] - Detector	HA725364H	Protein A affinity purified.	
Biotin Conjugated Human Cystatin SN Recombinant Rabbit mAb [PSH27-68] - Detector	HA724713B	Protein A affinity purified.	
Human Cystatin SN Recombinant Rabbit mAb [PSH27-68] - BSA and Azide free (Detector)	HA724713	Protein A affinity purified.	
Human Cystatin SN Recombinant Rabbit mAb [PSH27-67] - BSA and Azide free (Capture)	HA724712	Protein A affinity purified.	
Human Desmoglein 3 Recombinant Rabbit mAb [PSH27-30] - BSA and Azide free (Capture)	HA724686	Protein A affinity purified.	
Human Desmoglein 3 Recombinant Rabbit mAb [PSH27-31] - BSA and Azide free (Detector)	HA724687	Protein A affinity purified.	

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Biotin Conjugated Human Desmoglein 3 Recombinant Rabbit mAb [PSH27-31] - Detector	HA724687B	Protein A affinity purified.	
Human DNAM-1/CD226, Tag Free Protein	HA211629	>95% as determined by SDS-PAGE.	HEK293
Human DR6 Recombinant Rabbit mAb [PSH27-28] - BSA and Azide free (Capture)	HA724684	Protein A affinity purified.	
Human DR6 Recombinant Rabbit mAb [PSH27-29] - BSA and Azide free (Detector)	HA724685	Protein A affinity purified.	
Biotin Conjugated Human DR6 Recombinant Rabbit mAb [PSH27-29] - Detector	HA724685B	Protein A affinity purified.	
Human Dystrophin/DMD (CTD), C-His Tag Protein	HA211606	>95% as determined by SDS-PAGE.	HEK293
Human EpCAM / TROP1, C-His Tag Protein	HA211540	>95% as determined by SDS-PAGE.	HEK293
HRP Conjugated Human FABP4 Recombinant Rabbit mAb [PSH25-38] - Detector	HA725473H	Protein A affinity purified.	
Human FABP5/E-FABP, Tag Free Protein	HA211542	>95% as determined by SDS-PAGE.	HEK293
Human Factor XI, C-His Tag Protein	HA211612	>95% as determined by SDS-PAGE.	HEK293
Human Fas Ligand/TNFSF6 enzyme-linked immunoassay kit	EH0291		
HRP Conjugated Human Fetuin A Recombinant Rabbit mAb [PSH19-13] - Detector	HA724047H	Protein A affinity purified.	
Biotin Conjugated Human FGF7 Recombinant Rabbit mAb [PSH27-20] - Detector	HA725528B	Protein A affinity purified.	
Human FGF7 Recombinant Rabbit mAb [PSH27-20] - BSA and Azide free (Detector)	HA725528	Protein A affinity purified.	
Human FGF7 Recombinant Rabbit mAb [PSH27-19] - BSA and Azide free (Capture)	HA725527	Protein A affinity purified.	
Human FUS, C-His Tag Protein	HA211549	>95% as determined by SDS-PAGE.	HEK293
Human Galectin3 enzyme-linked immunoassay kit	EH0286		
HRP Conjugated Human GC / Vitamin D binding protein Recombinant Rabbit mAb [PSH26-30] - Detector	HA725374H	Protein A affinity purified.	
Human GLP1R Recombinant Rabbit mAb [PSH25-53] - BSA and Azide free (Capture)	HA724544	Protein A affinity purified.	
Human GLP1R Recombinant Rabbit mAb [PSH25-55] - BSA and Azide free (Capture)	HA724546	Protein A affinity purified.	
Biotin Conjugated Human GLP1R Recombinant Rabbit mAb [PSH25-54] - Detector	HA724545B	Protein A affinity purified.	
Human GLP1R Recombinant Rabbit mAb [PSH25-56] - BSA and Azide free (Detector)	HA724547	Protein A affinity purified.	
Biotin Conjugated Human GLP1R Recombinant Rabbit mAb [PSH25-56] - Detector	HA724547B	Protein A affinity purified.	

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Human GLP1R Recombinant Rabbit mAb [PSH25-54] - BSA and Azide free (Detector)	HA724545	Protein A affinity purified.	
HRP Conjugated Human Hemopexin Recombinant Rabbit mAb [PSH22-71] - Detector	HA724315H	Protein A affinity purified.	
Human HIF-1 alpha/HIF1A, Tag Free Protein	HA211604	>95% as determined by SDS-PAGE.	HEK293
Human HRAS, Tag Free Protein	HA211537	>95% as determined by SDS-PAGE.	E.coli
Human I-309/CCL1 enzyme-linked immunoassay kit	EH0287		
HRP Conjugated Human ICAM2 Recombinant Rabbit mAb [PSH25-48] - Detector	HA724539H	Protein A affinity purified.	
Biotin Conjugated Human IgA Recombinant Mouse mAb [5-C3-R] - Detector	HA601667B	Protein A affinity purified.	
Human IgA Recombinant Rabbit mAb [JM10-42] - BSA and Azide free (Capture)	HA724694	Protein A affinity purified.	
Human IgA Recombinant Mouse mAb [5-C3-R] - BSA and Azide free (Detector)	HA601667	Protein A affinity purified.	
HRP Conjugated Human IGFBP6 Recombinant Rabbit mAb [PSH22-19] - Detector	HA724280H	Protein A affinity purified.	
Human IL-19 enzyme-linked immunoassay kit	EH0285		
Human IL-34 enzyme-linked immunoassay kit	EH0284		
HRP Conjugated Human LDL Receptor Recombinant Rabbit mAb [PSH24-60] - Detector	HA724456H	Protein A affinity purified.	
Human MCP-2/CCL8 enzyme-linked immunoassay kit	EH0277		
Human MIP-5/CCL15 enzyme-linked immunoassay kit	EH0268		
HRP Conjugated Human/Monkey Alpha-1-Microglobulin Recombinant Rabbit mAb [PSH25-82] - Detector	HA724574H	Protein A affinity purified.	
HRP Conjugated Human/Monkey Apo-M Recombinant Rabbit mAb [PSH24-29] - Detector	HA724436H	Protein A affinity purified.	
Human/Monkey FGFR2 Recombinant Rabbit mAb [PSH27-34] - BSA and Azide free (Capture)	HA724690	Protein A affinity purified.	
Biotin Conjugated Human/Monkey FGFR2 Recombinant Rabbit mAb [PSH27-35] - Detector	HA724691B	Protein A affinity purified.	
Human/Monkey FGFR2 Recombinant Rabbit mAb [PSH27-35] - BSA and Azide free (Detector)	HA724691	Protein A affinity purified.	
HRP Conjugated Human/Monkey IGFBP2 Recombinant Rabbit mAb [PSH25-34] - Detector	HA724525H	Protein A affinity purified.	
HRP Conjugated Human/Monkey Mannose Receptor(CD206) Recombinant Rabbit mAb [PSH23-47] - Detector	HA724378H	Protein A affinity purified.	
HRP Conjugated Human/Monkey Pro-COL1A1 Recombinant Rabbit mAb [PSH25-16] - Detector	HA724510H	Protein A affinity purified.	
Human/Monkey SPINT1 / HAI-1 Recombinant Rabbit mAb [PSH27-69] - BSA and Azide free (Capture)	HA724714	Protein A affinity purified.	

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Biotin Conjugated Human/Monkey SPINT1 / HAI-1 Recombinant Rabbit mAb [PSH27-70] - Detector	HA724715B	Protein A affinity purified.	
Human/Monkey SPINT1 / HAI-1 Recombinant Rabbit mAb [PSH27-70] - BSA and Azide free (Detector)	HA724715	Protein A affinity purified.	
Biotin Conjugated Human/Mouse APLP1 Recombinant Rabbit mAb [PSH27-41] - Detector	HA725514B	Protein A affinity purified.	
Human/Mouse APLP1 Recombinant Rabbit mAb [PSH27-41] - BSA and Azide free (Detector)	HA725514	Protein A affinity purified.	
Human/Mouse APLP1 Recombinant Rabbit mAb [PSH27-40] - BSA and Azide free (Capture)	HA725513	Protein A affinity purified.	
Human Neuropeptide Y enzyme-linked immunoassay kit	EH0292		
Human Oncostatin M / OSM Recombinant Rabbit mAb [PSH27-60] - BSA and Azide free (Detector)	HA725530	Protein A affinity purified.	
Human Oncostatin M / OSM Recombinant Rabbit mAb [PSH27-59] - BSA and Azide free (Capture)	HA725529	Protein A affinity purified.	
Biotin Conjugated Human Oncostatin M / OSM Recombinant Rabbit mAb [PSH27-60] - Detector	HA725530B	Protein A affinity purified.	
Human OSMR Recombinant Rabbit mAb [PSH27-42] - BSA and Azide free (Capture)	HA725531	Protein A affinity purified.	
Biotin Conjugated Human OSMR Recombinant Rabbit mAb [PSH27-43] - Detector	HA725532B	Protein A affinity purified.	
Human OSMR Recombinant Rabbit mAb [PSH27-43] - BSA and Azide free (Detector)	HA725532	Protein A affinity purified.	
Human Peptide YY (PYY) (1-36)	HA211627		
HRP Conjugated Human Plasminogen Recombinant Rabbit mAb [PSH20-22] - Detector	HA724133H	Protein A affinity purified.	
Human PLGF enzyme-linked immunoassay kit	EH0288		
Human/Rat/Monkey BDNF Recombinant Rabbit mAb [PSH27-10] - BSA and Azide free (Detector)	HA724671	Protein A affinity purified.	
Human/Rat/Monkey BDNF Recombinant Rabbit mAb [PSH27-09] - BSA and Azide free (Detector)	HA724670	Protein A affinity purified.	
Biotin Conjugated Human/Rat/Monkey BDNF Recombinant Rabbit mAb [PSH27-09] - Detector	HA724670B	Protein A affinity purified.	
Human/Rat/Monkey BDNF Recombinant Rabbit mAb [PSH27-08] - BSA and Azide free (Capture)	HA724669	Protein A affinity purified.	
Biotin Conjugated Human/Rat/Monkey BDNF Recombinant Rabbit mAb [PSH27-10] - Detector	HA724671B	Protein A affinity purified.	
HRP Conjugated Human SCGB1A1 Recombinant Rabbit mAb [PSH23-84] - Detector	HA725471H	Protein A affinity purified.	
Human TACE/ADAM17, Tag Free Protein	HA211625	>95% as determined by SDS-PAGE.	HEK293
Human TIA1 Recombinant Rabbit mAb [PSH26-97] - BSA and Azide free (Detector)	HA725520	Protein A affinity purified.	

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Biotin Conjugated Human TIA1 Recombinant Rabbit mAb [PSH26-97] - Detector	HA725520B	Protein A affinity purified.	
Human TIA1 Recombinant Rabbit mAb [PSH26-96] - BSA and Azide free (Capture)	HA725519	Protein A affinity purified.	
Human TNFSF14 Recombinant Rabbit mAb [PSH27-26] - BSA and Azide free (Capture)	HA725533	Protein A affinity purified.	
Biotin Conjugated Human TNFSF14 Recombinant Rabbit mAb [PSH27-27] - Detector	HA725534B	Protein A affinity purified.	
Human TNFSF14 Recombinant Rabbit mAb [PSH27-27] - BSA and Azide free (Detector)	HA725534	Protein A affinity purified.	
Human TSH alpha/beta Heterodimer, Tag Free Protein	HA211624	>95% as determined by SDS-PAGE.	HEK293
IGF1R Recombinant Rabbit mAb [PSH26-79]	HA724644	Protein A affinity purified.	
IGF1R Recombinant Rabbit mAb [PSH26-79] - BSA and Azide free	HA752046	Protein A affinity purified.	
INSM1 Recombinant Rabbit mAb [PSH27-54]	HA724703	Protein A affinity purified.	
LINE-1 ORF1p Recombinant Rabbit mAb [PSH27-18]	HA724679	Protein A affinity purified.	
LINE-1 ORF1p Recombinant Rabbit mAb [PSH27-18] - BSA and Azide free	HA752071	Protein A affinity purified.	
L-Lactyl-Histone H3 (K18) Recombinant Rabbit mAb [PSH26-77] - BSA and Azide free	HA752044	Protein A affinity purified.	
L-Lactyl-Histone H3 (K18) Recombinant Rabbit mAb [PSH26-77]	HA724642	Protein A affinity purified.	
M13 Recombinant Mouse mAb [A5B3-R]	HA601665	Protein A affinity purified.	
MED1 Recombinant Rabbit mAb [PSH26-99]	HA724660	Protein A affinity purified.	
MED1 Recombinant Rabbit mAb [PSH26-99] - BSA and Azide free	HA752056	Protein A affinity purified.	
MORC2 Recombinant Rabbit mAb [PSH27-51] - BSA and Azide free	HA752076	Protein A affinity purified.	
MORC2 Recombinant Rabbit mAb [PSH27-51]	HA724700	Protein A affinity purified.	
Mouse Albumin/ALB, Tag Free Protein	HA211622	>95% as determined by SDS-PAGE.	HEK293
HRP Conjugated Mouse CD163 Recombinant Rabbit mAb [PSH22-48] - Detector	HA724298H	Protein A affinity purified.	
Mouse Chitinase-like protein 3 / YM1, C-His Tag	HA211603	>95% as determined by SDS-PAGE.	HEK293
HRP Conjugated Mouse CTGF Recombinant Rabbit mAb [PSH25-95] - Detector	HA725153H	Protein A affinity purified.	

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Mouse CXCL16 enzyme-linked immunoassay kit	EM0078		
Mouse DLL1, Tag Free Protein	HA211623	>95% as determined by SDS-PAGE.	HEK293
Mouse Dopamine beta Hydroxylase Recombinant Rabbit mAb [PSH27-37] - BSA and Azide free (Detector)	HA724693	Protein A affinity purified.	
Mouse Dopamine beta Hydroxylase Recombinant Rabbit mAb [PSH27-36] - BSA and Azide free (Capture)	HA724692	Protein A affinity purified.	
Biotin Conjugated Mouse Dopamine beta Hydroxylase Recombinant Rabbit mAb [PSH27-37] - Detector	HA724693B	Protein A affinity purified.	
Mouse IGFBP7, C-His Tag Protein	HA211539	>95% as determined by SDS-PAGE.	HEK293
Mouse IL-13, Tag Free Protein	HA211619	>95% as determined by SDS-PAGE.	HEK293
HRP Conjugated Mouse Lipocalin-2 Recombinant Rabbit mAb [PSH21-12] - Detector	HA725400H	Protein A affinity purified.	
Mouse M-CSF R / CSF1R / CD115, C-His Tag Protein	HA211533	>95% as determined by SDS-PAGE.	HEK293
HRP Conjugated Mouse Osteocalcin Recombinant Rabbit mAb [PSH22-74] - Detector	HA724318H	Protein A affinity purified.	
Mouse PAI-1/Serpin E1, Tag Free Protein	HA211562	>95% as determined by SDS-PAGE.	HEK293
HRP Conjugated Mouse/Rat IGFBP1 Recombinant Rabbit mAb [PSH26-22] - Detector	HA724607H	Protein A affinity purified.	
HRP Conjugated Mouse/Rat IGFBP6 Recombinant Rabbit mAb [PSH22-23] - Detector	HA724284H	Protein A affinity purified.	
Mouse/Rat Von Willebrand Factor/vWF enzyme-linked immunoassay kit (one step)	EMY0079		
HRP Conjugated Mouse SAA Recombinant Rabbit mAb [PSH25-07] - Detector	HA724501H	Protein A affinity purified.	
Mouse SP-D Recombinant Rabbit mAb [PSH27-33] - BSA and Azide free (Detector)	HA724689	Protein A affinity purified.	
Biotin Conjugated Mouse SP-D Recombinant Rabbit mAb [PSH27-33] - Detector	HA724689B	Protein A affinity purified.	
Mouse SP-D Recombinant Rabbit mAb [PSH27-32] - BSA and Azide free (Capture)	HA724688	Protein A affinity purified.	
Mouse SP-D/SFTPD, C-His Tag Protein	HA211630	>95% as determined by SDS-PAGE.	HEK293
MUC6 Recombinant Rabbit mAb [PSH26-89]	HA724654	Protein A affinity purified.	
MUC6 Recombinant Rabbit mAb [PSH26-89] - BSA and Azide free	HA752051	Protein A affinity purified.	
NAMPT Recombinant Rabbit mAb [PSH27-56] - BSA and Azide free (Capture)	HA725415	Protein A affinity purified.	

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Biotin Conjugated NAMPT Recombinant Rabbit mAb [PSH27-57] - Detector	HA725416B	Protein A affinity purified.	
NAMPT Recombinant Rabbit mAb [PSH27-57] - BSA and Azide free (Detector)	HA725416	Protein A affinity purified.	
NGFR Recombinant Rabbit mAb [PSH27-65] - BSA and Azide free (Capture)	HA724710	Protein A affinity purified.	
Biotin Conjugated NGFR Recombinant Rabbit mAb [PSH27-66] - Detector	HA724711B	Protein A affinity purified.	
NGFR Recombinant Rabbit mAb [PSH27-66] - BSA and Azide free (Detector)	HA724711	Protein A affinity purified.	
NPTXR Recombinant Rabbit mAb [PSH26-98] - BSA and Azide free	HA752055	Protein A affinity purified.	
NPTXR Recombinant Rabbit mAb [PSH26-98]	HA724659	Protein A affinity purified.	
pan-RFP Rabbit pAb	HA500643	Immunogen affinity purified.	
Parkin Recombinant Antibody [PSH08-11] - Guinea pig IgG2 (Chimeric)	HA601663	Protein A affinity purified.	
Parkin Recombinant Antibody [PSH08-11] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610449	Protein A affinity purified.	
Parkin Recombinant Antibody [PSH08-11] - Rat IgG1 (Chimeric)	HA601664	Protein A affinity purified.	
Parkin Recombinant Antibody [PSH08-11] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610448	Protein A affinity purified.	
PDGFR alpha Recombinant Rabbit mAb [PSH27-00]	HA724661	Protein A affinity purified.	
PDGFR alpha Recombinant Rabbit mAb [PSH27-00] - BSA and Azide free	HA752057	Protein A affinity purified.	
Perlecan Recombinant Rabbit mAb [PSH26-88] - BSA and Azide free	HA752050	Protein A affinity purified.	
Perlecan Recombinant Rabbit mAb [PSH26-88]	HA724653	Protein A affinity purified.	
Phospho-CBL (S669) Recombinant Rabbit mAb [PSH26-91]	HA724656	Protein A affinity purified.	
Phospho-EGFR (Y1110) Recombinant Rabbit mAb [PSH27-58]	HA724705	Protein A affinity purified.	
Phospho-IRF3 (S396) Recombinant Rabbit mAb [PSH27-05]	HA724666	Protein A affinity purified.	
Phospho-IRF3 (S396) Recombinant Rabbit mAb [PSH27-05] - BSA and Azide free	HA752062	Protein A affinity purified.	
PIKFYVE Recombinant Rabbit mAb [PSH27-04]	HA724665	Protein A affinity purified.	
PIKFYVE Recombinant Rabbit mAb [PSH27-04] - BSA and Azide free	HA752061	Protein A affinity purified.	

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

PIWIL2 Recombinant Rabbit mAb [PSH27-16] - BSA and Azide free	HA752069	Protein A affinity purified.	
PIWIL2 Recombinant Rabbit mAb [PSH27-16]	HA724677	Protein A affinity purified.	
Puromycin Recombinant Rabbit mAb [PSH27-46]	HA724697	Protein A affinity purified.	
Puromycin Recombinant Rabbit mAb [PSH27-46] - BSA and Azide free	HA752073	Protein A affinity purified.	
Rabbit IL-1 beta/IL1B, Tag Free Protein	HA211538	>95% as determined by SDS-PAGE.	HEK293
Rat CD200, C-His Tag Protein	HA211548	>95% as determined by SDS-PAGE.	HEK293
RBMX Recombinant Rabbit mAb [PSH27-53]	HA724702	Protein A affinity purified.	
RBMX Recombinant Rabbit mAb [PSH27-53] - BSA and Azide free	HA752077	Protein A affinity purified.	
SOX17 Recombinant Rabbit mAb [PSH27-45] - BSA and Azide free	HA752072	Protein A affinity purified.	
SOX17 Recombinant Rabbit mAb [PSH27-45]	HA724696	Protein A affinity purified.	
STMN2 Recombinant Rabbit mAb [PSH27-12] - BSA and Azide free	HA752065	Protein A affinity purified.	
STMN2 Recombinant Rabbit mAb [PSH27-12]	HA724673	Protein A affinity purified.	
Symmetric Di-Methyl-Histone H3 (R17) Recombinant Rabbit mAb [PSH27-44]	HA724695	Protein A affinity purified.	
Synapsin 1 Recombinant Rabbit mAb [PSH27-62] - BSA and Azide free	HA752080	Protein A affinity purified.	
Synapsin 1 Recombinant Rabbit mAb [PSH27-62]	HA724707	Protein A affinity purified.	
Testosterone Recombinant Rabbit mAb [PSH27-25] - BSA and Azide free	HA724683	Protein A affinity purified.	
Biotin Conjugated Testosterone Recombinant Rabbit mAb [PSH27-25]	HA724683B	Protein A affinity purified.	
Thrombospondin 1 Recombinant Rabbit mAb [PSH27-63] - BSA and Azide free (Capture)	HA724708	Protein A affinity purified.	
Thrombospondin 1 Recombinant Rabbit mAb [PSH27-64] - BSA and Azide free (Detector)	HA724709	Protein A affinity purified.	
Biotin Conjugated Thrombospondin 1 Recombinant Rabbit mAb [PSH27-64] - Detector	HA724709B	Protein A affinity purified.	
TOX Recombinant Rabbit mAb [PSH27-48] - BSA and Azide free	HA752075	Protein A affinity purified.	
TOX Recombinant Rabbit mAb [PSH27-48]	HA724699	Protein A affinity purified.	

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Vimentin Recombinant Antibody - Mouse IgG1 (Chimeric)	HA601662	Protein A affinity purified.	
Ym-1+Ym-2 Recombinant Rabbit mAb [PSH27-03]	HA724664	Protein A affinity purified.	
Ym-1+Ym-2 Recombinant Rabbit mAb [PSH27-03] - BSA and Azide free	HA752060	Protein A affinity purified.	
人CD14磁珠分选试剂(阳选)	HAK21014		