

干细胞通路研究-产品清单

高级验证	产品名	货号	反应种属	验证应用
	Activin A Receptor Type IB Recombinant Rabbit mAb [JE38-48]	HA721550	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	Activin A Receptor Type IC Rabbit pAb	ER1903-22	Human, Mouse, Rat	WB,FC
	Activin Receptor Type IA Rabbit pAb	HA500117	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Activin Receptor Type IIB / ACVR2B Recombinant Rabbit mAb [JE33-52]	HA722725	Human, Mouse, Rat	WB
	ACTR-IB Rabbit pAb	ER62592	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	ACVR2A Mouse mAb [A7C8]	HA601091	Human	WB
	ACVR2A Rabbit pAb	HA500084	Human, Mouse, Rat	WB
	ACVR2B Mouse mAb [A5H6]	HA600066	Human, Mouse, Rat, Zebrafish	WB
	ACVR2B Mouse mAb [A5H7]	HA600067	Human, Mouse, Rat, Zebrafish	WB
KD	AKT1 Recombinant Rabbit mAb [ST05-09] - BSA and Azide free	HA750177	Human, Mouse, Rat	WB,IHC-P,IP,FC,IHC-Fr,IF-Cell
KD	AKT1 Recombinant Rabbit mAb [ST05-09]	ET1609-47	Human, Mouse, Rat	WB,IHC-P,IP,FC,IHC-Fr,IF-Cell
	AKT1 Mouse mAb [D9-9-C9]	EM40507	Human, Mouse, Rat	WB,IHC-P
	AKT1 Recombinant Rabbit mAb [PSH0-50] - BSA and Azide free	HA750627	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC
	AKT1 Recombinant Rabbit mAb [PSH0-50]	HA721364	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC
	AKT1/2/3 Recombinant Rabbit mAb [ST48-09] - BSA and Azide free	HA750179	Human, Mouse, Rat, Monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
	AKT1/2/3 Recombinant Rabbit mAb [ST48-09]	ET1609-51	Human, Mouse, Rat, Monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
KD	AKT1/2/3 Recombinant Rabbit mAb [JE75-09]	HA721870	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IP
	AKT1/2/3 Rabbit pAb	ER1901-33	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	AKT1/2/3 Recombinant Rabbit mAb [JE75-09] - BSA and Azide free	HA750840	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IP
	Akt1/3 Rabbit pAb	ER62638	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Akt2 Mouse mAb [2A1]	RT1027	Human, Mouse, Rat	WB,IP,IF,IHC-P
	AKT2 Rabbit pAb	HA500091	Human, Mouse, Rat	WB,IHC-P
	APC Rabbit pAb	ER1802-39	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue,IHC-Fr
	APC Recombinant Rabbit mAb [SC61-03]	ET1610-80	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IP,IHC-P,IHC-Fr
	APC2 Mouse mAb [3A2G2]	EM1710-75	Human	IHC-P,IF-Cell,FC
	ATBF1 Rabbit pAb	ER2001-49	Human, Mouse, Rat	WB,FC
	ATBF1 Rabbit pAb	HA500181	Human, Rat	WB,IHC-P
KD	Axin1 Recombinant Rabbit mAb [PSH03-00] - BSA and Azide free	HA750855	Human	WB,IHC-P
KD	Axin1 Recombinant Rabbit mAb [PSH03-00]	HA721919	Human	WB,IHC-P
	Axin2 Recombinant Rabbit mAb [JM11-30]	ET1703-96	Human, Mouse, Rat	WB,IHC-P
	Beta Catenin Mouse mAb [A6-F8]	EM0306	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue,ChIP
KD	Beta Catenin Recombinant Rabbit mAb [SA30-04]	ET1601-5	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IP,mIHC,IF-Cell,IHC-Fr,FC
	iFluor™ 488 Conjugated Beta Catenin Recombinant Rabbit mAb [SA30-04]	HA720182F	Human, Mouse	IF,ICC
	Beta Catenin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS278RT	Human	mIHC
	Beta Catenin Recombinant Mouse mAb [A6-F8-R] - BSA and Azide free	HA610062	Human, Mouse, Rat	WB,IF-Cell,FC
KD	Beta Catenin Recombinant Rabbit mAb [SA30-04] - BSA and Azide free	HA750023	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IP,IF-Cell,IHC-Fr,FC
	Beta Catenin Mouse mAb [10-C0-B7]	M1405-6	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Beta Catenin Rabbit pAb	ER0805	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell,IF-Tissue
	Beta Catenin Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS277MS	Human	mIHC
	Beta Catenin Recombinant Mouse mAb [A6-F8-R]	HA601179	Human, Mouse, Rat	WB,IF-Cell,FC
	Beta Catenin Rabbit pAb	0407-16	Human, Mouse, Rat, Zebrafish	WB,IF-Cell
	Beta Catenin Recombinant Rabbit mAb	IRS043	Human	mIHC
	Beta Catenin Recombinant Rabbit mAb	IRS061	Mouse	mIHC
	bFGF Rabbit pAb	ER1904-22	Human, Mouse	WB,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Bmi1 Mouse mAb [A11-C11]	EM1701-55	Human, Mouse	WB,IHC-P,IF-Cell
	Bmi1 Recombinant Rabbit mAb [JJ093-3]	ET1701-89	Human, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	Bmi1 Mouse mAb [B3-G5]	EM20602	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	BMP4 Mouse mAb [1C-17]	EM1706-23	Human, Rat	WB,IF-Cell,FC
	BMP4 Recombinant Rabbit mAb [JE42-44]	HA721953	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue
	BMPR1A Rabbit pAb	R1510-1	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	BMPR1B Rabbit pAb	R1311-4	Human, Mouse	WB,IF-Cell,IHC-P
	BMPR2 Mouse mAb [F8-G6]	EM1709-92	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P
	CD130 / GP130 Mouse mAb [B11-B9-A8]	EM1701-43	Human, Mouse	WB,IHC-P,IF-Cell
	c-Myc Recombinant Rabbit mAb [PS00-11]	HA721182	Human, Mouse, Rat	WB,IHC-P,IP,ChIP
	c-Myc Recombinant Rabbit mAb [JE01-20]	HA722895	Human, Mouse, Rat, Monkey	WB,IF-Cell
	c-Myc Mouse mAb [6G6]	RT1149	Human, Mouse, Rat, Monkey	WB,IP,IF,IHC-P
	c-Myc Rabbit pAb	0912-2	Human	WB
	HRP Conjugated c-Myc Rabbit pAb	0912-4	Human	WB
	c-Myc Recombinant Rabbit mAb [PS00-11] - BSA and Azide free	HA750519	Human, Mouse, Rat	WB,IHC-P,IP,ChIP
	Dishevelled 2 Recombinant Rabbit mAb [PSH04-07]	HA722088	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Dishevelled 2 Recombinant Rabbit mAb [PSH04-07] - BSA and Azide free	HA750930	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Dishevelled 3 Recombinant Rabbit mAb [PSH15-93] - BSA and Azide free	HA751584	Human, Mouse, Rat, Monkey	WB,IP
	Dishevelled 3 Recombinant Rabbit mAb [PSH15-93]	HA723798	Human, Mouse, Rat, Monkey	WB,IP
	Dlx-5 Rabbit pAb	ER63256	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	DVL1 Rabbit pAb	ER1907-85	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Dvl-2 Rabbit pAb	ER63295	Human, Mouse	WB,IHC-P
	Dvl-3 Rabbit pAb	ER63296	Human, Mouse	WB,IHC-P
KD	ERK1 Recombinant Rabbit mAb [SP00-74]	ET1604-32	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	ERK1 Mouse mAb [7-D6-E5]	M1505-9	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	ERK1 Recombinant Rabbit mAb [SP05-09]	ET1604-16	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Erk1/2 Mouse mAb [2G1]	RT1484	Human, Mouse, Rat	WB,IP,IF
KD	ERK1/2 Recombinant Rabbit mAb [SA43-03] - BSA and Azide free	HA750019	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF- Tissue,IP,FC,IHC- P,IHC-Fr
	ERK1/2 Rabbit pAb	ER131011	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	ERK1/2 Recombinant Rabbit mAb [SA43-03]	ET1601-29	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF- Tissue,IP,FC,IHC- P,IHC-Fr
KD	ERK2 Recombinant Mouse mAb [5-D2-R]	HA601202	Human, Mouse, Rat	WB
KD	ERK2 Mouse mAb [B10-8]	EM1901-53	Human, Mouse, Rat	WB,IHC-P,FC
KD	ERK2 Recombinant Rabbit mAb [SZ25-01]	ET1603-23	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC- P,IP,FC,IHC-Fr
KD	ERK2 Mouse mAb [A10D2]	EM1901-54	Human, Mouse, Rat	WB,IHC-P,FC
	ERK2 Mouse mAb [5-D2]	M1407-3	Human, Mouse, Rat	WB,IF-Cell,IHC- P,FC
KD	ERK2 Rabbit pAb	ER131218	Human, Mouse, Rat	WB,IF-Cell,IHC- P,FC
	ERK2 Mouse mAb [C3B12]	EM1901-55	Human, Mouse, Rat	WB,IHC-P,FC
KD	ERK2 Recombinant Mouse mAb [5-D2-R] - BSA and Azide free	HA610080	Human, Mouse, Rat	WB
KD	ERK2 Recombinant Mouse mAb [C3B12-R]	HA601201	Human, Mouse, Rat	WB
KD	ERK2 Recombinant Mouse mAb [C3B12-R] - BSA and Azide free	HA610079	Human, Mouse, Rat	WB
	ESRRB Rabbit pAb	R1308-2	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	FGF2 Recombinant Rabbit mAb [JM28-10]	ET1703-18	Human, Mouse	WB,IP,FC,IHC-P
Relative expression, KD	FGFR1 Recombinant Rabbit mAb [SD08-25]	ET1612-62	Human	WB,IF-Cell
	FGFR2 Recombinant Rabbit mAb [JM10-60]	ET1703-78	Human, Mouse, Rat	WB,IP
	FGFR2/CD332 Rabbit pAb	R1311-7	Human	WB,IHC-P
	FGFR2/CD332 Rabbit pAb	R1310-10	Human	WB,IF-Cell,IHC- P,FC
	FGFR2/CD332 Recombinant Mouse mAb [A5-E10-G7-R]	EM50103	Human	WB
	FGFR2/CD332 Recombinant Mouse mAb [A5-E10-G7-R] - BSA and Azide free	HA610003	Human	WB
KD	FGFR3 Recombinant Rabbit mAb [JM110-33] - BSA and Azide free	HA750388	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC-P
KD	FGFR3 Recombinant Rabbit mAb [JM110-33]	ET1703-65	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	FGFR3 Recombinant Rabbit mAb [JM81-10]	ET1703-36	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
KD	FGFR3 Recombinant Rabbit mAb [JM81-10] - BSA and Azide free	HA750379	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	FGFR-4 Rabbit pAb	ER65688	Human	WB
	FGFR4 Recombinant Rabbit mAb [JE59-96]	HA721578	Human	WB,FC
	Frizzled-1 Rabbit pAb	ER63469	Human	WB,IF-Cell
	Frizzled 10/CD350 Rabbit pAb	ER1909-02	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Frizzled-2 Rabbit pAb	ER63470	Human, Mouse, Rat	WB,IF-Cell
	Frizzled-4 Rabbit pAb	ER65833	Human, Mouse, Rat	IHC-P
	Frizzled-5 Rabbit pAb	0905-3	Human	WB,IF-Cell,FC
	Frizzled-6 Rabbit pAb	ER63471	Human	WB,IHC-P,IF-Cell
	Frizzled 7 Rabbit pAb	ER1909-03	Human, Mouse, Rat	WB,FC
	Frizzled 7 Rabbit pAb	HA500433	Human	WB
	Frizzled 8 Recombinant Rabbit mAb [JB40-16]	ET7107-68	Human, Mouse	WB,IHC-P
	Frizzled-9 Rabbit pAb	ER63472	Human, Mouse, Monkey	WB,IF-Cell
	FZD2 Rabbit pAb	HA500406	Human, Mouse	WB,IHC-P
	FZD3 Rabbit pAb	ER1909-16	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	GRB2 Recombinant Rabbit mAb [JE49-73]	HA721923	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC
	GRB2 Rabbit pAb	ER31203	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	GRB2 Mouse mAb [11-B11]	EM1701-41	Human, Mouse, Rat	WB,IHC-P,IF-Cell
KD	GSK3 alpha Rabbit pAb	R1312-1	Human	WB,IF-Cell,IHC-P,FC
	GSK3 alpha Recombinant Rabbit mAb [JE33-27]	HA722217	Human, Mouse, Rat, Monkey	WB,IHC-P
	GSK3 alpha/beta Rabbit pAb	ER63648	Human, Mouse, Rat	WB,IHC-P
	GSK3 beta Rabbit pAb	R1510-4	Human, Mouse, Rat	IF-Cell,IHC-P,FC
KD	GSK3 beta Recombinant Rabbit mAb [SY28-03] - BSA and Azide free	HA750131	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IHC-Fr
KD	GSK3 beta Recombinant Rabbit mAb [SY28-03]	ET1607-71	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IHC-Fr
	GTPase HRAS Recombinant Rabbit mAb [JE42-15]	HA721355	Human, Mouse, Rat	WB,IHC-P
	HAND1 Rabbit pAb	ER63325	Human, Mouse, Rat	IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	HESX1 Rabbit pAb	ER60755	Human, Mouse	WB
	Hox-A1 Rabbit pAb	ER63703	Human, Mouse, Rat	WB,IHC-P
	Hox-A1/B1/D1 Rabbit pAb	ER63704	Human, Mouse	IHC-P
	Human FGF-2 Recombinant Rabbit mAb [PSH02-09] - BSA and Azide free (Detector)	HA721781	Human	ELISA(Det)
	Human FGF-2 Recombinant Rabbit mAb [PSH02-08] - BSA and Azide free (Capture)	HA721780	Human	ELISA(Cap)
	Biotin Conjugated Human FGF-2 Recombinant Rabbit mAb [PSH02-09]	HA722639B	Human	ELISA(Det)
	Human LIF Recombinant Rabbit mAb [PSH05-97] - BSA and Azide free (Detector)	HA722496	Human	ELISA(Det)
	Biotin Conjugated Human LIF Recombinant Rabbit mAb [PSH05-97]	HA722586B	Human	ELISA(Det)
	Human LIF Recombinant Rabbit mAb [PSH05-96] - BSA and Azide free (Capture)	HA722495	Human	ELISA(Cap)
	Biotin Conjugated Human WNT7A Recombinant Rabbit mAb [PSH19-62] - Detector	HA725314B	Human	ELISA(Det),ELISA
	Human WNT7A Recombinant Rabbit mAb [PSH19-62] - BSA and Azide free (Detector)	HA725314	Human	ELISA(Det)
	Human WNT7A Recombinant Rabbit mAb [PSH19-61] - BSA and Azide free (Capture)	HA725313	Human	ELISA(Cap)
	HXB1 Rabbit pAb	ER60763	Human	WB
	Id1 Recombinant Rabbit mAb	IRS159RB	Human	mIHC
	Id1 Recombinant Rabbit mAb [JM92-13]	ET1705-49	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Id1 Recombinant Rabbit mAb [JM92-13] - BSA and Azide free	HA750431	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
	ID2 Mouse mAb [A4-D4]	M1301-2	Human, Mouse	WB,IHC-P
	ID3 Rabbit pAb	ER62074	Human, Mouse, Rat	WB
	Id4 Rabbit pAb	ER63727	Human, Mouse, Rat	WB,IHC-P
Cell treatment	IGF1 Recombinant Rabbit mAb [JE30-01]	HA722527	Human	WB
	IGF1 (E peptide) Rabbit pAb	HA500107	Human	WB,IHC-P,FC
	IGF1R Rabbit pAb	R1310-6	Human	WB,IHC-P
	IGF1R Mouse mAb [D5-G2]	HA600033	Human	WB,FC
Relative expression	IGF1 Receptor beta Recombinant Rabbit mAb [PSH11-72]	HA723368	Human, Mouse, Rat	WB,IF-Cell
Relative expression	IGF1 Receptor beta Recombinant Rabbit mAb [PSH11-72] - BSA and Azide free	HA751427	Human, Mouse, Rat	WB,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	IGF I Rabbit pAb	ER1911-10	Human, Mouse, Rat, Dog	WB,IHC-P,IF-Tissue
	IGF-IR Rabbit pAb	ER63734	Human, Mouse, Rat	IF-Cell,WB,IHC-P
	Inhibin beta A Rabbit pAb	HA500271	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Inhibin beta A Rabbit pAb	ER1911-46	Human, Rat	WB,IF-Cell,FC
	Inhibin beta B Rabbit pAb	ER1911-47	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Inhibin beta B Rabbit pAb	HA500086	Human, Mouse	WB,IHC-P
	Inhibin beta-C Rabbit pAb	ER63761	Human	WB,IHC-P
	Inhibin beta-E Rabbit pAb	ER65538	Human, Mouse, Rat	WB,IHC-P
KD, Relative expression	Insulin Receptor Beta Recombinant Rabbit mAb [SN20-13]	ET1611-90	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC
KD, Relative expression	Insulin Receptor Beta Recombinant Rabbit mAb [SN20-13] - BSA and Azide free	HA750276	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC
	Islet 1 Recombinant Rabbit mAb [SC05-64] - BSA and Azide free	HA750211	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IHC-P,FC,IP,IHC-Fr
	Islet 1 Recombinant Rabbit mAb [SC05-64]	ET1610-33	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IHC-P,FC,IP,IHC-Fr
	JAK1 Recombinant Rabbit mAb [JM75-03]	ET1705-84	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell,IF-Tissue
	JAK1 Recombinant Rabbit mAb [JM75-03] - BSA and Azide free	HA750442	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell,IF-Tissue
	JAK1 Mouse mAb [1G1]	RT1341	Human, Mouse, Rat	WB,IP,IF,IHC-P
KD	JAK2 Mouse mAb [6-D3]	M1501-8	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	JAK2 Recombinant Rabbit mAb [SY0245] - BSA and Azide free	HA750118	Human, Mouse, Rat	WB,IF-Cell
KD	JAK2 Recombinant Rabbit mAb [SY0245]	ET1607-35	Human, Mouse, Rat	WB,IF-Cell
	JAK2 Rabbit pAb	ER1706-58	Human, Mouse, Rat	IF-Cell,IHC-P
	JAK2 Rabbit pAb	RT1343	Human, Mouse, Rat	WB,IP,IF,IHC-P
	JAK3 Recombinant Rabbit mAb [JM66-31]	ET1705-1	Human	WB,IHC-P
	Jarid2 Recombinant Rabbit mAb [PSH01-32] - BSA and Azide free	HA750710	Human, Mouse	IHC-P,IF-Cell
	Jarid2 Recombinant Rabbit mAb [PSH01-32]	HA721656	Human, Mouse	IHC-P,IF-Cell
	KAT6A Recombinant Rabbit mAb [PSH09-88] - BSA and Azide free	HA751322	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,ChIP
	KAT6A Recombinant Rabbit mAb [PSH09-88]	HA723157	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,ChIP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	KLF4 Recombinant Rabbit mAb [JF98-08]	ET1702-71	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	KLF4 Recombinant Rabbit mAb [JF98-08] - BSA and Azide free	HA750356	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	KLF4 Rabbit pAb	R1308-1	Human, Mouse, Rat	WB,IHC-P,FC
	KMT1E / SETDB1 Recombinant Rabbit mAb [PSH20-39]	HA724150	Human, Mouse, Rat, Monkey	WB,IHC-P,IP
	KMT1E / SETDB1 Recombinant Rabbit mAb [PSH20-40]	HA724151	Human, Mouse, Rat, Monkey	WB,IHC-P,IP
	K-Ras Rabbit pAb	ER65706	Human, Mouse, Rat	WB,IHC-P
	KRAS Recombinant Mouse mAb [A8E6] - BSA and Azide free	HA610008	Human, Mouse, Rat	WB,IHC-P
	KRAS Recombinant Mouse mAb [A8E6]	HA601060	Human, Mouse, Rat	WB,IHC-P
	KRAS Mouse mAb [A8E4]	HA601058	Human, Mouse, Rat	WB,IF-Cell,FC
	KRAS Mouse mAb [A8E5]	HA601059	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	LEFTY1+LEFTY2 Recombinant Rabbit mAb [JE38-88]	HA722177	Human, Mouse, Rat	WB
	LHX5 Rabbit pAb	ER60800	Human, Rat, Mouse	WB
	LIF Rabbit pAb	ER1912-32	Human, Mouse	WB
	LIFR Rabbit pAb	ER65635	Human	WB
	LIFR Rabbit pAb	ER1912-33	Human, Mouse, Rat	IHC-P,IF-Tissue
Relative expression	LIFR Recombinant Rabbit mAb [PSH08-76]	HA723028	Human	WB
Relative expression	LIFR Recombinant Rabbit mAb [PSH08-76] - BSA and Azide free	HA751241	Human	WB
	MEIS1 Recombinant Rabbit mAb [JA13-26]	ET1704-84	Human, Mouse, Rat	WB,IP
KD	MEK1 Recombinant Rabbit mAb [SZ22-01]	ET1603-20	Human, Mouse, Rat, Cow, Dog	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	MEK1/2 Recombinant Rabbit mAb [SR13-07] - BSA and Azide free	HA750042	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	MEK1/2 Recombinant Rabbit mAb [SR13-07]	ET1602-3	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	MEK2 Recombinant Rabbit mAb [SD207-08]	ET1612-6	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	MK13 Rabbit pAb	ER60146	Human, Mouse, Rat	WB
	MKP-4 Rabbit pAb	ER63965	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Biotin Conjugated Mouse GP130 / IL6ST Recombinant Rabbit mAb [PSH19-83] - Detector	HA724107B	Mouse	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse GP130 / IL6ST Recombinant Rabbit mAb [PSH19-83] - BSA and Azide free (Detector)	HA724107	Mouse	ELISA(Det)
	Mouse GP130 / IL6ST Recombinant Rabbit mAb [PSH19-82] - BSA and Azide free (Capture)	HA724106	Mouse	ELISA(Cap)
	Mouse LIF Recombinant Rabbit mAb [PSH02-00] - BSA and Azide free (Capture)	HA721758	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse LIF Recombinant Rabbit mAb [PSH02-01]	HA722643B	Mouse	ELISA(Det)
	Mouse LIF Recombinant Rabbit mAb [PSH02-01] - BSA and Azide free (Detector)	HA721759	Mouse	ELISA(Det)
	Myc tag Rabbit pAb	R1208-1	Species independent	WB,IP,IF-Cell,FC
	HRP Conjugated Myc tag Rabbit pAb	R1208-2	Species independent	WB
	Myf-5 Rabbit pAb	ER64072	Human, Mouse, Monkey	WB,IHC-P,IF-Cell
	Myf5 Recombinant Rabbit mAb [JE38-56]	HA722400	Mouse, Rat	WB
Relative expression	Nanog Recombinant Rabbit mAb [SC05-70]	ET1610-2	Human, Mouse	WB,IF-Tissue,FC,IHC-P,ChIP
Relative expression	Nanog Recombinant Rabbit mAb [SC05-70] - BSA and Azide free	HA750204	Human, Mouse	WB,IF-Tissue,FC,IHC-P,ChIP
	Nanog Rabbit pAb	ER1913-50	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Nanog Rabbit pAb	R1108-2	Human, Mouse	WB
	NGN1 Rabbit pAb	ER60918	Human	WB
	N/H/K-Ras Rabbit pAb	ER64095	Human, Mouse, Rat	WB,IHC-P
	Nodal Rabbit pAb	HA500041	Human, Rat	WB,IHC-P
	NRSF Rabbit pAb	ER1914-14	Human, Mouse	WB,IHC-P,IF-Tissue
KD	Oct4 Mouse mAb [A3-C11]	EM100306	Human, Mouse	WB,IF-Cell,FC
KO	Oct4 Recombinant Rabbit mAb [SD0750]	ET1612-20	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,ChIP
KD	Oct4 Recombinant Mouse mAb [A3-C11-R] - BSA and Azide free	HA610085	Human	WB,IF-Cell
KD	Oct4 Recombinant Mouse mAb [A3-C11-R]	HA601207	Human	WB,IF-Cell
	Oct4 Rabbit pAb	R1107-3	Human, Mouse	WB
	OCT-4 Recombinant Rabbit mAb	RE6048	Human	WB,IHC-P,FC,IP
	Otx1 Recombinant Rabbit mAb [JE42-10]	HA721881	Mouse, Rat, Human	WB
Cell treatment	Otx1 + Otx2 Recombinant Rabbit mAb [PSH01-62]	HA721705	Human, Rat, Mouse	WB,IHC-P,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Otx1 + Otx2 Recombinant Rabbit mAb [PSH01-62] - BSA and Azide free	HA750738	Human, Rat, Mouse	WB,IHC-P,IF-Cell,FC
Cell treatment	Otx2 Recombinant Rabbit mAb [PSH01-61] - BSA and Azide free	HA750737	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue,IF-Cell,FC
	Otx2 Recombinant Antibody [PSH01-61] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610239	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P
Cell treatment	Otx2 Recombinant Rabbit mAb [PSH01-61]	HA721704	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue,IF-Cell,FC
	Otx2 Recombinant Antibody [PSH01-61] - Rat IgG1 (Chimeric)	HA601399	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P
KO	p38 alpha / MAPK14 Recombinant Rabbit mAb [JF55-07]	ET1702-65	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	p38 alpha / MAPK14 Mouse mAb [9-B9]	EM40307	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Cell treatment	p38 alpha / MAPK14 Rabbit pAb	R1308-3	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD, Cell treatment	p38 alpha / MAPK14 Recombinant Rabbit mAb [SR43-04] - BSA and Azide free	HA750040	Human, Mouse, Rat	WB,IF-Cell,FC
KD, Cell treatment	p38 alpha / MAPK14 Recombinant Rabbit mAb [SR43-04]	ET1602-26	Human, Mouse, Rat	WB,IF-Cell,FC
	p38 alpha/MAPK14 + p38 beta/MAPK11 + p38 gamma/MAPK12 Recombinant Rabbit mAb [PSH03-60]	HA722019	Human, Mouse, Rat, Monkey	WB
	p38 alpha/MAPK14 + p38 beta/MAPK11 + p38 gamma/MAPK12 Recombinant Rabbit mAb [PSH03-60] - BSA and Azide free	HA750894	Human, Mouse, Rat, Monkey	WB
	p38 beta/MAPK11 + p38 gamma/MAPK12 + p38 delta/MAPK13 Recombinant Rabbit mAb [PSH03-59] - BSA and Azide free	HA750893	Human, Mouse, Rat, Monkey	WB
	p38 beta/MAPK11 + p38 gamma/MAPK12 + p38 delta/MAPK13 Recombinant Rabbit mAb [PSH03-59]	HA722018	Human, Mouse, Rat, Monkey	WB
	p38 delta / MAPK13 Rabbit pAb	HA500582	Human	WB
	p38 gamma/MAPK12 Recombinant Rabbit mAb [PSH03-58] - BSA and Azide free	HA750892	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	p38 gamma/MAPK12 Recombinant Rabbit mAb [PSH03-58]	HA722017	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	p38 MAPK Recombinant Rabbit Multiclonal Antibody [PSH03-72] - BSA and Azide free	HA750906	Human, Mouse, Rat, Monkey	WB
	p38 MAPK Recombinant Rabbit Multiclonal Antibody [PSH03-72]	HA722036	Human, Mouse, Rat, Monkey	WB
	PAX6 Recombinant Rabbit mAb	IRS066	Mouse	mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	PAX6 Recombinant Rabbit mAb [SD081-03]	ET1612-58	Human, Mouse, Rat	WB,IHC-P,FC,mIHC,IF-Cell,IF-Tissue,IHC-Fr
Relative expression	PAX6 Recombinant Rabbit mAb [SD08-31]	ET1612-89	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
Cell treatment	Phospho-AKT1 (S124) Recombinant Rabbit mAb [JJ08-46]	ET1701-36	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
Cell treatment	Phospho-AKT1 (S129) Recombinant Rabbit mAb [JE59-50]	HA721677	Human, Mouse, Rat	WB
Cell treatment	Phospho-AKT1 (T450) Recombinant Rabbit mAb [SD08-12]	ET1612-73	Human, Mouse, Rat	WB,IHC-P,IP,IF-Cell,FC
Cell treatment	Phospho-AKT2 (S474) Recombinant Rabbit mAb [PSH09-46]	HA723108	Human, Mouse, Rat	WB,FC
Cell treatment	Phospho-AKT2 (S474) Recombinant Rabbit mAb [PSH09-46] - BSA and Azide free	HA751293	Human, Mouse, Rat	WB,FC
Cell treatment	Phospho-AKT2 (S474) Rabbit pAb	HA500116	Human, Mouse	WB,IHC-P
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [PSH04-44]	HA722129	Human, Mouse, Rat	WB,IHC-P,FC
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [PSH04-44] - BSA and Azide free	HA750953	Human, Mouse, Rat	WB,IHC-P,FC
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [SY28-05]	ET1607-73	Human, Mouse, Rat, Dog	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [SY28-05] - BSA and Azide free	HA750132	Human, Mouse, Rat, Dog	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
Cell treatment	Phospho-AKT (T308) Recombinant Rabbit mAb [PSH08-10]	HA722951	Human	WB,IF-Cell,FC
Cell treatment	Phospho-AKT (T308) Recombinant Rabbit mAb [PSH08-10] - BSA and Azide free	HA751205	Human	WB,IF-Cell,FC
Cell treatment	Phospho-Beta Catenin (S33 + S37) Recombinant Rabbit mAb [JE59-59]	HA721580	Human, Rat, Mouse	WB
Cell treatment	Phospho-Beta Catenin (S552) Recombinant Rabbit mAb [PSH08-72]	HA723024	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Cell treatment	Phospho-Beta Catenin (S552) Recombinant Rabbit mAb [PSH08-72] - BSA and Azide free	HA751237	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Cell treatment	Phospho-Beta Catenin (T41/S45) Recombinant Rabbit mAb [JE54-02]	HA722316	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-c-Myc (S62) Recombinant Rabbit mAb [ST49-08] - BSA and Azide free	HA750186	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
Cell treatment	Phospho-c-Myc (S62) Recombinant Rabbit mAb [ST49-08]	ET1609-64	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
Cell treatment	Phospho-c-Myc (T58) Recombinant Rabbit mAb [SN60-01]	ET1611-24	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,FC,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-c-Myc (T58 + S62) Recombinant Rabbit mAb [SZ02-06]	ET1603-24	Human, Mouse	WB,IP,FC
Cell treatment	Phospho-Erk1 (T202)+Erk2 (T185) Recombinant Rabbit mAb [SZ2-4]	ET1603-22	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
Cell treatment	Phospho-Erk1 (T202)+Erk2 (T185) Recombinant Rabbit mAb [SZ2-4] - BSA and Azide free	HA750061	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
Cell treatment	Phospho-Erk1 (T202 + Y204) + Erk2 (T185 + Y187) Recombinant Rabbit mAb [SC58-01] - BSA and Azide free	HA750199	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP
Cell treatment	Phospho-Erk1 (T202 + Y204) + Erk2 (T185 + Y187) Recombinant Rabbit mAb [SC58-01]	ET1610-13	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP
Cell treatment	Phospho-GSK3 (alpha + beta) (Y216 + Y279) Recombinant Rabbit mAb [SY26-05] - BSA and Azide free	HA750127	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,IF-Tissue,FC
Cell treatment	Phospho-GSK3 (alpha + beta) (Y216 + Y279) Recombinant Rabbit mAb [SY26-05]	ET1607-54	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,IF-Tissue,FC
Cell treatment	Phospho-GSK3 alpha (S21) Recombinant Rabbit mAb [PSH09-11]	HA723069	Human, Mouse	WB,IHC-P,Dot Blot
Cell treatment	Phospho-GSK3 alpha (S21) Recombinant Rabbit mAb [PSH09-11] - BSA and Azide free	HA751265	Human, Mouse	WB,IHC-P,Dot Blot
Cell treatment	Phospho-GSK3 beta (S9) Recombinant Rabbit mAb [SY02-71] - BSA and Azide free	HA750128	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
Cell treatment	Phospho-GSK3 beta (S9) Recombinant Rabbit mAb [SY02-71]	ET1607-60	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
Cell treatment	Phospho-JAK2 (Y1007 + Y1008) Recombinant Rabbit mAb [SY24-03] - BSA and Azide free	HA750117	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
Cell treatment	Phospho-JAK2 (Y1007 + Y1008) Recombinant Rabbit mAb [SY24-03]	ET1607-34	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
Cell treatment	Phospho-MEK1/2 (S218 + S222) Recombinant Rabbit mAb [ST0490]	ET1609-50	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Phospho-MEK1 (S298) Recombinant Rabbit mAb [SD206-7]	ET1612-40	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	Phospho-MEK1 (T292) Recombinant Rabbit mAb [SD2088]	ET1612-42	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-p38 alpha (T180) Recombinant Rabbit mAb [PSH06-25]	HA722657	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-p38 alpha (T180) Recombinant Rabbit mAb [PSH06-25] - BSA and Azide free	HA751062	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-p38 alpha (T180 + Y182) Recombinant Rabbit mAb [PSH06-35]	HA722669	Human, Mouse, Rat	WB,IF-Cell,FC
Cell treatment	Phospho-p38 alpha (T180 + Y182) Recombinant Rabbit mAb [PSH06-35] - BSA and Azide free	HA751074	Human, Mouse, Rat	WB,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-p38 (T180) Rabbit pAb	ER2001-51	Human, Mouse, Rat	WB
Cell treatment	Phospho-p38 (T180) Recombinant Rabbit mAb [JE77-37]	HA722150	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-p38 (T180 + Y182) Rabbit pAb	ER2001-52	Human, Mouse, Rat	WB
	Phospho-p38 (T180 + Y182) Rabbit pAb	ER1903-01	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
Cell treatment	Phospho-PI3K p85 (Y467) + PI3K p55 (Y199) Recombinant Rabbit mAb [PSH01-39] - BSA and Azide free	HA750717	Human, Mouse	WB,IF-Cell
Cell treatment	Phospho-PI3K p85 (Y467) + PI3K p55 (Y199) Recombinant Rabbit mAb [PSH01-38]	HA721672	Human, Mouse	WB,IF-Cell
Cell treatment	Phospho-PI3K p85 (Y467) + PI3K p55 (Y199) Recombinant Rabbit mAb [PSH01-38] - BSA and Azide free	HA750716	Human, Mouse	WB,IF-Cell
Cell treatment	Phospho-PI3K p85 (Y467) + PI3K p55 (Y199) Recombinant Rabbit mAb [PSH01-39]	HA721673	Human, Mouse	WB,IF-Cell
Cell treatment	Phospho-Raf1 (S259) Recombinant Rabbit mAb [SD85-07] - BSA and Azide free	HA750304	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-Raf1 (S259) Recombinant Rabbit mAb [SD85-07]	ET1612-87	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	Phospho-Raf1 (S43) Recombinant Rabbit mAb [SD087-05]	ET1612-75	Human, Mouse, Rat	WB
	Phospho-Raf1 (S621) Recombinant Rabbit mAb [JJ085-05]	ET1701-3	Human, Rat, Mouse	WB,IHC-P
Cell treatment	Phospho-Smad1/5/9 (S463/S465/S467) Recombinant Rabbit mAb [JE59-46]	HA722566	Human, Mouse, Rat	WB,IHC,IF-Cell
Cell treatment	Phospho-Smad2 (S250) Recombinant Rabbit mAb [SD207-1]	ET1612-32	Human, Mouse, Rat	WB,FC,IP,IHC-P
Cell treatment	Phospho-Smad2 (S250) Recombinant Rabbit mAb [SD207-1] - BSA and Azide free	HA750291	Human, Mouse, Rat	WB,FC,IP,IHC-P
Cell treatment	Phospho-Smad2 (S255) Recombinant Rabbit mAb [JF0882]	ET1702-34	Human, Mouse, Rat	WB,IHC-P,IP
Cell treatment	Phospho-SMAD2 (S465+467) Recombinant Rabbit mAb [PSH05-76] - BSA and Azide free	HA751026	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-SMAD2 (S465+467) Recombinant Rabbit mAb [PSH05-76]	HA722443	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-Smad3 (S423 + S425) Recombinant Rabbit mAb [ST0493]	ET1609-41	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	Phospho-Smad3 (S423 + S425) Recombinant Rabbit mAb [ST0493] - BSA and Azide free	HA750172	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	Phospho-Smad5 (S463 + S465) Recombinant Rabbit mAb [SY09-03] - BSA and Azide free	HA750089	Human, Mouse, Rat	WB,IF-Cell,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-Smad5 (S463 + S465) Recombinant Rabbit mAb [SY09-03]	ET1605-5	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment, KD	Phospho-STAT3 (S727) Recombinant Rabbit mAb [SY24-09]	ET1607-39	Human, Rat, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP
Cell treatment, KD	Phospho-STAT3 (S727) Recombinant Rabbit mAb [SY24-09] - BSA and Azide free	HA750121	Human, Rat, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP
Cell treatment	Phospho-STAT3 (Y705) Recombinant Rabbit mAb [SZ43-01] - BSA and Azide free	HA750067	Human, Mouse, Rat	WB,IHC-P,IP
Cell treatment	Phospho-STAT3 (Y705) Recombinant Rabbit mAb [SZ43-01]	ET1603-40	Human, Mouse, Rat	WB,IHC-P,IP
	PI3 Kinase p110 beta Recombinant Rabbit mAb [JE65-38]	HA722474	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	PI 3 Kinase p110 delta Recombinant Rabbit mAb [JM83-24]	ET1705-46	Human, Mouse, Rat	WB,IHC-P
KD	PI3-kinase p110 subunit beta Recombinant Rabbit mAb [SC0580]	ET1610-36	Human	WB,IP,IHC-P,FC
	PI3 Kinase p110 α Recombinant Rabbit mAb [SJ0186]	ET1606-36	Human, Mouse	WB,IF-Cell,IF-Tissue,IP,IHC-P
	PI 3 Kinase p55 gamma Recombinant Rabbit mAb [SR44-06]	ET1602-27	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	PI 3 Kinase p85 alpha Recombinant Rabbit mAb [SU04-07] - BSA and Azide free	HA750162	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	PI 3 Kinase p85 alpha Recombinant Rabbit mAb [SU04-07]	ET1608-70	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	PI 3 Kinase p85 alpha Recombinant Mouse mAb [A3-D0-R] - BSA and Azide free	HA610084	Human	WB,IF-Cell
	PI 3 Kinase p85 alpha Rabbit pAb	R1312-6	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	PI 3 Kinase p85 alpha Recombinant Mouse mAb [A3-D0-R]	HA601206	Human	WB,IF-Cell
KD	PI 3 Kinase p85 alpha Mouse mAb [L3-D6]	EM1701-62	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue
KD	PI 3 Kinase p85 alpha Mouse mAb [A3-D0]	M1510-2	Human	WB,IF-Cell
	PI3 Kinase p85 alpha Recombinant Rabbit mAb [JE04-11]	HA722522	Human, Mouse, Rat	WB
	PI 3 Kinase p85 alpha/gamma Rabbit pAb	ER64588	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Tissue
	PI 3 Kinase p85 beta Recombinant Rabbit mAb [ST04-77]	ET1609-30	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
KD	Raf1 Recombinant Rabbit mAb [JJ086-06]	ET1701-21	Human, Mouse, Rat	WB
	Ras Recombinant Rabbit mAb [SA39-05]	ET1601-16	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,FC,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Ras Recombinant Rabbit mAb [JF10-11]	ET1702-94	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IP,FC
	Ras Recombinant Rabbit mAb [JE44-13]	HA721883	Human, Mouse, Rat	WB,IF-Cell,FC
	Ras Rabbit pAb	ER40115	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	REST Recombinant Rabbit mAb [JE32-25]	HA720081	Human	WB,IF-Cell,IHC-P
	RNF3 Rabbit pAb	ER1915-99	Human	WB
	SETDB1 Mouse mAb [A12-H6]	EM1709-76	Human, Mouse, Monkey	WB,IF-Cell
KO	Smad1 Recombinant Rabbit mAb [SY0254]	ET1607-42	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P
KO	Smad1 Recombinant Rabbit mAb [SY0254] - BSA and Azide free	HA750124	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P
	Smad1 Recombinant Rabbit mAb [SD2090]	ET1612-39	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Smad1/5/9 Rabbit pAb	ER64980	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Smad2 Mouse mAb [C9-B0]	EM1701-44	Human, Mouse, Rat	WB,IHC-P,IF-Cell
KO	Smad2 Recombinant Rabbit mAb [SP06-05] - BSA and Azide free	HA750073	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC
	Smad2 Rabbit pAb	R1511-29	Human, Mouse, Rat	IF-Cell,IHC-P
KO	Smad2 Recombinant Rabbit mAb [SP06-05]	ET1604-22	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC
	Smad2/3 Mouse mAb [4D3]	RT1566	Human, Mouse, Rat	WB,IP,IF,IHC-P
	Smad2/3 Rabbit pAb	ER64981	Human, Mouse, Rat, Monkey	WB,IHC-P
	Smad2/3 (Acetyl K19) Rabbit pAb	ER60061	Human, Mouse, Rat	IHC-P,WB
KO	Smad3 Recombinant Rabbit mAb [SY25-01] - BSA and Azide free	HA750123	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KO	Smad3 Recombinant Rabbit mAb [SY25-01]	ET1607-41	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Smad3 Rabbit pAb	R1510-24	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC
	Smad3 Mouse mAb [8E1]	RT1567	Human, Mouse, Rat	WB,IP,IF,IHC-P
	Smad3 Rabbit pAb	R1510-25	Human	WB,IF-Cell,IHC-P,FC
	Smad4 Recombinant Rabbit mAb [PD01-37]	HA721234	Human, Mouse, Rat	WB,IHC-P,IF-Cell
KO	Smad4 Recombinant Rabbit mAb [SP05-05] - BSA and Azide free	HA750070	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KO	Smad4 Recombinant Rabbit mAb [SP05-05]	ET1604-12	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Smad4 Rabbit pAb	ER1706-90	Human, Mouse, Rat	IF-Cell,IHC-P,FC
	Smad5 Recombinant Rabbit mAb [SI17-06]	ET1606-26	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	SnoN Rabbit pAb	ER64992	Human, Mouse	WB,IHC-P
	SOX2 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS274RT	Mouse	mIHC
	SOX2 Mouse mAb [C9-F10]	EM1708-84	Human	WB,IF-Cell,IHC-P
	SOX2 Rabbit pAb	R1106-1	Human, Mouse, Rat	WB
	SOX2 Recombinant Antibody [PO00-28] - Chicken IgY (Chimeric)	HA601499	Human, Mouse	IHC-Fr,WB
	SOX2 Recombinant Antibody [PO00-28] - Guinea pig IgG2 (Chimeric)	HA601500	Human, Mouse	IHC-Fr,WB
	SOX2 Recombinant Antibody [PO00-28] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610236	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB
	SOX2 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS273RT	Human	mIHC
	SOX2 Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS281MS	Human	mIHC
	SOX2 Recombinant Antibody [PO00-28] - Chicken IgY (Chimeric) - BSA and Azide free	HA610329	Human, Mouse	IHC-Fr,WB
Relative expression	SOX2 Recombinant Rabbit mAb [PO00-28]	HA721155	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-P,WB,IF-Cell,ChIP,IHC-Fr,IF-Tissue
	iFluor™ 488 Conjugated SOX2 Recombinant Rabbit mAb [PO00-28]	HA720179F	Human, Mouse	IF-Cell,FC
	SOX2 Rabbit pAb	IRS198RB	Human	mIHC
	SOX2 Recombinant Antibody [PO00-28] - Rat IgG1 (Chimeric)	HA601396	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB
	iFluor™ 594 Conjugated SOX2 Recombinant Rabbit mAb [PO00-28]	HA720180F	Human, Mouse	IF-Cell,FC
	SOX2 Recombinant Rabbit mAb	IRS240RB	Mouse	mIHC
	SOX2 Recombinant Antibody [PO00-28] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610330	Human, Mouse	IHC-Fr,WB
	SOX2 Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS282MS	Mouse	mIHC
KD	STAT3 Recombinant Rabbit mAb [SY34-01]	ET1605-45	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	STAT3 Recombinant Rabbit mAb [SY24-08]	ET1607-38	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC
KD	STAT3 Recombinant Rabbit mAb [SY24-08] - BSA and Azide free	HA750120	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	STAT3 Rabbit pAb	RT1591	Human, Mouse, Rat	WB,IP,IF,IHC-P
KD	STAT3 Recombinant Rabbit mAb [SY34-01] - BSA and Azide free	HA750087	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Tbx3 Rabbit pAb	R1306-4	Human, Mouse, Rat	WB
	Tbx3 Recombinant Rabbit mAb [JE33-38]	HA721184	Human	WB,FC
	TCF7 Recombinant Rabbit mAb [PSH10-10]	HA723179	Human	WB,IHC-P,IP
Relative expression	TCF7 Recombinant Rabbit mAb [PSH13-89]	HA723582	Human, Mouse, Rat	WB,IHC-P,IP,mIHC
	TCF7 Recombinant Rabbit mAb [PSH10-10] - BSA and Azide free	HA751340	Human	WB,IHC-P,IP
	TCF7 Rabbit pAb	ER65115	Human, Mouse	WB,IHC-P
Relative expression	TCF7 Recombinant Rabbit mAb [PSH13-89] - BSA and Azide free	HA751491	Human, Mouse, Rat	WB,IHC-P,IP
	TCF7 Recombinant Rabbit mAb	IRS204RB	Human	mIHC
	TGFB4/LEFTY2 Rabbit pAb	ER1917-67	Human	WB
	Wnt-1 Rabbit pAb	ER65317	Human, Mouse	WB,IHC-P,IF-Cell
	Wnt10a Rabbit pAb	HA500356	Human, Mouse	WB
	Wnt11 Rabbit pAb	ER1918-68	Human, Mouse	WB,IHC-P,IF-Tissue
	Wnt16 Rabbit pAb	HA500145	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	WNT2 Rabbit pAb	ER1511-4	Human, Mouse	WB,FC
	Wnt2b Recombinant Rabbit mAb [JB48-34]	ET7107-20	Human, Mouse, Rat	WB,IF-Cell,FC
Cell treatment	Wnt3a Rabbit pAb	HA500193	Human	WB,IF-Cell,IHC-P
	Wnt3a Recombinant Rabbit mAb [PSH19-96]	HA724116	Human, Mouse	WB,IHC-P,IP
	WNT4 Mouse mAb [4G1]	RT1662	Human, Mouse, Rat	WB,IP,IF,IHC-P
	Wnt5a Recombinant Rabbit mAb [JU30-30]	ET1706-33	Human	WB,IHC-P
	WNT5B Recombinant Rabbit mAb [JE53-67]	ET7110-34	Human	WB,FC
	Wnt6 Rabbit pAb	ER1918-70	Human, Mouse	WB,IHC-P,IF-Tissue
	WNT7A Rabbit pAb	ER1918-71	Human, Mouse, Rat	WB
	WNT7B Rabbit pAb	ER60344	Human, Mouse	WB
	WNT8A Recombinant Rabbit mAb [PSH0-39]	HA721341	Human, Mouse, Rat	WB,IF-Cell,FC
	WNT8A Recombinant Rabbit mAb [PSH0-39] - BSA and Azide free	HA750617	Human, Mouse, Rat	WB,IF-Cell,FC
	WNT8B Rabbit pAb	ER61525	Human, Mouse	WB
	WNT9A Rabbit pAb	ER60345	Human, Mouse	WB
	WNT9B Rabbit pAb	HA500412	Human, Mouse	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	ZNF144 Rabbit pAb	ER1918-94	Human, Rat	WB,IHC-P,IF-Tissue
--	-------------------	-----------	------------	--------------------