

常见PTM修饰-抗体清单

高级验证	产品名	货号	反应种属	验证应用
Cell treatment	Acetylated-Lysine Recombinant Rabbit Multiclonal Antibody [PSH10-04]	HA723173	Species independent	WB,IHC-P,ChIP,IF-Tissue
Cell treatment	Acetylated-Lysine Recombinant Rabbit mAb [PSH09-15]	HA723073	Species independent	WB,IHC-P,IF-Tissue,ChIP,IP
Cell treatment	Acetylated-Lysine Recombinant Rabbit mAb [PSH09-16]	HA723074	Species independent	WB,IHC-P,IF-Tissue
	ACSS1 (Acetyl K642) Rabbit pAb	ER60060	Human, Mouse, Rat	WB
	alpha Tubulin (Acetyl K112) Rabbit pAb	ER60054	Human, Mouse, Rat	WB
Cell treatment	alpha Tubulin (acetyl K40) Recombinant Rabbit mAb [PSH0-13]	HA721271	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Asymmetric Di-Methyl Arginine (ADMA) Recombinant Rabbit mAb [PSH09-55]	HA723117	Species independent	ELISA,WB
	E2F-4 (Acetyl K96) Rabbit pAb	ER60051	Human, Mouse, Rat	WB
	eIF5A/eIF5A2 (Acetyl K47) Rabbit pAb	ER60050	Human, Mouse, Rat	WB
	Histone H1 (Acetyl K25) Rabbit pAb	ER60044	Human, Monkey	WB,IHC-P
	Histone H1 (di methyl K25) Rabbit pAb	ER60071	Human, Mouse, Rat	WB
	Histone H1 (tri methyl K25) Rabbit pAb	ER60076	Human	WB
	Histone H2A (Acetyl K15) Rabbit pAb	ER60073	Human, Mouse, Rat	WB
	Histone H2A (Acetyl K5) Rabbit pAb	ER60045	Human, Mouse, Rat	IHC-P,WB
	Histone H2A (Acetyl K5) Rabbit pAb	ER60072	Human, Mouse, Rat	WB
Cell treatment	Histone H2A (acetyl K9) Recombinant Rabbit mAb [SR4-15]	ET1602-34	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,ChIP,CUT&Tag-seq
Cell treatment	Histone H2A (acetyl K9) Recombinant Rabbit mAb [PS01-45]	HA721284	Human, Mouse, Rat	WB,IHC-P,IF-Cell,CUT&Tag-seq,ChIP
	Histone H2A.Z (Acetyl K4) Rabbit pAb	ER60075	Human, Mouse, Rat	WB
	Histone H2A.Z (Acetyl K7) Rabbit pAb	ER60074	Human, Mouse, Rat	WB
	Histone H2B (Acetyl K12) Rabbit pAb	ER60047	Human, Mouse, Monkey	WB,IHC-P,IF-Cell
	Histone H2B (Acetyl K15) Rabbit pAb	ER60048	Human, Mouse	WB,IHC-P,IF-Cell
Cell treatment	Histone H2B (acetyl K20) Recombinant Rabbit mAb [SN72-09]	ET1611-48	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP,ChIP,CUT&Tag-seq
Cell treatment	Histone H2B (acetyl K20) Recombinant Rabbit mAb [PS01-46]	HA721285	Human, Mouse, Rat	WB,IHC-P,IF-Cell,ChIP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Histone H2B (Acetyl K5) Rabbit pAb	ER60046	Human, Mouse	WB,IHC-P,IF-Cell
	Histone H2B (mono methyl R79) Recombinant Rabbit mAb [ST0430]	ET1609-9	Human, Mouse, Rat	IF-Cell,IHC-P
Cell treatment	Histone H3 (acetyl K14) Recombinant Rabbit mAb [JU43-26]	ET1706-28	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,SNAP-ChIP,CUT&Tag-seq
Cell treatment	Histone H3 (acetyl K18) Mouse mAb [A7B6]	HA600091	Human	WB,IF-Cell
Cell treatment	Histone H3 (acetyl K18) Mouse mAb [A7B5]	HA600090	Human, Mouse, Rat	WB,IHC-P,ChIP
Cell treatment	Histone H3 (acetyl K18) Rabbit pAb	HA500047	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Cell treatment	Histone H3 (acetyl K23) Recombinant Rabbit mAb [JE46-39]	HA721840	Human, Mouse	WB,IF-Cell
Cell treatment	Histone H3 (acetyl K27) Mouse mAb [A6D7]	HA600047	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,ChIP
Cell treatment	Histone H3 (acetyl K27) Rabbit pAb	HA500046	Human, Mouse, Rat	WB,IHC-P,ChIP
Cell treatment	Histone H3 (acetyl K27) Mouse mAb [A6D6]	HA600048	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,ChIP
Cell treatment	Histone H3 (acetyl K27) Recombinant Rabbit mAb [PSH13-93]	HA723586	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Cell treatment	Histone H3 (acetyl K36) Recombinant Rabbit mAb [JE77-32]	HA722415	Human, Mouse, Rat	WB,IHC-P,IF-Cell,ChIP
Cell treatment	Histone H3 (acetyl K4) Recombinant Rabbit mAb [PSH03-71]	HA722035	Human, Mouse, Rat	WB,IF-Cell,ChIP,Dot Blot
Cell treatment	Histone H3 (acetyl K56) Recombinant Rabbit mAb [SU30-10]	ET1608-9	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,ChIP,CUT&Tag-seq
Cell treatment	Histone H3 (acetyl K9) Recombinant Rabbit mAb [PSH04-47] - ChIP Grade	HA722132	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,FC,ChIP,Dot Blot,IP
Cell treatment	Histone H3 (di methyl K27) Recombinant Rabbit mAb [PSH05-88]	HA722486	Human, Mouse, Rat	WB,IF-Cell,ChIP
Cell treatment	Histone H3 (di methyl K4) Rabbit pAb	R1110-3	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,Dot Blot
	Histone H3 (di methyl K4) Recombinant Rabbit mAb [JE00-98]	HA722135	Human, Mouse, Rat, Monkey	WB,IF-Cell,ChIP,Dot Blot,IHC-P
	Histone H3 (di methyl K9) Recombinant Rabbit mAb [SN07-30]	ET1611-51	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	Histone H3 (mono+di+tri methyl K36) Recombinant Rabbit mAb [PSH07-14]	HA722802	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,ChIP
	Histone H3 (mono+di+tri methyl K4) Recombinant Rabbit mAb [PSH06-30]	HA722662	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,ChIP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Histone H3 (mono+di+tri methyl K79) Recombinant Rabbit mAb [SR42-06]	ET1602-41	Human, Mouse, Rat	WB,IHC-P,ChIP
	Histone H3 (mono+di+tri methyl K79) Recombinant Rabbit mAb [PSH09-10]	HA723068	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,ChIP
	Histone H3 (mono methyl K14) Recombinant Rabbit mAb [JE43-29]	HA722232	Human, Mouse, Rat	WB,IF-Cell,ChIP,Dot Blot
	Histone H3 (mono methyl K18) Recombinant Rabbit mAb [SA42-07]	ET1601-14	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC-P
	Histone H3 (mono methyl K27) Recombinant Rabbit mAb [PSH05-67]	HA722376	Human, Mouse, Rat	WB,IF-Cell,FC,ChIP
	Histone H3 (mono methyl K36) Recombinant Rabbit mAb [SR04-20]	ET1602-48	Human, Mouse	WB,IF-Cell,IF- Tissue,ChIP
	Histone H3 (mono methyl K4) Rabbit pAb	R1110-1	Human, Mouse, Rat	WB,IF-Cell,IHC- P,FC,Dot Blot
	Histone H3 (mono methyl K4) Recombinant Rabbit mAb [PSH07-15]	HA722803	Human, Mouse, Rat	WB,IF-Cell,ChIP,Dot Blot
KD	Histone H3 (mono methyl K79) Recombinant Rabbit mAb [PSH07-66]	HA722848	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC- P,FC,ChIP
	Histone H3 (mono methyl K9) Recombinant Rabbit mAb [JE43-27]	HA722350	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC- P,FC,ChIP,IF-Tissue,Dot Blot
	Histone H3 (mono methyl R2) Recombinant Rabbit mAb [ST0427]	ET1609-8	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC-P,ChIP
Cell treatment	Histone H3 (tri methyl K27) Recombinant Rabbit mAb [PSH05-15]	HA722231	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF- Tissue,ChIP
	Histone H3 (tri methyl K36) Recombinant Rabbit mAb [PSH15-01]	HA723701	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,IP,Dot Blot,ChIP
	Histone H3 (tri methyl K9) Mouse mAb [1-6]	M1112-3	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF- Tissue
	Histone H3 (tri methyl K9) Mouse mAb [2G1]	RT1287	Human	WB,IF
Cell treatment	Histone H4 (acetyl K16) Recombinant Rabbit mAb [JB21-44]	ET7107-89	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell,IF- Tissue,ChIP
Cell treatment	Histone H4 (acetyl K5) Recombinant Rabbit mAb [SR31-07]	ET1602-40	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC- P,IP,ChIP,CUT&Tag-seq
	Histone H4 (acetyl K8) Recombinant Rabbit mAb [PS01-25]	HA721249	Human, Mouse, Rat	WB,IF-Cell,IHC-P,ChIP
	Histone H4 (tri methyl K20) Recombinant Rabbit mAb [ST044-07]	ET1609-16	Human, Mouse, Rat	IHC-P
	HSP70 (Acetyl K246) Rabbit pAb	ER60056	Human	WB
	L-Lactyl-alpha Tubulin (K40) Rabbit pAb	HA500554	Human, Mouse, Rat, Zebrafish	WB,IHC-P
Cell treatment	Mono-Methyl Arginine (R*GG) Recombinant Rabbit mAb [PSH07-52]	HA722840	Species independent	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	NF-E4 (Acetyl K43) Rabbit pAb	ER60055	Human	WB
Cell treatment	p53 (acetyl K370) Recombinant Rabbit mAb [SR40-09]	ET1602-38	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,FC
Cell treatment	p53 (acetyl K382) Recombinant Rabbit mAb [SD0801]	ET1612-71	Human, Mouse	WB,IF-Cell,IHC-P
Cell treatment	p53 (acetyl K382) Recombinant Rabbit mAb [PSH06-87]	HA722762	Human, Mouse	WB,IF-Cell
Cell treatment	p53 (acetyl K382) Recombinant Rabbit mAb [PSH06-88]	HA722763	Human	WB,IF-Cell
Cell treatment	Pan-Lactyl-lysine Recombinant Rabbit mAb [PSH03-74]	HA722038	Species independent	WB,IHC-P,IF-Cell,Dot Blot
Cell treatment	Pan-Lactyl-lysine Recombinant Rabbit mAb [PSH03-73]	HA722037	Species independent	WB,IHC-P,IF-Cell,Dot Blot
	Phospho-4E-BP1/2/3 (T45) Mouse mAb [3G2]	RT1004	Human, Mouse, Rat	WB,IP,IF
Cell treatment	Phospho-Acetyl Coenzyme A Carboxylase (S79) Recombinant Rabbit mAb [PSH07-13]	HA722801	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-Acetyl Coenzyme A Carboxylase (S79) Recombinant Rabbit mAb [JE63-95]	HA721714	Human, Mouse, Rat	WB
Cell treatment	Phospho-Acetyl Coenzyme A Carboxylase (S79) Recombinant Rabbit mAb [JE55-54]	HA721516	Human, Mouse, Rat	WB
Cell treatment	Phospho-AKT1 (S124) Recombinant Rabbit mAb [JJ08-46]	ET1701-36	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
Cell treatment	Phospho-AKT1 (S129) Recombinant Rabbit mAb [JE59-50]	HA721677	Human, Mouse, Rat	WB
Cell treatment	Phospho-AKT1 (T450) Recombinant Rabbit mAb [SD08-12]	ET1612-73	Human, Mouse, Rat	WB,IHC-P,IP,IF-Cell,FC
Cell treatment	Phospho-AKT2 (S474) Recombinant Rabbit mAb [PSH09-46]	HA723108	Human, Mouse, Rat	WB,FC
Cell treatment	Phospho-AKT2 (S474) Rabbit pAb	HA500116	Human, Mouse	WB,IHC-P
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [PSH04-44]	HA722129	Human, Mouse, Rat	WB,IHC-P,FC
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [SY28-05]	ET1607-73	Human, Mouse, Rat, Dog	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
Cell treatment	Phospho-AKT (T308) Recombinant Rabbit mAb [PSH08-10]	HA722951	Human	WB,IF-Cell,FC
	Phospho-alpha Synuclein (S129) Recombinant Rabbit mAb [JB22-44]	ET7107-30	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Phospho-AMPK alpha 1 (S485) Recombinant Rabbit mAb [PSH18-01]	HA723966	Human, Mouse, Rat	WB,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-AMPK alpha 1 (S496) Recombinant Rabbit mAb [SD0810]	ET1612-72	Human, Mouse	WB,IP,IHC-Fr
Cell treatment	Phospho-AMPK alpha 1 (T183) +AMPK alpha 2 (T172) Recombinant Rabbit mAb [PSH14-13]	HA723603	Human, Mouse, Rat	WB,IHC-P
	Phospho-AMPK alpha 2 (S345) Recombinant Rabbit mAb [JJ08-19]	ET1701-37	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,FC
Cell treatment	Phospho-AMPK alpha 2 (S491) Recombinant Rabbit mAb [JE59-49]	HA722567	Human	WB
	Phospho-AP2M1 (Thr156) Rabbit pAb	ER65961	Human, Mouse, Rat	WB
Cell treatment	Phospho-AS160 (T642) Recombinant Rabbit mAb [JE60-97]	HA721489	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-ATF2 (T71) Recombinant Rabbit mAb [SC05-90]	ET1610-30	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IP,IHC-P
Cell treatment	Phospho-ATM (S1981) Recombinant Rabbit mAb [JM93-25]	ET1705-50	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Cell treatment	Phospho-ATR (T1989) Rabbit pAb	HA500405	Human	WB
Cell treatment	Phospho-ATR (T1989) Recombinant Rabbit mAb [PS01-18]	HA721190	Human	WB,IF-Cell
Cell treatment	Phospho-Axl (Y702) Recombinant Rabbit mAb [PSH17-98]	HA723963	Human	WB
Cell treatment	Phospho-BCAR1 (Y410) Recombinant Rabbit mAb [JE47-97]	HA722124	Human	WB
Cell treatment	Phospho-BCKDHA (S293) Recombinant Rabbit mAb [PSH14-02]	HA723592	Human, Mouse, Rat, Monkey	WB
Cell treatment	Phospho-Bcl-2 (S70) Recombinant Rabbit mAb [PSH11-02]	HA723269	Human	WB,IF-Cell,FC
Cell treatment	Phospho-Bcln 1 (S30) Recombinant Rabbit mAb [PSH16-40]	HA723840	Human, Mouse, Rat	WB
Cell treatment	Phospho-Beta Catenin (S33 + S37) Recombinant Rabbit mAb [JE59-59]	HA721580	Human, Rat, Mouse	WB
Cell treatment	Phospho-Beta Catenin (S552) Recombinant Rabbit mAb [PSH08-72]	HA723024	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Cell treatment	Phospho-Beta Catenin (T41/S45) Recombinant Rabbit mAb [JE54-02]	HA722316	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-BRAF (T401) Recombinant Rabbit mAb [JJ08-72]	ET1701-20	Human, Mouse, Rat	WB,IF-Cell,IP,IHC-P
Cell treatment	Phospho-BTK (Y223) Recombinant Rabbit mAb [ST51-05]	ET1609-55	Human, Mouse	WB,IP
Cell treatment	Phospho-CaMKII (T286) Recombinant Rabbit mAb [JE51-81]	HA721794	Human, Mouse, Rat	WB
	Phospho-CaMKII (T286) Mouse mAb [2G1]	RT1092	Human, Mouse, Rat	WB,IP,IF

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-CBL (Y774) Recombinant Rabbit mAb [JE62-93]	HA722786	Human	WB
Cell treatment	Phospho-CDC37 (S13) Recombinant Rabbit mAb [SD08-50]	ET1612-74	Human, Mouse, Rat	WB,IP
Cell treatment	Phospho-Cdc6 (S54) Recombinant Rabbit mAb [SD08-47]	ET1612-96	Human, Mouse, Rat	WB,IF-Tissue,IHC-P
Cell treatment	Phospho-Cdk1/2/3 (T14) Recombinant Rabbit mAb [JE04-37]	HA722272	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-CDK1 (T161) Recombinant Rabbit mAb [PSH03-42]	HA721987	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Cell treatment	Phospho-CDK1 (T161) + CDK2/CDK3 (T160) Recombinant Rabbit mAb [JE75-37]	HA722582	Human	WB,Dot Blot
Cell treatment	Phospho-Cdk2 (Y15) Recombinant Rabbit mAb [SN72-04]	ET1611-52	Human, Mouse, Rat	WB,IHC-P,IP,IF-Tissue
Cell treatment	Phospho-CDK (Y15) Recombinant Rabbit mAb [PSH09-45]	HA723107	Human, Mouse	WB,IHC-P
Cell treatment	Phospho-Chk1 (S280) Recombinant Rabbit mAb [JE47-66]	HA721471	Human, Mouse, Rat	WB
Cell treatment	Phospho-Chk1 (S296) Recombinant Rabbit mAb [SN06-50]	ET1611-76	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Cell treatment	Phospho-Chk1 (S317) Recombinant Rabbit mAb [PSH03-84]	HA722050	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-Chk1 (S345) Recombinant Rabbit mAb [PS01-17]	HA721189	Human, Mouse	WB,IF-Cell
Cell treatment	Phospho-Chk2 (T68) Recombinant Rabbit mAb [JE43-12]	HA721633	Human, Mouse, Rat	WB
Cell treatment	Phospho-c-Jun (S63) Recombinant Rabbit mAb [SY0297]	ET1608-4	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-c-Jun (S73)+JunD (S100) Recombinant Rabbit mAb [JE65-46]	HA722475	Human, Mouse, Rat	WB,IF-Cell,FC
Cell treatment	Phospho-c-Jun (S73)+JunD (S100) Recombinant Rabbit mAb [PSH07-12]	HA722800	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-c-Jun (T91) Recombinant Rabbit mAb [JJ080-9]	ET1701-32	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
Cell treatment	Phospho-c-Myc (S62) Recombinant Rabbit mAb [ST49-08]	ET1609-64	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
Cell treatment	Phospho-c-Myc (T58) Recombinant Rabbit mAb [SN60-01]	ET1611-24	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,FC,IHC-P
Cell treatment	Phospho-c-Myc (T58 + S62) Recombinant Rabbit mAb [SZ02-06]	ET1603-24	Human, Mouse	WB,IP,FC
Cell treatment	Phospho-Cofilin (S3) Recombinant Rabbit mAb [PSH07-18]	HA722806	Human, Rat, Monkey	WB
Cell treatment	Phospho-Creb (S133) Recombinant Rabbit mAb [JB25-40]	ET7107-93	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-CrkL (Y207) Recombinant Rabbit mAb [JE55-58]	HA721675	Human, Mouse, Rat	WB
Cell treatment	Phospho-CRMP2 (T514) Recombinant Rabbit mAb [JE57-25]	HA721674	Human, Mouse, Rat	WB
	Phospho-Cyclin E1 (T77) Recombinant Rabbit mAb [SD2025]	ET1612-31	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	Phospho-Cytosolic Phospholipase A2 (S505) Recombinant Rabbit mAb [PSH07-56]	HA722849	Human, Mouse	WB,Dot Blot
Cell treatment	Phospho-Cytosolic Phospholipase A2 (S505) Recombinant Rabbit mAb [PSH07-16]	HA722804	Human, Mouse	WB
Cell treatment	Phospho-DNA PKcs (S2056) Recombinant Rabbit mAb [JE41-92]	HA722066	Human	WB,IF-Cell
Cell treatment	Phospho-DRP1 (S616) Recombinant Rabbit mAb [PSH06-77]	HA722760	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell,FC
Cell treatment	Phospho-EEF2 (T56) Recombinant Rabbit mAb [PSH06-36]	HA722670	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Phospho-EGFR (S695) Recombinant Rabbit mAb [SN07-36]	ET1611-50	Human	WB,IF-Cell,IF-Tissue
Cell treatment	Phospho-EGFR (T669) Recombinant Rabbit mAb [PSH10-03]	HA723172	Human, Mouse, Rat	WB,IP
Cell treatment	Phospho-EGFR (Y1045) Recombinant Rabbit mAb [PSH03-28]	HA721970	Human	WB,IF-Cell,FC
Cell treatment	Phospho-EGFR (Y1068) Recombinant Rabbit mAb [SD2055]	ET1612-30	Human	WB,IF-Cell,IF-Tissue,IP,FC
Cell treatment	Phospho-EGFR (Y1092) Recombinant Rabbit mAb [SJ0194]	ET1606-44	Human, Mouse	WB,IF-Tissue,IHC-P
	Phospho-EGFR (Y1148) Rabbit pAb	HA500580	Human	WB,IF-Cell
Cell treatment	Phospho-EGFR (Y1173) Recombinant Rabbit mAb [SC57-04]	ET1610-4	Human, Mouse	WB,IF-Cell,IHC-P,IP
Cell treatment	Phospho-EGFR (Y845) Recombinant Rabbit mAb [JE46-77]	HA721488	Human	WB
Cell treatment, Relative expression	Phospho-EIF2S1 (S51) Recombinant Rabbit mAb [SZ01-06]	ET1603-14	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Phospho-EIF2S1 (S52) Rabbit pAb	RT1198	Human, Mouse, Rat	WB,IP,IF,IHC-P
Cell treatment	Phospho-eIF4B (S406) Recombinant Rabbit mAb [JE63-80]	HA722149	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Cell treatment	Phospho-eIF4B (S406) Recombinant Rabbit mAb [JE58-32]	HA721676	Human, Mouse, Rat	WB
Cell treatment	Phospho-eIF4EBP1 (S65) Recombinant Rabbit mAb [PSH09-44]	HA723106	Human	WB,IF-Cell,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-eIF4EBP1 (T70) Recombinant Rabbit mAb [PSH09-91]	HA723222	Human, Mouse	WB,IHC-P,IP
Cell treatment	Phospho-eIF4E (S209) Recombinant Rabbit mAb [SU0396]	ET1608-66	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC
Cell treatment	Phospho-eIF4G (S1108) Recombinant Rabbit mAb [PS01-08]	HA721276	Human	WB,IHC-P
	Phospho-ELK1 (S383) Rabbit pAb	ER65967	Human, Mouse, Rat	WB,IHC
Cell treatment	Phospho-ELK1 (S383) Recombinant Rabbit mAb [PSH14-65]	HA723663	Human, Mouse	WB,IHC-P,IP
Cell treatment	Phospho-eNOS (S1177) Recombinant Rabbit mAb [PS01-10]	HA721188	Human, Dog, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-ErbB3/HER3 (Y1289) Recombinant Rabbit mAb [JE51-25]	HA722178	Human	WB,IF-Cell
Cell treatment	Phospho-Erk1 (T202)+Erk2 (T185) Recombinant Rabbit mAb [SZ2-4]	ET1603-22	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
Cell treatment	Phospho-Erk1 (T202 + Y204) + Erk2 (T185 + Y187) Recombinant Rabbit mAb [SC58-01]	ET1610-13	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP
Cell treatment	Phospho-Estrogen Receptor alpha (S118) Recombinant Rabbit mAb [SC05-87]	ET1610-32	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	Phospho-Ezrin (T567)/Radixin (T564)/Moesin (T558) Recombinant Rabbit mAb [PSH05-73]	HA722440	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Phospho-FAK (S722) Mouse mAb [4G1]	RT1216	Human, Mouse, Rat	WB,IP,IF,IHC-P
Cell treatment	Phospho-FAK (Y397) Recombinant Rabbit mAb [SC54-07]	ET1610-34	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	Phospho-FOXO1A (T24)+ FOXO3A (T32)+ FOXO4 (T28)+ FOXO6 (T26) Recombinant Rabbit mAb [PSH06-00]	HA722499	Human, Mouse	WB,IF-Cell,FC
Cell treatment	Phospho-FOXO3a (S253) Recombinant Rabbit mAb [ST49-01]	ET1609-49	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-GATA3 (S308) Recombinant Rabbit mAb [ST44-09]	ET1609-17	Human, Mouse, Rat	WB,IHC-P,IP,IF-Tissue
Cell treatment	Phospho-GCN2 (T899) Recombinant Rabbit mAb [JE00-51]	HA722249	Human, Mouse, Rat	WB
	Phospho-Glucose Transporter GLUT4 (S488) Recombinant Rabbit mAb [JE43-23]	HA721866	Mouse, Human, Rat	WB
	Phospho-GluR1 (AMPA subtype) (S831) Recombinant Rabbit mAb [JE44-96]	HA721865	Human, Mouse, Rat	WB
Cell treatment	Phospho-GluR1 (S845) Recombinant Rabbit mAb [JJ2009]	ET1701-28	Mouse, Rat, Human	WB,FC
Cell treatment	Phospho-Glycogen synthase 1 (S641) Recombinant Rabbit mAb [SR46-06]	ET1602-13	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-GSK3 (alpha + beta) (Y216 + Y279) Recombinant Rabbit mAb [SY26-05]	ET1607-54	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,IF-Tissue,FC
Cell treatment	Phospho-GSK3 alpha (S21) Recombinant Rabbit mAb [PSH09-11]	HA723069	Human, Mouse	WB,IHC-P,Dot Blot
Cell treatment	Phospho-GSK3 beta (S9) Recombinant Rabbit mAb [SY02-71]	ET1607-60	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Phospho-HER2 / ErbB2 (T686) Mouse mAb [4G3]	RT1275	Human, Mouse, Rat	WB,IP
Cell treatment	Phospho-HER2 / ErbB2 (Y1139) Recombinant Rabbit mAb [JE60-90]	HA722201	Human	WB,IF-Cell,FC
Cell treatment	Phospho-HER2 / ErbB2 (Y1221 + Y1222) Recombinant Rabbit mAb [JE44-12]	HA721433	Human	WB
Cell treatment	Phospho-HER2 / ErbB2 (Y1248) Recombinant Rabbit mAb [PSH04-03]	HA722084	Human	WB,IF-Cell
Cell treatment	Phospho-HER3/ErbB3 (Y1222) Recombinant Rabbit mAb [JE51-89]	ET7110-10	Human	WB,IHC-P
Cell treatment	Phospho-Histone H1.3 (T17)+Histone H1.4 (T17) Recombinant Rabbit mAb [SR38-03]	ET1602-11	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	Phospho-Histone H2A.X (S139) Recombinant Rabbit mAb [SR33-09]	ET1602-2	Human, Mouse, Rat, Cynomolgus monkey	WB,IHC-P,IF-Cell,IF-Tissue
Cell treatment	Phospho-Histone H3.3 (S31) Recombinant Rabbit mAb [JE62-92]	HA721034	Human, Rat	WB,IF-Cell
Cell treatment	Phospho-Histone H3 (S10) Recombinant Rabbit mAb [SA31-01]	ET1601-30	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC,ChIP
Cell treatment	Phospho-Histone H3 (S28) Recombinant Rabbit mAb [JE56-06]	HA721795	Human, Rat, Mouse	WB,IF-Cell,IHC-P,FC,ChIP
Cell treatment	Phospho-Histone H3 (T3) Recombinant Rabbit mAb [JE42-48]	HA721796	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-Hormone sensitive lipase (S660) Recombinant Rabbit mAb [PSH10-28]	HA723195	Mouse, Rat, Human	WB
Cell treatment	Phospho-Hormone sensitive lipase (S853) Recombinant Rabbit mAb [SN06-39]	ET1611-19	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-HSF1 (S326) Recombinant Rabbit mAb [SU31-03]	ET1608-11	Human	WB,IHC-P,IP
Cell treatment	Phospho-Hsp27 (S78) Recombinant Rabbit mAb [JJ08-70]	ET1701-19	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Phospho-Hsp27 (S82) Recombinant Rabbit mAb [SN06-73]	ET1611-16	Human, Mouse, Rat	WB,IP,IF-Cell,FC
Cell treatment	Phospho-IKB alpha (S32) Recombinant Rabbit mAb [ST53-05]	ET1609-78	Human, Mouse	WB,IF-Cell,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-IKB alpha (S32/S36) Recombinant Rabbit mAb [PSH06-76]	HA722770	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Cell treatment	Phospho-IKB alpha (S36) Recombinant Rabbit mAb [JE56-65]	HA721802	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-IRE1 (S724) Recombinant Rabbit mAb [PSH03-35]	HA721980	Human, Mouse	WB,IF-Cell,IHC-P,FC
Cell treatment	Phospho-IRF3 (S386) Recombinant Rabbit mAb [SU03-28]	ET1608-22	Human	WB,IF-Cell,IF-Tissue
Cell treatment	Phospho-IRF3 (S396) Recombinant Rabbit mAb [PSH06-75]	HA722772	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-JAK2 (Y1007 + Y1008) Recombinant Rabbit mAb [SY24-03]	ET1607-34	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
Cell treatment	Phospho-JNK1/2/3 (T183 + T183 + T221) Recombinant Rabbit mAb [ST500]	ET1609-42	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
Cell treatment	Phospho-JunD (S255) Recombinant Rabbit mAb [JJ08-21]	ET1701-35	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	Phospho-KAP1 (S824) Recombinant Rabbit mAb [JE50-99]	ET7110-11	Human, Mouse, Rat	WB,IHC-P,IP,Dot Blot,IF-Cell
Cell treatment	Phospho-Lamin A + Lamin C (S22) Recombinant Rabbit mAb [PSH11-47]	HA723339	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-Lamin A + Lamin C (S22) Recombinant Rabbit mAb [PSH11-48]	HA723340	Human, Mouse, Rat	WB
	Phospho-LATS1/2 (Ser909) Rabbit pAb	ER65962	Human, Mouse, Rat	WB,IHC
Cell treatment	Phospho-LRRK2 (S935) Recombinant Rabbit mAb [JE42-11]	HA722315	Mouse, Human	WB,IF-Cell,FC
Cell treatment, Relative expression	Phospho-LYN (Y397)/LCK (Y394)/HCK (Y411)/BLK (Y389) Recombinant Rabbit mAb [PSH09-43]	HA723105	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-MAP2 (S136) Recombinant Rabbit mAb [JE44-25]	HA722351	Mouse, Rat	WB,IHC-P,IHC-Fr,IF-Tissue
	Phospho-MCM2 (S108) Recombinant Rabbit mAb [JE63-52]	HA721021	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-MCM2 (S27) Recombinant Rabbit mAb [JE63-40]	HA722574	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
Cell treatment	Phospho-MCM2 (S40) Recombinant Rabbit mAb [JE35-01]	HA722216	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
Cell treatment	Phospho-MCM2 (S41) Recombinant Rabbit mAb [JE50-04]	HA722279	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
Cell treatment	Phospho-MDM2 (S166) Recombinant Rabbit mAb [JE64-20]	HA721791	Human, Mouse, Rat	WB
Cell treatment	Phospho-MEK1/2 (S218 + S222) Recombinant Rabbit mAb [ST0490]	ET1609-50	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Phospho-MEK1 (S298) Recombinant Rabbit mAb [SD206-7]	ET1612-40	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	Phospho-MEK1 (T292) Recombinant Rabbit mAb [SD2088]	ET1612-42	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-MEK3 (S189)+MEK6 (S207) Recombinant Rabbit mAb [PSH11-17]	HA723299	Human, Mouse	WB
Cell treatment	Phospho-Met (Y1234/Y1235) Recombinant Rabbit mAb [PSH06-90]	HA722765	Human	WB
Cell treatment	Phospho-Met (Y1349) Recombinant Rabbit mAb [JE51-68]	ET7110-14	Human, Rat	WB,IHC-P
Cell treatment	Phospho-MLKL (S345) Recombinant Rabbit mAb [JM92-37]	ET1705-51	Mouse	WB,IHC-P,IHC-Fr
Cell treatment	Phospho-MNK1 (T250 + T255) Recombinant Rabbit mAb [PSH17-95]	HA723960	Human, Mouse	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-MSK1 (S376) Recombinant Rabbit mAb [SN63-05]	ET1611-17	Human, Rat	WB,IF-Cell,IHC-P,IP
Cell treatment	Phospho-MST1/MST2 (T180 + T183) Recombinant Rabbit mAb [JE46-64]	HA721737	Human, Mouse, Rat	WB
Cell treatment	Phospho-mTOR (S2448) Mouse mAb [A5D5]	HA600094	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-mTOR (S2481) Recombinant Rabbit mAb [JE59-62]	HA721632	Human	WB
	Phospho-MYL9 (T18 + S19) Recombinant Rabbit mAb [PSH12-71]	HA723482	Mouse, Rat	WB
	Phospho-Myosin light chain (S20) Rabbit pAb	ER1913-46	Human, Mouse	WB,IHC-P,IF-Tissue
Cell treatment	Phospho-MYPT1 (T853) Recombinant Rabbit mAb [PSH08-71]	HA723023	Human, Mouse, Rat	WB
Cell treatment	Phospho-NDRG1 (T346) Recombinant Rabbit mAb [PSH07-17]	HA722805	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IF-Tissue
Cell treatment	Phospho-NF-kB p65 (S529) Recombinant Rabbit mAb [SP07-00]	ET1604-27	Human	WB,IHC-P,FC,IP
Cell treatment	Phospho-NF-kB p65 (S536) Recombinant Rabbit mAb [PSH10-58]	HA723223	Human, Mouse, Rat	WB,IF-Cell
	Phospho-NLRC4 (S533) Rat mAb [B7-B7]	EM1707-24	Human, Mouse	WB,IF-Cell
Cell treatment, Relative expression	Phospho-Nrf2 (S40) Recombinant Rabbit mAb [SU0334]	ET1608-28	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
Cell treatment	Phospho-Nucleolin (T76) Recombinant Rabbit mAb [JE64-03]	HA722942	Human	WB,IHC-P
Cell treatment	Phospho-Nucleolin (T84) Recombinant Rabbit mAb [JE64-18]	HA722920	Human	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Phospho-Nucleophosmin (S125) Recombinant Rabbit mAb [JE51-14]	HA722559	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-Nucleophosmin (T199) Recombinant Rabbit mAb [JE63-16]	HA722125	Human	WB
	Phospho-p130 (T986) Recombinant Rabbit mAb [JE64-95]	HA722472	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-p38 alpha (T180) Recombinant Rabbit mAb [PSH06-25]	HA722657	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-p38 alpha (T180 + Y182) Recombinant Rabbit mAb [PSH06-35]	HA722669	Human, Mouse, Rat	WB,IF-Cell,FC
Cell treatment	Phospho-p38 (T180) Rabbit pAb	ER2001-51	Human, Mouse, Rat	WB
Cell treatment	Phospho-p38 (T180) Recombinant Rabbit mAb [JE77-37]	HA722150	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-p38 (T180 + Y182) Rabbit pAb	ER2001-52	Human, Mouse, Rat	WB
	Phospho-p38 (T180 + Y182) Rabbit pAb	ER1903-01	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
Cell treatment	Phospho-p53 (S15) Recombinant Rabbit mAb [PSH01-98]	HA721756	Human	WB,IF-Cell
Cell treatment	Phospho-p53 (S20) Recombinant Rabbit mAb [PSH06-86]	HA722761	Human	WB
Cell treatment	Phospho-p53 (S33) Recombinant Rabbit mAb [JE41-97]	HA722703	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-p53 (S376) Recombinant Rabbit mAb [ST0440]	ET1609-14	Human, Mouse	WB
Cell treatment	Phospho-p53 (S392) Recombinant Rabbit mAb [SI17-04]	ET1606-24	Human, Mouse, Rat	WB,IHC-P,IP
Cell treatment	Phospho-p53 (T55) Recombinant Rabbit mAb [ST0439]	ET1609-13	Human	WB,IHC-P
	Phospho-p70 S6 Kinase (S434) Mouse mAb [2G1]	RT1456	Human, Mouse, Rat	WB,IP,IF,IHC-P
Cell treatment	Phospho-p70 S6 Kinase (T389) Recombinant Rabbit mAb [PSH02-23]	HA721803	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-p95/NBS1 (S343) Recombinant Rabbit mAb [SY0215]	ET1607-5	Human, Mouse	WB,IF-Tissue,IHC-P,IP
	Phospho-PAK1 (S144)+PAK2 (S141)+PAK3 (S139) Recombinant Rabbit mAb [SU03-70]	ET1608-58	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
Cell treatment	Phospho-PAK4+PAK5+PAK6 (S474 + S560 + S602) Recombinant Rabbit mAb [JE57-42]	HA722146	Human	WB
Cell treatment	Phospho-Paxillin (Y118) Recombinant Rabbit mAb [JE09-10]	HA722526	Human	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-PDGFR beta (Y1009) Recombinant Rabbit mAb [JE46-72]	HA721736	Human, Mouse, Rat	WB
Cell treatment	Phospho-PDGFR beta (Y1021) Recombinant Rabbit mAb [JE48-54]	HA721833	Mouse, Rat, Human	WB
Cell treatment	Phospho-PDGFR beta (Y740) Recombinant Rabbit mAb [SN0646]	ET1611-29	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue
Cell treatment	Phospho-PDGFR beta (Y751) Recombinant Rabbit mAb [PSH13-61]	HA723550	Mouse, Rat, Human	WB
Cell treatment	Phospho-PDHA1 (S293) Recombinant Rabbit mAb [JE58-27]	HA721443	Human, Mouse, Rat	WB,ICC,IHC-P
Cell treatment	Phospho-Phospholamban (S16) Recombinant Rabbit mAb [JE47-45]	HA722291	Mouse, Rat	WB
Cell treatment	Phospho-PI3K p85 (Y467) + PI3K p55 (Y199) Recombinant Rabbit mAb [PSH01-39]	HA721673	Human, Mouse	WB,IF-Cell
Cell treatment	Phospho-PI3K p85 (Y467) + PI3K p55 (Y199) Recombinant Rabbit mAb [PSH01-38]	HA721672	Human, Mouse	WB,IF-Cell
Cell treatment	Phospho-PKA alpha/beta/gamma (catalytic subunit) (T197) Recombinant Rabbit mAb [JE45-17]	HA721864	Human, Mouse, Rat	WB
Cell treatment	Phospho-PKA R2 (S99) Recombinant Rabbit mAb [SC54-04]	ET1610-29	Human, Mouse, Rat, Pig	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Cell treatment	Phospho-PKC alpha/beta II (T638/T641) Recombinant Rabbit mAb [PSH15-91]	HA723796	Human, Mouse, Rat	WB,IF-Cell,Dot Blot
Cell treatment	Phospho-PKC alpha/beta II (T638/T641) Recombinant Rabbit mAb [PSH16-27]	HA723825	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,Dot Blot
Cell treatment	Phospho-PKC alpha (T497) Recombinant Rabbit mAb [JE44-81]	HA722278	Human, Mouse, Rat	WB,IF-Cell,FC
Cell treatment	Phospho-PKC alpha (T638) Recombinant Rabbit mAb [JF0964]	ET1702-17	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,Dot Blot
Cell treatment	Phospho-PKC delta (Y311) Recombinant Rabbit mAb [JE44-20]	HA722458	Human, Mouse, Rat	WB
Cell treatment	Phospho-PKC (S729) Recombinant Rabbit mAb [JE77-54]	HA722846	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-PKC zeta (T560) Recombinant Rabbit mAb [SN0753]	ET1611-78	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-PKR (T446) Recombinant Rabbit mAb [SY230]	ET1607-20	Human	WB,IHC-P,IP,IF-Cell
Cell treatment	Phospho-PKR (T451) Recombinant Rabbit mAb [JE46-70]	HA722067	Human	WB,FC,IF-Cell
Cell treatment	Phospho-PLC-gamma-1 (S1248) Recombinant Rabbit mAb [PSH13-62]	HA723551	Human	WB,IP
Cell treatment	Phospho-PLC-gamma-1 (Y783) Recombinant Rabbit mAb [PSH06-29]	HA722661	Human, Mouse, Rat	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Phospho-PLK1 (T210) Recombinant Rabbit mAb [JJ20-58]	ET1701-33	Human, Mouse	WB,IF-Cell,IF-Tissue,FC,IHC-P
Cell treatment	Phospho-POLR2A (S2) Recombinant Rabbit mAb [JM11-51]	ET1703-86	Human, Mouse, Rat	WB,IP,IF-Cell,IF-Tissue,IHC-P,FC
Cell treatment	Phospho-POLR2A (S5) Recombinant Rabbit mAb [JM51-21]	ET1703-87	Human, Mouse, Rat	WB,IP,IF-Cell,IF-Tissue,IHC-P,FC
	Phospho-PP2A (Y307) Recombinant Rabbit mAb [ST49-05]	ET1609-40	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
Cell treatment	Phospho-PPP1CA (T320) Recombinant Rabbit mAb [PSH10-27]	HA723194	Human, Mouse, Rat	WB,IHC-P,IP
Cell treatment	Phospho-PRAS40 (T246) Recombinant Rabbit mAb [JE53-94]	HA721790	Human, Mouse, Rat	WB
	Phospho-PTEN (S380) Recombinant Rabbit mAb [JJ08-81]	ET1701-46	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-Rac1+Cdc42 (S71) Recombinant Rabbit mAb [JJ086-9]	ET1701-40	Human	WB,IF-Cell,IHC-P
Cell treatment	Phospho-Rad17 (S656) Recombinant Rabbit mAb [JE63-30]	HA722148	Human, Mouse	WB,IF-Cell,FC
Cell treatment	Phospho-Raf1 (S259) Recombinant Rabbit mAb [SD85-07]	ET1612-87	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	Phospho-Raf1 (S43) Recombinant Rabbit mAb [SD087-05]	ET1612-75	Human, Mouse, Rat	WB
	Phospho-Raf1 (S621) Recombinant Rabbit mAb [JJ085-05]	ET1701-3	Human, Rat, Mouse	WB,IHC-P
Cell treatment	Phospho-Rb (S608) Recombinant Rabbit mAb [PSH11-04]	HA723271	Human, Monkey	WB
Cell treatment	Phospho-Rb (S780) Recombinant Rabbit mAb [JE44-61]	HA721513	Human, Mouse, Rat	WB
Cell treatment	Phospho-Rb (S795) Recombinant Rabbit mAb [PSH06-89]	HA722764	Human	WB,IF-Cell
Cell treatment	Phospho-Rb (S807) Recombinant Rabbit mAb [SR3-82]	ET1602-36	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	Phospho-Rb (S807/S811) Recombinant Rabbit mAb [PSH11-01]	HA723268	Human, Mouse	WB,IHC-P,IF-Cell,FC,IP
Cell treatment	Phospho-RICTOR (T1135) Recombinant Rabbit mAb [PSH09-92]	HA723161	Human, Mouse	WB
Cell treatment	Phospho-RIP3 (S232) Recombinant Rabbit mAb [JE44-62]	HA721428	Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-RPA32/RPA2 (T21) Recombinant Rabbit mAb [SN06-36]	ET1611-18	Human, Mouse, Rat	WB,IP,IF-Cell
	Phospho-RPS6 (S235 + S236) Recombinant Rabbit mAb [PS01-04]	HA721275	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Phospho-RPS6 (S240 + S244) Recombinant Rabbit mAb [PS00-65]	HA721589	Human, Mouse, Rat	WB,IHC-P,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Phospho-RSK1 p90 (T359 + S363) Recombinant Rabbit mAb [SU03-65]	ET1608-53	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
Cell treatment	Phospho-RSK1 p90 (T573) Recombinant Rabbit mAb [JE50-12]	HA721793	Human, Mouse, Rat	WB
Cell treatment	Phospho-RSK1 (S380) Recombinant Rabbit mAb [SC05-32]	ET1610-28	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IP
	Phospho-SHIP (Y1020) Recombinant Rabbit mAb [JJ082-03]	ET1701-5	Human	WB,IF-Cell,IP
Cell treatment	Phospho-SHP2 (Y542) Recombinant Rabbit mAb [SN61-01]	ET1611-22	Mouse, Human	WB,IP
Cell treatment	Phospho-SIRT1 (S47) Recombinant Rabbit mAb [SN06-40]	ET1611-31	Human	WB,IF-Cell,IHC-P
	Phospho-SIRT1 (T530) Recombinant Rabbit mAb [JJ206-6]	ET1701-27	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
Cell treatment	Phospho-Smad1/5/9 (S463/S465/S467) Recombinant Rabbit mAb [JE59-46]	HA722566	Human, Mouse, Rat	WB,IHC,IF-Cell
Cell treatment	Phospho-Smad2 (S250) Recombinant Rabbit mAb [SD207-1]	ET1612-32	Human, Mouse, Rat	WB,FC,IP,IHC-P
Cell treatment	Phospho-Smad2 (S255) Recombinant Rabbit mAb [JF0882]	ET1702-34	Human, Mouse, Rat	WB,IHC-P,IP
Cell treatment	Phospho-SMAD2 (S465+467) Recombinant Rabbit mAb [PSH05-76]	HA722443	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-Smad3 (S423 + S425) Recombinant Rabbit mAb [ST0493]	ET1609-41	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	Phospho-Smad5 (S463 + S465) Recombinant Rabbit mAb [SY09-03]	ET1605-5	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-SMC1 (S957) Recombinant Rabbit mAb [SC65-03]	ET1610-55	Human	WB,IHC-P,IF-Cell
Cell treatment	Phospho-Src (S17) Recombinant Rabbit mAb [PSH11-03]	HA723270	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC,IP
Cell treatment	Phospho-Src (Y419) Recombinant Rabbit mAb [ST0800]	ET1609-15	Human, Mouse, Rat	WB,FC,IF-Cell
Cell treatment, Relative expression	Phospho-STAT1 (S727) Recombinant Rabbit mAb [SN67-04]	ET1611-20	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-STAT1 (Y701) Recombinant Rabbit mAb [PSH04-02]	HA722083	Human	WB,IF-Cell
Cell treatment, KD	Phospho-STAT3 (S727) Recombinant Rabbit mAb [SY24-09]	ET1607-39	Human, Rat, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP
Cell treatment	Phospho-STAT3 (Y705) Recombinant Rabbit mAb [SZ43-01]	ET1603-40	Human, Mouse, Rat	WB,IHC-P,IP
	Phospho-STAT4 (Y693) Rabbit pAb	ER1919-55	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-STAT5 (Y694) Recombinant Rabbit mAb [SC05-31]	ET1610-48	Human, Rat, Mouse	WB,IHC-P,IF-Cell
Cell treatment	Phospho-STAT6 (Y641) Recombinant Rabbit mAb [PSH06-98]	HA722776	Human, Mouse	WB
Cell treatment	Phospho-STING (S366) Recombinant Rabbit mAb [PSH09-72]	HA723137	Human	WB,Dot Blot,IP
Cell treatment	Phospho-Syk (Y352) Recombinant Rabbit mAb [PSH11-66]	HA723362	Human	WB,IF-Cell,FC
Cell treatment	Phospho-Syk (Y525 + Y526) Recombinant Rabbit mAb [PSH03-70]	HA722034	Human	WB,IF-Cell
	Phospho-Synapsin I (S9) Recombinant Rabbit mAb [SN63-04]	ET1611-26	Mouse, Rat, Human	WB,IHC-P
Cell treatment	Phospho-TAK1 (S412) Recombinant Rabbit mAb [PSH14-64]	HA723662	Human, Mouse, Rat	WB,FC
Cell treatment	Phospho-TAK1 (S439) Recombinant Rabbit mAb [JE44-93]	HA722242	Human, Mouse, Rat	WB
Cell treatment	Phospho-Tau (S198) Recombinant Rabbit mAb [JE66-79]	HA721668	Human, Mouse, Rat	WB
Cell treatment	Phospho-Tau (S199) Recombinant Rabbit mAb [JE46-95]	HA721800	Human, Mouse, Rat	WB
Cell treatment	Phospho-Tau (S199) Recombinant Rabbit mAb [JE66-80]	HA721667	Human, Mouse, Rat	WB
Cell treatment	Phospho-Tau (S202) Recombinant Rabbit mAb [JE66-81]	HA721666	Human, Mouse, Rat	WB
Cell treatment	Phospho-Tau (S214) Recombinant Rabbit mAb [JE66-83]	HA721797	Human, Mouse, Rat	WB
Cell treatment	Phospho-Tau (S324) Recombinant Rabbit mAb [JE30-73]	HA721798	Human, Mouse, Rat	WB
Cell treatment	Phospho-Tau (S356) Recombinant Rabbit mAb [JE46-86]	HA721799	Human, Mouse, Rat	WB
Cell treatment	Phospho-Tau (S396) Recombinant Rabbit mAb [SN62-09]	ET1611-68	Mouse, Rat, Human, Cynomolgus monkey, Pig	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
Cell treatment	Phospho-Tau (S404) Recombinant Rabbit mAb [JE66-82]	HA721665	Human, Mouse, Rat	WB,IHC
Cell treatment	Phospho-Tau (S404) Recombinant Rabbit mAb [JE47-88]	HA721801	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-Tau (T181) Recombinant Rabbit mAb [PSH05-42]	HA722271	Human, Mouse, Cynomolgus monkey	WB,IHC-P
Cell treatment	Phospho-Tau (T217) Recombinant Rabbit mAb [PSH09-33]	HA723091	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IHC-Fr
Cell treatment	Phospho-Tau (T231) Recombinant Rabbit mAb [SC58-08]	ET1610-31	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IP,IF-Tissue,IHC-Fr

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Phospho-Tau (T231) Recombinant Rabbit mAb [PSH01-04]	HA721570	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Phospho-Tau (T231) Recombinant Rabbit mAb [PSH01-05]	HA721569	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phosphothreonine Recombinant Rabbit mAb [PSH14-34]	HA723628	Human, Mouse, Rat, Species independent	WB,IHC-P,IF-Cell,IP,ELISA
Relative expression	Phospho-TrkB (Y817) Recombinant Rabbit mAb [SC0556]	ET1610-35	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
Cell treatment	Phospho-Tuberin (S939) Recombinant Rabbit mAb [JE54-06]	HA721762	Human, Mouse	WB
Cell treatment	Phospho-Tuberin (T1462) Recombinant Rabbit mAb [JE43-81]	HA721792	Human, Mouse, Rat	WB
Cell treatment	Phospho-UBF1 (S484) Recombinant Rabbit mAb [JE77-43]	HA722477	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-ULK1 (S556) Recombinant Rabbit mAb [JE59-28]	HA721445	Human, Mouse	WB
Cell treatment	Phospho-VASP (S156) Recombinant Rabbit mAb [JE45-91]	HA722145	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-Vimentin (S39) Recombinant Rabbit mAb [JE43-26]	HA721442	Human, Mouse, Rat	WB,IF-Cell
	Phospho-V-Myb+C-Myb (S11) Recombinant Rabbit mAb [SZ04-81]	ET1603-41	Human	IF-Cell,IF-Tissue,IHC-P
Cell treatment	Phospho-YAP1 (S127) Recombinant Rabbit mAb [SN0718]	ET1611-69	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	RIP140 (Acetyl K158) Rabbit pAb	ER60057	Human, Mouse, Rat	WB
	Smad2/3 (Acetyl K19) Rabbit pAb	ER60061	Human, Mouse, Rat	IHC-P,WB
	SOD2 (acetyl K68) Recombinant Rabbit mAb [JE62-28]	HA722463	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层