

配对抗体-产品清单

高级验证	产品名	货号	反应种属	验证应用
	Biotin Conjugated beta Amyloid 1-40 Recombinant Rabbit mAb [PSH03-98]	HA722618B	Human, Mouse	ELISA(Det)
	beta Amyloid 1-40 Recombinant Rabbit mAb [PSH02-82] - BSA and Azide free (Capture)	HA722077	Human, Mouse	ELISA(Cap)
	beta Amyloid 1-40 Recombinant Rabbit mAb [PSH03-98] - BSA and Azide free (Detector)	HA722078	Human, Mouse	ELISA(Det)
	beta Amyloid 1-42 Recombinant Rabbit mAb [PSH02-83] - BSA and Azide free (Capture)	HA722079	Human	ELISA(Cap)
	beta Amyloid 1-42 Recombinant Rabbit mAb [PSH03-99] - BSA and Azide free (Detector)	HA722080	Human	ELISA(Det)
	Biotin Conjugated beta Amyloid 1-42 Recombinant Rabbit mAb [PSH03-99]	HA722293	Human	ELISA(Det)
	Biotin Conjugated Beta Nerve Growth Factor Recombinant Rabbit mAb [PSH07-81]	HA722905B	Human	ELISA(Det)
	Biotin Conjugated Brain-derived tau/BD-tau Recombinant Rabbit mAb [PSH16-47] - Detector	HA723849B	Human	ELISA(Det),ELISA
	Brain-derived tau/BD-tau Recombinant Rabbit mAb [PSH16-46] - BSA and Azide free (Capture)	HA723847	Human, Mouse, Rat	ELISA(Cap)
	Brain-derived tau/BD-tau Recombinant Rabbit mAb [PSH16-47] - BSA and Azide free (Detector)	HA723848	Human	ELISA(Det)
	cIAP1 Recombinant Rabbit mAb [PSH02-10] - BSA and Azide free (Capture)	HA723281	Human	ELISA(Cap)
	Biotin Conjugated cIAP1 Recombinant Rabbit mAb [PSH06-71] - Detector	HA723283B	Human, Mouse, Rat	ELISA(Det),ELISA
	cIAP1 Recombinant Rabbit mAb [PSH06-72] - BSA and Azide free (Capture)	HA723280	Human, Mouse, Rat	ELISA(Cap)
	cIAP1 Recombinant Rabbit mAb [PSH06-71] - BSA and Azide free (Detector)	HA723282	Human, Mouse, Rat	ELISA(Det)
	GFAP Recombinant Rabbit mAb [PSH05-16] - BSA and Azide free (Capture)	HA722234	Human, Mouse, Rat	ELISA(Cap)
	GFAP Recombinant Mouse mAb [PSH05-17] - BSA and Azide free (Detector)	HA601342	Human, Mouse, Rat	ELISA(Det)
	Biotin Conjugated GFAP Mouse mAb [PSH05-17]	HA601355B	Human, Mouse, Rat	ELISA(Det)
	Biotin Conjugated H-FABP Recombinant Rabbit mAb [PSH16-99] - Detector	HA725178B	Human, Mouse	ELISA(Det),ELISA

	H-FABP Recombinant Rabbit mAb [PSH16-98] - BSA and Azide free (Capture)	HA725177	Human, Mouse	ELISA(Cap)
	H-FABP Recombinant Rabbit mAb [PSH16-99] - BSA and Azide free (Detector)	HA725178	Human, Mouse	ELISA(Det)
	HIV1 p24 Recombinant Rabbit mAb [PS01-31] - BSA and Azide free (Detector)	HA721528	Hiv1	ELISA(Det)
	HIV1 p24 Recombinant Mouse mAb [12G4-R] - BSA and Azide free (Capture)	HA601104	Hiv1	ELISA(Cap)
	HIV1 p24 Mouse mAb [12G10] - BSA and Azide free (Capture)	HA601159	Hiv1	ELISA(Cap)
	Biotin Conjugated HIV1 p24 Recombinant Mouse mAb [12G2-R]	HA601360B	Hiv1	ELISA(Det)
	HIV1 p24 Recombinant Rabbit mAb [PS01-35] - BSA and Azide free (Capture)	HA721527	Hiv1	ELISA(Cap)
	HIV1 p24 Recombinant Mouse mAb [12G2-R] - BSA and Azide free (Detector)	HA601105	Hiv1	ELISA(Det)
	HIV1 p24 Mouse mAb [12G2] - BSA and Azide free (Detector)	HA600009	Hiv1	ELISA(Det),WB
	Biotin Conjugated HIV1 p24 Recombinant Rabbit mAb [PS01-31] - Detector	HA721529	Hiv1	ELISA(Det)
	HRP Conjugated HIV1 p24 Mouse mAb [12G2]	HA601160	Hiv1	ELISA(Det),WB
	Biotin Conjugated HIV1 p24 Mouse mAb [12G2]	HA601363B	Hiv1	ELISA(Det),WB
	Biotin Conjugated Human ACE2 Recombinant Rabbit mAb [PSH17-18] - Detector	HA723905B	Human	ELISA(Det),ELISA
	Human ACE2 Recombinant Rabbit mAb [PSH17-18] - BSA and Azide free (Detector)	HA723905	Human	ELISA(Det)
	Human ACE2 Recombinant Rabbit mAb [PSH17-17] - BSA and Azide free (Capture)	HA723904	Human	ELISA(Cap)
	Human Acetylcholinesterase Recombinant Rabbit mAb [PSH17-02] - BSA and Azide free (Capture)	HA725223	Human	ELISA(Cap)
	Human Acetylcholinesterase Recombinant Rabbit mAb [PSH17-03] - BSA and Azide free (Detector)	HA725224	Human	ELISA(Det)
	Biotin Conjugated Human Acetylcholinesterase Recombinant Rabbit mAb [PSH17-03] - Detector	HA725224B	Human	ELISA(Det),ELISA
	Human Adiponectin Recombinant Rabbit mAb [PSH05-31] - BSA and Azide free (Capture)	HA722260	Human	ELISA(Cap)
	Biotin Conjugated Human Adiponectin Recombinant Rabbit mAb [PSH05-32] - Detector	HA722595B	Human	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human Adiponectin Recombinant Rabbit mAb [PSH05-32] - BSA and Azide free (Detector)	HA722261	Human	ELISA(Det)
	Human alpha Fetoprotein Recombinant Rabbit mAb [PSH14-54] - BSA and Azide free (Detector)	HA723645	Human	ELISA(Det)
	Human alpha Fetoprotein Recombinant Rabbit mAb [PSH14-53] - BSA and Azide free (Capture)	HA723644	Human	ELISA(Cap)
	Biotin Conjugated Human alpha Fetoprotein Recombinant Rabbit mAb [PSH14-55] - Detector	HA723648B	Human	ELISA(Det)
	Biotin Conjugated Human alpha Fetoprotein Recombinant Rabbit mAb [PSH14-54] - Detector	HA723646B	Human	ELISA(Det)
	Human alpha Fetoprotein Recombinant Rabbit mAb [PSH14-55] - BSA and Azide free (Detector)	HA723647	Human	ELISA(Det)
	Human Alpha-Synuclein Recombinant Rabbit mAb [JB42-33] - BSA and Azide free (Capture)	HA723467	Human	ELISA(Cap)
	Biotin Conjugated Human Alpha-Synuclein Recombinant Rabbit mAb [PSH12-56] - Detector	HA723469B	Human	ELISA(Det),ELISA
	Human Alpha-Synuclein Recombinant Rabbit mAb [PSH12-56] - BSA and Azide free (Detector)	HA723468	Human	ELISA(Det)
	Human Angiogenin Recombinant Rabbit mAb [PSH12-65] - BSA and Azide free (Detector)	HA725077	Human	ELISA(Det)
	Human Angiogenin Recombinant Rabbit mAb [PSH12-64] - BSA and Azide free (Capture)	HA725076	Human	ELISA(Cap)
	HRP Conjugated Human Angiogenin Recombinant Rabbit mAb [PSH12-65] - Detector	HA725078H	Human	ELISA(Det),ELISA
	Biotin Conjugated Human Angiogenin Recombinant Rabbit mAb [PSH12-65] - Detector	HA725078B	Human	ELISA(Det),ELISA
	Human ANGPTL3 Recombinant Rabbit mAb [PSH17-39] - BSA and Azide free (Capture)	HA725079	Human	ELISA(Cap)
	Biotin Conjugated Human ANGPTL3 Recombinant Rabbit mAb [PSH17-40] - Detector	HA725081B	Human	ELISA(Det),ELISA
	Human ANGPTL3 Recombinant Rabbit mAb [PSH17-40] - BSA and Azide free (Detector)	HA725080	Human	ELISA(Det)
	Biotin Conjugated Human Antithrombin III / ATIII Recombinant Rabbit mAb [PSH14-24] - Detector	HA723614B	Human	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human Antithrombin III / ATIII Recombinant Rabbit mAb [PSH14-26] - BSA and Azide free (Detector)	HA723616	Human	ELISA(Det)
	Biotin Conjugated Human Antithrombin III / ATIII Recombinant Rabbit mAb [PSH14-26] - Detector	HA723617B	Human	ELISA(Det),ELISA
	HRP Conjugated Human Antithrombin III / ATIII Recombinant Rabbit mAb [PSH14-26] - Detector	HA723617H	Human	ELISA(Det),ELISA
	Human Antithrombin III / ATIII Recombinant Rabbit mAb [PSH14-25] - BSA and Azide free (Capture)	HA723615	Human	ELISA(Cap)
	Human Antithrombin III / ATIII Recombinant Rabbit mAb [PSH14-23] - BSA and Azide free (Capture)	HA723612	Human	ELISA(Cap)
	Human Antithrombin III / ATIII Recombinant Rabbit mAb [PSH14-24] - BSA and Azide free (Detector)	HA723613	Human	ELISA(Det)
	Biotin Conjugated Human B7-H6 Recombinant Rabbit mAb [PSH10-20] - Detector	HA723192B	Human	ELISA(Det),ELISA
	Human B7-H6 Recombinant Rabbit mAb [PSH10-19] - BSA and Azide free (Capture)	HA723190	Human	ELISA(Cap)
	Human B7-H6 Recombinant Rabbit mAb [PSH10-20] - BSA and Azide free (Detector)	HA723191	Human	ELISA(Det)
	Human BACE1 Recombinant Rabbit mAb [PSH16-22] - BSA and Azide free (Capture)	HA725171	Human	ELISA(Cap)
	Biotin Conjugated Human BACE1 Recombinant Rabbit mAb [PSH16-23] - Detector	HA725172B	Human	ELISA(Det),ELISA
	Human BACE1 Recombinant Rabbit mAb [PSH16-23] - BSA and Azide free (Detector)	HA725172	Human	ELISA(Det)
	Biotin Conjugated Human BCMA Recombinant Rabbit mAb [PSH13-19] - Detector	HA725084B	Human	ELISA(Det),ELISA
	Human BCMA Recombinant Rabbit mAb [PSH13-19] - BSA and Azide free (Detector)	HA725083	Human	ELISA(Det)
	Human BCMA Recombinant Rabbit mAb [PSH13-18] - BSA and Azide free (Capture)	HA725082	Human	ELISA(Cap)
	Human Beta Nerve Growth Factor Recombinant Rabbit mAb [PSH07-81] - BSA and Azide free (Detector)	HA722904	Human	ELISA(Det)
	Human Beta Nerve Growth Factor Recombinant Rabbit mAb [PSH07-80] - BSA and Azide free (Capture)	HA722903	Human	ELISA(Cap)
	Human Beta Nerve Growth Factor Recombinant Rabbit mAb [PSH07-79] - BSA and Azide free (Capture)	HA722902	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human Caspase-9 Recombinant Rabbit mAb [PSH17-20] - Detector	HA723907B	Human	ELISA(Det),ELISA
	Human Caspase-9 Recombinant Rabbit mAb [PSH17-19] - BSA and Azide free (Capture)	HA723906	Human	ELISA(Cap)
	Human Caspase-9 Recombinant Rabbit mAb [PSH17-20] - BSA and Azide free (Detector)	HA723907	Human	ELISA(Det)
	Biotin Conjugated Human CCL11/Eotaxin Recombinant Rabbit mAb [PSH16-25] - Detector	HA725198B	Human	ELISA(Det),ELISA
	Human CCL11/Eotaxin Recombinant Rabbit mAb [PSH16-24] - BSA and Azide free (Capture)	HA725197	Human	ELISA(Cap)
	Human CCL11/Eotaxin Recombinant Rabbit mAb [PSH16-25] - BSA and Azide free (Detector)	HA725198	Human	ELISA(Det)
	Human CCL17 / TARC Recombinant Rabbit mAb [PSH13-35] - BSA and Azide free (Detector)	HA723516	Human	ELISA(Det)
	Human CCL17 / TARC Recombinant Rabbit mAb [PSH13-33] - BSA and Azide free (Capture)	HA723513	Human	ELISA(Cap)
	Biotin Conjugated Human CCL17 / TARC Recombinant Rabbit mAb [PSH13-34] - Detector	HA723515B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human CCL17 / TARC Recombinant Rabbit mAb [PSH13-35] - Detector	HA723517B	Human	ELISA(Det),ELISA
	Human CCL17 / TARC Recombinant Rabbit mAb [PSH13-34] - BSA and Azide free (Detector)	HA723514	Human	ELISA(Det)
	Human CCL20/MIP-3 alpha Recombinant Rabbit mAb [PSH12-32] - BSA and Azide free (Capture)	HA725010	Human	ELISA(Cap)
	Human CCL20/MIP-3 alpha Recombinant Rabbit mAb [PSH12-33] - BSA and Azide free (Detector)	HA725011	Human	ELISA(Det)
	Biotin Conjugated Human CCL20/MIP-3 alpha Recombinant Rabbit mAb [PSH12-33] - Detector	HA725012B	Human	ELISA(Det),ELISA
	Human CCL21 Recombinant Rabbit mAb [PSH13-36] - BSA and Azide free (Capture)	HA723518	Human	ELISA(Cap)
	Biotin Conjugated Human CCL21 Recombinant Rabbit mAb [PSH13-37] - Detector	HA723520B	Human	ELISA(Det),ELISA
	Human CCL21 Recombinant Rabbit mAb [PSH13-37] - BSA and Azide free (Detector)	HA723519	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human CCL23 Recombinant Rabbit mAb [PSH09-05] - BSA and Azide free (Capture)	HA723061	Human	ELISA(Cap)
	Human CCL23 Recombinant Rabbit mAb [PSH09-07] - BSA and Azide free (Detector)	HA723063	Human	ELISA(Det)
	Human CCL23 Recombinant Rabbit mAb [PSH09-06] - BSA and Azide free (Capture)	HA723062	Human	ELISA(Cap)
	Biotin Conjugated Human CCL23 Recombinant Rabbit mAb [PSH09-07] - BSA and Azide free (Detector)	HA723064B	Human	ELISA(Det),ELISA
	Human CCL25 Recombinant Rabbit mAb [PSH08-67] - BSA and Azide free (Detector)	HA723017	Human	ELISA(Det)
	Human CCL25 Recombinant Rabbit mAb [PSH08-66] - BSA and Azide free (Detector)	HA723016	Human	ELISA(Det)
	Human CCL25 Recombinant Rabbit mAb [PSH08-65] - BSA and Azide free (Capture)	HA723015	Human	ELISA(Cap)
	Biotin Conjugated Human CCL25 Recombinant Rabbit mAb [PSH08-66] - Detector	HA723018B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human CCL25 Recombinant Rabbit mAb [PSH08-67] - Detector	HA723019B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human CCL2 / MCP-1 Recombinant Rabbit mAb [PSH18-22] - Detector	HA725196B	Human	ELISA(Det)
	Human CCL2 / MCP-1 Recombinant Rabbit mAb [PSH18-22] - BSA and Azide free (Detector)	HA725196	Human	ELISA(Det)
	Human CCL2 / MCP-1 Recombinant Rabbit mAb [PSH18-21] - BSA and Azide free (Capture)	HA725195	Human	ELISA(Cap)
	Human CCL3 / MIP-1 alpha Recombinant Rabbit mAb [PSH11-78] - BSA and Azide free (Capture)	HA723380	Human	ELISA(Cap)
	Human CCL3 / MIP-1 alpha Recombinant Rabbit mAb [PSH11-77] - BSA and Azide free (Detector)	HA723378	Human	ELISA(Det)
	Biotin Conjugated Human CCL3 / MIP-1 alpha Recombinant Rabbit mAb [PSH11-77] - Detector	HA723379B	Human	ELISA(Det)
	Human CCL3 / MIP-1 alpha Recombinant Rabbit mAb [PSH11-76] - BSA and Azide free (Capture)	HA723377	Human	ELISA(Cap)
	Human CCL5 Recombinant Rabbit mAb [PSH13-16] - BSA and Azide free (Capture)	HA725055	Human	ELISA(Cap)
	HRP Conjugated Human CCL5 Recombinant Rabbit mAb [PSH13-17] - Detector	HA725057H	Human	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human CCL5 Recombinant Rabbit mAb [PSH13-17] - Detector	HA725057B	Human	ELISA(Det),ELISA
	Human CCL5 Recombinant Rabbit mAb [PSH13-17] - BSA and Azide free (Detector)	HA725056	Human	ELISA(Det)
	Human CD137 Recombinant Rabbit mAb [PSH08-63] - BSA and Azide free (Capture)	HA723012	Human	ELISA(Cap)
	Biotin Conjugated Human CD137 Recombinant Rabbit mAb [PSH08-64] - Detector	HA723014B	Human	ELISA(Det),ELISA
	Human CD137 Recombinant Rabbit mAb [PSH08-62] - BSA and Azide free (Detector)	HA723010	Human	ELISA(Det)
	Human CD137 Recombinant Rabbit mAb [PSH08-61] - BSA and Azide free (Capture)	HA723009	Human	ELISA(Cap)
	Biotin Conjugated Human CD137 Recombinant Rabbit mAb [PSH08-62] - Detector	HA723011B	Human	ELISA(Det),ELISA
	Human CD137 Recombinant Rabbit mAb [PSH08-64] - BSA and Azide free (Detector)	HA723013	Human	ELISA(Det)
	Human CD147 Recombinant Rabbit mAb [PSH13-07] - BSA and Azide free (Detector)	HA725086	Human	ELISA(Det)
	Biotin Conjugated Human CD147 Recombinant Rabbit mAb [PSH13-07] - Detector	HA725087B	Human	ELISA(Det),ELISA
	Human CD147 Recombinant Rabbit mAb [PSH13-06] - BSA and Azide free (Capture)	HA725085	Human	ELISA(Cap)
	Human CD163 Recombinant Rabbit mAb [PSH13-78] - BSA and Azide free (Capture)	HA723566	Human	ELISA(Cap)
	Biotin Conjugated Human CD163 Recombinant Rabbit mAb [PSH13-79] - Detector	HA723568B	Human	ELISA(Det),ELISA
	Human CD163 Recombinant Rabbit mAb [PSH13-80] - BSA and Azide free (Capture)	HA723569	Human	ELISA(Cap)
	Human CD163 Recombinant Rabbit mAb [PSH13-79] - BSA and Azide free (Detector)	HA723567	Human	ELISA(Det)
	Human CD200 Recombinant Rabbit mAb [PSH05-60] - BSA and Azide free (Detector)	HA722380	Human	ELISA(Det)
	Human CD200 Recombinant Rabbit mAb [PSH05-59] - BSA and Azide free (Capture)	HA722379	Human	ELISA(Cap)
	Biotin Conjugated Human CD200 Recombinant Rabbit mAb [PSH05-60]	HA722591B	Human	ELISA(Det)
	Human CD200R1 Recombinant Rabbit mAb [PSH05-27] - BSA and Azide free (Capture)	HA722256	Human	ELISA(Cap)
	Biotin Conjugated Human CD200R1 Recombinant Rabbit mAb [PSH05-28]	HA722598B	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human CD200R1 Recombinant Rabbit mAb [PSH05-28] - BSA and Azide free (Detector)	HA722257	Human	ELISA(Det)
	Biotin Conjugated Human CD276 Recombinant Rabbit mAb [PSH05-93] - Detector	HA722588B	Human	ELISA(Det),ELISA
	Human CD276 Recombinant Rabbit mAb [PSH05-93] - BSA and Azide free (Detector)	HA722492	Human	ELISA(Det)
	Human CD276 Recombinant Rabbit mAb [PSH05-92] - BSA and Azide free (Capture)	HA722491	Human	ELISA(Cap)
	Biotin Conjugated Human CD27/TNFRSF7 Recombinant Rabbit mAb [PSH08-33] - Detector	HA722976B	Human	ELISA(Det),ELISA
	Human CD27/TNFRSF7 Recombinant Rabbit mAb [PSH08-31] - BSA and Azide free (Capture)	HA722973	Human	ELISA(Cap)
	Human CD27/TNFRSF7 Recombinant Rabbit mAb [PSH08-33] - BSA and Azide free (Detector)	HA722975	Human	ELISA(Det)
	Human CD27/TNFRSF7 Recombinant Rabbit mAb [PSH08-32] - BSA and Azide free (Capture)	HA722974	Human	ELISA(Cap)
	Biotin Conjugated Human CD28 Recombinant Rabbit mAb [PSH12-51] - Detector	HA723459B	Human	ELISA(Det),ELISA
	Human CD28 Recombinant Rabbit mAb [PSH12-50] - BSA and Azide free (Capture)	HA723457	Human	ELISA(Cap)
	Human CD28 Recombinant Rabbit mAb [PSH12-51] - BSA and Azide free (Detector)	HA723458	Human	ELISA(Det)
	Human CD30 Recombinant Rabbit mAb [PSH04-81] - BSA and Azide free (Capture)	HA722187	Human	ELISA(Cap)
	Biotin Conjugated Human CD30 Recombinant Rabbit mAb [PSH04-82]	HA722601B	Human	ELISA(Det)
	Human CD30 Recombinant Rabbit mAb [PSH04-82] - BSA and Azide free (Detector)	HA722188	Human	ELISA(Det)
	Biotin Conjugated Human CD31 Recombinant Rabbit mAb [PSH07-83] - Detector	HA722908B	Human	ELISA(Det)
	Human CD31 Recombinant Rabbit mAb [PSH07-83] - BSA and Azide free (Detector)	HA722907	Human	ELISA(Det)
	Human CD31 Recombinant Rabbit mAb [PSH07-82] - BSA and Azide free (Capture)	HA722906	Human	ELISA(Cap)
	Human CD33 Recombinant Rabbit mAb [PSH05-26] - BSA and Azide free (Detector)	HA722255	Human	ELISA(Det)
	Human CD33 Recombinant Rabbit mAb [PSH05-25] - BSA and Azide free (Capture)	HA722254	Human	ELISA(Cap)
	Biotin Conjugated Human CD33 Recombinant Rabbit mAb [PSH05-26]	HA722599B	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human CD47 Recombinant Rabbit mAb [PSH11-06] - Detector	HA723274B	Human	ELISA(Det),ELISA
	Human CD47 Recombinant Rabbit mAb [PSH11-06] - BSA and Azide free (Detector)	HA723273	Human	ELISA(Det)
	Human CD47 Recombinant Rabbit mAb [PSH11-05] - BSA and Azide free (Capture)	HA723272	Human	ELISA(Cap)
	Biotin Conjugated Human CD48 Recombinant Rabbit mAb [PSH08-38] - Detector	HA722981B	Human	ELISA(Det),ELISA
	Human CD48 Recombinant Rabbit mAb [PSH08-37] - BSA and Azide free (Capture/Detector)	HA722978	Human	ELISA(Cap),ELISA(Det)
	Human CD48 Recombinant Rabbit mAb [PSH08-38] - BSA and Azide free (Detector)	HA722979	Human	ELISA(Det)
	Biotin Conjugated Human CD48 Recombinant Rabbit mAb [PSH08-37] - Detector	HA722980B	Human	ELISA(Det),ELISA
	Human CD48 Recombinant Rabbit mAb [PSH08-36] - BSA and Azide free (Capture)	HA722977	Human	ELISA(Cap)
	Human CD62E Recombinant Rabbit mAb [PSH14-43] - BSA and Azide free (Detector)	HA723638	Human	ELISA(Det)
	Human CD62E Recombinant Rabbit mAb [PSH14-41] - BSA and Azide free (Capture)	HA723635	Human	ELISA(Cap)
	HRP Conjugated Human CD62E Recombinant Rabbit mAb [PSH14-42] - Detector	HA723637H	Human	ELISA(Det),ELISA
	Biotin Conjugated Human CD62E Recombinant Rabbit mAb [PSH14-43] - Detector	HA723639B	Human	ELISA(Det),ELISA
	Human CD62E Recombinant Rabbit mAb [PSH14-42] - BSA and Azide free (Capture/Detector)	HA723636	Human	ELISA(Cap),ELISA(Det)
	Biotin Conjugated Human CD62E Recombinant Rabbit mAb [PSH14-42] - Detector	HA723637B	Human	ELISA(Det),ELISA
	Human CD62L Recombinant Rabbit mAb [PSH13-85] - BSA and Azide free (Capture/Detector)	HA723576	Human	ELISA(Det),ELISA(Cap)
	Biotin Conjugated Human CD62L Recombinant Rabbit mAb [PSH13-85] - Detector	HA723577B	Human	ELISA(Det),ELISA
	Human CD62L Recombinant Rabbit mAb [PSH13-86] - BSA and Azide free (Detector)	HA723578	Human	ELISA(Det)
	Human CD62L Recombinant Rabbit mAb [PSH13-84] - BSA and Azide free (Capture)	HA723575	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human CD62L Recombinant Rabbit mAb [PSH13-86] - Detector	HA723579B	Human	ELISA(Det),ELISA
	Human CD63 Recombinant Rabbit mAb [PSH08-43] - BSA and Azide free (Capture)	HA722988	Human	ELISA(Cap)
	Human CD63 Recombinant Rabbit mAb [PSH08-44] - BSA and Azide free (Detector)	HA722989	Human	ELISA(Det)
	Biotin Conjugated Human CD63 Recombinant Rabbit mAb [PSH08-44] - Detector	HA722990B	Human	ELISA(Det),ELISA
	Human CD73 Recombinant Rabbit mAb [PSH08-15] - BSA and Azide free (Capture)	HA722957	Human	ELISA(Cap)
	Human CD73 Recombinant Rabbit mAb [PSH08-16] - BSA and Azide free (Detector)	HA722958	Human	ELISA(Det)
	Biotin Conjugated Human CD73 Recombinant Rabbit mAb [PSH08-16] - Detector	HA722959B	Human	ELISA(Det)
	Human CD80 Recombinant Rabbit mAb [PSH09-61] - BSA and Azide free (Capture)	HA723123	Human	ELISA(Cap)
	Biotin Conjugated Human CD80 Recombinant Rabbit mAb [PSH09-63] - Detector	HA723127B	Human	ELISA(Det),ELISA
	Human CD80 Recombinant Rabbit mAb [PSH09-62] - BSA and Azide free (Detector)	HA723124	Human	ELISA(Det)
	Human CD80 Recombinant Rabbit mAb [PSH09-63] - BSA and Azide free (Detector)	HA723125	Human	ELISA(Det)
	Biotin Conjugated Human CD80 Recombinant Rabbit mAb [PSH09-62] - Detector	HA723126B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human CD85d Recombinant Rabbit mAb [PSH05-62]	HA722592B	Human	ELISA(Det)
	Human CD85d Recombinant Rabbit mAb [PSH05-62] - BSA and Azide free (Detector)	HA722382	Human	ELISA(Det)
	Human CD85d Recombinant Rabbit mAb [PSH05-61] - BSA and Azide free (Capture)	HA722381	Human	ELISA(Cap)
	Human CHIT1 Recombinant Rabbit mAb [PSH15-89] - BSA and Azide free (Capture)	HA725205	Human	ELISA(Cap)
	Human CHIT1 Recombinant Rabbit mAb [PSH15-90] - BSA and Azide free (Detector)	HA725206	Human	ELISA(Det)
	Biotin Conjugated Human CHIT1 Recombinant Rabbit mAb [PSH15-90] - Detector	HA725206B	Human	ELISA(Det),ELISA
	Human cIAP2 Recombinant Rabbit mAb [PSH09-38] - BSA and Azide free (Capture)	HA723098	Human	ELISA(Cap)
	Biotin Conjugated Human cIAP2 Recombinant Rabbit mAb [PSH09-39] - Detector	HA723101B	Human	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human cIAP2 Recombinant Rabbit mAb [PSH09-40] - BSA and Azide free (Detector)	HA723100	Human	ELISA(Det)
	Biotin Conjugated Human cIAP2 Recombinant Rabbit mAb [PSH09-40] - Detector	HA723102B	Human	ELISA(Det),ELISA
	Human cIAP2 Recombinant Rabbit mAb [PSH09-39] - BSA and Azide free (Capture/Detector)	HA723099	Human	ELISA(Det),ELISA(Cap)
	Human CKMB Recombinant Rabbit mAb [PSH11-56] - BSA and Azide free (Detector)	HA723349	Human, Rat, Mouse, Monkey	ELISA(Det)
	Human CKMB Recombinant Rabbit mAb [PSH11-55] - BSA and Azide free (Capture)	HA723348	Human, Rat, Mouse, Monkey	ELISA(Cap)
	Biotin Conjugated Human CKMB Recombinant Rabbit mAb [PSH11-56] - Detector	HA723350B	Human, Rat, Mouse, Monkey	ELISA(Det),ELISA
	HRP Conjugated Human Clusterin Recombinant Rabbit mAb [PSH15-56] - Detector	HA725093H	Human	ELISA(Det),ELISA
	Human Clusterin Recombinant Rabbit mAb [PSH15-56] - BSA and Azide free (Detector)	HA725092	Human	ELISA(Det)
	Human Clusterin Recombinant Rabbit mAb [PSH15-55] - BSA and Azide free (Capture)	HA725091	Human	ELISA(Cap)
	Biotin Conjugated Human Clusterin Recombinant Rabbit mAb [PSH15-56] - Detector	HA725093B	Human	ELISA(Det),ELISA
	Human CNTF Recombinant Rabbit mAb [PSH11-82] - BSA and Azide free (Detector)	HA723385	Human, Rat	ELISA(Det)
	Human CNTF Recombinant Rabbit mAb [PSH11-81] - BSA and Azide free (Capture)	HA723384	Human, Rat	ELISA(Cap)
	Biotin Conjugated Human CNTF Recombinant Rabbit mAb [PSH11-80] - Detector	HA723383B	Human, Rat	ELISA(Det),ELISA
	Biotin Conjugated Human CNTF Recombinant Rabbit mAb [PSH11-82] - Detector	HA723386B	Human, Rat	ELISA(Det),ELISA
	Human CNTF Recombinant Rabbit mAb [PSH11-80] - BSA and Azide free (Detector)	HA723382	Human, Rat	ELISA(Det)
	Human CNTF Recombinant Rabbit mAb [PSH11-79] - BSA and Azide free (Capture)	HA723381	Human, Rat	ELISA(Cap)
	Human Complement C3 Recombinant Rabbit mAb [PSH15-07] - BSA and Azide free (Capture)	HA725140	Human	ELISA(Cap)
	HRP Conjugated Human Complement C3 Recombinant Rabbit mAb [PSH15-08] - Detector	HA725142H	Human	ELISA(Det),ELISA
	Biotin Conjugated Human Complement C3 Recombinant Rabbit mAb [PSH15-08] - Detector	HA725142B	Human	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human Complement C3 Recombinant Rabbit mAb [PSH15-08] - BSA and Azide free (Detector)	HA725141	Human	ELISA(Det)
	HRP Conjugated Human Complement C5 Recombinant Rabbit mAb [PSH06-46] - Detector	HA725127H	Human	ELISA(Det)
	Human Complement C5 Recombinant Rabbit mAb [PSH06-45] - BSA and Azide free (Capture)	HA722682	Human	ELISA(Cap)
	Biotin Conjugated Human Complement C5 Recombinant Rabbit mAb [PSH06-46]	HA722684B	Human	ELISA(Det)
	Human Complement C5 Recombinant Rabbit mAb [PSH06-46] - BSA and Azide free (Detector)	HA722683	Human	ELISA(Det)
	Human CRP Recombinant Rabbit mAb [PSH04-26] - BSA and Azide free (Detector)	HA722109	Human	ELISA(Det)
	Biotin Conjugated Human CRP Recombinant Rabbit mAb [PSH04-26] - Detector	HA722616B	Human	ELISA(Det)
	Human CRP Recombinant Rabbit mAb [PSH04-25] - BSA and Azide free (Capture)	HA722108	Human	ELISA(Cap)
	Human CSF-1-R Recombinant Rabbit mAb [PSH12-68] - BSA and Azide free (Capture)	HA723477	Human	ELISA(Cap)
	Biotin Conjugated Human CSF-1-R Recombinant Rabbit mAb [PSH12-66] - Detector	HA723478B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human CSF-1-R Recombinant Rabbit mAb [PSH12-67] - Detector	HA723476B	Human	ELISA(Det),ELISA
	HRP Conjugated Human CSF-1-R Recombinant Rabbit mAb [PSH12-66] - Detector	HA723478H	Human	ELISA(Det),ELISA
	Human CSF-1-R Recombinant Rabbit mAb [PSH12-66] - BSA and Azide free (Capture/Detector)	HA723474	Human	ELISA(Cap),ELISA(Det)
	Human CSF-1-R Recombinant Rabbit mAb [PSH12-67] - BSA and Azide free (Detector)	HA723475	Human	ELISA(Det)
	Human CSPG4 Recombinant Rabbit mAb [PSH10-77] - BSA and Azide free (Capture)	HA723244	Human	ELISA(Cap)
	Human CSPG4 Recombinant Rabbit mAb [PSH10-79] - BSA and Azide free (Detector)	HA723246	Human	ELISA(Det)
	Human CSPG4 Recombinant Rabbit mAb [PSH10-78] - BSA and Azide free (Capture)	HA723245	Human	ELISA(Cap)
	Biotin Conjugated Human CSPG4 Recombinant Rabbit mAb [PSH10-79] - Detector	HA723247B	Human	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	HRP Conjugated Human CSPG4 Recombinant Rabbit mAb [PSH10-79] - Detector	HA723376H	Human	ELISA(Det)
	Biotin Conjugated Human CTGF Recombinant Rabbit mAb [PSH11-16] - Detector	HA723298B	Human	ELISA(Det),ELISA
	Human CTGF Recombinant Rabbit mAb [PSH11-15] - BSA and Azide free (Capture)	HA723296	Human	ELISA(Cap)
	Human CTGF Recombinant Rabbit mAb [PSH11-16] - BSA and Azide free (Detector)	HA723297	Human	ELISA(Det)
	HRP Conjugated Human CTGF Recombinant Rabbit mAb [PSH11-16] - Detector	HA723375H	Human	ELISA(Det)
	Biotin Conjugated Human CTLA-4 Recombinant Rabbit mAb [PSH08-14]	HA722956B	Human	ELISA(Det)
	Human CTLA-4 Recombinant Rabbit mAb [PSH08-13] - BSA and Azide free (Capture)	HA722954	Human	ELISA(Cap)
	Human CTLA-4 Recombinant Rabbit mAb [PSH08-14] - BSA and Azide free (Detector)	HA722955	Human	ELISA(Det)
	Human CX3CL1 Recombinant Rabbit mAb [PSH16-09] - BSA and Azide free (Detector)	HA723814	Human	ELISA(Det)
	Biotin Conjugated Human CX3CL1 Recombinant Rabbit mAb [PSH16-09] - Detector	HA723814B	Human	ELISA(Det),ELISA
	Human CX3CL1 Recombinant Rabbit mAb [PSH16-08] - BSA and Azide free (Capture)	HA723813	Human	ELISA(Cap)
	Biotin Conjugated Human CXADR Recombinant Rabbit mAb [PSH17-47] - Detector	HA725294B	Human	ELISA(Det),ELISA
	Human CXADR Recombinant Rabbit mAb [PSH17-47] - BSA and Azide free (Detector)	HA725294	Human	ELISA(Det)
	HRP Conjugated Human CXADR Recombinant Rabbit mAb [PSH17-47] - Detector	HA725294H	Human	ELISA(Det),ELISA
	Human CXADR Recombinant Rabbit mAb [PSH17-46] - BSA and Azide free (Capture)	HA725293	Human	ELISA(Cap)
	Biotin Conjugated Human CXCL10 / IP-10 Recombinant Rabbit mAb [PSH12-37] - Detector	HA725009B	Human	ELISA(Det),ELISA
	Human CXCL10 / IP-10 Recombinant Rabbit mAb [PSH12-37] - BSA and Azide free (Detector)	HA725008	Human	ELISA(Det)
	Human CXCL10 / IP-10 Recombinant Rabbit mAb [PSH12-36] - BSA and Azide free (Capture)	HA725007	Human	ELISA(Cap)
Cell treatment	Human CXCL13 Recombinant Rabbit mAb [PSH07-06] - BSA and Azide free (Detector)	HA722792	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Human CXCL13 Recombinant Rabbit mAb [PSH07-05] - BSA and Azide free (Capture)	HA722791	Human	ELISA(Cap)
Cell treatment	Biotin Conjugated Human CXCL13 Recombinant Rabbit mAb [PSH07-07] - Detector	HA722795B	Human	ELISA(Det),ELISA
Cell treatment	Biotin Conjugated Human CXCL13 Recombinant Rabbit mAb [PSH07-06] - Detector	HA722793B	Human	ELISA(Det),ELISA
Cell treatment	Human CXCL13 Recombinant Rabbit mAb [PSH07-07] - BSA and Azide free (Detector)	HA722794	Human	ELISA(Det)
	Human CXCL1/GRO alpha Recombinant Rabbit mAb [PSH05-94] - BSA and Azide free (Capture)	HA722493	Human	ELISA(Cap)
	Biotin Conjugated Human CXCL1/GRO alpha Recombinant Rabbit mAb [PSH05-95]	HA722587B	Human	ELISA(Det)
	Human CXCL1/GRO alpha Recombinant Rabbit mAb [PSH05-95] - BSA and Azide free (Detector)	HA722494	Human	ELISA(Det)
	Human CXCL2/MIP-2 Recombinant Rabbit mAb [PSH17-89] - BSA and Azide free (Detector)	HA725248	Human	ELISA(Det)
	Biotin Conjugated Human CXCL2/MIP-2 Recombinant Rabbit mAb [PSH17-89] - Detector	HA725248B	Human	ELISA(Det),ELISA
	Human CXCL2/MIP-2 Recombinant Rabbit mAb [PSH17-88] - BSA and Azide free (Capture)	HA725247	Human	ELISA(Cap)
	Human CXCL6/GCP-2 Recombinant Rabbit mAb [PSH11-31] - BSA and Azide free (Detector)	HA723315	Human	ELISA(Det)
	Human CXCL6/GCP-2 Recombinant Rabbit mAb [PSH11-30] - BSA and Azide free (Capture)	HA723314	Human	ELISA(Cap)
	Biotin Conjugated Human CXCL6/GCP-2 Recombinant Rabbit mAb [PSH11-31] - Detector	HA723316B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human CXCL9 Recombinant Rabbit mAb [PSH12-35] - Detector	HA725006B	Human	ELISA(Det),ELISA
	Human CXCL9 Recombinant Rabbit mAb [PSH12-35] - BSA and Azide free (Detector)	HA725005	Human	ELISA(Det)
	Human CXCL9 Recombinant Rabbit mAb [PSH12-34] - BSA and Azide free (Capture)	HA725004	Human	ELISA(Cap)
	Biotin Conjugated Human DR3/TNFRSF25 Recombinant Rabbit mAb [PSH16-68] - Detector	HA723870B	Human	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human DR3/TNFRSF25 Recombinant Rabbit mAb [PSH16-69] - BSA and Azide free (Detector)	HA723871	Human	ELISA(Det)
	Biotin Conjugated Human DR3/TNFRSF25 Recombinant Rabbit mAb [PSH16-69] - Detector	HA723871B	Human	ELISA(Det),ELISA
	Human DR3/TNFRSF25 Recombinant Rabbit mAb [PSH16-67] - BSA and Azide free (Capture)	HA723869	Human	ELISA(Cap)
	Human DR3/TNFRSF25 Recombinant Rabbit mAb [PSH16-68] - BSA and Azide free (Detector)	HA723870	Human	ELISA(Det)
	Human E-Cadherin Recombinant Rabbit mAb [SY0287] - BSA and Azide free (Detector)	HA723056	Human	ELISA(Det)
	Human E-Cadherin Recombinant Rabbit mAb [PSH09-02] - BSA and Azide free (Capture)	HA723055	Human	ELISA(Cap)
	Biotin Conjugated Human E-Cadherin Recombinant Rabbit mAb [SY0287] - BSA and Azide free (Detector)	HA723057B	Human	ELISA(Det),ELISA
	Human EGF Recombinant Rabbit mAb [PSH02-19] - BSA and Azide free (Capture)	HA721782	Human	ELISA(Cap)
	Human EGF Recombinant Rabbit mAb [PSH02-20] - BSA and Azide free (Detector)	HA721783	Human	ELISA(Det)
	Biotin Conjugated Human EGF Recombinant Rabbit mAb [PSH02-20] - Detector	HA722638B	Human	ELISA(Det)
	Biotin Conjugated Human Endoglin / CD105 Recombinant Rabbit mAb [PSH14-59] - Detector	HA723656B	Human	ELISA(Det),ELISA
	Human Endoglin / CD105 Recombinant Rabbit mAb [PSH14-59] - BSA and Azide free (Detector)	HA723655	Human	ELISA(Det)
	Human Endoglin / CD105 Recombinant Rabbit mAb [PSH14-58] - BSA and Azide free (Capture)	HA723654	Human	ELISA(Cap)
	Human Endoglin / CD105 Recombinant Rabbit mAb [PSH14-60] - BSA and Azide free (Capture)	HA723657	Human	ELISA(Cap)
	Human Endostatin Recombinant Rabbit mAb [PSH15-83] - BSA and Azide free (Detector)	HA723790	Human	ELISA(Det)
	Biotin Conjugated Human Endostatin Recombinant Rabbit mAb [PSH15-83] - Detector	HA723791B	Human	ELISA(Det),ELISA
	Human Endostatin Recombinant Rabbit mAb [PSH15-82] - BSA and Azide free (Capture)	HA723789	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	HRP Conjugated Human Endostatin Recombinant Rabbit mAb [PSH15-83] - Detector	HA723791H	Human	ELISA(Det),ELISA
	Human EphA2 Recombinant Rabbit mAb [PSH15-13] - BSA and Azide free (Detector)	HA723715	Human	ELISA(Det)
	Human EphA2 Recombinant Rabbit mAb [PSH15-12] - BSA and Azide free (Capture)	HA723714	Human	ELISA(Cap)
	Biotin Conjugated Human EphA2 Recombinant Rabbit mAb [PSH15-13] - Detector	HA723716B	Human	ELISA(Det),ELISA
	Human Fas / CD95 Recombinant Rabbit mAb [PSH16-37] - BSA and Azide free (Capture)	HA723835	Human	ELISA(Cap)
	Human Fas / CD95 Recombinant Rabbit mAb [PSH16-39] - BSA and Azide free (Capture)	HA723837	Human	ELISA(Cap)
	Biotin Conjugated Human Fas / CD95 Recombinant Rabbit mAb [PSH16-38] - Detector	HA723836B	Human	ELISA(Det),ELISA
	Human Fas / CD95 Recombinant Rabbit mAb [PSH16-38] - BSA and Azide free (Detector)	HA723836	Human	ELISA(Det)
	Human FGF-19 Recombinant Rabbit mAb [PSH13-29] - BSA and Azide free (Capture)	HA725100	Human	ELISA(Cap)
	Biotin Conjugated Human FGF-19 Recombinant Rabbit mAb [PSH13-30] - Detector	HA725102B	Human	ELISA(Det),ELISA
	Human FGF-19 Recombinant Rabbit mAb [PSH13-30] - BSA and Azide free (Detector)	HA725101	Human	ELISA(Det)
	Human FGF-2 Recombinant Rabbit mAb [PSH02-08] - BSA and Azide free (Capture)	HA721780	Human	ELISA(Cap)
	Biotin Conjugated Human FGF-2 Recombinant Rabbit mAb [PSH02-09]	HA722639B	Human	ELISA(Det)
	Human FGF-2 Recombinant Rabbit mAb [PSH02-09] - BSA and Azide free (Detector)	HA721781	Human	ELISA(Det)
	Human FOLR1 Recombinant Rabbit mAb [PSH11-59] - BSA and Azide free (Capture)	HA723354	Human	ELISA(Cap)
	Human FOLR1 Recombinant Rabbit mAb [PSH11-60] - BSA and Azide free (Detector)	HA723355	Human	ELISA(Det)
	Biotin Conjugated Human FOLR1 Recombinant Rabbit mAb [PSH11-60] - Detector	HA723356B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human G-CSF Recombinant Rabbit mAb [PSH13-13] - Detector	HA725003B	Human	ELISA(Det),ELISA
	Human G-CSF Recombinant Rabbit mAb [PSH13-12] - BSA and Azide free (Capture)	HA725001	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human G-CSF Recombinant Rabbit mAb [PSH13-13] - BSA and Azide free (Detector)	HA725002	Human	ELISA(Det)
	Human GDF-15 Recombinant Rabbit mAb [PSH16-02] - BSA and Azide free (Detector)	HA723807	Human	ELISA(Det)
	Biotin Conjugated Human GDF-15 Recombinant Rabbit mAb [PSH16-01] - Detector	HA723806B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human GDF-15 Recombinant Rabbit mAb [PSH16-02] - Detector	HA723807B	Human	ELISA(Det),ELISA
	Human GDF-15 Recombinant Rabbit mAb [PSH16-01] - BSA and Azide free (Detector)	HA723806	Human	ELISA(Det)
	Human GDF-15 Recombinant Rabbit mAb [PSH16-00] - BSA and Azide free (Capture)	HA723805	Human	ELISA(Cap)
	Biotin Conjugated Human GITR Recombinant Rabbit mAb [PSH05-64]	HA722593B	Human	ELISA(Det)
	Human GITR Recombinant Rabbit mAb [PSH05-63] - BSA and Azide free (Capture)	HA722383	Human	ELISA(Cap)
	Human GITR Recombinant Rabbit mAb [PSH05-64] - BSA and Azide free (Detector)	HA722384	Human	ELISA(Det)
	Human GM-CSF Recombinant Rabbit mAb [PS01-38] - BSA and Azide free (Detector)	HA721289	Human	ELISA(Det)
	Biotin Conjugated Human GM-CSF Recombinant Rabbit mAb [PS01-38]	HA722648B	Human	ELISA(Det)
	Human GM-CSF Recombinant Rabbit mAb [PS01-52] - BSA and Azide free (Capture)	HA721288	Human	ELISA(Cap)
	Human Granzyme A Recombinant Rabbit mAb [PSH09-34] - BSA and Azide free (Capture)	HA723092	Human	ELISA(Cap)
	Biotin Conjugated Human Granzyme A Recombinant Rabbit mAb [PSH09-35] - Detector	HA723094B	Human	ELISA(Det),ELISA
	Human Granzyme A Recombinant Rabbit mAb [PSH09-35] - BSA and Azide free (Detector)	HA723093	Human	ELISA(Det)
	Human Granzyme A Recombinant Rabbit mAb [PSH09-36] - BSA and Azide free (Capture)	HA723095	Human	ELISA(Cap)
	Human Granzyme A Recombinant Rabbit mAb [PSH09-37] - BSA and Azide free (Detector)	HA723096	Human	ELISA(Det)
	Biotin Conjugated Human Granzyme A Recombinant Rabbit mAb [PSH09-37] - Detector	HA723097B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human Granzyme B Recombinant Rabbit mAb [PSH15-99] - Detector	HA723804B	Human	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human Granzyme B Recombinant Rabbit mAb [PSH15-98] - BSA and Azide free (Capture)	HA723803	Human	ELISA(Cap)
	Human Granzyme B Recombinant Rabbit mAb [PSH15-99] - BSA and Azide free (Detector)	HA723804	Human	ELISA(Det)
	Biotin Conjugated Human Growth Hormone Recombinant Rabbit mAb [PSH12-82] - Detector	HA725108B	Human	ELISA(Det),ELISA
	Human Growth Hormone Recombinant Rabbit mAb [PSH12-81] - BSA and Azide free (Capture)	HA725106	Human	ELISA(Cap)
	Human Growth Hormone Recombinant Rabbit mAb [PSH12-82] - BSA and Azide free (Detector)	HA725107	Human	ELISA(Det)
	Human HGF Recombinant Rabbit mAb [PSH03-78] - BSA and Azide free (Capture)	HA722042	Human	ELISA(Cap)
	Human HGF Recombinant Rabbit mAb [PSH02-06] - BSA and Azide free (Capture)	HA721778	Human	ELISA(Cap)
	Human HGF Recombinant Rabbit mAb [PSH03-79] - BSA and Azide free (Detector)	HA722043	Human	ELISA(Det)
	Human HGF Recombinant Rabbit mAb [PSH02-07] - BSA and Azide free (Detector)	HA721779	Human	ELISA(Det)
	Biotin Conjugated Human HGF Recombinant Rabbit mAb [PSH02-07]	HA722640B	Human	ELISA(Det)
	Biotin Conjugated Human HGF Recombinant Rabbit mAb [PSH03-79]	HA722619B	Human	ELISA(Det)
	Human ICAM1 Recombinant Rabbit mAb [PSH10-15] - BSA and Azide free (Capture)	HA723184	Human	ELISA(Cap)
	Human ICAM1 Recombinant Rabbit mAb [PSH10-16] - BSA and Azide free (Detector)	HA723185	Human	ELISA(Det)
	Biotin Conjugated Human ICAM1 Recombinant Rabbit mAb [PSH10-16] - Detector	HA723186B	Human	ELISA(Det),ELISA
	HRP Conjugated Human ICAM1 Recombinant Rabbit mAb [PSH10-16] - Detector	HA723474H	Human	ELISA(Det),ELISA
	Human ICOSLG/B7-H2/CD275 Recombinant Rabbit mAb [PSH11-33] - BSA and Azide free (Detector)	HA723318	Human	ELISA(Det)
	Biotin Conjugated Human ICOSLG/B7-H2/CD275 Recombinant Rabbit mAb [PSH11-33] - Detector	HA723319B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human ICOSLG/B7-H2/CD275 Recombinant Rabbit mAb [PSH11-35] - Detector	HA723322B	Human	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human ICOSLG/B7-H2/CD275 Recombinant Rabbit mAb [PSH11-35] - BSA and Azide free (Detector)	HA723321	Human	ELISA(Det)
	Human ICOSLG/B7-H2/CD275 Recombinant Rabbit mAb [PSH11-32] - BSA and Azide free (Capture)	HA723317	Human	ELISA(Cap)
	Human ICOSLG/B7-H2/CD275 Recombinant Rabbit mAb [PSH11-34] - BSA and Azide free (Capture)	HA723320	Human	ELISA(Cap)
	Biotin Conjugated Human IDO1 Recombinant Rabbit mAb [PSH09-87] - Detector	HA723156B	Human	ELISA(Det),ELISA
	Human IDO1 Recombinant Rabbit mAb [PSH09-86] - BSA and Azide free (Capture)	HA723154	Human	ELISA(Cap)
	Human IDO1 Recombinant Rabbit mAb [PSH09-87] - BSA and Azide free (Detector)	HA723155	Human	ELISA(Det)
	Biotin Conjugated Human IFN-alpha Recombinant Rabbit mAb [PS00-79] - BSA and Azide free (Detector)	HA721369	Human	ELISA(Det)
	Human IFN-alpha Recombinant Rabbit mAb [PS00-79] - BSA and Azide free (Detector)	HA721291	Human	ELISA(Det)
	Human IFN-alpha Recombinant Rabbit mAb [PS00-80] - BSA and Azide free (Capture)	HA721290	Human	ELISA(Cap)
	Human IFN beta Recombinant Rabbit mAb [PSH17-29] - BSA and Azide free (Capture)	HA725052	Human	ELISA(Cap)
	Biotin Conjugated Human IFN beta Recombinant Rabbit mAb [PSH17-30] - Detector	HA725054B	Human	ELISA(Det),ELISA
	Human IFN beta Recombinant Rabbit mAb [PSH17-30] - BSA and Azide free (Detector)	HA725053	Human	ELISA(Det)
	Biotin Conjugated Human IFNGR1 Recombinant Rabbit mAb [PSH06-44]	HA722681B	Human	ELISA(Det)
	Human IFNGR1 Recombinant Rabbit mAb [PSH06-44] - BSA and Azide free (Detector)	HA722680	Human	ELISA(Det)
	Human IFNGR1 Recombinant Rabbit mAb [PSH06-43] - BSA and Azide free (Capture)	HA722679	Human	ELISA(Cap)
	Biotin Conjugated Human IGFBP1 Recombinant Rabbit mAb [PSH17-57] - Detector	HA723932B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human IGFBP1 Recombinant Rabbit mAb [PSH17-58] - Detector	HA723933B	Human	ELISA(Det),ELISA
	Human IGFBP1 Recombinant Rabbit mAb [PSH17-56] - BSA and Azide free (Capture)	HA723931	Human	ELISA(Cap)
	Human IGFBP1 Recombinant Rabbit mAb [PSH17-57] - BSA and Azide free (Detector)	HA723932	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human IGFBP1 Recombinant Rabbit mAb [PSH17-58] - BSA and Azide free (Detector)	HA723933	Human	ELISA(Det)
	Biotin Conjugated Human IGFBP3 Recombinant Rabbit mAb [PSH12-17] - Detector	HA723434B	Human	ELISA(Det),ELISA
	HRP Conjugated Human IGFBP3 Recombinant Rabbit mAb [PSH12-17] - Detector	HA723565H	Human	ELISA(Det)
	Human IGFBP3 Recombinant Rabbit mAb [PSH12-17] - BSA and Azide free (Detector)	HA723433	Human	ELISA(Det)
	Human IGFBP3 Recombinant Rabbit mAb [PSH12-16] - BSA and Azide free (Capture)	HA723432	Human	ELISA(Cap)
	Human IL-10 Recombinant Rabbit mAb [PS00-81] - BSA and Azide free (Detector)	HA721237	Human	ELISA(Det)
	Biotin Conjugated Human IL-10 Recombinant Rabbit mAb [PS00-81]	HA722652B	Human	ELISA(Det)
	Human IL-10 Recombinant Rabbit mAb [PS00-98] - BSA and Azide free (Capture)	HA721236	Human	ELISA(Cap)
	Human IL-10RB Recombinant Rabbit mAb [PSH17-44] - BSA and Azide free (Capture)	HA725295	Human	ELISA(Cap)
	Biotin Conjugated Human IL-10RB Recombinant Rabbit mAb [PSH17-45] - Detector	HA725296B	Human	ELISA(Det),ELISA
	Human IL-10RB Recombinant Rabbit mAb [PSH17-45] - BSA and Azide free (Detector)	HA725296	Human	ELISA(Det)
	Biotin Conjugated Human IL-12B p40 Recombinant Rabbit mAb [PSH03-04]	HA722629B	Human	ELISA(Det)
	Biotin Conjugated Human IL-12B p40 Recombinant Rabbit mAb [PSH02-58]	HA722633B	Human	ELISA(Det)
	Human IL-12B p40 Recombinant Rabbit mAb [PSH02-58] - BSA and Azide free (Detector)	HA721835	Human	ELISA(Det)
	Human IL-12B p40 Recombinant Rabbit mAb [PSH03-04] - BSA and Azide free (Detector)	HA721934	Human	ELISA(Det)
	Human IL-12B p40 Recombinant Rabbit mAb [PSH02-57] - BSA and Azide free (Capture)	HA721834	Human	ELISA(Cap)
	Human IL-13 Recombinant Rabbit mAb [PSH12-38] - BSA and Azide free (Capture)	HA725013	Human	ELISA(Cap)
	Biotin Conjugated Human IL-13 Recombinant Rabbit mAb [PSH12-39] - Detector	HA725015B	Human	ELISA(Det),ELISA
	Human IL-13 Recombinant Rabbit mAb [PSH12-39] - BSA and Azide free (Detector)	HA725014	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human IL-13R alpha 2 Recombinant Rabbit mAb [PSH07-64] - BSA and Azide free (Detector)	HA722881	Human	ELISA(Det)
	Human IL-13R alpha 2 Recombinant Rabbit mAb [PSH07-63] - BSA and Azide free (Capture)	HA722880	Human	ELISA(Cap)
	Biotin Conjugated Human IL-13R alpha 2 Recombinant Rabbit mAb [PSH07-64]	HA722882B	Human	ELISA(Det)
	Biotin Conjugated Human IL-15 Recombinant Rabbit mAb [PSH04-20]	HA722612B	Human	ELISA(Det)
	Human IL-15 Recombinant Rabbit mAb [PSH04-19] - BSA and Azide free (Capture)	HA722102	Human	ELISA(Cap)
	Human IL-15 Recombinant Rabbit mAb [PSH04-20] - BSA and Azide free (Detector)	HA722103	Human	ELISA(Det)
	Biotin Conjugated Human IL-17A Recombinant Rabbit mAb [PS01-24] - Detector	HA721631	Human	ELISA(Det)
	Human IL-17A Recombinant Rabbit mAb [PS01-24] - BSA and Azide free (Detector)	HA721233	Human	ELISA(Det)
	Human IL-17A Recombinant Rabbit mAb [PS01-20] - BSA and Azide free (Capture)	HA721232	Human	ELISA(Cap)
	HRP Conjugated Human IL-17A Recombinant Rabbit mAb [PS01-24]	HA721296	Human	ELISA(Det)
	Human IL-18 Recombinant Rabbit mAb [PSH01-71] - BSA and Azide free (Capture)	HA721722	Human	ELISA(Cap)
	Biotin Conjugated Human IL-18 Recombinant Rabbit mAb [PSH01-70]	HA722644B	Human	ELISA(Det)
	Human IL-18 Recombinant Rabbit mAb [PSH01-70] - BSA and Azide free (Detector)	HA721723	Human	ELISA(Det)
	Human IL-18BP Recombinant Rabbit mAb [PSH10-71] - BSA and Azide free (Capture)	HA723235	Human	ELISA(Cap)
	Biotin Conjugated Human IL-18BP Recombinant Rabbit mAb [PSH10-72] - Detector	HA723237B	Human	ELISA(Det),ELISA
	Human IL-18BP Recombinant Rabbit mAb [PSH10-72] - BSA and Azide free (Detector)	HA723236	Human	ELISA(Det)
	Human IL-1 beta Recombinant Rabbit mAb [PS01-39] - BSA and Azide free (Capture)	HA721278	Human	ELISA(Cap)
	Human IL-1 beta Recombinant Rabbit mAb [PS00-82] - BSA and Azide free (Detector)	HA721279	Human	ELISA(Det)
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-86] - BSA and Azide free (Detector)	HA723882	Human	ELISA(Det)
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-85] - BSA and Azide free (Capture)	HA723881	Human	ELISA(Cap)
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-87] - BSA and Azide free (Capture)	HA723883	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human IL-1 beta Recombinant Rabbit mAb [PS00-82]	HA722649B	Human	ELISA(Det)
	Biotin Conjugated Human IL-1 beta Recombinant Rabbit mAb [PSH16-86] - Detector	HA723882B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human IL-1RA Recombinant Rabbit mAb [PSH07-04]	HA722790B	Human	ELISA(Det)
	Human IL-1RA Recombinant Rabbit mAb [PSH07-03] - BSA and Azide free (Capture)	HA722788	Human	ELISA(Cap)
	Human IL-1RA Recombinant Rabbit mAb [PSH07-04] - BSA and Azide free (Detector)	HA722789	Human	ELISA(Det)
	Human IL-2 Mouse mAb [A8B7] - BSA and Azide free (Detector)	HA601006	Human	WB,ELISA(Det)
	Biotin Conjugated Human IL-2 Mouse mAb [A8B7]	HA601361B	Human	ELISA(Det)
	Human IL-2 Recombinant Rabbit mAb [PSH01-42] - BSA and Azide free (Capture)	HA721680	Human	ELISA(Cap)
	Human IL-21 Recombinant Mouse mAb [PSH02-49] - BSA and Azide free (Detector)	HA601196	Human	ELISA(Det)
	Biotin Conjugated Human IL-21 Recombinant Mouse mAb [PSH02-49]	HA601358B	Human	ELISA(Det)
	Human IL-21 Recombinant Mouse mAb [PSH02-48] - BSA and Azide free (Capture)	HA601195	Human	ELISA(Cap)
	Biotin Conjugated Human IL-22 Recombinant Rabbit mAb [PSH03-12] - Detector	HA722630B	Human	ELISA(Det),ELISA
	Human IL-22 Recombinant Rabbit mAb [PSH03-11] - BSA and Azide free (Capture)	HA721939	Human	ELISA(Cap)
	Human IL-22 Recombinant Rabbit mAb [PSH03-12] - BSA and Azide free (Detector)	HA721940	Human	ELISA(Det)
	Human IL-23A Recombinant Rabbit mAb [PSH03-03] - BSA and Azide free (Capture)	HA721933	Human	ELISA(Cap)
	Human IL-25 Recombinant Rabbit mAb [PSH15-23] - BSA and Azide free (Capture)	HA725016	Human	ELISA(Cap)
	Human IL-25 Recombinant Rabbit mAb [PSH15-24] - BSA and Azide free (Detector)	HA725017	Human	ELISA(Det)
	Biotin Conjugated Human IL-25 Recombinant Rabbit mAb [PSH15-24] - Detector	HA725018B	Human	ELISA(Det),ELISA
	Human IL-27 Recombinant Rabbit mAb [PSH03-82] - BSA and Azide free (Capture)	HA722046	Human	ELISA(Cap)
	Biotin Conjugated Human IL-27 Recombinant Rabbit mAb [PSH03-83] - Detector	HA722621B	Human	ELISA(Det)
	Human IL-27 Recombinant Rabbit mAb [PSH03-83] - BSA and Azide free (Detector)	HA722047	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human IL-27 Recombinant Rabbit mAb [PSH03-82] - Detector	HA722622B	Human	ELISA(Det)
	Human IL-27 Recombinant Rabbit mAb [PSH03-81] - BSA and Azide free (Capture)	HA722048	Human	ELISA(Cap)
	Human IL-27 Recombinant Rabbit mAb [PSH03-82] - BSA and Azide free (Detector)	HA722049	Human	ELISA(Det)
	Biotin Conjugated Human IL-27 p28 Recombinant Rabbit mAb [PSH03-81] - Detector	HA722620B	Human	ELISA(Det)
	Human IL-27 p28 Recombinant Rabbit mAb [PSH03-81] - BSA and Azide free (Detector)	HA722045	Human	ELISA(Det)
	Human IL-27 p28 Recombinant Rabbit mAb [PSH03-80] - BSA and Azide free (Capture)	HA722044	Human	ELISA(Cap)
	Human IL-29 Recombinant Rabbit mAb [PSH15-10] - BSA and Azide free (Detector)	HA723710	Human	ELISA(Det)
	Biotin Conjugated Human IL-29 Recombinant Rabbit mAb [PSH15-10] - Detector	HA723711B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human IL-29 Recombinant Rabbit mAb [PSH15-11] - Detector	HA723713B	Human	ELISA(Det),ELISA
	Human IL-29 Recombinant Rabbit mAb [PSH15-09] - BSA and Azide free (Capture)	HA723709	Human	ELISA(Cap)
	Human IL-29 Recombinant Rabbit mAb [PSH15-11] - BSA and Azide free (Detector)	HA723712	Human	ELISA(Det)
	Human IL-2 Receptor alpha Recombinant Rabbit mAb [PS00-52] - BSA and Azide free (Capture)	HA721280	Human	ELISA(Cap)
	Human IL-2 Receptor alpha Recombinant Rabbit mAb [PSH05-66] - BSA and Azide free (Detector)	HA722386	Human	ELISA(Det)
	Biotin Conjugated Human IL-2 Receptor alpha Recombinant Rabbit mAb [PSH05-66]	HA722594B	Human	ELISA(Det)
	Human IL-2 Receptor alpha Recombinant Rabbit mAb [PSH05-65] - BSA and Azide free (Capture)	HA722385	Human	ELISA(Cap)
	Human IL-2 Receptor alpha Recombinant Rabbit mAb [PS01-03] - BSA and Azide free (Detector)	HA721281	Human	ELISA(Det)
	Biotin Conjugated Human IL-2 Receptor alpha Recombinant Rabbit mAb [PS01-03]	HA722650B	Human	ELISA(Det)
	Biotin Conjugated Human IL-3 Mouse mAb [PSH0-98]	HA601359B	Human	ELISA(Det)
	Human IL-3 Mouse mAb [PSH0-96] - BSA and Azide free (Capture)	HA601153	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human IL-3 Mouse mAb [PSH0-98] - BSA and Azide free (Detector)	HA601155	Human	ELISA(Det)
	Human IL-3 Mouse mAb [PSH0-97] - BSA and Azide free (Capture)	HA601154	Human	ELISA(Cap)
	Human IL-31 Recombinant Rabbit mAb [PSH03-02] - BSA and Azide free (Detector)	HA721932	Human	ELISA(Det)
	Human IL-31 Recombinant Rabbit mAb [PSH03-01] - BSA and Azide free (Capture)	HA721931	Human	ELISA(Cap)
	Biotin Conjugated Human IL-31 Recombinant Rabbit mAb [PSH03-02] - Detector	HA722631B	Human	ELISA(Det),ELISA
	Human IL-33 Recombinant Rabbit mAb [PSH03-06] - BSA and Azide free (Detector)	HA721936	Human	ELISA(Det)
	Human IL-33 Recombinant Rabbit mAb [PSH03-05] - BSA and Azide free (Capture)	HA721935	Human	ELISA(Cap)
	Biotin Conjugated Human IL-33 Recombinant Rabbit mAb [PSH03-06]	HA722628B	Human	ELISA(Det)
	Human IL3RA/CD123 Recombinant Rabbit mAb [PSH11-57] - BSA and Azide free (Capture)	HA723351	Human	ELISA(Cap)
	Biotin Conjugated Human IL3RA/CD123 Recombinant Rabbit mAb [PSH11-58] - Detector	HA723353B	Human	ELISA(Det),ELISA
	Human IL3RA/CD123 Recombinant Rabbit mAb [PSH11-58] - BSA and Azide free (Detector)	HA723352	Human	ELISA(Det)
	Biotin Conjugated Human IL-4 Recombinant Rabbit mAb [PS00-66] - Detector	HA721368	Human	ELISA(Det)
	Human IL-4 Recombinant Rabbit mAb [PS00-67] - BSA and Azide free (Capture)	HA721230	Human	ELISA(Cap)
	HRP Conjugated Human IL-4 Recombinant Rabbit mAb [PS00-66] - Detector	HA721295	Human	ELISA(Det)
	Human IL-4 Recombinant Rabbit mAb [PS00-66] - BSA and Azide free (Detector)	HA721231	Human	ELISA(Det)
	Biotin Conjugated Human IL-4R Recombinant Rabbit mAb [PSH11-46] - Detector	HA723338B	Human	ELISA(Det),ELISA
	Human IL-4R Recombinant Rabbit mAb [PSH11-45] - BSA and Azide free (Capture)	HA723336	Human	ELISA(Cap)
	Human IL-4R Recombinant Rabbit mAb [PSH11-46] - BSA and Azide free (Detector)	HA723337	Human	ELISA(Det)
	Human IL-5 Recombinant Rabbit mAb [PSH04-65] - BSA and Azide free (Detector)	HA722165	Human	ELISA(Det)
	Human IL-5 Recombinant Rabbit mAb [PSH04-64] - BSA and Azide free (Capture)	HA722164	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human IL-5 Recombinant Rabbit mAb [PSH04-65] - Detector	HA722609B	Human	ELISA(Det)
	Human IL-6 Recombinant Rabbit mAb [PS00-21] - BSA and Azide free (Capture)	HA721167	Human	ELISA(Cap)
	Human IL-6 Recombinant Rabbit mAb [PS00-15] - BSA and Azide free (Detector)	HA721168	Human	ELISA(Det)
	Biotin Conjugated Human IL-6 Recombinant Rabbit mAb [PS00-15] - Detector	HA722233B	Human	ELISA(Det)
	HRP Conjugated Human IL-6 Recombinant Rabbit mAb [PS00-15] - Detector	HA721293	Human	ELISA(Det)
	Human IL-6R Recombinant Rabbit mAb [PSH15-39] - BSA and Azide free (Detector)	HA723744	Human	ELISA(Det)
	HRP Conjugated Human IL-6R Recombinant Rabbit mAb [PSH15-39] - Detector	HA723745H	Human	ELISA(Det),ELISA
	Human IL-6R Recombinant Rabbit mAb [PSH15-38] - BSA and Azide free (Capture)	HA723743	Human	ELISA(Cap)
	Biotin Conjugated Human IL-6R Recombinant Rabbit mAb [PSH15-39] - Detector	HA723745B	Human	ELISA(Det),ELISA
	Human IL-7 Recombinant Rabbit mAb [PSH03-09] - BSA and Azide free (Capture)	HA721937	Human	ELISA(Cap)
	Human IL-7 Recombinant Rabbit mAb [PSH03-10] - BSA and Azide free (Detector)	HA721938	Human	ELISA(Det)
	Biotin Conjugated Human IL-7 Recombinant Rabbit mAb [PSH03-10]	HA722632B	Human	ELISA(Det)
	Human IL-8 Recombinant Rabbit mAb [PS01-22] - BSA and Azide free (Detector)	HA722099	Human	ELISA(Det)
	Human IL-8 Recombinant Rabbit mAb [PS00-83] - BSA and Azide free (Capture)	HA722098	Human	ELISA(Cap)
	Biotin Conjugated Human IL-8 Recombinant Rabbit mAb [PS01-22]	HA722611B	Human	ELISA(Det)
	Human Interferon gamma Recombinant Rabbit mAb [PSH01-45] - BSA and Azide free (Capture)	HA721683	Human	ELISA(Cap)
	Biotin Conjugated Human Interferon gamma Recombinant Rabbit mAb [PSH01-47]	HA722645B	Human	ELISA(Det)
	Human Interferon gamma Recombinant Rabbit mAb [PS01-76] - BSA and Azide free (Detector)	HA721575	Human	ELISA(Det)
	Human Interferon gamma Recombinant Rabbit mAb [PSH01-47] - BSA and Azide free (Detector)	HA721685	Human	ELISA(Det)
	Biotin Conjugated Human Interferon gamma Recombinant Rabbit mAb [PS01-76] - Detector	HA721576	Human	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human Interferon gamma Recombinant Rabbit mAb [PSH01-46] - BSA and Azide free (Capture)	HA721684	Human	ELISA(Cap)
	Human Interferon gamma Recombinant Rabbit mAb [PS01-75] - BSA and Azide free (Capture)	HA721574	Human	ELISA(Cap)
	Human LAG-3 Recombinant Rabbit mAb [PSH05-29] - BSA and Azide free (Capture)	HA722258	Human	ELISA(Cap)
	Human LAG-3 Recombinant Rabbit mAb [PSH05-30] - BSA and Azide free (Detector)	HA722259	Human	ELISA(Det)
	Biotin Conjugated Human LAG-3 Recombinant Rabbit mAb [PSH05-30]	HA722596B	Human	ELISA(Det)
	Human LAMC2 Recombinant Rabbit mAb [PSH15-14] - BSA and Azide free (Capture)	HA723717	Human	ELISA(Cap)
	Human LAMC2 Recombinant Rabbit mAb [PSH15-15] - BSA and Azide free (Detector)	HA723718	Human	ELISA(Det)
	Biotin Conjugated Human LAMC2 Recombinant Rabbit mAb [PSH15-15] - Detector	HA723719B	Human	ELISA(Det),ELISA
	Human LAP3 Recombinant Rabbit mAb [PSH01-26] - BSA and Azide free (Capture)	HA723284	Human, Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Human LAP3 Recombinant Rabbit mAb [PSH01-34] - Detector	HA723286B	Human, Mouse, Rat	ELISA(Det),ELISA
	Human LAP3 Recombinant Rabbit mAb [PSH01-34] - BSA and Azide free (Detector)	HA723285	Human, Mouse, Rat	ELISA(Det)
	Human Leptin Recombinant Rabbit mAb [PSH12-07] - BSA and Azide free (Detector)	HA723422	Human	ELISA(Det)
	Biotin Conjugated Human Leptin Recombinant Rabbit mAb [PSH12-07] - Detector	HA723423B	Human	ELISA(Det),ELISA
	Human Leptin Recombinant Rabbit mAb [PSH12-06] - BSA and Azide free (Capture)	HA723421	Human	ELISA(Cap)
	Human Leptin Receptor Recombinant Rabbit mAb [PSH14-71] - BSA and Azide free (Capture)	HA725058	Human	ELISA(Cap)
	Human Leptin Receptor Recombinant Rabbit mAb [PSH14-72] - BSA and Azide free (Detector)	HA725059	Human	ELISA(Det)
	Biotin Conjugated Human Leptin Receptor Recombinant Rabbit mAb [PSH14-72] - Detector	HA725060B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human LIF Recombinant Rabbit mAb [PSH05-97]	HA722586B	Human	ELISA(Det)
	Human LIF Recombinant Rabbit mAb [PSH05-96] - BSA and Azide free (Capture)	HA722495	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human LIF Recombinant Rabbit mAb [PSH05-97] - BSA and Azide free (Detector)	HA722496	Human	ELISA(Det)
	HRP Conjugated Human Lipocalin-2 Recombinant Rabbit mAb [PSH13-68] - Detector	HA725111H	Human	ELISA(Det),ELISA
	Human Lipocalin-2 Recombinant Rabbit mAb [PSH13-68] - BSA and Azide free (Detector)	HA725110	Human	ELISA(Det)
	Human Lipocalin-2 Recombinant Rabbit mAb [PSH13-67] - BSA and Azide free (Capture)	HA725109	Human	ELISA(Cap)
	Biotin Conjugated Human Lipocalin-2 Recombinant Rabbit mAb [PSH13-68] - Detector	HA725111B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human Lipoprotein A Recombinant Rabbit mAb [PSH13-82] - Detector	HA723572B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human Lipoprotein A Recombinant Rabbit mAb [PSH13-83] - Detector	HA723574B	Human	ELISA(Det),ELISA
	Human Lipoprotein A Recombinant Rabbit mAb [PSH13-81] - BSA and Azide free (Capture)	HA723570	Human	ELISA(Cap)
	Human Lipoprotein A Recombinant Rabbit mAb [PSH13-82] - BSA and Azide free (Detector)	HA723571	Human	ELISA(Det)
	Human Lipoprotein A Recombinant Rabbit mAb [PSH13-83] - BSA and Azide free (Detector)	HA723573	Human	ELISA(Det)
	HRP Conjugated Human Lipoprotein A Recombinant Rabbit mAb [PSH13-83] - Detector	HA723574H	Human	ELISA(Det),ELISA
	Human Livin Recombinant Rabbit mAb [PSH12-96] - BSA and Azide free (Capture)	HA723501	Human	ELISA(Cap)
	Human Livin Recombinant Rabbit mAb [PSH12-94] - BSA and Azide free (Capture)	HA723498	Human	ELISA(Cap)
	Human Livin Recombinant Rabbit mAb [PSH12-95] - BSA and Azide free (Detector)	HA723499	Human	ELISA(Det)
	Biotin Conjugated Human Livin Recombinant Rabbit mAb [PSH12-95] - Detector	HA723500B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human MBL2 Recombinant Rabbit mAb [PSH04-86]	HA722604B	Human	ELISA(Det)
	Human MBL2 Recombinant Rabbit mAb [PSH04-86] - BSA and Azide free (Detector)	HA722192	Human	ELISA(Det)
	Human MBL2 Recombinant Rabbit mAb [PSH04-85] - BSA and Azide free (Capture)	HA722191	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human MCP4/CCL13 Recombinant Rabbit mAb [PSH15-69] - BSA and Azide free (Capture)	HA725199	Human	ELISA(Cap)
	Biotin Conjugated Human MCP4/CCL13 Recombinant Rabbit mAb [PSH15-70] - Detector	HA725200B	Human	ELISA(Det),ELISA
	Human MCP4/CCL13 Recombinant Rabbit mAb [PSH15-70] - BSA and Azide free (Detector)	HA725200	Human	ELISA(Det)
	Human M-CSF Recombinant Rabbit mAb [PSH11-36] - BSA and Azide free (Capture)	HA723323	Human	ELISA(Cap)
	Human M-CSF Recombinant Rabbit mAb [PSH11-38] - BSA and Azide free (Detector)	HA723325	Human	ELISA(Det)
	Human M-CSF Recombinant Rabbit mAb [PSH11-37] - BSA and Azide free (Capture)	HA723324	Human	ELISA(Cap)
	Biotin Conjugated Human M-CSF Recombinant Rabbit mAb [PSH11-38] - Detector	HA723326B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human MERTK Recombinant Rabbit mAb [PSH12-70] - Detector	HA723481B	Human	ELISA(Det),ELISA
	Human MERTK Recombinant Rabbit mAb [PSH12-70] - BSA and Azide free (Detector)	HA723480	Human	ELISA(Det)
	HRP Conjugated Human MERTK Recombinant Rabbit mAb [PSH12-70] - Detector	HA723481H	Human	ELISA(Det),ELISA
	Human MERTK Recombinant Rabbit mAb [PSH12-69] - BSA and Azide free (Capture)	HA723479	Human	ELISA(Cap)
	HRP Conjugated Human Mesothelin Recombinant Rabbit mAb [PSH0-57] - Detector	HA723289H	Human	ELISA(Det),ELISA
	Human Mesothelin Recombinant Rabbit mAb [PSH02-30] - BSA and Azide free (Capture)	HA723287	Human	ELISA(Cap)
	Human Mesothelin Recombinant Rabbit mAb [PSH0-57] - BSA and Azide free (Detector)	HA723288	Human	ELISA(Det)
	Biotin Conjugated Human Mesothelin Recombinant Rabbit mAb [PSH0-57] - Detector	HA723289B	Human	ELISA(Det),ELISA
	Human MICA Recombinant Rabbit mAb [PSH11-39] - BSA and Azide free (Capture)	HA723327	Human	ELISA(Cap)
	Human MICA Recombinant Rabbit mAb [PSH11-40] - BSA and Azide free (Detector)	HA723328	Human	ELISA(Det)
	Biotin Conjugated Human MICA Recombinant Rabbit mAb [PSH11-40] - Detector	HA723329B	Human	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human MICB Recombinant Rabbit mAb [PSH10-74] - BSA and Azide free (Detector)	HA723239	Human	ELISA(Det)
	Biotin Conjugated Human MICB Recombinant Rabbit mAb [PSH10-74]	HA723240B	Human	ELISA(Det),ELISA
	Human MICB Recombinant Rabbit mAb [PSH10-73] - BSA and Azide free (Capture)	HA723238	Human	ELISA(Cap)
	Biotin Conjugated Human MMP-1 Recombinant Rabbit mAb [PSH07-60]	HA722876B	Human	ELISA(Det)
	Human MMP-1 Recombinant Rabbit mAb [PSH07-59] - BSA and Azide free (Capture)	HA722874	Human	ELISA(Cap)
	Human MMP-1 Recombinant Rabbit mAb [PSH07-60] - BSA and Azide free (Detector)	HA722875	Human	ELISA(Det)
	Human MMP13 Recombinant Rabbit mAb [PSH11-13] - BSA and Azide free (Capture)	HA723293	Human	ELISA(Cap)
	Human MMP13 Recombinant Rabbit mAb [PSH11-14] - BSA and Azide free (Detector)	HA723294	Human	ELISA(Det)
	Biotin Conjugated Human MMP13 Recombinant Rabbit mAb [PSH11-14] - Detector	HA723295B	Human	ELISA(Det),ELISA
	Human MMP-3 Recombinant Rabbit mAb [PSH07-62] - BSA and Azide free (Detector)	HA722878	Human	ELISA(Det)
	Human MMP-3 Recombinant Rabbit mAb [PSH07-61] - BSA and Azide free (Capture)	HA722877	Human	ELISA(Cap)
	Biotin Conjugated Human MMP-3 Recombinant Rabbit mAb [PSH07-62]	HA722879B	Human	ELISA(Det)
	Human MMP-7 Recombinant Rabbit mAb [PSH07-94] - BSA and Azide free (Detector)	HA722927	Human	ELISA(Det)
	Biotin Conjugated Human MMP-7 Recombinant Rabbit mAb [PSH07-94] - Detector	HA722928B	Human	ELISA(Det)
	Human MMP-7 Recombinant Rabbit mAb [PSH07-93] - BSA and Azide free (Capture)	HA722926	Human	ELISA(Cap)
	Human MMP-8 Recombinant Rabbit mAb [PSH03-33] - BSA and Azide free (Capture)	HA721978	Human	ELISA(Cap)
	Human MMP-8 Recombinant Rabbit mAb [PSH03-34] - BSA and Azide free (Detector)	HA721979	Human	ELISA(Det)
	Biotin Conjugated Human MMP-8 Recombinant Rabbit mAb [PSH03-34]	HA722624B	Human	ELISA(Det)
	Human MMP-9 Recombinant Rabbit mAb [PSH04-23] - BSA and Azide free (Capture)	HA722106	Human	ELISA(Cap)
	Human MMP-9 Recombinant Rabbit mAb [PSH04-63] - BSA and Azide free (Detector)	HA722163	Human	ELISA(Det)
	Biotin Conjugated Human MMP-9 Recombinant Rabbit mAb [PSH04-24]	HA722615B	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human MMP-9 Recombinant Rabbit mAb [PSH04-62] - BSA and Azide free (Capture)	HA722162	Human	ELISA(Cap)
	Human MMP-9 Recombinant Rabbit mAb [PSH04-24] - BSA and Azide free (Detector)	HA722107	Human	ELISA(Det)
	Human/Mouse/Rat NF-L Recombinant Rabbit mAb [PSH11-42] - BSA and Azide free (Detector)	HA723331	Human, Mouse, Rat	ELISA(Det)
	Human/Mouse/Rat NF-L Recombinant Rabbit mAb [PSH11-41] - BSA and Azide free (Capture)	HA723330	Human, Mouse, Rat	ELISA(Cap)
	Human/Mouse/Rat NF-L Recombinant Rabbit mAb [PSH11-44] - BSA and Azide free (Detector)	HA723334	Human, Mouse, Rat	ELISA(Det)
	Biotin Conjugated Human/Mouse/Rat NF-L Recombinant Rabbit mAb [PSH11-42] - Detector	HA723332B	Human, Mouse, Rat	ELISA(Det),ELISA
	Biotin Conjugated Human/Mouse/Rat NF-L Recombinant Rabbit mAb [PSH11-44] - Detector	HA723335B	Human, Mouse, Rat	ELISA(Det),ELISA
	Human/Mouse/Rat NF-L Recombinant Rabbit mAb [PSH11-43] - BSA and Azide free (Capture)	HA723333	Human, Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Human/Mouse VEGFC Recombinant Rabbit mAb [PSH10-76] - Detector	HA723243B	Human, Mouse	ELISA(Det),ELISA
	Human/Mouse VEGFC Recombinant Rabbit mAb [PSH10-75] - BSA and Azide free (Capture)	HA723241	Human, Mouse	ELISA(Cap)
	Human/Mouse VEGFC Recombinant Rabbit mAb [PSH10-76] - BSA and Azide free (Detector)	HA723242	Human, Mouse	ELISA(Det)
	HRP Conjugated Human Myeloperoxidase Recombinant Rabbit mAb [PSH14-79] - Detector	HA723675H	Human	ELISA(Det),ELISA
	Human Myeloperoxidase Recombinant Rabbit mAb [PSH14-78] - BSA and Azide free (Capture)	HA723673	Human	ELISA(Cap)
	Human Myeloperoxidase Recombinant Rabbit mAb [PSH14-79] - BSA and Azide free (Detector)	HA723674	Human	ELISA(Det)
	Biotin Conjugated Human Myeloperoxidase Recombinant Rabbit mAb [PSH14-79] - Detector	HA723675B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human NCAM1 Recombinant Rabbit mAb [PSH08-09] - Detector	HA722950B	Human	ELISA(Det)
	Human NCAM1 Recombinant Rabbit mAb [PSH08-08] - BSA and Azide free (Capture/Detector)	HA722947	Human	ELISA(Det),ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human NCAM1 Recombinant Rabbit mAb [PSH08-08] - Detector	HA722948B	Human	ELISA(Det)
	Human NCAM1 Recombinant Rabbit mAb [PSH08-07] - BSA and Azide free (Capture)	HA722946	Human	ELISA(Cap)
	Human NCAM1 Recombinant Rabbit mAb [PSH08-09] - BSA and Azide free (Detector)	HA722949	Human	ELISA(Det)
	Human NCR1/NKp46 Recombinant Rabbit mAb [PSH17-31] - BSA and Azide free (Capture)	HA725321	Human	ELISA(Cap)
	Human NCR1/NKp46 Recombinant Rabbit mAb [PSH17-32] - BSA and Azide free (Detector)	HA725322	Human	ELISA(Det)
	Biotin Conjugated Human NCR1/NKp46 Recombinant Rabbit mAb [PSH17-32] - Detector	HA725322B	Human	ELISA(Det),ELISA
	Human Nectin-2 Recombinant Rabbit mAb [PSH12-53] - BSA and Azide free (Capture)	HA723463	Human	ELISA(Cap)
	Biotin Conjugated Human Nectin-2 Recombinant Rabbit mAb [PSH12-52] - Detector	HA723462B	Human	ELISA(Det),ELISA
	Human Nectin-2 Recombinant Rabbit mAb [PSH12-52] - BSA and Azide free (Detector)	HA723461	Human	ELISA(Det)
	Human Nectin-2 Recombinant Rabbit mAb [PSH11-54] - BSA and Azide free (Capture)	HA723460	Human	ELISA(Cap)
	Biotin Conjugated Human Nectin-4 Recombinant Rabbit mAb [PSH12-58] - Detector	HA723472B	Human	ELISA(Det),ELISA
	Human Nectin-4 Recombinant Rabbit mAb [PSH12-58] - BSA and Azide free (Detector)	HA723471	Human	ELISA(Det)
	Human Nectin-4 Recombinant Rabbit mAb [PSH12-57] - BSA and Azide free (Capture)	HA723470	Human	ELISA(Cap)
	Human Neurogranin Recombinant Rabbit mAb [PSH16-57] - BSA and Azide free (Detector)	HA725192	Human	ELISA(Det)
	Human Neurogranin Recombinant Rabbit mAb [PSH16-56] - BSA and Azide free (Capture)	HA725191	Human	ELISA(Cap)
	Biotin Conjugated Human Neurogranin Recombinant Rabbit mAb [PSH16-57] - Detector	HA725192B	Human	ELISA(Det),ELISA
	HRP Conjugated Human Neuropilin-1 Recombinant Rabbit mAb [PSH14-74] - Detector	HA725114H	Human	ELISA(Det),ELISA
	Human Neuropilin-1 Recombinant Rabbit mAb [PSH14-74] - BSA and Azide free (Detector)	HA725113	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human Neuropilin-1 Recombinant Rabbit mAb [PSH14-73] - BSA and Azide free (Capture)	HA725112	Human	ELISA(Cap)
	Human NKG2D Recombinant Rabbit mAb [PSH17-04] - BSA and Azide free (Capture)	HA725299	Human	ELISA(Cap)
	Human NKG2D Recombinant Rabbit mAb [PSH17-05] - BSA and Azide free (Detector)	HA725300	Human	ELISA(Det)
	Biotin Conjugated Human NKG2D Recombinant Rabbit mAb [PSH17-05] - Detector	HA725300B	Human	ELISA(Det),ELISA
	Human NLRP3 Recombinant Rabbit mAb [PSH11-09] - BSA and Azide free (Detector)	HA723278	Human	ELISA(Det)
	Human NLRP3 Recombinant Rabbit mAb [PSH11-08] - BSA and Azide free (Detector)	HA723276	Human	ELISA(Det)
	Human NLRP3 Recombinant Rabbit mAb [PSH11-07] - BSA and Azide free (Capture)	HA723275	Human	ELISA(Cap)
	Biotin Conjugated Human NLRP3 Recombinant Rabbit mAb [PSH11-09] - Detector	HA723279B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human NLRP3 Recombinant Rabbit mAb [PSH11-08] - Detector	HA723277B	Human	ELISA(Det),ELISA
	Human Osteopontin Recombinant Rabbit mAb [PSH08-60] - BSA and Azide free (Detector)	HA723006	Human	ELISA(Det)
	Biotin Conjugated Human Osteopontin Recombinant Rabbit mAb [PSH08-60] - Detector	HA723008B	Human	ELISA(Det),ELISA
	Human Osteopontin Recombinant Rabbit mAb [PSH08-59] - BSA and Azide free (Detector)	HA723005	Human	ELISA(Det)
	Human Osteopontin Recombinant Rabbit mAb [PSH08-58] - BSA and Azide free (Capture)	HA723004	Human	ELISA(Cap)
	Biotin Conjugated Human Osteopontin Recombinant Rabbit mAb [PSH08-59] - Detector	HA723007B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human Osteoprotegerin Recombinant Rabbit mAb [PSH13-32] - Detector	HA725069B	Human	ELISA(Det),ELISA
	Human Osteoprotegerin Recombinant Rabbit mAb [PSH13-32] - BSA and Azide free (Detector)	HA725068	Human	ELISA(Det)
	Human Osteoprotegerin Recombinant Rabbit mAb [PSH13-31] - BSA and Azide free (Capture)	HA725067	Human	ELISA(Cap)
	Human OX40/TNFRSF4 Recombinant Rabbit mAb [PSH10-18] - BSA and Azide free (Detector)	HA723188	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human OX40/TNFRSF4 Recombinant Rabbit mAb [PSH10-17] - BSA and Azide free (Capture)	HA723187	Human	ELISA(Cap)
	Biotin Conjugated Human OX40/TNFRSF4 Recombinant Rabbit mAb [PSH10-18] - Detector	HA723189B	Human	ELISA(Det),ELISA
	Human Oxytocin-neurophysin 1 Recombinant Rabbit mAb [PSH12-29] - BSA and Azide free (Detector)	HA723453	Human	ELISA(Det)
	Biotin Conjugated Human Oxytocin-neurophysin 1 Recombinant Rabbit mAb [PSH12-29] - Detector	HA723454B	Human	ELISA(Det),ELISA
	Human Oxytocin-neurophysin 1 Recombinant Rabbit mAb [PSH12-28] - BSA and Azide free (Capture)	HA723452	Human	ELISA(Cap)
	Biotin Conjugated Human PAI1 Recombinant Rabbit mAb [PSH12-18] - Detector	HA723439B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human PAI1 Recombinant Rabbit mAb [PSH12-19] - Detector	HA723437B	Human	ELISA(Det),ELISA
	Human PAI1 Recombinant Rabbit mAb [PSH12-18] - BSA and Azide free (Capture/Detector)	HA723435	Human	ELISA(Cap),ELISA(Det)
	Human PAI1 Recombinant Rabbit mAb [PSH12-20] - BSA and Azide free (Capture)	HA723438	Human	ELISA(Cap)
	HRP Conjugated Human PAI1 Recombinant Rabbit mAb [PSH12-18] - Detector	HA723439H	Human	ELISA(Det),ELISA
	Human PAI1 Recombinant Rabbit mAb [PSH12-19] - BSA and Azide free (Detector)	HA723436	Human	ELISA(Det)
	Human PARK7 Recombinant Mouse mAb [PSH15-67] - BSA and Azide free (Capture)	HA620001	Human	ELISA(Cap)
	HRP Conjugated Human PARK7 Recombinant Mouse mAb [PSH15-68] - Detector	HA620002H	Human	ELISA(Det),ELISA
	Human PARK7 Recombinant Mouse mAb [PSH15-68] - BSA and Azide free (Detector)	HA620002	Human	ELISA(Det)
	Biotin Conjugated Human PARK7 Recombinant Mouse mAb [PSH15-68] - Detector	HA620002B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human PCSK9 Recombinant Rabbit mAb [PSH12-86] - Detector	HA725126B	Human	ELISA(Det),ELISA
	Human PCSK9 Recombinant Rabbit mAb [PSH12-86] - BSA and Azide free (Detector)	HA725125	Human	ELISA(Det)
	HRP Conjugated Human PCSK9 Recombinant Rabbit mAb [PSH12-86] - Detector	HA725126H	Human	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human PCSK9 Recombinant Rabbit mAb [PSH12-85] - BSA and Azide free (Capture)	HA725124	Human	ELISA(Cap)
	Human PD1 Recombinant Rabbit mAb [SJ01-91] - BSA and Azide free (Capture)	HA723521	Human	ELISA(Cap)
	Human PD1 Recombinant Rabbit mAb [PSH13-38] - BSA and Azide free (Detector)	HA723522	Human	ELISA(Det)
	Biotin Conjugated Human PD1 Recombinant Rabbit mAb [PSH13-38] - Detector	HA723523B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human PD1 Recombinant Rabbit mAb [PSH13-39] - Detector	HA723525B	Human	ELISA(Det),ELISA
	Human PD1 Recombinant Rabbit mAb [PSH13-39] - BSA and Azide free (Detector)	HA723524	Human	ELISA(Det)
	Human PDGF-AA Recombinant Rabbit mAb [PSH13-28] - BSA and Azide free (Detector)	HA725062	Human	ELISA(Det)
	Human PDGF-AA Recombinant Rabbit mAb [PSH13-27] - BSA and Azide free (Capture)	HA725061	Human	ELISA(Cap)
	Biotin Conjugated Human PDGF-AA Recombinant Rabbit mAb [PSH13-28] - Detector	HA725063B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human PDGF-BB Recombinant Rabbit mAb [PSH05-99] - Detector	HA722585B	Human	ELISA(Det)
	Human PDGF-BB Recombinant Rabbit mAb [PSH05-99] - BSA and Azide free (Detector)	HA722498	Human	ELISA(Det)
	Human PDGF-BB Recombinant Rabbit mAb [PSH05-98] - BSA and Azide free (Capture)	HA722497	Human	ELISA(Cap)
	Human PD-L1 Recombinant Rabbit mAb [PSH15-37] - BSA and Azide free (Detector)	HA723741	Human	ELISA(Det)
	Human PD-L1 Recombinant Rabbit mAb [PSH15-36] - BSA and Azide free (Capture)	HA723740	Human	ELISA(Cap)
	Biotin Conjugated Human PD-L1 Recombinant Rabbit mAb [PSH15-37] - Detector	HA723742B	Human	ELISA(Det),ELISA
	Human PD-L2 Recombinant Rabbit mAb [PSH09-80] - BSA and Azide free (Capture)	HA723145	Human	ELISA(Cap)
	HRP Conjugated Human PD-L2 Recombinant Rabbit mAb [PSH09-82] - Detector	HA723373H	Human	ELISA(Det)
	Human PD-L2 Recombinant Rabbit mAb [PSH09-81] - BSA and Azide free (Capture)	HA723146	Human	ELISA(Cap)
	Biotin Conjugated Human PD-L2 Recombinant Rabbit mAb [PSH09-82] - Detector	HA723148B	Human	ELISA(Det),ELISA
	Human PD-L2 Recombinant Rabbit mAb [PSH09-82] - BSA and Azide free (Detector)	HA723147	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human Pentraxin 3 / PTX3 Recombinant Rabbit mAb [PSH11-85] - Detector	HA723391B	Human	ELISA(Det),ELISA
	Human Pentraxin 3 / PTX3 Recombinant Rabbit mAb [PSH11-85] - BSA and Azide free (Detector)	HA723390	Human	ELISA(Det)
	Human Pentraxin 3 / PTX3 Recombinant Rabbit mAb [PSH11-84] - BSA and Azide free (Detector)	HA723388	Human	ELISA(Det)
	Biotin Conjugated Human Pentraxin 3 / PTX3 Recombinant Rabbit mAb [PSH11-84] - Detector	HA723389B	Human	ELISA(Det),ELISA
	Human Pentraxin 3 / PTX3 Recombinant Rabbit mAb [PSH11-83] - BSA and Azide free (Capture)	HA723387	Human	ELISA(Cap)
	HRP Conjugated Human Poliovirus Receptor/PVR Recombinant Rabbit mAb [PSH12-55] - Detector	HA723466H	Human	ELISA(Det),ELISA
	Human Poliovirus Receptor/PVR Recombinant Rabbit mAb [PSH12-54] - BSA and Azide free (Capture)	HA723464	Human	ELISA(Cap)
	Biotin Conjugated Human Poliovirus Receptor/PVR Recombinant Rabbit mAb [PSH12-55] - Detector	HA723466B	Human	ELISA(Det),ELISA
	Human Poliovirus Receptor/PVR Recombinant Rabbit mAb [PSH12-55] - BSA and Azide free (Detector)	HA723465	Human	ELISA(Det)
	Human Prealbumin Recombinant Rabbit mAb [PSH17-42] - BSA and Azide free (Capture)	HA725259	Human	ELISA(Cap)
	Biotin Conjugated Human Prealbumin Recombinant Rabbit mAb [PSH17-43] - Detector	HA725260B	Human	ELISA(Det),ELISA
	Human Prealbumin Recombinant Rabbit mAb [PSH17-43] - BSA and Azide free (Detector)	HA725260	Human	ELISA(Det)
	HRP Conjugated Human Prealbumin Recombinant Rabbit mAb [PSH17-43] - Detector	HA725260H	Human	ELISA(Det),ELISA
	Biotin Conjugated Human Prolactin Recombinant Rabbit mAb [PSH03-14] - Detector	HA722627B	Human	ELISA(Det)
	Human Prolactin Recombinant Rabbit mAb [PSH03-14] - BSA and Azide free (Detector)	HA721942	Human	ELISA(Det)
	Human Prolactin Recombinant Rabbit mAb [PSH03-13] - BSA and Azide free (Capture)	HA721941	Human	ELISA(Cap)
	Biotin Conjugated Human RANKL Recombinant Rabbit mAb [PSH16-59] - Detector	HA723861B	Human	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human RANKL Recombinant Rabbit mAb [PSH16-59] - BSA and Azide free (Detector)	HA723861	Human	ELISA(Det)
	Human RANKL Recombinant Rabbit mAb [PSH16-60] - BSA and Azide free (Capture)	HA723862	Human	ELISA(Cap)
	Human RANKL Recombinant Rabbit mAb [PSH16-58] - BSA and Azide free (Capture)	HA723860	Human	ELISA(Cap)
	HRP Conjugated Human RBP4 Recombinant Rabbit mAb [PSH14-89] - Detector	HA723690H	Human	ELISA(Det),ELISA
	Human RBP4 Recombinant Rabbit mAb [PSH14-89] - BSA and Azide free (Detector)	HA723689	Human	ELISA(Det)
	Human RBP4 Recombinant Rabbit mAb [PSH14-88] - BSA and Azide free (Capture)	HA723688	Human	ELISA(Cap)
	Biotin Conjugated Human RBP4 Recombinant Rabbit mAb [PSH14-89] - Detector	HA723690B	Human	ELISA(Det),ELISA
	Human Renin Recombinant Rabbit mAb [PSH17-00] - BSA and Azide free (Capture)	HA725169	Human	ELISA(Cap)
	Biotin Conjugated Human Renin Recombinant Rabbit mAb [PSH17-01] - Detector	HA725170B	Human	ELISA(Det),ELISA
	Human Renin Recombinant Rabbit mAb [PSH17-01] - BSA and Azide free (Detector)	HA725170	Human	ELISA(Det)
	Human Resistin Recombinant Rabbit mAb [PSH14-61] - BSA and Azide free (Capture)	HA723658	Human	ELISA(Cap)
	Human Resistin Recombinant Rabbit mAb [PSH14-63] - BSA and Azide free (Capture)	HA723661	Human	ELISA(Cap)
	Biotin Conjugated Human Resistin Recombinant Rabbit mAb [PSH14-62] - Detector	HA723660B	Human	ELISA(Det),ELISA
	Human Resistin Recombinant Rabbit mAb [PSH14-62] - BSA and Azide free (Detector)	HA723659	Human	ELISA(Det)
	Biotin Conjugated Human S100A12 Recombinant Rabbit mAb [PSH17-68] - Detector	HA725216B	Human	ELISA(Det),ELISA
	Human S100A12 Recombinant Rabbit mAb [PSH17-67] - BSA and Azide free (Capture)	HA725215	Human	ELISA(Cap)
	Human S100A12 Recombinant Rabbit mAb [PSH17-68] - BSA and Azide free (Detector)	HA725216	Human	ELISA(Det)
	HRP Conjugated Human S100A12 Recombinant Rabbit mAb [PSH17-68] - Detector	HA725216H	Human	ELISA(Det),ELISA
	Human SAA1 Mouse mAb [A10E7] - BSA and Azide free (Capture)	HA601333	Human	ELISA(Cap)
	Human SAA1 Mouse mAb [A10E6] - BSA and Azide free (Detector)	HA601334	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human SAA1 Mouse mAb [A10E6]	HA601356B	Human	ELISA(Det)
	Human SCF Recombinant Rabbit mAb [PSH07-57] - BSA and Azide free (Capture)	HA722871	Human	ELISA(Cap)
	Human SCF Recombinant Rabbit mAb [PSH07-58] - BSA and Azide free (Detector)	HA722872	Human	ELISA(Det)
	Biotin Conjugated Human SCF Recombinant Rabbit mAb [PSH07-58]	HA722873B	Human	ELISA(Det)
	Human Siglec 7 Recombinant Rabbit mAb [PSH11-87] - BSA and Azide free (Capture/Detector)	HA723393	Human	ELISA(Cap),ELISA(Det)
	Human Siglec 7 Recombinant Rabbit mAb [PSH11-88] - BSA and Azide free (Detector)	HA723395	Human	ELISA(Det)
	Biotin Conjugated Human Siglec 7 Recombinant Rabbit mAb [PSH11-88] - Detector	HA723396B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human Siglec 7 Recombinant Rabbit mAb [PSH11-87] - Detector	HA723394B	Human	ELISA(Det),ELISA
	Human Siglec 7 Recombinant Rabbit mAb [PSH11-86] - BSA and Azide free (Capture)	HA723392	Human	ELISA(Cap)
	Biotin Conjugated Human Siglec 9 Recombinant Rabbit mAb [PSH11-91] - Detector	HA723401B	Human	ELISA(Det),ELISA
	Human Siglec 9 Recombinant Rabbit mAb [PSH11-89] - BSA and Azide free (Capture)	HA723397	Human	ELISA(Cap)
	Human Siglec 9 Recombinant Rabbit mAb [PSH11-91] - BSA and Azide free (Detector)	HA723400	Human	ELISA(Det)
	Biotin Conjugated Human Siglec 9 Recombinant Rabbit mAb [PSH11-90] - Detector	HA723399B	Human	ELISA(Det),ELISA
	Human Siglec 9 Recombinant Rabbit mAb [PSH11-90] - BSA and Azide free (Detector)	HA723398	Human	ELISA(Det)
	Biotin Conjugated Human SIRP alpha Recombinant Rabbit mAb [PSH17-64] - Detector	HA725348B	Human	ELISA(Det),ELISA
	Human SIRP alpha Recombinant Rabbit mAb [PSH17-64] - BSA and Azide free (Detector)	HA725348	Human	ELISA(Det)
	Human SIRP alpha Recombinant Rabbit mAb [PSH17-63] - BSA and Azide free (Capture)	HA725347	Human	ELISA(Cap)
	HRP Conjugated Human sTNF RI Recombinant Rabbit mAb [PSH09-90] - Detector	HA723374H	Human	ELISA(Det)
	Human sTNF RI Recombinant Rabbit mAb [PSH09-89] - BSA and Azide free (Capture)	HA723158	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human sTNF RI Recombinant Rabbit mAb [PSH09-90] - BSA and Azide free (Detector)	HA723160B	Human	ELISA(Det),ELISA
	Human sTNF RI Recombinant Rabbit mAb [PSH09-90] - BSA and Azide free (Detector)	HA723159	Human	ELISA(Det)
	Human Thrombopoietin Recombinant Rabbit mAb [PSH13-97] - BSA and Azide free (Detector)	HA725119	Human	ELISA(Det)
	Human Thrombopoietin Recombinant Rabbit mAb [PSH13-96] - BSA and Azide free (Capture)	HA725118	Human	ELISA(Cap)
	Biotin Conjugated Human Thrombopoietin Recombinant Rabbit mAb [PSH13-97] - Detector	HA725120B	Human	ELISA(Det),ELISA
	HRP Conjugated Human TIM 4 Recombinant Rabbit mAb [PSH14-01] - Detector	HA723591H	Human	ELISA(Det),ELISA
	Human TIM 4 Recombinant Rabbit mAb [PSH14-00] - BSA and Azide free (Capture)	HA723589	Human	ELISA(Cap)
	Biotin Conjugated Human TIM 4 Recombinant Rabbit mAb [PSH14-01] - Detector	HA723591B	Human	ELISA(Det),ELISA
	Human TIM 4 Recombinant Rabbit mAb [PSH14-01] - BSA and Azide free (Detector)	HA723590	Human	ELISA(Det)
	Biotin Conjugated Human TNF alpha Recombinant Rabbit mAb [PS01-07] - Detector	HA721559	Human	ELISA(Det)
	Human TNF alpha Recombinant Rabbit mAb [PS01-07] - BSA and Azide free (Detector)	HA721266	Human	ELISA(Det)
	HRP Conjugated Human TNF alpha Recombinant Rabbit mAb [PS01-07] - Detector	HA721294	Human	ELISA(Det)
	Human TNF alpha Recombinant Rabbit mAb [PS00-40] - BSA and Azide free (Capture)	HA721265	Human	ELISA(Cap)
	Human TNF-beta Recombinant Rabbit mAb [PS01-05] - BSA and Azide free (Capture)	HA721242	Human	ELISA(Cap)
	Biotin Conjugated Human TNF-beta Recombinant Rabbit mAb [PS01-19]	HA722651B	Human	ELISA(Det)
	Human TNF-beta Recombinant Rabbit mAb [PS01-19] - BSA and Azide free (Detector)	HA721243	Human	ELISA(Det)
	Biotin Conjugated Human TNF Receptor II Recombinant Rabbit mAb [PSH11-93] - Detector	HA723404B	Human	ELISA(Det),ELISA
	Human TNF Receptor II Recombinant Rabbit mAb [PSH11-92] - BSA and Azide free (Capture)	HA723402	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human TNF Receptor II Recombinant Rabbit mAb [PSH11-93] - BSA and Azide free (Detector)	HA723403	Human	ELISA(Det)
	Human TNF Receptor II Recombinant Rabbit mAb [PSH11-94] - BSA and Azide free (Capture)	HA723405	Human	ELISA(Cap)
	Biotin Conjugated Human tPA Recombinant Rabbit mAb [PSH13-11] - Detector	HA725117B	Human	ELISA(Det),ELISA
	Human tPA Recombinant Rabbit mAb [PSH13-10] - BSA and Azide free (Capture)	HA725115	Human	ELISA(Cap)
	Human tPA Recombinant Rabbit mAb [PSH13-11] - BSA and Azide free (Detector)	HA725116	Human	ELISA(Det)
	HRP Conjugated Human tPA Recombinant Rabbit mAb [PSH13-11] - Detector	HA725117H	Human	ELISA(Det),ELISA
	Human TRAF1 Recombinant Rabbit mAb [PSH14-28] - BSA and Azide free (Detector)	HA723619	Human	ELISA(Det)
	Biotin Conjugated Human TRAF1 Recombinant Rabbit mAb [PSH14-28] - Detector	HA723620B	Human	ELISA(Det),ELISA
	Human TRAF1 Recombinant Rabbit mAb [PSH14-27] - BSA and Azide free (Capture)	HA723618	Human	ELISA(Cap)
	Biotin Conjugated Human TRAIL Recombinant Rabbit mAb [PSH08-06]	HA722945B	Human	ELISA(Det)
	Human TRAIL Recombinant Rabbit mAb [PSH08-05] - BSA and Azide free (Capture)	HA722943	Human	ELISA(Cap)
	Human TRAIL Recombinant Rabbit mAb [PSH08-06] - BSA and Azide free (Detector)	HA722944	Human	ELISA(Det)
	Human TREM1 Recombinant Rabbit mAb [PSH17-85] - BSA and Azide free (Capture)	HA725233	Human	ELISA(Cap)
	Biotin Conjugated Human TREM1 Recombinant Rabbit mAb [PSH17-86] - Detector	HA725234B	Human	ELISA(Det),ELISA
	Human TREM1 Recombinant Rabbit mAb [PSH17-86] - BSA and Azide free (Detector)	HA725234	Human	ELISA(Det)
	Biotin Conjugated Human TREM2 Recombinant Rabbit mAb [PSH09-85] - Detector	HA723152B	Human	ELISA(Det),ELISA
	HRP Conjugated Human TREM2 Recombinant Rabbit mAb [PSH09-83] - Detector	HA723372	Human	ELISA(Det)
	Human TREM2 Recombinant Rabbit mAb [PSH09-84] - BSA and Azide free (Capture)	HA723150	Human	ELISA(Cap)
	Human TREM2 Recombinant Rabbit mAb [PSH09-83] - BSA and Azide free (Capture/Detector)	HA723149	Human	ELISA(Cap),ELISA(Det)
	Human TREM2 Recombinant Rabbit mAb [PSH09-85] - BSA and Azide free (Detector)	HA723151	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human TREM2 Recombinant Rabbit mAb [PSH09-83] - Detector	HA723153B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human TSLP Recombinant Rabbit mAb [PSH12-90] - Detector	HA723492B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human TSLP Recombinant Rabbit mAb [PSH12-88] - Detector	HA723489B	Human	ELISA(Det),ELISA
	Human TSLP Recombinant Rabbit mAb [PSH12-89] - BSA and Azide free (Capture)	HA723490	Human	ELISA(Cap)
	Human TSLP Recombinant Rabbit mAb [PSH12-87] - BSA and Azide free (Capture)	HA723487	Human	ELISA(Cap)
	Human TSLP Recombinant Rabbit mAb [PSH12-90] - BSA and Azide free (Detector)	HA723491	Human	ELISA(Det)
	Human TSLP Recombinant Rabbit mAb [PSH12-88] - BSA and Azide free (Detector)	HA723488	Human	ELISA(Det)
	Human ULBP1 Recombinant Rabbit mAb [PSH11-95] - BSA and Azide free (Capture)	HA723406	Human	ELISA(Cap)
	Human ULBP1 Recombinant Rabbit mAb [PSH11-97] - BSA and Azide free (Detector)	HA723409	Human	ELISA(Det)
	Biotin Conjugated Human ULBP1 Recombinant Rabbit mAb [PSH11-97] - Detector	HA723410B	Human	ELISA(Det),ELISA
	Human ULBP1 Recombinant Rabbit mAb [PSH11-96] - BSA and Azide free (Detector)	HA723407	Human	ELISA(Det)
	Biotin Conjugated Human ULBP1 Recombinant Rabbit mAb [PSH11-96] - Detector	HA723408B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human uPAR Mouse mAb [PSH13-77] - Detector	HA601439B	Human	ELISA(Det),ELISA
	Human uPAR Mouse mAb [PSH13-77] - BSA and Azide free (Detector)	HA601438	Human	ELISA(Det)
	Human uPAR Mouse mAb [PSH13-76] - BSA and Azide free (Capture)	HA601437	Human	ELISA(Cap)
	Human VCAM-1 Recombinant Rabbit mAb [PSH12-46] - BSA and Azide free (Capture)	HA725064	Human	ELISA(Cap)
	Human VCAM-1 Recombinant Rabbit mAb [PSH12-47] - BSA and Azide free (Detector)	HA725065	Human	ELISA(Det)
	HRP Conjugated Human VCAM-1 Recombinant Rabbit mAb [PSH12-47] - Detector	HA725066H	Human	ELISA(Det),ELISA
	Biotin Conjugated Human VCAM-1 Recombinant Rabbit mAb [PSH12-47] - Detector	HA725066B	Human	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human VE Cadherin Recombinant Rabbit mAb [PSH02-62] - Detector	HA722635B	Human	ELISA(Det)
	Human VE Cadherin Recombinant Rabbit mAb [PSH02-61] - BSA and Azide free (Capture)	HA721838	Human	ELISA(Cap)
	Human VE Cadherin Recombinant Rabbit mAb [PSH02-62] - BSA and Azide free (Detector)	HA721839	Human	ELISA(Det)
	Human VEGFA Recombinant Rabbit mAb [PSH12-92] - BSA and Azide free (Detector)	HA723494	Human	ELISA(Det)
	Human VEGFA Recombinant Rabbit mAb [PSH12-93] - BSA and Azide free (Detector)	HA723496	Human	ELISA(Det)
	Biotin Conjugated Human VEGFA Recombinant Rabbit mAb [PSH12-92] - Detector	HA723495B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human VEGFA Recombinant Rabbit mAb [PSH12-93] - Detector	HA723497B	Human	ELISA(Det),ELISA
	Human VEGFA Recombinant Rabbit mAb [PSH12-91] - BSA and Azide free (Capture)	HA723493	Human	ELISA(Cap)
	Biotin Conjugated Human VEGF Receptor 2 Recombinant Rabbit mAb [PSH12-84] - Detector	HA725123B	Human	ELISA(Det),ELISA
	Human VEGF Receptor 2 Recombinant Rabbit mAb [PSH12-83] - BSA and Azide free (Capture)	HA725121	Human	ELISA(Cap)
	Human VEGF Receptor 2 Recombinant Rabbit mAb [PSH12-84] - BSA and Azide free (Detector)	HA725122	Human	ELISA(Det)
	Human Von Willebrand Factor Recombinant Rabbit mAb [PSH09-03] - BSA and Azide free (Capture)	HA723058	Human	ELISA(Cap)
	Human Von Willebrand Factor Recombinant Rabbit mAb [PSH09-04] - BSA and Azide free (Detector)	HA723059	Human	ELISA(Det)
	HRP Conjugated Human Von Willebrand Factor Recombinant Rabbit mAb [PSH09-04] - Detector	HA723431H	Human	ELISA(Det)
	Biotin Conjugated Human Von Willebrand Factor Recombinant Rabbit mAb [PSH09-04] - BSA and Azide free (Detector)	HA723060B	Human	ELISA(Det),ELISA
	Human YKL-40 Recombinant Rabbit mAb [PSH13-08] - BSA and Azide free (Capture)	HA725088	Human	ELISA(Cap)
	HRP Conjugated Human YKL-40 Recombinant Rabbit mAb [PSH13-09] - Detector	HA725090H	Human	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human YKL-40 Recombinant Rabbit mAb [PSH13-09] - Detector	HA725090B	Human	ELISA(Det),ELISA
	Human YKL-40 Recombinant Rabbit mAb [PSH13-09] - BSA and Azide free (Detector)	HA725089	Human	ELISA(Det)
	Insulin Recombinant Rabbit mAb [PSH04-15] - BSA and Azide free (Capture)	HA722096	Human, Mouse, Rat	ELISA(Cap)
	Insulin Recombinant Rabbit mAb [PSH04-16] - BSA and Azide free (Detector)	HA722097	Human, Mouse, Rat	ELISA(Det)
	Biotin Conjugated Insulin Recombinant Rabbit mAb [PSH04-16] - Detector	HA722610B	Human, Mouse, Rat	ELISA(Det),ELISA
	Mouse Adiponectin Recombinant Rabbit mAb [PSH14-84] - BSA and Azide free (Capture)	HA723682	Mouse	ELISA(Cap)
	Mouse Adiponectin Recombinant Rabbit mAb [PSH14-83] - BSA and Azide free (Detector)	HA723680	Mouse	ELISA(Det)
	Biotin Conjugated Mouse Adiponectin Recombinant Rabbit mAb [PSH14-83] - Detector	HA723681B	Mouse	ELISA(Det),ELISA
	Mouse Adiponectin Recombinant Rabbit mAb [PSH14-82] - BSA and Azide free (Capture)	HA723679	Mouse	ELISA(Cap)
	Mouse Alpha 1 antitrypsin Recombinant Rabbit mAb [PSH16-20] - BSA and Azide free (Capture)	HA725146	Mouse	ELISA(Cap)
	Mouse Alpha 1 antitrypsin Recombinant Rabbit mAb [PSH16-21] - BSA and Azide free (Detector)	HA725147	Mouse	ELISA(Det)
	Biotin Conjugated Mouse Alpha 1 antitrypsin Recombinant Rabbit mAb [PSH16-21] - Detector	HA725148B	Mouse	ELISA(Det),ELISA
	HRP Conjugated Mouse Alpha 1 antitrypsin Recombinant Rabbit mAb [PSH16-21] - Detector	HA725148H	Mouse	ELISA(Det),ELISA
	Biotin Conjugated Mouse Angiotensinogen Recombinant Rabbit mAb [PSH17-10] - Detector	HA725133B	Mouse	ELISA(Det),ELISA
	HRP Conjugated Mouse Angiotensinogen Recombinant Rabbit mAb [PSH17-10] - Detector	HA725133H	Mouse	ELISA(Det),ELISA
	Mouse Angiotensinogen/AGT Recombinant Rabbit mAb [PSH17-09] - BSA and Azide free (Capture)	HA725131	Mouse	ELISA(Cap)
	Mouse Angiotensinogen/AGT Recombinant Rabbit mAb [PSH17-10] - BSA and Azide free (Detector)	HA725132	Mouse	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse Apolipoprotein B/Apo B Recombinant Rabbit mAb [PSH17-77] - BSA and Azide free (Capture)	HA723946	Mouse	ELISA(Cap)
	Mouse Apolipoprotein B/Apo B Recombinant Rabbit mAb [PSH17-76] - BSA and Azide free (Detector)	HA723945	Mouse	ELISA(Det)
	Biotin Conjugated Mouse Apolipoprotein B/Apo B Recombinant Rabbit mAb [PSH17-76] - Detector	HA723945B	Mouse	ELISA(Det),ELISA
	Mouse Apolipoprotein B/Apo B Recombinant Rabbit mAb [PSH17-75] - BSA and Azide free (Capture)	HA723944	Mouse	ELISA(Cap)
	Mouse Apolipoprotein E Recombinant Rabbit mAb [PSH16-84] - BSA and Azide free (Detector)	HA725222	Mouse	ELISA(Det)
	Biotin Conjugated Mouse Apolipoprotein E Recombinant Rabbit mAb [PSH16-84] - Detector	HA725222B	Mouse	ELISA(Det),ELISA
	Mouse Apolipoprotein E Recombinant Rabbit mAb [PSH16-83] - BSA and Azide free (Capture)	HA725221	Mouse	ELISA(Cap)
	Mouse CCL2 / MCP-1 Recombinant Rabbit mAb [PSH12-45] - BSA and Azide free (Detector)	HA725041	Mouse	ELISA(Det)
	Mouse CCL2 / MCP-1 Recombinant Rabbit mAb [PSH12-44] - BSA and Azide free (Capture)	HA725040	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse CCL2 / MCP-1 Recombinant Rabbit mAb [PSH12-45] - Detector	HA725042B	Mouse	ELISA(Det),ELISA
	Biotin Conjugated Mouse CCL3 Recombinant Rabbit mAb [PSH08-40] - Detector	HA722984B	Mouse	ELISA(Det),ELISA
	Mouse CCL3 Recombinant Rabbit mAb [PSH08-39] - BSA and Azide free (Capture)	HA722982	Mouse	ELISA(Cap)
	Mouse CCL3 Recombinant Rabbit mAb [PSH08-40] - BSA and Azide free (Detector)	HA722983	Mouse	ELISA(Det)
	Mouse CCL5/RANTES Recombinant Rabbit mAb [PSH14-15] - BSA and Azide free (Capture)	HA725046	Mouse	ELISA(Cap)
	Mouse CCL5/RANTES Recombinant Rabbit mAb [PSH14-16] - BSA and Azide free (Detector)	HA725047	Mouse	ELISA(Det)
	Biotin Conjugated Mouse CCL5/RANTES Recombinant Rabbit mAb [PSH14-16] - Detector	HA725048B	Mouse	ELISA(Det),ELISA
	Mouse CD30 Recombinant Rabbit mAb [PSH17-55] - BSA and Azide free (Detector)	HA723930	Mouse	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Mouse CD30 Recombinant Rabbit mAb [PSH17-55] - Detector	HA723930B	Mouse	ELISA(Det),ELISA
	Mouse CD30 Recombinant Rabbit mAb [PSH17-54] - BSA and Azide free (Capture)	HA723929	Mouse	ELISA(Cap)
	Mouse CD68 Recombinant Rabbit mAb [PSH12-24] - BSA and Azide free (Capture)	HA723446	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse CD68 Recombinant Rabbit mAb [PSH12-27] - Detector	HA723451B	Mouse	ELISA(Det),ELISA
	Mouse CD68 Recombinant Rabbit mAb [PSH12-27] - BSA and Azide free (Detector)	HA723450	Mouse	ELISA(Det)
	Biotin Conjugated Mouse CD68 Recombinant Rabbit mAb [PSH12-25] - Detector	HA723448B	Mouse	ELISA(Det),ELISA
	Mouse CD68 Recombinant Rabbit mAb [PSH12-25] - BSA and Azide free (Detector)	HA723447	Mouse	ELISA(Det)
	Mouse CD68 Recombinant Rabbit mAb [PSH12-26] - BSA and Azide free (Capture)	HA723449	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse CD86 Recombinant Rabbit mAb [PSH13-64] - Detector	HA723554B	Mouse	ELISA(Det),ELISA
	Mouse CD86 Recombinant Rabbit mAb [PSH13-63] - BSA and Azide free (Capture)	HA723552	Mouse	ELISA(Cap)
	Mouse CD86 Recombinant Rabbit mAb [PSH13-64] - BSA and Azide free (Detector)	HA723553	Mouse	ELISA(Det)
	Mouse CD86 Recombinant Rabbit mAb [PSH13-66] - BSA and Azide free (Detector)	HA723556	Mouse	ELISA(Det)
	Biotin Conjugated Mouse CD86 Recombinant Rabbit mAb [PSH13-66] - Detector	HA723557B	Mouse	ELISA(Det),ELISA
	Mouse CD86 Recombinant Rabbit mAb [PSH13-65] - BSA and Azide free (Capture)	HA723555	Mouse	ELISA(Cap)
	Mouse CRP Recombinant Rabbit mAb [PSH12-49] - BSA and Azide free (Detector)	HA725074	Mouse	ELISA(Det)
	Biotin Conjugated Mouse CRP Recombinant Rabbit mAb [PSH12-49] - Detector	HA725075B	Mouse	ELISA(Det),ELISA
	Mouse CRP Recombinant Rabbit mAb [PSH12-48] - BSA and Azide free (Capture)	HA725073	Mouse	ELISA(Cap)
	HRP Conjugated Mouse CRP Recombinant Rabbit mAb [PSH12-49] - Detector	HA725075H	Mouse	ELISA(Det),ELISA
	Mouse CXCL16 Recombinant Rabbit mAb [PSH17-73] - BSA and Azide free (Capture)	HA723942	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse CXCL16 Recombinant Rabbit mAb [PSH17-74] - Detector	HA723943B	Mouse	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse CXCL16 Recombinant Rabbit mAb [PSH17-71] - BSA and Azide free (Capture)	HA723940	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse CXCL16 Recombinant Rabbit mAb [PSH17-72] - Detector	HA723941B	Mouse	ELISA(Det),ELISA
	Mouse CXCL16 Recombinant Rabbit mAb [PSH17-74] - BSA and Azide free (Detector)	HA723943	Mouse	ELISA(Det)
	Mouse CXCL16 Recombinant Rabbit mAb [PSH17-72] - BSA and Azide free (Detector)	HA723941	Mouse	ELISA(Det)
	Biotin Conjugated Mouse CXCL2 Recombinant Rabbit mAb [PSH12-61] - Detector	HA725072B	Mouse	ELISA(Det),ELISA
	Mouse CXCL2 Recombinant Rabbit mAb [PSH12-60] - BSA and Azide free (Capture)	HA725070	Mouse	ELISA(Cap)
	Mouse CXCL2 Recombinant Rabbit mAb [PSH12-61] - BSA and Azide free (Detector)	HA725071	Mouse	ELISA(Det)
	Mouse CXCL9 Recombinant Rabbit mAb [PSH12-62] - BSA and Azide free (Capture)	HA725022	Mouse	ELISA(Cap)
	Mouse CXCL9 Recombinant Rabbit mAb [PSH12-63] - BSA and Azide free (Detector)	HA725023	Mouse	ELISA(Det)
	Biotin Conjugated Mouse CXCL9 Recombinant Rabbit mAb [PSH12-63]	HA725024B	Mouse	ELISA(Det),ELISA
	Biotin Conjugated Mouse EPO/Erythropoietin Recombinant Rabbit mAb [PSH17-22] - Detector	HA723909B	Mouse	ELISA(Det)
	Mouse EPO/Erythropoietin Recombinant Rabbit mAb [PSH17-21] - BSA and Azide free (Capture)	HA723908	Mouse	ELISA(Cap)
	Mouse EPO/Erythropoietin Recombinant Rabbit mAb [PSH17-23] - BSA and Azide free (Capture)	HA723910	Mouse	ELISA(Cap)
	Mouse EPO/Erythropoietin Recombinant Rabbit mAb [PSH17-22] - BSA and Azide free (Detector)	HA723909	Mouse	ELISA(Det)
	Mouse Factor XII Recombinant Rabbit mAb [PSH12-21] - BSA and Azide free (Capture)	HA723440	Mouse	ELISA(Cap)
	HRP Conjugated Mouse Factor XII Recombinant Rabbit mAb [PSH12-22] - Detector	HA723509H	Mouse	ELISA(Det),ELISA
	Biotin Conjugated Mouse Factor XII Recombinant Rabbit mAb [PSH12-22] - Detector	HA723442B	Mouse	ELISA(Det),ELISA
	Mouse Factor XII Recombinant Rabbit mAb [PSH12-22] - BSA and Azide free (Detector)	HA723441	Mouse	ELISA(Det)
	Mouse FGF-21 Recombinant Rabbit mAb [PSH15-27] - BSA and Azide free (Detector)	HA723732	Mouse	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Mouse FGF-21 Recombinant Rabbit mAb [PSH15-27] - Detector	HA723733B	Mouse	ELISA(Det),ELISA
	Mouse FGF-21 Recombinant Rabbit mAb [PSH15-26] - BSA and Azide free (Capture)	HA723731	Mouse	ELISA(Cap)
	Mouse G-CSF Recombinant Rabbit mAb [PSH15-74] - BSA and Azide free (Capture)	HA723777	Mouse	ELISA(Cap)
	Mouse G-CSF Recombinant Rabbit mAb [PSH15-75] - BSA and Azide free (Detector)	HA723778	Mouse	ELISA(Det)
	Biotin Conjugated Mouse G-CSF Recombinant Rabbit mAb [PSH15-75] - Detector	HA723779B	Mouse	ELISA(Det),ELISA
	Mouse GM-CSF Recombinant Rabbit mAb [PSH12-40] - BSA and Azide free (Capture)	HA725019	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse GM-CSF Recombinant Rabbit mAb [PSH12-41] - Detector	HA725021B	Mouse	ELISA(Det),ELISA
	Mouse GM-CSF Recombinant Rabbit mAb [PSH12-41] - BSA and Azide free (Detector)	HA725020	Mouse	ELISA(Det)
	Mouse ICAM1 Recombinant Rabbit mAb [PSH16-36] - BSA and Azide free (Capture)	HA723834	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse ICAM1 Recombinant Rabbit mAb [PSH16-35] - Detector	HA723833B	Mouse	ELISA(Det),ELISA
	HRP Conjugated Mouse ICAM1 Recombinant Rabbit mAb [PSH16-35] - Detector	HA723833H	Mouse	ELISA(Det),ELISA
	Mouse ICAM1 Recombinant Rabbit mAb [PSH16-35] - BSA and Azide free (Detector)	HA723833	Mouse	ELISA(Det)
	Mouse ICAM1 Recombinant Rabbit mAb [PSH16-34] - BSA and Azide free (Capture)	HA723832	Mouse	ELISA(Cap)
	Mouse IFN beta Recombinant Rabbit mAb [PSH04-72] - BSA and Azide free (Capture)	HA722172	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IFN beta Recombinant Rabbit mAb [PSH04-73]	HA722605B	Mouse	ELISA(Det),ELISA
	Mouse IFN beta Recombinant Rabbit mAb [PSH04-73] - BSA and Azide free (Detector)	HA722173	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IFN gamma Recombinant Rabbit mAb [PSH14-81] - Detector	HA723678B	Mouse	ELISA(Det),ELISA
	Mouse IFN gamma Recombinant Rabbit mAb [PSH14-80] - BSA and Azide free (Capture)	HA723676	Mouse	ELISA(Cap)
	Mouse IFN gamma Recombinant Rabbit mAb [PSH03-20] - BSA and Azide free (Capture)	HA721956	Mouse	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse IFN gamma Recombinant Rabbit mAb [PSH14-81] - BSA and Azide free (Detector)	HA723677	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-10 Recombinant Rabbit mAb [PS01-93] - Detector	HA721519	Mouse	ELISA(Det)
	Mouse IL-10 Recombinant Rabbit mAb [PS01-93] - BSA and Azide free (Detector)	HA721518	Mouse	ELISA(Det)
	Mouse IL-10 Recombinant Rabbit mAb [PS01-85] - BSA and Azide free (Capture)	HA721517	Mouse	ELISA(Cap)
	Mouse IL12 p40 Recombinant Rabbit mAb [PSH04-87] - BSA and Azide free (Capture)	HA722193	Mouse	ELISA(Cap)
	Mouse IL12 p40 Recombinant Rabbit mAb [PSH04-88] - BSA and Azide free (Detector)	HA722194	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL12 p40 Recombinant Rabbit mAb [PSH04-88]	HA722603B	Mouse	ELISA(Det)
	Mouse IL-13 Recombinant Rabbit mAb [PSH02-05] - BSA and Azide free (Detector)	HA721785	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-13 Recombinant Rabbit mAb [PSH02-05]	HA722641B	Mouse	ELISA(Det)
	Mouse IL-13 Recombinant Rabbit mAb [PSH02-04] - BSA and Azide free (Capture)	HA721784	Mouse	ELISA(Cap)
	Mouse IL-17A Recombinant Rabbit mAb [PSH05-36] - BSA and Azide free (Detector)	HA722265	Mouse, Rat	ELISA(Det)
	Biotin Conjugated Mouse IL-17A Recombinant Rabbit mAb [PS01-92] - Detector	HA721608	Mouse	ELISA(Det),ELISA
	Biotin Conjugated Mouse IL-17A Recombinant Rabbit mAb [PSH05-36]	HA722600B	Mouse, Rat	ELISA(Det),ELISA
	Mouse IL-17A Recombinant Rabbit mAb [PSH05-35] - BSA and Azide free (Capture)	HA722264	Mouse, Rat	ELISA(Cap)
	Mouse IL-17A Recombinant Rabbit mAb [PS01-91] - BSA and Azide free (Capture)	HA721606	Mouse	ELISA(Cap)
	Mouse IL-17A Recombinant Rabbit mAb [PS01-92] - BSA and Azide free (Detector)	HA721607	Mouse	ELISA(Det)
	Mouse IL-17E/IL-25 Recombinant Rabbit mAb [PSH15-33] - BSA and Azide free (Capture)	HA725043	Mouse	ELISA(Cap)
	Mouse IL-17E/IL-25 Recombinant Rabbit mAb [PSH15-34] - BSA and Azide free (Detector)	HA725044	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-17E/IL-25 Recombinant Rabbit mAb [PSH15-34] - Detector	HA725045B	Mouse	ELISA(Det),ELISA
	Mouse IL-17F Recombinant Rabbit mAb [PSH17-15] - BSA and Azide free (Capture/Detector)	HA723902	Mouse	ELISA(Det),ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse IL-17F Recombinant Rabbit mAb [PSH17-16] - BSA and Azide free (Detector)	HA723903	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-17F Recombinant Rabbit mAb [PSH17-16] - Detector	HA723903B	Mouse	ELISA(Det)
	Mouse IL-17F Recombinant Rabbit mAb [PSH17-14] - BSA and Azide free (Capture)	HA723901	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-17F Recombinant Rabbit mAb [PSH17-15] - Detector	HA723902B	Mouse	ELISA(Det)
	Mouse IL-18 Recombinant Rabbit mAb [PSH12-04] - BSA and Azide free (Detector)	HA723417	Mouse	ELISA(Det)
	Mouse IL-18 Recombinant Rabbit mAb [PSH12-03] - BSA and Azide free (Capture)	HA723416	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-18 Recombinant Rabbit mAb [PSH12-04] - Detector	HA723418B	Mouse	ELISA(Det),ELISA
	Mouse IL-1 alpha Recombinant Rabbit mAb [PSH13-42] - BSA and Azide free (Capture)	HA723529	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-1 alpha Recombinant Rabbit mAb [PSH13-43] - Detector	HA723531B	Mouse	ELISA(Det),ELISA
	Mouse IL-1 alpha Recombinant Rabbit mAb [PSH13-40] - BSA and Azide free (Capture)	HA723526	Mouse	ELISA(Cap)
	Mouse IL-1 alpha Recombinant Rabbit mAb [PSH13-41] - BSA and Azide free (Detector)	HA723527	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-1 alpha Recombinant Rabbit mAb [PSH13-41] - Detector	HA723528B	Mouse	ELISA(Det),ELISA
	Mouse IL-1 alpha Recombinant Rabbit mAb [PSH13-43] - BSA and Azide free (Detector)	HA723530	Mouse	ELISA(Det)
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH03-32] - BSA and Azide free (Detector)	HA721977	Mouse	ELISA(Det)
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH04-66] - BSA and Azide free (Capture)	HA722166	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-1 beta Recombinant Rabbit mAb [PSH03-32] - Detector	HA722625B	Mouse	ELISA(Det),ELISA
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH04-67] - BSA and Azide free (Detector)	HA722167	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-1 beta Recombinant Rabbit mAb [PSH04-67]	HA722608B	Mouse	ELISA(Det)
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH03-31] - BSA and Azide free (Capture)	HA721976	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-1RA Recombinant Rabbit mAb [PSH14-99] - Detector	HA723699B	Mouse	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse IL-1RA Recombinant Rabbit mAb [PSH14-99] - BSA and Azide free (Detector)	HA723698	Mouse	ELISA(Det)
	Mouse IL-1RA Recombinant Rabbit mAb [PSH14-98] - BSA and Azide free (Capture)	HA723697	Mouse	ELISA(Cap)
	Mouse IL-1RA Recombinant Rabbit mAb [PSH15-00] - BSA and Azide free (Capture)	HA723700	Mouse	ELISA(Cap)
	Mouse IL-2 Recombinant Rabbit mAb [PSH04-28] - BSA and Azide free (Detector)	HA722111	Mouse	ELISA(Det)
	Mouse IL-2 Recombinant Rabbit mAb [PSH04-27] - BSA and Azide free (Capture)	HA722110	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-2 Recombinant Rabbit mAb [PSH04-28]	HA722617B	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-21 Recombinant Rabbit mAb [PSH14-57] - Detector	HA723653B	Mouse	ELISA(Det),ELISA
	Biotin Conjugated Mouse IL-21 Recombinant Rabbit mAb [PSH14-56] - Detector	HA723651B	Mouse	ELISA(Det),ELISA
	Mouse IL-21 Recombinant Rabbit mAb [PSH14-57] - BSA and Azide free (Detector)	HA723652	Mouse	ELISA(Det)
	Mouse IL-21 Recombinant Rabbit mAb [PSH14-56] - BSA and Azide free (Detector)	HA723650	Mouse	ELISA(Det)
	Mouse IL-21 Recombinant Rabbit mAb [PSH04-74] - BSA and Azide free (Capture)	HA722174	Mouse	ELISA(Cap)
	Mouse IL-21 Recombinant Rabbit mAb [PSH04-75] - BSA and Azide free (Detector)	HA722175	Mouse	ELISA(Det)
	Mouse IL-21 Recombinant Rabbit mAb [PSH04-75] - BSA and Azide free (Capture)	HA723649	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-21 Recombinant Rabbit mAb [PSH04-75] - Detector	HA723420B	Mouse	ELISA(Det)
	Mouse IL-22 Recombinant Rabbit mAb [PSH13-04] - BSA and Azide free (Capture)	HA725031	Mouse	ELISA(Cap)
	Mouse IL-22 Recombinant Rabbit mAb [PSH13-05] - BSA and Azide free (Detector)	HA725032	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-22 Recombinant Rabbit mAb [PSH13-05] - Detector	HA725033B	Mouse	ELISA(Det),ELISA
	Mouse IL-27A p28 Recombinant Rabbit mAb [PSH14-32] - BSA and Azide free (Detector)	HA723625	Mouse	ELISA(Det)
	Mouse IL-27A p28 Recombinant Rabbit mAb [PSH14-33] - BSA and Azide free (Capture)	HA723627	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-27A p28 Recombinant Rabbit mAb [PSH14-32] - Detector	HA723626B	Mouse	ELISA(Det),ELISA
	Mouse IL-27A p28 Recombinant Rabbit mAb [PSH14-31] - BSA and Azide free (Capture)	HA723624	Mouse	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse IL-3 Recombinant Rabbit mAb [PSH14-85] - BSA and Azide free (Capture)	HA723683	Mouse	ELISA(Cap)
	Mouse IL-3 Recombinant Rabbit mAb [PSH14-87] - BSA and Azide free (Detector)	HA723686	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-3 Recombinant Rabbit mAb [PSH14-86] - Detector	HA723685B	Mouse	ELISA(Det),ELISA
	Mouse IL-3 Recombinant Rabbit mAb [PSH14-86] - BSA and Azide free (Detector)	HA723684	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-3 Recombinant Rabbit mAb [PSH14-87] - Detector	HA723687B	Mouse	ELISA(Det),ELISA
	Mouse IL-33 Recombinant Rabbit mAb [PSH08-42] - BSA and Azide free (Detector)	HA722986	Mouse	ELISA(Det)
	Mouse IL-33 Recombinant Rabbit mAb [PSH16-63] - BSA and Azide free (Capture)	HA723865	Mouse	ELISA(Cap)
	Mouse IL-33 Recombinant Rabbit mAb [PSH16-62] - BSA and Azide free (Detector)	HA723864	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-33 Recombinant Rabbit mAb [PSH08-42] - Detector	HA722987B	Mouse	ELISA(Det),ELISA
	Mouse IL-33 Recombinant Rabbit mAb [PSH08-41] - BSA and Azide free (Capture)	HA722985	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-33 Recombinant Rabbit mAb [PSH16-62] - Detector	HA723864B	Mouse	ELISA(Det),ELISA
	Mouse IL-33 Recombinant Rabbit mAb [PSH16-61] - BSA and Azide free (Capture)	HA723863	Mouse	ELISA(Cap)
	Mouse IL-4 Recombinant Rabbit mAb [PSH04-68] - BSA and Azide free (Capture)	HA722168	Mouse	ELISA(Cap)
	Mouse IL-4 Recombinant Rabbit mAb [PSH04-69] - BSA and Azide free (Detector)	HA722169	Mouse	ELISA(Det)
	Mouse IL-4 Recombinant Rabbit mAb [PSH05-34] - BSA and Azide free (Detector)	HA722263	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-4 Recombinant Rabbit mAb [PSH04-69]	HA722607B	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-4 Recombinant Rabbit mAb [PSH05-34]	HA722597B	Mouse	ELISA(Det)
	Mouse IL-4 Recombinant Rabbit mAb [PSH05-33] - BSA and Azide free (Capture)	HA722262	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-5 Recombinant Rabbit mAb [PSH12-78] - Detector	HA725036B	Mouse	ELISA(Det),ELISA
	Mouse IL-5 Recombinant Rabbit mAb [PSH12-77] - BSA and Azide free (Capture)	HA725034	Mouse	ELISA(Cap)
	Mouse IL-5 Recombinant Rabbit mAb [PSH12-78] - BSA and Azide free (Detector)	HA725035	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PS01-83] - BSA and Azide free (Detector)	HA721602	Mouse	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PSH0-60] - Detector	HA722647B	Mouse	ELISA(Det)
	Mouse IL-6 Recombinant Rabbit mAb [PSH0-59] - BSA and Azide free (Detector)	HA721412	Mouse	ELISA(Det)
	Mouse IL-6 Recombinant Rabbit mAb [PS01-74] - BSA and Azide free (Capture/Detector)	HA721600	Mouse	ELISA(Cap),ELISA(Det)
	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PSH0-59] - Detector	HA722646B	Mouse	ELISA(Det)
	Mouse IL-6 Recombinant Rabbit mAb [PSH14-47] - BSA and Azide free (Capture)	HA723642	Mouse	ELISA(Cap)
	Mouse IL-6 Recombinant Rabbit mAb [PS01-83] - BSA and Azide free (Detector)	HA721601	Mouse	ELISA(Det)
	Mouse IL-6 Recombinant Rabbit mAb [PSH0-60] - BSA and Azide free (Detector)	HA721413	Mouse	ELISA(Det)
	Mouse IL-6 Recombinant Rabbit mAb [PSH0-58] - BSA and Azide free (Capture)	HA721411	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PS01-74] - Detector	HA723643B	Mouse	ELISA(Det)
	Mouse IL-7 Recombinant Rabbit mAb [PSH14-29] - BSA and Azide free (Capture)	HA723621	Mouse	ELISA(Cap)
	Mouse IL-7 Recombinant Rabbit mAb [PSH14-30] - BSA and Azide free (Detector)	HA723622	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-7 Recombinant Rabbit mAb [PSH14-30] - Detector	HA723623B	Mouse	ELISA(Det),ELISA
	Mouse IL-9 Recombinant Rabbit mAb [PS01-82] - BSA and Azide free (Detector)	HA721604	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-9 Recombinant Rabbit mAb [PS01-82] - Detector	HA721605	Mouse	ELISA(Det),ELISA
	Mouse IL-9 Recombinant Rabbit mAb [PS01-81] - BSA and Azide free (Capture)	HA721603	Mouse	ELISA(Cap)
	Mouse Interferon alpha 1 Recombinant Rabbit mAb [PSH14-76] - BSA and Azide free (Detector)	HA723670	Mouse	ELISA(Det)
	Biotin Conjugated Mouse Interferon alpha 1 Recombinant Rabbit mAb [PSH14-76] - Detector	HA723671B	Mouse	ELISA(Det),ELISA
	Mouse Interferon alpha 1 Recombinant Rabbit mAb [PSH14-77] - BSA and Azide free (Capture)	HA723672	Mouse	ELISA(Cap)
	Mouse Interferon alpha 1 Recombinant Rabbit mAb [PSH14-75] - BSA and Azide free (Capture)	HA723669	Mouse	ELISA(Cap)
	Mouse LIF Recombinant Rabbit mAb [PSH02-00] - BSA and Azide free (Capture)	HA721758	Mouse	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse LIF Recombinant Rabbit mAb [PSH02-01] - BSA and Azide free (Detector)	HA721759	Mouse	ELISA(Det)
	Biotin Conjugated Mouse LIF Recombinant Rabbit mAb [PSH02-01]	HA722643B	Mouse	ELISA(Det)
	Mouse L-Selectin / CD62L Recombinant Rabbit mAb [PSH16-72] - BSA and Azide free (Capture)	HA723874	Mouse	ELISA(Cap)
	Mouse L-Selectin / CD62L Recombinant Rabbit mAb [PSH16-73] - BSA and Azide free (Detector)	HA723875	Mouse	ELISA(Det)
	Biotin Conjugated Mouse L-Selectin / CD62L Recombinant Rabbit mAb [PSH16-73] - Detector	HA723875B	Mouse	ELISA(Det),ELISA
	Mouse L-Selectin / CD62L Recombinant Rabbit mAb [PSH16-74] - BSA and Azide free (Capture)	HA723876	Mouse	ELISA(Cap)
	Mouse M-CSF/CSF1 Recombinant Rabbit mAb [PSH15-85] - BSA and Azide free (Detector)	HA723793	Mouse	ELISA(Det)
	Mouse M-CSF/CSF1 Recombinant Rabbit mAb [PSH15-84] - BSA and Azide free (Capture)	HA723792	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse M-CSF/CSF1 Recombinant Rabbit mAb [PSH15-85] - Detector	HA723794B	Mouse	ELISA(Det),ELISA
	Biotin Conjugated Mouse MDC/CCL22 Recombinant Rabbit mAb [PSH15-77] - Detector	HA723782B	Mouse	ELISA(Det),ELISA
	Mouse MDC/CCL22 Recombinant Rabbit mAb [PSH15-77] - BSA and Azide free (Detector)	HA723781	Mouse	ELISA(Det)
	Mouse MDC/CCL22 Recombinant Rabbit mAb [PSH15-78] - BSA and Azide free (Capture)	HA723783	Mouse	ELISA(Cap)
	Mouse MDC/CCL22 Recombinant Rabbit mAb [PSH15-76] - BSA and Azide free (Capture)	HA723780	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse MMP-9 Recombinant Rabbit mAb [PSH15-81] - Detector	HA723788B	Mouse	ELISA(Det),ELISA
	Biotin Conjugated Mouse MMP-9 Recombinant Rabbit mAb [PSH15-80] - Detector	HA723786B	Mouse	ELISA(Det),ELISA
	HRP Conjugated Mouse MMP-9 Recombinant Rabbit mAb [PSH15-81] - Detector	HA723788H	Mouse	ELISA(Det),ELISA
	Mouse MMP-9 Recombinant Rabbit mAb [PSH15-81] - BSA and Azide free (Detector)	HA723787	Mouse	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse MMP-9 Recombinant Rabbit mAb [PSH15-80] - BSA and Azide free (Capture/Detector)	HA723785	Mouse	ELISA(Det),ELISA(Cap)
	Mouse MMP-9 Recombinant Rabbit mAb [PSH15-79] - BSA and Azide free (Capture)	HA723784	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse Osteopontin Recombinant Rabbit mAb [PSH17-50] - Detector	HA723925B	Mouse	ELISA(Det),ELISA
	Mouse Osteopontin Recombinant Rabbit mAb [PSH17-49] - BSA and Azide free (Capture)	HA723924	Mouse	ELISA(Cap)
	Mouse Osteopontin Recombinant Rabbit mAb [PSH17-50] - BSA and Azide free (Detector)	HA723925	Mouse	ELISA(Det)
	Mouse PD1 Recombinant Rabbit mAb [PSH04-17] - BSA and Azide free (Capture)	HA722100	Mouse	ELISA(Cap)
	Mouse PD1 Recombinant Rabbit mAb [PSH04-18] - BSA and Azide free (Detector)	HA722101	Mouse	ELISA(Det)
	Biotin Conjugated Mouse PD1 Recombinant Rabbit mAb [PSH04-18]	HA722613B	Mouse	ELISA(Det)
	Mouse PD-L1 Recombinant Rabbit mAb [PSH04-83] - BSA and Azide free (Capture)	HA722189	Mouse	ELISA(Cap)
	Mouse PD-L1 Recombinant Rabbit mAb [PSH04-84] - BSA and Azide free (Detector)	HA722190	Mouse	ELISA(Det)
	Biotin Conjugated Mouse PD-L1 Recombinant Rabbit mAb [PSH04-84]	HA722602B	Mouse	ELISA(Det)
	Biotin Conjugated Mouse Pro-Collagen I alpha 1 Recombinant Rabbit mAb [PSH07-99]	HA722936B	Mouse	ELISA(Det)
	Mouse Pro-Collagen I alpha 1 Recombinant Rabbit mAb [PSH07-99] - BSA and Azide free (Detector)	HA722935	Mouse	ELISA(Det)
	HRP Conjugated Mouse Pro-Collagen I alpha 1 Recombinant Rabbit mAb [PSH07-99] - Detector	HA722936H	Mouse	ELISA(Det)
	Mouse Pro-Collagen I alpha 1 Recombinant Rabbit mAb [PSH07-98] - BSA and Azide free (Capture)	HA722934	Mouse	ELISA(Cap)
Cell treatment	Mouse RANKL Recombinant Rabbit mAb [PSH07-92] - BSA and Azide free (Detector)	HA722924	Mouse	ELISA(Det)
Cell treatment	Mouse RANKL Recombinant Rabbit mAb [PSH07-90] - BSA and Azide free (Capture)	HA722922	Mouse	ELISA(Cap)
Cell treatment	Mouse RANKL Recombinant Rabbit mAb [PSH07-91] - BSA and Azide free (Capture)	HA722923	Mouse	ELISA(Cap)
Cell treatment	Biotin Conjugated Mouse RANKL Recombinant Rabbit mAb [PSH07-92]	HA722925B	Mouse	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	HRP Conjugated Mouse/Rat Complement C3 Recombinant Rabbit mAb [PSH16-82] - Detector	HA725145H	Mouse, Rat	ELISA(Det),ELISA
	Biotin Conjugated Mouse/Rat Complement C3 Recombinant Rabbit mAb [PSH16-82] - Detector	HA725145B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat Complement C3 Recombinant Rabbit mAb [PSH16-82] - BSA and Azide free (Detector)	HA725144	Mouse, Rat	ELISA(Det)
	Mouse/Rat Complement C3 Recombinant Rabbit mAb [PSH16-81] - BSA and Azide free (Capture)	HA725143	Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Mouse/Rat CXCL1/GRO alpha Recombinant Rabbit mAb [PSH13-15] - Detector	HA725030B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat CXCL1/GRO alpha Recombinant Rabbit mAb [PSH13-15] - BSA and Azide free (Detector)	HA725029	Mouse, Rat	ELISA(Det)
	Mouse/Rat CXCL1/GRO alpha Recombinant Rabbit mAb [PSH13-14] - BSA and Azide free (Capture)	HA725028	Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Mouse/Rat E-Cadherin Recombinant Rabbit mAb [PSH17-25] - Detector	HA723912B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat E-Cadherin Recombinant Rabbit mAb [PSH17-24] - BSA and Azide free (Capture)	HA723911	Mouse, Rat	ELISA(Cap)
	Mouse/Rat E-Cadherin Recombinant Rabbit mAb [PSH17-25] - BSA and Azide free (Detector)	HA723912	Mouse, Rat	ELISA(Det)
	Mouse/Rat Factor XII Recombinant Rabbit mAb [PSH12-23] - BSA and Azide free (Detector)	HA723444	Mouse, Rat	ELISA(Det)
	Biotin Conjugated Mouse/Rat Factor XII Recombinant Rabbit mAb [PSH12-23] - Detector	HA723445B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat Factor XII Recombinant Rabbit mAb [PSH12-21] - BSA and Azide free (Capture)	HA723443	Mouse, Rat	ELISA(Cap)
	Mouse/Rat GDF15 Recombinant Rabbit mAb [PSH16-66] - BSA and Azide free (Capture)	HA723868	Mouse, Rat	ELISA(Cap)
	Mouse/Rat GDF15 Recombinant Rabbit mAb [PSH16-64] - BSA and Azide free (Capture)	HA723866	Mouse, Rat	ELISA(Cap)
	Mouse/Rat GDF15 Recombinant Rabbit mAb [PSH16-65] - BSA and Azide free (Detector)	HA723867	Mouse, Rat	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Mouse/Rat GDF15 Recombinant Rabbit mAb [PSH16-65] - Detector	HA723867B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat IGF1 Recombinant Rabbit mAb [PSH12-42] - BSA and Azide free (Capture)	HA725037	Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Mouse/Rat IGF1 Recombinant Rabbit mAb [PSH12-43] - Detector	HA725039B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat IGF1 Recombinant Rabbit mAb [PSH12-43] - BSA and Azide free (Detector)	HA725038	Mouse, Rat	ELISA(Det)
	HRP Conjugated Mouse/Rat IGF1 Recombinant Rabbit mAb [PSH12-43] - Detector	HA725039H	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat Leptin Recombinant Rabbit mAb [PSH08-95] - BSA and Azide free (Detector)	HA723047	Mouse, Rat	ELISA(Det)
	Biotin Conjugated Mouse/Rat Leptin Recombinant Rabbit mAb [PSH08-95] - Detector	HA723048B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat Leptin Recombinant Rabbit mAb [PSH08-94] - BSA and Azide free (Capture)	HA723046	Mouse, Rat	ELISA(Cap)
	Mouse/Rat PDGF-BB Recombinant Rabbit mAb [PSH16-05] - BSA and Azide free (Capture)	HA723810	Mouse, Rat	ELISA(Cap)
	Mouse/Rat PDGF-BB Recombinant Rabbit mAb [PSH16-07] - BSA and Azide free (Capture)	HA723812	Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Mouse/Rat PDGF-BB Recombinant Rabbit mAb [PSH16-06] - Detector	HA723811B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat PDGF-BB Recombinant Rabbit mAb [PSH16-06] - BSA and Azide free (Detector)	HA723811	Mouse, Rat	ELISA(Det)
	Mouse/Rat Prealbumin Recombinant Rabbit mAb [PSH15-88] - BSA and Azide free (Detector)	HA725129	Mouse, Rat	ELISA(Det)
	Mouse/Rat Prealbumin Recombinant Rabbit mAb [PSH15-87] - BSA and Azide free (Capture)	HA725128	Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Mouse/Rat Prealbumin Recombinant Rabbit mAb [PSH15-88] - Detector	HA725129B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat SCF Recombinant Rabbit mAb [PSH16-91] - BSA and Azide free (Capture)	HA723887	Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Mouse/Rat SCF Recombinant Rabbit mAb [PSH16-92] - Detector	HA723888B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat SCF Recombinant Rabbit mAb [PSH16-93] - BSA and Azide free (Capture)	HA723889	Mouse, Rat	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse/Rat SCF Recombinant Rabbit mAb [PSH16-92] - BSA and Azide free (Detector)	HA723888	Mouse, Rat	ELISA(Det)
	Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-30] - BSA and Azide free (Detector)	HA721975	Mouse, Rat	ELISA(Det)
	Biotin Conjugated Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-30]	HA722626B	Mouse, Rat	ELISA(Det)
	Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-29] - BSA and Azide free (Capture)	HA721974	Mouse, Rat	ELISA(Cap)
	Mouse Resistin Recombinant Rabbit mAb [PSH16-04] - BSA and Azide free (Detector)	HA723809	Mouse	ELISA(Det)
	HRP Conjugated Mouse Resistin Recombinant Rabbit mAb [PSH16-04] - Detector	HA723809H	Mouse	ELISA(Det),ELISA
	Biotin Conjugated Mouse Resistin Recombinant Rabbit mAb [PSH16-04] - Detector	HA723809B	Mouse	ELISA(Det),ELISA
	Mouse Resistin Recombinant Rabbit mAb [PSH16-03] - BSA and Azide free (Capture)	HA723808	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse TNF alpha Recombinant Rabbit mAb [PSH03-77]	HA722623B	Mouse	ELISA(Det)
	Mouse TNF alpha Recombinant Rabbit mAb [PSH03-77] - BSA and Azide free (Detector)	HA722041	Mouse	ELISA(Det)
	Mouse TNF alpha Recombinant Rabbit mAb [PSH03-76] - BSA and Azide free (Capture)	HA722040	Mouse	ELISA(Cap)
	Mouse TREM2 Recombinant Rabbit mAb [PSH02-03] - BSA and Azide free (Detector)	HA721761	Mouse	ELISA(Det)
	Biotin Conjugated Mouse TREM2 Recombinant Rabbit mAb [PSH02-03] - Detector	HA722642B	Mouse	ELISA(Det),ELISA
	Mouse TREM2 Recombinant Rabbit mAb [PSH02-02] - BSA and Azide free (Capture)	HA721760	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse VCAM1 Recombinant Rabbit mAb [PSH16-32] - Detector	HA723830B	Mouse	ELISA(Det),ELISA
	Mouse VCAM1 Recombinant Rabbit mAb [PSH16-33] - BSA and Azide free (Capture)	HA723831	Mouse	ELISA(Cap)
	HRP Conjugated Mouse VCAM1 Recombinant Rabbit mAb [PSH16-32] - Detector	HA723830H	Mouse	ELISA(Det),ELISA
	Mouse VCAM1 Recombinant Rabbit mAb [PSH16-32] - BSA and Azide free (Detector)	HA723830	Mouse	ELISA(Det)
	Mouse VCAM1 Recombinant Rabbit mAb [PSH16-31] - BSA and Azide free (Capture)	HA723829	Mouse	ELISA(Cap)
	Mouse VEGF Recombinant Rabbit mAb [PSH02-73] - BSA and Azide free (Detector)	HA721851	Mouse	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse VEGF Recombinant Rabbit mAb [PSH02-72] - BSA and Azide free (Capture)	HA721850	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse VEGF Recombinant Rabbit mAb [PSH02-73]	HA722636B	Mouse	ELISA(Det)
	Neuropeptide Y Recombinant Rabbit mAb [PSH16-90] - BSA and Azide free (Detector)	HA723886	Human, Mouse, Rat	ELISA(Det)
	Neuropeptide Y Recombinant Rabbit mAb [PSH16-89] - BSA and Azide free (Detector)	HA723885	Human, Mouse, Rat	ELISA(Det)
	Neuropeptide Y Recombinant Rabbit mAb [PSH16-88] - BSA and Azide free (Capture)	HA723884	Human, Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Neuropeptide Y Recombinant Rabbit mAb [PSH16-90] - Detector	HA723886B	Human, Mouse, Rat	ELISA(Det),ELISA
	Biotin Conjugated Neuropeptide Y Recombinant Rabbit mAb [PSH16-89] - Detector	HA723885B	Human, Mouse, Rat	ELISA(Det),ELISA
	Non-Phospho-Tau (T217) Recombinant Rabbit mAb [PSH16-45] - BSA and Azide free (Detector)	HA723845	Human	ELISA(Det)
	Non-Phospho-Tau (T217) Recombinant Rabbit mAb [PSH16-43] - BSA and Azide free (Capture)	HA723843	Human, Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Non-Phospho-Tau (T217) Recombinant Rabbit mAb [PSH16-45] - Detector	HA723846B	Human	ELISA(Det),ELISA
	Non-Phospho-Tau (T217) Recombinant Rabbit mAb [PSH16-44] - BSA and Azide free (Capture)	HA723844	Human	ELISA(Cap)
	NuMA Mouse mAb [PSH03-07] - BSA and Azide free (Capture)	HA601223	Human	ELISA(Cap)
	NuMA Mouse mAb [PSH03-08] - BSA and Azide free (Detector)	HA601224	Human	ELISA(Det)
	Biotin Conjugated NuMA Mouse mAb [PSH03-08] - Detector	HA601357B	Human	ELISA(Det),ELISA
	Phospho-Tau (T217) Recombinant Rabbit mAb [PSH17-41] - BSA and Azide free (Capture)	HA723922	Human, Mouse	ELISA(Cap)
	Phospho-Tau (T217) Recombinant Rabbit mAb [PSH16-42] - BSA and Azide free (Capture)	HA723842	Human	ELISA(Cap)
	Phospho-Tau (T217) Recombinant Rabbit mAb [PSH16-75] - BSA and Azide free (Capture)	HA723877	Human	ELISA(Cap)
	Phospho-Tau (T217) Recombinant Rabbit mAb [PSH16-41] - BSA and Azide free (Capture)	HA723841	Human, Mouse	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	P-Selectin / CD62P Recombinant Rabbit mAb [PSH16-71] - BSA and Azide free (Detector)	HA723873	Human, Mouse, Rat	ELISA(Det)
	P-Selectin / CD62P Recombinant Rabbit mAb [PSH16-70] - BSA and Azide free (Capture)	HA723872	Human, Mouse, Rat	ELISA(Cap)
	Biotin Conjugated P-Selectin / CD62P Recombinant Rabbit mAb [PSH16-71] - Detector	HA723873B	Human, Mouse, Rat	ELISA(Det),ELISA
	Biotin Conjugated Rat CRP Recombinant Rabbit mAb [PSH12-80] - Detector	HA725051B	Rat	ELISA(Det),ELISA
	HRP Conjugated Rat CRP Recombinant Rabbit mAb [PSH12-80] - Detector	HA725051H	Rat	ELISA(Det),ELISA
	Rat CRP Recombinant Rabbit mAb [PSH12-79] - BSA and Azide free (Capture)	HA725049	Rat	ELISA(Cap)
	Rat CRP Recombinant Rabbit mAb [PSH12-80] - BSA and Azide free (Detector)	HA725050	Rat	ELISA(Det)
	Biotin Conjugated Rat ICAM1 Recombinant Rabbit mAb [PSH17-52] - Detector	HA723927B	Rat	ELISA(Det),ELISA
	Rat ICAM1 Recombinant Rabbit mAb [PSH17-52] - BSA and Azide free (Detector)	HA723927	Rat	ELISA(Det)
	Rat ICAM1 Recombinant Rabbit mAb [PSH17-53] - BSA and Azide free (Capture)	HA723928	Rat	ELISA(Cap)
	Rat ICAM1 Recombinant Rabbit mAb [PSH17-51] - BSA and Azide free (Capture)	HA723926	Rat	ELISA(Cap)
	Rat IFN gamma Recombinant Rabbit mAb [PSH05-69] - BSA and Azide free (Capture)	HA722421	Rat	ELISA(Cap)
	Rat IFN gamma Recombinant Rabbit mAb [PS01-97] - BSA and Azide free (Detector)	HA721619	Rat	ELISA(Det)
	Rat IFN gamma Recombinant Rabbit mAb [PS01-96] - BSA and Azide free (Capture)	HA721618	Rat	ELISA(Cap)
	Biotin Conjugated Rat IFN gamma Recombinant Rabbit mAb [PSH05-70]	HA722589B	Rat	ELISA(Det)
	Biotin Conjugated Rat IFN gamma Recombinant Rabbit mAb [PS01-97] - Detector	HA721620	Rat	ELISA(Det)
	Rat IFN gamma Recombinant Rabbit mAb [PSH05-70] - BSA and Azide free (Detector)	HA722422	Rat	ELISA(Det)
	Rat IL-10 Recombinant Rabbit mAb [PSH02-21] - BSA and Azide free (Capture)	HA721786	Rat	ELISA(Cap)
	Rat IL-10 Recombinant Rabbit mAb [PSH02-22] - BSA and Azide free (Detector)	HA721787	Rat	ELISA(Det)
	Biotin Conjugated Rat IL-10 Recombinant Rabbit mAb [PSH02-22]	HA722637B	Rat	ELISA(Det)
	Rat IL-1 beta Recombinant Rabbit mAb [PSH02-60] - BSA and Azide free (Detector)	HA721837	Rat	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Rat IL-1 beta Recombinant Rabbit mAb [PS01-86] - BSA and Azide free (Capture)	HA721609	Rat	ELISA(Cap)
	Biotin Conjugated Rat IL-1 beta Recombinant Rabbit mAb [PSH02-60]	HA722634B	Rat	ELISA(Det)
	Biotin Conjugated Rat IL-1 beta Recombinant Rabbit mAb [PS01-87] - Detector	HA721611	Rat	ELISA(Det),ELISA
	Rat IL-1 beta Recombinant Rabbit mAb [PSH02-59] - BSA and Azide free (Capture)	HA721836	Rat	ELISA(Cap)
	Rat IL-1 beta Recombinant Rabbit mAb [PS01-87] - BSA and Azide free (Detector)	HA721610	Rat	ELISA(Det)
	Rat IL-2 Recombinant Rabbit mAb [PS01-88] - BSA and Azide free (Capture)	HA721612	Rat	ELISA(Cap)
	Biotin Conjugated Rat IL-2 Recombinant Rabbit mAb [PS01-89] - Detector	HA721614	Rat	ELISA(Det)
	Rat IL-2 Recombinant Rabbit mAb [PSH05-72] - BSA and Azide free (Detector)	HA722424	Rat	ELISA(Det)
	Rat IL-2 Recombinant Rabbit mAb [PSH05-71] - BSA and Azide free (Capture)	HA722423	Rat	ELISA(Cap)
	Biotin Conjugated Rat IL-2 Recombinant Rabbit mAb [PSH05-72]	HA722590B	Rat	ELISA(Det)
	Rat IL-2 Recombinant Rabbit mAb [PS01-89] - BSA and Azide free (Detector)	HA721613	Rat	ELISA(Det)
	Biotin Conjugated Rat IL-4 Recombinant Rabbit mAb [PSH04-22]	HA722614B	Rat	ELISA(Det)
	Rat IL-4 Recombinant Rabbit mAb [PSH04-21] - BSA and Azide free (Capture)	HA722104	Rat	ELISA(Cap)
	Rat IL-4 Recombinant Rabbit mAb [PSH04-22] - BSA and Azide free (Detector)	HA722105	Rat	ELISA(Det)
	Biotin Conjugated Rat IL-6 Recombinant Rabbit mAb [PS01-94] - BSA and Azide free (Detector)	HA721617	Rat	ELISA(Det)
	Rat IL-6 Recombinant Rabbit mAb [PS01-94] - BSA and Azide free (Detector)	HA721616	Rat	ELISA(Det)
	Rat IL-6 Recombinant Rabbit mAb [PS01-90] - BSA and Azide free (Capture)	HA721615	Rat	ELISA(Cap)
	Rat TNF alpha Recombinant Rabbit mAb [PSH04-71] - BSA and Azide free (Detector)	HA722171	Rat	ELISA(Det)
	Rat TNF alpha Recombinant Rabbit mAb [PSH04-70] - BSA and Azide free (Capture)	HA722170	Rat	ELISA(Cap)
	Biotin Conjugated Rat TNF alpha Recombinant Rabbit mAb [PSH04-71]	HA722606B	Rat	ELISA(Det)
	S100B Recombinant Rabbit mAb [SC57-02] - BSA and Azide free (Capture)	HA723728	Human, Mouse, Rat	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated S100B Recombinant Rabbit mAb [PSH15-25] - Detector	HA723730B	Human, Mouse, Rat	ELISA(Det),ELISA
	S100B Recombinant Rabbit mAb [PSH15-25] - BSA and Azide free (Detector)	HA723729	Human, Mouse, Rat	ELISA(Det)
	SARS-CoV-2 Nucleocapsid protein Mouse mAb [A5B8] - BSA and Azide free (Capture)	HA600008	Sars-cov-2	WB,ELISA(Cap)
	SARS-CoV-2 Nucleocapsid protein Mouse mAb [H4G1] - BSA and Azide free (Capture)	HA600006	Sars-cov-2	WB,ELISA(Cap)
	SARS-CoV-2 Nucleocapsid protein Mouse mAb [A5B10] - BSA and Azide free (Capture)	HA600007	Sars-cov-2	WB,ELISA(Cap)
	SARS-CoV-2 Nucleocapsid protein Mouse mAb [A4G9] - BSA and Azide free (Detector)	HA600005	Sars-cov-2	WB,ELISA(Det)
	SNAP25 Recombinant Rabbit mAb [PSH17-66] - BSA and Azide free (Detector)	HA725228	Human, Mouse	ELISA(Det)
	Biotin Conjugated SNAP25 Recombinant Rabbit mAb [PSH17-66] - Detector	HA725228B	Human, Mouse	ELISA(Det),ELISA
	SNAP25 Recombinant Rabbit mAb [PSH17-65] - BSA and Azide free (Capture)	HA725227	Human, Mouse	ELISA(Cap)
	Tau Recombinant Rabbit mAb [PSH16-49] - BSA and Azide free (Detector)	HA723851	Human	ELISA(Det)
	Tau Recombinant Rabbit mAb [PSH16-48] - BSA and Azide free (Capture)	HA723850	Human, Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Tau Recombinant Rabbit mAb [PSH16-49] - Detector	HA723852B	Human	ELISA(Det),ELISA
	Biotin Conjugated Tau Recombinant Rabbit mAb [PSH17-11] - Detector	HA723898B	Human	ELISA(Det),ELISA
	Tau Recombinant Rabbit mAb [PSH10-42] - BSA and Azide free (Detector)	HA723209	Human, Mouse, Rat	ELISA(Det)
	Biotin Conjugated Tau Recombinant Rabbit mAb [PSH10-42] - Detector	HA723209B	Human, Mouse, Rat	ELISA(Det),ELISA
	Tau Recombinant Rabbit mAb [PSH17-11] - BSA and Azide free (Detector)	HA723898	Human	ELISA(Det)
	TGF beta 1 Recombinant Rabbit mAb [PSH15-21] - BSA and Azide free (Capture)	HA723725	Human, Mouse, Rat	ELISA(Cap)
	TGF beta 1 Recombinant Rabbit mAb [PSH15-22] - BSA and Azide free (Detector)	HA723726	Human, Mouse, Rat	ELISA(Det)
	Biotin Conjugated TGF beta 1 Recombinant Rabbit mAb [PSH15-22] - Detector	HA723727B	Human, Mouse, Rat	ELISA(Det),ELISA
	Biotin Conjugated TNFRSF14/HVEM Recombinant Rabbit mAb [PSH06-42]	HA722678B	Human	ELISA(Det)
	TNFRSF14/HVEM Recombinant Rabbit mAb [PSH06-42] - BSA and Azide free (Detector)	HA722677	Human	ELISA(Det)
	TNFRSF14/HVEM Recombinant Rabbit mAb [PSH06-41] - BSA and Azide free (Capture)	HA722676	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

