

外泌体研究-产品清单

高级验证	产品名	货号	反应种属	验证应用
	A2M Recombinant Rabbit mAb [JM81-41]	ET1703-69	Human	WB,IHC-P
	AChE Rabbit pAb	ER62571	Human, Mouse, Rat	WB,IHC-P
	ADAM10 Rabbit pAb	ER1706-97	Human	WB,IHC-P
KD	ADAM10 Recombinant Rabbit mAb [JM32-11]	ET1703-60	Human, Mouse, Rat	WB,IHC-P,IP
	ALCAM Mouse mAb [C7-E12]	EM1708-92	Human	WB,IF-Cell,IHC-P,FC
	ALCAM/CD166 Mouse mAb [7C2]	RT1029	Human, Mouse, Rat	WB,IP,IF,IHC-P
KD	ALIX Recombinant Rabbit mAb [JM85-31]	ET1705-74	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell
KD	ALIX Recombinant Rabbit mAb [JM85-31] - BSA and Azide free	HA750439	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell
	alpha 1 Antitrypsin Rabbit pAb	R1103-1	Human	WB
	alpha 1 Antitrypsin Recombinant Rabbit mAb [JF10-03]	ET1702-87	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Alpha-2-macroglobulin Mouse mAb [8-B10]	M1510-23	Human, Mouse, Rat	ELISA,IF-Cell,IHC-P
	Alpha-2-macroglobulin Mouse mAb [B1-D7]	M1510-24	Human	WB,IHC-P,ELISA
	Alpha-2-macroglobulin Mouse mAb [A5-A9]	M1510-22	Human	WB,IF-Cell,IHC-P,ELISA
	Annexin A1 Rabbit pAb	ER1803-70	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Annexin A1 Recombinant Rabbit mAb [PSH04-13]	HA722094	Human	IHC-P,IF-Tissue
KD	Annexin A2 Recombinant Rabbit mAb [JA42-30]	ET1704-49	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC
	Annexin IV Recombinant Rabbit mAb [JG37-83]	ET7108-56	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	Annexin V Recombinant Rabbit mAb [JF50-11]	ET1702-62	Human, Mouse, Rat	WB,FC
	ANXA6 Recombinant Rabbit mAb [PSH05-82]	HA722481	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IP
	ANXA6 Recombinant Rabbit mAb [PSH05-82] - BSA and Azide free	HA751030	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IP
KD	ANXA7 Recombinant Rabbit mAb [JE65-36]	HA721078	Human, Rat	WB
KD	ANXA7 Recombinant Rabbit mAb [JE65-60]	HA721068	Human	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Aquaporin 2 Mouse mAb [7H8B5]	EM1710-44	Human	WB
	Aquaporin 2 Rabbit pAb	ER1903-86	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Aquaporin 2 Rabbit pAb	ER1903-88	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Aquaporin 2 Recombinant Rabbit mAb [PSH09-19]	HA723077	Human, Mouse, Rat	IHC-P,IHC-Fr,IF-Tissue
	Aquaporin 2 Recombinant Rabbit mAb [PSH09-19] - BSA and Azide free	HA751273	Human, Mouse, Rat	IHC-P,IHC-Fr,IF-Tissue
KD	ARF1 Recombinant Rabbit mAb [SD2002]	ET1612-12	Human, Mouse, Rat	WB,IHC-P,FC
	ARF1 Rabbit pAb	0807-7	Human	WB,IHC-P,IF-Cell,FC
Relative expression	beta 2 Microglobulin Recombinant Rabbit mAb [PSH12-14]	HA723430	Human	WB,IF-Cell,IHC-P,FC,IP
	beta 2 Microglobulin Rabbit pAb	0808-2	Human	WB,IHC-P,FC
Relative expression	beta 2 Microglobulin Recombinant Rabbit mAb [PSH12-14] - BSA and Azide free	HA751443	Human	WB,IF-Cell,IHC-P,FC,IP
KD	beta 2 Microglobulin Recombinant Rabbit mAb [SU35-05]	ET1608-45	Human, Mouse, Rat	WB,IF-Cell
KD	Calnexin Rabbit pAb	ER1803-42	Human, Mouse, Rat	WB,IHC-P,FC,IP
KD	Calnexin Recombinant Rabbit mAb [SN20-54]	ET1611-86	Human, Mouse, Rat	WB,IP,IHC-P,IF-Tissue
KD	Calnexin Recombinant Rabbit mAb [SN20-54] - BSA and Azide free	HA750275	Human, Mouse, Rat	WB,IP,IHC-P,IF-Tissue
	CCT2 Mouse mAb [C4-B5]	EM1709-68	Human, Mouse, Rat, Monkey	WB,FC
	CCT2 Recombinant Rabbit mAb [JE37-96]	HA722240	Human, Mouse, Rat	WB,IP
	CD105 Rabbit pAb	0902-1	Mouse	WB,IF-Cell,IHC-P,FC
	CD105 Mouse mAb [A5E8]	HA600077	Human	WB
	CD105 Rabbit pAb	HA500266	Human	WB,IHC-P
Cell treatment	CD105 Recombinant Rabbit mAb [JE60-59]	HA720072	Human	WB,IHC-P,IF-Tissue
Relative expression	CD105 Recombinant Rabbit mAb [PSH15-97]	HA723802	Human	FC
Relative expression	CD105 Recombinant Rabbit mAb [PSH15-96]	HA723801	Human	FC
Relative expression	CD105 Recombinant Rabbit mAb [PSH15-96] - BSA and Azide free	HA751587	Human	FC
Relative expression	CD105 Recombinant Rabbit mAb [PSH15-97] - BSA and Azide free	HA751588	Human	FC
	CD105 Mouse mAb [H4-A9]	EM1709-99	Human	WB,IF-Cell,IHC-P,FC
	CD146 Rabbit pAb	0806-5	Human, Mouse, Rat	WB,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	CD146 Recombinant Rabbit mAb [JB42-35]	ET7107-91	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	CD147 Rabbit pAb	1003-10	Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	CD147 Recombinant Rabbit mAb [JF1-045]	ET1702-58	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
Relative expression	CD151 Mouse mAb [A6G11]	HA600059	Human, Mouse	WB,IF-Cell,IHC-P
	CD151 Mouse mAb [3G1]	RT1111	Human	WB,IP,IF,IHC-P
	CD16 Recombinant Rabbit mAb	IRS017	Human	mIHC
	CD16 Mouse mAb [Q32]	EM60067	Human	WB,IHC-P
Relative expression	CD16 Recombinant Rabbit mAb [PD00-12]	HA721149	Human	IHC-P,IF-Tissue,IF-Cell,mIHC,WB
	CD16 Recombinant Rabbit mAb [JE49-79]	ET7109-97	Human, Mouse, Rat	WB
	CD164 Rabbit pAb	ER62936	Human, Mouse, Rat	IHC-P
KD	CD166 Mouse mAb [A6A9]	HA600082	Human	IF-Cell,IHC-P,WB
KD	CD166 Mouse mAb [A6A8]	HA600083	Human	IF-Cell,IHC-P,WB
KD	CD166 Mouse mAb [5-5]	EM1902-36	Human	WB,IHC-P,IF-Cell,FC
	CD166 Recombinant Rabbit mAb [JE57-47]	HA722705	Human, Mouse, Rat	WB,IHC-P,FC,IHC-Fr
KD	CD166 Mouse mAb [11-6]	M1012-7	Human, Mouse, Rat	WB,IHC-P,FC
KD	CD63 Recombinant Rabbit mAb [SY21-02] - BSA and Azide free	HA750112	Human	WB,IHC-P
Relative expression	CD63 Recombinant Rabbit mAb [PSH06-64]	HA722731	Human	WB,IF-Cell,FC,IP,IHC-P
Relative expression	CD63 Recombinant Rabbit mAb [PSH06-64] - BSA and Azide free	HA751088	Human	WB,IF-Cell,FC,IP,IHC-P
KD	CD63 Recombinant Rabbit mAb [SY21-02]	ET1607-2	Human	WB,IHC-P
Relative expression	CD81 Recombinant Rabbit mAb [PSH10-26]	HA723193	Human, Rat	WB,IHC-P,FC
Cell treatment	CD82 Recombinant Rabbit mAb [JE32-17]	HA722533	Human	WB
	CD86 Recombinant Rabbit mAb [SJ20-00] - BSA and Azide free	HA750109	Human, Mouse, Rat	WB,IF-Cell,IP
	CD86 Recombinant Rabbit mAb [SJ20-00]	ET1606-50	Human, Mouse, Rat	WB,IF-Cell,IP
	CD86 Rabbit pAb	ER1906-01	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
Relative expression, KD	CD9 Recombinant Rabbit mAb [PSH0-95] - BSA and Azide free	HA750675	Human	WB,IHC-P
	CD9 Recombinant Rabbit mAb [SA35-08]	ET1601-9	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC
	CD9 Rabbit pAb	ER80402	Human, Mouse	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression, KD	CD9 Recombinant Rabbit mAb [PSH0-95]	HA721533	Human	WB,IHC-P
	CD9 Recombinant Rabbit mAb [SA35-08] - BSA and Azide free	HA750027	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC
	CDC42 Recombinant Rabbit mAb [JJ086-04] - BSA and Azide free	HA750320	Human, Mouse, Rat	WB,FC,IF-Cell
	CDC42 Recombinant Rabbit mAb [JJ086-04]	ET1701-7	Human, Mouse, Rat	WB,FC,IF-Cell
Relative expression	CNPase Recombinant Rabbit mAb [JF10-25]	ET1702-46	Human, Mouse, Rat, Cynomolgus monkey	WB,IHC-P,IF-Tissue
KO	DPP4 / CD26 Recombinant Rabbit mAb [JM11-42]	ET1703-93	Human, Mouse, Rat	WB
KO	DPP4 / CD26 Recombinant Rabbit mAb [JM11-42] - BSA and Azide free	HA750396	Human, Mouse, Rat	WB
Cell treatment	DPP4 / CD26 Recombinant Rabbit mAb [PSH12-59]	HA723473	Human	WB,IHC-P,IF-Cell,FC,IP
Cell treatment	DPP4 / CD26 Recombinant Rabbit mAb [PSH12-59] - BSA and Azide free	HA751446	Human	WB,IHC-P,IF-Cell,FC,IP
	Dynamin 1 Mouse mAb [A1A12]	EM1901-44	Human	WB,IF-Cell,FC
	Dynamin 1 Recombinant Mouse mAb [A1B1-R] - BSA and Azide free	HA610105	Mouse, Rat	WB
	Dynamin 1 Recombinant Mouse mAb [A1B1-R]	HA601230	Mouse, Rat	WB
	Dynamin 1 Recombinant Mouse mAb [A1B2-R] - BSA and Azide free	HA610197	Human, Mouse, Rat	WB,IHC-P
	Dynamin 1 Rabbit pAb	ER1802-77	Mouse, Rat, Human	WB,IHC-P
	Dynamin 1 Mouse mAb [A1B2]	EM1901-42	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Dynamin 1 Recombinant Mouse mAb [A1B2-R]	HA601325	Human, Mouse, Rat	WB,IHC-P
	Dynamin 1 Mouse mAb [A1B1]	EM1901-43	Human, Mouse, Rat	IHC-P,WB,IF-Cell
	EBP50 Recombinant Rabbit mAb [JJ0946]	ET1701-96	Human, Rat	WB,IHC-P
Relative expression, KD	EGFR Recombinant Rabbit mAb [PD00-75] - BSA and Azide free	HA750550	Human, Monkey	WB,IHC-P,IF-Cell,FC
Relative expression, KO	EGFR Recombinant Rabbit mAb [SP00-86]	ET1604-44	Human, Mouse, Rat, Monkey	WB,IF-Cell,IP,FC
	EGFR Rabbit pAb	0407-21	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P
	EGFR Rabbit pAb	R1608-3	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression, KD	EGFR Recombinant Rabbit mAb [SZ40-19]	ET1603-37	Human, Mouse, Monkey, Rat	WB,IF-Cell,IP,FC,mIHC
Relative expression, KD	EGFR Mouse mAb [A8E12]	HA601073	Human, Monkey	WB,IHC-P
	EGFR Rabbit pAb	ER1512-6	Human, Mouse, Monkey	WB,IF-Cell,IHC-P,FC
Relative expression, KD	EGFR Recombinant Rabbit mAb [PD00-75]	HA721220	Human, Monkey	WB,IHC-P,IF-Cell,FC
	EGFR Rabbit pAb	R1511-18	Human, Monkey	WB,IF-Cell,IHC-P
Relative expression, KO	EGFR Mouse mAb [7-F2-F8]	EM1901-67	Human, Monkey	WB,IHC-P,FC
	EGFR Recombinant Rabbit mAb	IRS046	Human	mIHC
	EGFR Mouse mAb [B6-E5-D9]	EM1701-22	Human, Mouse, Monkey, Rat	WB,IHC-P,IF-Cell,FC
Relative expression, KO	EGFR Recombinant Rabbit mAb [SP00-86] - BSA and Azide free	HA750078	Human, Mouse, Rat, Monkey	WB,IF-Cell,IP,FC
	ENO1 Recombinant Rabbit mAb [JM63-53]	ET1705-56	Human, Mouse, Rat, Zebrafish	WB,FC,IF-Tissue,IHC-P,IF-Cell
	EpCAM Recombinant Rabbit mAb	IRS254RB	Human	mIHC
	iFluor™ 488 Conjugated EpCAM Recombinant Rabbit mAb [PS01-69]	HA720217F	Human	IF-Cell
	EpCAM Recombinant Rabbit mAb [PS01-69] - BSA and Azide free	HA750703	Human	IHC-P,IF-Cell,FC
	iFluor™ 647 Conjugated EpCAM Recombinant Rabbit mAb [PS01-69]	HA720185F	Human	FC,IF-Cell
	iFluor™ 555 Conjugated EpCAM Recombinant Rabbit mAb [PS01-69]	HA720218F	Human	FC
	EpCAM Recombinant Rabbit mAb [SU03-32]	ET1608-26	Human, Mouse, Rat	WB,IHC-P
	EpCAM Mouse mAb [3-E1-E2]	EM1111	Human, Mouse	WB,IF-Cell
	EpCAM Rabbit pAb	ER1802-80	Human, Mouse, Rat	WB,IHC-P
	EpCAM Recombinant Rabbit mAb [PS01-69]	HA721644	Human	IHC-P,IF-Cell,FC
	EpCAM Mouse mAb [7E11]	EM1712-43	Human	WB,IHC-P,ELISA
	EZR Mouse mAb [F1-9A]	EM1706-9	Human, Mouse, Monkey	WB,IF-Cell,FC
	Ezrin Recombinant Rabbit mAb [PSH19-97]	HA724117	Human, Mouse, Rat	WB,IHC-P
	Ezrin Recombinant Rabbit mAb [JM10-65]	ET1703-25	Human	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	Fatty Acid Synthase Rabbit pAb	R1706-8	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KO	Fatty Acid Synthase Recombinant Rabbit mAb [JJ0939] - BSA and Azide free	HA750327	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
KO	Fatty Acid Synthase Recombinant Rabbit mAb [JJ0939]	ET1701-91	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	FCG3B Rabbit pAb	ER62049	Human	WB
KD	Filamin A Recombinant Rabbit mAb [SA30-08]	ET1601-3	Human, Mouse, Rat, Monkey	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Filamin A Rabbit pAb	R1601-3	Human, Mouse, Rat	WB
	Flotillin 1 Recombinant Rabbit mAb [JB19-45]	ET7107-82	Human, Mouse, Rat	WB,IHC-P
	Flotillin-2 Rabbit pAb	ER60092	Human, Rat, Mouse	WB,IHC-P
KD	Flotillin-2 Recombinant Rabbit mAb [PSH03-41]	HA721986	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
KD	Flotillin-2 Recombinant Rabbit mAb [PSH03-41] - BSA and Azide free	HA750875	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
KD	Gelsolin Recombinant Rabbit mAb [JB36-68]	ET7107-45	Human, Mouse, Monkey	WB,IF-Tissue,IHC-P
	GM130 Recombinant Rabbit mAb [JE42-53]	HA721282	Human, Mouse	WB,IHC-P,IF-Cell
	GM130 (cis-Golgi Marker) Rabbit pAb	R1608-7	Human, Mouse, Rat	WB,IHC-P
	GPC1 Rabbit pAb	R1512-15	Human	WB,IF-Cell,IHC-P,FC
Cell treatment	GRP78 / BIP Recombinant Mouse mAb [C9-9-R]	HA601076	Human, Mouse, Rat, Zebrafish	WB,IHC-P,IF-Cell
	GRP78 / BIP Recombinant Mouse mAb [A5-C9-G6-R]	HA601318	Human, Mouse, Rat	WB
	GRP78 / BIP Recombinant Mouse mAb [A5-C9-G6-R] - BSA and Azide free	HA610190	Human, Mouse, Rat	WB
	GRP78 / BIP Recombinant Rabbit mAb [JE01-35]	HA722202	Human, Mouse, Rat, Chicken	WB,IHC-P,IF-Tissue
Cell treatment	GRP78 / BIP Mouse mAb [2-7]	M1505-13	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,IF-Tissue,FC
Cell treatment	GRP78 / BIP Recombinant Mouse mAb [C9-9-R] - BSA and Azide free	HA610013	Human, Mouse, Rat, Zebrafish	WB,IHC-P,IF-Cell
Cell treatment	GRP78 / BIP Mouse mAb [A5-C9-G6]	M1701-1	Human, Mouse, Rat	WB
Cell treatment	GRP78 / BIP Rabbit pAb	ER1706-50	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Cell treatment	GRP78 / BIP Mouse mAb [C9-9]	M1506-2	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P
KD	GST3 Recombinant Mouse mAb [11D5-R]	HA601272	Human	WB,IF-Cell,IHC-P
KD	GST3 Mouse mAb [11D5]	EM1801-18	Human	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	GST3 Rabbit pAb	ER1802-63	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	GST3 Recombinant Rabbit mAb [JB40-79]	ET7107-71	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue
KD	GST3 Recombinant Mouse mAb [11D5-R] - BSA and Azide free	HA610147	Human	WB,IF-Cell,IHC-P
	HER2 / ErbB2 Recombinant Rabbit mAb	IRS034	Human	mIHC
	HER2 / ErbB2 Recombinant Rabbit mAb [PD00-53]	HA721178	Human	IHC-P,WB,IF-Cell,FC,mIHC
	HER2 / ErbB2 Rabbit pAb	R1511-33	Human	WB,IF-Cell,IHC-P,FC
	HER2 / ErbB2 Rabbit pAb	ER0106	Human	WB,IHC-P
	HER2 / ErbB2 Recombinant Rabbit mAb [PD00-99]	HA721210	Human	WB,IHC-P,IF-Tissue
	Hsc70 Recombinant Rabbit mAb [SR39-04] - BSA and Azide free	HA750044	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	Hsc70 Recombinant Rabbit mAb [SR39-04]	ET1602-33	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	Hsc70 Rabbit pAb	R1511-6	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Hsp20 Recombinant Rabbit mAb [SD086-03]	ET1612-81	Human, Mouse, Rat	WB,IHC-P,FC
	Hsp27 Mouse mAb [10-8]	M1505-3	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Hsp27 Recombinant Rabbit mAb [JJ09-13]	ET1701-70	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Hsp27 Recombinant Rabbit mAb [JJ09-13] - BSA and Azide free	HA750321	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Hsp27 Rabbit pAb	ER31216	Human	WB,IF-Cell,IHC-P
	HSP60 Recombinant Rabbit mAb [ST48-04] - BSA and Azide free	HA750175	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IP
	HSP60 Mouse mAb [9-49]	M1007-9	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	HSP60 Recombinant Rabbit mAb [ST48-04]	ET1609-45	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IP
	HSP60 Rabbit pAb	1007-1	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	HSP60 Mouse mAb [1-80]	EM00704	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	HSP60 Rabbit pAb	1007-3	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	Hsp70 Recombinant Rabbit mAb [SA0379]	ET1601-11	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Hsp70 Recombinant Rabbit mAb [PD01-65]	HA721497	Human	WB,IHC-P,IF-Tissue
	Hsp70 Rabbit pAb	ER50802	Human, Mouse, Rat, Cow	WB,IF-Cell,IHC-P,FC
	HSP70 Mouse mAb [3G10]	EM60062	Human, Mouse, Rat	WB,IHC-P,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	HSP70 (Acetyl K246) Rabbit pAb	ER60056	Human	WB
	Hsp90 Recombinant Rabbit mAb [JE05-06]	HA722689	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC,IP
	Hsp90 alpha Recombinant Rabbit mAb [SY14-06]	ET1605-57	Human, Mouse, Rat, Monkey	WB,IF-Cell,IF-Tissue,IHC-P
	Hsp90 alpha Mouse mAb [A10-A5-D2]	M1603-3	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Hsp90 alpha Rabbit pAb	R1510-29	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Hsp90 alpha Recombinant Rabbit mAb [SY14-06] - BSA and Azide free	HA750091	Human, Mouse, Rat, Monkey	WB,IF-Cell,IF-Tissue,IHC-P
KD	Hsp90 beta Recombinant Rabbit mAb [SY46-01] - BSA and Azide free	HA750090	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	Hsp90 beta Recombinant Rabbit mAb [SY46-01]	ET1605-56	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	HSP90 Beta Recombinant Mouse mAb [2-1-G3-R] - BSA and Azide free	HA610060	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	HSP90 Beta Recombinant Mouse mAb [2-1-G3-R]	HA601176	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	HSP90 Beta Mouse mAb [2-1-G3]	EM21103	Human, Mouse, Rat	WB,IHC-P
KD	HSPA4 Recombinant Rabbit mAb [SD0861]	ET1612-83	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Human Acetylcholinesterase Recombinant Rabbit mAb [PSH17-03] - BSA and Azide free (Detector)	HA725224	Human	ELISA(Det)
	Biotin Conjugated Human Acetylcholinesterase Recombinant Rabbit mAb [PSH17-03] - Detector	HA725224B	Human	ELISA(Det),ELISA
	Human Acetylcholinesterase Recombinant Rabbit mAb [PSH17-02] - BSA and Azide free (Capture)	HA725223	Human	ELISA(Cap)
	Biotin Conjugated Human CD147 Recombinant Rabbit mAb [PSH13-07] - Detector	HA725087B	Human	ELISA(Det),ELISA
	Human CD147 Recombinant Rabbit mAb [PSH13-07] - BSA and Azide free (Detector)	HA725086	Human	ELISA(Det)
	Human CD147 Recombinant Rabbit mAb [PSH13-06] - BSA and Azide free (Capture)	HA725085	Human	ELISA(Cap)
	Human CD16 Recombinant Mouse mAb [PSH04-95]	HA601340	Human	FC
	PE Conjugated Human CD16 Recombinant Mouse mAb [PSH04-95]	HA600129F	Human	FC
	Biotin Conjugated Human CD63 Recombinant Rabbit mAb [PSH08-44] - Detector	HA722990B	Human	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human CD63 Recombinant Rabbit mAb [PSH08-44] - BSA and Azide free (Detector)	HA722989	Human	ELISA(Det)
	Human CD63 Recombinant Rabbit mAb [PSH08-43] - BSA and Azide free (Capture)	HA722988	Human	ELISA(Cap)
	Human Endoglin / CD105 Recombinant Rabbit mAb [PSH14-58] - BSA and Azide free (Capture)	HA723654	Human	ELISA(Cap)
	Biotin Conjugated Human Endoglin / CD105 Recombinant Rabbit mAb [PSH14-59] - Detector	HA723656B	Human	ELISA(Det),ELISA
	Human Endoglin / CD105 Recombinant Rabbit mAb [PSH14-59] - BSA and Azide free (Detector)	HA723655	Human	ELISA(Det)
	Human Endoglin / CD105 Recombinant Rabbit mAb [PSH14-60] - BSA and Azide free (Capture)	HA723657	Human	ELISA(Cap)
	Biotin Conjugated Human ICAM1 Recombinant Rabbit mAb [PSH10-16] - Detector	HA723186B	Human	ELISA(Det),ELISA
	Human ICAM1 Recombinant Rabbit mAb [PSH10-15] - BSA and Azide free (Capture)	HA723184	Human	ELISA(Cap)
	Human ICAM1 Recombinant Rabbit mAb [PSH10-16] - BSA and Azide free (Detector)	HA723185	Human	ELISA(Det)
	HRP Conjugated Human ICAM1 Recombinant Rabbit mAb [PSH10-16] - Detector	HA723474H	Human	ELISA(Det),ELISA
	Human MMP-9 Recombinant Rabbit mAb [PSH04-23] - BSA and Azide free (Capture)	HA722106	Human	ELISA(Cap)
	Human MMP-9 Recombinant Rabbit mAb [PSH04-63] - BSA and Azide free (Detector)	HA722163	Human	ELISA(Det)
	Human MMP-9 Recombinant Rabbit mAb [PSH04-24] - BSA and Azide free (Detector)	HA722107	Human	ELISA(Det)
	Biotin Conjugated Human MMP-9 Recombinant Rabbit mAb [PSH04-24]	HA722615B	Human	ELISA(Det)
	Human MMP-9 Recombinant Rabbit mAb [PSH04-62] - BSA and Azide free (Capture)	HA722162	Human	ELISA(Cap)
	Human NCAM1 Recombinant Rabbit mAb [PSH08-07] - BSA and Azide free (Capture)	HA722946	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human NCAM1 Recombinant Rabbit mAb [PSH08-09] - BSA and Azide free (Detector)	HA722949	Human	ELISA(Det)
	Human NCAM1 Recombinant Rabbit mAb [PSH08-08] - BSA and Azide free (Capture/Detector)	HA722947	Human	ELISA(Det),ELISA(Cap)
	Biotin Conjugated Human NCAM1 Recombinant Rabbit mAb [PSH08-08] - Detector	HA722948B	Human	ELISA(Det)
	Biotin Conjugated Human NCAM1 Recombinant Rabbit mAb [PSH08-09] - Detector	HA722950B	Human	ELISA(Det)
	PE Conjugated Human NCAM1 / CD56 Recombinant Rabbit mAb [PSH06-53]	HA720215F	Human	FC
	APC Conjugated Human NCAM1 / CD56 Recombinant Rabbit mAb [PSH06-53]	HA723757F	Human	FC
Relative expression	ICAM1 Recombinant Rabbit mAb [ST0487] - BSA and Azide free	HA750176	Human	WB,IF-Tissue,IHC-P
Relative expression	ICAM1 Recombinant Rabbit mAb [ST0487]	ET1609-46	Human	WB,IF-Tissue,IHC-P
KD	ICAM1 Mouse mAb [F11-A4]	M1505-7	Human	WB,IF-Cell,FC
	ICAM1 Rabbit pAb	ER131207	Human	WB,IF-Cell,IHC-P,FC
	ICAM1 Mouse mAb [J0-F7]	M1511-6	Human	WB,IF-Cell,FC
Relative expression, KD	ICAM1 Recombinant Rabbit mAb [PSH09-31]	HA723089	Human	WB,IHC-P,IF-Tissue
KD	Lactate Dehydrogenase Recombinant Rabbit mAb [SU39-06] - BSA and Azide free	HA750155	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
KD	Lactate Dehydrogenase Recombinant Rabbit mAb [SU39-06]	ET1608-57	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	LDHA Recombinant Rabbit mAb [PSH13-58]	HA723547	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
KD	LDHA Rabbit pAb	ER00702	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell
	LDHA Mouse mAb [20-179]	M1007-7	Human, Mouse	IHC-P
	LDHA Recombinant Rabbit mAb [PSH13-58] - BSA and Azide free	HA751477	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	LDHA + LDHB Recombinant Rabbit mAb [PSH12-75] - BSA and Azide free	HA751451	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
	LDHA + LDHB Recombinant Rabbit mAb [PSH12-75]	HA723486	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
KD	LDHB Rabbit pAb	0807-1	Human, Mouse, Rat	WB,IF-Cell,FC,IHC-P
	LDHC Rabbit pAb	ER1912-25	Human	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	LRP-1 Mouse mAb [A4-F5]	M1211-4	Mouse	WB,IHC-P
KD	LRP1 Recombinant Rabbit mAb [SA0290]	ET1601-1	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,IP,IHC-Fr
	MDH1 Rabbit pAb	ER1912-90	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
	MDH1 Recombinant Rabbit mAb [JE35-03]	HA722196	Human, Mouse, Rat	WB,IF-Cell,FC
	MFGE8 Recombinant Rabbit mAb [PSH05-75] - BSA and Azide free	HA751025	Human	WB
	MFGE8 Recombinant Rabbit mAb [PSH05-75]	HA722442	Human	WB
	MMP-2 Mouse mAb [B11-E6]	EM1706-99	Human, Mouse	WB,FC
Relative expression	MMP-2 Recombinant Rabbit mAb [PSH11-26] - BSA and Azide free	HA751405	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
Relative expression	MMP-2 Recombinant Rabbit mAb [PSH11-26]	HA723308	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
	MMP-2 Recombinant Rabbit mAb [SI15-04]	ET1606-4	Human, Mouse, Rat	WB,IHC-P,IP,IF-Tissue
	MMP-9 Recombinant Mouse mAb [10A1]	EM1801-22	Human	ELISA
	MMP-9 Mouse mAb [2G3]	RT1401	Human	WB,IP,IF,IHC-P
Cell treatment, KD	MMP-9 Recombinant Rabbit mAb [JA80-73]	ET1704-69	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,FC
Cell treatment	MMP-9 Recombinant Rabbit mAb [PSH06-33]	HA722665	Human	WB,IHC-P
Cell treatment, KD	MMP-9 Recombinant Rabbit mAb [JA80-73] - BSA and Azide free	HA750420	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,FC
Cell treatment	MMP-9 Recombinant Rabbit mAb [PSH06-33] - BSA and Azide free	HA751070	Human	WB,IHC-P
	MMP-9 Recombinant Rabbit mAb	IRS172RB	Human	mIHC
KD	Moesin Recombinant Rabbit mAb [SC69-01]	ET1610-65	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	PE Conjugated Mouse CD16/32 Recombinant Rabbit mAb [PSH06-60]	HA720216F	Mouse	FC
	Biotin Conjugated Mouse CD86 Recombinant Rabbit mAb [PSH13-64] - Detector	HA723554B	Mouse	ELISA(Det),ELISA
	Mouse CD86 Recombinant Rabbit mAb [PSH13-64] - BSA and Azide free (Detector)	HA723553	Mouse	ELISA(Det)
	Mouse CD86 Recombinant Rabbit mAb [PSH13-63] - BSA and Azide free (Capture)	HA723552	Mouse	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Mouse CD86 Recombinant Rabbit mAb [PSH13-66] - Detector	HA723557B	Mouse	ELISA(Det),ELISA
	Mouse CD86 Recombinant Rabbit mAb [PSH13-66] - BSA and Azide free (Detector)	HA723556	Mouse	ELISA(Det)
	Mouse CD86 Recombinant Rabbit mAb [PSH13-65] - BSA and Azide free (Capture)	HA723555	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse MMP-9 Recombinant Rabbit mAb [PSH15-81] - Detector	HA723788B	Mouse	ELISA(Det),ELISA
	Mouse MMP-9 Recombinant Rabbit mAb [PSH15-81] - BSA and Azide free (Detector)	HA723787	Mouse	ELISA(Det)
	HRP Conjugated Mouse MMP-9 Recombinant Rabbit mAb [PSH15-81] - Detector	HA723788H	Mouse	ELISA(Det),ELISA
	Mouse MMP-9 Recombinant Rabbit mAb [PSH15-79] - BSA and Azide free (Capture)	HA723784	Mouse	ELISA(Cap)
	Mouse MMP-9 Recombinant Rabbit mAb [PSH15-80] - BSA and Azide free (Capture/Detector)	HA723785	Mouse	ELISA(Det),ELISA(Cap)
	Biotin Conjugated Mouse MMP-9 Recombinant Rabbit mAb [PSH15-80] - Detector	HA723786B	Mouse	ELISA(Det),ELISA
Relative expression	NCAM1 / CD56 Recombinant Rabbit mAb [PSH06-53]	HA722756	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
Relative expression	NCAM1 / CD56 Recombinant Rabbit mAb [PSH06-53] - BSA and Azide free	HA751105	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
	NCAM1 / CD56 Rabbit pAb	R1204-1	Human, Mouse, Rat	WB,IF-Cell
	NCAM1 / CD56 Recombinant Rabbit mAb	IRS026	Human	mIHC
Relative expression	NCAM1 / CD56 Recombinant Rabbit mAb [PSH06-82]	HA722755	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Cell,FC,IF-Tissue
Relative expression	NCAM1 / CD56 Recombinant Mouse mAb [PDH0-07]	HA601126	Human	IHC-P,WB
	NCAM1 / CD56 Recombinant Rabbit mAb [JF1021] - BSA and Azide free	HA750346	Human, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
Relative expression	NCAM1 / CD56 Recombinant Rabbit mAb [PSH06-82] - BSA and Azide free	HA751104	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Cell,FC,IF-Tissue
	NCAM1 / CD56 Recombinant Rabbit mAb [JF1021]	ET1702-43	Human, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP,mIHC
Relative expression	NCAM1 / CD56 Recombinant Rabbit mAb [PSH08-89]	HA723041	Human	WB,IHC-P,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	Nicastrin Recombinant Rabbit mAb [JG68-33]	ET7108-38	Human, Mouse, Rat	WB
KD	P4HB Recombinant Rabbit mAb [JE54-98]	ET7110-92	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC
KD	PDIA3 / ERp57 Recombinant Rabbit mAb [JM11-44]	ET1703-68	Human	WB,IP,IHC-P
KD	PDIA3 / ERp57 Recombinant Rabbit mAb [PSH08-79] - BSA and Azide free	HA751244	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell
KD	PDIA3 / ERp57 Recombinant Rabbit mAb [PSH08-79]	HA723031	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell
	Peroxiredoxin 2 Mouse mAb [7F4]	EM1701-71	Human	WB,IF-Cell,FC,IHC-P
KD	Peroxiredoxin 2 Rabbit pAb	R1510-37	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	Peroxiredoxin 2 Recombinant Mouse mAb [7F2-R] - BSA and Azide free	HA610123	Human	WB
KD	Peroxiredoxin 2 Recombinant Rabbit mAb [JJ090-3]	ET1701-61	Human, Mouse, Rat	WB,IF-Tissue,IHC-P
	Peroxiredoxin 2 Recombinant Mouse mAb [7F5-R] - BSA and Azide free	HA610082	Human	WB,IF-Cell,IHC-P
KD	Peroxiredoxin 2 Mouse mAb [7F5]	EM1701-72	Human, Mouse, Rat	WB,IHC-P,FC
	Peroxiredoxin 2 Mouse mAb [7F2]	EM1701-70	Human	WB,IHC-P,FC,IF-Cell
KD	Peroxiredoxin 2 Recombinant Mouse mAb [7F2-R]	HA601248	Human	WB
	Peroxiredoxin 2 Recombinant Mouse mAb [7F5-R]	HA601204	Human	WB,IF-Cell,IHC-P
	PGK1 Recombinant Rabbit mAb [ST49-07]	ET1609-63	Human, Mouse, Rat	WB
	Phospho-EGFR (S695) Recombinant Rabbit mAb [SN07-36]	ET1611-50	Human	WB,IF-Cell,IF-Tissue
Cell treatment	Phospho-EGFR (T669) Recombinant Rabbit mAb [PSH10-03]	HA723172	Human, Mouse, Rat	WB,IP
Cell treatment	Phospho-EGFR (T669) Recombinant Rabbit mAb [PSH10-03] - BSA and Azide free	HA751333	Human, Mouse, Rat	WB,IP
Cell treatment	Phospho-EGFR (Y1045) Recombinant Rabbit mAb [PSH03-28]	HA721970	Human	WB,IF-Cell,FC
Cell treatment	Phospho-EGFR (Y1045) Recombinant Rabbit mAb [PSH03-28] - BSA and Azide free	HA750868	Human	WB,IF-Cell,FC
Cell treatment	Phospho-EGFR (Y1068) Recombinant Rabbit mAb [SD2055] - BSA and Azide free	HA750290	Human	WB,IF-Cell,IF-Tissue,IP,FC
Cell treatment	Phospho-EGFR (Y1068) Recombinant Rabbit mAb [SD2055]	ET1612-30	Human	WB,IF-Cell,IF-Tissue,IP,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-EGFR (Y1092) Recombinant Rabbit mAb [SJ0194]	ET1606-44	Human, Mouse	WB,IF-Tissue,IHC-P
	Phospho-EGFR (Y1148) Rabbit pAb	HA500580	Human	WB,IF-Cell
Cell treatment	Phospho-EGFR (Y1173) Recombinant Rabbit mAb [SC57-04]	ET1610-4	Human, Mouse	WB,IF-Cell,IHC-P,IP
Cell treatment	Phospho-EGFR (Y845) Recombinant Rabbit mAb [JE46-77]	HA721488	Human	WB
Cell treatment	Phospho-Ezrin (T567)/Radixin (T564)/Moesin (T558) Recombinant Rabbit mAb [PSH05-73]	HA722440	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-Ezrin (T567)/Radixin (T564)/Moesin (T558) Recombinant Rabbit mAb [PSH05-73] - BSA and Azide free	HA751023	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Phospho-HER2 / ErbB2 (T686) Mouse mAb [4G3]	RT1275	Human, Mouse, Rat	WB,IP
Cell treatment	Phospho-HER2 / ErbB2 (Y1139) Recombinant Rabbit mAb [JE60-90]	HA722201	Human	WB,IF-Cell,FC
Cell treatment	Phospho-HER2 / ErbB2 (Y1221 + Y1222) Recombinant Rabbit mAb [JE44-12]	HA721433	Human	WB
Cell treatment	Phospho-HER2 / ErbB2 (Y1248) Recombinant Rabbit mAb [PSH04-03]	HA722084	Human	WB,IF-Cell
Cell treatment	Phospho-HER2 / ErbB2 (Y1248) Recombinant Rabbit mAb [PSH04-03] - BSA and Azide free	HA750926	Human	WB,IF-Cell
Cell treatment	Phospho-Hsp27 (S78) Recombinant Rabbit mAb [JJ08-70]	ET1701-19	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Phospho-Hsp27 (S82) Recombinant Rabbit mAb [SN06-73]	ET1611-16	Human, Mouse, Rat	WB,IP,IF-Cell,FC
Cell treatment	Phospho-Rac1+Cdc42 (S71) Recombinant Rabbit mAb [JJ086-9]	ET1701-40	Human	WB,IF-Cell,IHC-P
Relative expression	PKM1 Recombinant Rabbit mAb [PSH04-37] - BSA and Azide free	HA750946	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Relative expression	PKM1 Recombinant Rabbit mAb [PSH04-37]	HA722120	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	PKM1/2 Recombinant Rabbit mAb [JE07-34]	HA722524	Human, Mouse, Rat	WB,IP
KO	PKM2 Rabbit pAb	ER1901-90	Human, Mouse	WB,FC
KO	PKM2 Rabbit pAb	ER1802-70	Human	WB,IF-Cell,IHC-P
Relative expression	PKM2 Recombinant Rabbit mAb [PSH11-74]	HA723370	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC,IP
Relative expression	PKM2 Recombinant Rabbit mAb [PSH11-74] - BSA and Azide free	HA751429	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC,IP
Relative expression	PKM2 Rabbit pAb	R1603-5	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	PSMA Mouse mAb [A5G5]	HA600045	Human	WB,IHC-P
Relative expression	PSMA Mouse mAb [A5F12]	HA600046	Human	WB,IHC-P
Relative expression	PSMA Mouse mAb [A9H10]	HA601163	Human	WB,IHC-P,IF-Cell
Relative expression	PSMA Recombinant Rabbit mAb [PD00-72]	HA721721	Human	WB,IHC,IF-Tissue
Relative expression, KD	PSMA Mouse mAb [A10A2]	HA601169	Human	WB,IHC-P,IF-Cell
Relative expression	PSMA Recombinant Rabbit mAb [JE42-45]	HA721593	Human	WB,IHC-P
Relative expression	PSMA Mouse mAb [A9H9]	HA601164	Human	WB,IHC-P,IF-Cell
Relative expression	PSMA Rabbit pAb	ER1802-10	Human, Mouse, Rat	WB
	Rab10 Mouse mAb [F8-E2]	EM1709-43	Human, Mouse	WB,IF-Cell
KD	RAB10 Recombinant Rabbit mAb [PSH03-24]	HA721966	Human, Mouse, Rat	WB
KD	RAB10 Recombinant Rabbit mAb [PSH03-24] - BSA and Azide free	HA750864	Human, Mouse, Rat	WB
	Rab 11B Rabbit pAb	ER64710	Human, Mouse, Rat	WB
	Rab13 Mouse mAb	EM1707-79	Human	IF-Cell,IHC-P,FC
	RAB14 Rabbit pAb	ER62485	Human, Mouse, Rat	WB
	RAB15 Rabbit pAb	ER61095	Human, Mouse, Rat	WB
	Rab 18 Rabbit pAb	ER64711	Human, Mouse, Rat	WB
	RAB19 Rabbit pAb	ER61097	Human, Rat, Mouse	WB
Relative expression, KD	Rab1A Recombinant Rabbit mAb [PSH03-62]	HA722021	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	Rab1A Rabbit pAb	ER61098	Human, Rat, Mouse	WB,IF-Cell,IHC-P,FC
Relative expression, KD	Rab1A Recombinant Rabbit mAb [PSH03-62] - BSA and Azide free	HA750896	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	Rab1b Mouse mAb [G2-A12]	EM1708-26	Human	WB,IHC-P,FC
	Rab 20 Rabbit pAb	ER64712	Human	WB,IHC-P
	RAB21 Rabbit pAb	ER61099	Human, Rat, Mouse	WB
KD	RAB22A Recombinant Rabbit mAb [JE66-08]	HA721387	Human, Mouse, Rat	WB,IF-Cell
	RAB23 Rabbit pAb	ER61101	Human	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Rab 24 Rabbit pAb	ER64713	Human, Mouse, Rat	IHC-P
	Rab 26 Rabbit pAb	ER65425	Human, Mouse, Rat	WB,IHC-P
	RAB27A Mouse mAb [9D-7C]	EM1706-32	Human	IHC-P
	RAB28 Rabbit pAb	ER61102	Human, Rat, Mouse	WB
	Rab 2A Rabbit pAb	ER64714	Human, Mouse, Rat	WB,IHC-P
	Rab 2B Rabbit pAb	ER64715	Human, Mouse, Rat	WB
	RAB30 Rabbit pAb	ER61103	Human, Mouse	WB
	Rab 31 Rabbit pAb	ER64718	Human, Mouse, Rat	WB
	RAB32 Rabbit pAb	ER61104	Human, Mouse	WB
	Rab 34 Rabbit pAb	ER64720	Human	WB,IHC-P
	Rab 35 Rabbit pAb	ER64721	Human, Mouse, Rat	WB
	RAB36 Rabbit pAb	ER61106	Human, Mouse	WB
	Rab 37 Rabbit pAb	ER64722	Human	WB,IHC-P
	Rab 38 Rabbit pAb	ER64723	Human, Mouse, Rat	WB
	Rab3D Rabbit pAb	HA500165	Human, Mouse, Rat	WB,IHC-P,FC
	Rab 40B Rabbit pAb	ER64724	Human	WB,IHC-P
	Rab4A Recombinant Rabbit mAb [SR27-08]	ET1602-15	Human, Rat, Monkey, Mouse	WB,IF-Cell,IF-Tissue,IP
KD	Rab5 Recombinant Rabbit mAb [ST47-04] - BSA and Azide free	HA750167	Human, Mouse, Rat	WB,IHC-P,IF-Cell
KD	Rab5 Recombinant Rabbit mAb [ST47-04]	ET1609-27	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Rab5a Mouse mAb [C11-H5]	EM1708-08	Human	WB,IF-Cell,IHC-P,FC
	RAB5B Rabbit pAb	ER61113	Human, Mouse	WB
KD	RAB5C Recombinant Rabbit mAb [JE63-14]	HA721012	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Rab 6A Rabbit pAb	ER64727	Human, Mouse, Rat	WB,IHC-P
	Rab6b Mouse mAb [E4-D12]	EM1707-48	Human, Rat	WB,IF-Cell,IHC-P,FC
	Rab 6C Rabbit pAb	ER64728	Human, Mouse, Rat	WB
KD	RAB7 Recombinant Rabbit mAb [SN202-03] - BSA and Azide free	HA750279	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	RAB7 Recombinant Rabbit mAb [SN202-03]	ET1611-96	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	RAB7 Rabbit pAb	ER1802-53	Human, Mouse, Rat	WB,IHC-P
	Rab 7L1 Rabbit pAb	ER64729	Human, Mouse, Rat	WB
KD	RAB8A Recombinant Rabbit mAb [JE64-79]	HA721127	Human, Mouse, Rat	WB,FC,IF-Cell,IP
	RAB8B Rabbit pAb	ER61114	Human, Mouse, Rat	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	Rab9 Recombinant Rabbit mAb [SR45-05]	ET1602-29	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	RAB9B Rabbit pAb	ER61115	Human, Mouse	WB
	RAC1 Recombinant Rabbit mAb [PSH05-09] - BSA and Azide free	HA750982	Human, Rat	WB,IF-Cell
	RAC1 Recombinant Rabbit mAb [PSH05-09]	HA722225	Human, Rat	WB,IF-Cell
	Rac1-2-3 Recombinant Rabbit mAb [JM11-29]	ET1703-80	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	RAC1/2/3 Rabbit pAb	HA500424	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	Radixin Recombinant Rabbit mAb [SC06-22]	ET1610-41	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	RAP1A + RAP1B Recombinant Rabbit mAb [JE35-61]	HA721335	Human, Mouse	WB,FC
	RAP1B Rabbit pAb	ER61929	Human, Rat, Mouse	WB
	RB33B Rabbit pAb	ER61105	Human, Mouse	WB
	RB39A Rabbit pAb	ER61107	Human, Mouse	WB
	RB39B Rabbit pAb	ER61108	Human, Mouse	WB
	RB40A Rabbit pAb	ER61109	Human	WB
	RB40C Rabbit pAb	ER61111	Human, Mouse	WB
	RB40L Rabbit pAb	ER61110	Human	WB
	RIM 2 Mouse mAb [3G2]	RT1539	Human, Rat	WB,IP
	RPS3 Recombinant Rabbit mAb [SA46-08]	ET1601-27	Human, Mouse, Rat	WB,IHC-P
KD	Syntenin Recombinant Rabbit mAb [JE40-72]	ET7109-14	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
KD	TAPA1/CD81 Recombinant Rabbit mAb [SN206-01]	ET1611-87	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Transferrin Receptor 2 Recombinant Rabbit mAb [JG58-34]	ET7108-21	Human	WB,IHC-P
Relative expression	Transferrin Receptor (CD71) Recombinant Rabbit mAb [PSH13-56]	HA723545	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP,IF-Tissue
KD	Transferrin Receptor (CD71) Recombinant Rabbit mAb [JF0956]	ET1702-06	Human, Mouse, Rat	WB,IP,FC
KD	Transferrin Receptor (CD71) Recombinant Rabbit mAb [JF0956] - BSA and Azide free	HA750333	Human, Mouse, Rat	WB,IP,FC
	Transferrin Receptor (CD71) Rabbit pAb	ER1706-35	Human	IHC-P
	Triosephosphate isomerase Recombinant Rabbit mAb [PSH14-38]	HA723632	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
	Triosephosphate isomerase Recombinant Rabbit mAb [PSH14-38] - BSA and Azide free	HA751520	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Triosephosphate isomerase Rabbit pAb	HA500283	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
KD	TSG101 Recombinant Rabbit mAb [JJ0900]	ET1701-59	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	TSG101 Recombinant Rabbit mAb [JJ0900] - BSA and Azide free	HA750316	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC