

荧光直标抗体-产品清单

高级验证	产品名	货号	反应种属	验证应用
	iFluor™ 488 Conjugated 14-3-3 sigma Recombinant Rabbit mAb [SD2070]	HA723707F	Human	IF-Tissue
	FITC Conjugated 6X His tag Recombinant Rabbit mAb [PSH07-10]	HA720222F	Species independent	IF-Cell,FC
	iFluor™ 594 Conjugated 6X His tag Recombinant Rabbit mAb [PSH07-10]	HA720221F	Species independent	IF-Cell,FC
	iFluor™ 488 Conjugated Actin Recombinant Rabbit mAb [JJ09-29]	HA720160F	Human, Mouse, Rat	IF-Cell,IF-Tissue,FC,WB
	FITC Conjugated Bcl-2 Recombinant Rabbit mAb [JF104-8]	HA720210F	Human, Mouse	FC
	iFluor™ 488 Conjugated Beta Catenin Recombinant Rabbit mAb [SA30-04]	HA720182F	Human, Mouse	IF,ICC
	iFluor™ 488 Conjugated Beta III Tubulin Mouse mAb [A8-D10]	HA600106F	Human, Mouse, Rat	IF-Cell,FC,IF-Tissue
	iFluor™ 594 Conjugated Beta III Tubulin Mouse mAb [A8-D10]	HA600133F	Human, Rat	IF,ICC,FC
	iFluor™ 594 Conjugated beta Tubulin Recombinant Mouse mAb [A1-A4-R]	HA600110F	Human, Mouse	IF-Cell,IF-Tissue,FC
	iFluor™ 647 Conjugated BrdU Recombinant Rabbit mAb [PSH0-18]	HA720188F	Species independent	IF-Tissue,FC
	iFluor™ 488 Conjugated BrdU Recombinant Rabbit mAb [PSH0-18]	HA720186F	Species independent	IF-Cell,IF-Tissue,FC
	iFluor™ 594 Conjugated BrdU Recombinant Rabbit mAb [PSH0-18]	HA720187F	Species independent	IF-Tissue,FC
	iFluor™ 488 Conjugated Calretinin Recombinant Rabbit mAb [JM12-93]	HA720153F	Human	IF-Cell,FC
	iFluor™ 488 Conjugated CCR7 Recombinant Rabbit mAb [SR36-04]	HA720157F	Human, Mouse, Rat	IF-Cell,IF-Tissue
	APC Conjugated CD127 Recombinant Mouse mAb [PSH12-76]	HA601536F	Human	FC
	FITC Conjugated CD15 Recombinant Rabbit mAb [PD00-42]	HA720195F	Human	FC
	FITC Conjugated CD19 Recombinant Rabbit mAb [PSH13-44]	HA723762F	Human	FC
	FITC Conjugated CD19 Recombinant Mouse mAb [PSH14-97]	HA601518F	Human	FC
	PE Conjugated CD19 Recombinant Rabbit mAb [PSH13-44]	HA723763F	Human	FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	APC Conjugated CD19 Recombinant Rabbit mAb [PSH13-44]	HA723761F	Human	FC
	APC Conjugated CD20 Recombinant Rabbit mAb [PSH08-12]	HA723838F	Human	FC
	PE Conjugated CD20 Recombinant Rabbit mAb [PSH08-12]	HA723839F	Human	FC
	APC Conjugated CD25 / IL-2 Receptor alpha Recombinant Mouse mAb [PSH14-46]	HA601530F	Human	FC
	APC Conjugated CD27 Recombinant Rabbit mAb [PSH08-90]	HA723828F	Human	FC
	PE Conjugated CD45RO Recombinant Mouse mAb [PSH13-94]	HA601531F	Human	FC
	iFluor™ 488 Conjugated CD47 Mouse mAb [1B1]	HA600104F	Human	FC
	APC Conjugated CD94 Recombinant Mouse mAb [PSH14-52]	HA601529F	Human	FC
	PE Conjugated CD94 Recombinant Mouse mAb [PSH14-52]	HA601539F	Human	FC
	iFluor™ 488 Conjugated Collagen VI Recombinant Rabbit mAb [SD83-03]	HA720155F	Human, Mouse, Rat	IF-Cell,IF-Tissue
	iFluor™ 488 Conjugated Cytokeratin 1 Recombinant Rabbit mAb [SN72-08]	HA720127F	Human, Mouse, Rat	IF-Tissue
	iFluor™ 594 Conjugated Cytokeratin 1 Recombinant Rabbit mAb [SN72-08]	HA720124F	Human, Mouse, Rat	IF-Tissue
	iFluor™ 647 Conjugated Cytokeratin 1 Recombinant Rabbit mAb [SN72-08]	HA720141F	Human, Mouse, Rat	IF-Tissue
	iFluor™ 488 Conjugated Cytokeratin 10 Recombinant Rabbit mAb [ST05-43]	HA720134F	Human, Rat	IF-Tissue
	iFluor™ 594 Conjugated Cytokeratin 10 Recombinant Rabbit mAb [ST05-43]	HA720119F	Human, Rat	IF-Tissue
	iFluor™ 647 Conjugated Cytokeratin 10 Recombinant Rabbit mAb [ST05-43]	HA720146F	Human, Rat, Mouse	IF-Tissue
	iFluor™ 594 Conjugated Cytokeratin 14 Recombinant Rabbit mAb [SC65-06]	HA720117F	Rat, Human, Mouse	IF-Tissue
	iFluor™ 647 Conjugated Cytokeratin 14 Recombinant Rabbit mAb [SC65-06]	HA720147F	Human, Mouse, Rat	IF-Cell,IF-Tissue
	iFluor™ 488 Conjugated Cytokeratin 14 Recombinant Rabbit mAb [SC65-06]	HA720135F	Mouse, Rat, Human	IF-Tissue,IF-Cell
	iFluor™ 488 Conjugated Cytokeratin 15 Recombinant Rabbit mAb [ST04-85]	HA720136F	Human, Mouse, Rat	IF-Tissue,IF-Cell,FC
	iFluor™ 594 Conjugated Cytokeratin 15 Recombinant Rabbit mAb [ST04-85]	HA720122F	Human, Mouse, Rat	IF-Tissue
	iFluor™ 647 Conjugated Cytokeratin 15 Recombinant Rabbit mAb [ST04-85]	HA720148F	Human, Mouse, Rat	IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	iFluor™ 647 Conjugated Cytokeratin 16 Recombinant Rabbit mAb [SC52-09]	HA720149F	Human	IF-Tissue
	iFluor™ 488 Conjugated Cytokeratin 16 Recombinant Rabbit mAb [SC52-09]	HA720137F	Human	IF-Tissue,IF-Cell,FC
	iFluor™ 594 Conjugated Cytokeratin 16 Recombinant Rabbit mAb [SC52-09]	HA720114F	Human	IF-Cell,IF-Tissue
	iFluor™ 594 Conjugated Cytokeratin 17 Recombinant Rabbit mAb [SR45-06]	HA720116F	Human, Mouse	IF-Tissue,IF-Cell
	iFluor™ 647 Conjugated Cytokeratin 17 Recombinant Rabbit mAb [SR45-06]	HA720150F	Human, Mouse, Rat	IF-Cell,IF-Tissue
	iFluor™ 488 Conjugated Cytokeratin 17 Recombinant Rabbit mAb [SR45-06]	HA720138F	Human, Mouse	IF-Tissue,IF-Cell
	iFluor™ 488 Conjugated Cytokeratin 18 Recombinant Rabbit mAb [SZ80-07]	HA720139F	Human, Mouse, Rat	IF-Tissue,IF-Cell
	iFluor™ 647 Conjugated Cytokeratin 18 Recombinant Rabbit mAb [SZ80-07]	HA720151F	Human, Mouse, Rat	IF-Tissue
	iFluor™ 594 Conjugated Cytokeratin 19 Recombinant Rabbit mAb [SA30-06]	HA720115F	Human, Mouse	IF-Cell,IF-Tissue,FC
	iFluor™ 647 Conjugated Cytokeratin 19 Recombinant Rabbit mAb [SA30-06]	HA720152F	Human, Mouse	IF-Cell,IF-Tissue,FC
	iFluor™ 488 Conjugated Cytokeratin 19 Recombinant Rabbit mAb [SA30-06]	HA720140F	Human, Mouse	IF-Cell,IF-Tissue,FC
	iFluor™ 488 Conjugated Cytokeratin 20 Recombinant Rabbit mAb [SA35-03]	HA720189F	Human	IF-Cell,IF-Tissue
	iFluor™ 647 Conjugated Cytokeratin 4 Recombinant Rabbit mAb [SN74-03]	HA720142F	Human	IF-Cell,IF-Tissue
	iFluor™ 488 Conjugated Cytokeratin 4 Recombinant Rabbit mAb [SN74-03]	HA720128F	Human	IF-Tissue,IF-Cell
	iFluor™ 594 Conjugated Cytokeratin 4 Recombinant Rabbit mAb [SN74-03]	HA720123F	Human	IF-Tissue
	iFluor™ 647 Conjugated Cytokeratin 5 Recombinant Rabbit mAb [SC62-04]	HA720126F	Human	IF-Cell,IF-Tissue
	iFluor™ 594 Conjugated Cytokeratin 5 Recombinant Rabbit mAb [SC62-04]	HA720121F	Human	IF-Tissue
	iFluor™ 488 Conjugated Cytokeratin 5 Recombinant Rabbit mAb [SC62-04]	HA720129F	Human	IF-Tissue,IF-Cell
	iFluor™ 647 Conjugated Cytokeratin 6 Recombinant Rabbit mAb [SN71-07]	HA720143F	Human, Rat	IF-Cell,IF-Tissue
	iFluor™ 488 Conjugated Cytokeratin 6 Recombinant Rabbit mAb [SN71-07]	HA720130F	Human, Rat	IF-Cell
	iFluor™ 488 Conjugated Cytokeratin 7 Recombinant Rabbit mAb [ST50-05]	HA720131F	Human, Mouse	IF-Tissue,IF-Cell,FC
	iFluor™ 594 Conjugated Cytokeratin 7 Recombinant Rabbit mAb [ST50-05]	HA720120F	Human, Mouse	IF-Cell,IF-Tissue,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	iFluor™ 647 Conjugated Cytokeratin 7 Recombinant Rabbit mAb [ST50-05]	HA720144F	Human, Mouse	IF-Cell,IF-Tissue,FC
	iFluor™ 594 Conjugated Cytokeratin 8 Recombinant Rabbit mAb [SU0338]	HA720118F	Human, Mouse	IF-Cell,IF-Tissue,FC
	iFluor™ 647 Conjugated Cytokeratin 8 Recombinant Rabbit mAb [SU0338]	HA720145F	Human, Mouse	IF-Cell,IF-Tissue
	iFluor™ 488 Conjugated Cytokeratin 8 Recombinant Rabbit mAb [SU0338]	HA720132F	Human	IF-Tissue,IF-Cell
	iFluor™ 488 Conjugated Cytokeratin 9 Recombinant Rabbit mAb [SD201-09]	HA720133F	Human, Mouse	IF-Tissue
	iFluor™ 594 Conjugated DYKDDDDK Tag (FLAG) Recombinant Mouse mAb [A2-A4-R]	HA600138F	Species independent	IF-Cell
	FITC Conjugated DYKDDDDK Tag (FLAG) Recombinant Mouse mAb [A2-A4-R]	HA600139F	Species independent	IF-Cell
	iFluor™ 594 Conjugated E-Cadherin Recombinant Rabbit mAb [SY0287]	HA720164F	Human	IF-Tissue
	iFluor™ 488 Conjugated E-Cadherin Recombinant Rabbit mAb [SY0287]	HA720159F	Human	IF-Cell,FC
	iFluor™ 647 Conjugated EpCAM Recombinant Rabbit mAb [PS01-69]	HA720185F	Human	FC,IF-Cell
	iFluor™ 488 Conjugated FOXO1A Recombinant Rabbit mAb [SU33-01]	HA720161F	Human, Mouse, Rat	IF-Cell,FC
	iFluor™ 647 Conjugated GFAP Mouse mAb [1-D4]	HA600103F	Human, Rat	IF-Tissue
	iFluor™ 594 Conjugated GFAP Mouse mAb [1-D4]	HA600102F	Human, Rat	IF-Tissue
	iFluor™ 488 Conjugated GFAP Mouse mAb [1-D4]	HA600101F	Human, Rat	IF-Tissue,IHC-Fr
	iFluor™ 647 Conjugated Glucose Transporter GLUT1 Recombinant Rabbit mAb [SA0377]	HA720170F	Human, Mouse	IF-Tissue
	iFluor™ 488 Conjugated Glucose Transporter GLUT1 Recombinant Rabbit mAb [SA0377]	HA720154F	Human, Mouse, Rat	IF-Cell,IF-Tissue
	FITC Conjugated GST Recombinant Rabbit mAb [SN70-05]	HA720225F	Species independent	IF-Cell
	FITC Conjugated HA tag Recombinant Rabbit mAb [PSH10-66]	HA720223F	Species independent	IF-Cell,FC
	iFluor™ 594 Conjugated HA tag Recombinant Rabbit mAb [PSH10-66]	HA720224F	Species independent	IF-Cell,FC
	iFluor™ 594 Conjugated Heme Oxygenase 1 (HO-1) Recombinant Rabbit mAb [SP08-07]	HA720169F	Human	IF-Tissue
	iFluor™ 647 Conjugated HER3 Recombinant Rabbit mAb [R5]	HA720181F	Human	FC
	FITC Conjugated Human CD10 Recombinant Rabbit mAb [SN75-07]	HA720209F	Human	FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	FITC Conjugated Human CD11b Recombinant Mouse mAb [A2-9-R]	HA600130F	Human	FC
	FITC Conjugated Human CD14 Recombinant Rabbit mAb [SC69-02]	HA720194F	Human	FC
	iFluor™ 647 Conjugated Human CD14 Recombinant Rabbit mAb [SC69-02]	HA720184F	Human	FC
	PE Conjugated Human CD16 Recombinant Mouse mAb [PSH04-95]	HA600129F	Human	FC
	FITC Conjugated Human CD200 Recombinant Rabbit mAb [PSH01-48]	HA720208F	Human	FC
	FITC Conjugated Human CD200R1 Recombinant Rabbit mAb [PSH05-02]	HA720212F	Human	FC
	FITC Conjugated Human CD3 Recombinant Mouse mAb [PSH04-92]	HA600116F	Human	FC
	PE Conjugated Human CD3 Recombinant Mouse mAb [PSH04-92]	HA600117F	Human	FC
	Percp-Cy5.5 Conjugated Human CD3 Recombinant Mouse mAb [PSH04-92]	HA601548F	Human	FC
	APC Conjugated Human CD3 Recombinant Mouse mAb [PSH04-92]	HA600118F	Human	FC
	FITC Conjugated Human CD33 Recombinant Rabbit mAb [PSH04-96]	HA720197F	Human	FC
	FITC Conjugated Human CD38 Recombinant Rabbit mAb [SN07-20]	HA720196F	Human	FC
	PE-Cy7 Conjugated Human CD4 Recombinant Mouse mAb [PSH04-93]	HA601542F	Human	FC
	APC Conjugated Human CD4 Recombinant Mouse mAb [PSH04-93]	HA600121F	Human	FC
	PE Conjugated Human CD4 Recombinant Mouse mAb [PSH04-93]	HA600120F	Human	FC
	FITC Conjugated Human CD4 Recombinant Mouse mAb [PSH04-93]	HA600119F	Human	FC
	PE Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600126F	Human	FC
	APC Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600127F	Human	FC
	FITC Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600125F	Human	FC
	Percp-Cy5.5 Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600128F	Human	FC
	FITC Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600112F	Human	FC
	iFluor™ 555 Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA601556F	Human	FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	iFluor™ 488 Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600109F	Human	FC
	FITC Conjugated Human CD45RA Recombinant Mouse mAb [PSH04-97]	HA600131F	Human	FC
	PE Conjugated Human CD45RA Recombinant Mouse mAb [PSH04-97]	HA601525F	Human	FC
	FITC Conjugated Human CD5 Recombinant Rabbit mAb [PSH05-01]	HA720211F	Human	FC
	FITC Conjugated Human CD58 Recombinant Rabbit mAb [PSH02-74]	HA720198F	Human	FC
	FITC Conjugated Human CD68 Recombinant Mouse mAb [A3C3-R]	HA600132F	Human	FC
	APC Conjugated Human CD8 alpha Recombinant Mouse mAb [PSH04-94]	HA600124F	Human	FC
	FITC Conjugated Human CD8 alpha Recombinant Mouse mAb [PSH04-94]	HA600122F	Human	FC
	APC-Cy7 Conjugated Human CD8 alpha Recombinant Mouse mAb [PSH04-94]	HA601541F	Human	FC
	PE Conjugated Human CD8 alpha Recombinant Mouse mAb [PSH04-94]	HA600123F	Human	FC
	FITC Conjugated Human IgD Recombinant Mouse mAb [PSH14-48]	HA601526F	Human	FC
	PE Conjugated Human IgE Recombinant Mouse mAb [PSH14-50]	HA601527F	Human	FC
	PE Conjugated Human NCAM1 / CD56 Recombinant Rabbit mAb [PSH06-53]	HA720215F	Human	FC
	APC Conjugated Human NCAM1 / CD56 Recombinant Rabbit mAb [PSH06-53]	HA723757F	Human	FC
	iFluor™ 594 Conjugated Iba1 Recombinant Rabbit mAb [JM36-62]	HA723821F	Mouse	IHC-Fr
	iFluor™ 488 Conjugated Iba1 Recombinant Rabbit mAb [JM36-62]	HA720158F	Human, Mouse, Rat	IF-Cell,FC
	iFluor™ 488 Conjugated Ki67 Recombinant Rabbit mAb [SR00-02]	HA720162F	Human	IF-Cell,IF-Tissue,FC
	iFluor™ 647 Conjugated Ki67 Recombinant Rabbit mAb [SR00-02]	HA720163F	Human, Mouse, Rat	IF-Cell,IF-Tissue,FC,IHC-Fr
	iFluor™ 488 Conjugated Lamin B1 Recombinant Rabbit mAb [SI17-07]	HA720171F	Human, Mouse, Rat	IF-Tissue
	iFluor™ 647 Conjugated Lamin B1 Recombinant Rabbit mAb [SI17-07]	HA720172F	Human, Mouse, Rat	IF-Tissue
	PE Conjugated Ly6g Recombinant Rabbit mAb [PSH14-95]	HA723853F	Mouse	FC
	FITC Conjugated Mannose Receptor(CD206) Recombinant Rabbit mAb [JF0953]	HA720214F	Human	FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	PE Conjugated Mouse CD16/32 Recombinant Rabbit mAb [PSH06-60]	HA720216F	Mouse	FC
	APC Conjugated Anti-Mouse CD3 Recombinant Antibody [PSH04-98] - Rabbit IgG (Chimeric)	HA720201F	Mouse	FC
	PE Conjugated Anti-Mouse CD3 Recombinant Antibody [PSH04-98] - Rabbit IgG (Chimeric)	HA720200F	Mouse	FC
	FITC Conjugated Anti-Mouse CD3 Recombinant Antibody [PSH04-98] - Rabbit IgG (Chimeric)	HA720199F	Mouse	FC
	APC Conjugated Anti-Mouse CD4 Recombinant Antibody [PSH04-99] - Rabbit IgG (Chimeric)	HA720204F	Mouse	FC
	FITC Conjugated Anti-Mouse CD4 Recombinant Antibody [PSH04-99] - Rabbit IgG (Chimeric)	HA720202F	Mouse	FC
	PE Conjugated Anti-Mouse CD4 Recombinant Antibody [PSH04-99] - Rabbit IgG (Chimeric)	HA720203F	Mouse	FC
	FITC Conjugated Anti-Mouse CD8 alpha Recombinant Antibody [PSH05-00] - Rabbit IgG (Chimeric)	HA720205F	Mouse	FC
	PE Conjugated Anti-Mouse CD8 alpha Recombinant Antibody [PSH05-00] - Rabbit IgG (Chimeric)	HA720206F	Mouse	FC
	APC Conjugated Anti-Mouse CD8 alpha Recombinant Antibody [PSH05-00] - Rabbit IgG (Chimeric)	HA720207F	Mouse	FC
Cell treatment	FITC Conjugated Mouse PD-L1 Recombinant Rabbit mAb [PSH05-04]	HA720213F	Mouse	FC
	iFluor™ 647 Conjugated MUC1 Mouse mAb [PSH0-48]	HA600107F	Human	FC
	iFluor™ 647 Conjugated MUC1 Mouse mAb [PSH0-49]	HA600108F	Human	FC
	FITC Conjugated Myc tag Recombinant Mouse mAb [A3C8-R]	HA600136F	Species independent	IF-Cell
	iFluor™ 594 Conjugated Myc tag Recombinant Mouse mAb [A3C8-R]	HA600137F	Species independent	IF-Cell
	iFluor™ 488 Conjugated NeuN Recombinant Rabbit mAb [SR45-07]	HA720168F	Human, Mouse, Rat	IF-Tissue
	iFluor™ 594 Conjugated NeuN Recombinant Rabbit mAb [SR45-07]	HA720167F	Human, Mouse, Rat	IF-Tissue,IHC-Fr
	FITC Conjugated NK1.1 Recombinant Rabbit mAb [PSH08-46]	HA720219F	Mouse	FC
	iFluor™ 488 Conjugated OGT Recombinant Rabbit mAb [JB44-39]	HA720156F	Human, Mouse, Rat	IF-Cell
	iFluor™ 488 Conjugated p75 NGF Receptor Recombinant Rabbit mAb [SA39-02]	HA720190F	Human, Mouse, Rat	IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	iFluor™ 555 Conjugated pan Cytokeratin Recombinant Mouse mAb [PDH09-10]	HA600134F	Human	IF-Cell
	iFluor™ 488 Conjugated pan Cytokeratin Recombinant Mouse mAb [PDH0-10]	HA600111F	Human	IF-Tissue
	iFluor™ 488 Conjugated pan Cytokeratin Recombinant Mouse mAb [PDH09-10]	HA600135F	Human	IF-Cell,IF-Tissue
	iFluor™ 647 Conjugated pan Cytokeratin Recombinant Mouse mAb [PD00-15]	HA600105F	Human	IF-Tissue,FC
	APC Conjugated Recombinant Mouse IgG1 [PSH04-91] - Isotype Control	HA600115F		FC
	PE Conjugated Recombinant Mouse IgG1 [PSH04-91] - Isotype Control	HA600114F		FC
	FITC Conjugated Recombinant Mouse IgG1 [PSH04-91] - Isotype Control	HA600113F		FC
	PE Conjugated Recombinant Rabbit IgG [PSH04-42] - Isotype control	HA720192F		FC
	APC Conjugated Recombinant Rabbit IgG [PSH04-42] - Isotype control	HA720193F		FC
	FITC Conjugated Recombinant Rabbit IgG [PSH04-42] - Isotype control	HA720191F		FC
	iFluor™ 594 Conjugated S100 beta Recombinant Rabbit mAb [SC57-02]	HA720125F	Human	IF-Tissue
	iFluor™ 647 Conjugated SATB2 Recombinant Rabbit mAb [JM52-44]	HA720173F	Human, Mouse, Rat	IF-Tissue,FC
	FITC Conjugated Siglec 8 Recombinant Mouse mAb [PSH14-96]	HA601517F	Human	FC
	iFluor™ 488 Conjugated Sodium Potassium ATPase Recombinant Rabbit mAb [ST0533]	HA720174F	Human, Mouse, Rat	IF-Tissue
	iFluor™ 594 Conjugated Sodium Potassium ATPase Recombinant Rabbit mAb [ST0533]	HA720175F	Human, Mouse, Rat	IF-Tissue
	iFluor™ 647 Conjugated Sodium Potassium ATPase Recombinant Rabbit mAb [ST0533]	HA720176F	Human, Mouse, Rat	IF-Tissue,FC
	iFluor™ 594 Conjugated SOX2 Recombinant Rabbit mAb [PO00-28]	HA720180F	Human, Mouse	IF-Cell,FC
	iFluor™ 488 Conjugated SOX2 Recombinant Rabbit mAb [PO00-28]	HA720179F	Human, Mouse	IF-Cell,FC
	iFluor™ 594 Conjugated Strep-Tag II Recombinant Rabbit mAb [PSH05-13]	HA720228F	Species independent	IF-Cell
	FITC Conjugated Strep-Tag II Recombinant Rabbit mAb [PSH05-13]	HA720229F	Species independent	IF-Cell
	iFluor™ 594 Conjugated V5 tag Recombinant Rabbit mAb [SY30-01]	HA720226F	Species independent	IF-Cell
	FITC Conjugated V5 tag Recombinant Rabbit mAb [SY30-01]	HA720227F	Species independent	IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	iFluor™ 488 Conjugated Vimentin Recombinant Rabbit mAb [SC60-05]	HA720166F	Human, Mouse, Rat	IF-Cell,IF-Tissue,FC
	iFluor™ 647 Conjugated Vimentin Recombinant Rabbit mAb [SC60-05]	HA720165F	Human, Mouse, Rat	IF-Cell,IF-Tissue,FC
	iFluor™ 488 Conjugated VISTA Recombinant Rabbit mAb [PSH0-65]	HA720183F	Human	FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层