

AD & 金属离子稳态失衡研究-产品清单

高级验证	产品名	货号	反应种属	验证应用
KD	ACO1 Recombinant Rabbit mAb [JE63-98]	HA722071	Human, Mouse, Rat	WB,FC
	ACO1 Rabbit pAb	R1512-16	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	ADAM10 Rabbit pAb	ER1706-97	Human	WB,IHC-P
KD	ADAM10 Recombinant Rabbit mAb [JM32-11]	ET1703-60	Human, Mouse, Rat	WB,IHC-P,IP
	ALOX12 Rabbit pAb	ER1903-58	Human, Mouse	WB
	AMID Recombinant Rabbit mAb [PSH15-60]	HA723767	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	AMID Rabbit pAb	ER62655	Human, Mouse, Monkey	WB,IHC-P,IF-Cell
	AMID Recombinant Rabbit mAb [PSH15-60] - BSA and Azide free	HA751571	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	ATP7A Rabbit pAb	ER62772	Human, Mouse, Rat	IHC-P
	BACE1 Recombinant Rabbit mAb [PSH05-19] - BSA and Azide free	HA750990	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue
	BACE1 Recombinant Rabbit mAb [PSH05-19]	HA722244	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue
	BACE1 Mouse mAb [7H4]	RT1066	Human, Mouse, Rat	WB,IP,IF,IHC-P
	BACE1 Recombinant Rabbit mAb [JE30-40]	HA721421	Human, Mouse, Rat	WB
KD	Cdk5 Recombinant Rabbit mAb [SR39-01]	ET1602-17	Human, Mouse, Rat	WB,IHC-P
	DMT1 Recombinant Rabbit mAb [PSH09-20] - BSA and Azide free	HA751274	Human	WB,IHC-P,FC
	DMT1 Rabbit pAb	ER1907-55	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	DMT1 Recombinant Rabbit mAb [PSH09-20]	HA723078	Human	WB,IHC-P,FC
KD	ERK1 Recombinant Rabbit mAb [SP00-74]	ET1604-32	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	ERK1 Mouse mAb [7-D6-E5]	M1505-9	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	ERK1 Recombinant Rabbit mAb [SP05-09]	ET1604-16	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Erk1/2 Mouse mAb [2G1]	RT1484	Human, Mouse, Rat	WB,IP,IF
	ERK1/2 Rabbit pAb	ER131011	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	ERK1/2 Recombinant Rabbit mAb [SA43-03] - BSA and Azide free	HA750019	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IP,FC,IHC-P,IHC-Fr

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	ERK1/2 Recombinant Rabbit mAb [SA43-03]	ET1601-29	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IP,FC,IHC-P,IHC-Fr
KD	ERK2 Mouse mAb [B10-8]	EM1901-53	Human, Mouse, Rat	WB,IHC-P,FC
KD	ERK2 Recombinant Mouse mAb [C3B12-R]	HA601201	Human, Mouse, Rat	WB
KD	ERK2 Mouse mAb [A10D2]	EM1901-54	Human, Mouse, Rat	WB,IHC-P,FC
KD	ERK2 Recombinant Mouse mAb [5-D2-R] - BSA and Azide free	HA610080	Human, Mouse, Rat	WB
KD	ERK2 Rabbit pAb	ER131218	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	ERK2 Recombinant Mouse mAb [C3B12-R] - BSA and Azide free	HA610079	Human, Mouse, Rat	WB
KD	ERK2 Recombinant Mouse mAb [5-D2-R]	HA601202	Human, Mouse, Rat	WB
	ERK2 Mouse mAb [5-D2]	M1407-3	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	ERK2 Mouse mAb [C3B12]	EM1901-55	Human, Mouse, Rat	WB,IHC-P,FC
KD	ERK2 Recombinant Rabbit mAb [SZ25-01]	ET1603-23	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
KD	GPX4 Recombinant Rabbit mAb [JU11-31] - BSA and Azide free	HA750455	Human, Mouse, Rat, Zebrafish, Pig, Cynomolgus monkey	WB,IHC-P,IF-Cell,FC
KD	GPX4 Recombinant Rabbit mAb [JU11-31]	ET1706-45	Human, Mouse, Rat, Zebrafish, Pig, Cynomolgus monkey	WB,IHC-P,IF-Cell,FC
	GPX4 Rabbit pAb	ER1803-15	Mouse, Rat	WB
	GSK3 alpha Recombinant Rabbit mAb [JE33-27]	HA722217	Human, Mouse, Rat, Monkey	WB,IHC-P
KD	GSK3 alpha Rabbit pAb	R1312-1	Human	WB,IF-Cell,IHC-P,FC
	GSK3 alpha/beta Rabbit pAb	ER63648	Human, Mouse, Rat	WB,IHC-P
KD	GSK3 beta Recombinant Rabbit mAb [SY28-03] - BSA and Azide free	HA750131	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IHC-Fr
	GSK3 beta Rabbit pAb	R1510-4	Human, Mouse, Rat	IF-Cell,IHC-P,FC
KD	GSK3 beta Recombinant Rabbit mAb [SY28-03]	ET1607-71	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IHC-Fr
	Human BACE1 Recombinant Rabbit mAb [PSH16-23] - BSA and Azide free (Detector)	HA725172	Human	ELISA(Det)
	Biotin Conjugated Human BACE1 Recombinant Rabbit mAb [PSH16-23] - Detector	HA725172B	Human	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human BACE1 Recombinant Rabbit mAb [PSH16-22] - BSA and Azide free (Capture)	HA725171	Human	ELISA(Cap)
	JNK1 Recombinant Rabbit mAb [JE63-04]	HA721899	Human, Mouse, Rat	WB,IHC-P
	JNK1 Rabbit pAb	R1309-1	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC
	JNK1+JNK2+JNK3 Recombinant Rabbit mAb [SA43-06] - BSA and Azide free	HA750018	Human, Mouse, Rat, Cow, Monkey, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP
	JNK1+JNK2+JNK3 Recombinant Rabbit mAb [SA43-06]	ET1601-28	Human, Mouse, Rat, Cow, Monkey, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP
	JNK1+JNK3 Recombinant Rabbit mAb [JJ08-86]	ET1701-47	Human, Mouse, Rat	WB,IP,IHC-P
	JNK2 Rabbit pAb	HA500110	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	JNK2 Recombinant Rabbit mAb [SC05-85]	ET1610-11	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
	JNK3 Recombinant Rabbit mAb [SD082-09]	ET1612-68	Human, Mouse, Rat	WB,IF-Tissue,IHC-P
	LRP-1 Mouse mAb [A4-F5]	M1211-4	Mouse	WB,IHC-P
KD	LRP1 Recombinant Rabbit mAb [SA0290]	ET1601-1	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,IP,IHC-Fr
Cell treatment	Phospho-Erk1 (T202)+Erk2 (T185) Recombinant Rabbit mAb [SZ2-4] - BSA and Azide free	HA750061	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
Cell treatment	Phospho-Erk1 (T202)+Erk2 (T185) Recombinant Rabbit mAb [SZ2-4]	ET1603-22	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
Cell treatment	Phospho-Erk1 (T202 + Y204) + Erk2 (T185 + Y187) Recombinant Rabbit mAb [SC58-01] - BSA and Azide free	HA750199	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP
Cell treatment	Phospho-Erk1 (T202 + Y204) + Erk2 (T185 + Y187) Recombinant Rabbit mAb [SC58-01]	ET1610-13	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP
Cell treatment	Phospho-GSK3 (alpha + beta) (Y216 + Y279) Recombinant Rabbit mAb [SY26-05] - BSA and Azide free	HA750127	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,IF-Tissue,FC
Cell treatment	Phospho-GSK3 (alpha + beta) (Y216 + Y279) Recombinant Rabbit mAb [SY26-05]	ET1607-54	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,IF-Tissue,FC
Cell treatment	Phospho-GSK3 alpha (S21) Recombinant Rabbit mAb [PSH09-11] - BSA and Azide free	HA751265	Human, Mouse	WB,IHC-P,Dot Blot
Cell treatment	Phospho-GSK3 alpha (S21) Recombinant Rabbit mAb [PSH09-11]	HA723069	Human, Mouse	WB,IHC-P,Dot Blot
Cell treatment	Phospho-GSK3 beta (S9) Recombinant Rabbit mAb [SY02-71]	ET1607-60	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
Cell treatment	Phospho-GSK3 beta (S9) Recombinant Rabbit mAb [SY02-71] - BSA and Azide free	HA750128	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-JNK1/2/3 (T183 + T183 + T221) Recombinant Rabbit mAb [ST500] - BSA and Azide free	HA750173	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
Cell treatment	Phospho-JNK1/2/3 (T183 + T183 + T221) Recombinant Rabbit mAb [ST500]	ET1609-42	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
	SLC31A1 / CTR1 Recombinant Rabbit mAb [JE56-93]	HA720005	Human, Mouse, Rat	WB
	SLC31A1 / CTR1 Recombinant Rabbit mAb [PSH09-67]	HA723132	Human, Mouse, Rat	WB
	SLC40A1 Mouse mAb [A10A1]	HA601178	Human	WB
	SLC40A1 Rabbit pAb	ER1916-80	Human	WB
	SP1 Recombinant Rabbit mAb [JF0950]	ET1702-02	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,ChIP,IP
	SP1 Recombinant Rabbit mAb [JF0950] - BSA and Azide free	HA750330	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,ChIP,IP
	STEAP1 Rabbit pAb	ER1917-24	Human	WB,IHC-P,IF-Tissue
	Transferrin Mouse mAb [D7-D11]	M0809-5	Human	WB
	Transferrin Rabbit pAb	ER2001-53	Human	WB,IHC-P,FC
	Transferrin Rabbit pAb	R1212-1	Human	WB,IF-Cell
	Transferrin Mouse mAb [E4-E6]	M1510-17	Human	IF-Cell,IHC-P,FC,WB
KD	Transferrin Receptor (CD71) Recombinant Rabbit mAb [JF0956]	ET1702-06	Human, Mouse, Rat	WB,IP,FC
Relative expression	Transferrin Receptor (CD71) Recombinant Rabbit mAb [PSH13-56]	HA723545	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP,IF-Tissue
	Transferrin Receptor (CD71) Rabbit pAb	ER1706-35	Human	IHC-P
KD	Transferrin Receptor (CD71) Recombinant Rabbit mAb [JF0956] - BSA and Azide free	HA750333	Human, Mouse, Rat	WB,IP,FC
	Vitamin D3 Receptor Rabbit pAb	ER1803-18	Human, Rat	WB,IF-Cell,FC
KO	Vitamin D Receptor Recombinant Rabbit mAb [JA11-16]	ET1704-09	Human, Mouse, Rat	WB,IHC-P,FC
	Vitamin D Receptor Rabbit pAb	ER1802-40	Human	WB,FC
KO	Vitamin D Receptor Recombinant Rabbit mAb [JA11-16] - BSA and Azide free	HA750402	Human, Mouse, Rat	WB,IHC-P,FC
	xCT / SLC7A11 Mouse mAb [A7C5]	HA600097	Human, Mouse	WB,IHC-P,FC
Relative expression	xCT / SLC7A11 Recombinant Rabbit mAb [PSH19-94]	HA724114	Human, Mouse, Rat	WB,IHC-P
	xCT / SLC7A11 Recombinant Mouse mAb [A7C6-R] - BSA and Azide free	HA610011	Human, Mouse	WB,IHC-P,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	xCT / SLC7A11 Mouse mAb [A7C6]	HA600098	Human, Mouse	WB,IF-Cell,IHC-P,FC
	xCT / SLC7A11 Recombinant Rabbit mAb [JE31-93]	HA721868	Human, Mouse, Rat	WB,IP
	xCT / SLC7A11 Recombinant Mouse mAb [A7C6-R]	HA601071	Human, Mouse	WB,IHC-P,IF-Cell,FC
	ZNT3 Rabbit pAb	ER61422	Human, Rat, Mouse	WB