

重组蛋白-产品清单

| 产品名 | 货号 | Purity | Source |
|---|----------|--|--------------------|
| beta Amyloid 1-40 Protein | HA210498 | 87.46% | Synthetic peptide. |
| beta Amyloid 1-42 Protein | HA210499 | 85.18% | Synthetic peptide. |
| Biotinylated Mouse PD-L1/B7-H1 Protein (Primary Amine Labeling) | HA21401 | > 95% as determined by Tris-Bis PAGE | HEK293 |
| Biotinylated SARS-COV-2 Spike RBD Protein | HA210049 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Cynomolgus ACE2/ACEH Protein | HA210048 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Cynomolgus CD3E&CD3D/CD3 epsilon&CD3 delta Protein | HA210154 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Cynomolgus CD3E/CD3 epsilon Protein | HA210059 | > 95% as determined by Tris-Bis PAGE | HEK293 |
| Cynomolgus CD3E/CD3 epsilon Protein | HA210060 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Cynomolgus CD3E&CD3G/CD3 epsilon&CD3 gamma Protein | HA210153 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Cynomolgus FGL1 Protein | HA210194 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Cynomolgus IL-17A/CTLA-8 Protein | HA210122 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Cynomolgus IL-2 R beta/CD122 Protein | HA210089 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Cynomolgus IL-2 R gamma/CD132 Protein | HA210088 | > 95% as determined by Tris-Bis PAGE | HEK293 |
| Cynomolgus Nectin-4 Protein | HA210009 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Cynomolgus PD-1/PDCD1 Protein | HA210227 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Cynomolgus/Rhesus macaque BCMA/TNFRSF17 Protein | HA210198 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|---|----------|--|--------|
| Cynomolgus/Rhesus macaque PD-L1/B7-H1 Protein | HA210159 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Cynomolgus Siglec-15/CD33L3 Protein | HA210188 | > 95% as determined by Tris-Bis PAGE> 95% as determined by HPLC | HEK293 |
| Cynomolgus TNFR2/CD120b/TNFRSF1B Protein | HA210077 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Cynomolgus TROP-2/TACSTD2 Protein | HA210034 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Escherichia phage P1 Recombinase cre, C-His Tag Protein | HA210699 | >95% as determined by SDS-PAGE. | E.coli |
| Glu98-Leu231(NKG2C) & Ser34-Ile179(CD94) | HA210039 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| HIV-1 Gag-p24, C-His Tag Protein | HA210738 | >95% as determined by SDS-PAGE. | E.coli |
| Human 14-3-3 zeta, Tag Free Protein | HA211154 | >95% as determined by SDS-PAGE. | HEK293 |
| Human 2B4/CD244/SLAMF4 Protein | HA210218 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human 4-1BB Ligand/TNFSF9 Trimer Protein | HA210273 | > 95% as determined by Tris-Bis PAGE | HEK293 |
| Human 4-1BB Ligand/TNFSF9 Trimer Protein | HA210151 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human 4-1BB/TNFRSF9, C-His Tag Protein | HA210839 | >95% as determined by SDS-PAGE. | HEK293 |
| Human 5'-nucleotidase/CD73/NT5E, C-His Tag Protein | HA210915 | >90% as determined by SDS-PAGE. | HEK293 |
| Human A2M/alpha-2-macroglobulin, C-His Tag Protein | HA210693 | >95% as determined by SDS-PAGE. | HEK293 |
| Human ABL1/JTK7/p150, C-His, Flag Tag Protein | HA210682 | >95% as determined by SDS-PAGE. | HEK293 |
| Human ACE2, C-His Tag (ECD) Protein | HA21008 | >95% as determined by SDS-PAGE. | HEK293 |
| Human ACE, Tag Free Protein | HA211209 | >95% as determined by SDS-PAGE. | HEK293 |
| Human ACHE/Acetylcholinesterase, C-His Tag Protein | HA211178 | >90% as determined by SDS-PAGE. | HEK293 |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|---|----------|--------------------------------------|--------|
| Human ACP1, C-His Tag Protein | HA210510 | >95% as determined by SDS-PAGE. | E.coli |
| Human ACP2, C-His Tag Protein | HA210851 | >95% as determined by SDS-PAGE. | HEK293 |
| Human ACP3/ACPP, C His, Flag Tag Protein | HA210858 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Adiponectin/Acrp30, Tag Free Protein | HA210935 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Agrin, C-His Tag Protein | HA211173 | >90% as determined by SDS-PAGE. | HEK293 |
| Human AGRP, Tag Free Protein | HA211302 | >90% as determined by SDS-PAGE. | HEK293 |
| Human AHNAK, N-His Tag, C-His Tag Protein | HA210552 | >90% as determined by SDS-PAGE. | E.coli |
| Human AIDA, N-His, C-His Tag Protein | HA210586 | >95% as determined by SDS-PAGE. | E.coli |
| Human Alkaline phosphatase/APLI, C-His Tag Protein | HA210960 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Alpha-fetoprotein, C-His Tag | HA210967 | >95% as determined by SDS-PAGE. | HEK293 |
| Human alpha-Synuclein, C-His Tag Protein | HA211150 | >95% as determined by SDS-PAGE. | E.coli |
| Human Angiogenin, Tag Free Protein | HA210986 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Angiopoietin-1/ANGPT1, Tag Free Protein | HA210884 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Angiotensinogen, C-His Tag Protein | HA211165 | >95% as determined by SDS-PAGE. | HEK293 |
| Human ANGPT2/Angiopoietin-2 Protein | HA210249 | > 95% as determined by Tris-Bis PAGE | HEK293 |
| Human ANGPTL3/Angiopoietin-like 3, Tag Free Protein | HA211029 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Annexin A6/ANXA6, N-Twin Strep, C-His, Flag Tag Protein | HA210671 | >95% as determined by SDS-PAGE. | E.coli |
| Human ANPEP/CD13, Tag Free (ECD) Protein | HA210821 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Antithrombin III, C-His Tag Protein | HA210985 | >95% as determined by SDS-PAGE. | HEK293 |
| Human ANXA6, C-His Tag Protein | HA210830 | >95% as determined by SDS-PAGE. | E.coli |
| Human APOC3, Tag Free Protein | HA211274 | >95% as determined by SDS-PAGE. | HEK293 |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|---|----------|--|--------|
| Human Apolipoprotein A-I/APOA1, Tag Free | HA211090 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Apolipoprotein E/APOE4 (C130R), Tag Free | HA210965 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Apolipoprotein E, Tag Free | HA211078 | >95% as determined by SDS-PAGE. | HEK293 |
| Human APP, C-His Tag (ECD) Protein | HA211180 | >95% as determined by SDS-PAGE. | HEK293 |
| Human APRIL/TNFSF13 Trimer Protein | HA210271 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human Arginase 1/ARG1, C-His Tag Protein | HA210914 | >95% as determined by SDS-PAGE. | HEK293 |
| Human ASAH1, C-His, Flag Tag Protein | HA210848 | >95% as determined by SDS-PAGE. | HEK293 |
| Human ASGR1, Tag Free (ECD) Protein | HA210807 | >95% as determined by SDS-PAGE. | HEK293 |
| Human ASGR2, Tag Free (ECD) Protein | HA210842 | >95% as determined by SDS-PAGE. | HEK293 |
| Human ATF6, Tag Free Protein | HA210806 | >95% as determined by SDS-PAGE. | HEK293 |
| Human ATL3, C-His Tag Protein | HA210690 | >95% as determined by SDS-PAGE. | E.coli |
| Human Azurocidin/CAP37/AZU1/HBP Protein | HA210050 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human B7-1/CD80 Protein | HA210152 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human B7-1/CD80 Protein | HA210269 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human B7-1/CD80, Tag Free Protein | HA210994 | >95% as determined by SDS-PAGE. | HEK293 |
| Human B7-2/CD86 Protein | HA210270 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human B7-2/CD86 Protein | HA210008 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human B7-DC/PD-L2/CD273, C-Flag Tag (ECD) Protein | HA210599 | >95% as determined by SDS-PAGE. | HEK293 |
| Human B7-H3/CD276 Protein | HA210268 | > 95% as determined by Tris-Bis PAGE | HEK293 |
| Human B7-H3, C-His Tag (ECD) Protein | HA210964 | >95% as determined by SDS-PAGE. | HEK293 |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|---|----------|--|--------|
| Human B7-H4 Protein | HA210267 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human B7-H6/NCR3LG1, C-His Tag Protein | HA210889 | >95% as determined by SDS-PAGE. | HEK293 |
| Human B7-H6/NCR3LG1 Protein | HA210265 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human B7-H7/HLA2 Protein | HA210199 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human BACE1, C-His Tag Protein | HA211177 | >95% as determined by SDS-PAGE. | HEK293 |
| Human BAFF/BLyS/TNFSF13B, Tag Free (ECD) Protein | HA210874 | >95% as determined by SDS-PAGE. | HEK293 |
| Human BAFFR/TNFRSF13C Protein | HA210160 | > 95% as determined by Tris-Bis PAGE | HEK293 |
| Human BAFFR/TNFRSF13C Protein | HA210167 | > 95% as determined by Tris-Bis PAGE | HEK293 |
| Human BAFF/TNFSF13B/CD257 Protein | HA210264 | > 90% as determined by Tris-Bis PAGE > 90% as determined by HPLC | HEK293 |
| Human BAFF/TNFSF13B/CD257 Trimer Protein | HA210272 | > 95% as determined by Tris-Bis PAGE > 94% as determined by HPLC | HEK293 |
| Human BCMA/TNFRSF17 (1-54), Tag Free Protein | HA210968 | >95% as determined by SDS-PAGE. | HEK293 |
| Human BCMA/TNFRSF17, C-His, Flag Tag Protein | HA210776 | >95% as determined by SDS-PAGE. | E.coli |
| Human Beta-2-microglobulin/B2M, C-His Tag Protein | HA210978 | >95% as determined by SDS-PAGE. | HEK293 |
| Human beta-NGF, Tag Free Protein | HA211073 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Beta-synuclein/SNCB, Tag Free Protein | HA211157 | >95% as determined by SDS-PAGE. | HEK293 |
| Human BHLHE40, N-Strep II, C-His Tag Protein | HA210574 | >95% as determined by SDS-PAGE. | E.coli |
| Human BRCA2, C-Flag Protein | HA210793 | >95% as determined by SDS-PAGE. | E.coli |
| Human BST2/Tetherin, Tag Free Protein | HA211225 | >95% as determined by SDS-PAGE. | HEK293 |
| Human BTLA Protein | HA210195 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|---|----------|--------------------------------------|--------|
| Human C1q R1 / CD93, Tag Free Protein | HA211294 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CA125/MUC16, C-His Tag Protein | HA211195 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CA125/MUC16 Protein | HA210124 | > 95% as determined by Tris-Bis PAGE | HEK293 |
| Human Cadherin-3, C-His Tag Protein | HA211162 | >90% as determined by SDS-PAGE. | HEK293 |
| Human CALCA, C-His Tag Protein | HA210832 | >95% as determined by SDS-PAGE. | E.coli |
| Human CALCOCO2/NDP52, Tag Free Protein | HA210782 | >95% as determined by SDS-PAGE. | E.coli |
| Human Calprotectin (S100A8&S100A9) Protein | HA210041 | > 95% as determined by Tris-Bis PAGE | E.coli |
| Human Carbonic Anhydrase XII/CA12, Tag Free (ECD) Protein | HA210702 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Cardiotrophin-1/CT-1, Tag Free Protein | HA211250 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Caspase-9/CASP9, Tag Free Protein | HA211015 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CASQ1, Tag Free Protein | HA210799 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CASQ2, Tag Free Protein | HA210559 | >90% as determined by SDS-PAGE. | HEK293 |
| Human Cathepsin D, C-His Tag Protein | HA211077 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CBX8, C-His Tag Protein | HA210826 | >95% as determined by SDS-PAGE. | E.coli |
| Human CCL11/Eotaxin, Tag Free Protein | HA210865 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CCL16/LEC, Tag Free Protein | HA211292 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CCL17/TARC, Tag Free Protein | HA210895 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CCL18/PARC/MIP-4 Protein | HA211080 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CCL19 / MIP-3 beta, Tag Free Protein | HA211211 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CCL20, Tag Free Protein | HA211012 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CCL21/6CKine, Tag Free Protein | HA211144 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CCL22 / MDC, C-hFc Tag Protein | HA211312 | >95% as determined by SDS-PAGE. | HEK293 |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|---|----------|--|--------|
| Human CCL23/MPIF-1, Tag Free Protein | HA210894 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CCL25/TECK, Tag Free Protein | HA211000 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CCL26/Eotaxin-3, Tag Free Protein | HA211152 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CCL27, Tag Free Protein | HA211133 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CCL2/MCP1, Tag Free Protein | HA210787 | >95% as determined by SDS-PAGE. | E.coli |
| Human CCL3L1, N-Twin Strep, C-His, Flag Tag Protein | HA210749 | >95% as determined by SDS-PAGE. | E.coli |
| Human CCL3/MIP-1, Tag Free Protein | HA210901 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CCL4/MIP-1, Tag Free Protein | HA211213 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CCL5/RANTES, Tag Free Protein | HA210998 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CD147/Basigin/EMMPRIN, C-His Tag Protein | HA210959 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CD14, C-His Tag Protein | HA211004 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CD155/PVR, C-His Tag (ECD) Protein | HA210957 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CD155/PVR Protein | HA210149 | > 95% as determined by Tris-Bis PAGE > 95% as determined by SEC-HPLC | HEK293 |
| Human CD155/PVR Protein | HA210243 | > 95% as determined by Tris-Bis PAGE | HEK293 |
| Human CD155/PVR, Tag Free (ECD) Protein | HA211135 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CD160 Protein | HA210117 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human CD163, Tag Free (ECD) Protein | HA211143 | >90% as determined by SDS-PAGE. | HEK293 |
| Human CD200, C-His Tag (ECD) Protein | HA210665 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CD200R1, C-His, Flag (ECD) Protein | HA210572 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CD207, C-Flag Tag Protein | HA210802 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CD24 Protein | HA210170 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|---|----------|---|--------|
| Human CD24 Protein | HA210171 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human CD24 Protein | HA210123 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human CD272/BTLA, C-His, Flag Tag (ECD) Protein | HA210546 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CD27 Ligand/CD70 Protein | HA210201 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human CD27/TNFRSF7 Protein | HA210247 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human CD27/TNFRSF7, Tag Free Protein | HA210880 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CD28, C-His Tag Protein | HA210920 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CD2, C-His, Flag Tag (ECD) Protein | HA210554 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CD2, N-Twin Strep, C-His, Flag Tag Protein | HA210512 | >95% as determined by SDS-PAGE. | E.coli |
| Human CD30 Ligand/TNFSF8 Protein | HA210191 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human CD30 Ligand / TNFSF8, Tag Free Protein | HA211293 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CD30/TNFRSF8, C-His, Flag Tag (ECD) Protein | HA210649 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CD31/PECAM-1, C-His Tag Protein | HA210891 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CD38 Protein | HA210246 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human CD3E&CD3D/CD3 epsilon&CD3 delta Protein | HA210260 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human CD3E/CD3 epsilon Protein | HA210146 | > 95% as determined by Tris-Bis PAGE | HEK293 |
| Human CD3E/CD3 epsilon Protein | HA210065 | > 95% as determined by Tris-Bis PAGE | HEK293 |
| Human CD3E&CD3G/CD3 epsilon&CD3 gamma Protein | HA210259 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human CD40 Ligand/TNFSF5, Tag Free (ECD) Protein | HA210939 | >90% as determined by SDS-PAGE. | HEK293 |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|---|----------|--|--------|
| Human CD40 Ligand/TNFSF5 Trimer Protein | HA210233 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human CD40 Ligand/TNFSF5 Trimer Protein | HA210234 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human CD40/TNFRSF5 Protein | HA210204 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human CD46, Tag Free Protein | HA211227 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CD47, C-His Tag (ECD) Protein | HA211050 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CD47 Protein | HA210208 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human CD47 Protein | HA210207 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human CD48, C-His Tag Protein | HA210906 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CD48/SLAMF2 Protein | HA210245 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human CD4/LEU3 Protein | HA210102 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human CD52 Protein | HA210114 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human CD58/LFA-3, C-His Tag (ECD) Protein | HA210605 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CD5, Tag Free (ECD) Protein | HA210822 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CD63, Tag Free (ECD) Protein | HA210764 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CD73/NT5E Protein | HA210244 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human CD79B Protein | HA210169 | > 95% as determined by Tris-Bis PAGE | HEK293 |
| Human CD83, C-His Tag Protein | HA211203 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CD94 Protein | HA210162 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human CDC42, Tag Free Protein | HA210767 | >95% as determined by SDS-PAGE. | E.coli |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|---|----------|--|--------|
| Human CDCP1/CD318, N-Twin Strep, C-His, Flag Tag Protein | HA210664 | >95% as determined by SDS-PAGE. | E.coli |
| Human CDNF, C-His Tag Protein | HA210616 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CEACAM-1/CD66a, C-His Tag (ECD) Protein | HA210526 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CEACAM-5/CD66e, Tag Free Protein | HA211260 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CEACAM-6/CD66c Protein | HA210013 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human CENPF, C-His Tag Protein | HA210523 | | E.coli |
| Human CFH, Tag Free Protein | HA211200 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CGAS, C-His Tag Protein | HA210827 | >95% as determined by SDS-PAGE. | E.coli |
| Human Chitinase 3-like 1/CHI3L1, C-His Tag Protein | HA211030 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Chitotriosidase/CHIT1, C-His Tag Protein | HA211215 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CHT, C-His Tag Protein | HA210645 | >95% as determined by SDS-PAGE. | E.coli |
| Human cIAP1/BIRC2, Tag Free Protein | HA211301 | >90% as determined by SDS-PAGE. | HEK293 |
| Human cIAP-2/BIRC3, N-Twin Strep, C-His, Flag Tag Protein | HA210619 | >90% as determined by SDS-PAGE. | E.coli |
| Human Cingulin/CGN, N-Twin Strep, C-His, Flag Tag Protein | HA210584 | >95% as determined by SDS-PAGE. | E.coli |
| Human CKB, C-His, Flag Tag Protein | HA210777 | >90% as determined by SDS-PAGE. | E.coli |
| Human CKMB heterodimer, C-Flag, StrepII Tag&C-His Tag Protein | HA210947 | >95% as determined by SDS-PAGE. | E.coli |
| Human CKM, C-His, Flag Tag Protein | HA210778 | >95% as determined by SDS-PAGE. | E.coli |
| Human Clusterin, C-His, Flag Tag Protein | HA210955 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CNOT7, C-His Tag Protein | HA210537 | >95% as determined by SDS-PAGE. | E.coli |
| Human CNTF, N-His Tag, C-His Tag Protein | HA210597 | >95% as determined by SDS-PAGE. | E.coli |
| Human CNTF, Tag Free Protein | HA210944 | >95% as determined by SDS-PAGE. | HEK293 |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|--|----------|--|--------|
| Human Coagulation factor VIII, Tag Free Protein | HA211190 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Complement C3, C-His Tag Protein | HA211092 | >90% as determined by SDS-PAGE. | HEK293 |
| Human Complement Component C5a, C-His Tag Protein | HA210636 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Complement Factor D/CFD, C-His Tag Protein | HA211128 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Corticoliberin/CRH, C-hFc Tag Protein | HA211316 | >95% as determined by SDS-PAGE. | HEK293 |
| Human COX5A, C-His Tag Protein | HA210534 | >90% as determined by SDS-PAGE. | E.coli |
| Human C-Reactive Protein/CRP, C-His, Flag Tag | HA210648 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Cripto-1/TDGF1, C-His, Flag Tag Protein | HA210543 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CRKL, C-His Tag Protein | HA210551 | >95% as determined by SDS-PAGE. | E.coli |
| Human CRLF3, N-Twin Strep, C-His, Flag Tag Protein | HA210724 | >95% as determined by SDS-PAGE. | E.coli |
| Human CSF3R/G-CSFR/CD114, C-His Tag Protein | HA211262 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CSK, C-His Tag Protein | HA210519 | >95% as determined by SDS-PAGE. | E.coli |
| Human CSPP1, C-His, Flag Tag Protein | HA210757 | >95% as determined by SDS-PAGE. | E.coli |
| Human CTGF/CCN2, C-Flag Tag Protein | HA211051 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CTLA-4/CD152 Protein | HA210241 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human CTLA-4, C-His Tag Protein | HA210885 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CXADR/CAR, C-His Tag Protein | HA211206 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CXCL10, C-Flag Tag Protein | HA210740 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CXCL11/I-TAC, Tag Free Protein | HA211110 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CXCL12/SDF-1, Tag Free Protein | HA210900 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CXCL13 / BCA-1, Tag Free Protein | HA211065 | >95% as determined by SDS-PAGE. | HEK293 |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|---|----------|--|--------|
| Human CXCL16, Tag Free Protein | HA211138 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CXCL1, Tag Free Protein | HA210838 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CXCL2/GRO beta, Tag Free Protein | HA210904 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CXCL5, Tag Free Protein | HA210887 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CXCL6/GCP-2, Tag Free Protein | HA211111 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CXCL9, Tag Free Protein | HA210774 | >95% as determined by SDS-PAGE. | HEK293 |
| Human CYLD, C-His Tag Protein | HA210704 | >95% as determined by SDS-PAGE. | E.coli |
| Human/Cynomolgus/Rhesus macaque CD28 Protein | HA210209 | > 95% as determined by Tris-Bis PAGE | HEK293 |
| Human/Cynomolgus/Rhesus macaque ROR1 Protein | HA210115 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human/Cynomolgus/Rhesus macaque ROR1 Protein | HA210116 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human Cystatin C, C-His Tag Protein | HA210995 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Cytokeratin 8, Tag Free Protein | HA211101 | >95% as determined by SDS-PAGE. | HEK293 |
| Human DDX1 (436-740), N-Twin Strep, C-His, Flag Tag Protein | HA210753 | >95% as determined by SDS-PAGE. | E.coli |
| Human DDX1, N-Twin Strep, C-His, Flag Tag Protein | HA210517 | >95% as determined by SDS-PAGE. | E.coli |
| Human DICER1, C-His Tag Protein | HA210673 | >95% as determined by SDS-PAGE. | E.coli |
| Human DICR/CLEC4A, Tag Free Protein | HA211226 | >95% as determined by SDS-PAGE. | HEK293 |
| Human DKK1, C-His Tag Protein | HA211040 | >95% as determined by SDS-PAGE. | HEK293 |
| Human DLD, N-Twin Strep, C-His Tag Protein | HA210504 | >95% as determined by SDS-PAGE. | E.coli |
| Human DLK1, C-His Tag (ECD) Protein | HA210634 | >95% as determined by SDS-PAGE. | HEK293 |
| Human DLL3 Protein | HA210242 | > 95% as determined by Tris-Bis PAGE | HEK293 |
| Human DNAM-1/CD226 Protein | HA210240 | > 95% as determined by Tris-Bis PAGE > 95% as determined by SEC-HPLC | HEK293 |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|---|----------|--|--------|
| Human DOT1L, N-Strep II, C-His Tag Protein | HA210700 | >95% as determined by SDS-PAGE. | E.coli |
| Human DPP-4/CD26/DPPIV, C-His Tag (ECD) Protein | HA211034 | >95% as determined by SDS-PAGE. | HEK293 |
| Human DR3/TNFRSF25, Tag Free Protein | HA211322 | >95% as determined by SDS-PAGE. | HEK293 |
| Human DVL2, C Twin Strep, Flag Protein | HA210781 | >95% as determined by SDS-PAGE. | E.coli |
| Human E-Cadherin, C-His Tag (ECD) Protein | HA210886 | >95% as determined by SDS-PAGE. | HEK293 |
| Human ECM1, C-His Tag Protein | HA210933 | >95% as determined by SDS-PAGE. | HEK293 |
| Human EGFR, C-His Tag Protein | HA211104 | >95% as determined by SDS-PAGE. | HEK293 |
| Human EGFR/HER1 Protein | HA210147 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human EGF, Tag Free Protein | HA210717 | >95% as determined by SDS-PAGE. | HEK293 |
| Human EIF2AK3/PERK, C-His Tag Protein | HA210525 | >95% as determined by SDS-PAGE. | E.coli |
| Human EIF3K, N-Twin Strep, C-His, Flag Tag Protein | HA210535 | >95% as determined by SDS-PAGE. | E.coli |
| Human EIF4B, C-His Tag Protein | HA210835 | >95% as determined by SDS-PAGE. | E.coli |
| Human EIF4G1, N-Twin Strep, C-His, Flag Tag Protein | HA210818 | >95% as determined by SDS-PAGE. | E.coli |
| Human EN2, Tag Free Protein | HA210803 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Endoglin/CD105, C-His Tag (ECD) Protein | HA211006 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Endoglin/CD105 Protein | HA210069 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human Endorepellin/Perlecan, C-His Tag Protein | HA211127 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Endostatin/COL18A1, Tag Free Protein | HA211001 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Enolase 2/ENO2, C-His Tag Protein | HA211021 | >95% as determined by SDS-PAGE. | E.coli |
| Human EpCAM/TROP1 Protein | HA210205 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human EPCR/PROCR/CD201, C-His, Flag Tag (ECD) Protein | HA210662 | >95% as determined by SDS-PAGE. | HEK293 |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|---|----------|--|--------|
| Human EPHA2, C-His Tag (ECD) Protein | HA211075 | >95% as determined by SDS-PAGE. | HEK293 |
| Human ERO1A, C-His, Flag Tag Protein | HA210868 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Erythropoietin/EPO, C-His, Flag Tag Protein | HA210852 | >95% as determined by SDS-PAGE. | HEK293 |
| Human E-Selectin/CD62E, Tag Free (ECD) Protein | HA210996 | >95% as determined by SDS-PAGE. | HEK293 |
| Human FABP3 / H-FABP, Tag Free Protein | HA211159 | >95% as determined by SDS-PAGE. | HEK293 |
| Human FABP7/B-FABP, C-His Tag Protein | HA210677 | >95% as determined by SDS-PAGE. | E.coli |
| Human Factor B, C-His Tag Protein | HA211099 | >95% as determined by SDS-PAGE. | HEK293 |
| Human FAM171A2, Tag Free Protein | HA211284 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Fas/TNFRSF6/CD95, C-His Tag (ECD) Protein | HA211064 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Fc gamma RI/CD64 Protein | HA210142 | > 95% as determined by Tris-Bis PAGE | HEK293 |
| Human Fc gamma RIIA/CD32a (H167) Protein | HA210144 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human Fc gamma RIIA/CD32a (R167) Protein | HA210180 | > 95% as determined by Tris-Bis PAGE> 95% as determined by HPLC | HEK293 |
| Human Fc gamma RIIB/CD32b Protein | HA210192 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human Fc gamma RIIB/CD16a (F176) Protein | HA210143 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human Fc gamma RIIB/CD16a (V176) Protein | HA210178 | > 95% as determined by Tris-Bis PAGE> 95% as determined by HPLC | HEK293 |
| Human Fc gamma RIIB/CD16a (V176) Protein | HA210281 | > 95% as determined by Bis-Tris PAGE > 95% as determined by HPLC | HEK293 |
| Human Fc gamma RIIB/CD16b (NA2) Protein | HA210179 | > 95% as determined by Tris-Bis PAGE> 95% as determined by HPLC | HEK293 |
| Human FcRn Protein | HA210239 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human Fetuin A/AHSG, Tag Free Protein | HA211146 | >95% as determined by SDS-PAGE. | HEK293 |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|---|----------|--|--------|
| Human FGF10, C-Flag Protein | HA210794 | >95% as determined by SDS-PAGE. | E.coli |
| Human FGF19, Tag Free Protein | HA210981 | >95% as determined by SDS-PAGE. | HEK293 |
| Human FGF1, C-His Tag Protein | HA211079 | >95% as determined by SDS-PAGE. | E.coli |
| Human FGF-21, C-His Tag Protein | HA210921 | >95% as determined by SDS-PAGE. | HEK293 |
| Human FGF21 Protein | HA210001 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human FGF2/b-FGF, C-His, Flag Tag Protein | HA210582 | >90% as determined by SDS-PAGE. | E.coli |
| Human FGF-4, C-His Tag Protein | HA210946 | >95% as determined by SDS-PAGE. | E.coli |
| Human FGF7, C-Flag Tag Protein | HA210769 | >90% as determined by SDS-PAGE. | E.coli |
| Human FGF7, Tag Free Protein | HA210829 | >95% as determined by SDS-PAGE. | E.coli |
| Human FGF Basic/FGF2, Tag Free Protein | HA210834 | >95% as determined by SDS-PAGE. | E.coli |
| Human FGF basic Protein | HA210238 | > 95% as determined by Tris-Bis PAGE | E.coli |
| Human FGFR3 alpha (IIIb) Protein | HA210086 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human FGFR3 alpha (IIIc) Protein | HA210085 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human FGFR4 Protein | HA210140 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human FGFR4 Protein | HA210139 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human FGL1 Protein | HA210275 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human FGL1 Protein | HA210255 | > 95% as determined by Tris-Bis PAGE > 90% as determined by HPLC | HEK293 |
| Human Fibrillin-1/FBN1, Tag Free Protein | HA211184 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Fibronectin/FN1, Tag Free Protein | HA211024 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Ficolin-2, C-His Tag Protein | HA211174 | >90% as determined by SDS-PAGE. | HEK293 |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|--|----------|---|--------|
| Human FLI1, Tag Free Protein | HA210798 | >95% as determined by SDS-PAGE. | E.coli |
| Human Flt3 Ligand/FLT3LG, Tag Free Protein | HA211222 | >95% as determined by SDS-PAGE. | HEK293 |
| Human FOLR1 Protein | HA210193 | > 95% as determined by Tris-Bis PAGE; > 95% as determined by HPLC | HEK293 |
| Human FOLR1, Tag Free Protein | HA211055 | >95% as determined by SDS-PAGE. | HEK293 |
| Human FOXP2. C-His Tag Protein | HA210654 | >95% as determined by SDS-PAGE. | E.coli |
| Human Fractalkine/CX3CL1, C-His Tag (ECD) Protein | HA211057 | >95% as determined by SDS-PAGE. | HEK293 |
| Human FSH Beta, C-His Tag Protein | HA210932 | >95% as determined by SDS-PAGE. | HEK293 |
| Human GADD45A, N-Twin Strep, C-His, Flag Tag Protein | HA210743 | >95% as determined by SDS-PAGE. | E.coli |
| Human GADD45B, N-Twin Strep, C-His, Flag Tag Protein | HA210533 | >90% as determined by SDS-PAGE. | E.coli |
| Human Galectin 9/LGALS9 Protein | HA210237 | > 95% as determined by Tris-Bis PAGE | HEK293 |
| Human Galectin-9 / LGALS9, Tag Free Protein | HA211310 | >95% as determined by SDS-PAGE. | HEK293 |
| Human GANAB, C-His Tag Protein | HA210610 | >95% as determined by SDS-PAGE. | HEK293 |
| Human G-CSF, C-His Tag Protein | HA210840 | >95% as determined by SDS-PAGE. | HEK293 |
| Human GDNF, Tag Free Protein | HA211323 | >90% as determined by SDS-PAGE. | HEK293 |
| Human GFAP, Tag Free Protein | HA210990 | >95% as determined by SDS-PAGE. | E.coli |
| Human GFRAL/GFR alpha-like Protein | HA210111 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human GGA1, C-His Tag Protein | HA210748 | >95% as determined by SDS-PAGE. | E.coli |
| Human GH1, C-His Tag Protein | HA211027 | >95% as determined by SDS-PAGE. | HEK293 |
| Human GHRL, C-His, Flag Tag Protein | HA210855 | >95% as determined by SDS-PAGE. | HEK293 |
| Human GITR Ligand/TNFSF18, Tag Free Protein | HA211265 | >95% as determined by SDS-PAGE. | HEK293 |
| Human GITR Ligand/TNFSF18 Trimer Protein | HA210032 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|---|----------|--|--------|
| Human GITR/TNFRSF18/AITR, C-Flag Tag (ECD) Protein | HA210945 | >95% as determined by SDS-PAGE. | HEK293 |
| Human GM-CSF Protein | HA210236 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | E.coli |
| Human GM-CSF, Tag Free Protein | HA210712 | >95% as determined by SDS-PAGE. | HEK293 |
| Human GPA33/A33 Protein | HA210174 | > 95% as determined by Tris-Bis PAGE> 95% as determined by HPLC | HEK293 |
| Human GPC3/Glypican 3 Protein | HA210235 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human Grancalcin/GCA, N-Twin Strep, C-His, Flag Tag Protein | HA210563 | >95% as determined by SDS-PAGE. | E.coli |
| Human GranzymeA/GZMA, Tag Free Protein | HA210916 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Granzyme B, Tag Free Protein | HA211076 | >95% as determined by SDS-PAGE. | HEK293 |
| Human GRO gamma/CXCL3, Tag Free Protein | HA211249 | >95% as determined by SDS-PAGE. | HEK293 |
| Human GRP56, C-His, Flag Tag (ECD) Protein | HA210812 | >95% as determined by SDS-PAGE. | HEK293 |
| Human GSTO1, C-His Tag Protein | HA210509 | >95% as determined by SDS-PAGE. | E.coli |
| Human GTPase HRAS, C-His Tag Protein | HA211122 | >95% as determined by SDS-PAGE. | E.coli |
| Human GUCY2C/Guanylyl cyclase C Protein | HA210037 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human Haptoglobin/HP, C-His Tag Protein | HA211100 | >95% as determined by SDS-PAGE. | HEK293 |
| Human HBEGF, Tag Free Protein | HA211320 | >95% as determined by SDS-PAGE. | HEK293 |
| Human HCLS1, N-Strep II, C-His Tag Protein | HA210624 | >90% as determined by SDS-PAGE. | E.coli |
| Human Her3/ErbB3 Protein | HA210105 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human Her3/ErbB3 Protein | HA210203 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human HEXB, C-His Tag Protein | HA210540 | >95% as determined by SDS-PAGE. | HEK293 |
| Human HGF, C-Flag Tag Protein | HA210739 | >90% as determined by SDS-PAGE. | HEK293 |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|---|----------|---|--------|
| Human HGF R/c-MET Protein | HA210187 | > 95% as determined by Tris-Bis PAGE; > 95% as determined by HPLC | HEK293 |
| Human HGF, Tag Free Protein | HA210913 | >95% as determined by SDS-PAGE. | HEK293 |
| Human HLA-A*02:01&B2M&GP100 (YLEPGPVTA) Monomer Protein | HA210173 | > 95% as determined by Tris-Bis PAGE> 95% as determined by HPLC | HEK293 |
| Human HLA-A*02:01&B2M&WT-1 (RMFPNAPYL) Monomer Protein | HA210202 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human HSPA1A, Tag Free Protein | HA210756 | >90% as determined by SDS-PAGE. | E.coli |
| Human Human LOX-1/OLR1, C-His Tag (ECD) Protein | HA210941 | >95% as determined by SDS-PAGE. | HEK293 |
| Human HVEM/TNFRSF14, C-Flag Tag Protein | HA210556 | >95% as determined by SDS-PAGE. | HEK293 |
| Human HVEM/TNFRSF14, Tag Free (ECD) Protein | HA210849 | >95% as determined by SDS-PAGE. | HEK293 |
| Human I-309/CCL1, Tag Free Protein | HA211081 | >90% as determined by SDS-PAGE. | HEK293 |
| Human ICAM-1/CD54, C-His Tag (ECD) Protein | HA210857 | >95% as determined by SDS-PAGE. | HEK293 |
| Human ICAM-2/CD102, C-His Tag Protein | HA211314 | >95% as determined by SDS-PAGE. | HEK293 |
| Human ICAM-3/CD50, C-His Tag (ECD) Protein | HA211145 | >95% as determined by SDS-PAGE. | HEK293 |
| Human ICOS Ligand/B7-H2/ICOSLG, C-His Tag (ECD) Protein | HA210593 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IDO, C-His Protein | HA210938 | >90% as determined by SDS-PAGE. | E.coli |
| Human IFNA1/IFNA13, C-Tag Protein | HA210732 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IFN-alpha/beta R1, C-His (ECD) Protein | HA210626 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IFN-alpha/beta R2, C-His Tag (ECD) Protein | HA210608 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IFN-gamma, C-His Tag Protein | HA210762 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IFN gamma/IFNG Protein | HA210005 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human IFN-gamma R1, C-His Tag (ECD) Protein | HA210514 | >95% as determined by SDS-PAGE. | HEK293 |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|---|----------|---|--------|
| Human IFN-gamma, Tag Free Protein | HA211129 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IGF1R/CD221 Protein | HA210184 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human IGF1, Tag Free Protein | HA210866 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IGFBP-1, Tag Free Protein | HA211115 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IGFBP2, C-His Tag Protein | HA211103 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IGFBP3/IBP3, C-His, Flag Tag Protein | HA210585 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IGFBP4, Tag Free Protein | HA211096 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IGFBP-6, Tag Free Protein | HA211086 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IGFBP7, C-His Tag Protein | HA211108 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IGF-II/IGF2, Tag Free Protein | HA211319 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL-10, C-His Tag Protein | HA210729 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL-10 Protein | HA210040 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human IL-10RB, C-His Tag Protein | HA211205 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL-10, Tag Free Protein | HA210988 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL12B/P40, C-His Tag Protein | HA210655 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL-12 (linked heterodimer), C-His Tag Protein | HA210910 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL-12 Protein | HA210019 | > 95% as determined by Tris-Bis PAGE > 95% as determined by SEC-HPLC | HEK293 |
| Human IL12RB1, C-His Tag (ECD) Protein | HA210728 | >90% as determined by SDS-PAGE. | HEK293 |
| Human IL-13Ra2 Protein | HA210134 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human IL-13 R alpha 2/IL13RA2, Tag Free (ECD) Protein | HA210969 | >95% as determined by SDS-PAGE. | HEK293 |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|---|----------|--|--------|
| Human IL-13, Tag Free Protein | HA210908 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL-15, C-Flag Tag Protein | HA210773 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL-15 Protein | HA210010 | > 95% as determined by Tris-Bis PAGE > 95% as determined by SEC-HPLC | E.coli |
| Human IL-16, C-His Tag Protein | HA210917 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL-17A, Tag Free Protein | HA210731 | >90% as determined by SDS-PAGE. | HEK293 |
| Human IL17B, Tag Free Protein | HA211270 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL-17C, C-hFc Tag Protein | HA211315 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL-17F, Tag Free Protein | HA210716 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL17RA/CD217, C-His Tag (ECD) Protein | HA210581 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL-18 BPα, C-His Tag Protein | HA210686 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL-18, Tag Free Protein | HA210779 | >95% as determined by SDS-PAGE. | E.coli |
| Human IL-18, Tag Free (HEK293) Protein | HA211285 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL-19, Tag Free Protein | HA211269 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL-1 α/IL1A, C-His, Flag Tag Protein | HA210602 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL-1 β/IL1B, C-His Tag Protein | HA210576 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL1R1, C-His Tag Protein | HA211261 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL-1R3/IL-1 RAcP Protein | HA210101 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human IL1RA, C-His Tag Protein | HA210837 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL-1 RII, C-His Tag (ECD) Protein | HA210684 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL1RL1/DER4/ST2, C-hFc, His (ECD) Protein | HA210589 | >90% as determined by SDS-PAGE. | HEK293 |
| Human IL1RL1/DER4/ST2, C-His (ECD) Protein | HA210713 | >95% as determined by SDS-PAGE. | HEK293 |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|---|----------|--|--------|
| Human IL-1 Rrp2/IL-1 R6 Protein | HA210100 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human IL-20, C-hFc Tag Protein | HA211268 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL-21, C-His, Flag Tag Protein | HA210578 | >90% as determined by SDS-PAGE. | HEK293 |
| Human IL-21, Tag Free Protein | HA211070 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL-22, C-His, Flag Tag Protein | HA210691 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL-23, C-His Tag Protein | HA210909 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL-23R Protein | HA210175 | > 95% as determined by Tris-Bis PAGE> 95% as determined by HPLC | HEK293 |
| Human IL-25/IL-17E, Tag Free Protein | HA211098 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL-27, C-Flag Tag Protein | HA210912 | >90% as determined by SDS-PAGE. | HEK293 |
| Human IL-28A/IFN-lambda 2, Tag Free Protein | HA210956 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL-29/IFN-lambda 1, Tag Free Protein | HA210970 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL-2, C-His Tag Protein | HA210670 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL-2 Protein | HA210230 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | E.coli |
| Human IL-2R alpha, Tag Free Protein | HA210771 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL-2RB/CD122, C-His Tag (ECD) Protein | HA210544 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL2RG/CD132, C-His Tag (ECD) Protein | HA210567 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL-2, Tag Free Protein | HA210795 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL-31, C-Flag Tag Protein | HA210923 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL-33, C-His Tag Protein | HA210539 | >95% as determined by SDS-PAGE. | E.coli |
| Human IL-34, Tag Free Protein | HA211267 | >90% as determined by SDS-PAGE. | HEK293 |
| Human IL-35 (linked heterodimer), C-hFc Tag Protein | HA211212 | >90% as determined by SDS-PAGE. | HEK293 |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|--|----------|--|--------|
| Human IL-3R alpha/CD123, C-His Tag Protein | HA211048 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL-3 R alpha/CD123 Protein | HA210224 | > 95% as determined by Tris-Bis PAGE | HEK293 |
| Human IL-3, Tag Free Protein | HA210928 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL-4 Protein | HA210002 | > 95% as determined by Tris-Bis PAGE | HEK293 |
| Human IL-4R, C-His, Tag (ECD) Protein | HA210720 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL-4, Tag Free Protein | HA210763 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL-5, C-His Tag Protein | HA210594 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL5RA, C-His Tag (ECD) Protein | HA210727 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL-6, C-His Tag Protein | HA210639 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL-6 Protein | HA210229 | > 95% as determined by Tris-Bis PAGE > 90% as determined by HPLC | E.coli |
| Human IL-6 R alpha/CD126 Protein | HA210164 | > 95% as determined by Tris-Bis PAGE> 95% as determined by HPLC | HEK293 |
| Human IL-6R, C-His (ECD) Protein | HA210640 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL-6, Tag Free Protein | HA210797 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL-7, C-His, Flag Tag Protein | HA210595 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL-7 R alpha / CD127, Tag Free Protein | HA211291 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL-7, Tag Free Protein | HA210815 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL-8/CXCL8 Tag Free Protein | HA210642 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IL-9, C-Flag Tag Protein | HA210911 | >95% as determined by SDS-PAGE. | HEK293 |
| Human ILKAP, N-His, C-His Tag Protein | HA210744 | >95% as determined by SDS-PAGE. | E.coli |
| Human IMMT/Mitofilin, C-His Tag Protein | HA210532 | >95% as determined by SDS-PAGE. | E.coli |
| Human Insulin, Tag Free Protein | HA211183 | >95% as determined by SDS-PAGE. | E.coli |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|---|----------|--|--------|
| Human Integrin alpha 5 beta 1 (ITGA5&ITGB1) Heterodimer Protein | HA210128 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human Integrin alpha V beta 3 (ITGAV&ITGB3) Heterodimer Protein | HA210126 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human Integrin alpha V beta 5 (ITGAV&ITGB5) Heterodimer Protein | HA210125 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human Interferon alpha-2/IFNA2, C-His Tag Protein | HA210644 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Interferon omega-1, Tag Free Protein | HA211271 | >95% as determined by SDS-PAGE. | HEK293 |
| Human IP-10/CXCL10, Tag Free Protein | HA210896 | >95% as determined by SDS-PAGE. | HEK293 |
| Human KAT2A/GCN5, N-His, C-His Tag Protein | HA210633 | >95% as determined by SDS-PAGE. | E.coli |
| Human KDM3A, C-His Tag Protein | HA210614 | >90% as determined by SDS-PAGE. | E.coli |
| Human KEAP1, C-His Tag Protein | HA210752 | >95% as determined by SDS-PAGE. | E.coli |
| Human KHK, C-His Tag Protein | HA210502 | >95% as determined by SDS-PAGE. | E.coli |
| Human KIF5B, N-Twin Strep, C-His Tag Protein | HA210680 | >95% as determined by SDS-PAGE. | E.coli |
| Human KIR2DL1 Protein | HA210070 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human KIR2DL2 Protein | HA210016 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human KIR2DL3 Protein | HA210072 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human KIR3DL3 Protein | HA210021 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human KIRREL1, C-His, Flag Tag (ECD) Protein | HA210661 | >95% as determined by SDS-PAGE. | HEK293 |
| Human KLK3, Tag Free Protein | HA211208 | >95% as determined by SDS-PAGE. | HEK293 |
| Human KSR1, C-His Tag Protein | HA210518 | >95% as determined by SDS-PAGE. | E.coli |
| Human LAG3/CD223 Protein | HA210225 | > 95% as determined by Tris-Bis PAGE | HEK293 |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|--|----------|--|--------|
| Human LAG3/CD223 Protein | HA210226 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human LAG-3, Tag Free Protein | HA210934 | >95% as determined by SDS-PAGE. | HEK293 |
| Human LAIR1, N-Twin Strep, C-His, Flag Tag Protein | HA210566 | >95% as determined by SDS-PAGE. | E.coli |
| Human LAMC2, Tag Free Protein | HA211219 | >95% as determined by SDS-PAGE. | HEK293 |
| Human LAMP3/CD208, Tag Free Protein | HA211259 | >95% as determined by SDS-PAGE. | HEK293 |
| Human LAP3, C-His Tag (Isoform 2) Protein | HA211317 | >95% as determined by SDS-PAGE. | E.coli |
| Human LDHB, C-His Tag Protein | HA211198 | >95% as determined by SDS-PAGE. | E.coli |
| Human LDLR/LDL Receptor, Tag Free Protein | HA211287 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Leptin Receptor, C-His Tag (ECD) Protein | HA211148 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Leptin, Tag Free Protein | HA210961 | >95% as determined by SDS-PAGE. | HEK293 |
| Human LIAS, N-Twin Strep, C-His, Flag Tag Protein | HA210500 | >95% as determined by SDS-PAGE. | E.coli |
| Human LIF, C-His Tag Protein | HA210736 | >95% as determined by SDS-PAGE. | HEK293 |
| Human LIGHT/TNFSF14 Protein | HA210217 | > 95% as determined by Tris-Bis PAGE | HEK293 |
| Human LIGHT/TNFSF14 Trimer Protein | HA210223 | > 90% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human LILRA1/CD85i/LIR-6 Protein | HA210109 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human LILRA3/CD85e Protein | HA210110 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human LILRA5/CD85f/ILT11 Protein | HA210112 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human LILRA6/CD85b/ILT8 Protein | HA210113 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human LILRB1/CD85j/ILT2, C-His, Flag Tag (ECD) Protein | HA210718 | >95% as determined by SDS-PAGE. | HEK293 |
| Human LILRB2/CD85d/ILT4, C-His, Flag Tag (ECD) Protein | HA210643 | >95% as determined by SDS-PAGE. | HEK293 |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|--|----------|--|--------|
| Human LILRB3/CD85a/ILT5 Protein | HA210106 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human LILRB4/ILT3, C-His, Flag Tag (ECD) Protein | HA210719 | >95% as determined by SDS-PAGE. | HEK293 |
| Human LILRB5/CD85c/LIR-8 Protein | HA210108 | > 95% as determined by Tris-Bis PAGE | HEK293 |
| Human Lipocalin-2/NGAL, C-His Tag Protein | HA211028 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Livin/BIRC7, C-His Tag Protein | HA211325 | >90% as determined by SDS-PAGE. | E.coli |
| Human LPA, Tag Free Protein | HA210979 | >95% as determined by SDS-PAGE. | HEK293 |
| Human LRRC15/LIB Protein | HA210017 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | CHO |
| Human L-selectin/CD62L, C-His Tag Protein | HA211181 | >95% as determined by SDS-PAGE. | HEK293 |
| Human LTA, Tag Free Protein | HA210742 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Lumican, C-His Tag Protein | HA211171 | >95% as determined by SDS-PAGE. | HEK293 |
| Human MAG/Siglec-4a, C-His, Flag Tag (ECD) Protein | HA210613 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Mammaglobin-A, C-His, Flag Tag Protein | HA210786 | >95% as determined by SDS-PAGE. | HEK293 |
| Human MAP4K1/HPK1, N-Twin Strep, C-His, Flag Tag Protein | HA210723 | >95% as determined by SDS-PAGE. | E.coli |
| Human MAPK12, C-His Tag Protein | HA210513 | >95% as determined by SDS-PAGE. | E.coli |
| Human MAPK15, N-Twin Strep, C-His, Flag Tag Protein | HA210715 | >90% as determined by SDS-PAGE. | E.coli |
| Human MAPK6, C-His Tag Protein | HA210591 | >95% as determined by SDS-PAGE. | E.coli |
| Human Mature TGF beta 1 Protein | HA210007 | > 95% as determined by Tris-Bis PAGE | HEK293 |
| Human MBL2, Tag Free Protein | HA210628 | >95% as determined by SDS-PAGE. | HEK293 |
| Human MCP-2/CCL8, Tag Free Protein | HA211228 | >95% as determined by SDS-PAGE. | HEK293 |
| Human MCP-4/CCL13, Tag Free Protein | HA211153 | >95% as determined by SDS-PAGE. | HEK293 |
| Human M-CSF/CSF1 (33-190), Tag Free Protein | HA210919 | >95% as determined by SDS-PAGE. | HEK293 |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|--|----------|--|--------|
| Human M-CSF/CSF1 (33-255), C His Tag Protein | HA210918 | >95% as determined by SDS-PAGE. | HEK293 |
| Human M-CSF/CSF-1 Protein | HA210090 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human M-CSF R/CSF1R/CD115 Protein | HA210098 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human M-CSFR/CSF1R/CD115, Tag Free (ECD) Protein | HA211149 | >90% as determined by SDS-PAGE. | HEK293 |
| Human MDH1, Tag Free Protein | HA211217 | >95% as determined by SDS-PAGE. | HEK293 |
| Human MERTK/Mer, C-His Tag (ECD) Protein | HA211172 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Mesothelin, C-His Tag Protein | HA210521 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Mesothelin/MSLN, Tag Free Protein | HA210542 | >95% as determined by SDS-PAGE. | HEK293 |
| Human MICA, C-His (ECD) Protein | HA210951 | >95% as determined by SDS-PAGE. | HEK293 |
| Human MICA Protein | HA210092 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human MICB, C-His Tag Protein | HA210907 | >95% as determined by SDS-PAGE. | HEK293 |
| Human MICB Protein | HA210091 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human MIF, C-His Tag Protein | HA210825 | >95% as determined by SDS-PAGE. | E.coli |
| Human MIP-5/CCL15, C-His Tag Protein | HA211207 | >95% as determined by SDS-PAGE. | HEK293 |
| Human MLH1, C-His Tag Protein | HA210652 | >95% as determined by SDS-PAGE. | E.coli |
| Human MMP12, C-Flag Tag Protein | HA210856 | >95% as determined by SDS-PAGE. | HEK293 |
| Human MMP13, C-His, Flag Tag Protein | HA210854 | >95% as determined by SDS-PAGE. | HEK293 |
| Human MMP1, C-Flag Protein | HA210903 | >95% as determined by SDS-PAGE. | HEK293 |
| Human MMP2/CLG4A, C-His, Flag Tag Protein | HA210660 | >95% as determined by SDS-PAGE. | HEK293 |
| Human MMP3, C-His, Flag Tag Protein | HA210853 | >95% as determined by SDS-PAGE. | HEK293 |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|--|----------|---|--------|
| Human MMP-7, Tag Free Protein | HA210989 | >90% as determined by SDS-PAGE. | HEK293 |
| Human MMP8/CLG1, C-His, Flag Tag Protein | HA210646 | >95% as determined by SDS-PAGE. | HEK293 |
| Human MMP-9, C-His, Flag Tag Protein | HA210698 | >95% as determined by SDS-PAGE. | HEK293 |
| Human MMP-9 Protein | HA210015 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human MPG, C-His Tag Protein | HA210765 | >95% as determined by SDS-PAGE. | E.coli |
| Human MPZL1, C-His Tag (ECD) Protein | HA210522 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Mucin-6/MUC6, Tag Free Protein | HA211305 | >95% as determined by SDS-PAGE. | HEK293 |
| Human MUSK, C-His, Flag Tag (ECD) Protein | HA210565 | >90% as determined by SDS-PAGE. | HEK293 |
| Human Myeloperoxidase/MPO, Tag Free Protein | HA211022 | >95% as determined by SDS-PAGE. | HEK293 |
| Human N-Cadherin, C-His Tag Protein | HA211124 | >95% as determined by SDS-PAGE. | HEK293 |
| Human NCAM1/CD56, Tag Free (ECD) Protein | HA210820 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Nectin-1/PVRL1/CD111 Protein | HA210118 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human Nectin-2/CD112 Protein | HA210119 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human Nectin-2, C-His Tag (ECD) Protein | HA210905 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Nectin-3/CD113 Protein | HA210097 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human Nectin-4, C-His,Flag Tag (ECD) Protein | HA210883 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Nectin-4 Protein | HA210163 | > 95% as determined by Tris-Bis PAGE > 95% as determined by SEC-HPLC | HEK293 |
| Human Nectin-4 Protein | HA210104 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human Neurogranin, Tag Free Protein | HA211156 | >90% as determined by SDS-PAGE. | HEK293 |
| Human Neuropeptide Y, Tag Free Protein | HA211142 | >95% as determined by SDS-PAGE. | HEK293 |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|---|----------|--|--------|
| Human Neuropilin-1/NRP1, C-His Tag (ECD) Protein | HA211017 | >95% as determined by SDS-PAGE. | HEK293 |
| Human NFL, Tag Free Protein | HA211087 | >95% as determined by SDS-PAGE. | HEK293 |
| Human NG2/MCSP/CSPG4, Tag Free (ECD) Protein | HA211047 | >95% as determined by SDS-PAGE. | HEK293 |
| Human NKG2D, Tag Free Protein | HA211290 | >95% as determined by SDS-PAGE. | HEK293 |
| Human NKp30/NCR3/CD337 Protein | HA210068 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human NKp46/NCR1/CD335 Protein | HA210183 | > 95% as determined by Tris-Bis PAGE> 95% as determined by HPLC | HEK293 |
| Human NKp46/NCR1, Tag Free Protein | HA211272 | >95% as determined by SDS-PAGE. | HEK293 |
| Human NKX2-2, N-Twin Strep, C-His, Flag Tag Protein | HA210755 | >90% as determined by SDS-PAGE. | E.coli |
| Human NKX3-1, N-Twin Strep, C-His, Flag Tag Protein | HA210568 | >95% as determined by SDS-PAGE. | E.coli |
| Human Noggin/NOG, Tag Free Protein | HA210726 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Notch 1 Protein | HA210074 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human NPTX1, Tag Free Protein | HA211176 | >95% as determined by SDS-PAGE. | HEK293 |
| Human NPTX2, Tag Free Protein | HA211175 | >95% as determined by SDS-PAGE. | HEK293 |
| Human NPTXR, C-His Tag Protein | HA211216 | >90% as determined by SDS-PAGE. | HEK293 |
| Human NR1H2, C-Strep II Tag Protein | HA210766 | >90% as determined by SDS-PAGE. | E.coli |
| Human NR4A2/NURR1, C-His Tag Protein | HA210687 | >95% as determined by SDS-PAGE. | E.coli |
| Human NRG1, C-Flag (176-246) Protein | HA210750 | >95% as determined by SDS-PAGE. | E.coli |
| Human NRG1, Tag Free Protein | HA210772 | >95% as determined by SDS-PAGE. | HEK293 |
| Human NSUN2, N-Strep II, C-His Tag Protein | HA210650 | >95% as determined by SDS-PAGE. | E.coli |
| Human NUDT21, N-Twin Strep, C-His, Flag Tag Protein | HA210632 | >95% as determined by SDS-PAGE. | E.coli |
| Human NUMA1 (1-153), C-His Tag Protein | HA210958 | >95% as determined by SDS-PAGE. | E.coli |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|--|----------|--|--------|
| Human NUMA1, C-His Tag Protein | HA210925 | >95% as determined by SDS-PAGE. | E.coli |
| Human Oncostatin M/OSM, C-His Tag Protein | HA211164 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Oncostatin M/OSM Protein | HA210087 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human Osteoactivin/GPNMB, C-His Tag (ECD) Protein | HA210623 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Osteopontin/OPN, C-His Tag Protein | HA210993 | >90% as determined by SDS-PAGE. | HEK293 |
| Human Osteoprotegerin / TNFRSF11B, C-His Tag Protein | HA210977 | >95% as determined by SDS-PAGE. | HEK293 |
| Human OTUB1, Tag Free Protein | HA210630 | >95% as determined by SDS-PAGE. | E.coli |
| Human OX40 Ligand/TNFSF4 Trimer Protein | HA210231 | > 90% as determined by Tris-Bis PAGE | HEK293 |
| Human OX40/TNFRSF4/CD134 Protein | HA210232 | > 95% as determined by Tris-Bis PAGE | HEK293 |
| Human OX40/TNFRSF4/CD134 Protein | HA210157 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human OX40/TNFRSF4, C-His Tag Protein | HA210892 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Oxytocin/OXT, Tag Free Protein | HA211134 | >95% as determined by SDS-PAGE. | HEK293 |
| Human PABPC1, N-Twin Strep, C-His, Flag Tag Protein | HA210588 | >95% as determined by SDS-PAGE. | E.coli |
| Human PAI-1/Serpin E1, C-His Tag Protein | HA210982 | >95% as determined by SDS-PAGE. | HEK293 |
| Human PARK7/DJ-1, Tag Free Protein | HA211158 | >95% as determined by SDS-PAGE. | HEK293 |
| Human PAX5, N-Twin Strep, C-His, Flag Tag Protein | HA210824 | >95% as determined by SDS-PAGE. | E.coli |
| Human PCNT, C-His Tag Protein | HA210575 | >95% as determined by SDS-PAGE. | HEK293 |
| Human PCSK9, C-His Tag Protein | HA210991 | >95% as determined by SDS-PAGE. | HEK293 |
| Human PCSK9 Protein | HA210215 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human PD-1, C-His Tag Protein | HA210922 | >95% as determined by SDS-PAGE. | HEK293 |
| Human PD-1/PDCD1 Protein | HA210141 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|---|----------|--|--------|
| Human PD-1/PDCD1 Protein | HA210219 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human PDE1B, N-Twin Strep, C-His, Flag Tag Protein | HA210579 | >95% as determined by SDS-PAGE. | E.coli |
| Human PDGF-AA, Tag Free Protein | HA211046 | >95% as determined by SDS-PAGE. | HEK293 |
| Human PDGF-B, Tag Free Protein | HA210859 | >95% as determined by SDS-PAGE. | HEK293 |
| Human PDGF R alpha/PDGFR A Protein | HA210024 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human PDGF R beta/CD140b Protein | HA210075 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human PD-L1/B7-H1, C-hFc Tag Protein | HA210873 | >95% as determined by SDS-PAGE. | HEK293 |
| Human PD-L1/B7-H1, C-His Tag Protein | HA210954 | >95% as determined by SDS-PAGE. | HEK293 |
| Human PD-L1/B7-H1 Protein | HA210257 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human PD-L1/B7-H1 Protein | HA210228 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human PDLIM3/ALP, N-Twin Strep, C-His, Flag Tag Protein | HA210621 | >90% as determined by SDS-PAGE. | E.coli |
| Human Pentraxin 3/PTX3, Tag Free Protein | HA210943 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Periostin/OSF-2, C-His, Flag Tag Protein | HA210681 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Peroxiredoxin 6/PRDX6, C-His Tag Protein | HA211191 | >95% as determined by SDS-PAGE. | HEK293 |
| Human PHGDH, N-Twin Strep, C-His, Flag Tag Protein | HA210678 | >95% as determined by SDS-PAGE. | E.coli |
| Human Plexin B2, C-His Tag (ECD) Protein | HA210550 | >95% as determined by SDS-PAGE. | HEK293 |
| Human PLK4, C-His Tag Protein | HA210831 | >95% as determined by SDS-PAGE. | E.coli |
| Human PLK4, N-Twin Strep, C-His, Flag Tag Protein | HA210672 | >95% as determined by SDS-PAGE. | E.coli |
| Human PMCA1/ATP2B1, C-His Tag Protein | HA210711 | >95% as determined by SDS-PAGE. | E.coli |
| Human PPID, C-His Tag Protein | HA210669 | >90% as determined by SDS-PAGE. | E.coli |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|---|----------|--|--------|
| Human PRDM16, C-His Tag Protein | HA210615 | >95% as determined by SDS-PAGE. | E.coli |
| Human PRDM1, N-Twin Strep, C-His, Flag Tag Protein | HA210703 | >95% as determined by SDS-PAGE. | E.coli |
| Human PRDX5, C-His Tag Protein | HA210506 | >95% as determined by SDS-PAGE. | E.coli |
| Human Progranulin/GRN, Tag Free Protein | HA211273 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Prolactin/PRL, C-His, Flag Tag Protein | HA210676 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Pro-neuropeptide Y (29-64 peptide) Protein | HA211303 | | |
| Human P-Selectin/CD62P, C-His Tag (ECD) Protein | HA210931 | >95% as determined by SDS-PAGE. | HEK293 |
| Human PSMA/FOLH1 Protein | HA210254 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human PTGR1, C-His Tag Protein | HA210657 | >95% as determined by SDS-PAGE. | E.coli |
| Human RAC1, C-Flag Tag Protein | HA210758 | >95% as determined by SDS-PAGE. | E.coli |
| Human RAGE, C-His Tag Protein | HA211263 | >95% as determined by SDS-PAGE. | HEK293 |
| Human RBM15B, N-Twin Strep, C-His, Flag Tag Protein | HA210553 | >90% as determined by SDS-PAGE. | E.coli |
| Human RBP4, C-His Tag Protein | HA210999 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Renin, C-His Tag Protein | HA211179 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Resistin, Tag Free Protein | HA211002 | >95% as determined by SDS-PAGE. | HEK293 |
| Human REST, Tag Free Protein | HA211258 | >95% as determined by SDS-PAGE. | HEK293 |
| Human RGMa Protein | HA210030 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human ROR1, Tag Free Protein | HA211231 | >95% as determined by SDS-PAGE. | HEK293 |
| Human ROR2/NTRKR2 Protein | HA210103 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human RORC, N-Twin Strep, C-His, Flag Tag Protein | HA210606 | >95% as determined by SDS-PAGE. | E.coli |
| Human R-spondin-1/RSPO1, Tag Free Protein | HA210730 | >95% as determined by SDS-PAGE. | HEK293 |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|--|----------|--|--------|
| Human S100A12,Tag Free Protein | HA211137 | >95% as determined by SDS-PAGE. | HEK293 |
| Human S100B, C-His Tag Protein | HA211097 | >95% as determined by SDS-PAGE. | E.coli |
| Human S100B, Tag Free Protein | HA211252 | >95% as determined by SDS-PAGE. | HEK293 |
| Human SAA1, C-Flag Tag Protein | HA210927 | >95% as determined by SDS-PAGE. | HEK293 |
| Human SCF/KITLG, C-His Tag (ECD) Protein | HA210618 | >95% as determined by SDS-PAGE. | HEK293 |
| Human SCF/KITLG, Tag Free Protein | HA211223 | >95% as determined by SDS-PAGE. | HEK293 |
| Human SCF/KITLG, Tag Free (26-189) Protein | HA211286 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Serpin A1/Alpha-1-antitrypsin, C-His Tag Protein | HA211106 | >95% as determined by SDS-PAGE. | HEK293 |
| Human SETD7/SET7/SET9, C-His Tag Protein | HA210701 | >95% as determined by SDS-PAGE. | E.coli |
| Human sFRP-1, Tag Free Protein | HA211218 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Siglec-10 Protein | HA210172 | > 95% as determined by Tris-Bis PAGE> 95% as determined by HPLC | HEK293 |
| Human Siglec-10 Protein | HA210093 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human Siglec-15/CD33L3 Protein | HA210220 | > 95% as determined by Tris-Bis PAGE | HEK293 |
| Human Siglec-15 Protein | HA210221 | > 95% as determined by SDS-PAGE. > 95% as determined by HPLC. | HEK293 |
| Human Siglec-2/CD22 Protein | HA210248 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human Siglec-3/CD33 Protein | HA210182 | > 95% as determined by Tris-Bis PAGE> 95% as determined by HPLC | HEK293 |
| Human Siglec-3/CD33 Protein | HA210129 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human Siglec-3/CD33, Tag Free Protein | HA210949 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Siglec-6/CD327 Protein | HA210130 | > 95% as determined by Tris-Bis PAGE | HEK293 |
| Human Siglec-7/CD328, C-His (ECD) Protein | HA210942 | >95% as determined by SDS-PAGE. | HEK293 |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|---|----------|--|--------|
| Human Siglec-9 Protein | HA210136 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human Siglec-9 Protein | HA210135 | > 95% as determined by Tris-Bis PAGE | HEK293 |
| Human SIGLEC9, Tag Free (ECD) Protein | HA211074 | >95% as determined by SDS-PAGE. | HEK293 |
| Human SIRP alpha/CD172a Protein | HA210251 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human SIRP alpha/CD172a Protein | HA210252 | > 95% as determined by Tris-Bis PAGE > 95% as determined by SEC-HPLC | HEK293 |
| Human SIRP alpha V8 Protein | HA210038 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human SIRPA, Tag Free Protein | HA211204 | >95% as determined by SDS-PAGE. | HEK293 |
| Human SIRP Beta/CD172b Protein | HA210081 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human SIRP Gamma/CD172g Protein | HA210082 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human SLAMF7/CRACC/CD319 Protein | HA210181 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human SLAMF7/CRACC/CD319 Protein | HA210096 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human SLC16A3/MCT4, N-Twin Strep, C-His, Flag Tag Protein | HA210557 | >95% as determined by SDS-PAGE. | E.coli |
| Human SLC1A4, N-Twin Strep, C-His, Flag Tag Protein | HA210658 | >95% as determined by SDS-PAGE. | E.coli |
| Human SLIT2, Tag Free Protein | HA211214 | >95% as determined by SDS-PAGE. | HEK293 |
| Human SMC2, Tag Free Protein | HA210590 | >95% as determined by SDS-PAGE. | E.coli |
| Human SMN1/SMN2, Tag Free Protein | HA211276 | >95% as determined by SDS-PAGE. | HEK293 |
| Human SMOC-1, C-His Tag Protein | HA210549 | >95% as determined by SDS-PAGE. | HEK293 |
| Human SNAP25, Tag Free Protein | HA211155 | >95% as determined by SDS-PAGE. | HEK293 |
| Human SOD1/Cu-Zn SOD, C-His Tag Protein | HA211121 | >95% as determined by SDS-PAGE. | E.coli |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|---|----------|--|--------|
| Human Sonic Hedgehog/SHH (C24II), Tag Free Protein | HA211279 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Sonic Hedgehog/SHH, Tag Free Protein | HA211278 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Sorcin/SRI, N-Twin Strep, C-His, Flag Tag Protein | HA210501 | >90% as determined by SDS-PAGE. | E.coli |
| Human SP110, N-Twin Strep, C-His, Flag Tag Protein | HA210674 | >95% as determined by SDS-PAGE. | E.coli |
| Human SPARC, C-His Tag Protein | HA210516 | >95% as determined by SDS-PAGE. | HEK293 |
| Human SP-D/SFTPD, C-His Tag Protein | HA211160 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Spondin-1/SPON1, C-His Tag Protein | HA210596 | >95% as determined by SDS-PAGE. | HEK293 |
| Human SPR, C-His Tag Protein | HA210511 | >95% as determined by SDS-PAGE. | E.coli |
| Human Syndecan-1/CD138, C-His Tag (ECD) Protein | HA210953 | >95% as determined by SDS-PAGE. | HEK293 |
| Human TACI/TNFRSF13B, Tag Free Protein | HA211248 | >95% as determined by SDS-PAGE. | HEK293 |
| Human TACSTD2/TROP2, C-His Tag (ECD) Protein | HA210737 | >95% as determined by SDS-PAGE. | HEK293 |
| Human TAF5/FAM19A5, Tag Free Protein | HA211192 | >95% as determined by SDS-PAGE. | HEK293 |
| Human TAG1/Contactin 2/CNTN2, Tag Free Protein | HA211194 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Tau/Tau441, C-His (isoform 2N4R) Protein | HA210937 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Tau/Tau758, C-His (isoform PNS-tau) Protein | HA211246 | >95% as determined by SDS-PAGE. | HEK293 |
| Human TBC1D4/AS160, N-Twin Strep, C-His, Flag Tag Protein | HA210612 | >90% as determined by SDS-PAGE. | E.coli |
| Human TFAP2C, N-Twin Strep, C-His, Flag Tag Protein | HA210722 | >95% as determined by SDS-PAGE. | E.coli |
| Human TFPI Protein | HA210003 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human TGF-alpha, Tag Free Protein | HA211318 | >90% as determined by SDS-PAGE. | HEK293 |
| Human TGF-beta RII, Tag Free (ECD) Protein | HA210555 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Thioredoxin/TXN, C-His Tag Protein | HA211189 | >95% as determined by SDS-PAGE. | HEK293 |
| Human TIE2/TEK, C-His, Flag Tag (ECD) Protein | HA210845 | >95% as determined by SDS-PAGE. | HEK293 |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|---|----------|---|--------|
| Human TIGIT Protein | HA210250 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human TIM-1/KIM-1/HAVCR , C-Flag Tag Protein | HA211313 | >95% as determined by SDS-PAGE. | HEK293 |
| Human TIM-3/HAVCR2, C-His, Flag Tag (ECD) Protein | HA210850 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Tim-3/HAVCR2 Protein | HA210253 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human Tim4/TIMD4, Tag Free (ECD) Protein | HA211014 | >95% as determined by SDS-PAGE. | HEK293 |
| Human TIMP1, Tag Free Protein | HA211003 | >95% as determined by SDS-PAGE. | HEK293 |
| Human TL1A/TNFSF15, Tag Free Protein | HA211266 | >95% as determined by SDS-PAGE. | HEK293 |
| Human TNF-alpha, C-His Tag Protein | HA210761 | >95% as determined by SDS-PAGE. | HEK293 |
| Human TNF-alpha, Tag Free Protein | HA210760 | >95% as determined by SDS-PAGE. | HEK293 |
| Human TNFR1/CD120a/TNFRSF1A Protein | HA210177 | > 95% as determined by Tris-Bis PAGE> 95% as determined by HPLC | HEK293 |
| Human TNFR2/CD120b/TNFRSF1B Protein | HA210176 | > 95% as determined by Tris-Bis PAGE> 95% as determined by HPLC | HEK293 |
| Human TNFR2/CD120b/TNFRSF1B Protein | HA210057 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human TNF RII/TNFRSF1B, C-His (ECD) Protein | HA210950 | >95% as determined by SDS-PAGE. | HEK293 |
| Human TNF RI/TNFRSF1A, Tag Free (ECD) Protein | HA210562 | >95% as determined by SDS-PAGE. | HEK293 |
| Human TNFSF11/RANKL/CD254, Tag Free Protein | HA211170 | >95% as determined by SDS-PAGE. | HEK293 |
| Human TNFSF15 Trimer Protein | HA210023 | > 95% as determined by Tris-Bis PAGE; > 95% as determined by HPLC | HEK293 |
| Human TNKS2, C-His Tag Protein | HA210560 | >95% as determined by SDS-PAGE. | E.coli |
| Human TOPBP1, C-His Tag Protein | HA210666 | >95% as determined by SDS-PAGE. | E.coli |
| Human tPA, C-His Tag Protein | HA211026 | >90% as determined by SDS-PAGE. | HEK293 |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|---|----------|--|--------|
| Human TPO/Thrombopoietin, C-His Tag Protein | HA211023 | >95% as determined by SDS-PAGE. | HEK293 |
| Human TRAF1, Tag Free Protein | HA211289 | >95% as determined by SDS-PAGE. | HEK293 |
| Human TRAF3, C-His Tag Protein | HA210601 | >95% as determined by SDS-PAGE. | E.coli |
| Human TRAIL/TNFSF10, C-His Tag Protein | HA210893 | >95% as determined by SDS-PAGE. | E.coli |
| Human Transthyretin/TTR, Tag Free Protein | HA211085 | >95% as determined by SDS-PAGE. | HEK293 |
| Human TREM1, C-His, Flag Tag Protein | HA210846 | >95% as determined by SDS-PAGE. | HEK293 |
| Human TREM2, N-Strep II, C-His Tag Protein | HA210647 | >90% as determined by SDS-PAGE. | HEK293 |
| Human TREM2, Tag Free Protein | HA211052 | >95% as determined by SDS-PAGE. | HEK293 |
| Human TrkC, C-His, Flag Tag (ECD) Protein | HA210620 | >95% as determined by SDS-PAGE. | HEK293 |
| Human TROP-2/TACSTD2 Protein | HA210256 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human TSLP, C-His Tag Protein | HA210952 | >95% as determined by SDS-PAGE. | HEK293 |
| Human TSLP Protein | HA210156 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human UBE2K, C-His Tag Protein | HA210508 | >95% as determined by SDS-PAGE. | E.coli |
| Human UCH-L1/PGP9.5, Tag Free Protein | HA211112 | >95% as determined by SDS-PAGE. | HEK293 |
| Human UFC1, His Tag Protein | HA210507 | >95% as determined by SDS-PAGE. | E.coli |
| Human ULBP1/NKG2DL1, C-His Tag Protein | HA210963 | >95% as determined by SDS-PAGE. | HEK293 |
| Human ULBP-1 Protein | HA210022 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Human ULBP3, C-His Tag Protein | HA210983 | >95% as determined by SDS-PAGE. | HEK293 |
| Human uPAR/PLAUR, Tag Free Protein | HA211102 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Uromodulin/UMOD, Tag Free Protein | HA211187 | >95% as determined by SDS-PAGE. | HEK293 |
| Human VCAM-1 / CD106, C-His Tag Protein | HA210976 | >95% as determined by SDS-PAGE. | HEK293 |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|---|----------|---|--------|
| Human VE-cadherin, Tag Free (ECD) Protein | HA210635 | >95% as determined by SDS-PAGE. | HEK293 |
| Human VEGF121, Tag Free Protein | HA210929 | >95% as determined by SDS-PAGE. | HEK293 |
| Human VEGF165, Tag Free Protein | HA210784 | >95% as determined by SDS-PAGE. | HEK293 |
| Human VEGF-C, C-His, Flag Tag Protein | HA210869 | >95% as determined by SDS-PAGE. | HEK293 |
| Human VEGFR1/FLT1 C-His Tag (ECD) Protein | HA211083 | >90% as determined by SDS-PAGE. | HEK293 |
| Human VEGF R1/FLT-1 Protein | HA210132 | > 95% as determined by Tris-Bis PAGE | HEK293 |
| Human VEGF R2/KDR Protein | HA210186 | > 95% as determined by Tris-Bis PAGE; > 95% as determined by HPLC | HEK293 |
| Human VEGFR2/KDR, Tag Free (ECD) Protein | HA211025 | >95% as determined by SDS-PAGE. | HEK293 |
| Human VEGFR3/FLT4, C-His Tag Protein | HA211107 | >90% as determined by SDS-PAGE. | HEK293 |
| Human VIPR2, C-hFc tag Protein | HA211224 | >95% as determined by SDS-PAGE. | HEK293 |
| Human VISTA/B7-H5/PD-1H, C-His Tag (ECD) Protein | HA210663 | >95% as determined by SDS-PAGE. | HEK293 |
| Human Von Willebrand Factor/vWF, Tag Free Protein | HA210992 | >90% as determined by SDS-PAGE. | HEK293 |
| Human VSIG4, C-His, Flag Tag (ECD) Protein | HA210847 | >95% as determined by SDS-PAGE. | HEK293 |
| Human VSNL1, C-His Tag Protein | HA210561 | >95% as determined by SDS-PAGE. | E.coli |
| Human YTHDC2, C-His Tag Protein | HA210524 | >95% as determined by SDS-PAGE. | E.coli |
| Human ZO1, N-Twin Strep, C-His, Flag Tag Protein | HA210528 | >95% as determined by SDS-PAGE. | E.coli |
| Monkey Albumin/ALB, C-His Tag Protein | HA211229 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse ACE, Tag Free Protein | HA211125 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse Adiponectin, C-His Tag Protein | HA210973 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse alpha-Synuclein, C-His Tag Protein | HA211105 | >95% as determined by SDS-PAGE. | E.coli |
| Mouse Angiotensinogen, C-His Tag Protein | HA211120 | >95% as determined by SDS-PAGE. | HEK293 |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|---|----------|--|--------|
| Mouse APOC3, Tag Free Protein | HA211230 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse ApoE, Tag Free Protein | HA211210 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse B7-2/CD86, C-His, Flag Tag Protein | HA210879 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse B7-H5/Gi24/VISTA Protein | HA210206 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Mouse BCL2L1, C-His Tag Protein | HA210828 | >95% as determined by SDS-PAGE. | E.coli |
| Mouse BCMA/TNFRSF17 Protein | HA210196 | > 95% as determined by Tris-Bis PAGE | HEK293 |
| Mouse BCMA/TNFRSF17 Protein | HA210197 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Mouse Cadherin-3/CDH3, Tag Free Protein | HA211188 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse CCL19/MIP-3b, Tag Free Protein | HA211277 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse CCL20/MIP-3, Tag Free Protein | HA211071 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse CCL21a. Tag free Protein | HA210813 | >95% as determined by SDS-PAGE. | E.coli |
| Mouse CCL2, Tag Free Protein | HA210775 | >90% as determined by SDS-PAGE. | HEK293 |
| Mouse CCL3/MIP-1 alpha, Tag Free Protein | HA210790 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse CCL5/RANTES, C-hFc Tag Protein | HA211095 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse CD24 Protein | HA210078 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Mouse CD25/IL-2R alpha, Tag Free Protein | HA211117 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse CD30/TNFRSF8, C-His Tag Protein | HA211062 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse CD36/SR-B3 Protein | HA210213 | > 95% as determined by Tris-Bis PAGE | HEK293 |
| Mouse CD3E&CD3D/CD3 epsilon&CD3 delta Protein | HA210145 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Mouse CD40 Ligand/TNFSF5 Trimer Protein | HA210210 | > 95% as determined by Tris-Bis PAGE | HEK293 |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|---|----------|--|--------|
| Mouse CD40/TNFRSF5 Protein | HA210212 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Mouse CD4, Tag Free (ECD) Protein | HA210863 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse CD68/SR-D1, C-His Tag (ECD) Protein | HA211131 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse Coagulation factor XII/F12, C-His Tag Protein | HA210984 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse Complement C3, C-His Tag Protein | HA211093 | >90% as determined by SDS-PAGE. | HEK293 |
| Mouse Complement C5, C-His Tag Protein | HA211091 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse Complement Factor D, C-His Tag Protein | HA211116 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse C-Reactive Protein/CRP, Tag Free | HA211020 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse CTGF/CCN2, Tag Free Protein | HA211094 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse CXCL10, Tag Free Protein | HA210805 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse CXCL12 / SDF-1, Tag Free Protein | HA211308 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse CXCL16, C-His Tag Protein | HA211059 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse CXCL1, Tag Free Protein | HA210792 | >95% as determined by SDS-PAGE. | E.coli |
| Mouse CXCL1, Tag Free (HEK293) Protein | HA211049 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse CXCL9/MIG, Tag Free Protein | HA211031 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse DKK1, C-His Tag Protein | HA211082 | >90% as determined by SDS-PAGE. | HEK293 |
| Mouse EGF, Tag Free Protein | HA211220 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse Erythropoietin/EPO, C-His Tag Protein | HA211068 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse Fas/CD95/TNFRSF6, Tag Free (ECD) Protein | HA211151 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse FGF-21, Tag Free Protein | HA211166 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse FGF basic/FGF2/bFGF, C-His Tag Protein | HA211283 | >95% as determined by SDS-PAGE. | E.coli |
| Mouse FGFR3, N Strep II, C-His Tag (ECD) Protein | HA210541 | >95% as determined by SDS-PAGE. | HEK293 |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|---|----------|--|--------|
| Mouse Flt-3 Ligand/FLT3LG, Tag Free Protein | HA211282 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse G-CSF, C-His Tag Protein | HA210972 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse GDF-15, Tag Free Protein | HA211089 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse GFRAL/GFR alpha-like Protein | HA210029 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Mouse GLP1R, C-His Tag Protein | HA211109 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse GM-CSF/CSF2, Tag Free Protein | HA210930 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse gp130/IL6ST, C-His Tag Protein | HA211196 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse Granzyme B, Tag Free Protein | HA211011 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse GRO Beta/CXCL2, Tag Free Protein | HA211013 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse HGF, Tag Free Protein | HA211169 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse HMGB1, Tag Free Protein | HA210800 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse ICAM-1/CD54, C-His Tag (ECD) Protein | HA211066 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse IFN-beta, Tag Free Protein | HA210789 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse IFN-gamma, C-His, Flag Tag Protein | HA210609 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse IFN-gamma, Tag Free Protein | HA210974 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse IGF1, Tag Free Protein | HA210936 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse IGFBP-2, Tag Free Protein | HA211197 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse IGFBP-4, Tag Free Protein | HA211240 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse IGFBP5, Tag Free Protein | HA211201 | >90% as determined by SDS-PAGE. | HEK293 |
| Mouse IgG Protein | HA1027 | Protein G affinity purified. | |
| Mouse IL-10, C-His, Flag Tag Protein | HA210569 | >90% as determined by SDS-PAGE. | HEK293 |
| Mouse IL-12/IL-23 p40, C-His Tag Protein | HA210890 | >95% as determined by SDS-PAGE. | HEK293 |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|--|----------|---------------------------------|--------|
| Mouse IL-12 (linked heterodimer), C-Flag Tag Protein | HA210785 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse IL-13, C-His, Flag Tag Protein | HA210571 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse IL-15, Tag Free Protein | HA211321 | >90% as determined by SDS-PAGE. | HEK293 |
| Mouse IL-17A, C-Flag Protein | HA210791 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse IL-17F, Tag Free Protein | HA211035 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse IL-18, Tag Free Protein | HA210940 | >95% as determined by SDS-PAGE. | E.coli |
| Mouse IL-1 alpha/IL1A, C-His Tag Protein | HA210971 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse IL-1 beta/IL1B, C-His, Flag Tag Protein | HA210706 | >95% as determined by SDS-PAGE. | E.coli |
| Mouse IL-1 beta/IL1B, C-His Tag Protein | HA210962 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse IL-1ra, C-His Tag Protein | HA210520 | >95% as determined by SDS-PAGE. | E.coli |
| Mouse IL-1ra/IL-1F3, Tag Free Protein | HA211041 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse IL-20, Tag Free Protein | HA211139 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse IL-21, Tag Free Protein | HA210801 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse IL-22, C-His Tag Protein | HA211016 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse IL-23, C-His Tag Protein | HA211060 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse IL-23 P19/IL23A, Tag Free Protein | HA211007 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse IL-25/IL-17E, Tag Free Protein | HA211130 | >90% as determined by SDS-PAGE. | HEK293 |
| Mouse IL-27 P28/ IL27A, Tag Free Protein | HA211005 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse IL-2, C-His, Flag Tag Protein | HA210653 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse IL-2, Tag Free Protein | HA210816 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse IL-33, Tag Free Protein | HA210768 | >95% as determined by SDS-PAGE. | E.coli |
| Mouse IL-3, Tag Free Protein | HA211168 | >95% as determined by SDS-PAGE. | HEK293 |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|---|----------|--|--------|
| Mouse IL-4, C-His, Flag Tag Protein | HA210741 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse IL-4, Tag Free Protein | HA210902 | >90% as determined by SDS-PAGE. | HEK293 |
| Mouse IL-5, Tag Free Protein | HA211018 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse IL-6, C-His Tag Protein | HA210667 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse IL-6, Tag Free Protein | HA210796 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse IL-7, Tag Free Protein | HA211019 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse IL-9, C-His, Flag Tag Protein | HA210631 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse Interferon alpha-1/IFNA1, C-His Tag Protein | HA210987 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse LAG3/CD223 Protein | HA210200 | > 95% as determined by Tris-Bis PAGE | HEK293 |
| Mouse Leptin, Tag Free Protein | HA210898 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse LIF, C-His Tag Protein | HA210735 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse L-selectin/CD62L, C-His Tag Protein | HA211126 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse LYVE1, C-His, Flag Tag (ECD) Protein | HA210844 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse MCSF/CSF-1, C-His Tag (33-187) Protein | HA211182 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse MCSF/CSF-1, C-His Tag (ECD) Protein | HA210725 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse MDC/CCL22, Tag Free Protein | HA211072 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse MMP-3, Tag Free Protein | HA211280 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse MMP-9, C-His Tag Protein | HA211045 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse Myeloperoxidase/MPO, C-His Tag Protein | HA211038 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse Nectin-4 Protein | HA210006 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Mouse Ninjurin-1, C-His, Flag Tag (ECD) Protein | HA210788 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse Noggin/NOG, Tag Free Protein | HA211281 | >95% as determined by SDS-PAGE. | HEK293 |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|---|----------|--|--------|
| Mouse Osteopontin/OPN, Tag Free Protein | HA210997 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse PCSK9, Tag Free Protein | HA211202 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse PD1/PDCD1/CD279, C-His, Flag Tag (ECD) Protein | HA210625 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse PD-1/PDCD1 Protein | HA21400 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Mouse PDGF-AA, Tag Free Protein | HA211136 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse PDGFBB, Tag Free Protein | HA211037 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse PD-L1, C-His, Flag Tag Protein | HA210689 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse Pro-Collagen I alpha 1/COL1A1, Tag Free Protein | HA211113 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse Prolactin, Tag Free Protein | HA211036 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse P-Selectin/CD62P, Tag Free (ECD) Protein | HA211140 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse Resistin, Tag Free Protein | HA211042 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse R-spondin-1/RSPO1, Tag Free Protein | HA211306 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse SCF, C-His Tag Protein | HA211067 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse Serpin A1a, C-His Tag Protein | HA211161 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse Serpin C1/Antithrombin-III, C-His Tag Protein | HA211123 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse ST2/IL-1 RL1, C-His Tag Protein | HA211193 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse Thrombopoietin / TPO, C-His Tag Protein | HA211311 | >90% as determined by SDS-PAGE. | HEK293 |
| Mouse Tim-3/HAVCR2 Protein | HA210211 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Mouse TNF-alpha, C-His, Flag Tag Protein | HA210527 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse TNF-alpha, Tag Free Protein | HA211054 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse TRANCE/RANK L/TNFSF11, N-GST, His Tag Protein | HA210876 | >90% as determined by SDS-PAGE. | E.coli |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|---|----------|--|--------|
| Mouse TRANCE/RANK L/TNFSF11, Tag Free Protein | HA210558 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse Transthyretin/TTR, C-His Tag Protein | HA211118 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse TREM2, C-Strep II Tag (ECD) Protein | HA210924 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse TREM2, Tag Free (ECD) Protein | HA211061 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse TROP-2/TACSTD2 Protein | HA210033 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Mouse TSLP, Tag Free Protein | HA211088 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse VCAM-1/CD106, C-His Tag (ECD) Protein | HA211043 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse VEGF164, Tag Free Protein | HA210780 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse VEGF188, Tag Free Protein | HA210770 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse VEGFR1/FLT1, C-His Tag Protein | HA211053 | >95% as determined by SDS-PAGE. | HEK293 |
| Mouse Von Willebrand Factor/vWF Protein | HA211245 | >95% as determined by SDS-PAGE. | HEK293 |
| Oryctolagus cuniculus SDC1, C-His, Flag Tag Protein | HA210638 | >95% as determined by SDS-PAGE. | HEK293 |
| Rabbit IgG Protein | HA1002 | Protein A affinity purified. | |
| Rat CCL2/MCP-1, Tag Free Protein | HA211114 | >95% as determined by SDS-PAGE. | HEK293 |
| Rat CCL3/MIP-1 alpha, Tag Free Protein | HA211056 | >95% as determined by SDS-PAGE. | HEK293 |
| Rat C-Reactive Protein/CRP, C-His Tag | HA211008 | >95% as determined by SDS-PAGE. | HEK293 |
| Rat CXCL1, C-His Tag Protein | HA211010 | >90% as determined by SDS-PAGE. | HEK293 |
| Rat CXCL2/MIP-2, Tag Free Protein | HA211132 | >95% as determined by SDS-PAGE. | |
| Rat GM-CSF/CSF2, Tag Free Protein | HA211147 | >95% as determined by SDS-PAGE. | HEK293 |
| Rat ICAM-1/CD54, C-His Tag (ECD) Protein | HA211063 | >95% as determined by SDS-PAGE. | HEK293 |
| Rat IFN-gamma, C-His, Flag Tag Protein | HA210577 | >95% as determined by SDS-PAGE. | HEK293 |
| Rat IGF1, Tag Free Protein | HA211009 | >90% as determined by SDS-PAGE. | HEK293 |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|--|----------|--|---|
| Rat IL-10, C-His, Flag Tag Protein | HA210714 | >90% as determined by SDS-PAGE. | HEK293 |
| Rat IL-18, Tag Free Protein | HA211199 | >95% as determined by SDS-PAGE. | HEK293 |
| Rat IL-1 alpha/IL1A, Tag Free Protein | HA211058 | >95% as determined by SDS-PAGE. | HEK293 |
| Rat IL-1 beta/IL1B, C-Flag Tag Protein | HA210708 | >95% as determined by SDS-PAGE. | E.coli |
| Rat IL-2, C-His, Flag Tag Protein | HA210683 | >95% as determined by SDS-PAGE. | HEK293 |
| Rat IL-4, C-His, Flag Tag Protein | HA210759 | >95% as determined by SDS-PAGE. | HEK293 |
| Rat IL-6, C-His, Flag Tag Protein | HA210629 | >95% as determined by SDS-PAGE. | HEK293 |
| Rat MMP-9, Tag Free Protein | HA211141 | >95% as determined by SDS-PAGE. | HEK293 |
| Rat TIMP1, Tag Free Protein | HA211163 | >95% as determined by SDS-PAGE. | HEK293 |
| Rat TNF-alpha, C-His, Flag Tag Protein | HA210538 | >95% as determined by SDS-PAGE. | HEK293 |
| Recombinant HIV-1 p24 protein | HA210278 | | Recombinant p24 Protein is expressed from E.coli with a TRX and His tag aa 133-363. |
| Recombinant HIV-1 p24 protein | HA210279 | >90% as determined by SDS-PAGE. | E.coli |
| SARS-COV-2 Nucleocapsid Protein | HA210045 | > 95% as determined by Tris-Bis PAGE | E.coli |
| SARS-COV-2 Nucleocapsid Protein CTD Domain Protein | HA210043 | > 95% as determined by Tris-Bis PAGE | E.coli |
| SARS-CoV-2 Nucleocapsid protein, N-His Tag | HAg2002 | >95% as determined by SDS-PAGE. | E.coli |
| SARS-COV-2 Nucleocapsid Protein NTD Domain Protein | HA210044 | > 95% as determined by Tris-Bis PAGE | E.coli |
| SARS-COV-2 Spike RBD Protein | HA210064 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| SARS-COV-2 Spike S1 Protein | HA210046 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| SARS-COV-2 Spike S1 Protein | HA210061 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| SARS-COV-2 Spike S Trimer Protein | HA210047 | > 95% as determined by Tris-Bis PAGE | HEK293 |

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

| | | | |
|--|----------|--|--------|
| SARS Spike RBD Protein | HA210053 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| SARS Spike S1 Protein | HA210055 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| SARS Spike S1 Protein | HA210052 | > 95% as determined by Tris-Bis PAGE > 95% as determined by HPLC | HEK293 |
| Streptomyces avidinii Streptavidin, C-His Tag Protein | HA210570 | >95% as determined by SDS-PAGE. | E.coli |
| Wnt surrogate, Tag Free Protein | HA210783 | >95% as determined by SDS-PAGE. | HEK293 |