

B细胞研究-产品清单

高级验证	产品名	货号	反应种属	验证应用
	ALCAM Mouse mAb [C7-E12]	EM1708-92	Human	WB,IF-Cell,IHC-P,FC
	ALCAM/CD166 Mouse mAb [7C2]	RT1029	Human, Mouse, Rat	WB,IP,IF,IHC-P
	AP-A Recombinant Rabbit mAb [JE46-03]	ET7109-45	Human	WB,IHC-P
	B7-1/CD80 Rabbit pAb	ER1904-08	Mouse, Rat	WB
	BLNK Recombinant Rabbit mAb [JE63-50]	HA721028	Human, Mouse, Rat	WB
	BLNK Mouse mAb [G1-F9]	EM1709-94	Human, Mouse	WB,IF-Cell,IHC-P,FC
	BOB1 Recombinant Rabbit mAb [JE47-78]	ET7109-68	Human	WB,IHC-P,FC
	BOB1 Rabbit pAb	ER1803-97	Human, Mouse, Rat	WB,FC
	BST2 Recombinant Rabbit mAb [JU11-93]	ET1706-46	Human	WB,IHC-P
	Calcineurin A Recombinant Rabbit mAb [JA64-11]	ET1704-02	Human, Mouse, Rat	WB,FC
	Calcineurin B Rabbit pAb	HA500346	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	CARD11 Recombinant Rabbit mAb [JE37-28]	HA722309	Human, Mouse, Rat	WB,IHC-P
Relative expression	CD10 Recombinant Rabbit mAb [JE60-45]	HA722864	Human	WB,IHC-P
Relative expression	CD10 Recombinant Rabbit mAb [SN75-07]	ET1611-82	Human, Mouse, Rat	WB,FC,IHC-P,IF-Cell,IF-Tissue
Relative expression	CD10 Recombinant Mouse mAb [A1G3-R] - BSA and Azide free	HA610121	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	CD10 Mouse mAb [A1G5]	EM1901-26	Human, Mouse	WB,IHC-P,IF-Tissue,FC
Relative expression	CD10 Recombinant Rabbit mAb [SN75-07] - BSA and Azide free	HA750273	Human, Mouse, Rat	WB,FC,IHC-P,IF-Cell,IF-Tissue
Relative expression	CD10 Recombinant Mouse mAb [A1G3-R]	HA601246	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	CD10 Mouse mAb	IRS101MS	Human	mIHC
	CD10 Mouse mAb [A1G4]	EM1901-25	Human	WB,IHC-P,FC
Relative expression	CD10 Rabbit pAb	0805-4	Human, Mouse, Rat	WB,IHC-P
	CD10 Mouse mAb	IRS222MS	Mouse	mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	CD10 Mouse mAb [A1G3]	EM1901-24	Human, Mouse	WB,IHC-P,FC,IF-Tissue
Relative expression	CD134 Recombinant Rabbit mAb [JE36-48]	HA722700	Human	WB
	CD134 Mouse mAb [2H1]	EM1701-24	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Relative expression	CD134 Recombinant Rabbit mAb [PSH10-14] - BSA and Azide free	HA751343	Human	WB,IF-Cell,FC,IP
Relative expression	CD134 Recombinant Rabbit mAb [PSH10-14]	HA723183	Human	WB,IF-Cell,FC,IP
KD	CD166 Mouse mAb [5-5]	EM1902-36	Human	WB,IHC-P,IF-Cell,FC
	CD166 Recombinant Rabbit mAb [JE57-47]	HA722705	Human, Mouse, Rat	WB,IHC-P,FC,IHC-Fr
KD	CD166 Mouse mAb [A6A9]	HA600082	Human	IF-Cell,IHC-P,WB
KD	CD166 Mouse mAb [A6A8]	HA600083	Human	IF-Cell,IHC-P,WB
KD	CD166 Mouse mAb [11-6]	M1012-7	Human, Mouse, Rat	WB,IHC-P,FC
Cell treatment	CD169 Recombinant Rabbit mAb [PSH05-21]	HA722246	Human, Mouse	WB,FC
Cell treatment	CD169 Recombinant Rabbit mAb [PSH05-21] - BSA and Azide free	HA750992	Human, Mouse	WB,FC
	CD180 Rabbit pAb	ER2001-06	Human, Mouse, Rat	WB,IHC-P,FC
	CD19 Recombinant mAb	IRS110RB	Mouse	mIHC
Relative expression	CD19 Recombinant Rabbit mAb [PSH03-94] - BSA and Azide free	HA750919	Mouse	WB,IHC-P,IF-Tissue,IP,IHC-Fr,FC,IF-Cell
	FITC Conjugated CD19 Recombinant Mouse mAb [PSH14-97]	HA601518F	Human	FC
Relative expression	CD19 Recombinant Rabbit mAb [PSH03-94]	HA722073	Mouse	WB,IHC-P,IF-Tissue,IP,IHC-Fr,FC,IF-Cell,mIHC
	FITC Conjugated CD19 Recombinant Rabbit mAb [PSH13-44]	HA723762F	Human	FC
	CD19 Recombinant Rabbit mAb	IRS018	Human	mIHC
	CD19 Recombinant Rabbit mAb [JF099-9]	ET1702-74	Human, Mouse, Rat	WB,FC
Relative expression	CD19 Recombinant Rabbit mAb [JF100-06]	ET1702-93	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC,mIHC
Relative expression	CD19 Recombinant Rabbit mAb [PSH13-44] - BSA and Azide free	HA751463	Human	WB,IF-Cell,FC
Relative expression	APC Conjugated CD19 Recombinant Rabbit mAb [PSH13-44]	HA723761F	Human	FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	CD19 Recombinant Rabbit mAb [PSH13-44]	HA723533	Human	WB,IF-Cell,FC
	PE Conjugated CD19 Recombinant Rabbit mAb [PSH13-44]	HA723763F	Human	FC
	CD19 Mouse mAb [A3-A8]	EM40308	Human, Mouse	WB,IHC-P,FC
Relative expression	CD19 Recombinant Rabbit mAb [JF100-06] - BSA and Azide free	HA750365	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
	CD19 Rabbit pAb	R1107-7	Human, Mouse	WB
Relative expression	CD20 Recombinant Rabbit mAb [PD00-02]	HA721138	Human, Mouse	IHC-P,mIHC,WB,IF-Cell,IP,IHC-Fr
Relative expression	CD20 Recombinant Rabbit mAb [PSH08-12] - BSA and Azide free	HA751207	Human	WB,IF-Cell,IHC-P,FC,IP
Relative expression	CD20 Recombinant Rabbit mAb [PSH08-12]	HA722953	Human	WB,IF-Cell,IHC-P,FC,IP
Relative expression	CD20 Recombinant Antibody [PD00-02] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610360	Human, Mouse	WB,IHC-P,IHC-Fr
	CD20 Rabbit pAb	IRS114RB	Mouse	mIHC
	APC Conjugated CD20 Recombinant Rabbit mAb [PSH08-12]	HA723838F	Human	FC
	CD20 Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS112MS	Human	mIHC
Relative expression	CD20 Recombinant Antibody [PD00-02] - Mouse IgG1 (Chimeric)	HA601544	Human, Mouse	WB,IHC-P,IHC-Fr
	PE Conjugated CD20 Recombinant Rabbit mAb [PSH08-12]	HA723839F	Human	FC
	CD20 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS113RT	Human	mIHC
Relative expression	CD20 Recombinant Rabbit mAb [SY12-01]	ET1605-33	Human, Monkey	WB,IHC-P,IF-Tissue,IF-Cell,FC,IP
	CD20 Recombinant Rabbit mAb	IRS003	Human	mIHC
	CD21 Recombinant Rabbit mAb	IRS019	Human	mIHC
	CD21 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS265RT	Human	mIHC
	CD21 Rabbit pAb	ER1706-56	Human	IF-Cell,IHC-P,FC
	CD21 Recombinant Rabbit mAb [PD00-23]	HA721163	Human	WB,IHC-P,IF-Cell,FC,mIHC,IF-Tissue
Relative expression	CD21 Mouse mAb [A4E5]	HA600001	Human	WB,IF-Cell,IHC-P
Relative expression	CD21 Recombinant Rabbit mAb [SC0681]	ET1610-61	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CD22 Recombinant Rabbit mAb [PO00-31]	HA721157	Human	WB,IHC-P,IF-Tissue
	CD22 Mouse mAb [A2G6]	EM1902-13	Human	WB,IHC-P,FC
	CD23 Rabbit pAb	HA500386	Human	IHC-P,FC
	CD23 Recombinant Rabbit mAb	IRS021	Human	mIHC
	CD23 Recombinant Rabbit mAb [PD00-03]	HA721139	Human	IHC-P,mIHC,IF-Tissue
	CD24 Rabbit pAb	R1303-1	Human	WB,IHC-P
Cell treatment, Relative expression	CD24 Recombinant Rabbit mAb [PSH10-85] - BSA and Azide free	HA751377	Human	IF-Cell,FC,WB
Cell treatment, Relative expression	CD24 Recombinant Rabbit mAb [PSH10-85]	HA723252	Human	IF-Cell,FC,WB
	CD24 Rabbit pAb	0804-4	Human	WB,IF-Cell,FC
Relative expression	CD27 Recombinant Rabbit mAb [JB40-98]	ET7107-73	Human, Mouse, Rat	WB,IHC-P
	CD27 Recombinant Rabbit mAb	IRS115RB	Human	mIHC
Relative expression	CD27 Recombinant Rabbit mAb [PSH08-90] - BSA and Azide free	HA751255	Human	WB,IF-Cell,IHC-P,FC,IP
	APC Conjugated CD27 Recombinant Rabbit mAb [PSH08-90]	HA723828F	Human	FC
Relative expression	CD27 Recombinant Rabbit mAb [PSH08-90]	HA723042	Human	WB,IF-Cell,IHC-P,FC,IP
	CD28 Rabbit pAb	ER2001-42	Human, Mouse, Rat	WB,IHC-P,FC
	CD3 Recombinant Rabbit mAb	IRS022	Human	mIHC
Relative expression	CD3 Recombinant Rabbit mAb [PSH19-70]	HA724094	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	CD3 Recombinant Rabbit mAb	IRS074	Mouse	mIHC
	CD3 Recombinant Mouse mAb [PD01-43]	HA601110	Human	WB,IHC-P,IF-Tissue
	CD3 Recombinant Rabbit mAb [JE80-02]	HA720082	Human	WB,IHC-P,IF-Tissue,IP,mIHC,IF-Cell
Relative expression	CD30 Recombinant Rabbit mAb [PSH04-12]	HA722093	Human	WB,IF-Cell,IHC-P,IF-Tissue
Relative expression	CD30 Recombinant Rabbit mAb [PSH04-11]	HA722092	Human	WB,IF-Cell,IHC-P
Relative expression	CD30 Recombinant Rabbit mAb [JE57-71]	HA721579	Human	WB,IHC-P
Relative expression	CD30 Recombinant Rabbit mAb [PSH04-10]	HA722091	Human	WB,IHC-P,FC,IF-Tissue
	CD32 Rabbit pAb	ER62945	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	CD32 Recombinant Rabbit mAb [JM10-70] - BSA and Azide free	HA750382	Human	WB,IHC-P,IF-Cell,FC,IP
Cell treatment	CD32 Recombinant Rabbit mAb [JM10-70]	ET1703-43	Human	WB,IHC-P,IF-Cell,FC,IP
	CD32A + CD32B + CD32C Recombinant Rabbit mAb [JE64-17]	HA721930	Human	WB,IHC-P
	CD32B Mouse mAb [A9C11]	HA601121	Human	WB,IHC-P
	CD32-C Rabbit pAb	ER62946	Human	WB
Relative expression	CD38 Recombinant Antibody [PD01-49] - Rat IgG1 (Chimeric)	HA601417		
	CD38 Recombinant Antibody [PD01-49] - Mouse IgG1 (Chimeric)	HA601554	Human	WB,IHC-P,FC
	CD38 Recombinant Rabbit mAb [JE29-35]	ET7110-53	Human	WB,IHC-P,FC,IF-Tissue
	CD38 Recombinant Rabbit mAb [SN07-20]	ET1611-66	Human	WB,IHC-P,FC,IF-Tissue
	CD38 Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS120MS	Human	mIHC
	CD38 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS119RT	Human	mIHC
Relative expression	CD38 Recombinant Antibody [PD01-49] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610257		
	CD38 Recombinant Antibody [PD01-49] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610369	Human	WB,IHC-P,FC
	CD38 Recombinant Rabbit mAb	IRS006	Human	mIHC
Relative expression	CD38 Recombinant Rabbit mAb [PD01-49]	HA721268	Human	WB,IHC-P,IF-Cell,FC,mIHC,IF-Tissue
	CD39 Recombinant Rabbit mAb [JA90-36] - BSA and Azide free	HA750421	Human, Mouse, Rat	WB,IHC-P
	CD39 Recombinant Rabbit mAb [JA90-36]	ET1704-74	Human, Mouse, Rat	WB,IHC-P
	CD3D Recombinant Rabbit mAb [JJ08-97]	ET1701-57	Human	WB,IF-Cell,IHC-P,FC,IP,IF-Tissue
	CD3E Rabbit pAb	HA500344	Human, Mouse, Rat	WB
	CD3 epsilon Rabbit pAb	ER80501	Human, Mouse, Rat	WB
	CD3 epsilon Recombinant Rabbit mAb [SY0239]	ET1607-29	Human	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP
	CD3 epsilon Recombinant Rabbit mAb [JE17-74]	ET7110-63	Human	WB,IHC-P,FC
	CD3 epsilon Recombinant Rabbit mAb [SY0239] - BSA and Azide free	HA750115	Human	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CD3G Recombinant Rabbit mAb [JB38-29]	ET7107-55	Human, Mouse	WB,IHC-P,FC,IP,IF-Cell,IF-Tissue
	CD3 zeta Rabbit pAb	0805-7	Human	WB,IF-Cell,FC
Relative expression	CD3 zeta Recombinant Rabbit mAb [JE09-37]	HA722692	Human, Rat	WB,IHC-P,IF-Tissue
	CD4 Rabbit pAb	R1211-1	Human	WB,IF-Cell
Relative expression	CD4 Recombinant Rabbit mAb [PSH08-22]	HA722966	Mouse	WB,IF-Cell,IHC-P,mIHC,IHC-Fr
	CD4 Recombinant Rabbit Multiclonal Antibody [PSH13-24]	HA723511	Human, Mouse	WB,IF-Cell,IHC-P
Relative expression	CD4 Recombinant Rabbit mAb [PSH08-22] - BSA and Azide free	HA751213	Mouse	WB,IF-Cell,IHC-P,IHC-Fr
	CD4 Recombinant Rabbit mAb [JE56-36]	ET7110-64	Human	WB,IHC-P
Relative expression	CD4 Recombinant Antibody [PSH08-22] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610353	Mouse	WB,IF-Tissue,IHC-P,IF-Cell
KD	CD4 Recombinant Rabbit mAb [ST0488]	ET1609-52	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC,mIHC
KD	CD4 Rabbit pAb	0108-6	Human, Mouse, Rat	WB,IF-Cell,FC
	CD4 Recombinant Rabbit Multiclonal Antibody [PSH13-24] - BSA and Azide free	HA751460	Human, Mouse	WB,IF-Cell,IHC-P
	CD4 Recombinant Antibody [ST0488] - Mouse IgG1 (Chimeric)	HA601557	Human	WB,IF-Cell,IHC-P
	CD4 Recombinant Rabbit mAb	IRS005	Human	mIHC
	CD4 Rabbit pAb	ER1803-84	Human	WB,FC
	CD4 Recombinant Antibody [PSH08-22] - Rat IgG1 (Chimeric)	HA601534	Mouse	WB,IF-Tissue,IHC-P,IHC-Fr
Relative expression	CD4 Recombinant Antibody [PSH08-22] - Mouse IgG1 (Chimeric)	HA601533	Mouse	WB,IF-Tissue,IHC-P,IF-Cell
	CD4 Recombinant Antibody [PSH08-22] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610354	Mouse	WB,IF-Tissue,IHC-P,IHC-Fr
	CD4 Recombinant Rabbit mAb	IRS051	Mouse	mIHC
	CD40 Recombinant Rabbit mAb [PSH10-01] - BSA and Azide free	HA751331	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-P,IHC-Fr,IF-Cell,FC,IF-Tissue
Relative expression	CD40 Recombinant Rabbit mAb [PSH05-90]	HA722488	Human, Cynomolgus monkey, Pig	WB,IF-Cell,FC
	CD40 Recombinant Rabbit mAb [PSH10-01]	HA723170	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-P,IHC-Fr,IF-Cell,FC,IF-Tissue
	CD40 Mouse mAb [A1E5]	EM1901-41	Human	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	CD40 Recombinant Rabbit mAb [PSH05-90] - BSA and Azide free	HA751037	Human, Cynomolgus monkey, Pig	WB,IF-Cell,FC
Relative expression	CD40 Mouse mAb [A1E6]	EM1901-40	Human	WB,IF-Cell,IHC-P
	CD40 Rabbit pAb	ER1803-54	Human, Rat, Mouse	WB,IF-Cell,IHC-P,FC
	CD45 Recombinant Rabbit mAb	IRS025	Human	mIHC
Relative expression	CD45 Recombinant Rabbit Multiclonal Antibody [PSH14-19] - BSA and Azide free	HA751510	Human, Mouse	WB,IHC-P,IF-Cell
	CD45 Rabbit pAb	ER1905-93	Human	FC
	CD45 Rabbit pAb	ER1901-29	Human	WB,IHC-P,FC
Relative expression	CD45 Recombinant Rabbit Multiclonal Antibody [PSH14-19]	HA723606	Human, Mouse	WB,IHC-P,IF-Cell
	CD45 Rabbit pAb	ER1902-02	Human	WB,FC
	CD45 Recombinant Rabbit mAb [JE03-05]	ET7111-03	Human	WB,IHC-P,FC,IF-Tissue,mIHC
	CD45RO Recombinant Mouse mAb	IRS124MS	Human	mIHC
	PE Conjugated CD45RO Recombinant Mouse mAb [PSH13-94]	HA601531F	Human	FC
	CD45RO Recombinant Mouse mAb [PSH13-94]	HA601441	Human	IHC-P,FC
	CD5 Recombinant Rabbit mAb [PD00-01]	HA721137	Human	IHC-P,IF-Tissue,WB
	CD5 Recombinant Rabbit mAb [JE25-41]	ET7110-56	Human	WB,IHC-P
	CD74 Recombinant Rabbit mAb [JF40-10]	ET1702-51	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
Relative expression	CD74 Recombinant Mouse mAb [PDM0-22]	HA601117	Human	WB,IHC-P,IF-Cell,IF-Tissue
	CD79a Recombinant Rabbit mAb [JE32-14]	HA722896	Human	WB,IF-Cell,IHC-P,IF-Tissue
	CD79a Rabbit pAb	R1204-4	Human, Mouse	WB,IF-Cell
	CD79a Mouse mAb [C1-E8]	EM1902-29	Human, Rat, Mouse	WB,FC
	CD79a Recombinant Rabbit mAb [JE43-50]	ET7109-33	Human	WB,IF-Tissue,IHC-P,FC
	CD79a Recombinant Mouse mAb [PDM0-23]	HA601118	Human	WB,IHC-P
	CD79b Recombinant Rabbit mAb [JE58-34]	HA720043	Human	WB,IHC-P,FC
Relative expression	CD80 Recombinant Rabbit mAb [PSH09-79]	HA723144	Human	WB,IF-Cell,FC
Relative expression	CD80 Recombinant Rabbit mAb [PSH09-79] - BSA and Azide free	HA751321	Human	WB,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CD80 Recombinant Mouse mAb [8-E5-R] - BSA and Azide free	HA610122	Human	WB
	CD80 Rabbit pAb	1007-8	Human	WB,IF-Cell
Relative expression	CD80 Recombinant Rabbit mAb [JF100-4]	ET1702-95	Human	WB,IHC-P
Relative expression, Cell treatment	CD80 Recombinant Rabbit mAb [PSH18-13]	HA723976	Mouse	WB,IF-Cell,FC,IP
Relative expression	CD80 Mouse mAb [8-E5]	M1007-10	Human, Mouse	WB,IHC-P,FC
	CD80 Recombinant Mouse mAb [8-E5-R]	HA601247	Human	WB
	CD80 Rabbit pAb	ER1905-99	Human, Mouse, Rat	WB
Relative expression	CD81 Recombinant Rabbit mAb [PSH10-26]	HA723193	Human, Rat	WB,IHC-P,FC
	CD83 Rabbit pAb	ER62949	Human	WB,IHC-P
	CD83 Recombinant Rabbit mAb [JE33-61]	HA721556	Human	WB
	CD84 Recombinant Rabbit mAb [JE59-27]	HA720037	Human	WB,IF-Cell,IHC-P,FC
	CXCR5 Recombinant Rabbit mAb [JB11-40] - BSA and Azide free	HA750463	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell
	CXCR5 Recombinant Rabbit mAb [JB11-40]	ET7107-74	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell
Cell treatment	DPP4 / CD26 Recombinant Rabbit mAb [PSH12-59]	HA723473	Human	WB,IHC-P,IF-Cell,FC,IP
KO	DPP4 / CD26 Recombinant Rabbit mAb [JM11-42] - BSA and Azide free	HA750396	Human, Mouse, Rat	WB
Cell treatment	DPP4 / CD26 Recombinant Rabbit mAb [PSH12-59] - BSA and Azide free	HA751446	Human	WB,IHC-P,IF-Cell,FC,IP
KO	DPP4 / CD26 Recombinant Rabbit mAb [JM11-42]	ET1703-93	Human, Mouse, Rat	WB
	Flt3 Rabbit pAb	ER2001-25	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell
	FLT3 Mouse mAb [C3-B7]	EM1709-22	Human	WB,IHC-P
	FITC Conjugated Human CD10 Recombinant Rabbit mAb [SN75-07]	HA720209F	Human	FC
	Human CD27/TNFRSF7 Recombinant Rabbit mAb [PSH08-32] - BSA and Azide free (Capture)	HA722974	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human CD27/TNFRSF7 Recombinant Rabbit mAb [PSH08-33] - Detector	HA722976B	Human	ELISA(Det),ELISA
	Human CD27/TNFRSF7 Recombinant Rabbit mAb [PSH08-33] - BSA and Azide free (Detector)	HA722975	Human	ELISA(Det)
	Human CD27/TNFRSF7 Recombinant Rabbit mAb [PSH08-31] - BSA and Azide free (Capture)	HA722973	Human	ELISA(Cap)
	Biotin Conjugated Human CD28 Recombinant Rabbit mAb [PSH12-51] - Detector	HA723459B	Human	ELISA(Det),ELISA
	Human CD28 Recombinant Rabbit mAb [PSH12-51] - BSA and Azide free (Detector)	HA723458	Human	ELISA(Det)
	Human CD28 Recombinant Rabbit mAb [PSH12-50] - BSA and Azide free (Capture)	HA723457	Human	ELISA(Cap)
	PE Conjugated Human CD3 Recombinant Mouse mAb [PSH04-92]	HA600117F	Human	FC
	APC Conjugated Human CD3 Recombinant Mouse mAb [PSH04-92]	HA600118F	Human	FC
	Human CD3 Recombinant Mouse mAb [PSH04-92]	HA601337	Human	FC
	Percp-Cy5.5 Conjugated Human CD3 Recombinant Mouse mAb [PSH04-92]	HA601548F	Human	FC
	FITC Conjugated Human CD3 Recombinant Mouse mAb [PSH04-92]	HA600116F	Human	FC
	Human CD30 Recombinant Rabbit mAb [PSH04-82] - BSA and Azide free (Detector)	HA722188	Human	ELISA(Det)
	Human CD30 Recombinant Rabbit mAb [PSH04-81] - BSA and Azide free (Capture)	HA722187	Human	ELISA(Cap)
	Biotin Conjugated Human CD30 Recombinant Rabbit mAb [PSH04-82]	HA722601B	Human	ELISA(Det)
	FITC Conjugated Human CD38 Recombinant Rabbit mAb [SN07-20]	HA720196F	Human	FC
	APC Conjugated Human CD4 Recombinant Mouse mAb [PSH04-93]	HA600121F	Human	FC
	FITC Conjugated Human CD4 Recombinant Mouse mAb [PSH04-93]	HA600119F	Human	FC
	PE Conjugated Human CD4 Recombinant Mouse mAb [PSH04-93]	HA600120F	Human	FC
	PE-Cy7 Conjugated Human CD4 Recombinant Mouse mAb [PSH04-93]	HA601542F	Human	FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human CD4 Recombinant Mouse mAb [PSH04-93]	HA601338	Human	FC
	iFluor™ 488 Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600109F	Human	FC
	Human CD45 Recombinant Mouse mAb [PSH03-64]	HA601285	Human	IF-Cell,FC
	FITC Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600112F	Human	FC
	FITC Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600125F	Human	FC
	APC Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600127F	Human	FC
	PE Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600126F	Human	FC
	iFluor™ 555 Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA601556F	Human	FC
	Human CD45 Recombinant Mouse mAb [PSH03-64] - BSA and Azide free	HA610158	Human	IF-Cell,FC
	Percp-Cy5.5 Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600128F	Human	FC
	Human CD45RA Recombinant Mouse mAb [PSH04-97]	HA601341	Human	FC,IF-Cell,IHC-P
	FITC Conjugated Human CD45RA Recombinant Mouse mAb [PSH04-97]	HA600131F	Human	FC
	PE Conjugated Human CD45RA Recombinant Mouse mAb [PSH04-97]	HA601525F	Human	FC
	FITC Conjugated Human CD5 Recombinant Rabbit mAb [PSH05-01]	HA720211F	Human	FC
Relative expression	Human CD5 Recombinant Rabbit mAb [PSH05-01]	HA722210	Human	IF-Cell,FC
	Human CD80 Recombinant Rabbit mAb [PSH09-62] - BSA and Azide free (Detector)	HA723124	Human	ELISA(Det)
	Human CD80 Recombinant Rabbit mAb [PSH09-61] - BSA and Azide free (Capture)	HA723123	Human	ELISA(Cap)
	Biotin Conjugated Human CD80 Recombinant Rabbit mAb [PSH09-63] - Detector	HA723127B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human CD80 Recombinant Rabbit mAb [PSH09-62] - Detector	HA723126B	Human	ELISA(Det),ELISA
	Human CD80 Recombinant Rabbit mAb [PSH09-63] - BSA and Azide free (Detector)	HA723125	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human CD83 Recombinant Rabbit mAb [PSH20-08] - BSA and Azide free (Capture)	HA725305	Human	ELISA(Cap)
	Biotin Conjugated Human CD83 Recombinant Rabbit mAb [PSH20-09] - Detector	HA725306B	Human	ELISA(Det),ELISA
	Human CD83 Recombinant Rabbit mAb [PSH20-09] - BSA and Azide free (Detector)	HA725306	Human	ELISA(Det)
	Human IgG4 Mouse mAb [2G1]	RT1315	Human	ELISA
	HRP Conjugated Mouse Anti-Human IgG heavy chain [54-76-R] monoclonal Antibody	HA601284	Human	WB
	Human IgG heavy chain Recombinant Mouse mAb [54-76-R]	HA601283	Human	WB,ELISA,IHC-P
	Human IgG heavy chain Mouse mAb [54-76]	M0806-1	Human	WB,ELISA
	human IgM Mouse mAb [2-B5]	M1510-7	Human	WB,ELISA
	human IgM Mouse mAb [B7-3]	M1510-5	Human	WB,IHC-P,ELISA
	human IgM Mouse mAb [B6-8]	M1510-6	Human	WB,ELISA
	human IgM Mouse mAb [9-C1]	M1510-9	Human	WB,IHC-P,ELISA
	Human IgM Rabbit pAb	R1511-23	Human	WB,IHC-P
Cell treatment	Human IgM Recombinant Rabbit mAb [JF84-09]	ET1702-81	Human	WB,IHC-P,IP,FC,IF-Cell
	Human IgM Mouse mAb [4-B2]	M1510-8	Human	WB,IHC-P,ELISA
	Human IL-5 Recombinant Rabbit mAb [PSH04-65] - BSA and Azide free (Detector)	HA722165	Human	ELISA(Det)
	Human IL-5 Recombinant Rabbit mAb [PSH04-64] - BSA and Azide free (Capture)	HA722164	Human	ELISA(Cap)
	Biotin Conjugated Human IL-5 Recombinant Rabbit mAb [PSH04-65] - Detector	HA722609B	Human	ELISA(Det)
	Human OX40/TNFRSF4 Recombinant Rabbit mAb [PSH10-17] - BSA and Azide free (Capture)	HA723187	Human	ELISA(Cap)
	Human OX40/TNFRSF4 Recombinant Rabbit mAb [PSH10-18] - BSA and Azide free (Detector)	HA723188	Human	ELISA(Det)
	Biotin Conjugated Human OX40/TNFRSF4 Recombinant Rabbit mAb [PSH10-18] - Detector	HA723189B	Human	ELISA(Det),ELISA
	IgA Recombinant Rabbit mAb [JM10-42]	ET1703-10	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	IgA Mouse mAb [3-B6]	M1510-31	Human	WB,ELISA
	IgA Mouse mAb [5-C3]	M1510-30	Human	WB,ELISA
	IgA Mouse mAb [3-F5]	HA601137	Human	ELISA,WB
	IgG4 Recombinant Rabbit mAb [ST05-12]	ET1609-56	Human	WB,IHC-P
	IL-5 Rabbit pAb	HA500081	Human	WB,IF-Cell,ELISA
	intestinal alkaline phosphatase Recombinant Rabbit mAb [PSH17-83]	HA723952	Human, Mouse, Rat	IHC-P
	LY75 Recombinant Rabbit mAb [JB87-35]	ET7107-58	Human, Mouse	WB,IHC-P,IP
	Mouse anti-Human IgG Mouse mAb [A5D12]	HA600076	Human	WB,ELISA
	Mouse anti-Human IgG Recombinant Mouse mAb [A5D12-R]	HA601220	Human	WB,ELISA
	PE Conjugated Anti-Mouse CD3 Recombinant Antibody [PSH04-98] - Rabbit IgG (Chimeric)	HA720200F	Mouse	FC
	APC Conjugated Anti-Mouse CD3 Recombinant Antibody [PSH04-98] - Rabbit IgG (Chimeric)	HA720201F	Mouse	FC
	FITC Conjugated Anti-Mouse CD3 Recombinant Antibody [PSH04-98] - Rabbit IgG (Chimeric)	HA720199F	Mouse	FC
	Anti-Mouse CD3 Recombinant Antibody [PSH04-98] - Rabbit IgG (Chimeric)	HA723347	Mouse	FC
	PE Conjugated Anti-Mouse CD4 Recombinant Antibody [PSH04-99] - Rabbit IgG (Chimeric)	HA720203F	Mouse	FC
	APC Conjugated Anti-Mouse CD4 Recombinant Antibody [PSH04-99] - Rabbit IgG (Chimeric)	HA720204F	Mouse	FC
	FITC Conjugated Anti-Mouse CD4 Recombinant Antibody [PSH04-99] - Rabbit IgG (Chimeric)	HA720202F	Mouse	FC
	PTPRC Mouse mAb [C11-C10]	EM1708-69	Human	WB,IHC-P
	PU.1/Spi1 Recombinant Rabbit mAb [JB72-35]	ET7107-49	Human	WB,IHC-P,FC
	PU.1 / Spi1 Recombinant Rabbit mAb [PD00-14]	HA721162	Human	WB,IHC-P
	HRP Conjugated Rabbit anti-Human IgG polyclonal Antibody	HA1024	Human	WB,ELISA
KD	TAPA1/CD81 Recombinant Rabbit mAb [SN206-01]	ET1611-87	Human, Mouse, Rat	WB,IF-Cell,IHC- P,FC
Relative expression	T-bet / Tbx21 Recombinant Rabbit mAb [JE60-20]	HA721901	Human	WB,IHC-P,IF- Cell,FC
	T-bet / Tbx21 Rabbit pAb	ER1802-65	Human	WB,FC,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	TRK fused gene Recombinant Rabbit mAb [SD082-0]	ET1612-52	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
----	---	-----------	-------------------	----------------------------