

NK细胞研究-产品清单

高级验证	产品名	货号	反应种属	验证应用
	2B4 Rabbit pAb	HA500263	Human, Rat	WB
	2B4 Rabbit pAb	HA500180	Human	IHC-P
	B3GAT1 Rabbit pAb	ER1901-83	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	B3GAT1 Rabbit pAb	ER1901-15	Human, Mouse, Rat	WB,IHC-P,FC
	B3GAT1 Rabbit pAb	HA500443	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	CD10 Mouse mAb	IRS101MS	Human	mIHC
	CD10 Mouse mAb [A1G5]	EM1901-26	Human, Mouse	WB,IHC-P,IF-Tissue,FC
Relative expression	CD10 Recombinant Rabbit mAb [SN75-07]	ET1611-82	Human, Mouse, Rat	WB,FC,IHC-P,IF-Cell,IF-Tissue
Relative expression	CD10 Mouse mAb [A1G3]	EM1901-24	Human, Mouse	WB,IHC-P,FC,IF-Tissue
Relative expression	CD10 Recombinant Mouse mAb [A1G3-R] - BSA and Azide free	HA610121	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	CD10 Mouse mAb	IRS222MS	Mouse	mIHC
	CD10 Mouse mAb [A1G4]	EM1901-25	Human	WB,IHC-P,FC
Relative expression	CD10 Recombinant Rabbit mAb [JE60-45]	HA722864	Human	WB,IHC-P
Relative expression	CD10 Recombinant Mouse mAb [A1G3-R]	HA601246	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Relative expression	CD10 Recombinant Rabbit mAb [SN75-07] - BSA and Azide free	HA750273	Human, Mouse, Rat	WB,FC,IHC-P,IF-Cell,IF-Tissue
Relative expression	CD10 Rabbit pAb	0805-4	Human, Mouse, Rat	WB,IHC-P
KD	CD100 Recombinant Rabbit mAb [JE63-21]	HA721013	Human	WB
	CD100 Mouse mAb [5H6E3]	EM1710-47	Human	WB,IHC-P,FC
Relative expression	CD11a Recombinant Rabbit mAb [SC56-01]	ET1610-1	Human	WB,IF-Tissue,IHC-P,IP
Relative expression	CD11b Recombinant Rabbit mAb [PSH03-96] - BSA and Azide free	HA750921	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	CD11b Recombinant Rabbit mAb	IRS014	Human	mIHC
	CD11b Mouse mAb [4-6]	EM1701-42	Human	IHC-P,IF-Cell,FC,WB,IF-Tissue
Relative expression	CD11b Recombinant Rabbit mAb [JE01-17]	HA722856	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	CD11b Mouse mAb [A2-9]	M1211-6	Human	WB,IHC-P,FC,IF-Tissue
Relative expression	CD11b Recombinant Rabbit mAb [PSH03-96]	HA722075	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	CD11b Mouse mAb [A2-B5-D2]	EM1902-35	Human, Mouse	WB
Relative expression	CD11b Recombinant Rabbit mAb [JU93-81]	ET1706-04	Human	WB,IP,IHC-P,IF-Cell,IF-Tissue,mIHC
Relative expression	CD11b Recombinant Rabbit mAb [JU93-81] - BSA and Azide free	HA750451	Human	WB,IP,IHC-P,IF-Cell,IF-Tissue
	CD11b/ITAM Mouse mAb [4G2]	RT1107	Human, Mouse	WB,IP,IF,FC
	CD11c Recombinant Rabbit mAb	IRS001	Human	mIHC
	CD11c Recombinant Antibody [SI19-06] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610258		
	CD11c Recombinant Antibody [SI19-06] - Mouse IgG1 (Chimeric)	HA601545	Human	IF-Tissue,IHC-P
Relative expression	CD11c Recombinant Rabbit mAb [PSH07-42]	HA722830	Mouse	WB,IHC-P,IHC-Fr,IF-Tissue
	CD11c Recombinant Antibody - Rat IgG1 (Chimeric)	IRS103RT	Human	mIHC
	CD11c Recombinant Antibody [SI19-06] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610361	Human	IF-Tissue,IHC-P
Relative expression	CD11c Recombinant Rabbit mAb [PSH07-42] - BSA and Azide free	HA751163	Mouse	WB,IHC-P,IHC-Fr,IF-Tissue
	CD11c Recombinant Rabbit mAb [SI19-06]	ET1606-19	Human	WB,IHC-P,IP,mIHC,IF-Tissue
	CD11c Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS102MS	Human	mIHC
	CD11c Recombinant Antibody [SI19-06] - Rat IgG1 (Chimeric)	HA601418		
	CD11c Recombinant Rabbit mAb	IRS084	Mouse	mIHC
	CD11c / ITGAX Mouse mAb [3G1]	RT1108	Human	WB,IP,IF,IHC-P
	CD127 Rabbit pAb	0804-9	Human	WB,IHC-P,IF-Cell,FC
	CD127 Recombinant Mouse mAb [PSH12-76]	HA601429	Human	FC
	APC Conjugated CD127 Recombinant Mouse mAb [PSH12-76]	HA601536F	Human	FC
KD	CD127 Recombinant Rabbit mAb [JE37-10]	HA721214	Human, Mouse	WB,IF-Cell,IHC-P,FC
	CD16 Mouse mAb [Q32]	EM60067	Human	WB,IHC-P
	CD16 Recombinant Rabbit mAb	IRS017	Human	mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CD16 Recombinant Rabbit mAb [JE49-79]	ET7109-97	Human, Mouse, Rat	WB
Relative expression	CD16 Recombinant Rabbit mAb [PD00-12]	HA721149	Human	IHC-P,IF-Tissue,IF-Cell,mIHC,WB
	CD161 Mouse mAb [A8E1]	HA601063	Human	WB,IHC-P,FC
	CD161c Rabbit pAb	ER1905-79	Mouse, Rat	WB
	CD18 Recombinant Rabbit mAb [JM71-31]	ET1703-94	Human	WB,IP
	CD223 Mouse mAb [4H6A3]	EM1710-97	Human	WB,FC
	CD27 Recombinant Rabbit mAb	IRS115RB	Human	mIHC
Relative expression	CD27 Recombinant Rabbit mAb [PSH08-90]	HA723042	Human	WB,IF-Cell,IHC-P,FC,IP
Relative expression	CD27 Recombinant Rabbit mAb [PSH08-90] - BSA and Azide free	HA751255	Human	WB,IF-Cell,IHC-P,FC,IP
Relative expression	CD27 Recombinant Rabbit mAb [JB40-98]	ET7107-73	Human, Mouse, Rat	WB,IHC-P
	APC Conjugated CD27 Recombinant Rabbit mAb [PSH08-90]	HA723828F	Human	FC
Relative expression	CD30 Recombinant Rabbit mAb [PSH04-12]	HA722093	Human	WB,IF-Cell,IHC-P,IF-Tissue
Relative expression	CD30 Recombinant Rabbit mAb [JE57-71]	HA721579	Human	WB,IHC-P
Relative expression	CD30 Recombinant Rabbit mAb [PSH04-11]	HA722092	Human	WB,IF-Cell,IHC-P
Relative expression	CD30 Recombinant Rabbit mAb [PSH04-10]	HA722091	Human	WB,IHC-P,FC,IF-Tissue
Relative expression, Cell treatment	CD33 Recombinant Rabbit mAb [PSH02-50]	HA721826	Human	WB,IF-Cell,FC
Relative expression, Cell treatment	CD33 Recombinant Rabbit mAb [PSH02-50] - BSA and Azide free	HA750813	Human	WB,IF-Cell,FC
	CD33 Mouse mAb [B7-C6]	EM1709-20	Human	IHC-P,FC
Relative expression	CD33 Recombinant Rabbit mAb [PSH02-51] - BSA and Azide free	HA750814	Human	WB,IF-Cell,FC
	CD33 Rabbit pAb	ER1905-89	Mouse	WB,IHC-P,IF-Tissue
Relative expression	CD33 Recombinant Rabbit mAb [PSH02-51]	HA721827	Human	WB,IF-Cell,FC
	CD34 Recombinant Rabbit mAb	IRS050	Mouse	mIHC
	CD34 Recombinant Rabbit mAb [SI16-01] - BSA and Azide free	HA750092	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	CD34 Recombinant Rabbit mAb	IRS024	Human	mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CD34 Recombinant Mouse mAb [PDM0-12]	HA601109	Human	IHC-P,IF-Tissue
	CD34 Rabbit pAb	0108-2	Human, Mouse, Zebrafish	WB,IF-Cell,IHC-P
	CD34 Rabbit pAb	ER1706-66	Human, Mouse	WB,IHC-P
	CD34 Mouse mAb [15H1]	EM1901-01	Human, Mouse, Rat	IHC-P,FC
	CD34 Recombinant Rabbit mAb [SI16-01]	ET1606-11	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,mIHC,FC
	CD34 Rabbit pAb	ER0802	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC
Relative expression	CD38 Recombinant Rabbit mAb [PD01-49]	HA721268	Human	WB,IHC-P,IF-Cell,FC,mIHC,IF-Tissue
	CD38 Recombinant Rabbit mAb [SN07-20]	ET1611-66	Human	WB,IHC-P,FC,IF-Tissue
Relative expression	CD38 Recombinant Antibody [PD01-49] - Rat IgG1 (Chimeric)	HA601417		
	CD38 Recombinant Rabbit mAb	IRS006	Human	mIHC
	CD38 Recombinant Rabbit mAb [JE29-35]	ET7110-53	Human	WB,IHC-P,FC,IF-Tissue
	CD38 Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS120MS	Human	mIHC
	CD38 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS119RT	Human	mIHC
	CD38 Recombinant Antibody [PD01-49] - Mouse IgG1 (Chimeric)	HA601554	Human	WB,IHC-P,FC
	CD38 Recombinant Antibody [PD01-49] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610369	Human	WB,IHC-P,FC
Relative expression	CD38 Recombinant Antibody [PD01-49] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610257		
	CD39 Recombinant Rabbit mAb [JA90-36] - BSA and Azide free	HA750421	Human, Mouse, Rat	WB,IHC-P
	CD39 Recombinant Rabbit mAb [JA90-36]	ET1704-74	Human, Mouse, Rat	WB,IHC-P
	CD45 Rabbit pAb	ER1905-93	Human	FC
	CD45 Rabbit pAb	ER1902-02	Human	WB,FC
Relative expression	CD45 Recombinant Rabbit Multiclonal Antibody [PSH14-19]	HA723606	Human, Mouse	WB,IHC-P,IF-Cell
Relative expression	CD45 Recombinant Rabbit Multiclonal Antibody [PSH14-19] - BSA and Azide free	HA751510	Human, Mouse	WB,IHC-P,IF-Cell
	CD45 Recombinant Rabbit mAb	IRS025	Human	mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CD45 Recombinant Rabbit mAb [JE03-05]	ET7111-03	Human	WB,IHC-P,FC,IF-Tissue,mIHC
	CD45 Rabbit pAb	ER1901-29	Human	WB,IHC-P,FC
	CD45RO Recombinant Mouse mAb	IRS124MS	Human	mIHC
	CD45RO Recombinant Mouse mAb [PSH13-94]	HA601441	Human	IHC-P,FC
	PE Conjugated CD45RO Recombinant Mouse mAb [PSH13-94]	HA601531F	Human	FC
	CD66b Rabbit pAb	HA500100	Human	WB,IHC-P,mIHC,IF-Tissue
	CD66b pAb	IRS127RB	Human	mIHC
	CD69 Rabbit pAb	HA500073	Human, Mouse, Rat	WB
	CD69 Mouse mAb [E6-B6]	EM1709-21	Human	WB,IHC-P,FC
Cell treatment	CD69 Mouse mAb [A8C4]	HA601013	Human	WB
Relative expression	CD7 Recombinant Rabbit mAb [JE57-72]	HA720023	Human	WB,IF-Cell,IHC-P
Cell treatment	CD82 Recombinant Rabbit mAb [JE32-17]	HA722533	Human	WB
	c-kit Rabbit pAb	0108-4	Human	WB,IF-Cell
Relative expression, KD	c-Kit Recombinant Rabbit mAb [ST04-99]	ET1609-60	Human	WB,IHC-P
Relative expression	c-Kit Recombinant Rabbit mAb [PSH17-96]	HA723961	Mouse, Rat	WB,IHC-Fr,IHC-P
Relative expression, KD	c-Kit Recombinant Rabbit mAb [PD00-24]	HA721154	Human	IHC-P,WB,mIHC
	c-Kit Recombinant Rabbit mAb	IRS013	Human	mIHC
Relative expression	c-Kit Recombinant Rabbit mAb [PSH17-79]	HA723948	Mouse	IHC-Fr,IF-Cell,FC
	CX3CR1 Mouse mAb [A5C7]	HA600074	Human	IF-Cell
	Histamine H2 receptor Rabbit pAb	HA500154	Human	WB,IHC-P
	FITC Conjugated Human CD10 Recombinant Rabbit mAb [SN75-07]	HA720209F	Human	FC
	FITC Conjugated Human CD11b Recombinant Mouse mAb [A2-9-R]	HA600130F	Human	FC
	Biotin Conjugated Human CD127 Recombinant Rabbit mAb [PSH18-28] - Detector	HA725346B	Human	ELISA(Det),ELISA
	Human CD127 Recombinant Rabbit mAb [PSH18-27] - BSA and Azide free (Capture)	HA725345	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human CD127 Recombinant Rabbit mAb [PSH18-28] - BSA and Azide free (Detector)	HA725346	Human	ELISA(Det)
	PE Conjugated Human CD16 Recombinant Mouse mAb [PSH04-95]	HA600129F	Human	FC
	Human CD16 Recombinant Mouse mAb [PSH04-95]	HA601340	Human	FC
	Human CD27/TNFRSF7 Recombinant Rabbit mAb [PSH08-33] - BSA and Azide free (Detector)	HA722975	Human	ELISA(Det)
	Human CD27/TNFRSF7 Recombinant Rabbit mAb [PSH08-31] - BSA and Azide free (Capture)	HA722973	Human	ELISA(Cap)
	Biotin Conjugated Human CD27/TNFRSF7 Recombinant Rabbit mAb [PSH08-33] - Detector	HA722976B	Human	ELISA(Det),ELISA
	Human CD27/TNFRSF7 Recombinant Rabbit mAb [PSH08-32] - BSA and Azide free (Capture)	HA722974	Human	ELISA(Cap)
	Human CD30 Recombinant Rabbit mAb [PSH04-82] - BSA and Azide free (Detector)	HA722188	Human	ELISA(Det)
	Human CD30 Recombinant Rabbit mAb [PSH04-81] - BSA and Azide free (Capture)	HA722187	Human	ELISA(Cap)
	Biotin Conjugated Human CD30 Recombinant Rabbit mAb [PSH04-82]	HA722601B	Human	ELISA(Det)
	Biotin Conjugated Human CD33 Recombinant Rabbit mAb [PSH05-26]	HA722599B	Human	ELISA(Det)
	Human CD33 Recombinant Rabbit mAb [PSH05-25] - BSA and Azide free (Capture)	HA722254	Human	ELISA(Cap)
	Human CD33 Recombinant Rabbit mAb [PSH05-26] - BSA and Azide free (Detector)	HA722255	Human	ELISA(Det)
	FITC Conjugated Human CD33 Recombinant Rabbit mAb [PSH04-96]	HA720197F	Human	FC
	FITC Conjugated Human CD38 Recombinant Rabbit mAb [SN07-20]	HA720196F	Human	FC
	iFluor™ 555 Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA601556F	Human	FC
	iFluor™ 488 Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600109F	Human	FC
	Percp-Cy5.5 Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600128F	Human	FC
	Human CD45 Recombinant Mouse mAb [PSH03-64]	HA601285	Human	IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	APC Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600127F	Human	FC
	FITC Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600112F	Human	FC
	FITC Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600125F	Human	FC
	PE Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600126F	Human	FC
	Human CD45 Recombinant Mouse mAb [PSH03-64] - BSA and Azide free	HA610158	Human	IF-Cell,FC
	PE Conjugated Human CD45RA Recombinant Mouse mAb [PSH04-97]	HA601525F	Human	FC
	FITC Conjugated Human CD45RA Recombinant Mouse mAb [PSH04-97]	HA600131F	Human	FC
	Human CD45RA Recombinant Mouse mAb [PSH04-97]	HA601341	Human	FC,IF-Cell,IHC-P
	Biotin Conjugated Human IL-12B p40 Recombinant Rabbit mAb [PSH03-04]	HA722629B	Human	ELISA(Det)
	Human IL-12B p40 Recombinant Rabbit mAb [PSH02-58] - BSA and Azide free (Detector)	HA721835	Human	ELISA(Det)
	Human IL-12B p40 Recombinant Rabbit mAb [PSH03-04] - BSA and Azide free (Detector)	HA721934	Human	ELISA(Det)
	Human IL-12B p40 Recombinant Rabbit mAb [PSH02-57] - BSA and Azide free (Capture)	HA721834	Human	ELISA(Cap)
	Biotin Conjugated Human IL-12B p40 Recombinant Rabbit mAb [PSH02-58]	HA722633B	Human	ELISA(Det)
	Human IL-18 Recombinant Rabbit mAb [PSH01-71] - BSA and Azide free (Capture)	HA721722	Human	ELISA(Cap)
Relative expression	Human IL-18 Recombinant Rabbit mAb [PSH02-09]	HA721768	Human	WB
	Human IL-18 Recombinant Rabbit mAb [PSH02-08]	HA721767	Human	WB,IHC-P,IF-Tissue,IP
	Human IL-18 Recombinant Rabbit mAb [PSH01-70] - BSA and Azide free (Detector)	HA721723	Human	ELISA(Det)
Relative expression	Human IL-18 Recombinant Rabbit mAb [PSH02-09] - BSA and Azide free	HA750779	Human	WB
	Human IL-18 Recombinant Rabbit mAb [PSH02-08] - BSA and Azide free	HA750778	Human	WB,IHC-P,IF-Tissue
	Biotin Conjugated Human IL-18 Recombinant Rabbit mAb [PSH01-70]	HA722644B	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human IL-2 Mouse mAb [A8B7] - BSA and Azide free (Detector)	HA601006	Human	WB,ELISA(Det)
	Biotin Conjugated Human IL-2 Mouse mAb [A8B7]	HA601361B	Human	ELISA(Det)
	Human IL-2 Recombinant Rabbit mAb [PSH01-42] - BSA and Azide free (Capture)	HA721680	Human	ELISA(Cap)
	Human IL-7 Recombinant Rabbit mAb [PSH03-09] - BSA and Azide free (Capture)	HA721937	Human	ELISA(Cap)
	Human IL-7 Recombinant Rabbit mAb [PSH03-10] - BSA and Azide free (Detector)	HA721938	Human	ELISA(Det)
	Biotin Conjugated Human IL-7 Recombinant Rabbit mAb [PSH03-10]	HA722632B	Human	ELISA(Det)
	Human Interferon gamma Recombinant Rabbit mAb [PSH01-46] - BSA and Azide free (Capture)	HA721684	Human	ELISA(Cap)
	Human Interferon gamma Recombinant Rabbit mAb [PSH01-47] - BSA and Azide free (Detector)	HA721685	Human	ELISA(Det)
	Human Interferon gamma Recombinant Rabbit mAb [PSH01-45] - BSA and Azide free (Capture)	HA721683	Human	ELISA(Cap)
	Human Interferon gamma Recombinant Rabbit mAb [PS01-76] - BSA and Azide free (Detector)	HA721575	Human	ELISA(Det)
	Biotin Conjugated Human Interferon gamma Recombinant Rabbit mAb [PS01-76] - Detector	HA721576	Human	ELISA(Det),ELISA
	Human Interferon gamma Recombinant Rabbit mAb [PS01-75] - BSA and Azide free (Capture)	HA721574	Human	ELISA(Cap)
	Biotin Conjugated Human Interferon gamma Recombinant Rabbit mAb [PSH01-47]	HA722645B	Human	ELISA(Det)
	Human LAG-3 Recombinant Rabbit mAb [PSH05-29] - BSA and Azide free (Capture)	HA722258	Human	ELISA(Cap)
	Biotin Conjugated Human LAG-3 Recombinant Rabbit mAb [PSH05-30]	HA722596B	Human	ELISA(Det)
	Human LAG-3 Recombinant Rabbit mAb [PSH05-30] - BSA and Azide free (Detector)	HA722259	Human	ELISA(Det)
	Human MICA Recombinant Rabbit mAb [PSH11-39] - BSA and Azide free (Capture)	HA723327	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human MICA Recombinant Rabbit mAb [PSH11-40] - BSA and Azide free (Detector)	HA723328	Human	ELISA(Det)
	Biotin Conjugated Human MICA Recombinant Rabbit mAb [PSH11-40] - Detector	HA723329B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human MICB Recombinant Rabbit mAb [PSH10-74]	HA723240B	Human	ELISA(Det),ELISA
	Human MICB Recombinant Rabbit mAb [PSH10-73] - BSA and Azide free (Capture)	HA723238	Human	ELISA(Cap)
	Human MICB Recombinant Rabbit mAb [PSH10-74] - BSA and Azide free (Detector)	HA723239	Human	ELISA(Det)
	Human NCAM1 Recombinant Rabbit mAb [PSH08-08] - BSA and Azide free (Capture/Detector)	HA722947	Human	ELISA(Det),ELISA(Cap)
	Human NCAM1 Recombinant Rabbit mAb [PSH08-07] - BSA and Azide free (Capture)	HA722946	Human	ELISA(Cap)
	Human NCAM1 Recombinant Rabbit mAb [PSH08-09] - BSA and Azide free (Detector)	HA722949	Human	ELISA(Det)
	Biotin Conjugated Human NCAM1 Recombinant Rabbit mAb [PSH08-09] - Detector	HA722950B	Human	ELISA(Det)
	Biotin Conjugated Human NCAM1 Recombinant Rabbit mAb [PSH08-08] - Detector	HA722948B	Human	ELISA(Det)
	APC Conjugated Human NCAM1 / CD56 Recombinant Rabbit mAb [PSH06-53]	HA723757F	Human	FC
	PE Conjugated Human NCAM1 / CD56 Recombinant Rabbit mAb [PSH06-53]	HA720215F	Human	FC
	Human TNF alpha Recombinant Rabbit mAb	IRS040	Human	mIHC
	Human TNF alpha Recombinant Rabbit mAb [PS00-40] - BSA and Azide free	HA750593	Human	IHC-P
	Human TNF alpha Recombinant Rabbit mAb [PS00-40]	HA721286	Human	IHC-P,mIHC
	HRP Conjugated Human TNF alpha Recombinant Rabbit mAb [PS01-07] - Detector	HA721294	Human	ELISA(Det)
	Biotin Conjugated Human TNF alpha Recombinant Rabbit mAb [PS01-07] - Detector	HA721559	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human TNF alpha Recombinant Rabbit mAb [PS01-07] - BSA and Azide free (Detector)	HA721266	Human	ELISA(Det)
	Human TNF alpha Recombinant Rabbit mAb [PS00-40] - BSA and Azide free (Capture)	HA721265	Human	ELISA(Cap)
	Biotin Conjugated Human ULBP1 Recombinant Rabbit mAb [PSH11-96] - Detector	HA723408B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human ULBP1 Recombinant Rabbit mAb [PSH11-97] - Detector	HA723410B	Human	ELISA(Det),ELISA
	Human ULBP1 Recombinant Rabbit mAb [PSH11-97] - BSA and Azide free (Detector)	HA723409	Human	ELISA(Det)
	Human ULBP1 Recombinant Rabbit mAb [PSH11-96] - BSA and Azide free (Detector)	HA723407	Human	ELISA(Det)
	Human ULBP1 Recombinant Rabbit mAb [PSH11-95] - BSA and Azide free (Capture)	HA723406	Human	ELISA(Cap)
	IL12B Rabbit pAb	ER1803-31	Human, Mouse	WB
	IL-18 Recombinant Rabbit mAb [PS00-41]	HA721536	Human	WB,IHC-P
Relative expression	IL-18 Recombinant Rabbit mAb [PSH18-15]	HA723978	Human, Mouse, Rat	WB,IHC-P,IP
	IL-18 Recombinant Rabbit mAb [PS01-06] - BSA and Azide free	HA750677	Human	WB,FC
	IL-18 Recombinant Rabbit mAb [PS00-41] - BSA and Azide free	HA750676	Human	WB,IHC-P
	IL-18 Recombinant Rabbit mAb [PS01-06]	HA721537	Human	WB,FC
	IL-2 Mouse mAb [A8B9]	HA601016	Human	WB,ELISA
	IL-2 Rabbit pAb	ER63748	Human, Rat	WB,IHC-P
	IL-2 Mouse mAb [A8B8]	HA601015	Human	WB,ELISA
	IL-7 Rabbit pAb	HA500171	Human	WB,IF-Cell,IHC-P
	IL-7 Mouse mAb [A8B11]	HA601048	Human	WB
	IL7 Mouse mAb [AH31-34]	EM170412	Human, Rat	ELISA,WB,ICC,IHC-P
	IL-7R alpha Rabbit pAb	ER1706-68	Human, Zebrafish	WB,IHC-P,FC
	Interferon gamma Recombinant Rabbit mAb [JM10-10]	ET1703-17	Human	WB,FC
	Interferon gamma Recombinant Rabbit mAb [JM10-10] - BSA and Azide free	HA750373	Human	WB,FC
	LAG-3 Recombinant Rabbit mAb	IRS020	Human	mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	LAG-3 Recombinant Rabbit mAb [PSH0-42] - BSA and Azide free	HA750620	Mouse	WB,IF-Cell
Cell treatment	LAG-3 Recombinant Rabbit mAb [PSH0-43]	HA721345	Human, Mouse	WB,IHC-P
Cell treatment	LAG-3 Recombinant Rabbit mAb [PS01-33]	HA721358	Human	IHC-P,IF-Cell,WB,mIHC
Cell treatment	LAG-3 Recombinant Rabbit mAb [PSH0-43] - BSA and Azide free	HA750621	Human, Mouse	WB,IHC-P
Cell treatment	LAG-3 Recombinant Rabbit mAb [PSH0-44] - BSA and Azide free	HA750622	Human, Mouse	WB,IHC-P,IF-Cell
Cell treatment	LAG-3 Recombinant Rabbit mAb [PSH0-44]	HA721346	Human, Mouse	WB,IHC-P,IF-Cell
Cell treatment	LAG-3 Recombinant Rabbit mAb [PSH0-42]	HA721344	Mouse	WB,IF-Cell
Relative expression	LAIR1 Recombinant Mouse mAb [PSH04-50] - BSA and Azide free	HA610201	Human	WB,IHC-P
Relative expression	LAIR1 Recombinant Mouse mAb [PSH04-50]	HA601329	Human	WB,IHC-P
Relative expression	MICA Recombinant Rabbit mAb [PSH13-60] - BSA and Azide free	HA751479	Human	WB,FC,IP
Relative expression	MICA Recombinant Rabbit mAb [PSH13-60]	HA723549	Human	WB,FC,IP
	MICA Mouse mAb [A4C5]	HA600002	Human	WB
Relative expression	MICA + MICB Recombinant Rabbit mAb [PSH10-88] - BSA and Azide free	HA751380	Human	WB,FC,IP
Relative expression	MICA + MICB Recombinant Rabbit mAb [PSH10-88]	HA723255	Human	WB,FC,IP
Relative expression	MICA + MICB Rabbit pAb	ER1803-95	Human	WB
	PE Conjugated Mouse CD16/32 Recombinant Rabbit mAb [PSH06-60]	HA720216F	Mouse	FC
	Mouse IL12 p40 Recombinant Rabbit mAb [PSH04-88] - BSA and Azide free (Detector)	HA722194	Mouse	ELISA(Det)
	Mouse IL12 p40 Recombinant Rabbit mAb [PSH04-87] - BSA and Azide free (Capture)	HA722193	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL12 p40 Recombinant Rabbit mAb [PSH04-88]	HA722603B	Mouse	ELISA(Det)
	Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-29] - BSA and Azide free (Capture)	HA721974	Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-30]	HA722626B	Mouse, Rat	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-30] - BSA and Azide free (Detector)	HA721975	Mouse, Rat	ELISA(Det)
	Biotin Conjugated Mouse TNF alpha Recombinant Rabbit mAb [PSH03-77]	HA722623B	Mouse	ELISA(Det)
	Mouse TNF alpha Recombinant Rabbit mAb [PSH03-77] - BSA and Azide free (Detector)	HA722041	Mouse	ELISA(Det)
	Mouse TNF alpha Recombinant Rabbit mAb [PSH03-76] - BSA and Azide free (Capture)	HA722040	Mouse	ELISA(Cap)
Relative expression	NCAM1 / CD56 Recombinant Rabbit mAb [PSH06-53]	HA722756	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
Relative expression	NCAM1 / CD56 Recombinant Rabbit mAb [PSH08-89]	HA723041	Human	WB,IHC-P,IF-Tissue
Relative expression	NCAM1 / CD56 Recombinant Rabbit mAb [PSH06-82]	HA722755	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Cell,FC,IF-Tissue
Relative expression	NCAM1 / CD56 Recombinant Rabbit mAb [PSH06-53] - BSA and Azide free	HA751105	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
Relative expression	NCAM1 / CD56 Recombinant Mouse mAb [PDH0-07]	HA601126	Human	IHC-P,WB
	NCAM1 / CD56 Rabbit pAb	R1204-1	Human, Mouse, Rat	WB,IF-Cell
Relative expression	NCAM1 / CD56 Recombinant Rabbit mAb [PSH06-82] - BSA and Azide free	HA751104	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Cell,FC,IF-Tissue
	NCAM1 / CD56 Recombinant Rabbit mAb [JF1021] - BSA and Azide free	HA750346	Human, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	NCAM1 / CD56 Recombinant Rabbit mAb [JF1021]	ET1702-43	Human, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP,mIHC
	NCAM1 / CD56 Recombinant Rabbit mAb	IRS026	Human	mIHC
	Neutrophil Elastase Recombinant Rabbit mAb [JF098-6]	ET1702-78	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P
	Neutrophil Elastase Recombinant Rabbit mAb [JF098-6] - BSA and Azide free	HA750360	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P
	NK1.1 Recombinant Rabbit mAb [PSH08-47] - BSA and Azide free	HA751222	Mouse	WB,IHC-P
	NK1.1 Recombinant Rabbit mAb [PSH08-46] - BSA and Azide free	HA751221	Mouse	FC
	FITC Conjugated NK1.1 Recombinant Rabbit mAb [PSH08-46]	HA720219F	Mouse	FC
	NK1.1 Recombinant Rabbit mAb [PSH08-46]	HA722992	Mouse	FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	NK1.1 Recombinant Rabbit mAb [PSH08-47]	HA722993	Mouse	WB,IHC-P
	NK-p44 Recombinant Rabbit mAb [JE61-04]	HA720070	Human	WB
	PTPRC Mouse mAb [C11-C10]	EM1708-69	Human	WB,IHC-P
	TNF alpha Rabbit pAb	ER65189	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Cell treatment	TNF alpha Rabbit pAb	R1203-1	Mouse	WB,IF-Cell,ELISA
	TNF alpha Rabbit pAb	ER1919-22	Human, Mouse	WB
Cell treatment	TNF alpha Recombinant Rabbit mAb [PS01-07] - BSA and Azide free	HA750897	Human, Mouse	WB
Cell treatment	TNF alpha Recombinant Rabbit mAb [PS01-07]	HA722022	Human, Mouse	WB
	ULBP1 Recombinant Rabbit mAb [JE65-75]	HA720092	Human	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层