

T细胞研究-产品清单

高级验证	产品名	货号	反应种属	验证应用
Cell treatment	4-1BBL Recombinant Rabbit mAb [JE33-91]	HA721357	Human, Mouse	WB
	ALCAM Mouse mAb [C7-E12]	EM1708-92	Human	WB,IF-Cell,IHC-P,FC
	ALCAM/CD166 Mouse mAb [7C2]	RT1029	Human, Mouse, Rat	WB,IP,IF,IHC-P
Relative expression	B7H4 Recombinant Rabbit mAb [PSH03-65] - BSA and Azide free	HA750899	Human, Mouse, Rat	WB,IF-Cell
Relative expression	B7H4 Recombinant Rabbit mAb [PSH03-65]	HA722024	Human, Mouse, Rat	WB,IF-Cell
	B7-H4(VTCN1) Rabbit pAb	R1512-1	Human, Mouse	WB,IF-Cell,IHC-P,FC
	B7-H4(VTCN1) Rabbit pAb	ER1802-56	Human, Rat	WB,IF-Cell,IHC-P,FC
Relative expression	CD134 Recombinant Rabbit mAb [PSH10-14] - BSA and Azide free	HA751343	Human	WB,IF-Cell,FC,IP
	CD134 Mouse mAb [2H1]	EM1701-24	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Relative expression	CD134 Recombinant Rabbit mAb [PSH10-14]	HA723183	Human	WB,IF-Cell,FC,IP
Relative expression	CD134 Recombinant Rabbit mAb [JE36-48]	HA722700	Human	WB
KD	CD166 Mouse mAb [5-5]	EM1902-36	Human	WB,IHC-P,IF-Cell,FC
KD	CD166 Mouse mAb [A6A9]	HA600082	Human	IF-Cell,IHC-P,WB
KD	CD166 Mouse mAb [A6A8]	HA600083	Human	IF-Cell,IHC-P,WB
	CD166 Recombinant Rabbit mAb [JE57-47]	HA722705	Human, Mouse, Rat	WB,IHC-P,FC,IHC-Fr
KD	CD166 Mouse mAb [11-6]	M1012-7	Human, Mouse, Rat	WB,IHC-P,FC
	CD1a Recombinant Rabbit mAb [JE28-74]	ET7110-48	Human	WB,IHC-P
Relative expression	CD1a Mouse mAb [A4B10]	HA600031	Human	IHC-P,FC
Relative expression	CD1a Recombinant Rabbit mAb [JM21-33]	ET1705-17	Human	WB,IHC-P,IF-Tissue
	CD1a Rabbit pAb	ER1902-73	Human	WB,FC
	CD1a Mouse mAb [A3F9]	EM1901-89	Human	WB,IHC-P
	CD1c Recombinant Rabbit mAb [JE55-00]	HA722562	Human	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CD1c Recombinant mAb	IRS111RB	Human	mIHC
	CD1C Mouse mAb [3G1B3]	EM1710-74	Human	WB,IHC-P,FC
Relative expression	CD2 Recombinant Rabbit mAb [PSH05-85] - BSA and Azide free	HA751033	Human	WB,IHC-P,IF-Cell,FC,IF-Tissue
Relative expression	CD2 Recombinant Rabbit mAb [PSH05-84]	HA722483	Human	WB,IHC-P
Relative expression	CD2 Recombinant Rabbit mAb [PSH05-85]	HA722484	Human	WB,IHC-P,IF-Cell,FC,IF-Tissue
	CD2 Mouse mAb [B4-D1]	EM1707-56	Human, Mouse	WB,IHC-P,FC
Relative expression	CD200 Recombinant Rabbit mAb [PSH01-48]	HA721691	Human	WB,IHC-P,IF-Tissue,FC
Relative expression	CD200 Recombinant Rabbit mAb [PSH01-49]	HA721692	Human	WB,IHC-P,FC,IF-Tissue
	CD223 Mouse mAb [4H6A3]	EM1710-97	Human	WB,FC
	APC Conjugated CD25 / IL-2 Receptor alpha Recombinant Mouse mAb [PSH14-46]	HA601530F	Human	FC
	CD25 / IL-2 Receptor alpha Recombinant Mouse mAb [PSH14-46]	HA601486	Human, Mouse	FC,IHC-P,IF-Tissue
	CD27 Recombinant Rabbit mAb	IRS115RB	Human	mIHC
Relative expression	CD27 Recombinant Rabbit mAb [PSH08-90]	HA723042	Human	WB,IF-Cell,IHC-P,FC,IP
Relative expression	CD27 Recombinant Rabbit mAb [JB40-98]	ET7107-73	Human, Mouse, Rat	WB,IHC-P
Relative expression	CD27 Recombinant Rabbit mAb [PSH08-90] - BSA and Azide free	HA751255	Human	WB,IF-Cell,IHC-P,FC,IP
	APC Conjugated CD27 Recombinant Rabbit mAb [PSH08-90]	HA723828F	Human	FC
	CD276 Mouse mAb [C10-A1]	EM1709-85	Human	WB,IHC-P,FC
Relative expression, KD	CD276 Recombinant Rabbit mAb [PSH02-46]	HA721823	Human	WB,IF-Cell,FC
Relative expression, KD	CD276 Recombinant Rabbit mAb [PD00-36]	HA721245	Human	IHC-P,FC,WB
KD	CD276 Rabbit pAb	ER1803-36	Human	WB,IF-Cell
Relative expression, KD	CD276 Recombinant Rabbit mAb [PSH02-46] - BSA and Azide free	HA750811	Human	WB,IF-Cell,FC
	CD28 Rabbit pAb	ER2001-42	Human, Mouse, Rat	WB,IHC-P,FC
	CD3 Recombinant Mouse mAb [PD01-43]	HA601110	Human	WB,IHC-P,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	CD3 Recombinant Rabbit mAb [PSH19-70]	HA724094	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	CD3 Recombinant Rabbit mAb [JE80-02]	HA720082	Human	WB,IHC-P,IF-Tissue,IP,mIHC,IF-Cell
	CD3 Recombinant Rabbit mAb	IRS022	Human	mIHC
	CD3 Recombinant Rabbit mAb	IRS074	Mouse	mIHC
	CD38 Recombinant Antibody [PD01-49] - Mouse IgG1 (Chimeric)	HA601554	Human	WB,IHC-P,FC
Relative expression	CD38 Recombinant Rabbit mAb [PD01-49]	HA721268	Human	WB,IHC-P,IF-Cell,FC,mIHC,IF-Tissue
Relative expression	CD38 Recombinant Antibody [PD01-49] - Rat IgG1 (Chimeric)	HA601417		
	CD38 Recombinant Antibody [PD01-49] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610369	Human	WB,IHC-P,FC
Relative expression	CD38 Recombinant Antibody [PD01-49] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610257		
	CD38 Recombinant Rabbit mAb [SN07-20]	ET1611-66	Human	WB,IHC-P,FC,IF-Tissue
	CD38 Recombinant Rabbit mAb [JE29-35]	ET7110-53	Human	WB,IHC-P,FC,IF-Tissue
	CD38 Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS120MS	Human	mIHC
	CD38 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS119RT	Human	mIHC
	CD38 Recombinant Rabbit mAb	IRS006	Human	mIHC
	CD39 Recombinant Rabbit mAb [JA90-36] - BSA and Azide free	HA750421	Human, Mouse, Rat	WB,IHC-P
	CD39 Recombinant Rabbit mAb [JA90-36]	ET1704-74	Human, Mouse, Rat	WB,IHC-P
	CD3D Recombinant Rabbit mAb [JJ08-97]	ET1701-57	Human	WB,IF-Cell,IHC-P,FC,IP,IF-Tissue
	CD3E Rabbit pAb	HA500344	Human, Mouse, Rat	WB
	CD3 epsilon Recombinant Rabbit mAb [SY0239]	ET1607-29	Human	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP
	CD3 epsilon Recombinant Rabbit mAb [SY0239] - BSA and Azide free	HA750115	Human	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP
	CD3 epsilon Recombinant Rabbit mAb [JE17-74]	ET7110-63	Human	WB,IHC-P,FC
	CD3 epsilon Rabbit pAb	ER80501	Human, Mouse, Rat	WB
	CD3G Recombinant Rabbit mAb [JB38-29]	ET7107-55	Human, Mouse	WB,IHC-P,FC,IP,IF-Cell,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	CD3 zeta Recombinant Rabbit mAb [JE09-37]	HA722692	Human, Rat	WB,IHC-P,IF-Tissue
	CD3 zeta Rabbit pAb	0805-7	Human	WB,IF-Cell,FC
	CD4 Recombinant Rabbit Multiclonal Antibody [PSH13-24] - BSA and Azide free	HA751460	Human, Mouse	WB,IF-Cell,IHC-P
	CD4 Rabbit pAb	R1211-1	Human	WB,IF-Cell
	CD4 Recombinant Rabbit mAb	IRS051	Mouse	mIHC
	CD4 Recombinant Antibody [PSH08-22] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610354	Mouse	WB,IF-Tissue,IHC-P,IHC-Fr
Relative expression	CD4 Recombinant Antibody [PSH08-22] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610353	Mouse	WB,IF-Tissue,IHC-P,IF-Cell
	CD4 Rabbit pAb	ER1803-84	Human	WB,FC
	CD4 Recombinant Rabbit Multiclonal Antibody [PSH13-24]	HA723511	Human, Mouse	WB,IF-Cell,IHC-P
	CD4 Recombinant Rabbit mAb	IRS005	Human	mIHC
	CD4 Recombinant Antibody [PSH08-22] - Rat IgG1 (Chimeric)	HA601534	Mouse	WB,IF-Tissue,IHC-P,IHC-Fr
Relative expression	CD4 Recombinant Antibody [PSH08-22] - Mouse IgG1 (Chimeric)	HA601533	Mouse	WB,IF-Tissue,IHC-P,IF-Cell
	CD4 Recombinant Antibody [ST0488] - Mouse IgG1 (Chimeric)	HA601557	Human	WB,IF-Cell,IHC-P
KD	CD4 Recombinant Rabbit mAb [ST0488]	ET1609-52	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC,mIHC
Relative expression	CD4 Recombinant Rabbit mAb [PSH08-22]	HA722966	Mouse	WB,IF-Cell,IHC-P,mIHC,IHC-Fr
Relative expression	CD4 Recombinant Rabbit mAb [PSH08-22] - BSA and Azide free	HA751213	Mouse	WB,IF-Cell,IHC-P,IHC-Fr
KD	CD4 Rabbit pAb	0108-6	Human, Mouse, Rat	WB,IF-Cell,FC
	CD4 Recombinant Rabbit mAb [JE56-36]	ET7110-64	Human	WB,IHC-P
	CD40L Recombinant Rabbit mAb [JM11-34]	ET1703-64	Human, Mouse, Rat	WB
	CD40L Recombinant Rabbit mAb [JM11-34] - BSA and Azide free	HA750387	Human, Mouse, Rat	WB
	CD45 Rabbit pAb	ER1905-93	Human	FC
Relative expression	CD45 Recombinant Rabbit Multiclonal Antibody [PSH14-19]	HA723606	Human, Mouse	WB,IHC-P,IF-Cell
	CD45 Recombinant Rabbit mAb	IRS025	Human	mIHC
	CD45 Rabbit pAb	ER1901-29	Human	WB,IHC-P,FC
	CD45 Rabbit pAb	ER1902-02	Human	WB,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	CD45 Recombinant Rabbit Multiclonal Antibody [PSH14-19] - BSA and Azide free	HA751510	Human, Mouse	WB,IHC-P,IF-Cell
	CD45 Recombinant Rabbit mAb [JE03-05]	ET7111-03	Human	WB,IHC-P,FC,IF-Tissue,mIHC
	CD45RO Recombinant Mouse mAb [PSH13-94]	HA601441	Human	IHC-P,FC
	PE Conjugated CD45RO Recombinant Mouse mAb [PSH13-94]	HA601531F	Human	FC
	CD45RO Recombinant Mouse mAb	IRS124MS	Human	mIHC
	CD5 Recombinant Rabbit mAb [JE25-41]	ET7110-56	Human	WB,IHC-P
	CD5 Recombinant Rabbit mAb [PD00-01]	HA721137	Human	IHC-P,IF-Tissue,WB
	CD62P Recombinant Rabbit mAb [JM10-95]	ET1703-49	Human	WB,IHC-P
	CD69 Rabbit pAb	HA500073	Human, Mouse, Rat	WB
	CD69 Mouse mAb [E6-B6]	EM1709-21	Human	WB,IHC-P,FC
Cell treatment	CD69 Mouse mAb [A8C4]	HA601013	Human	WB
Relative expression	CD7 Recombinant Rabbit mAb [JE57-72]	HA720023	Human	WB,IF-Cell,IHC-P
	CD8 Rabbit pAb	IRS083	Mouse	mIHC
	CD8 Rabbit pAb	0108-7	Mouse	WB,IF-Cell,FC
Cell treatment	CD82 Recombinant Rabbit mAb [JE32-17]	HA722533	Human	WB
	CD83 Recombinant Rabbit mAb [JE33-61]	HA721556	Human	WB
	CD83 Rabbit pAb	ER62949	Human	WB,IHC-P
	CD86 Recombinant Rabbit mAb [SJ20-00] - BSA and Azide free	HA750109	Human, Mouse, Rat	WB,IF-Cell,IP
	CD86 Rabbit pAb	ER1906-01	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
	CD86 Recombinant Rabbit mAb [SJ20-00]	ET1606-50	Human, Mouse, Rat	WB,IF-Cell,IP
	CD8a Mouse mAb	IRS007	Human	mIHC
	CD8 alpha Recombinant Rabbit mAb [SI18-01] - BSA and Azide free	HA750101	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
	CD8 alpha Recombinant Rabbit mAb [SI18-01]	ET1606-31	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
Relative expression	CD8 alpha Recombinant Rabbit mAb [PSH19-81]	HA724105	Human, Rat	WB,IHC-P,IF-Tissue,FC
	CD90 / THY1 Recombinant Rabbit mAb [SU35-07] - BSA and Azide free	HA750149	Human, Rat	WB,IHC-P
	CD90 / THY1 Rabbit pAb	R1107-9	Human, Mouse, Rat	WB
KD	CD90 / THY1 Recombinant Rabbit mAb [JF10-09]	ET1702-92	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	CD90 / THY1 Recombinant Rabbit mAb [JF10-09] - BSA and Azide free	HA750364	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
	CD90 / THY1 Recombinant Rabbit mAb [SU35-07]	ET1608-46	Human, Rat	WB,IHC-P
	CD99 Rabbit pAb	ER1803-81	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	CD99 Recombinant Rabbit mAb [PD00-13]	HA721150	Human	IHC-P,IF-Tissue
	CD99 Recombinant Rabbit mAb [JF0991]	ET1702-35	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	ChT1 Rabbit pAb	HA500237	Human, Mouse, Rat	WB,IHC-P
	ChT1 Rabbit pAb	HA500314	Human, Mouse	WB,IF-Cell,IHC-P,FC
	COPS8 / CSN8 Rabbit pAb	ER2001-29	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	CTLA-4 Recombinant Rabbit mAb [PS01-34] - BSA and Azide free	HA750581	Human, Mouse	IHC-P
	CTLA-4 Recombinant Rabbit mAb	IRS016	Human	mIHC
	CTLA-4 Recombinant Rabbit mAb [PS01-34]	HA721269	Human, Mouse	IHC-P,mIHC
	CTLA4 Rabbit pAb	ER1906-81	Human, Mouse, Rat	WB,FC
	CXCR4 Rabbit pAb	ER1802-28	Human, Rat	WB,IF-Cell,FC
	CXCR4 Recombinant Rabbit mAb [JE03-35]	HA722304	Human	WB,IF-Cell,IHC-P
	DC-SIGN Rabbit pAb	ER1907-23	Human	WB,IHC-P,IF-Tissue,FC
	DC-SIGN Recombinant Rabbit mAb [JE31-24]	HA721382	Human	WB
KO	DPP4 / CD26 Recombinant Rabbit mAb [JM11-42] - BSA and Azide free	HA750396	Human, Mouse, Rat	WB
Cell treatment	DPP4 / CD26 Recombinant Rabbit mAb [PSH12-59] - BSA and Azide free	HA751446	Human	WB,IHC-P,IF-Cell,FC,IP
Cell treatment	DPP4 / CD26 Recombinant Rabbit mAb [PSH12-59]	HA723473	Human	WB,IHC-P,IF-Cell,FC,IP
KO	DPP4 / CD26 Recombinant Rabbit mAb [JM11-42]	ET1703-93	Human, Mouse, Rat	WB
Relative expression	FOXP3 Recombinant Rabbit mAb [PSH07-47] - BSA and Azide free	HA751168	Human, Mouse, Rat	WB,IHC-P
Relative expression	FOXP3 Recombinant Rabbit mAb [PSH07-47]	HA722835	Human, Mouse, Rat	WB,IHC-P,mIHC
	FOXP3 Recombinant Rabbit mAb [JF0898] - BSA and Azide free	HA750335	Human, Mouse, Rat	WB,IF-Cell,FC,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	FOXP3 Recombinant Rabbit mAb	IRS053	Mouse	mIHC
	FOXP3 Rabbit pAb	IRS147RB	Mouse	mIHC
	FOXP3 Recombinant Rabbit mAb [JF0926]	ET1702-16	Human	WB,FC,IHC-P
	FOXP3 Recombinant Rabbit mAb	IRS008	Human	mIHC
	FOXP3 Recombinant Rabbit mAb [JF0898]	ET1702-12	Human, Mouse, Rat	WB,IF-Cell,FC,IHC-P
Relative expression	GATA2 Recombinant Rabbit mAb [PSH04-53]	HA722152	Human, Rat	WB,IF-Cell
Relative expression	GATA2 Recombinant Rabbit mAb [PSH04-53] - BSA and Azide free	HA750960	Human, Rat	WB,IF-Cell
	GATA-2/3 Rabbit pAb	ER63507	Human, Mouse	WB,IHC-P,IF-Cell
Relative expression	GATA2 + GATA3 Recombinant Rabbit mAb [JE04-46]	HA722713	Human, Mouse, Rat	WB,IHC-P
	GATA3 Rabbit pAb	ER1901-20	Human, Rat	WB,IF-Cell,IHC-P,FC
	GATA3 Rabbit pAb	ER1902-69	Human, Rat	WB,IHC-P,FC
	GATA3 Recombinant Mouse mAb [A3G3-R]	HA601189	Human	WB,IF-Cell,IHC-P,FC
	GATA3 Mouse mAb [A3G3]	EM1902-33	Human	WB,IF-Cell,IHC-P
	GATA3 Recombinant Mouse mAb [PD00-41]	HA601132	Human	WB,IHC-P,IF-Cell
	GATA3 Mouse mAb [A3G6]	EM1902-23	Human	WB
	GATA3 Recombinant Mouse mAb [A3G3-R] - BSA and Azide free	HA610069	Human	WB,IF-Cell,IHC-P,FC
	GITR Recombinant Rabbit mAb [PSH04-06] - BSA and Azide free	HA750929	Human	WB,IF-Cell,FC
	GITR Recombinant Rabbit mAb [PSH04-06]	HA722087	Human	WB,IF-Cell,FC
Relative expression	Granzyme B Recombinant Rabbit mAb [PSH14-66]	HA723664	Human, Mouse	WB,IHC-P
	Granzyme B Mouse mAb [4G9]	RT1259	Human, Mouse, Rat	WB,IP,IF,IHC-P
Relative expression	Granzyme B Recombinant Rabbit mAb [PSH14-66] - BSA and Azide free	HA751527	Human, Mouse	WB,IHC-P
	Granzyme B Rabbit pAb	HA500252	Human	IHC-P,mIHC
	Granzyme B Rabbit pAb	IRS032	Human	mIHC
	FITC Conjugated Human CD200 Recombinant Rabbit mAb [PSH01-48]	HA720208F	Human	FC
	Biotin Conjugated Human CD200 Recombinant Rabbit mAb [PSH05-60]	HA722591B	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human CD200 Recombinant Rabbit mAb [PSH05-59] - BSA and Azide free (Capture)	HA722379	Human	ELISA(Cap)
	Human CD200 Recombinant Rabbit mAb [PSH05-60] - BSA and Azide free (Detector)	HA722380	Human	ELISA(Det)
	Human CD200R1 Recombinant Rabbit mAb [PSH05-27] - BSA and Azide free (Capture)	HA722256	Human	ELISA(Cap)
	Human CD200R1 Recombinant Rabbit mAb [PSH05-28] - BSA and Azide free (Detector)	HA722257	Human	ELISA(Det)
	Biotin Conjugated Human CD200R1 Recombinant Rabbit mAb [PSH05-28]	HA722598B	Human	ELISA(Det)
Cell treatment	Human CD200R1 Recombinant Rabbit mAb [PSH05-02]	HA722211	Human	IF-Cell,FC
	FITC Conjugated Human CD200R1 Recombinant Rabbit mAb [PSH05-02]	HA720212F	Human	FC
	Biotin Conjugated Human CD276 Recombinant Rabbit mAb [PSH05-93] - Detector	HA722588B	Human	ELISA(Det),ELISA
	Human CD276 Recombinant Rabbit mAb [PSH05-93] - BSA and Azide free (Detector)	HA722492	Human	ELISA(Det)
	Human CD276 Recombinant Rabbit mAb [PSH05-92] - BSA and Azide free (Capture)	HA722491	Human	ELISA(Cap)
	Human CD27/TNFRSF7 Recombinant Rabbit mAb [PSH08-31] - BSA and Azide free (Capture)	HA722973	Human	ELISA(Cap)
	Human CD27/TNFRSF7 Recombinant Rabbit mAb [PSH08-32] - BSA and Azide free (Capture)	HA722974	Human	ELISA(Cap)
	Biotin Conjugated Human CD27/TNFRSF7 Recombinant Rabbit mAb [PSH08-33] - Detector	HA722976B	Human	ELISA(Det),ELISA
	Human CD27/TNFRSF7 Recombinant Rabbit mAb [PSH08-33] - BSA and Azide free (Detector)	HA722975	Human	ELISA(Det)
	Biotin Conjugated Human CD28 Recombinant Rabbit mAb [PSH12-51] - Detector	HA723459B	Human	ELISA(Det),ELISA
	Human CD28 Recombinant Rabbit mAb [PSH12-50] - BSA and Azide free (Capture)	HA723457	Human	ELISA(Cap)
	Human CD28 Recombinant Rabbit mAb [PSH12-51] - BSA and Azide free (Detector)	HA723458	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human CD3 Recombinant Mouse mAb [PSH04-92]	HA601337	Human	FC
	APC Conjugated Human CD3 Recombinant Mouse mAb [PSH04-92]	HA600118F	Human	FC
	FITC Conjugated Human CD3 Recombinant Mouse mAb [PSH04-92]	HA600116F	Human	FC
	Percp-Cy5.5 Conjugated Human CD3 Recombinant Mouse mAb [PSH04-92]	HA601548F	Human	FC
	PE Conjugated Human CD3 Recombinant Mouse mAb [PSH04-92]	HA600117F	Human	FC
	FITC Conjugated Human CD38 Recombinant Rabbit mAb [SN07-20]	HA720196F	Human	FC
	PE Conjugated Human CD4 Recombinant Mouse mAb [PSH04-93]	HA600120F	Human	FC
	PE-Cy7 Conjugated Human CD4 Recombinant Mouse mAb [PSH04-93]	HA601542F	Human	FC
	APC Conjugated Human CD4 Recombinant Mouse mAb [PSH04-93]	HA600121F	Human	FC
	FITC Conjugated Human CD4 Recombinant Mouse mAb [PSH04-93]	HA600119F	Human	FC
	Human CD4 Recombinant Mouse mAb [PSH04-93]	HA601338	Human	FC
	iFluor™ 555 Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA601556F	Human	FC
	iFluor™ 488 Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600109F	Human	FC
	APC Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600127F	Human	FC
	PE Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600126F	Human	FC
	FITC Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600112F	Human	FC
	Human CD45 Recombinant Mouse mAb [PSH03-64] - BSA and Azide free	HA610158	Human	IF-Cell,FC
	Human CD45 Recombinant Mouse mAb [PSH03-64]	HA601285	Human	IF-Cell,FC
	FITC Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600125F	Human	FC
	Percp-Cy5.5 Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600128F	Human	FC
	PE Conjugated Human CD45RA Recombinant Mouse mAb [PSH04-97]	HA601525F	Human	FC
	Human CD45RA Recombinant Mouse mAb [PSH04-97]	HA601341	Human	FC,IF-Cell,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	FITC Conjugated Human CD45RA Recombinant Mouse mAb [PSH04-97]	HA600131F	Human	FC
Relative expression	Human CD5 Recombinant Rabbit mAb [PSH05-01]	HA722210	Human	IF-Cell,FC
	FITC Conjugated Human CD5 Recombinant Rabbit mAb [PSH05-01]	HA720211F	Human	FC
	Human CD83 Recombinant Rabbit mAb [PSH20-09] - BSA and Azide free (Detector)	HA725306	Human	ELISA(Det)
	Biotin Conjugated Human CD83 Recombinant Rabbit mAb [PSH20-09] - Detector	HA725306B	Human	ELISA(Det),ELISA
	Human CD83 Recombinant Rabbit mAb [PSH20-08] - BSA and Azide free (Capture)	HA725305	Human	ELISA(Cap)
	APC-Cy7 Conjugated Human CD8 alpha Recombinant Mouse mAb [PSH04-94]	HA601541F	Human	FC
	PE Conjugated Human CD8 alpha Recombinant Mouse mAb [PSH04-94]	HA600123F	Human	FC
	FITC Conjugated Human CD8 alpha Recombinant Mouse mAb [PSH04-94]	HA600122F	Human	FC
	APC Conjugated Human CD8 alpha Recombinant Mouse mAb [PSH04-94]	HA600124F	Human	FC
	Human CD8 alpha Recombinant Mouse mAb [PSH04-94]	HA601339	Human	FC
	Human CTLA-4 Recombinant Rabbit mAb [PSH08-13] - BSA and Azide free (Capture)	HA722954	Human	ELISA(Cap)
	Biotin Conjugated Human CTLA-4 Recombinant Rabbit mAb [PSH08-14]	HA722956B	Human	ELISA(Det)
	Human CTLA-4 Recombinant Rabbit mAb [PSH08-14] - BSA and Azide free (Detector)	HA722955	Human	ELISA(Det)
	Human GITR Recombinant Rabbit mAb [PSH05-63] - BSA and Azide free (Capture)	HA722383	Human	ELISA(Cap)
	Human GITR Recombinant Rabbit mAb [PSH05-64] - BSA and Azide free (Detector)	HA722384	Human	ELISA(Det)
	Biotin Conjugated Human GITR Recombinant Rabbit mAb [PSH05-64]	HA722593B	Human	ELISA(Det)
	Human Granzyme B Recombinant Rabbit mAb [PSH15-99] - BSA and Azide free (Detector)	HA723804	Human	ELISA(Det)
	Biotin Conjugated Human Granzyme B Recombinant Rabbit mAb [PSH15-99] - Detector	HA723804B	Human	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human Granzyme B Recombinant Rabbit mAb [PSH15-98] - BSA and Azide free (Capture)	HA723803	Human	ELISA(Cap)
	Human ICAM1 Recombinant Rabbit mAb [PSH10-15] - BSA and Azide free (Capture)	HA723184	Human	ELISA(Cap)
	Biotin Conjugated Human ICAM1 Recombinant Rabbit mAb [PSH10-16] - Detector	HA723186B	Human	ELISA(Det),ELISA
	HRP Conjugated Human ICAM1 Recombinant Rabbit mAb [PSH10-16] - Detector	HA723474H	Human	ELISA(Det),ELISA
	Human ICAM1 Recombinant Rabbit mAb [PSH10-16] - BSA and Azide free (Detector)	HA723185	Human	ELISA(Det)
	Human ICOSLG/B7-H2/CD275 Recombinant Rabbit mAb [PSH11-32] - BSA and Azide free (Capture)	HA723317	Human	ELISA(Cap)
	Biotin Conjugated Human ICOSLG/B7-H2/CD275 Recombinant Rabbit mAb [PSH11-35] - Detector	HA723322B	Human	ELISA(Det),ELISA
	Human ICOSLG/B7-H2/CD275 Recombinant Rabbit mAb [PSH11-34] - BSA and Azide free (Capture)	HA723320	Human	ELISA(Cap)
	Human ICOSLG/B7-H2/CD275 Recombinant Rabbit mAb [PSH11-35] - BSA and Azide free (Detector)	HA723321	Human	ELISA(Det)
	Human ICOSLG/B7-H2/CD275 Recombinant Rabbit mAb [PSH11-33] - BSA and Azide free (Detector)	HA723318	Human	ELISA(Det)
	Biotin Conjugated Human ICOSLG/B7-H2/CD275 Recombinant Rabbit mAb [PSH11-33] - Detector	HA723319B	Human	ELISA(Det),ELISA
	Human IL-10 Recombinant Rabbit mAb [PS00-98] - BSA and Azide free (Capture)	HA721236	Human	ELISA(Cap)
	Biotin Conjugated Human IL-10 Recombinant Rabbit mAb [PS00-81]	HA722652B	Human	ELISA(Det)
	Human IL-10 Recombinant Rabbit mAb [PS00-81] - BSA and Azide free (Detector)	HA721237	Human	ELISA(Det)
	Human IL-2 Receptor alpha Recombinant Rabbit mAb [PSH05-66] - BSA and Azide free (Detector)	HA722386	Human	ELISA(Det)
	Biotin Conjugated Human IL-2 Receptor alpha Recombinant Rabbit mAb [PSH05-66]	HA722594B	Human	ELISA(Det)
	Human IL-2 Receptor alpha Recombinant Rabbit mAb [PS00-52] - BSA and Azide free (Capture)	HA721280	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human IL-2 Receptor alpha Recombinant Rabbit mAb [PSH05-65] - BSA and Azide free (Capture)	HA722385	Human	ELISA(Cap)
	Human IL-2 Receptor alpha Recombinant Rabbit mAb [PS01-03] - BSA and Azide free (Detector)	HA721281	Human	ELISA(Det)
	Biotin Conjugated Human IL-2 Receptor alpha Recombinant Rabbit mAb [PS01-03]	HA722650B	Human	ELISA(Det)
	Human IL-4 Recombinant Rabbit mAb [PS00-67] - BSA and Azide free (Capture)	HA721230	Human	ELISA(Cap)
	Human IL-4 Recombinant Rabbit mAb [PS00-66] - BSA and Azide free (Detector)	HA721231	Human	ELISA(Det)
	Biotin Conjugated Human IL-4 Recombinant Rabbit mAb [PS00-66] - Detector	HA721368	Human	ELISA(Det)
	HRP Conjugated Human IL-4 Recombinant Rabbit mAb [PS00-66] - Detector	HA721295	Human	ELISA(Det)
	Human Interferon gamma Recombinant Rabbit mAb [PS01-75] - BSA and Azide free (Capture)	HA721574	Human	ELISA(Cap)
	Human Interferon gamma Recombinant Rabbit mAb [PS01-76] - BSA and Azide free (Detector)	HA721575	Human	ELISA(Det)
	Human Interferon gamma Recombinant Rabbit mAb [PSH01-46] - BSA and Azide free (Capture)	HA721684	Human	ELISA(Cap)
	Human Interferon gamma Recombinant Rabbit mAb [PSH01-47] - BSA and Azide free (Detector)	HA721685	Human	ELISA(Det)
	Human Interferon gamma Recombinant Rabbit mAb [PSH01-45] - BSA and Azide free (Capture)	HA721683	Human	ELISA(Cap)
	Biotin Conjugated Human Interferon gamma Recombinant Rabbit mAb [PS01-76] - Detector	HA721576	Human	ELISA(Det),ELISA
	Biotin Conjugated Human Interferon gamma Recombinant Rabbit mAb [PSH01-47]	HA722645B	Human	ELISA(Det)
	Biotin Conjugated Human KIM-1/HAVCR1 Recombinant Rabbit mAb [PSH18-81]	HA725250B	Human	ELISA(Det),ELISA
	Human KIM-1/HAVCR1 Recombinant Rabbit mAb [PSH18-81] - BSA and Azide free (Detector)	HA725250	Human	ELISA(Det)
	Human KIM-1/HAVCR1 Recombinant Rabbit mAb [PSH18-80] - BSA and Azide free (Capture)	HA725249	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human LAG-3 Recombinant Rabbit mAb [PSH05-30]	HA722596B	Human	ELISA(Det)
	Human LAG-3 Recombinant Rabbit mAb [PSH05-30] - BSA and Azide free (Detector)	HA722259	Human	ELISA(Det)
	Human LAG-3 Recombinant Rabbit mAb [PSH05-29] - BSA and Azide free (Capture)	HA722258	Human	ELISA(Cap)
	Human Lipocalin-2 Recombinant Rabbit mAb [PSH13-67] - BSA and Azide free (Capture)	HA725109	Human	ELISA(Cap)
	HRP Conjugated Human Lipocalin-2 Recombinant Rabbit mAb [PSH13-68] - Detector	HA725111H	Human	ELISA(Det),ELISA
	Biotin Conjugated Human Lipocalin-2 Recombinant Rabbit mAb [PSH13-68] - Detector	HA725111B	Human	ELISA(Det),ELISA
	Human Lipocalin-2 Recombinant Rabbit mAb [PSH13-68] - BSA and Azide free (Detector)	HA725110	Human	ELISA(Det)
	Human Neuropilin-1 Recombinant Rabbit mAb [PSH14-73] - BSA and Azide free (Capture)	HA725112	Human	ELISA(Cap)
	HRP Conjugated Human Neuropilin-1 Recombinant Rabbit mAb [PSH14-74] - Detector	HA725114H	Human	ELISA(Det),ELISA
	Human Neuropilin-1 Recombinant Rabbit mAb [PSH14-74] - BSA and Azide free (Detector)	HA725113	Human	ELISA(Det)
	Biotin Conjugated Human OX40/TNFRSF4 Recombinant Rabbit mAb [PSH10-18] - Detector	HA723189B	Human	ELISA(Det),ELISA
	Human OX40/TNFRSF4 Recombinant Rabbit mAb [PSH10-17] - BSA and Azide free (Capture)	HA723187	Human	ELISA(Cap)
	Human OX40/TNFRSF4 Recombinant Rabbit mAb [PSH10-18] - BSA and Azide free (Detector)	HA723188	Human	ELISA(Det)
	Human PD1 Recombinant Rabbit mAb [PSH13-39] - BSA and Azide free (Detector)	HA723524	Human	ELISA(Det)
	Human PD1 Recombinant Rabbit mAb [PSH13-38] - BSA and Azide free (Detector)	HA723522	Human	ELISA(Det)
	Human PD1 Recombinant Rabbit mAb [SJ01-91] - BSA and Azide free (Capture)	HA723521	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human PD1 Recombinant Rabbit mAb [PSH13-38] - Detector	HA723523B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human PD1 Recombinant Rabbit mAb [PSH13-39] - Detector	HA723525B	Human	ELISA(Det),ELISA
	Human RANKL Recombinant Rabbit mAb [PSH16-59] - BSA and Azide free (Detector)	HA723861	Human	ELISA(Det)
	Human RANKL Recombinant Rabbit mAb [PSH16-58] - BSA and Azide free (Capture)	HA723860	Human	ELISA(Cap)
	Human RANKL Recombinant Rabbit mAb [PSH16-60] - BSA and Azide free (Capture)	HA723862	Human	ELISA(Cap)
	Biotin Conjugated Human RANKL Recombinant Rabbit mAb [PSH16-59] - Detector	HA723861B	Human	ELISA(Det),ELISA
	ICAM1 Mouse mAb [J0-F7]	M1511-6	Human	WB,IF-Cell,FC
Relative expression, KD	ICAM1 Recombinant Rabbit mAb [PSH09-31]	HA723089	Human	WB,IHC-P,IF-Tissue
	ICAM1 Rabbit pAb	ER131207	Human	WB,IF-Cell,IHC-P,FC
Relative expression	ICAM1 Recombinant Rabbit mAb [ST0487] - BSA and Azide free	HA750176	Human	WB,IF-Tissue,IHC-P
KD	ICAM1 Mouse mAb [F11-A4]	M1505-7	Human	WB,IF-Cell,FC
Relative expression	ICAM1 Recombinant Rabbit mAb [ST0487]	ET1609-46	Human	WB,IF-Tissue,IHC-P
Cell treatment, Relative expression	ICOSLG / B7-H2 / CD275 Recombinant Rabbit mAb [JE52-85]	HA721299	Human	WB
Cell treatment	IL-10 Recombinant Rabbit mAb [JE59-71]	HA722032	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	IL-10 Rabbit pAb	ER65464	Human, Mouse, Rat	WB,IHC-P
	IL-10 Rabbit pAb	ER1911-19	Human, Mouse	WB,IHC-P,IF-Tissue
	IL-2 receptor alpha Rabbit pAb	HA500076	Human, Rat	WB
	IL-4 Rabbit pAb	ER1706-55	Human	WB,IF-Cell,IHC-P
	IL-4 Recombinant Rabbit mAb [JA86-10]	ET1704-28	Human	WB,IF-Cell,IF-Tissue,IHC-P
	IL-4 Mouse mAb [A7D7]	HA601034	Human	ELISA
KD	Integrin alpha 2 Recombinant Rabbit mAb [SN0752]	ET1611-57	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
KD	Integrin alpha 2 Recombinant Rabbit mAb [SN0752] - BSA and Azide free	HA750263	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Integrin alpha E Recombinant Rabbit mAb [SN62-06]	ET1611-27	Human	IHC-P
	Interferon gamma Recombinant Rabbit mAb [JM10-10]	ET1703-17	Human	WB,FC
	Interferon gamma Recombinant Rabbit mAb [JM10-10] - BSA and Azide free	HA750373	Human	WB,FC
Cell treatment	LAG-3 Recombinant Rabbit mAb [PSH0-42] - BSA and Azide free	HA750620	Mouse	WB,IF-Cell
	LAG-3 Recombinant Rabbit mAb	IRS020	Human	mIHC
Cell treatment	LAG-3 Recombinant Rabbit mAb [PSH0-43] - BSA and Azide free	HA750621	Human, Mouse	WB,IHC-P
Cell treatment	LAG-3 Recombinant Rabbit mAb [PSH0-43]	HA721345	Human, Mouse	WB,IHC-P
Cell treatment	LAG-3 Recombinant Rabbit mAb [PSH0-44]	HA721346	Human, Mouse	WB,IHC-P,IF-Cell
Cell treatment	LAG-3 Recombinant Rabbit mAb [PS01-33]	HA721358	Human	IHC-P,IF-Cell,WB,mIHC
Cell treatment	LAG-3 Recombinant Rabbit mAb [PSH0-44] - BSA and Azide free	HA750622	Human, Mouse	WB,IHC-P,IF-Cell
Cell treatment	LAG-3 Recombinant Rabbit mAb [PSH0-42]	HA721344	Mouse	WB,IF-Cell
	Lck Rabbit pAb	ER63849	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Relative expression	Lck Recombinant Rabbit mAb [PSH09-25]	HA723083	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC
	Lck Rabbit pAb	HA500280	Human, Rat	WB,FC
Relative expression	Lck Recombinant Rabbit mAb [PSH09-25] - BSA and Azide free	HA751279	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC
	Lipocalin-2 Recombinant Rabbit mAb [JM10-67]	ET1703-39	Human	WB,IF-Tissue,IHC-P
	LSP1 Recombinant Rabbit mAb [SN06-25]	ET1611-9	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	LY6A Recombinant Rabbit mAb	IRS166RB	Mouse	mIHC
	LY6A Recombinant Rabbit mAb [JM53-11]	ET1703-67	Mouse	WB,IHC-P
	Lysozyme Recombinant Mouse mAb [A5A8-R]	HA601229	Human	WB,IHC-P
	Lysozyme Recombinant Mouse mAb [A5A8-R] - BSA and Azide free	HA610104	Human	WB,IHC-P
	Lysozyme Recombinant Rabbit mAb [ST50-02] - BSA and Azide free	HA750170	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Lysozyme Mouse mAb [A5A9]	EM1901-99	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Lysozyme Recombinant Rabbit mAb	IRS168RB	Mouse	mIHC
	Lysozyme Rabbit pAb	ER1803-88	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Lysozyme Mouse mAb [A5A8]	EM1901-98	Human	WB,IHC-P
	Lysozyme Recombinant Rabbit mAb [ST50-02]	ET1609-35	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP,mIHC
	PE Conjugated Anti-Mouse CD3 Recombinant Antibody [PSH04-98] - Rabbit IgG (Chimeric)	HA720200F	Mouse	FC
	Anti-Mouse CD3 Recombinant Antibody [PSH04-98] - Rabbit IgG (Chimeric)	HA723347	Mouse	FC
	APC Conjugated Anti-Mouse CD3 Recombinant Antibody [PSH04-98] - Rabbit IgG (Chimeric)	HA720201F	Mouse	FC
	FITC Conjugated Anti-Mouse CD3 Recombinant Antibody [PSH04-98] - Rabbit IgG (Chimeric)	HA720199F	Mouse	FC
	APC Conjugated Anti-Mouse CD4 Recombinant Antibody [PSH04-99] - Rabbit IgG (Chimeric)	HA720204F	Mouse	FC
	FITC Conjugated Anti-Mouse CD4 Recombinant Antibody [PSH04-99] - Rabbit IgG (Chimeric)	HA720202F	Mouse	FC
	PE Conjugated Anti-Mouse CD4 Recombinant Antibody [PSH04-99] - Rabbit IgG (Chimeric)	HA720203F	Mouse	FC
	Mouse CD86 Recombinant Rabbit mAb [PSH13-63] - BSA and Azide free (Capture)	HA723552	Mouse	ELISA(Cap)
	Mouse CD86 Recombinant Rabbit mAb [PSH13-64] - BSA and Azide free (Detector)	HA723553	Mouse	ELISA(Det)
	Mouse CD86 Recombinant Rabbit mAb [PSH13-65] - BSA and Azide free (Capture)	HA723555	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse CD86 Recombinant Rabbit mAb [PSH13-66] - Detector	HA723557B	Mouse	ELISA(Det),ELISA
	Mouse CD86 Recombinant Rabbit mAb [PSH13-66] - BSA and Azide free (Detector)	HA723556	Mouse	ELISA(Det)
	Biotin Conjugated Mouse CD86 Recombinant Rabbit mAb [PSH13-64] - Detector	HA723554B	Mouse	ELISA(Det),ELISA
	Anti-Mouse CD8 alpha Recombinant Antibody [PSH05-00] - Rabbit IgG (Chimeric) - BSA and Azide free	HA751593	Mouse	FC,IHC-Fr
	APC Conjugated Anti-Mouse CD8 alpha Recombinant Antibody [PSH05-00] - Rabbit IgG (Chimeric)	HA720207F	Mouse	FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	PE Conjugated Anti-Mouse CD8 alpha Recombinant Antibody [PSH05-00] - Rabbit IgG (Chimeric)	HA720206F	Mouse	FC
	Mouse CD8 alpha Recombinant Rat mAb [PSH16-78]	HA601543	Mouse	FC
	FITC Conjugated Anti-Mouse CD8 alpha Recombinant Antibody [PSH05-00] - Rabbit IgG (Chimeric)	HA720205F	Mouse	FC
	Anti-Mouse CD8 alpha Recombinant Antibody [PSH05-00] - Rabbit IgG (Chimeric)	HA723820	Mouse	FC,IHC-Fr
	Mouse CD8 alpha Recombinant Rat mAb [PSH16-78] - BSA and Azide free	HA610359	Mouse	FC
	Mouse IL-10 Recombinant Rabbit mAb [PS01-85] - BSA and Azide free (Capture)	HA721517	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-10 Recombinant Rabbit mAb [PS01-93] - Detector	HA721519	Mouse	ELISA(Det)
	Mouse IL-10 Recombinant Rabbit mAb [PS01-93] - BSA and Azide free (Detector)	HA721518	Mouse	ELISA(Det)
Cell treatment	Mouse RANKL Recombinant Rabbit mAb [PSH07-91] - BSA and Azide free (Capture)	HA722923	Mouse	ELISA(Cap)
Cell treatment	Biotin Conjugated Mouse RANKL Recombinant Rabbit mAb [PSH07-92]	HA722925B	Mouse	ELISA(Det)
Cell treatment	Mouse RANKL Recombinant Rabbit mAb [PSH07-92] - BSA and Azide free (Detector)	HA722924	Mouse	ELISA(Det)
Cell treatment	Mouse RANKL Recombinant Rabbit mAb [PSH07-90] - BSA and Azide free (Capture)	HA722922	Mouse	ELISA(Cap)
	Munc13-4 Rabbit pAb	HA500052	Human, Rat	WB,IF-Cell,IHC-P
Relative expression, KD	Neuropilin 1 Recombinant Rabbit mAb [ST05-30] - BSA and Azide free	HA750189	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
Relative expression, KD	Neuropilin 1 Recombinant Rabbit mAb [ST05-30]	ET1609-69	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Neuropilin 1 Mouse mAb [A3-D1]	M0910-5	Human, Mouse	WB,IF-Cell,IHC-P
Relative expression	NFAT1 Recombinant Rabbit mAb [JA11-08] - BSA and Azide free	HA750405	Human	WB,IF-Tissue,IHC-P
Relative expression	NFAT1 Recombinant Rabbit mAb [JA11-08]	ET1704-14	Human	WB,IF-Tissue,IHC-P
	NFAT2 Recombinant Rabbit mAb [JA81-03]	ET1704-45	Human	WB,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	NK-p44 Recombinant Rabbit mAb [JE61-04]	HA720070	Human	WB
	NLRC3 Rabbit pAb	ER1706-45	Human, Mouse	WB,IF-Cell,IHC-P,FC
	NLRC3 Rabbit pAb	ER1706-69	Human, Mouse	IF-Cell,IHC-P,FC
Relative expression	NOX2 / gp91phox Recombinant Rabbit mAb [PSH10-59]	HA723224	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IHC-Fr,IP
Relative expression	NOX2 / gp91phox Rabbit pAb	ER1511-34	Human, Mouse	IF-Cell,IHC-P,FC,WB
	NOX2 / gp91phox Recombinant Rabbit mAb [SN07-16] - BSA and Azide free	HA750255	Human, Mouse, Rat	WB
	NOX2 / gp91phox Recombinant Rabbit mAb [SN07-16]	ET1611-44	Human, Mouse, Rat	WB
Relative expression	NOX2 / gp91phox Recombinant Rabbit mAb [PSH10-59] - BSA and Azide free	HA751362	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IHC-Fr,IP
	NOXA2 Recombinant Rabbit mAb [JG86-31]	ET7108-72	Human, Mouse, Rat	WB,IHC-P,IP
	Orai1 Rabbit pAb	ER1803-11	Human, Mouse	WB,IHC-P
Cell treatment	PD1 Mouse mAb [1F2]	EM1707-60	Human	WB,IF-Cell,IHC-P,FC
	PD1 Mouse mAb	IRS190MS	Human	mIHC
	PD1 Recombinant Rabbit mAb [SJ01-91]	ET1606-41	Human	WB,IHC-P
	PD1 Recombinant Rabbit mAb [SJ01-91] - BSA and Azide free	HA750103	Human	WB,IHC-P
	PD1 Blocking Mouse mAb [3A11]	EM1701-17	Human	ELISA,FC
	Perforin Rabbit pAb	ER1803-77	Human	WB,IF-Cell,IHC-P,FC
Relative expression	Perforin-1 Recombinant Rabbit mAb [PSH15-92]	HA723797	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Relative expression	Perforin-1 Recombinant Rabbit mAb [PSH15-92] - BSA and Azide free	HA751583	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-GATA3 (S308) Recombinant Rabbit mAb [ST44-09]	ET1609-17	Human, Mouse, Rat	WB,IHC-P,IP,IF-Tissue
Cell treatment, Relative expression	Phospho-LYN (Y397)/LCK (Y394)/HCK (Y411)/BLK (Y389) Recombinant Rabbit mAb [PSH09-43]	HA723105	Human, Mouse, Rat	WB,IF-Cell
Cell treatment, Relative expression	Phospho-LYN (Y397)/LCK (Y394)/HCK (Y411)/BLK (Y389) Recombinant Rabbit mAb [PSH09-43] - BSA and Azide free	HA751290	Human, Mouse, Rat	WB,IF-Cell
	Phospho-SHIP (Y1020) Recombinant Rabbit mAb [JJ082-03]	ET1701-5	Human	WB,IF-Cell,IP
	P-Selectin / CD62P Recombinant Rabbit mAb [PSH16-70] - BSA and Azide free (Capture)	HA723872	Human, Mouse, Rat	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	P-Selectin / CD62P Recombinant Rabbit mAb [PSH19-03]	HA724039	Human, Mouse, Rat	WB,IHC-P,FC,IHC-Fr,IP
	HRP Conjugated P-Selectin / CD62P Recombinant Rabbit mAb [PSH16-71] - Detector	HA723873H	Human, Mouse, Rat	ELISA(Det),ELISA
	Biotin Conjugated P-Selectin / CD62P Recombinant Rabbit mAb [PSH16-71] - Detector	HA723873B	Human, Mouse, Rat	ELISA(Det),ELISA
	P-Selectin / CD62P Recombinant Rabbit mAb [PSH16-71] - BSA and Azide free (Detector)	HA723873	Human, Mouse, Rat	ELISA(Det)
	PSMB10 Rabbit pAb	HA500393	Human, Mouse	WB,IHC-P,IF-Cell
Cell treatment	PSMB10 Recombinant Rabbit mAb [PSH01-58] - BSA and Azide free	HA750734	Human, Mouse	WB,IHC-P
Cell treatment	PSMB10 Recombinant Rabbit mAb [PSH01-58]	HA721701	Human, Mouse	WB,IHC-P
	PTPRC Mouse mAb [C11-C10]	EM1708-69	Human	WB,IHC-P
	RANKL Rabbit pAb	HA500369	Human, Mouse, Rat	WB,IF-Cell
	Rat IL-10 Recombinant Rabbit mAb [PSH02-22] - BSA and Azide free (Detector)	HA721787	Rat	ELISA(Det)
	Rat IL-10 Recombinant Rabbit mAb [PSH02-21] - BSA and Azide free (Capture)	HA721786	Rat	ELISA(Cap)
	Biotin Conjugated Rat IL-10 Recombinant Rabbit mAb [PSH02-22]	HA722637B	Rat	ELISA(Det)
	Rat KIM-1 / HAVCR1 Recombinant Rabbit mAb [PSH20-49] - BSA and Azide free (Capture)	HA725381	Rat	ELISA(Cap)
	Biotin Conjugated Rat KIM-1 / HAVCR1 Recombinant Rabbit mAb [PSH20-50] - Detector	HA725382B	Rat	ELISA(Det),ELISA
	Rat KIM-1 / HAVCR1 Recombinant Rabbit mAb [PSH20-50] - BSA and Azide free (Detector)	HA725382	Rat	ELISA(Det)
	SCA-1 Rabbit pAb	0804-10	Mouse	WB,IF-Cell
	SET Recombinant Rabbit mAb [JB60-34]	ET7107-18	Human, Mouse	WB,IF-Cell,IHC-P,IP,IF-Tissue
	SHIP Recombinant Rabbit mAb [SY11-08]	ET1605-30	Human	WB,IHC-P,IP,FC
	T-bet / Tbx21 Rabbit pAb	ER1802-65	Human	WB,FC,IHC-P
Relative expression	T-bet / Tbx21 Recombinant Rabbit mAb [JE60-20]	HA721901	Human	WB,IHC-P,IF-Cell,FC
	Telomerase reverse transcriptase Recombinant Rabbit mAb [JF0586]	ET1702-22	Human	WB,FC
	TERT Mouse mAb [G9-H2]	EM1707-13	Human	IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	TFE3 Rabbit pAb	ER65143	Human, Mouse	WB,IHC-P
	TFE3 Recombinant Rabbit mAb [JE60-60]	HA721896	Human, Mouse, Rat	WB,IHC-P,FC
	Themis Recombinant Rabbit mAb [JE55-73]	ET7111-20	Human, Rat	WB,IF-Cell,IHC-P,FC
	Thy-1 (OX7) Mouse mAb [1G2]	RT1615	Mouse, Rat	WB,IP,IF,IHC-P,FC
	TIA1 Recombinant Rabbit mAb [JM42-11] - BSA and Azide free	HA750385	Human, Mouse	WB,IF-Tissue,IHC-P,IP
	TIA1 Recombinant Rabbit mAb [JM42-11]	ET1703-59	Human, Mouse	WB,IF-Tissue,IHC-P,IP
	TIM 1 Recombinant Rabbit mAb [JE39-41]	HA721535	Human, Mouse, Rat	WB
	TLR4 Rabbit pAb	ER1706-43	Human, Mouse, Rat	IF-Cell,IHC-P,FC
	TLR4 Mouse mAb	RT1666	Human, Mouse, Rat	WB
	TNFRSF18 Mouse mAb [4H2D6]	EM1711-16	Human	WB,IHC-P,FC
	TNFSF11 Mouse mAb [8A7B9]	EM1710-42	Human, Monkey	WB,IHC-P,IF-Cell,FC
	TNFSF14 Rabbit pAb	HA500384	Human, Mouse	WB,IHC-P
	TNR18 Rabbit pAb	ER62416	Human	WB
Relative expression	ZAP70 Mouse mAb [A1B7]	EM1902-28	Human, Mouse, Rat	WB,IF-Cell,FC
Relative expression	ZAP70 Recombinant Mouse mAb [A1B5-R] - BSA and Azide free	HA610096	Human, Mouse, Rat	WB,IHC-P
Relative expression	ZAP70 Recombinant Mouse mAb [A1B5-R]	HA601218	Human, Mouse, Rat	WB,IHC-P
Relative expression	ZAP70 Recombinant Rabbit mAb [JU08-39]	ET1706-42	Human	WB,IF-Cell,IHC-P,FC,IP
Relative expression	ZAP70 Mouse mAb [A1B5]	EM1901-27	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell
Relative expression	ZAP70 Recombinant Rabbit mAb [PSH11-52]	HA723344	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层