

肠道屏障功能研究-产品清单

高级验证	产品名	货号	反应种属	验证应用
	Actin Recombinant Rabbit mAb [JJ09-29]	ET1701-80	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P
	Actin Recombinant Rabbit mAb [JJ09-29] - BSA and Azide free	HA750323	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P
	HRP Conjugated Actin Recombinant Rabbit mAb [JJ09-29]	HA721186	Human, Mouse, Rat	WB
	Actin Recombinant Antibody [JM36-62] - Human IgG1 (Chimeric)	HA721878	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Actin Recombinant Antibody [JM36-62] - Human IgG1 (Chimeric) - BSA and Azide free	HA750844	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	iFluor™ 488 Conjugated Actin Recombinant Rabbit mAb [JJ09-29]	HA720160F	Human, Mouse, Rat	IF-Cell,IF-Tissue,FC,WB
	Actin Recombinant Rabbit mAb [JF47-01]	ET1702-52	Human, Mouse, Rat, Zebrafish	WB,IHC-P,IF-Cell,FC
	ADORA2B Rabbit pAb	ER1903-44	Human, Mouse, Rat	WB
	AFP Rabbit pAb	ER80601	Human	WB,IF-Cell,IHC-P,FC
	AFP Rabbit pAb	R1601-1	Human, Rat, Mouse	WB,IF-Cell,IHC-P,FC
KD	AFP Mouse mAb [A5-G3]	EM1701-31	Human	WB,IHC-P
Relative expression	AFP Recombinant Rabbit mAb [PSH11-75] - BSA and Azide free	HA751430	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue
Relative expression	AFP Recombinant Rabbit mAb [PSH11-75]	HA723371	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue
Relative expression	AFP Recombinant Rabbit mAb [PD00-71]	HA721238	Human	IHC-P,WB
	Alpha-cardiac actin Mouse mAb [33-32]	M1206-1	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,IF-Tissue,IHC-Fr
	Alpha-cardiac actin Recombinant Mouse mAb [33-32-R]	HA601208	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Alpha-cardiac actin Rabbit pAb	0407-3	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P
	Alpha-cardiac actin Recombinant Mouse mAb [33-32-R] - BSA and Azide free	HA610086	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
KD, Relative expression	alpha smooth muscle Actin Recombinant Antibody [SY02-64] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610267	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	alpha smooth muscle Actin Recombinant Antibody [SY02-64] - Mouse IgG1 (Chimeric)	HA601546	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	alpha smooth muscle Actin Recombinant Rabbit mAb	IRS012	Human	mIHC
	alpha smooth muscle Actin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS220RT	Mouse	mIHC
	alpha smooth muscle Actin Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS219MS	Mouse	mIHC
KD	alpha smooth muscle Actin Recombinant Rabbit mAb [SY25-03]	ET1607-43	Human, Mouse, Rat	WB,IHC-P,FC,IP
KD, Relative expression	alpha smooth muscle Actin Recombinant Rabbit mAb [SY02-64] - BSA and Azide free	HA750126	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	alpha smooth muscle Actin Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS093MS	Human	mIHC
KD, Relative expression	alpha smooth muscle Actin Recombinant Antibody [SY02-64] - Rat IgG1 (Chimeric)	HA601427	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	alpha smooth muscle Actin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS094RT	Human	mIHC
KD, Relative expression	alpha smooth muscle Actin Recombinant Rabbit mAb [SY02-64]	ET1607-53	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,mIHC
	alpha smooth muscle Actin Recombinant Antibody [SY02-64] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610362	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	alpha smooth muscle Actin Recombinant Rabbit mAb	IRS048	Mouse	mIHC
	Arginase1 Recombinant Rabbit mAb	IRS082RB	Mouse	mIHC
	Arginase1 Recombinant Rabbit mAb	IRS081RB	Human	mIHC
KD, Relative expression	CDX2 Recombinant Rabbit mAb [SY09-02]	ET1605-4	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,CUT&Tag-seq
	c-kit Rabbit pAb	0108-4	Human	WB,IF-Cell
Relative expression, KD	c-Kit Recombinant Rabbit mAb [PD00-24]	HA721154	Human	IHC-P,WB,mIHC
	c-Kit Recombinant Rabbit mAb	IRS013	Human	mIHC
Relative expression, KD	c-Kit Recombinant Rabbit mAb [ST04-99]	ET1609-60	Human	WB,IHC-P
Relative expression, KD	Claudin 1 Recombinant Rabbit mAb [JE30-33]	HA721999	Human, Mouse, Rat	WB,IHC-P
KD	Claudin 1 Recombinant Rabbit mAb [PS01-43]	HA721248	Human	WB,IHC-P,IF-Tissue,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	Claudin 1 Recombinant Rabbit mAb [PS01-43] - BSA and Azide free	HA750567	Human	WB,IHC-P,IF-Tissue,IF-Cell
	Claudin 1 Mouse mAb [8G1]	RT1141	Human, Mouse, Rat	WB,IP,IF,IHC-P
	Claudin 2 Mouse mAb [8H11]	EM60133	Human, Mouse, Rat	WB,IHC-P
	COL1A1 Recombinant Rabbit mAb [ST58-04]	ET1609-68	Human, Cow	WB,IHC-P
	COL1A1 Recombinant Rabbit mAb [ST58-04] - BSA and Azide free	HA750188	Human, Cow	WB,IHC-P
	COL1A1 Recombinant Rabbit mAb [PSH06-20] - BSA and Azide free	HA751058	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell,FC,IP,IHC-Fr
	COL1A1 Recombinant mAb	IRS137RB	Human	mIHC
	COL1A1 Recombinant Rabbit mAb [PSH06-20]	HA722517	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell,FC,IP,mIHC,IHC-Fr
	COL1A1 Recombinant Rabbit mAb	IRS228RB	Mouse	mIHC
Relative expression	COL1A2 Recombinant Rabbit mAb [PSH07-21] - BSA and Azide free	HA751142	Human, Rat	WB,IF-Cell,FC
Relative expression	COL1A2 Recombinant Rabbit mAb [PSH07-21]	HA722809	Human, Rat	WB,IF-Cell,FC
	Collagen III Recombinant Rabbit mAb [JE60-79]	HA720050	Human, Mouse, Rat	WB,IF-Cell,FC
	Collagen III Rabbit pAb	ER1906-50	Human	WB,FC,IF-Cell
	iFluor™ 594 Conjugated E-Cadherin Recombinant Rabbit mAb [SY0287]	HA720164F	Human	IF-Tissue
	iFluor™ 488 Conjugated E-Cadherin Recombinant Rabbit mAb [SY0287]	HA720159F	Human	IF-Cell,FC
Relative expression	E-Cadherin Recombinant Rabbit mAb [PSH13-75] - BSA and Azide free	HA751488	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Cell,FC,IP
	E-Cadherin Recombinant Rabbit mAb	IRS030	Human	mIHC
	E-Cadherin Mouse mAb [3-F9]	M1405-3	Human, Mouse	IHC-P
	E-Cadherin Recombinant Mouse mAb	IRS229MS	Mouse	mIHC
	E-Cadherin Rabbit pAb	0407-25	Human, Mouse	WB,IF-Cell,FC
KD	E-Cadherin Mouse mAb [A0-G11-2]	EM0502	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Relative expression	E-Cadherin Rabbit pAb	ER63312	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IHC-Fr,FC
	E-Cadherin Recombinant Mouse mAb [A0-G11-2-R]	HA601143	Human, Mouse, Rat	WB,IHC,mIHC
	E-Cadherin Recombinant Mouse mAb [A0-G11-2-R] - BSA and Azide free	HA610047	Human, Mouse, Rat	WB,IHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression, KD	E-Cadherin Recombinant Rabbit mAb [SY0287] - BSA and Azide free	HA750133	Human, Cynomolgus monkey	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
Relative expression, KD	E-Cadherin Recombinant Rabbit mAb [SY0287]	ET1607-75	Human, Cynomolgus monkey	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP,mIHC
Relative expression	E-Cadherin Recombinant Rabbit mAb [PSH13-75]	HA723564	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Cell,FC,IP
	E-Cadherin Recombinant Mouse mAb	IRS031	Human	mIHC
	E-Cadherin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS284RT	Mouse	mIHC
	EpCAM Recombinant Rabbit mAb [PS01-69]	HA721644	Human	IHC-P,IF-Cell,FC
	EpCAM Recombinant Rabbit mAb [SU03-32]	ET1608-26	Human, Mouse, Rat	WB,IHC-P
	iFluor™ 555 Conjugated EpCAM Recombinant Rabbit mAb [PS01-69]	HA720218F	Human	FC
	EpCAM Recombinant Rabbit mAb	IRS254RB	Human	mIHC
	EpCAM Rabbit pAb	ER1802-80	Human, Mouse, Rat	WB,IHC-P
	iFluor™ 647 Conjugated EpCAM Recombinant Rabbit mAb [PS01-69]	HA720185F	Human	FC,IF-Cell
	iFluor™ 488 Conjugated EpCAM Recombinant Rabbit mAb [PS01-69]	HA720217F	Human	IF-Cell
	EpCAM Mouse mAb [3-E1-E2]	EM1111	Human, Mouse	WB,IF-Cell
	EpCAM Mouse mAb [7E11]	EM1712-43	Human	WB,IHC-P,ELISA
	EpCAM Recombinant Rabbit mAb [PS01-69] - BSA and Azide free	HA750703	Human	IHC-P,IF-Cell,FC
	Fibronectin Recombinant mAb	IRS144RB	Mouse	mIHC
	Fibronectin Recombinant Rabbit mAb [PSH14-22]	HA723611	Human, Mouse	WB,IF-Cell,FC
KD	Fibronectin Recombinant Rabbit mAb [JF0582]	ET1702-25	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P
	Fibronectin Recombinant Rabbit mAb [PSH14-21]	HA723610	Human, Mouse	WB,IHC-P,IF-Cell,FC
	Fibronectin Recombinant Rabbit mAb [PSH14-22] - BSA and Azide free	HA751515	Human, Mouse	WB,IF-Cell,FC
KD	Fibronectin Recombinant Rabbit mAb [JF0582] - BSA and Azide free	HA750340	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P
	Fibronectin Recombinant Rabbit mAb [PSH14-21] - BSA and Azide free	HA751514	Human, Mouse	WB,IHC-P,IF-Cell,FC
	Fibronectin Mouse mAb [3G4]	RT1224	Human, Mouse, Rat	WB,IP,IF,IHC-P
	Biotin Conjugated Human alpha Fetoprotein Recombinant Rabbit mAb [PSH14-54] - Detector	HA723646B	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human alpha Fetoprotein Recombinant Rabbit mAb [PSH14-55] - BSA and Azide free (Detector)	HA723647	Human	ELISA(Det)
	Human alpha Fetoprotein Recombinant Rabbit mAb [PSH14-53] - BSA and Azide free (Capture)	HA723644	Human	ELISA(Cap)
	Human alpha Fetoprotein Recombinant Rabbit mAb [PSH14-54] - BSA and Azide free (Detector)	HA723645	Human	ELISA(Det)
	Biotin Conjugated Human alpha Fetoprotein Recombinant Rabbit mAb [PSH14-55] - Detector	HA723648B	Human	ELISA(Det)
	Human E-Cadherin Recombinant Rabbit mAb [SY0287] - BSA and Azide free (Detector)	HA723056	Human	ELISA(Det)
	Biotin Conjugated Human E-Cadherin Recombinant Rabbit mAb [SY0287] - BSA and Azide free (Detector)	HA723057B	Human	ELISA(Det),ELISA
	Human E-Cadherin Recombinant Rabbit mAb [PSH09-02] - BSA and Azide free (Capture)	HA723055	Human	ELISA(Cap)
	Human IL-12B p40 Recombinant Rabbit mAb [PSH03-04] - BSA and Azide free (Detector)	HA721934	Human	ELISA(Det)
	Human IL-12B p40 Recombinant Rabbit mAb [PSH02-57] - BSA and Azide free (Capture)	HA721834	Human	ELISA(Cap)
	Biotin Conjugated Human IL-12B p40 Recombinant Rabbit mAb [PSH03-04]	HA722629B	Human	ELISA(Det)
	Biotin Conjugated Human IL-12B p40 Recombinant Rabbit mAb [PSH02-58]	HA722633B	Human	ELISA(Det)
	Human IL-12B p40 Recombinant Rabbit mAb [PSH02-58] - BSA and Azide free (Detector)	HA721835	Human	ELISA(Det)
	Human IL-23A Recombinant Rabbit mAb [PSH03-03] - BSA and Azide free (Capture)	HA721933	Human	ELISA(Cap)
	Human VE Cadherin Recombinant Rabbit mAb [PSH02-62] - BSA and Azide free (Detector)	HA721839	Human	ELISA(Det)
	Human VE Cadherin Recombinant Rabbit mAb [PSH02-61] - BSA and Azide free (Capture)	HA721838	Human	ELISA(Cap)
	Biotin Conjugated Human VE Cadherin Recombinant Rabbit mAb [PSH02-62] - Detector	HA722635B	Human	ELISA(Det)
	IL-12 Rabbit pAb	ER1911-22	Rat	WB,IHC-P,IF-Tissue
	IL12 Rabbit pAb	ER1911-21	Mouse, Rat	WB,IHC-P,IF-Tissue,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	IL-12A Recombinant Rabbit mAb [PSH16-10] - BSA and Azide free	HA723815	Human	ELISA
	IL12B Rabbit pAb	ER1803-31	Human, Mouse	WB
	IL-23 Recombinant Rabbit mAb [JE38-99]	HA721037	Human	WB,IHC-P,IF-Tissue
Relative expression	Liver Arginase Recombinant Rabbit mAb [PSH09-78]	HA723143	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Relative expression	Liver Arginase Recombinant Rabbit mAb [JB21-49]	ET7107-90	Human	WB,IHC-P,IP
	Liver Arginase Recombinant Rabbit mAb [PD00-05]	HA721147	Human	WB,IHC-P,IF-Tissue
Cell treatment	Liver Arginase Recombinant Rabbit mAb [SY09-06]	ET1605-8	Human, Mouse, Rat	WB
	Lysozyme Rabbit pAb	ER1803-88	Human, Mouse, Rat	WB,IHC-P
	Lysozyme Recombinant Mouse mAb [A5A8-R]	HA601229	Human	WB,IHC-P
	Lysozyme Recombinant Mouse mAb [A5A8-R] - BSA and Azide free	HA610104	Human	WB,IHC-P
	Lysozyme Mouse mAb [A5A8]	EM1901-98	Human	WB,IHC-P
	Lysozyme Recombinant Rabbit mAb [ST50-02]	ET1609-35	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP,mIHC
	Lysozyme Recombinant Rabbit mAb	IRS168RB	Mouse	mIHC
	Lysozyme Mouse mAb [A5A9]	EM1901-99	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Lysozyme Recombinant Rabbit mAb [ST50-02] - BSA and Azide free	HA750170	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Mast Cell Tryptase Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS210MS	Human	mIHC
	Mast Cell Tryptase Recombinant Rabbit mAb	IRS041	Human	mIHC
	Mast Cell Tryptase Recombinant Antibody [SC68-07] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610259	0	
	Mast Cell Tryptase Recombinant Antibody - Rat IgG1 (Chimeric)	IRS209RT	Human	mIHC
	Mast Cell Tryptase Recombinant Rabbit mAb [SC68-07]	ET1610-64	Human, Mouse	WB,IHC-P,IP,mIHC
	Mast Cell Tryptase Recombinant Antibody [SC68-07] - Rat IgG1 (Chimeric)	HA601419		
	Mast Cell Tryptase Recombinant Rabbit mAb [SC68-07] - BSA and Azide free	HA750229	Human, Mouse	WB,IHC-P,IP
	Mouse IL12 p40 Recombinant Rabbit mAb [PSH04-87] - BSA and Azide free (Capture)	HA722193	Mouse	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse IL12 p40 Recombinant Rabbit mAb [PSH04-88] - BSA and Azide free (Detector)	HA722194	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL12 p40 Recombinant Rabbit mAb [PSH04-88]	HA722603B	Mouse	ELISA(Det)
	Biotin Conjugated Mouse Pro-Collagen I alpha 1 Recombinant Rabbit mAb [PSH07-99]	HA722936B	Mouse	ELISA(Det)
	Mouse Pro-Collagen I alpha 1 Recombinant Rabbit mAb [PSH07-99] - BSA and Azide free (Detector)	HA722935	Mouse	ELISA(Det)
	HRP Conjugated Mouse Pro-Collagen I alpha 1 Recombinant Rabbit mAb [PSH07-99] - Detector	HA722936H	Mouse	ELISA(Det)
	Mouse Pro-Collagen I alpha 1 Recombinant Rabbit mAb [PSH07-98] - BSA and Azide free (Capture)	HA722934	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse/Rat E-Cadherin Recombinant Rabbit mAb [PSH17-25] - Detector	HA723912B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat E-Cadherin Recombinant Rabbit mAb [PSH17-25] - BSA and Azide free (Detector)	HA723912	Mouse, Rat	ELISA(Det)
	Mouse/Rat E-Cadherin Recombinant Rabbit mAb [PSH17-24] - BSA and Azide free (Capture)	HA723911	Mouse, Rat	ELISA(Cap)
	HRP Conjugated Mouse/Rat E-Cadherin Recombinant Rabbit mAb [PSH17-25] - Detector	HA723912H	Mouse, Rat	ELISA(Det),ELISA
	MUC2 Recombinant Rabbit mAb [JA10-92] - BSA and Azide free	HA750401	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	MUC2 Recombinant Rabbit mAb [JA10-92]	ET1704-06	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	MUC2 Recombinant Rabbit mAb	IRS175RB	Human	mIHC
	Myosin light chain kinase / MLCK Rabbit pAb	R1511-26	Human, Mouse	IF-Cell,IHC-P,FC,WB
	Myosin light chain kinase / MLCK Rabbit pAb	R1510-18	Human, Mouse	IF-Cell,IHC-P,FC
	Myosin light chain kinase / MLCK Recombinant Rabbit mAb [SU40-06]	ET1608-68	Human, Mouse, Rat	WB,IF-Tissue,IHC-P
	N Cadherin Recombinant Antibody [SY02-46] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610298	Human, Mouse, Rat	IHC-Fr,IHC-P
Relative expression, KD	N Cadherin Recombinant Rabbit mAb [SY02-46] - BSA and Azide free	HA750119	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue
	N Cadherin Recombinant Antibody [SY02-46] - Mouse IgG1 (Chimeric)	HA601464	Human, Mouse, Rat	IHC-Fr,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression, KD	N Cadherin Recombinant Rabbit mAb [SY02-46]	ET1607-37	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue
	N-Cadherin Rabbit pAb	ER0503	Human, Mouse, Rat	WB,IHC,IF
	N-Cadherin Mouse mAb [A10-4]	M1304-1	Human, Mouse	WB,IHC-P
KD	Occludin Recombinant Rabbit mAb [JJ091-08]	ET1701-76	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IP,FC,IF-Cell,IHC-P,IHC-Fr,IF-Tissue
	Occludin Recombinant Rabbit mAb	IRS183RB	Mouse	mIHC
	Occludin Rabbit pAb	R1510-33	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Pan-Actin Mouse mAb [A2E1]	HA600032	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Pan-Actin Recombinant Mouse mAb [A2E1-R] - BSA and Azide free	HA610160	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue
	Pan-Actin Recombinant Mouse mAb [A2E1-R]	HA601287	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue
	Pan Cadherin Recombinant Rabbit mAb	RE6045	Human, Mouse, Rat	WB,ICC,IHC-P,FC,IP
	Pan-Cadherin Recombinant Rabbit mAb [ST54-01] - BSA and Azide free	HA750190	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,IHC-Fr
	Pan-Cadherin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS287RT	Human	mIHC
	Pan-Cadherin Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS286MS	Mouse	mIHC
	Pan-Cadherin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS288RT	Mouse	mIHC
	Pan-Cadherin Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS285MS	Human	mIHC
	Pan-Cadherin Recombinant Rabbit mAb [ST54-01]	ET1609-70	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,IHC-Fr
	Pan-Cadherin Recombinant Rabbit mAb	IRS185RB	Mouse	mIHC
	R-cadherin Rabbit pAb	ER64773	Human, Mouse	IHC-P,IF-Cell,WB
	smooth muscle Myosin heavy chain 11 Recombinant Rabbit mAb [JA03-35] - BSA and Azide free	HA750417	Human, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC
	smooth muscle Myosin heavy chain 11 Recombinant Rabbit mAb [JA03-35]	ET1704-61	Human, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC
	smooth muscle Myosin heavy chain 11 Recombinant Rabbit mAb [JE60-21]	HA722863	Human, Mouse, Rat	WB,IHC-P
	smooth muscle Myosin heavy chain 11 Recombinant Rabbit mAb	IRS176RB	Human	mIHC
KD	VE Cadherin Recombinant Rabbit mAb [PSH03-15] - BSA and Azide free	HA750856	Human	WB,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	VE Cadherin Recombinant Rabbit mAb [PSH03-15]	HA721943	Human	WB,IF-Cell,FC
	VPAC1 Recombinant Rabbit mAb [PD00-96]	HA721209	Human, Mouse, Rat, Monkey	WB,IHC-P
	VPAC1 Rabbit pAb	ER65303	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	ZO1 Recombinant Rabbit mAb [PSH07-09]	HA722797	Human, Mouse, Monkey, Cynomolgus monkey	WB,IF-Cell
	ZO1 Recombinant Rabbit mAb	IRS217RB	Mouse	mIHC
	ZO1 Recombinant Rabbit mAb [PSH07-09] - BSA and Azide free	HA751131	Human, Mouse, Monkey, Cynomolgus monkey	WB,IF-Cell
	ZO1 Rabbit pAb	ER41204	Human, Mouse, Rat, Cynomolgus monkey	IF-Cell,IHC-P