

巨噬细胞标志物研究-产品清单

高级验证	产品名	货号	反应种属	验证应用
	Arginase1 Recombinant Rabbit mAb	IRS082RB	Mouse	mIHC
	Arginase1 Recombinant Rabbit mAb	IRS081RB	Human	mIHC
Relative expression	Axl Recombinant Rabbit mAb [PSH06-68] - BSA and Azide free	HA751097	Human, Monkey	WB,IF-Cell,FC
	Axl Mouse mAb [C4-A8]	M1310-5	Human	WB,IF-Cell
Relative expression	Axl Recombinant Rabbit mAb [PSH06-68]	HA722740	Human, Monkey	WB,IF-Cell,FC
Relative expression	Axl Rabbit pAb	R1406-3	Human, Mouse	WB,IF-Cell,IHC-P,FC
	CD11b Mouse mAb [A2-B5-D2]	EM1902-35	Human, Mouse	WB
Relative expression	CD11b Recombinant Rabbit mAb [PSH03-96]	HA722075	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	CD11b Recombinant Rabbit mAb	IRS014	Human	mIHC
Relative expression	CD11b Recombinant Rabbit mAb [PSH03-96] - BSA and Azide free	HA750921	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Relative expression	CD11b Mouse mAb [A2-9]	M1211-6	Human	WB,IHC-P,FC,IF-Tissue
Relative expression	CD11b Recombinant Rabbit mAb [JU93-81] - BSA and Azide free	HA750451	Human	WB,IP,IHC-P,IF-Cell,IF-Tissue
Relative expression	CD11b Recombinant Rabbit mAb [JU93-81]	ET1706-04	Human	WB,IP,IHC-P,IF-Cell,IF-Tissue,mIHC
Relative expression	CD11b Recombinant Rabbit mAb [JE01-17]	HA722856	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Tissue
	CD11b Mouse mAb [4-6]	EM1701-42	Human	IHC-P,IF-Cell,FC,WB,IF-Tissue
	CD11b/ITAM Mouse mAb [4G2]	RT1107	Human, Mouse	WB,IP,IF,FC
KD	CD14 Recombinant Rabbit mAb [PD01-63]	HA721321	Human, Mouse	WB,IHC-P,IF-Tissue
	CD14 Rabbit pAb	R1107-2	Human, Mouse	WB,IF-Cell,FC
	CD14 Recombinant Rabbit mAb	IRS015	Human	mIHC
KD	CD14 Recombinant Antibody [SC69-02] - Rat IgG1 (Chimeric)	HA601420		
KD	CD14 Recombinant Antibody [SC69-02] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610260	0	

KD	CD14 Recombinant Rabbit mAb [SC69-02] - BSA and Azide free	HA750232	Human, Mouse	WB,IF-Cell,IHC-P,FC
	CD14 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS106RT	Human	mIHC
KD	CD14 Recombinant Rabbit mAb [SC69-02]	ET1610-85	Human, Mouse	WB,IF-Cell,IHC-P,mIHC,FC
Relative expression	CD16 Recombinant Rabbit mAb [PD00-12]	HA721149	Human	IHC-P,IF-Tissue,IF-Cell,mIHC,WB
	CD16 Recombinant Rabbit mAb	IRS017	Human	mIHC
	CD16 Recombinant Rabbit mAb [JE49-79]	ET7109-97	Human, Mouse, Rat	WB
	CD16 Mouse mAb [Q32]	EM60067	Human	WB,IHC-P
	CD163 Recombinant Rabbit mAb [JA51-30] - BSA and Azide free	HA750412	Human	WB,IHC-P,IF-Tissue,IP
	CD163 Recombinant Rabbit mAb [JA51-30]	ET1704-43	Human	WB,IHC-P,mIHC,IF-Tissue,IP
	CD163 Rabbit pAb	IRS109RB	Mouse	mIHC
Relative expression	CD163 Recombinant Rabbit mAb [PSH18-98]	HA724034	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	CD163 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS276RT	Human	mIHC
	CD163 Recombinant Antibody [JA51-30] - Rat IgG1 (Chimeric)	HA601549	Human	IHC-P
	CD163 Mouse mAb [A3B3]	EM1901-90	Human	WB,IHC-P,FC
	CD163 Rabbit pAb	ER1804-03	Human, Mouse, Rat	WB,IHC-P
Relative expression	CD163 Recombinant Rabbit mAb [PSH15-06]	HA723706	Human	IHC-P,WB
	CD163 Mouse mAb [A3B5]	EM1901-91	Human	IHC-P,FC
	CD163 Recombinant Rabbit mAb	IRS002	Human	mIHC
	CD163 Recombinant Antibody [JA51-30] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610364	Human	IHC-P
Cell treatment	CD204 Recombinant Mouse mAb [A8G2-R]	HA601165	Human, Rat	WB,IHC-P,IF-Cell
Cell treatment	CD204 Mouse mAb [A8G2]	HA601090	Human	WB,IHC-P
Cell treatment	CD204 Recombinant Mouse mAb [A8G2-R] - BSA and Azide free	HA610052	Human, Rat	WB,IHC-P,IF-Cell
	CD68 Recombinant Antibody - Rabbit IgG (Chimeric)	IRS266RB	Human	mIHC
Relative expression	CD68 Mouse mAb [A3C4]	EM1901-95	Human	WB,IHC-P,FC,mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	CD68 Recombinant Rabbit mAb [PSH05-47] - BSA and Azide free	HA751006	Mouse, Rat	WB,IF-Cell,IHC-P,IHC-Fr,IF-Tissue
	CD68 Recombinant Rabbit mAb [SN07-27]	ET1611-53	Human	WB
Relative expression	CD68 Recombinant Mouse mAb [A3C3-R]	HA601290	Human	WB,IF-Cell,IHC-P,FC
	CD68 Mouse mAb [F3-D7]	EM1706-11	Human	WB,IHC-P,IF-Cell,FC
	CD68 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS264RT	Human	mIHC
Relative expression	CD68 Mouse mAb [A3C2]	EM1901-93	Human	WB,IHC-P,FC
Relative expression, Cell treatment	CD68 Recombinant Rabbit mAb [PSH06-65]	HA722737	Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,IHC-Fr
	CD68 Recombinant Antibody [PDM0-13] - Rat IgG1 (Chimeric)	HA601520	Human	IF-Tissue,IHC-P
Cell treatment	CD68 Recombinant Rabbit mAb [PSH05-47]	HA722285	Mouse, Rat	WB,IF-Cell,IHC-P,IHC-Fr,IF-Tissue
	CD68 Recombinant Antibody [PDM0-13] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610347	Human	IF-Tissue,IHC-P
Relative expression	CD68 Mouse mAb [A3C3]	EM1901-94	Human	WB,IHC-P,FC
	CD68 Rabbit pAb	ER1901-32	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	CD68 Recombinant Mouse mAb	IRS004	Human	mIHC
Relative expression	CD68 Recombinant Mouse mAb [A3C4-R]	HA601292	Human	WB,IF-Cell,IHC-P
Relative expression	CD68 Recombinant Mouse mAb [A3C2-R]	HA601291	Human	WB,IF-Cell,IHC-P
	CD68 Recombinant Antibody [PDM0-13] - Rabbit IgG (Chimeric) - BSA and Azide free	HA751544	Human	IF-Tissue,IHC-P
Relative expression	CD68 Recombinant Mouse mAb [PDM0-13] - BSA and Azide free	HA610034	Human	IHC-P,IF-Cell,FC,IF-Tissue
Relative expression	CD68 Recombinant Mouse mAb [A3C4-R] - BSA and Azide free	HA610165	Human	WB,IF-Cell,IHC-P
Relative expression, Cell treatment	CD68 Recombinant Rabbit mAb [PSH06-65] - BSA and Azide free	HA751094	Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,IHC-Fr
Relative expression	CD68 Recombinant Mouse mAb [PDM0-13]	HA601115	Human	IHC-P,IF-Cell,FC,mIHC,IF-Tissue
Relative expression	CD68 Recombinant Mouse mAb [A3C2-R] - BSA and Azide free	HA610164	Human	WB,IF-Cell,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CD68 Recombinant Antibody [PDM0-13] - Rabbit IgG (Chimeric)	HA723708	Human	IF-Tissue,IHC-P
	CD74 Recombinant Rabbit mAb [JF40-10]	ET1702-51	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
Relative expression	CD74 Recombinant Mouse mAb [PDM0-22]	HA601117	Human	WB,IHC-P,IF-Cell,IF-Tissue
	F4/80 Recombinant Rabbit mAb	IRS052	Mouse	mIHC
	F4/80 Recombinant Antibody [PSH0-85] - Rat IgG1 (Chimeric)	HA601485	Mouse, Rat	IHC-P
Relative expression	F4/80 Recombinant Antibody [PSH01-87] - Mouse IgG1 (Chimeric)	HA601537	Mouse	WB,IHC-P
Relative expression	F4/80 Recombinant Antibody [PSH01-87] - Guinea pig IgG2 (Chimeric)	HA601567	Mouse	WB,IHC-Fr
	F4/80 Recombinant Antibody [PSH01-87] - Rat IgG1 (Chimeric)	HA601538	Mouse	WB,IHC-Fr,IHC-P
	F4/80 Recombinant Antibody [PSH0-85] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610319	Mouse, Rat	IHC-P
	F4/80 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS143RT	Mouse	mIHC
Relative expression	F4/80 Recombinant Rabbit mAb [PSH0-85]	HA721520	Mouse, Rat	WB,IHC-P,IF-Tissue,IHC-Fr,mIHC
Relative expression	F4/80 Recombinant Rabbit mAb [PSH0-85] - BSA and Azide free	HA750665	Mouse, Rat	WB,IHC-P,IF-Tissue,IHC-Fr
Relative expression	F4/80 Recombinant Antibody [PSH01-87] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610356	Mouse	WB,IHC-P
	F4/80 Recombinant Antibody [PSH01-87] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610357	Mouse	WB,IHC-Fr,IHC-P
Relative expression	F4/80 Recombinant Rabbit mAb [PSH01-87]	HA721745	Mouse	WB,IHC-P,IF-Tissue,IHC-Fr
Relative expression	F4/80 Recombinant Rabbit mAb [PSH01-87] - BSA and Azide free	HA750761	Mouse	WB,IHC-P,IF-Tissue,IHC-Fr
	GPNUMB Recombinant Rabbit mAb [PSH0-82] - BSA and Azide free	HA750662	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	GPNUMB Recombinant Rabbit mAb [PSH0-82]	HA721511	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	GPNUMB Mouse mAb [E5-C10]	EM1708-60	Human	WB,IHC-P
	Histamine H2 receptor Rabbit pAb	HA500154	Human	WB,IHC-P
	HLA-DR Recombinant Rabbit mAb [SC06-78]	ET1610-66	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	HLA-DR Recombinant Rabbit mAb [SC06-78] - BSA and Azide free	HA750230	Human, Mouse, Rat	WB,IF-Cell,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	HLA-DR Recombinant Mouse mAb [10-D8-R] - BSA and Azide free	HA610127	Human	WB,IF-Cell,IHC-P,IF-Tissue
	HLA-DR Recombinant Mouse mAb [10-D8-R]	HA601252	Human	WB,IF-Cell,IHC-P,IF-Tissue
	HLA-DR Recombinant Mouse mAb	IRS154MS	Human	mIHC
	HLA-DR Mouse mAb [10-D8]	M1701-5	Human	WB,IHC-P,IF-Cell
	HLA-DR Rabbit pAb	R1510-35	Human	WB,IF-Cell,IHC-P,FC
	HLA-DR Rabbit pAb	R1510-34	Human, Mouse	WB,IF-Cell,IHC-P,FC
	HLA-DRB4 Recombinant Rabbit mAb [JE55-85]	ET7111-17	Human	WB,IHC-P,FC
	FITC Conjugated Human CD11b Recombinant Mouse mAb [A2-9-R]	HA600130F	Human	FC
	iFluor™ 647 Conjugated Human CD14 Recombinant Rabbit mAb [SC69-02]	HA720184F	Human	FC
	FITC Conjugated Human CD14 Recombinant Rabbit mAb [SC69-02]	HA720194F	Human	FC
	PE Conjugated Human CD16 Recombinant Mouse mAb [PSH04-95]	HA600129F	Human	FC
	Human CD16 Recombinant Mouse mAb [PSH04-95]	HA601340	Human	FC
	Human CD163 Recombinant Rabbit mAb [PSH13-78] - BSA and Azide free (Capture)	HA723566	Human	ELISA(Cap)
	Human CD163 Recombinant Rabbit mAb [PSH13-80] - BSA and Azide free (Capture)	HA723569	Human	ELISA(Cap)
	Human CD163 Recombinant Rabbit mAb [PSH13-79] - BSA and Azide free (Detector)	HA723567	Human	ELISA(Det)
	Biotin Conjugated Human CD163 Recombinant Rabbit mAb [PSH13-79] - Detector	HA723568B	Human	ELISA(Det),ELISA
	FITC Conjugated Human CD68 Recombinant Mouse mAb [A3C3-R]	HA600132F	Human	FC
	Biotin Conjugated Human IDO1 Recombinant Rabbit mAb [PSH09-87] - Detector	HA723156B	Human	ELISA(Det),ELISA
	Human IDO1 Recombinant Rabbit mAb [PSH09-87] - BSA and Azide free (Detector)	HA723155	Human	ELISA(Det)
	Human IDO1 Recombinant Rabbit mAb [PSH09-86] - BSA and Azide free (Capture)	HA723154	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human MERTK Recombinant Rabbit mAb [PSH12-69] - BSA and Azide free (Capture)	HA723479	Human	ELISA(Cap)
	HRP Conjugated Human MERTK Recombinant Rabbit mAb [PSH12-70] - Detector	HA723481H	Human	ELISA(Det),ELISA
	Human MERTK Recombinant Rabbit mAb [PSH12-70] - BSA and Azide free (Detector)	HA723480	Human	ELISA(Det)
	Biotin Conjugated Human MERTK Recombinant Rabbit mAb [PSH12-70] - Detector	HA723481B	Human	ELISA(Det),ELISA
	Human Osteopontin Recombinant Rabbit mAb [PSH08-58] - BSA and Azide free (Capture)	HA723004	Human	ELISA(Cap)
	Biotin Conjugated Human Osteopontin Recombinant Rabbit mAb [PSH08-60] - Detector	HA723008B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human Osteopontin Recombinant Rabbit mAb [PSH08-59] - Detector	HA723007B	Human	ELISA(Det),ELISA
	Human Osteopontin Recombinant Rabbit mAb [PSH08-60] - BSA and Azide free (Detector)	HA723006	Human	ELISA(Det)
	Human Osteopontin Recombinant Rabbit mAb [PSH08-59] - BSA and Azide free (Detector)	HA723005	Human	ELISA(Det)
Cell treatment, Relative expression	IDO Recombinant Rabbit mAb [PSH12-31]	HA723456	Human	WB,IF-Cell,IHC-P,FC,IP
	IDO Recombinant Rabbit mAb	IRS045	Human	mIHC
	IDO Rabbit pAb	ER1706-94	Human, Mouse	WB,IF-Cell,IHC-P,FC
Cell treatment	IDO Recombinant Rabbit mAb [PD00-62]	HA721331	Human	IHC-P,IF-Cell,WB,mIHC
	IDO Rat mAb [1G1]	RT1318	Human, Mouse	WB,IP,IF,IHC-P
Cell treatment, Relative expression	IDO Recombinant Rabbit mAb [PSH12-31] - BSA and Azide free	HA751445	Human	WB,IF-Cell,IHC-P,FC,IP
Cell treatment	iNOS Rabbit pAb	ER1706-89	Human, Mouse, Rat, Chicken	IF-Cell,FC,WB,IHC-Fr
Cell treatment	iNOS Recombinant Rabbit mAb [PSH03-69]	HA722031	Human, Mouse	WB,IHC-P,IF-Cell
Cell treatment	iNOS Recombinant Rabbit mAb [JE59-61]	HA721583	Human, Mouse	WB,IF-Cell
Cell treatment	iNOS Recombinant Rabbit mAb [PSH03-69] - BSA and Azide free	HA750903	Human, Mouse	WB,IHC-P,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	Liver Arginase Recombinant Rabbit mAb [PSH09-78]	HA723143	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Liver Arginase Recombinant Rabbit mAb [PD00-05]	HA721147	Human	WB,IHC-P,IF-Tissue
Cell treatment	Liver Arginase Recombinant Rabbit mAb [SY09-06]	ET1605-8	Human, Mouse, Rat	WB
Relative expression	Liver Arginase Recombinant Rabbit mAb [JB21-49]	ET7107-90	Human	WB,IHC-P,IP
Relative expression	Ly6c Recombinant Rabbit mAb [PSH04-51]	HA722134	Mouse	WB
Relative expression	Ly6c Rabbit pAb	HA500088	Mouse	WB,IHC-P,FC
Relative expression	Ly6c Recombinant Rabbit mAb [PSH04-51] - BSA and Azide free	HA750958	Mouse	WB
	Ly6c Rabbit pAb	HA500087	Mouse	WB,IHC-P,FC
Cell treatment	Mannose Receptor(CD206) Recombinant Rabbit mAb [PSH07-76]	HA722892	Mouse, Rat	WB,IHC-P,IF-Tissue,IHC-Fr,IF-Cell,IP,mIHC
	Mannose Receptor(CD206) Recombinant Antibody [PSH07-76] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610287	Mouse, Rat	IHC-Fr,IF-Tissue,IHC-P,WB,IF-Cell
	Mannose Receptor(CD206) Recombinant Antibody [PSH07-76] - Guinea pig IgG2 (Chimeric)	HA601452	Mouse, Rat	IHC-Fr,WB
	Mannose Receptor(CD206) Recombinant Antibody [PSH07-76] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610286	Mouse, Rat	IHC-Fr,WB
Cell treatment	Mannose Receptor(CD206) Recombinant Rabbit mAb [PSH07-76] - BSA and Azide free	HA751188	Mouse, Rat	WB,IHC-P,IF-Tissue,IHC-Fr,IF-Cell,IP
KD	Mannose Receptor(CD206) Recombinant Rabbit mAb [JF0953]	ET1702-04	Human, Mouse	WB,IF-Cell,FC
KD	Mannose Receptor(CD206) Recombinant Rabbit mAb [JF0953] - BSA and Azide free	HA750332	Human, Mouse	WB,IF-Cell,FC
	Mannose Receptor(CD206) Recombinant Antibody [PSH07-76] - Mouse IgG1 (Chimeric)	HA601453	Mouse, Rat	IHC-Fr,IF-Tissue,IHC-P,WB,IF-Cell
	FITC Conjugated Mannose Receptor(CD206) Recombinant Rabbit mAb [JF0953]	HA720214F	Human	FC
	Mannose Receptor(CD206) Recombinant Rabbit mAb	IRS049	Mouse	mIHC
	MARCO Rabbit pAb	ER60815	Human, Mouse	WB
Cell treatment	MARCO Recombinant Rabbit mAb [PSH04-31] - BSA and Azide free	HA750940	Human, Mouse	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	MARCO Recombinant Rabbit mAb [PSH04-31]	HA722114	Human, Mouse	WB,IHC-P
KD	MERTK Recombinant Rabbit mAb [SR29-07] - BSA and Azide free	HA750036	Human	WB,IHC-P,IP,FC,IF-Tissue
KD	MERTK Recombinant Rabbit mAb [SR29-07]	ET1602-21	Human	WB,IHC-P,FC,IF-Tissue
	PE Conjugated Mouse CD16/32 Recombinant Rabbit mAb [PSH06-60]	HA720216F	Mouse	FC
	Biotin Conjugated Mouse CD68 Recombinant Rabbit mAb [PSH12-27] - Detector	HA723451B	Mouse	ELISA(Det),ELISA
	Mouse CD68 Recombinant Rabbit mAb [PSH12-27] - BSA and Azide free (Detector)	HA723450	Mouse	ELISA(Det)
	Mouse CD68 Recombinant Rabbit mAb [PSH12-25] - BSA and Azide free (Detector)	HA723447	Mouse	ELISA(Det)
	Biotin Conjugated Mouse CD68 Recombinant Rabbit mAb [PSH12-25] - Detector	HA723448B	Mouse	ELISA(Det),ELISA
	Mouse CD68 Recombinant Rabbit mAb [PSH12-26] - BSA and Azide free (Capture)	HA723449	Mouse	ELISA(Cap)
	Mouse CD68 Recombinant Rabbit mAb [PSH12-24] - BSA and Azide free (Capture)	HA723446	Mouse	ELISA(Cap)
	Mouse Osteopontin Recombinant Rabbit mAb [PSH17-50] - BSA and Azide free (Detector)	HA723925	Mouse	ELISA(Det)
	Biotin Conjugated Mouse Osteopontin Recombinant Rabbit mAb [PSH17-50] - Detector	HA723925B	Mouse	ELISA(Det),ELISA
	Mouse Osteopontin Recombinant Rabbit mAb [PSH17-49] - BSA and Azide free (Capture)	HA723924	Mouse	ELISA(Cap)
	HRP Conjugated Mouse Osteopontin Recombinant Rabbit mAb [PSH17-50] - Detector	HA723925H	Mouse	ELISA(Det),ELISA
	MSR1 Rabbit pAb	ER1913-21	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
	MSR1 / CD204 Rabbit pAb	ER1902-71	Human	WB,FC
	Osteopontin Rabbit pAb	0806-8	Human, Zebrafish	WB,IHC-P
Cell treatment	Osteopontin Rabbit pAb	ER1802-16	Human, Mouse	IF-Cell,IHC-P,ELISA
	Osteopontin Rabbit pAb	HA500464	Human, Mouse	WB,FC
	Osteopontin Rabbit pAb	ER1802-15	Human, Mouse	IHC-P,ELISA
Cell treatment	Osteopontin Rabbit pAb	0806-6	Mouse, Human, Rat	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Osteopontin Recombinant Rabbit mAb [PSH09-23] - BSA and Azide free	HA751277	Human	IHC-Fr,IHC-P,IF-Tissue
Cell treatment	Osteopontin Recombinant Rabbit mAb [PSH09-24] - BSA and Azide free	HA751278	Human, Mouse	WB,IF-Cell,FC
Cell treatment	Osteopontin Recombinant Rabbit mAb [PSH09-24]	HA723082	Human, Mouse	WB,IF-Cell,FC
	Osteopontin Recombinant Rabbit mAb [PSH09-23]	HA723081	Human	IHC-Fr,IHC-P,IF-Tissue
	Osteopontin Recombinant Rabbit mAb [PSH14-11]	HA723601	Mouse, Rat	IHC-P,IHC-Fr,WB,IF-Cell
	Osteopontin Recombinant Rabbit mAb	IRS200RB	Human	mIHC
	Osteopontin Recombinant Rabbit mAb [PSH14-11] - BSA and Azide free	HA751506	Mouse, Rat	IHC-P,IHC-Fr,WB,IF-Cell
Cell treatment	Phospho-Axl (Y702) Recombinant Rabbit mAb [PSH17-98]	HA723963	Human	WB
Relative expression, KD, Cell treatment	VISTA Recombinant Rabbit mAb [PSH0-64]	HA721452	Human	WB,IHC-P,IF-Tissue
	iFluor™ 488 Conjugated VISTA Recombinant Rabbit mAb [PSH0-65]	HA720183F	Human	FC
	VISTA Recombinant Rabbit mAb [PSH0-65] - BSA and Azide free	HA750642	Human	WB,IHC-P,IF-Cell,FC,IF-Tissue
Relative expression, KD, Cell treatment	VISTA Recombinant Rabbit mAb [PSH0-64] - BSA and Azide free	HA750641	Human	WB,IHC-P,IF-Tissue
	VISTA Recombinant Rabbit mAb [PSH0-65]	HA721453	Human	WB,IHC-P,IF-Cell,FC,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层