

巨噬细胞极化研究-产品清单

高级验证	产品名	货号	反应种属	验证应用
	Arginase1 Recombinant Rabbit mAb	IRS081RB	Human	mIHC
	Arginase1 Recombinant Rabbit mAb	IRS082RB	Mouse	mIHC
	B7-1/CD80 Rabbit pAb	ER1904-08	Mouse, Rat	WB
KD	CD14 Recombinant Rabbit mAb [PD01-63]	HA721321	Human, Mouse	WB,IHC-P,IF-Tissue
	CD14 Rabbit pAb	R1107-2	Human, Mouse	WB,IF-Cell,FC
KD	CD14 Recombinant Antibody [SC69-02] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610260	0	
KD	CD14 Recombinant Rabbit mAb [SC69-02] - BSA and Azide free	HA750232	Human, Mouse	WB,IF-Cell,IHC-P,FC
	CD14 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS106RT	Human	mIHC
KD	CD14 Recombinant Rabbit mAb [SC69-02]	ET1610-85	Human, Mouse	WB,IF-Cell,IHC-P,mIHC,FC
KD	CD14 Recombinant Antibody [SC69-02] - Rat IgG1 (Chimeric)	HA601420		
	CD14 Recombinant Rabbit mAb	IRS015	Human	mIHC
	CD16 Recombinant Rabbit mAb [JE49-79]	ET7109-97	Human, Mouse, Rat	WB
Relative expression	CD16 Recombinant Rabbit mAb [PD00-12]	HA721149	Human	IHC-P,IF-Tissue,IF-Cell,mIHC,WB
	CD16 Mouse mAb [Q32]	EM60067	Human	WB,IHC-P
	CD16 Recombinant Rabbit mAb	IRS017	Human	mIHC
	CD163 Mouse mAb [A3B3]	EM1901-90	Human	WB,IHC-P,FC
	CD163 Recombinant Rabbit mAb [JA51-30]	ET1704-43	Human	WB,IHC-P,mIHC,IF-Tissue,IP
	CD163 Mouse mAb [A3B5]	EM1901-91	Human	IHC-P,FC
	CD163 Rabbit pAb	ER1804-03	Human, Mouse, Rat	WB,IHC-P
	CD163 Recombinant Rabbit mAb	IRS002	Human	mIHC
	CD163 Recombinant Rabbit mAb [JA51-30] - BSA and Azide free	HA750412	Human	WB,IHC-P,IF-Tissue,IP
Relative expression	CD163 Recombinant Rabbit mAb [PSH15-06]	HA723706	Human	IHC-P,WB

Relative expression	CD163 Recombinant Rabbit mAb [PSH18-98]	HA724034	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	CD163 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS276RT	Human	mIHC
	CD163 Recombinant Antibody [JA51-30] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610364	Human	IHC-P
	CD163 Recombinant Antibody [JA51-30] - Rat IgG1 (Chimeric)	HA601549	Human	IHC-P
	CD163 Rabbit pAb	IRS109RB	Mouse	mIHC
Cell treatment	CD204 Recombinant Mouse mAb [A8G2-R] - BSA and Azide free	HA610052	Human, Rat	WB,IHC-P,IF-Cell
Cell treatment	CD204 Mouse mAb [A8G2]	HA601090	Human	WB,IHC-P
Cell treatment	CD204 Recombinant Mouse mAb [A8G2-R]	HA601165	Human, Rat	WB,IHC-P,IF-Cell
	CD32 Rabbit pAb	ER62945	Human, Mouse, Rat	WB,IHC-P
Cell treatment	CD32 Recombinant Rabbit mAb [JM10-70]	ET1703-43	Human	WB,IHC-P,IF-Cell,FC,IP
Cell treatment	CD32 Recombinant Rabbit mAb [JM10-70] - BSA and Azide free	HA750382	Human	WB,IHC-P,IF-Cell,FC,IP
	CD32A + CD32B + CD32C Recombinant Rabbit mAb [JE64-17]	HA721930	Human	WB,IHC-P
	CD32B Mouse mAb [A9C11]	HA601121	Human	WB,IHC-P
	CD32-C Rabbit pAb	ER62946	Human	WB
	CD64 Mouse mAb [A7E4]	HA601003	Human	WB
	CD64 Mouse mAb [5F5F8]	EM1710-99	Human	WB,IHC-P,FC
Relative expression	CD68 Mouse mAb [A3C3]	EM1901-94	Human	WB,IHC-P,FC
Relative expression	CD68 Mouse mAb [A3C2]	EM1901-93	Human	WB,IHC-P,FC
Relative expression, Cell treatment	CD68 Recombinant Rabbit mAb [PSH06-65] - BSA and Azide free	HA751094	Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,IHC-Fr
	CD68 Recombinant Antibody [PDM0-13] - Rabbit IgG (Chimeric) - BSA and Azide free	HA751544	Human	IF-Tissue,IHC-P
Relative expression	CD68 Recombinant Mouse mAb [A3C3-R]	HA601290	Human	WB,IF-Cell,IHC-P,FC
Cell treatment	CD68 Recombinant Rabbit mAb [PSH05-47] - BSA and Azide free	HA751006	Mouse, Rat	WB,IF-Cell,IHC-P,IHC-Fr,IF-Tissue
	CD68 Recombinant Antibody [PDM0-13] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610347	Human	IF-Tissue,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	CD68 Recombinant Mouse mAb [A3C4-R] - BSA and Azide free	HA610165	Human	WB,IF-Cell,IHC-P
Relative expression	CD68 Recombinant Mouse mAb [A3C4-R]	HA601292	Human	WB,IF-Cell,IHC-P
Relative expression	CD68 Recombinant Mouse mAb [A3C2-R]	HA601291	Human	WB,IF-Cell,IHC-P
	CD68 Recombinant Antibody - Rabbit IgG (Chimeric)	IRS266RB	Human	mIHC
Relative expression	CD68 Recombinant Mouse mAb [PDM0-13] - BSA and Azide free	HA610034	Human	IHC-P,IF-Cell,FC,IF-Tissue
	CD68 Recombinant Antibody [PDM0-13] - Rat IgG1 (Chimeric)	HA601520	Human	IF-Tissue,IHC-P
	CD68 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS264RT	Human	mIHC
Relative expression	CD68 Mouse mAb [A3C4]	EM1901-95	Human	WB,IHC-P,FC,mIHC
	CD68 Mouse mAb [F3-D7]	EM1706-11	Human	WB,IHC-P,IF-Cell,FC
Relative expression	CD68 Recombinant Mouse mAb [PDM0-13]	HA601115	Human	IHC-P,IF-Cell,FC,mIHC,IF-Tissue
	CD68 Recombinant Antibody [PDM0-13] - Rabbit IgG (Chimeric)	HA723708	Human	IF-Tissue,IHC-P
	CD68 Recombinant Mouse mAb	IRS004	Human	mIHC
	CD68 Rabbit pAb	ER1901-32	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Cell treatment	CD68 Recombinant Rabbit mAb [PSH05-47]	HA722285	Mouse, Rat	WB,IF-Cell,IHC-P,IHC-Fr,IF-Tissue
	CD68 Recombinant Rabbit mAb [SN07-27]	ET1611-53	Human	WB
Relative expression, Cell treatment	CD68 Recombinant Rabbit mAb [PSH06-65]	HA722737	Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,IHC-Fr
Relative expression	CD68 Recombinant Mouse mAb [A3C2-R] - BSA and Azide free	HA610164	Human	WB,IF-Cell,IHC-P
	CD80 Recombinant Mouse mAb [8-E5-R]	HA601247	Human	WB
Relative expression, Cell treatment	CD80 Recombinant Rabbit mAb [PSH18-13]	HA723976	Mouse	WB,IF-Cell,FC,IP
Relative expression	CD80 Mouse mAb [8-E5]	M1007-10	Human, Mouse	WB,IHC-P,FC
	CD80 Rabbit pAb	1007-8	Human	WB,IF-Cell
Relative expression	CD80 Recombinant Rabbit mAb [PSH09-79]	HA723144	Human	WB,IF-Cell,FC
	CD80 Rabbit pAb	ER1905-99	Human, Mouse, Rat	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CD80 Recombinant Mouse mAb [8-E5-R] - BSA and Azide free	HA610122	Human	WB
Relative expression	CD80 Recombinant Rabbit mAb [PSH09-79] - BSA and Azide free	HA751321	Human	WB,IF-Cell,FC
Relative expression	CD80 Recombinant Rabbit mAb [JF100-4]	ET1702-95	Human	WB,IHC-P
	CD86 Recombinant Rabbit mAb [SJ20-00]	ET1606-50	Human, Mouse, Rat	WB,IF-Cell,IP
	CD86 Recombinant Rabbit mAb [SJ20-00] - BSA and Azide free	HA750109	Human, Mouse, Rat	WB,IF-Cell,IP
	CD86 Rabbit pAb	ER1906-01	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
	CSF-1-R Recombinant Rabbit mAb [PD00-45]	HA721180	Human	IHC-P,FC
Relative expression	CSF-1-R Recombinant Rabbit mAb [PSH09-77]	HA723142	Human, Mouse, Rat	WB,IHC-P
Relative expression	CSF-1-R Recombinant Rabbit mAb [PSH09-77] - BSA and Azide free	HA751319	Human, Mouse, Rat	WB,IHC-P
	CSF1R Mouse mAb [B9-B9]	EM1708-56	Human	WB,IF-Cell,IHC-P
	DC-SIGN Rabbit pAb	ER1907-23	Human	WB,IHC-P,IF-Tissue,FC
	DC-SIGN Recombinant Rabbit mAb [JE31-24]	HA721382	Human	WB
	FCG3B Rabbit pAb	ER62049	Human	WB
	FCGR1A/CD64 Rabbit pAb	ER1908-78	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue
	FITC Conjugated Human CD14 Recombinant Rabbit mAb [SC69-02]	HA720194F	Human	FC
	iFluor™ 647 Conjugated Human CD14 Recombinant Rabbit mAb [SC69-02]	HA720184F	Human	FC
	PE Conjugated Human CD16 Recombinant Mouse mAb [PSH04-95]	HA600129F	Human	FC
	Human CD16 Recombinant Mouse mAb [PSH04-95]	HA601340	Human	FC
	Human CD163 Recombinant Rabbit mAb [PSH13-80] - BSA and Azide free (Capture)	HA723569	Human	ELISA(Cap)
	Human CD163 Recombinant Rabbit mAb [PSH13-78] - BSA and Azide free (Capture)	HA723566	Human	ELISA(Cap)
	Human CD163 Recombinant Rabbit mAb [PSH13-79] - BSA and Azide free (Detector)	HA723567	Human	ELISA(Det)
	Biotin Conjugated Human CD163 Recombinant Rabbit mAb [PSH13-79] - Detector	HA723568B	Human	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	FITC Conjugated Human CD68 Recombinant Mouse mAb [A3C3-R]	HA600132F	Human	FC
	Biotin Conjugated Human CD80 Recombinant Rabbit mAb [PSH09-63] - Detector	HA723127B	Human	ELISA(Det),ELISA
	Human CD80 Recombinant Rabbit mAb [PSH09-63] - BSA and Azide free (Detector)	HA723125	Human	ELISA(Det)
	Human CD80 Recombinant Rabbit mAb [PSH09-62] - BSA and Azide free (Detector)	HA723124	Human	ELISA(Det)
	Biotin Conjugated Human CD80 Recombinant Rabbit mAb [PSH09-62] - Detector	HA723126B	Human	ELISA(Det),ELISA
	Human CD80 Recombinant Rabbit mAb [PSH09-61] - BSA and Azide free (Capture)	HA723123	Human	ELISA(Cap)
	Human CSF-1-R Recombinant Rabbit mAb [PSH12-67] - BSA and Azide free (Detector)	HA723475	Human	ELISA(Det)
	Human CSF-1-R Recombinant Rabbit mAb [PSH12-66] - BSA and Azide free (Capture/Detector)	HA723474	Human	ELISA(Cap),ELISA(Det)
	Biotin Conjugated Human CSF-1-R Recombinant Rabbit mAb [PSH12-66] - Detector	HA723478B	Human	ELISA(Det),ELISA
	Human CSF-1-R Recombinant Rabbit mAb [PSH12-68] - BSA and Azide free (Capture)	HA723477	Human	ELISA(Cap)
	HRP Conjugated Human CSF-1-R Recombinant Rabbit mAb [PSH12-66] - Detector	HA723478H	Human	ELISA(Det),ELISA
	Biotin Conjugated Human CSF-1-R Recombinant Rabbit mAb [PSH12-67] - Detector	HA723476B	Human	ELISA(Det),ELISA
	Human IDO1 Recombinant Rabbit mAb [PSH09-87] - BSA and Azide free (Detector)	HA723155	Human	ELISA(Det)
	Human IDO1 Recombinant Rabbit mAb [PSH09-86] - BSA and Azide free (Capture)	HA723154	Human	ELISA(Cap)
	Biotin Conjugated Human IDO1 Recombinant Rabbit mAb [PSH09-87] - Detector	HA723156B	Human	ELISA(Det),ELISA
	Human IL-10 Recombinant Rabbit mAb [PS00-81] - BSA and Azide free (Detector)	HA721237	Human	ELISA(Det)
	Biotin Conjugated Human IL-10 Recombinant Rabbit mAb [PS00-81]	HA722652B	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human IL-10 Recombinant Rabbit mAb [PS00-98] - BSA and Azide free (Capture)	HA721236	Human	ELISA(Cap)
	Biotin Conjugated Human IL-12B p40 Recombinant Rabbit mAb [PSH02-58]	HA722633B	Human	ELISA(Det)
	Human IL-12B p40 Recombinant Rabbit mAb [PSH02-57] - BSA and Azide free (Capture)	HA721834	Human	ELISA(Cap)
	Biotin Conjugated Human IL-12B p40 Recombinant Rabbit mAb [PSH03-04]	HA722629B	Human	ELISA(Det)
	Human IL-12B p40 Recombinant Rabbit mAb [PSH03-04] - BSA and Azide free (Detector)	HA721934	Human	ELISA(Det)
	Human IL-12B p40 Recombinant Rabbit mAb [PSH02-58] - BSA and Azide free (Detector)	HA721835	Human	ELISA(Det)
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS00-82] - BSA and Azide free	HA750745	Human	WB
	Human IL-1 beta Recombinant Rabbit mAb [PS01-39] - BSA and Azide free (Capture)	HA721278	Human	ELISA(Cap)
	Biotin Conjugated Human IL-1 beta Recombinant Rabbit mAb [PS00-82]	HA722649B	Human	ELISA(Det)
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS01-39]	HA721712	Human	WB
	Human IL-1 beta Recombinant Rabbit mAb [PS00-82] - BSA and Azide free (Detector)	HA721279	Human	ELISA(Det)
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS00-82]	HA721713	Human	WB
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-85] - BSA and Azide free (Capture)	HA723881	Human	ELISA(Cap)
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-86] - BSA and Azide free (Detector)	HA723882	Human	ELISA(Det)
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS01-39] - BSA and Azide free	HA750744	Human	WB
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-87] - BSA and Azide free (Capture)	HA723883	Human	ELISA(Cap)
	Biotin Conjugated Human IL-1 beta Recombinant Rabbit mAb [PSH16-86] - Detector	HA723882B	Human	ELISA(Det),ELISA
	Human IL-23A Recombinant Rabbit mAb [PSH03-03] - BSA and Azide free (Capture)	HA721933	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human IL-6 Mouse mAb [A6D2]	HA601050	Human	WB,ELISA
	Human IL-6 Recombinant Rabbit mAb [PS00-21] - BSA and Azide free (Capture)	HA721167	Human	ELISA(Cap)
	HRP Conjugated Human IL-6 Recombinant Rabbit mAb [PS00-15] - Detector	HA721293	Human	ELISA(Det)
	Human IL-6 Mouse mAb [A6D1]	HA601051	Human	WB,ELISA
	Human IL-6 Recombinant Rabbit mAb [PS00-15] - BSA and Azide free (Detector)	HA721168	Human	ELISA(Det)
	Biotin Conjugated Human IL-6 Recombinant Rabbit mAb [PS00-15] - Detector	HA722233B	Human	ELISA(Det)
	Human Interferon gamma Recombinant Rabbit mAb [PSH01-45] - BSA and Azide free (Capture)	HA721683	Human	ELISA(Cap)
	Human Interferon gamma Recombinant Rabbit mAb [PS01-75] - BSA and Azide free (Capture)	HA721574	Human	ELISA(Cap)
	Human Interferon gamma Recombinant Rabbit mAb [PSH01-47] - BSA and Azide free (Detector)	HA721685	Human	ELISA(Det)
	Human Interferon gamma Recombinant Rabbit mAb [PSH01-46] - BSA and Azide free (Capture)	HA721684	Human	ELISA(Cap)
	Human Interferon gamma Recombinant Rabbit mAb [PS01-76] - BSA and Azide free (Detector)	HA721575	Human	ELISA(Det)
	Biotin Conjugated Human Interferon gamma Recombinant Rabbit mAb [PS01-76] - Detector	HA721576	Human	ELISA(Det),ELISA
	Biotin Conjugated Human Interferon gamma Recombinant Rabbit mAb [PSH01-47]	HA722645B	Human	ELISA(Det)
	Biotin Conjugated Human MERTK Recombinant Rabbit mAb [PSH12-70] - Detector	HA723481B	Human	ELISA(Det),ELISA
	Human MERTK Recombinant Rabbit mAb [PSH12-70] - BSA and Azide free (Detector)	HA723480	Human	ELISA(Det)
	Human MERTK Recombinant Rabbit mAb [PSH12-69] - BSA and Azide free (Capture)	HA723479	Human	ELISA(Cap)
	HRP Conjugated Human MERTK Recombinant Rabbit mAb [PSH12-70] - Detector	HA723481H	Human	ELISA(Det),ELISA
	Human TNF alpha Recombinant Rabbit mAb [PS00-40]	HA721286	Human	IHC-P,mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human TNF alpha Recombinant Rabbit mAb [PS00-40] - BSA and Azide free (Capture)	HA721265	Human	ELISA(Cap)
	HRP Conjugated Human TNF alpha Recombinant Rabbit mAb [PS01-07] - Detector	HA721294	Human	ELISA(Det)
	Biotin Conjugated Human TNF alpha Recombinant Rabbit mAb [PS01-07] - Detector	HA721559	Human	ELISA(Det)
	Human TNF alpha Recombinant Rabbit mAb [PS00-40] - BSA and Azide free	HA750593	Human	IHC-P
	Human TNF alpha Recombinant Rabbit mAb	IRS040	Human	mIHC
	Human TNF alpha Recombinant Rabbit mAb [PS01-07] - BSA and Azide free (Detector)	HA721266	Human	ELISA(Det)
	IDO Recombinant Rabbit mAb	IRS045	Human	mIHC
	IDO Rabbit pAb	ER1706-94	Human, Mouse	WB,IF-Cell,IHC-P,FC
	IDO Rat mAb [1G1]	RT1318	Human, Mouse	WB,IP,IF,IHC-P
Cell treatment	IDO Recombinant Rabbit mAb [PD00-62]	HA721331	Human	IHC-P,IF-Cell,WB,mIHC
Cell treatment, Relative expression	IDO Recombinant Rabbit mAb [PSH12-31]	HA723456	Human	WB,IF-Cell,IHC-P,FC,IP
Cell treatment, Relative expression	IDO Recombinant Rabbit mAb [PSH12-31] - BSA and Azide free	HA751445	Human	WB,IF-Cell,IHC-P,FC,IP
Cell treatment	IL-10 Recombinant Rabbit mAb [JE59-71]	HA722032	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	IL-10 Rabbit pAb	ER1911-19	Human, Mouse	WB,IHC-P,IF-Tissue
	IL-10 Rabbit pAb	ER65464	Human, Mouse, Rat	WB,IHC-P
	IL12B Rabbit pAb	ER1803-31	Human, Mouse	WB
Cell treatment	IL-1 alpha Recombinant Rabbit mAb [PSH15-65] - BSA and Azide free	HA751576	Mouse	WB
	IL-1 alpha Recombinant Rabbit mAb [PSH15-66]	HA723773	Mouse	IF-Cell,FC
	IL-1 alpha Rabbit pAb	ER1802-36	Human, Mouse	IHC-P,FC,WB
	IL-1 alpha Recombinant Rabbit mAb [PSH15-66] - BSA and Azide free	HA751577	Mouse	IF-Cell,FC
Cell treatment	IL-1 alpha Recombinant Rabbit mAb [PSH15-65]	HA723772	Mouse	WB
Cell treatment	IL-1 beta Recombinant Rabbit mAb [JJ087-3]	ET1701-39	Human	WB,IF-Cell,IF-Tissue,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	IL-1 beta Recombinant Rabbit mAb [PSH18-00]	HA723965	Mouse, Rat	WB
Cell treatment	IL-1 beta Recombinant Mouse mAb [A7F7-R] - BSA and Azide free	HA610004	Human, Mouse, Rat	WB,ELISA,IHC-P,IF-Tissue
Cell treatment	IL-1 beta Mouse mAb [A7F9]	HA601002	Human, Mouse, Rat	WB,ELISA,IF-Cell
Cell treatment	IL-1 beta Recombinant Mouse mAb [A7F7-R]	HA601036	Human, Mouse, Rat	WB,ELISA,IHC-P,IF-Tissue
Cell treatment	IL-1 beta Mouse mAb [A7A10]	HA601035	Human	WB,ELISA,IF-Cell
	IL-23 Recombinant Rabbit mAb [JE38-99]	HA721037	Human	WB,IHC-P,IF-Tissue
	IL-6 Mouse mAb [3-A4]	M1303-3	Mouse	WB,IF-Cell,FC,ELISA
	IL-6 Rabbit pAb	R1412-2	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	IL-6 Rabbit pAb	R1203-2	Mouse	WB,ELISA
	IL-6 Recombinant Rabbit mAb	IRS160RB	Human	mIHC
	IL-6 Mouse mAb [1-6]	EM1701-45	Mouse, Rat	WB,IF-Cell
	IL-6 Mouse mAb	EM30301	Human	WB
	Interferon gamma Recombinant Rabbit mAb [JM10-10]	ET1703-17	Human	WB,FC
	Interferon gamma Recombinant Rabbit mAb [JM10-10] - BSA and Azide free	HA750373	Human	WB,FC
	IRF5 Recombinant Rabbit mAb [SN201-05] - BSA and Azide free	HA750277	Human	WB,IHC-P,FC
	IRF5 Recombinant Rabbit mAb [SD203-07]	ET1612-15	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,IP
	IRF5 Recombinant Rabbit mAb [SN201-05]	ET1611-94	Human	WB,IHC-P,FC
	Liver Arginase Recombinant Rabbit mAb [PD00-05]	HA721147	Human	WB,IHC-P,IF-Tissue
Relative expression	Liver Arginase Recombinant Rabbit mAb [JB21-49]	ET7107-90	Human	WB,IHC-P,IP
Cell treatment	Liver Arginase Recombinant Rabbit mAb [SY09-06]	ET1605-8	Human, Mouse, Rat	WB
Relative expression	Liver Arginase Recombinant Rabbit mAb [PSH09-78]	HA723143	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Relative expression	Ly6c Recombinant Rabbit mAb [PSH04-51]	HA722134	Mouse	WB
Relative expression	Ly6c Recombinant Rabbit mAb [PSH04-51] - BSA and Azide free	HA750958	Mouse	WB
Relative expression	Ly6c Rabbit pAb	HA500088	Mouse	WB,IHC-P,FC
	Ly6c Rabbit pAb	HA500087	Mouse	WB,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mannose Receptor(CD206) Recombinant Antibody [PSH07-76] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610287	Mouse, Rat	IHC-Fr,IF-Tissue,IHC- P,WB,IF-Cell
	Mannose Receptor(CD206) Recombinant Rabbit mAb	IRS049	Mouse	mIHC
Cell treatment	Mannose Receptor(CD206) Recombinant Rabbit mAb [PSH07-76]	HA722892	Mouse, Rat	WB,IHC-P,IF- Tissue,IHC-Fr,IF- Cell,IP,mIHC
	Mannose Receptor(CD206) Recombinant Antibody [PSH07-76] - Mouse IgG1 (Chimeric)	HA601453	Mouse, Rat	IHC-Fr,IF-Tissue,IHC- P,WB,IF-Cell
	FITC Conjugated Mannose Receptor(CD206) Recombinant Rabbit mAb [JF0953]	HA720214F	Human	FC
KD	Mannose Receptor(CD206) Recombinant Rabbit mAb [JF0953]	ET1702-04	Human, Mouse	WB,IF-Cell,FC
	Mannose Receptor(CD206) Recombinant Antibody [PSH07-76] - Guinea pig IgG2 (Chimeric)	HA601452	Mouse, Rat	IHC-Fr,WB
Cell treatment	Mannose Receptor(CD206) Recombinant Rabbit mAb [PSH07-76] - BSA and Azide free	HA751188	Mouse, Rat	WB,IHC-P,IF- Tissue,IHC-Fr,IF-Cell,IP
KD	Mannose Receptor(CD206) Recombinant Rabbit mAb [JF0953] - BSA and Azide free	HA750332	Human, Mouse	WB,IF-Cell,FC
	Mannose Receptor(CD206) Recombinant Antibody [PSH07-76] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610286	Mouse, Rat	IHC-Fr,WB
KD	MERTK Recombinant Rabbit mAb [SR29-07]	ET1602-21	Human	WB,IHC-P,FC,IF-Tissue
KD	MERTK Recombinant Rabbit mAb [SR29-07] - BSA and Azide free	HA750036	Human	WB,IHC-P,IP,FC,IF- Tissue
	MHC Class II Recombinant Rabbit mAb [JA10-94] - BSA and Azide free	HA750404	Human	WB,IP,IF-Cell,IF- Tissue,IHC-P
	MHC Class II Recombinant Rabbit mAb [JA10-94]	ET1704-13	Human	WB,IP,IF-Cell,IF- Tissue,IHC-P,mIHC
	PE Conjugated Mouse CD16/32 Recombinant Rabbit mAb [PSH06-60]	HA720216F	Mouse	FC
	Biotin Conjugated Mouse CD68 Recombinant Rabbit mAb [PSH12-27] - Detector	HA723451B	Mouse	ELISA(Det),ELISA
	Biotin Conjugated Mouse CD68 Recombinant Rabbit mAb [PSH12-25] - Detector	HA723448B	Mouse	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse CD68 Recombinant Rabbit mAb [PSH12-25] - BSA and Azide free (Detector)	HA723447	Mouse	ELISA(Det)
	Mouse CD68 Recombinant Rabbit mAb [PSH12-27] - BSA and Azide free (Detector)	HA723450	Mouse	ELISA(Det)
	Mouse CD68 Recombinant Rabbit mAb [PSH12-26] - BSA and Azide free (Capture)	HA723449	Mouse	ELISA(Cap)
	Mouse CD68 Recombinant Rabbit mAb [PSH12-24] - BSA and Azide free (Capture)	HA723446	Mouse	ELISA(Cap)
	Mouse CD86 Recombinant Rabbit mAb [PSH13-65] - BSA and Azide free (Capture)	HA723555	Mouse	ELISA(Cap)
	Mouse CD86 Recombinant Rabbit mAb [PSH13-66] - BSA and Azide free (Detector)	HA723556	Mouse	ELISA(Det)
	Mouse CD86 Recombinant Rabbit mAb [PSH13-64] - BSA and Azide free (Detector)	HA723553	Mouse	ELISA(Det)
	Biotin Conjugated Mouse CD86 Recombinant Rabbit mAb [PSH13-64] - Detector	HA723554B	Mouse	ELISA(Det),ELISA
	Mouse CD86 Recombinant Rabbit mAb [PSH13-63] - BSA and Azide free (Capture)	HA723552	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse CD86 Recombinant Rabbit mAb [PSH13-66] - Detector	HA723557B	Mouse	ELISA(Det),ELISA
	Mouse IL-10 Recombinant Rabbit mAb [PS01-85] - BSA and Azide free (Capture)	HA721517	Mouse	ELISA(Cap)
	Mouse IL-10 Recombinant Rabbit mAb [PS01-93] - BSA and Azide free (Detector)	HA721518	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-10 Recombinant Rabbit mAb [PS01-93] - Detector	HA721519	Mouse	ELISA(Det)
	Mouse IL12 p40 Recombinant Rabbit mAb [PSH04-87] - BSA and Azide free (Capture)	HA722193	Mouse	ELISA(Cap)
	Mouse IL12 p40 Recombinant Rabbit mAb [PSH04-88] - BSA and Azide free (Detector)	HA722194	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL12 p40 Recombinant Rabbit mAb [PSH04-88]	HA722603B	Mouse	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Mouse IL-1 alpha Recombinant Rabbit mAb [PSH13-43] - Detector	HA723531B	Mouse	ELISA(Det),ELISA
	Biotin Conjugated Mouse IL-1 alpha Recombinant Rabbit mAb [PSH13-41] - Detector	HA723528B	Mouse	ELISA(Det),ELISA
	Mouse IL-1 alpha Recombinant Rabbit mAb [PSH13-40] - BSA and Azide free (Capture)	HA723526	Mouse	ELISA(Cap)
	Mouse IL-1 alpha Recombinant Rabbit mAb [PSH13-42] - BSA and Azide free (Capture)	HA723529	Mouse	ELISA(Cap)
	Mouse IL-1 alpha Recombinant Rabbit mAb [PSH13-43] - BSA and Azide free (Detector)	HA723530	Mouse	ELISA(Det)
	Mouse IL-1 alpha Recombinant Rabbit mAb [PSH13-41] - BSA and Azide free (Detector)	HA723527	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-1 beta Recombinant Rabbit mAb [PSH03-32] - Detector	HA722625B	Mouse	ELISA(Det),ELISA
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH09-59] - BSA and Azide free	HA723121	Mouse	ELISA
	Biotin Conjugated Mouse IL-1 beta Recombinant Rabbit mAb [PSH04-67]	HA722608B	Mouse	ELISA(Det)
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH04-67] - BSA and Azide free (Detector)	HA722167	Mouse	ELISA(Det)
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH09-58] - BSA and Azide free	HA723120	Mouse	ELISA
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH03-31] - BSA and Azide free (Capture)	HA721976	Mouse	ELISA(Cap)
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH04-66] - BSA and Azide free (Capture)	HA722166	Mouse	ELISA(Cap)
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH03-32] - BSA and Azide free (Detector)	HA721977	Mouse	ELISA(Det)
	Mouse IL-6 Recombinant Rabbit mAb [PS01-74] - BSA and Azide free (Capture/Detector)	HA721600	Mouse	ELISA(Cap),ELISA(Det)
	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PSH0-59] - Detector	HA722646B	Mouse	ELISA(Det)
	Mouse IL-6 Recombinant Rabbit mAb [PSH0-60] - BSA and Azide free (Detector)	HA721413	Mouse	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PSH0-60] - Detector	HA722647B	Mouse	ELISA(Det)
	HRP Conjugated Mouse IL-6 Mouse mAb [3-A4]	M1303-2	Mouse	WB,ELISA
	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PS01-83] - BSA and Azide free (Detector)	HA721602	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PS01-74] - Detector	HA723643B	Mouse	ELISA(Det)
	Mouse IL-6 Recombinant Rabbit mAb [PS01-83] - BSA and Azide free (Detector)	HA721601	Mouse	ELISA(Det)
	Mouse IL-6 Recombinant Rabbit mAb [PSH0-58] - BSA and Azide free (Capture)	HA721411	Mouse	ELISA(Cap)
	Mouse IL-6 Recombinant Rabbit mAb [PSH14-47] - BSA and Azide free (Capture)	HA723642	Mouse	ELISA(Cap)
	Mouse IL-6 Recombinant Rabbit mAb [PSH0-59] - BSA and Azide free (Detector)	HA721412	Mouse	ELISA(Det)
	Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-30] - BSA and Azide free (Detector)	HA721975	Mouse, Rat	ELISA(Det)
	Biotin Conjugated Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-30]	HA722626B	Mouse, Rat	ELISA(Det)
	Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-29] - BSA and Azide free (Capture)	HA721974	Mouse, Rat	ELISA(Cap)
	Mouse TNF alpha Recombinant Rabbit mAb [PSH03-77] - BSA and Azide free (Detector)	HA722041	Mouse	ELISA(Det)
	Biotin Conjugated Mouse TNF alpha Recombinant Rabbit mAb [PSH03-77]	HA722623B	Mouse	ELISA(Det)
	Mouse TNF alpha Recombinant Rabbit mAb [PSH03-76] - BSA and Azide free (Capture)	HA722040	Mouse	ELISA(Cap)
	MSR1 Rabbit pAb	ER1913-21	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
	MSR1 / CD204 Rabbit pAb	ER1902-71	Human	WB,FC
Relative expression	MUM1 Mouse mAb [A7D1]	HA601030	Human	IHC-P,WB
	MUM1 Mouse mAb [A7C12]	HA601031	Human	IHC-P
Relative expression	MUM1 Recombinant Rabbit mAb [JE52-60]	ET7110-44	Human	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment, Relative expression	Phospho-STAT1 (S727) Recombinant Rabbit mAb [SN67-04]	ET1611-20	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment, Relative expression	Phospho-STAT1 (S727) Recombinant Rabbit mAb [SN67-04] - BSA and Azide free	HA750247	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-STAT1 (Y701) Recombinant Rabbit mAb [PSH04-02]	HA722083	Human	WB,IF-Cell
Cell treatment	Phospho-STAT1 (Y701) Recombinant Rabbit mAb [PSH04-02] - BSA and Azide free	HA750925	Human	WB,IF-Cell
Cell treatment	Phospho-STAT6 (Y641) Recombinant Rabbit mAb [PSH06-98]	HA722776	Human, Mouse	WB
Cell treatment	Phospho-STAT6 (Y641) Recombinant Rabbit mAb [PSH06-98] - BSA and Azide free	HA751125	Human, Mouse	WB
	Rat IL-10 Recombinant Rabbit mAb [PSH02-21] - BSA and Azide free (Capture)	HA721786	Rat	ELISA(Cap)
	Biotin Conjugated Rat IL-10 Recombinant Rabbit mAb [PSH02-22]	HA722637B	Rat	ELISA(Det)
	Rat IL-10 Recombinant Rabbit mAb [PSH02-22] - BSA and Azide free (Detector)	HA721787	Rat	ELISA(Det)
	Rat IL-1 beta Recombinant Rabbit mAb [PSH02-60] - BSA and Azide free (Detector)	HA721837	Rat	ELISA(Det)
	Rat IL-1 beta Recombinant Rabbit mAb [PS01-87] - BSA and Azide free (Detector)	HA721610	Rat	ELISA(Det)
	Rat IL-1 beta Recombinant Rabbit mAb [PSH02-59] - BSA and Azide free (Capture)	HA721836	Rat	ELISA(Cap)
	Biotin Conjugated Rat IL-1 beta Recombinant Rabbit mAb [PSH02-60]	HA722634B	Rat	ELISA(Det)
	Rat IL-1 beta Recombinant Rabbit mAb [PS01-86] - BSA and Azide free (Capture)	HA721609	Rat	ELISA(Cap)
	Biotin Conjugated Rat IL-1 beta Recombinant Rabbit mAb [PS01-87] - Detector	HA721611	Rat	ELISA(Det),ELISA
	STAT1 Rabbit pAb	R1408-2	Human, Mouse, Rat, Zebrafish	WB,IF-Cell
KD	STAT1 Recombinant Rabbit mAb [SD20-75]	ET1612-22	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	STAT1 Mouse mAb [G3-B11]	M1407-1	Human	WB,IF-Cell,IHC-P,IF-Tissue
	STAT1 Recombinant Mouse mAb [G3-B11-R]	HA601309	Human	WB,IHC-P
	STAT1 Recombinant Mouse mAb [G3-B11-R] - BSA and Azide free	HA610181	Human	WB,IHC-P
KD	STAT1 alpha Recombinant Rabbit mAb [SJ01-89]	ET1606-39	Human, Mouse	WB,IP,FC,IHC-P,IF-Cell,IF-Tissue
KD	STAT1 alpha Recombinant Rabbit mAb [SJ01-89] - BSA and Azide free	HA750102	Human, Mouse	WB,IP,FC,IHC-P,IF-Cell,IF-Tissue
	STAT1 alpha/beta Rabbit pAb	ER31221	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	STAT6 Recombinant Rabbit mAb [SY13-09]	ET1605-49	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC,IF-Tissue
KD	STAT6 Recombinant Rabbit mAb [SY02-72]	ET1607-61	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	STAT6 Rabbit pAb	ER40121	Human	WB,IF-Cell,IHC-P,FC
	TGF beta 1 Recombinant Rabbit mAb [PSH18-59]	HA724012	Human, Mouse, Monkey	WB,IF-Cell,IHC-P,FC
Relative expression	TGF beta 1 Recombinant Rabbit mAb [PD00-17]	HA721143	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Tissue
	TGF beta 1 Recombinant Rabbit mAb [PSH15-22] - BSA and Azide free (Detector)	HA723726	Human, Mouse, Rat	ELISA(Det)
	TGF beta 1 Mouse mAb [B2-F6]	EM1709-93	Human	WB,IHC-P,FC
	TGF beta 1 Recombinant Rabbit mAb [PSH15-21] - BSA and Azide free (Capture)	HA723725	Human, Mouse, Rat	ELISA(Cap)
	TGF beta 1 Recombinant Rabbit mAb	IRS205RB	Human	mIHC
	Biotin Conjugated TGF beta 1 Recombinant Rabbit mAb [PSH15-22] - Detector	HA723727B	Human, Mouse, Rat	ELISA(Det),ELISA
	TGF beta 1 Rabbit pAb	HA500496	Human	WB
Relative expression	TGF beta 1 Recombinant Rabbit mAb [PD00-17] - BSA and Azide free	HA750485	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Tissue
Cell treatment	TNF alpha Rabbit pAb	R1203-1	Mouse	WB,IF-Cell,ELISA
	TNF alpha Rabbit pAb	ER65189	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Cell treatment	TNF alpha Recombinant Rabbit mAb [PS01-07]	HA722022	Human, Mouse	WB
Cell treatment	TNF alpha Recombinant Rabbit mAb [PS01-07] - BSA and Azide free	HA750897	Human, Mouse	WB
	TNF alpha Rabbit pAb	ER1919-22	Human, Mouse	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

