

AD & 神经炎症假说研究-产品清单

高级验证	产品名	货号	反应种属	验证应用
	ALDH1L1 Recombinant Rabbit mAb [JU53-54]	ET7106-64	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr
	ALDH1L1 Recombinant Rabbit mAb [JU53-54] - BSA and Azide free	HA750457	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr
	Arginase1 Recombinant Rabbit mAb	IRS081RB	Human	mIHC
	Arginase1 Recombinant Rabbit mAb	IRS082RB	Mouse	mIHC
	B7-1/CD80 Rabbit pAb	ER1904-08	Mouse, Rat	WB
	BACH1 Rabbit pAb	ER62786	Human	WB,IHC-P,IF-Cell
	BACH1 Recombinant Rabbit mAb [PSH18-52]	HA724005	Human	WB,FC,IP
	Beta TRCP Rabbit pAb	HA500358	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	BTRC Mouse mAb [D8-C5]	EM1706-79	Human	WB,IHC-P,FC
	CCL17 Rabbit pAb	ER65950	Human	WB
	CCL18 Rabbit pAb	ER65928	Human	IHC-P
	CCL2 Mouse mAb [E11-D8]	EM1710-22	Human, Mouse, Monkey	WB,IF-Cell,IHC-P,FC
	CCL2 / MCP-1 Rabbit pAb	HA500042	Human	WB,IHC-P,IF-Cell
	CCL2 / MCP-1 Rabbit pAb	HA500267	Mouse	WB,IF-Cell,FC
	CCL2 / MCP-1 Recombinant Rabbit mAb [PSH12-44] - BSA and Azide free	HA751533	Mouse	WB,IF-Cell,FC
Cell treatment	CCL2 / MCP-1 Recombinant Rabbit mAb [PSH18-16]	HA723979	Mouse, Rat	WB,IF-Cell,FC
	CCL2 / MCP-1 Recombinant Rabbit mAb [PSH12-44]	HA723691	Mouse	WB,IF-Cell,FC
	CCL3 / MIP-1 alpha Rabbit pAb	R1510-9	Human, Mouse, Rat	IF-Cell,IHC-P,FC,ELISA
Cell treatment	CCL3 / MIP-1 alpha Recombinant Rabbit mAb [PSH19-79]	HA724103	Mouse, Rat	WB,IF-Cell,FC,IP
Cell treatment, Relative expression	CCL3 / MIP-1 alpha Recombinant Rabbit mAb [PSH11-28] - BSA and Azide free	HA751407	Human, Mouse	WB,IF-Cell,IHC-P
Cell treatment, Relative expression	CCL3 / MIP-1 alpha Recombinant Rabbit mAb [PSH11-28]	HA723310	Human, Mouse	WB,IF-Cell,IHC-P

Cell treatment, Relative expression	CCL3 / MIP-1 alpha Recombinant Rabbit mAb [PSH11-27] - BSA and Azide free	HA751406	Human	WB,IF-Cell,FC
Cell treatment	CCL3 / MIP-1 alpha Rabbit pAb	R1510-10	Human	WB
Cell treatment, Relative expression	CCL3 / MIP-1 alpha Recombinant Rabbit mAb [PSH11-27]	HA723309	Human	WB,IF-Cell,FC
	CCR2 Recombinant Rabbit mAb [SN707]	ET1611-65	Human, Mouse	WB,IHC-P,IP,FC
	CCR2 Recombinant Rabbit mAb [SN707] - BSA and Azide free	HA750266	Human, Mouse	WB,IHC-P,IP,FC
Relative expression	CD11b Mouse mAb [A2-9]	M1211-6	Human	WB,IHC-P,FC,IF-Tissue
	CD11b Mouse mAb [4-6]	EM1701-42	Human	IHC-P,IF-Cell,FC,WB,IF-Tissue
Relative expression	CD11b Recombinant Rabbit mAb [JE01-17]	HA722856	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Tissue
Relative expression	CD11b Recombinant Rabbit mAb [JU93-81] - BSA and Azide free	HA750451	Human	WB,IP,IHC-P,IF-Cell,IF-Tissue
Relative expression	CD11b Recombinant Rabbit mAb [PSH03-96] - BSA and Azide free	HA750921	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	CD11b Recombinant Rabbit mAb	IRS014	Human	mIHC
Relative expression	CD11b Recombinant Rabbit mAb [JU93-81]	ET1706-04	Human	WB,IP,IHC-P,IF-Cell,IF-Tissue,mIHC
Relative expression	CD11b Recombinant Rabbit mAb [PSH03-96]	HA722075	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	CD11b Mouse mAb [A2-B5-D2]	EM1902-35	Human, Mouse	WB
	CD11b/ITAM Mouse mAb [4G2]	RT1107	Human, Mouse	WB,IP,IF,FC
	CD16 Recombinant Rabbit mAb	IRS017	Human	mIHC
	CD16 Mouse mAb [Q32]	EM60067	Human	WB,IHC-P
Relative expression	CD16 Recombinant Rabbit mAb [PD00-12]	HA721149	Human	IHC-P,IF-Tissue,IF-Cell,mIHC,WB
	CD16 Recombinant Rabbit mAb [JE49-79]	ET7109-97	Human, Mouse, Rat	WB
	CD163 Rabbit pAb	IRS109RB	Mouse	mIHC
	CD163 Recombinant Rabbit mAb [JA51-30] - BSA and Azide free	HA750412	Human	WB,IHC-P,IF-Tissue,IP
	CD163 Recombinant Rabbit mAb	IRS002	Human	mIHC
	CD163 Rabbit pAb	ER1804-03	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	CD163 Recombinant Rabbit mAb [PSH15-06]	HA723706	Human	IHC-P,WB
	CD163 Recombinant Antibody [JA51-30] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610364	Human	IHC-P
	CD163 Recombinant Antibody [JA51-30] - Rat IgG1 (Chimeric)	HA601549	Human	IHC-P
	CD163 Mouse mAb [A3B3]	EM1901-90	Human	WB,IHC-P,FC
	CD163 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS276RT	Human	mIHC
	CD163 Mouse mAb [A3B5]	EM1901-91	Human	IHC-P,FC
Relative expression	CD163 Recombinant Rabbit mAb [PSH18-98]	HA724034	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	CD163 Recombinant Rabbit mAb [JA51-30]	ET1704-43	Human	WB,IHC-P,mIHC,IF-Tissue,IP
	CD182 Mouse mAb [5F7A11]	EM1711-55	Human	WB,IHC-P,IF-Cell,FC
	CD23 Rabbit pAb	HA500386	Human	IHC-P,FC
	CD23 Recombinant Rabbit mAb [PD00-03]	HA721139	Human	IHC-P,mIHC,IF-Tissue
	CD23 Recombinant Rabbit mAb	IRS021	Human	mIHC
	APC Conjugated CD25 / IL-2 Receptor alpha Recombinant Mouse mAb [PSH14-46]	HA601530F	Human	FC
	CD25 / IL-2 Receptor alpha Recombinant Mouse mAb [PSH14-46]	HA601486	Human, Mouse	FC,IHC-P,IF-Tissue
	CD32 Rabbit pAb	ER62945	Human, Mouse, Rat	WB,IHC-P
Cell treatment	CD32 Recombinant Rabbit mAb [JM10-70]	ET1703-43	Human	WB,IHC-P,IF-Cell,FC,IP
Cell treatment	CD32 Recombinant Rabbit mAb [JM10-70] - BSA and Azide free	HA750382	Human	WB,IHC-P,IF-Cell,FC,IP
	CD32A + CD32B + CD32C Recombinant Rabbit mAb [JE64-17]	HA721930	Human	WB,IHC-P
	CD32B Mouse mAb [A9C11]	HA601121	Human	WB,IHC-P
	CD32-C Rabbit pAb	ER62946	Human	WB
Relative expression	CD33 Recombinant Rabbit mAb [PSH02-51] - BSA and Azide free	HA750814	Human	WB,IF-Cell,FC
Relative expression	CD33 Recombinant Rabbit mAb [PSH02-51]	HA721827	Human	WB,IF-Cell,FC
Relative expression, Cell treatment	CD33 Recombinant Rabbit mAb [PSH02-50] - BSA and Azide free	HA750813	Human	WB,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CD33 Mouse mAb [B7-C6]	EM1709-20	Human	IHC-P,FC
Relative expression, Cell treatment	CD33 Recombinant Rabbit mAb [PSH02-50]	HA721826	Human	WB,IF-Cell,FC
	CD33 Rabbit pAb	ER1905-89	Mouse	WB,IHC-P,IF-Tissue
Cell treatment	CD36 Recombinant Rabbit mAb [JJ2005]	ET1701-24	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
Cell treatment	CD36 Recombinant Rabbit mAb [JJ2005] - BSA and Azide free	HA750307	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	CD36 Rabbit pAb	ER1802-7	Human, Mouse, Rat	IF-Cell,IHC-P,FC
Relative expression	CD36 Recombinant Rabbit mAb [PSH18-47]	HA724000	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	CD36 Recombinant Rabbit mAb	IRS118RB	Human	mIHC
Relative expression	CD40 Recombinant Rabbit mAb [PSH05-90]	HA722488	Human, Cynomolgus monkey, Pig	WB,IF-Cell,FC
Relative expression	CD40 Recombinant Rabbit mAb [PSH05-90] - BSA and Azide free	HA751037	Human, Cynomolgus monkey, Pig	WB,IF-Cell,FC
Relative expression	CD40 Mouse mAb [A1E6]	EM1901-40	Human	WB,IF-Cell,IHC-P
	CD40 Recombinant Rabbit mAb [PSH10-01]	HA723170	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-P,IHC-Fr,IF-Cell,FC,IF-Tissue
	CD40 Recombinant Rabbit mAb [PSH10-01] - BSA and Azide free	HA751331	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-P,IHC-Fr,IF-Cell,FC,IF-Tissue
	CD40 Mouse mAb [A1E5]	EM1901-41	Human	WB,IHC-P
	CD40 Rabbit pAb	ER1803-54	Human, Rat, Mouse	WB,IF-Cell,IHC-P,FC
Relative expression	CD47 Recombinant Rabbit mAb [PSH10-31] - BSA and Azide free	HA751348	Human	WB,IF-Cell,FC
Relative expression	CD47 Recombinant Rabbit mAb [PSH10-31]	HA723198	Human	WB,IF-Cell,FC
	CD47 Mouse mAb [1A9]	EM1701-38	Human	IF-Cell,IHC-P,FC
	CD47 Recombinant Rabbit mAb [PD00-33]	HA721158	Human	IHC-P,FC
	iFluor™ 488 Conjugated CD47 Mouse mAb [1B1]	HA600104F	Human	FC
	CD59 Recombinant Rabbit mAb [JM10-71]	ET1703-28	Human	WB,IP
Relative expression	CD68 Recombinant Mouse mAb [A3C3-R]	HA601290	Human	WB,IF-Cell,IHC-P,FC
Relative expression	CD68 Recombinant Mouse mAb [A3C2-R]	HA601291	Human	WB,IF-Cell,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	CD68 Mouse mAb [A3C3]	EM1901-94	Human	WB,IHC-P,FC
	CD68 Recombinant Antibody [PDM0-13] - Rabbit IgG (Chimeric)	HA723708	Human	IF-Tissue,IHC-P
Relative expression, Cell treatment	CD68 Recombinant Rabbit mAb [PSH06-65]	HA722737	Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,IHC-Fr
Relative expression	CD68 Recombinant Mouse mAb [PDM0-13]	HA601115	Human	IHC-P,IF-Cell,FC,mIHC,IF-Tissue
Relative expression	CD68 Mouse mAb [A3C2]	EM1901-93	Human	WB,IHC-P,FC
Relative expression	CD68 Recombinant Mouse mAb [A3C4-R]	HA601292	Human	WB,IF-Cell,IHC-P
	CD68 Recombinant Rabbit mAb [SN07-27]	ET1611-53	Human	WB
	CD68 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS264RT	Human	mIHC
	CD68 Recombinant Antibody - Rabbit IgG (Chimeric)	IRS266RB	Human	mIHC
Cell treatment	CD68 Recombinant Rabbit mAb [PSH05-47] - BSA and Azide free	HA751006	Mouse, Rat	WB,IF-Cell,IHC-P,IHC-Fr,IF-Tissue
Relative expression, Cell treatment	CD68 Recombinant Rabbit mAb [PSH06-65] - BSA and Azide free	HA751094	Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,IHC-Fr
Relative expression	CD68 Recombinant Mouse mAb [A3C4-R] - BSA and Azide free	HA610165	Human	WB,IF-Cell,IHC-P
Relative expression	CD68 Recombinant Mouse mAb [PDM0-13] - BSA and Azide free	HA610034	Human	IHC-P,IF-Cell,FC,IF-Tissue
	CD68 Recombinant Antibody [PDM0-13] - Rabbit IgG (Chimeric) - BSA and Azide free	HA751544	Human	IF-Tissue,IHC-P
	CD68 Recombinant Antibody [PDM0-13] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610347	Human	IF-Tissue,IHC-P
Relative expression	CD68 Mouse mAb [A3C4]	EM1901-95	Human	WB,IHC-P,FC,mIHC
	CD68 Recombinant Mouse mAb	IRS004	Human	mIHC
	CD68 Recombinant Antibody [PDM0-13] - Rat IgG1 (Chimeric)	HA601520	Human	IF-Tissue,IHC-P
	CD68 Rabbit pAb	ER1901-32	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Relative expression	CD68 Recombinant Mouse mAb [A3C2-R] - BSA and Azide free	HA610164	Human	WB,IF-Cell,IHC-P
Cell treatment	CD68 Recombinant Rabbit mAb [PSH05-47]	HA722285	Mouse, Rat	WB,IF-Cell,IHC-P,IHC-Fr,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CD68 Mouse mAb [F3-D7]	EM1706-11	Human	WB,IHC-P,IF-Cell,FC
Relative expression, Cell treatment	CD80 Recombinant Rabbit mAb [PSH18-13]	HA723976	Mouse	WB,IF-Cell,FC,IP
	CD80 Rabbit pAb	ER1905-99	Human, Mouse, Rat	WB
Relative expression	CD80 Recombinant Rabbit mAb [PSH09-79]	HA723144	Human	WB,IF-Cell,FC
	CD80 Rabbit pAb	1007-8	Human	WB,IF-Cell
Relative expression	CD80 Mouse mAb [8-E5]	M1007-10	Human, Mouse	WB,IHC-P,FC
	CD80 Recombinant Mouse mAb [8-E5-R]	HA601247	Human	WB
	CD80 Recombinant Mouse mAb [8-E5-R] - BSA and Azide free	HA610122	Human	WB
Relative expression	CD80 Recombinant Rabbit mAb [PSH09-79] - BSA and Azide free	HA751321	Human	WB,IF-Cell,FC
Relative expression	CD80 Recombinant Rabbit mAb [JF100-4]	ET1702-95	Human	WB,IHC-P
	CD86 Rabbit pAb	ER1906-01	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
	CD86 Recombinant Rabbit mAb [SJ20-00]	ET1606-50	Human, Mouse, Rat	WB,IF-Cell,IP
	CD86 Recombinant Rabbit mAb [SJ20-00] - BSA and Azide free	HA750109	Human, Mouse, Rat	WB,IF-Cell,IP
	Cdk11A/B Rabbit pAb	ER62965	Human	WB
	Cullin 3 Recombinant Rabbit mAb [JE37-29]	HA721711	Human, Mouse, Rat	WB,IHC-P
	Cullin 3 Recombinant Rabbit mAb [JE57-80]	HA721882	Human, Mouse, Rat	WB,IHC-P
	Cullin-3 Rabbit pAb	ER1706-70	Human, Mouse	WB,IHC-P,FC
	CX3CR1 Mouse mAb [A5C7]	HA600074	Human	IF-Cell
	CXCL16 Rabbit pAb	ER1906-85	Human, Mouse	IHC-P,IF-Tissue,FC
Cell treatment	CXCL1 / GRO alpha Recombinant Rabbit mAb [PSH07-53]	HA722841	Human	WB
Cell treatment	CXCL1 / GRO alpha Recombinant Rabbit mAb [PSH07-53] - BSA and Azide free	HA751174	Human	WB
Cell treatment	CXCL9 Recombinant Rabbit mAb [PSH12-34]	HA723608	Human	WB,IF-Cell
Cell treatment	CXCL9 Recombinant Rabbit mAb [PSH12-34] - BSA and Azide free	HA751512	Human	WB,IF-Cell
	CXCR1 Rabbit pAb	ER1906-86	Human	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CXCR1 Rabbit pAb	HA500526	Human	WB,IF-Cell
	CXCR2 Rabbit pAb	ER1906-87	Human	WB
	DC-SIGNR1/CD209b Rabbit pAb	ER1907-24	Mouse	WB
Cell treatment	EAAT1 Recombinant Rabbit mAb [JA30-35]	ET1704-54	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr
Cell treatment	EAAT1 Recombinant Rabbit mAb [JA30-35] - BSA and Azide free	HA750415	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr
	Eotaxin-2 Rabbit pAb	ER65789	Human	IHC-P
	FCG3B Rabbit pAb	ER62049	Human	WB
KD	GCLC Recombinant Rabbit mAb [JA08-03]	ET1704-38	Human, Rat	WB,IHC-P,FC
	GCLM Recombinant Rabbit mAb [JM39-61]	ET1705-87	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,IP,IHC-Fr
	GFAP Rabbit pAb	R1308-9	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	GFAP Recombinant Rabbit mAb [PSH05-16] - BSA and Azide free (Capture)	HA722234	Human, Mouse, Rat	ELISA(Cap)
	GFAP Recombinant Mouse mAb [PSH05-17] - BSA and Azide free (Detector)	HA601342	Human, Mouse, Rat	ELISA(Det)
	GFAP Mouse mAb [1-D4]	EM140707	Human, Rat	WB,IHC-P,IF-Tissue,IHC-Fr
	iFluor™ 594 Conjugated GFAP Mouse mAb [1-D4]	HA600102F	Human, Rat	IF-Tissue
	iFluor™ 488 Conjugated GFAP Mouse mAb [1-D4]	HA600101F	Human, Rat	IF-Tissue,IHC-Fr
Relative expression	GFAP Recombinant Mouse mAb [PSH10-24] - BSA and Azide free	HA610225	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Tissue
	GFAP Rabbit pAb	R1308-10	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	GFAP Recombinant Antibody - Rat IgG1 (Chimeric)	IRS148RT	Mouse	mIHC
Relative expression	GFAP Recombinant Antibody [SA03-04] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610317	Mouse, Rat	IHC-Fr,WB
	Biotin Conjugated GFAP Mouse mAb [PSH05-17]	HA601355B	Human, Mouse, Rat	ELISA(Det)
Relative expression	GFAP Recombinant Antibody [SA03-04] - Rat IgG1 (Chimeric)	HA601483	Mouse, Rat	IHC-Fr,WB
	GFAP Recombinant Rabbit mAb [SA03-04]	ET1601-23	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr,mIHC,IP
	iFluor™ 647 Conjugated GFAP Mouse mAb [1-D4]	HA600103F	Human, Rat	IF-Tissue
	GFAP Recombinant Rabbit mAb	IRS069	Mouse	mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	GFAP Recombinant Mouse mAb [PSH10-24]	HA601380	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Tissue
	Glutamine Synthetase Recombinant Mouse mAb [A3G2-R] - BSA and Azide free	HA610100	Human, Mouse, Rat	WB,IHC-P
	Glutamine Synthetase Recombinant Mouse mAb [A3G2-R]	HA601225	Human, Mouse, Rat	WB,IHC-P
	Glutamine Synthetase Recombinant Antibody [PSH10-48] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610268	Human, Mouse, Rat	IHC-Fr,WB
	Glutamine Synthetase Mouse mAb [A3G2]	EM1902-39	Human, Mouse, Rat	WB,IHC-P,FC,mIHC,IF-Tissue
	Glutamine Synthetase Recombinant Mouse mAb	IRS231MS	Mouse	mIHC
	Glutamine Synthetase Mouse mAb	IRS033	Human	mIHC
	Glutamine Synthetase Recombinant Rabbit mAb [PSH10-48] - BSA and Azide free	HA751351	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-Fr
	Glutamine Synthetase Recombinant Antibody [PSH10-48] - Rat IgG1 (Chimeric)	HA601428	Human, Mouse, Rat	IHC-Fr,WB
	Glutamine Synthetase Recombinant Rabbit mAb [PSH10-48]	HA723212	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-Fr
	Glutamine Synthetase Recombinant Antibody [PSH10-48] - Chicken IgY (Chimeric)	HA601516	Mouse, Rat	IHC-Fr,WB
	Glutamine Synthetase Recombinant Antibody [PSH10-48] - Chicken IgY (Chimeric) - BSA and Azide free	HA610346	Mouse, Rat	IHC-Fr,WB
	GRO alpha Rabbit pAb	ER63645	Human, Mouse, Rat	IHC-P
KD	GSK3 alpha Rabbit pAb	R1312-1	Human	WB,IF-Cell,IHC-P,FC
	GSK3 alpha Recombinant Rabbit mAb [JE33-27]	HA722217	Human, Mouse, Rat, Monkey	WB,IHC-P
	GSK3 alpha/beta Rabbit pAb	ER63648	Human, Mouse, Rat	WB,IHC-P
KD	GSK3 beta Recombinant Rabbit mAb [SY28-03]	ET1607-71	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IHC-Fr
	GSK3 beta Rabbit pAb	R1510-4	Human, Mouse, Rat	IF-Cell,IHC-P,FC
KD	GSK3 beta Recombinant Rabbit mAb [SY28-03] - BSA and Azide free	HA750131	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IHC-Fr
	Heme Oxygenase 1 (HO-1) Recombinant Rabbit mAb	IRS156RB	Human	mIHC
KD	Heme Oxygenase 1 (HO-1) Recombinant Rabbit mAb [PSH02-76] - BSA and Azide free	HA750831	Human, Mouse	WB,IF-Cell,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	Heme Oxygenase 1 (HO-1) Recombinant Rabbit mAb [PSH02-76]	HA721854	Human, Mouse	WB,IF-Cell,IHC-P
	Heme Oxygenase 1 (HO-1) Recombinant Rabbit mAb [SP08-07] - BSA and Azide free	HA750079	Human	WB,IHC-P,IP,FC,IF-Tissue
	Heme Oxygenase 1 (HO-1) Recombinant Rabbit mAb [SP08-07]	ET1604-45	Human	WB,IHC-P,IP,FC,IF-Tissue
	iFluor™ 594 Conjugated Heme Oxygenase 1 (HO-1) Recombinant Rabbit mAb [SP08-07]	HA720169F	Human	IF-Tissue
	Heme Oxygenase 1 (HO-1) Recombinant Rabbit mAb	IRS232RB	Mouse	mIHC
KD	Heme Oxygenase 1 (HO-1) Rabbit pAb	ER1802-73	Human, Mouse	WB,IF-Cell,IHC-P,FC
	HLA-DR Rabbit pAb	R1510-35	Human	WB,IF-Cell,IHC-P,FC
	HLA-DR Rabbit pAb	R1510-34	Human, Mouse	WB,IF-Cell,IHC-P,FC
	HLA-DR Mouse mAb [10-D8]	M1701-5	Human	WB,IHC-P,IF-Cell
	HLA-DR Recombinant Rabbit mAb [SC06-78] - BSA and Azide free	HA750230	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	HLA-DR Recombinant Mouse mAb	IRS154MS	Human	mIHC
	HLA-DR Recombinant Rabbit mAb [SC06-78]	ET1610-66	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	HLA-DR Recombinant Mouse mAb [10-D8-R] - BSA and Azide free	HA610127	Human	WB,IF-Cell,IHC-P,IF-Tissue
	HLA-DR Recombinant Mouse mAb [10-D8-R]	HA601252	Human	WB,IF-Cell,IHC-P,IF-Tissue
	Human CCL17 / TARC Recombinant Rabbit mAb [PSH13-33] - BSA and Azide free (Capture)	HA723513	Human	ELISA(Cap)
	Biotin Conjugated Human CCL17 / TARC Recombinant Rabbit mAb [PSH13-35] - Detector	HA723517B	Human	ELISA(Det),ELISA
	Human CCL17 / TARC Recombinant Rabbit mAb [PSH13-35] - BSA and Azide free (Detector)	HA723516	Human	ELISA(Det)
	Human CCL17 / TARC Recombinant Rabbit mAb [PSH13-34] - BSA and Azide free (Detector)	HA723514	Human	ELISA(Det)
	Biotin Conjugated Human CCL17 / TARC Recombinant Rabbit mAb [PSH13-34] - Detector	HA723515B	Human	ELISA(Det),ELISA
	Human CCL2 / MCP-1 Recombinant Rabbit mAb [PSH18-21] - BSA and Azide free (Capture)	HA725195	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human CCL2 / MCP-1 Recombinant Rabbit mAb [PSH18-22] - BSA and Azide free (Detector)	HA725196	Human	ELISA(Det)
	Biotin Conjugated Human CCL2 / MCP-1 Recombinant Rabbit mAb [PSH18-22] - Detector	HA725196B	Human	ELISA(Det)
	Human CCL3 / MIP-1 alpha Recombinant Rabbit mAb [PSH11-77] - BSA and Azide free (Detector)	HA723378	Human	ELISA(Det)
	Biotin Conjugated Human CCL3 / MIP-1 alpha Recombinant Rabbit mAb [PSH11-77] - Detector	HA723379B	Human	ELISA(Det)
	Human CCL3 / MIP-1 alpha Recombinant Rabbit mAb [PSH11-78] - BSA and Azide free (Capture)	HA723380	Human	ELISA(Cap)
	Human CCL3 / MIP-1 alpha Recombinant Rabbit mAb [PSH11-76] - BSA and Azide free (Capture)	HA723377	Human	ELISA(Cap)
	Human CCL5 Recombinant Rabbit mAb [PSH13-17] - BSA and Azide free (Detector)	HA725056	Human	ELISA(Det)
	HRP Conjugated Human CCL5 Recombinant Rabbit mAb [PSH13-17] - Detector	HA725057H	Human	ELISA(Det),ELISA
	Human CCL5 Recombinant Rabbit mAb [PSH13-16] - BSA and Azide free (Capture)	HA725055	Human	ELISA(Cap)
	Biotin Conjugated Human CCL5 Recombinant Rabbit mAb [PSH13-17] - Detector	HA725057B	Human	ELISA(Det),ELISA
	FITC Conjugated Human CD11b Recombinant Mouse mAb [A2-9-R]	HA600130F	Human	FC
	Human CD16 Recombinant Mouse mAb [PSH04-95]	HA601340	Human	FC
	PE Conjugated Human CD16 Recombinant Mouse mAb [PSH04-95]	HA600129F	Human	FC
	Human CD163 Recombinant Rabbit mAb [PSH13-78] - BSA and Azide free (Capture)	HA723566	Human	ELISA(Cap)
	Biotin Conjugated Human CD163 Recombinant Rabbit mAb [PSH13-79] - Detector	HA723568B	Human	ELISA(Det),ELISA
	Human CD163 Recombinant Rabbit mAb [PSH13-80] - BSA and Azide free (Capture)	HA723569	Human	ELISA(Cap)
	Human CD163 Recombinant Rabbit mAb [PSH13-79] - BSA and Azide free (Detector)	HA723567	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human CD33 Recombinant Rabbit mAb [PSH05-26]	HA722599B	Human	ELISA(Det)
	FITC Conjugated Human CD33 Recombinant Rabbit mAb [PSH04-96]	HA720197F	Human	FC
	Human CD33 Recombinant Rabbit mAb [PSH05-25] - BSA and Azide free (Capture)	HA722254	Human	ELISA(Cap)
	Human CD33 Recombinant Rabbit mAb [PSH05-26] - BSA and Azide free (Detector)	HA722255	Human	ELISA(Det)
	Biotin Conjugated Human CD47 Recombinant Rabbit mAb [PSH11-06] - Detector	HA723274B	Human	ELISA(Det),ELISA
	Human CD47 Recombinant Rabbit mAb [PSH11-05] - BSA and Azide free (Capture)	HA723272	Human	ELISA(Cap)
	Human CD47 Recombinant Rabbit mAb [PSH11-06] - BSA and Azide free (Detector)	HA723273	Human	ELISA(Det)
	FITC Conjugated Human CD68 Recombinant Mouse mAb [A3C3-R]	HA600132F	Human	FC
	Biotin Conjugated Human CD80 Recombinant Rabbit mAb [PSH09-63] - Detector	HA723127B	Human	ELISA(Det),ELISA
	Human CD80 Recombinant Rabbit mAb [PSH09-62] - BSA and Azide free (Detector)	HA723124	Human	ELISA(Det)
	Human CD80 Recombinant Rabbit mAb [PSH09-61] - BSA and Azide free (Capture)	HA723123	Human	ELISA(Cap)
	Human CD80 Recombinant Rabbit mAb [PSH09-63] - BSA and Azide free (Detector)	HA723125	Human	ELISA(Det)
	Biotin Conjugated Human CD80 Recombinant Rabbit mAb [PSH09-62] - Detector	HA723126B	Human	ELISA(Det),ELISA
	Human CXCL10 / IP-10 Recombinant Rabbit mAb [PSH12-36] - BSA and Azide free (Capture)	HA725007	Human	ELISA(Cap)
	Biotin Conjugated Human CXCL10 / IP-10 Recombinant Rabbit mAb [PSH12-37] - Detector	HA725009B	Human	ELISA(Det),ELISA
	Human CXCL10 / IP-10 Recombinant Rabbit mAb [PSH12-37] - BSA and Azide free (Detector)	HA725008	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human CXCL1/GRO alpha Recombinant Rabbit mAb [PSH05-95]	HA722587B	Human	ELISA(Det)
	Human CXCL1/GRO alpha Recombinant Rabbit mAb [PSH05-95] - BSA and Azide free (Detector)	HA722494	Human	ELISA(Det)
	Human CXCL1/GRO alpha Recombinant Rabbit mAb [PSH05-94] - BSA and Azide free (Capture)	HA722493	Human	ELISA(Cap)
	Human CXCL9 Recombinant Rabbit mAb [PSH12-34] - BSA and Azide free (Capture)	HA725004	Human	ELISA(Cap)
	Human CXCL9 Recombinant Rabbit mAb [PSH12-35] - BSA and Azide free (Detector)	HA725005	Human	ELISA(Det)
	Biotin Conjugated Human CXCL9 Recombinant Rabbit mAb [PSH12-35] - Detector	HA725006B	Human	ELISA(Det),ELISA
	Human IL-10 Recombinant Rabbit mAb [PS00-98] - BSA and Azide free (Capture)	HA721236	Human	ELISA(Cap)
	Biotin Conjugated Human IL-10 Recombinant Rabbit mAb [PS00-81]	HA722652B	Human	ELISA(Det)
	Human IL-10 Recombinant Rabbit mAb [PS00-81] - BSA and Azide free (Detector)	HA721237	Human	ELISA(Det)
	Biotin Conjugated Human IL-12B p40 Recombinant Rabbit mAb [PSH02-58]	HA722633B	Human	ELISA(Det)
	Human IL-12B p40 Recombinant Rabbit mAb [PSH02-58] - BSA and Azide free (Detector)	HA721835	Human	ELISA(Det)
	Human IL-12B p40 Recombinant Rabbit mAb [PSH03-04] - BSA and Azide free (Detector)	HA721934	Human	ELISA(Det)
	Human IL-12B p40 Recombinant Rabbit mAb [PSH02-57] - BSA and Azide free (Capture)	HA721834	Human	ELISA(Cap)
	Biotin Conjugated Human IL-12B p40 Recombinant Rabbit mAb [PSH03-04]	HA722629B	Human	ELISA(Det)
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS00-82] - BSA and Azide free	HA750745	Human	WB
	Biotin Conjugated Human IL-1 beta Recombinant Rabbit mAb [PSH16-86] - Detector	HA723882B	Human	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human IL-1 beta Recombinant Rabbit mAb [PS00-82]	HA722649B	Human	ELISA(Det)
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS01-39]	HA721712	Human	WB
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS00-82]	HA721713	Human	WB
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-86] - BSA and Azide free (Detector)	HA723882	Human	ELISA(Det)
	Human IL-1 beta Recombinant Rabbit mAb [PS00-82] - BSA and Azide free (Detector)	HA721279	Human	ELISA(Det)
	Human IL-1 beta Recombinant Rabbit mAb [PS01-39] - BSA and Azide free (Capture)	HA721278	Human	ELISA(Cap)
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-85] - BSA and Azide free (Capture)	HA723881	Human	ELISA(Cap)
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS01-39] - BSA and Azide free	HA750744	Human	WB
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-87] - BSA and Azide free (Capture)	HA723883	Human	ELISA(Cap)
	Human IL-23A Recombinant Rabbit mAb [PSH03-03] - BSA and Azide free (Capture)	HA721933	Human	ELISA(Cap)
	Human IL-2 Receptor alpha Recombinant Rabbit mAb [PSH05-66] - BSA and Azide free (Detector)	HA722386	Human	ELISA(Det)
	Biotin Conjugated Human IL-2 Receptor alpha Recombinant Rabbit mAb [PS01-03]	HA722650B	Human	ELISA(Det)
	Biotin Conjugated Human IL-2 Receptor alpha Recombinant Rabbit mAb [PSH05-66]	HA722594B	Human	ELISA(Det)
	Human IL-2 Receptor alpha Recombinant Rabbit mAb [PS01-03] - BSA and Azide free (Detector)	HA721281	Human	ELISA(Det)
	Human IL-2 Receptor alpha Recombinant Rabbit mAb [PSH05-65] - BSA and Azide free (Capture)	HA722385	Human	ELISA(Cap)
	Human IL-2 Receptor alpha Recombinant Rabbit mAb [PS00-52] - BSA and Azide free (Capture)	HA721280	Human	ELISA(Cap)
	Human IL-33 Recombinant Rabbit mAb [PSH03-06] - BSA and Azide free (Detector)	HA721936	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human IL-33 Recombinant Rabbit mAb [PSH03-06]	HA722628B	Human	ELISA(Det)
	Human IL-33 Recombinant Rabbit mAb [PSH03-05] - BSA and Azide free (Capture)	HA721935	Human	ELISA(Cap)
	Biotin Conjugated Human IL-4 Recombinant Rabbit mAb [PS00-66] - Detector	HA721368	Human	ELISA(Det)
	Human IL-4 Recombinant Rabbit mAb [PS00-66] - BSA and Azide free (Detector)	HA721231	Human	ELISA(Det)
	Human IL-4 Recombinant Rabbit mAb [PS00-67] - BSA and Azide free (Capture)	HA721230	Human	ELISA(Cap)
	HRP Conjugated Human IL-4 Recombinant Rabbit mAb [PS00-66] - Detector	HA721295	Human	ELISA(Det)
	Biotin Conjugated Human IL-6 Recombinant Rabbit mAb [PS00-15] - Detector	HA722233B	Human	ELISA(Det)
	HRP Conjugated Human IL-6 Recombinant Rabbit mAb [PS00-15] - Detector	HA721293	Human	ELISA(Det)
	Human IL-6 Recombinant Rabbit mAb [PS00-21] - BSA and Azide free (Capture)	HA721167	Human	ELISA(Cap)
	Human IL-6 Mouse mAb [A6D2]	HA601050	Human	WB,ELISA
	Human IL-6 Recombinant Rabbit mAb [PS00-15] - BSA and Azide free (Detector)	HA721168	Human	ELISA(Det)
	Human IL-6 Mouse mAb [A6D1]	HA601051	Human	WB,ELISA
	Human IL-8 Recombinant Rabbit mAb [PS01-22] - BSA and Azide free (Detector)	HA722099	Human	ELISA(Det)
	Biotin Conjugated Human IL-8 Recombinant Rabbit mAb [PS01-22]	HA722611B	Human	ELISA(Det)
	Human IL-8 Recombinant Rabbit mAb [PS00-83] - BSA and Azide free (Capture)	HA722098	Human	ELISA(Cap)
	HRP Conjugated Human sTNF RI Recombinant Rabbit mAb [PSH09-90] - Detector	HA723374H	Human	ELISA(Det)
	Human sTNF RI Recombinant Rabbit mAb [PSH09-90] - BSA and Azide free (Detector)	HA723159	Human	ELISA(Det)
	Human sTNF RI Recombinant Rabbit mAb [PSH09-89] - BSA and Azide free (Capture)	HA723158	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human sTNF RI Recombinant Rabbit mAb [PSH09-90] - BSA and Azide free (Detector)	HA723160B	Human	ELISA(Det),ELISA
	Human TNF alpha Recombinant Rabbit mAb [PS00-40]	HA721286	Human	IHC-P,mIHC
	HRP Conjugated Human TNF alpha Recombinant Rabbit mAb [PS01-07] - Detector	HA721294	Human	ELISA(Det)
	Biotin Conjugated Human TNF alpha Recombinant Rabbit mAb [PS01-07] - Detector	HA721559	Human	ELISA(Det)
	Human TNF alpha Recombinant Rabbit mAb [PS00-40] - BSA and Azide free	HA750593	Human	IHC-P
	Human TNF alpha Recombinant Rabbit mAb [PS00-40] - BSA and Azide free (Capture)	HA721265	Human	ELISA(Cap)
	Human TNF alpha Recombinant Rabbit mAb [PS01-07] - BSA and Azide free (Detector)	HA721266	Human	ELISA(Det)
	Human TNF alpha Recombinant Rabbit mAb	IRS040	Human	mIHC
	HRP Conjugated Human TREM2 Recombinant Rabbit mAb [PSH09-83] - Detector	HA723372	Human	ELISA(Det)
	Biotin Conjugated Human TREM2 Recombinant Rabbit mAb [PSH09-85] - Detector	HA723152B	Human	ELISA(Det),ELISA
	Human TREM2 Recombinant Rabbit mAb [PSH09-84] - BSA and Azide free (Capture)	HA723150	Human	ELISA(Cap)
	Human TREM2 Recombinant Rabbit mAb [PSH09-83] - BSA and Azide free (Capture/Detector)	HA723149	Human	ELISA(Cap),ELISA(Det)
	Biotin Conjugated Human TREM2 Recombinant Rabbit mAb [PSH09-83] - Detector	HA723153B	Human	ELISA(Det),ELISA
	Human TREM2 Recombinant Rabbit mAb [PSH09-85] - BSA and Azide free (Detector)	HA723151	Human	ELISA(Det)
	Iba1 Recombinant Antibody [JM36-62] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610214	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,WB
	Iba1 Rabbit pAb	ER1802-20	Human, Mouse	IHC-P
	Iba1 Recombinant Mouse mAb [8G8]	RT1316	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,FC,IP,IHC-Fr,IF-Cell,IF-Tissue
	Iba1 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS158RT	Mouse	mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Iba1 Recombinant Antibody [JM36-62] - Rat IgG1 (Chimeric)	HA601367	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P
	Iba1 Recombinant Mouse mAb [8G8] - BSA and Azide free	HA610321	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,FC,IP,IHC-Fr,IF-Cell,IF-Tissue
	iFluor™ 594 Conjugated Iba1 Recombinant Rabbit mAb [JM36-62]	HA723821F	Mouse	IHC-Fr
	Iba1 Recombinant Antibody [JM36-62] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610213	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P
	Iba1 Recombinant Antibody [JM36-62] - Guinea pig IgG2 (Chimeric)	HA601368	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,WB
Relative expression	Iba1 Recombinant Antibody [JM36-62] - Chicken IgY (Chimeric)	HA601376	Human, Mouse, Rat, Pig, Cynomolgus monkey	IHC-Fr,WB,FC
	Iba1 Recombinant Rabbit mAb	IRS068	Mouse	mIHC
Relative expression	Iba1 Recombinant Rabbit mAb [JM36-62]	ET1705-78	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue,IF-Cell,FC,IP,mIHC
Relative expression	Iba1 Recombinant Rabbit mAb [JM36-62] - BSA and Azide free	HA750440	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue,IF-Cell,FC,IP
	iFluor™ 488 Conjugated Iba1 Recombinant Rabbit mAb [JM36-62]	HA720158F	Human, Mouse, Rat	IF-Cell,FC
Relative expression	Iba1 Recombinant Antibody [JM36-62] - Chicken IgY (Chimeric) - BSA and Azide free	HA610221	Human, Mouse, Rat, Pig, Cynomolgus monkey	IHC-Fr,WB,FC
Cell treatment	IL-10 Recombinant Rabbit mAb [JE59-71]	HA722032	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	IL-10 Rabbit pAb	ER65464	Human, Mouse, Rat	WB,IHC-P
	IL-10 Rabbit pAb	ER1911-19	Human, Mouse	WB,IHC-P,IF-Tissue
	IL-12 Rabbit pAb	ER1911-22	Rat	WB,IHC-P,IF-Tissue
	IL12 Rabbit pAb	ER1911-21	Mouse, Rat	WB,IHC-P,IF-Tissue,FC
	IL-12A Recombinant Rabbit mAb [PSH16-10] - BSA and Azide free	HA723815	Human	ELISA
	IL12B Rabbit pAb	ER1803-31	Human, Mouse	WB
	IL-1 alpha Recombinant Rabbit mAb [PSH15-66]	HA723773	Mouse	IF-Cell,FC
Cell treatment	IL-1 alpha Recombinant Rabbit mAb [PSH15-65] - BSA and Azide free	HA751576	Mouse	WB
Cell treatment	IL-1 alpha Recombinant Rabbit mAb [PSH15-65]	HA723772	Mouse	WB
	IL-1 alpha Rabbit pAb	ER1802-36	Human, Mouse	IHC-P,FC,WB
	IL-1 alpha Recombinant Rabbit mAb [PSH15-66] - BSA and Azide free	HA751577	Mouse	IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	IL-1 beta Recombinant Rabbit mAb [PSH18-00]	HA723965	Mouse, Rat	WB
Cell treatment	IL-1 beta Recombinant Mouse mAb [A7F7-R] - BSA and Azide free	HA610004	Human, Mouse, Rat	WB,ELISA,IHC-P,IF-Tissue
Cell treatment	IL-1 beta Recombinant Rabbit mAb [JJ087-3]	ET1701-39	Human	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	IL-1 beta Recombinant Mouse mAb [A7F7-R]	HA601036	Human, Mouse, Rat	WB,ELISA,IHC-P,IF-Tissue
Cell treatment	IL-1 beta Mouse mAb [A7F9]	HA601002	Human, Mouse, Rat	WB,ELISA,IF-Cell
Cell treatment	IL-1 beta Mouse mAb [A7A10]	HA601035	Human	WB,ELISA,IF-Cell
	IL-23 Recombinant Rabbit mAb [JE38-99]	HA721037	Human	WB,IHC-P,IF-Tissue
	IL2RA/CD25 Rabbit pAb	ER1911-34	Human, Mouse	WB,FC
	IL-2 receptor alpha Rabbit pAb	HA500076	Human, Rat	WB
	IL33 Rabbit pAb	ER1911-36	Human, Rat	WB,IHC-P,IF-Tissue
	IL-4 Rabbit pAb	ER1706-55	Human	WB,IF-Cell,IHC-P
	IL-4 Mouse mAb [A7D7]	HA601034	Human	ELISA
	IL-4 Recombinant Rabbit mAb [JA86-10]	ET1704-28	Human	WB,IF-Cell,IF-Tissue,IHC-P
	IL-6 Recombinant Rabbit mAb	IRS160RB	Human	mIHC
	IL-6 Rabbit pAb	R1412-2	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	IL-6 Mouse mAb	EM30301	Human	WB
	IL-6 Mouse mAb [1-6]	EM1701-45	Mouse, Rat	WB,IF-Cell
	IL-6 Rabbit pAb	R1203-2	Mouse	WB,ELISA
	IL-6 Mouse mAb [3-A4]	M1303-3	Mouse	WB,IF-Cell,FC,ELISA
Cell treatment	IL-8 Recombinant Rabbit mAb [PSH10-96]	HA723263	Human	WB,FC
Cell treatment	IL-8 Recombinant Rabbit mAb [PSH10-96] - BSA and Azide free	HA751386	Human	WB,FC
	IL-8 Rabbit pAb	ER1706-67	Human	ELISA,IF-Cell,IHC-P
Cell treatment	IL-8 Rabbit pAb	ER1901-61	Human	WB,IHC-P,IF-Cell,FC
Cell treatment	IL-8 Rabbit pAb	R1511-15	Human	WB,IHC-P,IF-Cell,FC
Cell treatment	IL-8 Recombinant Rabbit mAb [PSH10-97]	HA723264	Human	WB
Cell treatment	IL-8 Recombinant Rabbit mAb [PSH10-97] - BSA and Azide free	HA751387	Human	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	iNOS Recombinant Rabbit mAb [PSH03-69] - BSA and Azide free	HA750903	Human, Mouse	WB,IHC-P,IF-Cell
Cell treatment	iNOS Recombinant Rabbit mAb [PSH03-69]	HA722031	Human, Mouse	WB,IHC-P,IF-Cell
Cell treatment	iNOS Rabbit pAb	ER1706-89	Human, Mouse, Rat, Chicken	IF-Cell,FC,WB,IHC-Fr
Cell treatment	iNOS Recombinant Rabbit mAb [JE59-61]	HA721583	Human, Mouse	WB,IF-Cell
Cell treatment	IP10 Recombinant Rabbit mAb [JA10-82] - BSA and Azide free	HA750411	Human	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	IP10 Recombinant Rabbit mAb [JA10-82]	ET1704-27	Human	WB,IF-Cell,IF-Tissue,IHC-P
Relative expression	Liver Arginase Recombinant Rabbit mAb [JB21-49]	ET7107-90	Human	WB,IHC-P,IP
Cell treatment	Liver Arginase Recombinant Rabbit mAb [SY09-06]	ET1605-8	Human, Mouse, Rat	WB
Relative expression	Liver Arginase Recombinant Rabbit mAb [PSH09-78]	HA723143	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Liver Arginase Recombinant Rabbit mAb [PD00-05]	HA721147	Human	WB,IHC-P,IF-Tissue
Cell treatment	Mannose Receptor(CD206) Recombinant Rabbit mAb [PSH07-76]	HA722892	Mouse, Rat	WB,IHC-P,IF-Tissue,IHC-Fr,IF-Cell,IP,mIHC
Cell treatment	Mannose Receptor(CD206) Recombinant Rabbit mAb [PSH07-76] - BSA and Azide free	HA751188	Mouse, Rat	WB,IHC-P,IF-Tissue,IHC-Fr,IF-Cell,IP
	FITC Conjugated Mannose Receptor(CD206) Recombinant Rabbit mAb [JF0953]	HA720214F	Human	FC
	Mannose Receptor(CD206) Recombinant Antibody [PSH07-76] - Guinea pig IgG2 (Chimeric)	HA601452	Mouse, Rat	IHC-Fr,WB
KD	Mannose Receptor(CD206) Recombinant Rabbit mAb [JF0953] - BSA and Azide free	HA750332	Human, Mouse	WB,IF-Cell,FC
	Mannose Receptor(CD206) Recombinant Antibody [PSH07-76] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610287	Mouse, Rat	IHC-Fr,IF-Tissue,IHC-P,WB,IF-Cell
	Mannose Receptor(CD206) Recombinant Rabbit mAb	IRS049	Mouse	mIHC
KD	Mannose Receptor(CD206) Recombinant Rabbit mAb [JF0953]	ET1702-04	Human, Mouse	WB,IF-Cell,FC
	Mannose Receptor(CD206) Recombinant Antibody [PSH07-76] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610286	Mouse, Rat	IHC-Fr,WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mannose Receptor(CD206) Recombinant Antibody [PSH07-76] - Mouse IgG1 (Chimeric)	HA601453	Mouse, Rat	IHC-Fr,IF-Tissue,IHC- P,WB,IF-Cell
	Biotin Conjugated Mouse CCL2 / MCP-1 Recombinant Rabbit mAb [PSH12-45] - Detector	HA725042B	Mouse	ELISA(Det),ELISA
	Mouse CCL2 / MCP-1 Recombinant Rabbit mAb [PSH12-44] - BSA and Azide free (Capture)	HA725040	Mouse	ELISA(Cap)
	Mouse CCL2 / MCP-1 Recombinant Rabbit mAb [PSH12-45] - BSA and Azide free (Detector)	HA725041	Mouse	ELISA(Det)
	Mouse CCL3 Recombinant Rabbit mAb [PSH08-40] - BSA and Azide free (Detector)	HA722983	Mouse	ELISA(Det)
	Mouse CCL3 Recombinant Rabbit mAb [PSH08-39] - BSA and Azide free (Capture)	HA722982	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse CCL3 Recombinant Rabbit mAb [PSH08- 40] - Detector	HA722984B	Mouse	ELISA(Det),ELISA
	Mouse CCL5 Recombinant Rat mAb	HA601395	Human, Mouse	ELISA
	Mouse CCL5/RANTES Recombinant Rabbit mAb [PSH14-15] - BSA and Azide free (Capture)	HA725046	Mouse	ELISA(Cap)
	Mouse CCL5/RANTES Recombinant Rabbit mAb [PSH14-16] - BSA and Azide free (Detector)	HA725047	Mouse	ELISA(Det)
	Biotin Conjugated Mouse CCL5/RANTES Recombinant Rabbit mAb [PSH14-16] - Detector	HA725048B	Mouse	ELISA(Det),ELISA
	PE Conjugated Mouse CD16/32 Recombinant Rabbit mAb [PSH06- 60]	HA720216F	Mouse	FC
	Biotin Conjugated Mouse CD68 Recombinant Rabbit mAb [PSH12- 25] - Detector	HA723448B	Mouse	ELISA(Det),ELISA
	Mouse CD68 Recombinant Rabbit mAb [PSH12-27] - BSA and Azide free (Detector)	HA723450	Mouse	ELISA(Det)
	Mouse CD68 Recombinant Rabbit mAb [PSH12-26] - BSA and Azide free (Capture)	HA723449	Mouse	ELISA(Cap)
	Mouse CD68 Recombinant Rabbit mAb [PSH12-24] - BSA and Azide free (Capture)	HA723446	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse CD68 Recombinant Rabbit mAb [PSH12- 27] - Detector	HA723451B	Mouse	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse CD68 Recombinant Rabbit mAb [PSH12-25] - BSA and Azide free (Detector)	HA723447	Mouse	ELISA(Det)
	Mouse CD86 Recombinant Rabbit mAb [PSH13-65] - BSA and Azide free (Capture)	HA723555	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse CD86 Recombinant Rabbit mAb [PSH13-64] - Detector	HA723554B	Mouse	ELISA(Det),ELISA
	Mouse CD86 Recombinant Rabbit mAb [PSH13-64] - BSA and Azide free (Detector)	HA723553	Mouse	ELISA(Det)
	Biotin Conjugated Mouse CD86 Recombinant Rabbit mAb [PSH13-66] - Detector	HA723557B	Mouse	ELISA(Det),ELISA
	Mouse CD86 Recombinant Rabbit mAb [PSH13-66] - BSA and Azide free (Detector)	HA723556	Mouse	ELISA(Det)
	Mouse CD86 Recombinant Rabbit mAb [PSH13-63] - BSA and Azide free (Capture)	HA723552	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse CXCL16 Recombinant Rabbit mAb [PSH17-74] - Detector	HA723943B	Mouse	ELISA(Det),ELISA
	Mouse CXCL16 Recombinant Rabbit mAb [PSH17-72] - BSA and Azide free (Detector)	HA723941	Mouse	ELISA(Det)
	Mouse CXCL16 Recombinant Rabbit mAb [PSH17-71] - BSA and Azide free (Capture)	HA723940	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse CXCL16 Recombinant Rabbit mAb [PSH17-72] - Detector	HA723941B	Mouse	ELISA(Det),ELISA
	Mouse CXCL16 Recombinant Rabbit mAb [PSH17-74] - BSA and Azide free (Detector)	HA723943	Mouse	ELISA(Det)
	Mouse CXCL16 Recombinant Rabbit mAb [PSH17-73] - BSA and Azide free (Capture)	HA723942	Mouse	ELISA(Cap)
	Mouse IL-10 Recombinant Rabbit mAb [PS01-93] - BSA and Azide free (Detector)	HA721518	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-10 Recombinant Rabbit mAb [PS01-93] - Detector	HA721519	Mouse	ELISA(Det)
	Mouse IL-10 Recombinant Rabbit mAb [PS01-85] - BSA and Azide free (Capture)	HA721517	Mouse	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Mouse IL12 p40 Recombinant Rabbit mAb [PSH04-88]	HA722603B	Mouse	ELISA(Det)
	Mouse IL12 p40 Recombinant Rabbit mAb [PSH04-87] - BSA and Azide free (Capture)	HA722193	Mouse	ELISA(Cap)
	Mouse IL12 p40 Recombinant Rabbit mAb [PSH04-88] - BSA and Azide free (Detector)	HA722194	Mouse	ELISA(Det)
	Mouse IL-1 alpha Recombinant Rabbit mAb [PSH13-41] - BSA and Azide free (Detector)	HA723527	Mouse	ELISA(Det)
	Mouse IL-1 alpha Recombinant Rabbit mAb [PSH13-40] - BSA and Azide free (Capture)	HA723526	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-1 alpha Recombinant Rabbit mAb [PSH13-41] - Detector	HA723528B	Mouse	ELISA(Det),ELISA
	Biotin Conjugated Mouse IL-1 alpha Recombinant Rabbit mAb [PSH13-43] - Detector	HA723531B	Mouse	ELISA(Det),ELISA
	Mouse IL-1 alpha Recombinant Rabbit mAb [PSH13-42] - BSA and Azide free (Capture)	HA723529	Mouse	ELISA(Cap)
	Mouse IL-1 alpha Recombinant Rabbit mAb [PSH13-43] - BSA and Azide free (Detector)	HA723530	Mouse	ELISA(Det)
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH03-31] - BSA and Azide free (Capture)	HA721976	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-1 beta Recombinant Rabbit mAb [PSH04-67]	HA722608B	Mouse	ELISA(Det)
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH04-66] - BSA and Azide free (Capture)	HA722166	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-1 beta Recombinant Rabbit mAb [PSH03-32] - Detector	HA722625B	Mouse	ELISA(Det),ELISA
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH09-59] - BSA and Azide free	HA723121	Mouse	ELISA
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH04-67] - BSA and Azide free (Detector)	HA722167	Mouse	ELISA(Det)
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH09-58] - BSA and Azide free	HA723120	Mouse	ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse IL-1 beta Recombinant Rabbit mAb [PSH03-32] - BSA and Azide free (Detector)	HA721977	Mouse	ELISA(Det)
	Mouse IL-6 Recombinant Rabbit mAb [PS01-83] - BSA and Azide free (Detector)	HA721601	Mouse	ELISA(Det)
	Mouse IL-6 Recombinant Rabbit mAb [PSH0-58] - BSA and Azide free (Capture)	HA721411	Mouse	ELISA(Cap)
	Mouse IL-6 Recombinant Rabbit mAb [PSH0-59] - BSA and Azide free (Detector)	HA721412	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PS01-74] - Detector	HA723643B	Mouse	ELISA(Det)
	Mouse IL-6 Recombinant Rabbit mAb [PSH0-60] - BSA and Azide free (Detector)	HA721413	Mouse	ELISA(Det)
	Mouse IL-6 Recombinant Rabbit mAb [PS01-74] - BSA and Azide free (Capture/Detector)	HA721600	Mouse	ELISA(Cap),ELISA(Det)
	Mouse IL-6 Recombinant Rabbit mAb [PSH14-47] - BSA and Azide free (Capture)	HA723642	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PSH0-59] - Detector	HA722646B	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PSH0-60] - Detector	HA722647B	Mouse	ELISA(Det)
	HRP Conjugated Mouse IL-6 Mouse mAb [3-A4]	M1303-2	Mouse	WB,ELISA
	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PS01-83] - BSA and Azide free (Detector)	HA721602	Mouse	ELISA(Det)
	Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-30] - BSA and Azide free (Detector)	HA721975	Mouse, Rat	ELISA(Det)
	Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-29] - BSA and Azide free (Capture)	HA721974	Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-30]	HA722626B	Mouse, Rat	ELISA(Det)
	Biotin Conjugated Mouse TNF alpha Recombinant Rabbit mAb [PSH03-77]	HA722623B	Mouse	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse TNF alpha Recombinant Rabbit mAb [PSH03-77] - BSA and Azide free (Detector)	HA722041	Mouse	ELISA(Det)
	Mouse TNF alpha Recombinant Rabbit mAb [PSH03-76] - BSA and Azide free (Capture)	HA722040	Mouse	ELISA(Cap)
KD	MRP1 Recombinant Rabbit mAb [JA93-03]	ET1704-56	Human	WB,IHC-P,FC,IF-Cell
	MTCO2 Recombinant Rabbit mAb [SC06-88]	ET1610-72	Human	WB,IF-Cell,IF-Tissue,IHC-P
KD	NF-kB p100 / NFKB2 Mouse mAb [16B3]	EM1901-79	Human, Mouse, Rat	WB,IHC-P,FC
KD	NF-kB p100 / NFKB2 Mouse mAb [16B2]	EM1901-78	Human	WB,IHC-P,FC
KD	NF-kB p100 / NFKB2 Recombinant Rabbit mAb [JM82-03]	ET1705-45	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
	NF-kB p105 / p50 Rabbit pAb	R1309-9	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	NF-kB p105 / p50 Mouse mAb [A2-7]	M1505-8	Human, Mouse	WB,IF-Cell,IHC-P
KD	NF-kB p105 / p50 Recombinant Rabbit mAb [SZ20-01]	ET1603-18	Human, Mouse, Rat	WB,IHC-P,ICC,FC,IF-Tissue
KD	NF-kB p105 / p50 Recombinant Rabbit mAb [SZ20-01] - BSA and Azide free	HA750060	Human, Mouse, Rat	WB,IHC-P,ICC,FC,IF-Tissue
KO, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [PSH0-27]	HA721307	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
	NF-kB p65 Rabbit pAb	HA500402	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
KO, KD, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [SZ10-04] - BSA and Azide free	HA750058	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KO, KD, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [SZ10-04]	ET1603-12	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	NF-kB p65 Rabbit pAb	ER0815	Human, Mouse, Rat, Zebrafish	WB,IHC-P,FC
KO, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [PSH0-27] - BSA and Azide free	HA750601	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
KD	Notch 1 Recombinant Rabbit mAb [SJ205] - BSA and Azide free	HA750111	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	Notch 1 Recombinant Rabbit mAb [SJ205]	ET1606-55	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Relative expression, KD	NQO1 Recombinant Mouse mAb [13A1-R]	HA601209	Human	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression, KD	NQO1 Recombinant Mouse mAb [13A1-R] - BSA and Azide free	HA610087	Human	WB
	NQO1 Rabbit pAb	ER1802-85	Human, Mouse, Rat	WB,FC
Relative expression, KD	NQO1 Mouse mAb [13A1]	EM1801-21	Human, Rat	WB
KD	NQO1 Recombinant Rabbit mAb [JF440-1] - BSA and Azide free	HA750349	Human, Mouse, Rat	WB,IP,FC,IHC-P
KD	NQO1 Recombinant Rabbit mAb [JF440-1]	ET1702-50	Human, Mouse, Rat	WB,IP,FC,IHC-P
Cell treatment	Nrf2 Recombinant Rabbit mAb [PSH11-20]	HA723302	Human, Mouse, Rat	WB,IF-Cell,ChIP
Cell treatment	Nrf2 Rabbit pAb	ER1706-41	Human, Mouse, Rat	IHC-P,IF-Cell,WB
	Nrf2 Rabbit pAb	R1312-8	Human, Mouse	WB,IHC-P,IF-Tissue
Cell treatment	Nrf2 Recombinant Rabbit mAb [PSH11-20] - BSA and Azide free	HA751401	Human, Mouse, Rat	WB,IF-Cell,ChIP
Cell treatment	Nrf2 Recombinant Rabbit mAb [JE30-14]	HA721432	Human, Mouse	WB,IF-Cell
Cell treatment	Phospho-GSK3 (alpha + beta) (Y216 + Y279) Recombinant Rabbit mAb [SY26-05] - BSA and Azide free	HA750127	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,IF-Tissue,FC
Cell treatment	Phospho-GSK3 (alpha + beta) (Y216 + Y279) Recombinant Rabbit mAb [SY26-05]	ET1607-54	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,IF-Tissue,FC
Cell treatment	Phospho-GSK3 alpha (S21) Recombinant Rabbit mAb [PSH09-11] - BSA and Azide free	HA751265	Human, Mouse	WB,IHC-P,Dot Blot
Cell treatment	Phospho-GSK3 alpha (S21) Recombinant Rabbit mAb [PSH09-11]	HA723069	Human, Mouse	WB,IHC-P,Dot Blot
Cell treatment	Phospho-GSK3 beta (S9) Recombinant Rabbit mAb [SY02-71]	ET1607-60	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
Cell treatment	Phospho-GSK3 beta (S9) Recombinant Rabbit mAb [SY02-71] - BSA and Azide free	HA750128	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
Cell treatment	Phospho-NF-kB p65 (S529) Recombinant Rabbit mAb [SP07-00]	ET1604-27	Human	WB,IHC-P,FC,IP
Cell treatment	Phospho-NF-kB p65 (S536) Recombinant Rabbit mAb [PSH10-58]	HA723223	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-NF-kB p65 (S536) Recombinant Rabbit mAb [PSH10-58] - BSA and Azide free	HA751361	Human, Mouse, Rat	WB,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment, Relative expression	Phospho-Nrf2 (S40) Recombinant Rabbit mAb [SU0334]	ET1608-28	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
Cell treatment, Relative expression	Phospho-Nrf2 (S40) Recombinant Rabbit mAb [SU0334] - BSA and Azide free	HA750142	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
Cell treatment	Phospho-SQSTM1 / p62 (S349) Recombinant Rabbit mAb [PSH08-97] - BSA and Azide free	HA751258	Human, Mouse	WB,IF-Cell,IHC-P
Cell treatment	Phospho-SQSTM1 / p62 (S349) Recombinant Rabbit mAb [PSH08-97]	HA723050	Human, Mouse	WB,IF-Cell,IHC-P
	PITSLRE Rabbit pAb	ER64610	Human, Mouse	WB,IHC-P
	RANTES Recombinant Rabbit mAb [JM03-45]	ET1705-70	Human, Mouse	WB,IHC-P,IF-Cell,IF-Tissue
	Biotin Conjugated Rat IL-10 Recombinant Rabbit mAb [PSH02-22]	HA722637B	Rat	ELISA(Det)
	Rat IL-10 Recombinant Rabbit mAb [PSH02-22] - BSA and Azide free (Detector)	HA721787	Rat	ELISA(Det)
	Rat IL-10 Recombinant Rabbit mAb [PSH02-21] - BSA and Azide free (Capture)	HA721786	Rat	ELISA(Cap)
	Rat IL-1 beta Recombinant Rabbit mAb [PSH02-59] - BSA and Azide free (Capture)	HA721836	Rat	ELISA(Cap)
	Rat IL-1 beta Recombinant Rabbit mAb [PS01-87] - BSA and Azide free (Detector)	HA721610	Rat	ELISA(Det)
	Biotin Conjugated Rat IL-1 beta Recombinant Rabbit mAb [PS01-87] - Detector	HA721611	Rat	ELISA(Det),ELISA
	Biotin Conjugated Rat IL-1 beta Recombinant Rabbit mAb [PSH02-60]	HA722634B	Rat	ELISA(Det)
	Rat IL-1 beta Recombinant Rabbit mAb [PS01-86] - BSA and Azide free (Capture)	HA721609	Rat	ELISA(Cap)
	Rat IL-1 beta Recombinant Rabbit mAb [PSH02-60] - BSA and Azide free (Detector)	HA721837	Rat	ELISA(Det)
	RBX1 Recombinant Rabbit mAb [PSH07-01]	HA722779	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IP,IF-Tissue
	RBX1 Rabbit pAb	ER1802-9	Human, Mouse	WB,IF-Cell,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	RBX1 Recombinant Rabbit mAb [JE32-69]	HA722534	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
	RBX1 Recombinant Rabbit mAb [PSH07-01] - BSA and Azide free	HA751128	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IP,IF-Tissue
	ROS1 Recombinant Rabbit mAb [PD01-27]	HA721420	Human	IHC-P
	ROS1 Recombinant Rabbit mAb [PD01-27] - BSA and Azide free	HA750637	Human	IHC-P
	ROS1 Rabbit pAb	ER1916-10	Human, Rat	WB,FC
	S100B Recombinant Rabbit mAb [SC57-02] - BSA and Azide free (Capture)	HA723728	Human, Mouse, Rat	ELISA(Cap)
	S100B Recombinant Rabbit mAb [PSH15-25] - BSA and Azide free (Detector)	HA723729	Human, Mouse, Rat	ELISA(Det)
	Biotin Conjugated S100B Recombinant Rabbit mAb [PSH15-25] - Detector	HA723730B	Human, Mouse, Rat	ELISA(Det),ELISA
	S100 beta Recombinant Rabbit mAb [SC57-02]	ET1610-3	Human, Mouse, Rat, Zebrafish, Goat, Cynomolgus monkey, Green monkey, Pig, Rabbit	WB,IF-Tissue,IHC-P,IHC-Fr,IF-Cell,FC,IP
	S100 beta Recombinant Rabbit mAb [SC57-02] - BSA and Azide free	HA750209	Human, Mouse, Rat, Zebrafish, Goat, Cynomolgus monkey, Green monkey, Pig, Rabbit	WB,IF-Tissue,IHC-P,IHC-Fr,IF-Cell,FC,IP
	S100 beta Recombinant Antibody [SC57-02] - Rat IgG1 (Chimeric)	HA601552	Mouse, Rat	IHC-Fr,IHC-P
	S100 beta Recombinant Rabbit mAb [PSH17-13] - BSA and Azide free	HA751617	Human, Mouse, Rat	IHC-Fr,IHC-P,WB,IF-Tissue
	S100 beta Recombinant Antibody [SC57-02] - Guinea pig IgG2 (Chimeric)	HA601550	Mouse	IHC-Fr
Relative expression	S100 beta Recombinant Rabbit mAb [PSH17-13]	HA723900	Human, Mouse, Rat	IHC-Fr,IHC-P,WB,IF-Tissue
	S100 beta Rabbit pAb	RT1547	Human, Mouse, Rat	WB,IP,IF
Relative expression	S100 beta Recombinant Mouse mAb [PD00-11]	HA601098	Human, Mouse, Rat	IHC-P,WB,IF-Tissue
	iFluor™ 594 Conjugated S100 beta Recombinant Rabbit mAb [SC57-02]	HA720125F	Human	IF-Tissue
	S100 beta Recombinant Antibody [SC57-02] - Mouse IgG1 (Chimeric)	HA601551	Mouse, Rat	IHC-Fr,IHC-P
	S100 beta Recombinant Antibody [SC57-02] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610366	Mouse, Rat	IHC-Fr,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	S100 beta Recombinant Antibody [SC57-02] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610365	Mouse	IHC-Fr
	S100 beta Recombinant Antibody [SC57-02] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610367	Mouse, Rat	IHC-Fr,IHC-P
	SQSTM1 / p62 Rabbit pAb	R1309-8	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	SQSTM1 / p62 Mouse mAb [10-E10]	EM0704	Human, Mouse, Rat	IF-Cell,IHC-P,FC
Cell treatment, KD	SQSTM1 / p62 Recombinant Rabbit mAb [PS00-61]	HA721171	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Cell,FC,IF-Tissue,IHC-Fr
Cell treatment, KD	SQSTM1 / p62 Recombinant Rabbit mAb [PS00-61] - BSA and Azide free	HA750510	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Cell,FC,IF-Tissue,IHC-Fr
	Superoxide Dismutase 3 Recombinant Rabbit mAb [JE31-84]	HA721408	Human	WB
	SYVN1 Rabbit pAb	ER1917-45	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
	SYVN1 Recombinant Rabbit mAb [JE44-67]	HA721874	Human	WB,IF-Cell,FC
	TGF beta 1 Recombinant Rabbit mAb	IRS205RB	Human	mIHC
	Biotin Conjugated TGF beta 1 Recombinant Rabbit mAb [PSH15-22] - Detector	HA723727B	Human, Mouse, Rat	ELISA(Det),ELISA
	TGF beta 1 Mouse mAb [B2-F6]	EM1709-93	Human	WB,IHC-P,FC
Relative expression	TGF beta 1 Recombinant Rabbit mAb [PD00-17] - BSA and Azide free	HA750485	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Tissue
	TGF beta 1 Recombinant Rabbit mAb [PSH18-59]	HA724012	Human, Mouse, Monkey	WB,IF-Cell,IHC-P,FC
Relative expression	TGF beta 1 Recombinant Rabbit mAb [PD00-17]	HA721143	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Tissue
	TGF beta 1 Recombinant Rabbit mAb [PSH15-22] - BSA and Azide free (Detector)	HA723726	Human, Mouse, Rat	ELISA(Det)
	TGF beta 1 Recombinant Rabbit mAb [PSH15-21] - BSA and Azide free (Capture)	HA723725	Human, Mouse, Rat	ELISA(Cap)
	TGF beta 1 Rabbit pAb	HA500496	Human	WB
KO	TLR2 Recombinant Rabbit mAb [JM22-41] - BSA and Azide free	HA750447	Human, Mouse, Rat	WB,FC,IF-Cell
KO	TLR2 Recombinant Rabbit mAb [JM22-41]	ET1705-92	Human, Mouse, Rat	WB,FC,IF-Cell
	TLR4 Mouse mAb	RT1666	Human, Mouse, Rat	WB
	TLR4 Rabbit pAb	ER1706-43	Human, Mouse, Rat	IF-Cell,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	TMEM119 Recombinant Rabbit mAb [PSH12-98]	HA723503	Human, Mouse, Rat	WB,IF-Cell,FC,IP
	TMEM119 Rabbit pAb	HA500260	Human, Mouse	WB,IF-Cell,FC
	TMEM119 Recombinant Rabbit mAb [PSH12-98] - BSA and Azide free	HA751453	Human, Mouse, Rat	WB,IF-Cell,FC,IP
	TMEM119 Rabbit pAb	ER2001-37	Human	WB
Cell treatment	TNF alpha Recombinant Rabbit mAb [PS01-07] - BSA and Azide free	HA750897	Human, Mouse	WB
	TNF alpha Rabbit pAb	ER65189	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Cell treatment	TNF alpha Rabbit pAb	R1203-1	Mouse	WB,IF-Cell,ELISA
Cell treatment	TNF alpha Recombinant Rabbit mAb [PS01-07]	HA722022	Human, Mouse	WB
	TNF alpha Rabbit pAb	ER1919-22	Human, Mouse	WB
	TNFR1 Mouse mAb [2G2]	RT1624	Human, Mouse, Rat	WB,IP,IF,IHC-P
	TNF Receptor I Rabbit pAb	HA500140	Human	WB,IHC-P,FC
	Total OXPHOS Antibody Cocktail Recombinant Rabbit Multiclonal Antibody [PSH15-50]	HA723756	Human	WB,IHC-P
	TREM2 Rabbit pAb	ER2001-59	Human, Rat	WB,FC
Relative expression	TREM2 Rabbit pAb	ER1902-96	Human	WB,IHC-P,FC
	TREM2 Rabbit pAb	ER1918-04	Human, Mouse	WB
	TREM2 Rabbit pAb	ER2001-05	Human	WB,IHC-P,FC
	TXNRD1 Recombinant Rabbit mAb [JA11-32]	ET1704-81	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层