

流式抗体-产品清单

高级验证	产品名	货号	反应种属	验证应用
	FITC Conjugated Bcl-2 Recombinant Rabbit mAb [JF104-8]	HA720210F	Human, Mouse	FC
	BCMA Recombinant Rabbit mAb [PSH16-80]	HA723880	Human	FC
Relative expression	BCMA Recombinant Rabbit mAb [PSH17-94]	HA723959	Human	FC
Relative expression	CD105 Recombinant Rabbit mAb [PSH15-96]	HA723801	Human	FC
Relative expression	CD105 Recombinant Rabbit mAb [PSH15-97]	HA723802	Human	FC
	PE Conjugated CD123 / IL3RA Recombinant Mouse mAb [PSH14-51]	HA601528F	Human	FC
	CD123 / IL3RA Recombinant Mouse mAb [PSH14-51]	HA601490	Human	FC
	APC Conjugated CD127 Recombinant Mouse mAb [PSH12-76]	HA601536F	Human	FC
	CD127 Recombinant Mouse mAb [PSH12-76]	HA601429	Human	FC
	FITC Conjugated CD15 Recombinant Rabbit mAb [PD00-42]	HA720195F	Human	FC
	FITC Conjugated CD19 Recombinant Mouse mAb [PSH14-97]	HA601518F	Human	FC
	PE Conjugated CD19 Recombinant Rabbit mAb [PSH13-44]	HA723763F	Human	FC
Relative expression	APC Conjugated CD19 Recombinant Rabbit mAb [PSH13-44]	HA723761F	Human	FC
	FITC Conjugated CD19 Recombinant Rabbit mAb [PSH13-44]	HA723762F	Human	FC
	APC Conjugated CD20 Recombinant Rabbit mAb [PSH08-12]	HA723838F	Human	FC
	PE Conjugated CD20 Recombinant Rabbit mAb [PSH08-12]	HA723839F	Human	FC
	APC Conjugated CD25 / IL-2 Receptor alpha Recombinant Mouse mAb [PSH14-46]	HA601530F	Human	FC
	CD25 / IL-2 Receptor alpha Recombinant Mouse mAb [PSH14-46]	HA601486	Human, Mouse	FC,IHC-P,IF-Tissue
	APC Conjugated CD27 Recombinant Rabbit mAb [PSH08-90]	HA723828F	Human	FC
	PE Conjugated CD45RO Recombinant Mouse mAb [PSH13-94]	HA601531F	Human	FC

	iFluor™ 488 Conjugated CD47 Mouse mAb [1B1]	HA600104F	Human	FC
Relative expression	CD62L Recombinant Rabbit mAb [PSH14-92]	HA723694	Human	FC,IF-Cell
Relative expression	CD62L Recombinant Rabbit mAb [PSH15-59]	HA723766	Human	FC
Relative expression	CD81 Recombinant Rabbit mAb [PSH10-26]	HA723193	Human, Rat	WB,IHC-P,FC
	PE Conjugated CD94 Recombinant Mouse mAb [PSH14-52]	HA601539F	Human	FC
	CD94 Recombinant Mouse mAb [PSH14-52]	HA601491	Human	FC
	APC Conjugated CD94 Recombinant Mouse mAb [PSH14-52]	HA601529F	Human	FC
	Glycophorin A / CD235a Recombinant Mouse mAb [PSH17-08]	HA601553	Human	FC
	HER3 Recombinant Rabbit mAb [R5]	HA721126	Human	FC
	HER3 Recombinant Rabbit mAb [R4]	HA721135	Human	FC
	iFluor™ 647 Conjugated HER3 Recombinant Rabbit mAb [R5]	HA720181F	Human	FC
	FITC Conjugated Human CD10 Recombinant Rabbit mAb [SN75-07]	HA720209F	Human	FC
	FITC Conjugated Human CD11b Recombinant Mouse mAb [A2-9-R]	HA600130F	Human	FC
	FITC Conjugated Human CD14 Recombinant Rabbit mAb [SC69-02]	HA720194F	Human	FC
	iFluor™ 647 Conjugated Human CD14 Recombinant Rabbit mAb [SC69-02]	HA720184F	Human	FC
	Human CD16 Recombinant Mouse mAb [PSH04-95]	HA601340	Human	FC
	PE Conjugated Human CD16 Recombinant Mouse mAb [PSH04-95]	HA600129F	Human	FC
	FITC Conjugated Human CD200 Recombinant Rabbit mAb [PSH01-48]	HA720208F	Human	FC
	FITC Conjugated Human CD200R1 Recombinant Rabbit mAb [PSH05-02]	HA720212F	Human	FC
Cell treatment	Human CD200R1 Recombinant Rabbit mAb [PSH05-02]	HA722211	Human	IF-Cell,FC
	PE Conjugated Human CD3 Recombinant Mouse mAb [PSH04-92]	HA600117F	Human	FC
	Human CD3 Recombinant Mouse mAb [PSH04-92]	HA601337	Human	FC
	APC Conjugated Human CD3 Recombinant Mouse mAb [PSH04-92]	HA600118F	Human	FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Percp-Cy5.5 Conjugated Human CD3 Recombinant Mouse mAb [PSH04-92]	HA601548F	Human	FC
	FITC Conjugated Human CD3 Recombinant Mouse mAb [PSH04-92]	HA600116F	Human	FC
	FITC Conjugated Human CD33 Recombinant Rabbit mAb [PSH04-96]	HA720197F	Human	FC
	FITC Conjugated Human CD38 Recombinant Rabbit mAb [SN07-20]	HA720196F	Human	FC
	PE-Cy7 Conjugated Human CD4 Recombinant Mouse mAb [PSH04-93]	HA601542F	Human	FC
	Human CD4 Recombinant Mouse mAb [PSH04-93]	HA601338	Human	FC
	FITC Conjugated Human CD4 Recombinant Mouse mAb [PSH04-93]	HA600119F	Human	FC
	PE Conjugated Human CD4 Recombinant Mouse mAb [PSH04-93]	HA600120F	Human	FC
	APC Conjugated Human CD4 Recombinant Mouse mAb [PSH04-93]	HA600121F	Human	FC
	APC Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600127F	Human	FC
	iFluor™ 488 Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600109F	Human	FC
	iFluor™ 555 Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA601556F	Human	FC
	FITC Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600125F	Human	FC
	PE Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600126F	Human	FC
	Percp-Cy5.5 Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600128F	Human	FC
	FITC Conjugated Human CD45RA Recombinant Mouse mAb [PSH04-97]	HA600131F	Human	FC
	PE Conjugated Human CD45RA Recombinant Mouse mAb [PSH04-97]	HA601525F	Human	FC
	Human CD45RA Recombinant Mouse mAb [PSH04-97]	HA601341	Human	FC,IF-Cell,IHC-P
	FITC Conjugated Human CD5 Recombinant Rabbit mAb [PSH05-01]	HA720211F	Human	FC
Relative expression	Human CD5 Recombinant Rabbit mAb [PSH05-01]	HA722210	Human	IF-Cell,FC
	FITC Conjugated Human CD58 Recombinant Rabbit mAb [PSH02-74]	HA720198F	Human	FC
	FITC Conjugated Human CD68 Recombinant Mouse mAb [A3C3-R]	HA600132F	Human	FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human CD8 alpha Recombinant Mouse mAb [PSH04-94]	HA601339	Human	FC
	APC Conjugated Human CD8 alpha Recombinant Mouse mAb [PSH04-94]	HA600124F	Human	FC
	FITC Conjugated Human CD8 alpha Recombinant Mouse mAb [PSH04-94]	HA600122F	Human	FC
	APC-Cy7 Conjugated Human CD8 alpha Recombinant Mouse mAb [PSH04-94]	HA601541F	Human	FC
	PE Conjugated Human CD8 alpha Recombinant Mouse mAb [PSH04-94]	HA600123F	Human	FC
	FITC Conjugated Human IgD Recombinant Mouse mAb [PSH14-48]	HA601526F	Human	FC
	Human IgD Recombinant Mouse mAb [PSH14-48]	HA601487	Human	FC
	Human IgE Recombinant Mouse mAb [PSH14-50]	HA601489	Human	FC
	PE Conjugated Human IgE Recombinant Mouse mAb [PSH14-50]	HA601527F	Human	FC
	Human IgM Recombinant Mouse mAb [PSH14-49]	HA601488	Human	FC
	PE Conjugated Human NCAM1 / CD56 Recombinant Rabbit mAb [PSH06-53]	HA720215F	Human	FC
	APC Conjugated Human NCAM1 / CD56 Recombinant Rabbit mAb [PSH06-53]	HA723757F	Human	FC
	Integrin alpha 1 Recombinant Rabbit mAb [PSH14-04]	HA723594	Human	FC
	PE Conjugated Ly6g Recombinant Rabbit mAb [PSH14-95]	HA723853F	Mouse	FC
	Ly6g Recombinant Rabbit mAb [PSH14-95]	HA723696	Mouse	FC,IHC-P
	PE Conjugated Mouse CD16/32 Recombinant Rabbit mAb [PSH06-60]	HA720216F	Mouse	FC
	APC Conjugated Anti-Mouse CD3 Recombinant Antibody [PSH04-98] - Rabbit IgG (Chimeric)	HA720201F	Mouse	FC
	Anti-Mouse CD3 Recombinant Antibody [PSH04-98] - Rabbit IgG (Chimeric)	HA723347	Mouse	FC
	PE Conjugated Anti-Mouse CD3 Recombinant Antibody [PSH04-98] - Rabbit IgG (Chimeric)	HA720200F	Mouse	FC
	FITC Conjugated Anti-Mouse CD3 Recombinant Antibody [PSH04-98] - Rabbit IgG (Chimeric)	HA720199F	Mouse	FC
	APC Conjugated Anti-Mouse CD4 Recombinant Antibody [PSH04-99] - Rabbit IgG (Chimeric)	HA720204F	Mouse	FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	FITC Conjugated Anti-Mouse CD4 Recombinant Antibody [PSH04-99] - Rabbit IgG (Chimeric)	HA720202F	Mouse	FC
	PE Conjugated Anti-Mouse CD4 Recombinant Antibody [PSH04-99] - Rabbit IgG (Chimeric)	HA720203F	Mouse	FC
	FITC Conjugated Anti-Mouse CD8 alpha Recombinant Antibody [PSH05-00] - Rabbit IgG (Chimeric)	HA720205F	Mouse	FC
	APC Conjugated Anti-Mouse CD8 alpha Recombinant Antibody [PSH05-00] - Rabbit IgG (Chimeric)	HA720207F	Mouse	FC
	Mouse CD8 alpha Recombinant Rat mAb [PSH16-78]	HA601543	Mouse	FC
	PE Conjugated Anti-Mouse CD8 alpha Recombinant Antibody [PSH05-00] - Rabbit IgG (Chimeric)	HA720206F	Mouse	FC
Cell treatment	Mouse PD-L1 Recombinant Rabbit mAb [PSH05-03]	HA722208	Mouse	FC
Cell treatment	Mouse PD-L1 Recombinant Rabbit mAb [PSH05-04]	HA722209	Mouse	FC
Cell treatment	FITC Conjugated Mouse PD-L1 Recombinant Rabbit mAb [PSH05-04]	HA720213F	Mouse	FC
	iFluor™ 647 Conjugated MUC1 Mouse mAb [PSH0-48]	HA600107F	Human	FC
	iFluor™ 647 Conjugated MUC1 Mouse mAb [PSH0-49]	HA600108F	Human	FC
	FITC Conjugated NK1.1 Recombinant Rabbit mAb [PSH08-46]	HA720219F	Mouse	FC
	APC Conjugated Recombinant Mouse IgG1 [PSH04-91] - Isotype Control	HA600115F		FC
	PE Conjugated Recombinant Mouse IgG1 [PSH04-91] - Isotype Control	HA600114F		FC
	FITC Conjugated Recombinant Mouse IgG1 [PSH04-91] - Isotype Control	HA600113F		FC
	PE Conjugated Recombinant Rabbit IgG [PSH04-42] - Isotype control	HA720192F		FC
	FITC Conjugated Recombinant Rabbit IgG [PSH04-42] - Isotype control	HA720191F		FC
	APC Conjugated Recombinant Rabbit IgG [PSH04-42] - Isotype control	HA720193F		FC
	FITC Conjugated Siglec 8 Recombinant Mouse mAb [PSH14-96]	HA601517F	Human	FC
	TCR beta chain Recombinant Mouse mAb [PSH16-30]	HA601535	Mouse	FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	TIM 3 Recombinant Rabbit mAb [PSH10-00]	HA723169	Human	WB,IF-Cell,FC,IP
	TRBC1 Recombinant Mouse mAb [PSH15-51]	HA601524	Human	FC
	iFluor™ 488 Conjugated VISTA Recombinant Rabbit mAb [PSH0-65]	HA720183F	Human	FC