

免疫细胞类型标记物-产品清单

高级验证	产品名	货号	反应种属	验证应用
	2B4 Rabbit pAb	HA500180	Human	IHC-P
	2B4 Rabbit pAb	HA500263	Human, Rat	WB
Cell treatment	4-1BBL Recombinant Rabbit mAb [JE33-91]	HA721357	Human, Mouse	WB
	ALCAM Mouse mAb [C7-E12]	EM1708-92	Human	WB,IF-Cell,IHC-P,FC
	ALCAM/CD166 Mouse mAb [7C2]	RT1029	Human, Mouse, Rat	WB,IP,IF,IHC-P
	Angiotensinogen Recombinant Rabbit mAb [SD201-02]	ET1612-48	Human	WB,IF-Cell,IF-Tissue,IP
	AP-A Recombinant Rabbit mAb [JE46-03]	ET7109-45	Human	WB,IHC-P
	Arginase1 Recombinant Rabbit mAb	IRS082RB	Mouse	mIHC
	Arginase1 Recombinant Rabbit mAb	IRS081RB	Human	mIHC
	B3GAT1 Rabbit pAb	HA500443	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	B3GAT1 Rabbit pAb	ER1901-15	Human, Mouse, Rat	WB,IHC-P,FC
	B3GAT1 Rabbit pAb	ER1901-83	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	B7-1/CD80 Rabbit pAb	ER1904-08	Mouse, Rat	WB
Relative expression	B7H4 Recombinant Rabbit mAb [PSH03-65]	HA722024	Human, Mouse, Rat	WB,IF-Cell
Relative expression	B7H4 Recombinant Rabbit mAb [PSH03-65] - BSA and Azide free	HA750899	Human, Mouse, Rat	WB,IF-Cell
	B7-H4(VTCN1) Rabbit pAb	R1512-1	Human, Mouse	WB,IF-Cell,IHC-P,FC
	B7-H4(VTCN1) Rabbit pAb	ER1802-56	Human, Rat	WB,IF-Cell,IHC-P,FC
	BLNK Mouse mAb [G1-F9]	EM1709-94	Human, Mouse	WB,IF-Cell,IHC-P,FC
	BLNK Recombinant Rabbit mAb [JE63-50]	HA721028	Human, Mouse, Rat	WB
	BOB1 Recombinant Rabbit mAb [JE47-78]	ET7109-68	Human	WB,IHC-P,FC
	BOB1 Rabbit pAb	ER1803-97	Human, Mouse, Rat	WB,FC
	BST2 Recombinant Rabbit mAb [JU11-93]	ET1706-46	Human	WB,IHC-P
	Calcineurin A Recombinant Rabbit mAb [JA64-11]	ET1704-02	Human, Mouse, Rat	WB,FC
	Calcineurin B Rabbit pAb	HA500346	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	CARD11 Recombinant Rabbit mAb [JE37-28]	HA722309	Human, Mouse, Rat	WB,IHC-P
Relative expression	CD10 Recombinant Mouse mAb [A1G3-R]	HA601246	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	CD10 Mouse mAb [A1G5]	EM1901-26	Human, Mouse	WB,IHC-P,IF-Tissue,FC
	CD10 Mouse mAb [A1G4]	EM1901-25	Human	WB,IHC-P,FC
	CD10 Mouse mAb	IRS222MS	Mouse	mIHC
	CD10 Mouse mAb	IRS101MS	Human	mIHC
Relative expression	CD10 Recombinant Mouse mAb [A1G3-R] - BSA and Azide free	HA610121	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Relative expression	CD10 Rabbit pAb	0805-4	Human, Mouse, Rat	WB,IHC-P
Relative expression	CD10 Recombinant Rabbit mAb [SN75-07] - BSA and Azide free	HA750273	Human, Mouse, Rat	WB,FC,IHC-P,IF-Cell,IF-Tissue
Relative expression	CD10 Recombinant Rabbit mAb [SN75-07]	ET1611-82	Human, Mouse, Rat	WB,FC,IHC-P,IF-Cell,IF-Tissue
Relative expression	CD10 Recombinant Rabbit mAb [JE60-45]	HA722864	Human	WB,IHC-P
Relative expression	CD10 Mouse mAb [A1G3]	EM1901-24	Human, Mouse	WB,IHC-P,FC,IF-Tissue
	CD100 Mouse mAb [5H6E3]	EM1710-47	Human	WB,IHC-P,FC
KD	CD100 Recombinant Rabbit mAb [JE63-21]	HA721013	Human	WB
Relative expression	CD11a Recombinant Rabbit mAb [SC56-01]	ET1610-1	Human	WB,IF-Tissue,IHC-P,IP
	CD11b Mouse mAb [A2-B5-D2]	EM1902-35	Human, Mouse	WB
Relative expression	CD11b Mouse mAb [A2-9]	M1211-6	Human	WB,IHC-P,FC,IF-Tissue
	CD11b Recombinant Rabbit mAb	IRS014	Human	mIHC
Relative expression	CD11b Recombinant Rabbit mAb [JE01-17]	HA722856	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Tissue
Relative expression	CD11b Recombinant Rabbit mAb [JU93-81]	ET1706-04	Human	WB,IP,IHC-P,IF-Cell,IF-Tissue,mIHC
Relative expression	CD11b Recombinant Rabbit mAb [JU93-81] - BSA and Azide free	HA750451	Human	WB,IP,IHC-P,IF-Cell,IF-Tissue
Relative expression	CD11b Recombinant Rabbit mAb [PSH03-96] - BSA and Azide free	HA750921	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Relative expression	CD11b Recombinant Rabbit mAb [PSH03-96]	HA722075	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	CD11b Mouse mAb [4-6]	EM1701-42	Human	IHC-P,IF-Cell,FC,WB,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CD11b/ITAM Mouse mAb [4G2]	RT1107	Human, Mouse	WB,IP,IF,FC
	CD11c Recombinant Rabbit mAb	IRS001	Human	mIHC
	CD11c Recombinant Rabbit mAb [SI19-06]	ET1606-19	Human	WB,IHC-P,IP,mIHC,IF-Tissue
	CD11c Recombinant Antibody [SI19-06] - Mouse IgG1 (Chimeric)	HA601545	Human	IF-Tissue,IHC-P
	CD11c Recombinant Antibody [SI19-06] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610361	Human	IF-Tissue,IHC-P
	CD11c Recombinant Rabbit mAb	IRS084	Mouse	mIHC
	CD11c Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS102MS	Human	mIHC
	CD11c Recombinant Antibody [SI19-06] - Rat IgG1 (Chimeric)	HA601418		
Relative expression	CD11c Recombinant Rabbit mAb [PSH07-42] - BSA and Azide free	HA751163	Mouse	WB,IHC-P,IHC-Fr,IF-Tissue
	CD11c Recombinant Antibody [SI19-06] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610258		
Relative expression	CD11c Recombinant Rabbit mAb [PSH07-42]	HA722830	Mouse	WB,IHC-P,IHC-Fr,IF-Tissue
	CD11c Recombinant Antibody - Rat IgG1 (Chimeric)	IRS103RT	Human	mIHC
	CD11c / ITGAX Mouse mAb [3G1]	RT1108	Human	WB,IP,IF,IHC-P
KD	CD127 Recombinant Rabbit mAb [JE37-10]	HA721214	Human, Mouse	WB,IF-Cell,IHC-P,FC
	APC Conjugated CD127 Recombinant Mouse mAb [PSH12-76]	HA601536F	Human	FC
	CD127 Recombinant Mouse mAb [PSH12-76]	HA601429	Human	FC
	CD127 Rabbit pAb	0804-9	Human	WB,IHC-P,IF-Cell,FC
Relative expression	CD134 Recombinant Rabbit mAb [PSH10-14]	HA723183	Human	WB,IF-Cell,FC,IP
	CD134 Mouse mAb [2H1]	EM1701-24	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Relative expression	CD134 Recombinant Rabbit mAb [JE36-48]	HA722700	Human	WB
Relative expression	CD134 Recombinant Rabbit mAb [PSH10-14] - BSA and Azide free	HA751343	Human	WB,IF-Cell,FC,IP
	CD14 Rabbit pAb	R1107-2	Human, Mouse	WB,IF-Cell,FC
KD	CD14 Recombinant Rabbit mAb [SC69-02]	ET1610-85	Human, Mouse	WB,IF-Cell,IHC-P,mIHC,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	CD14 Recombinant Rabbit mAb [PD01-63]	HA721321	Human, Mouse	WB,IHC-P,IF-Tissue
	CD14 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS106RT	Human	mIHC
	CD14 Recombinant Rabbit mAb	IRS015	Human	mIHC
KD	CD14 Recombinant Rabbit mAb [SC69-02] - BSA and Azide free	HA750232	Human, Mouse	WB,IF-Cell,IHC-P,FC
KD	CD14 Recombinant Antibody [SC69-02] - Rat IgG1 (Chimeric)	HA601420		
KD	CD14 Recombinant Antibody [SC69-02] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610260	0	
	CD16 Recombinant Rabbit mAb	IRS017	Human	mIHC
	CD16 Mouse mAb [Q32]	EM60067	Human	WB,IHC-P
Relative expression	CD16 Recombinant Rabbit mAb [PD00-12]	HA721149	Human	IHC-P,IF-Tissue,IF-Cell,mIHC,WB
	CD16 Recombinant Rabbit mAb [JE49-79]	ET7109-97	Human, Mouse, Rat	WB
	CD161 Mouse mAb [A8E1]	HA601063	Human	WB,IHC-P,FC
	CD161c Rabbit pAb	ER1905-79	Mouse, Rat	WB
	CD163 Rabbit pAb	ER1804-03	Human, Mouse, Rat	WB,IHC-P
Relative expression	CD163 Recombinant Rabbit mAb [PSH15-06]	HA723706	Human	IHC-P,WB
	CD163 Recombinant Rabbit mAb [JA51-30] - BSA and Azide free	HA750412	Human	WB,IHC-P,IF-Tissue,IP
	CD163 Mouse mAb [A3B5]	EM1901-91	Human	IHC-P,FC
	CD163 Recombinant Antibody [JA51-30] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610364	Human	IHC-P
	CD163 Rabbit pAb	IRS109RB	Mouse	mIHC
Relative expression	CD163 Recombinant Rabbit mAb [PSH18-98]	HA724034	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	CD163 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS276RT	Human	mIHC
	CD163 Mouse mAb [A3B3]	EM1901-90	Human	WB,IHC-P,FC
	CD163 Recombinant Antibody [JA51-30] - Rat IgG1 (Chimeric)	HA601549	Human	IHC-P
	CD163 Recombinant Rabbit mAb [JA51-30]	ET1704-43	Human	WB,IHC-P,mIHC,IF-Tissue,IP
	CD163 Recombinant Rabbit mAb	IRS002	Human	mIHC
	CD166 Recombinant Rabbit mAb [JE57-47]	HA722705	Human, Mouse, Rat	WB,IHC-P,FC,IHC-Fr

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	CD166 Mouse mAb [A6A8]	HA600083	Human	IF-Cell,IHC-P,WB
KD	CD166 Mouse mAb [5-5]	EM1902-36	Human	WB,IHC-P,IF-Cell,FC
KD	CD166 Mouse mAb [11-6]	M1012-7	Human, Mouse, Rat	WB,IHC-P,FC
KD	CD166 Mouse mAb [A6A9]	HA600082	Human	IF-Cell,IHC-P,WB
Cell treatment	CD169 Recombinant Rabbit mAb [PSH05-21] - BSA and Azide free	HA750992	Human, Mouse	WB,FC
Cell treatment	CD169 Recombinant Rabbit mAb [PSH05-21]	HA722246	Human, Mouse	WB,FC
	CD18 Recombinant Rabbit mAb [JM71-31]	ET1703-94	Human	WB,IP
	CD180 Rabbit pAb	ER2001-06	Human, Mouse, Rat	WB,IHC-P,FC
	CD19 Recombinant mAb	IRS110RB	Mouse	mIHC
Relative expression	CD19 Recombinant Rabbit mAb [PSH13-44]	HA723533	Human	WB,IF-Cell,FC
Relative expression	CD19 Recombinant Rabbit mAb [JF100-06] - BSA and Azide free	HA750365	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
Relative expression	CD19 Recombinant Rabbit mAb [PSH13-44] - BSA and Azide free	HA751463	Human	WB,IF-Cell,FC
	CD19 Rabbit pAb	R1107-7	Human, Mouse	WB
Relative expression	CD19 Recombinant Rabbit mAb [PSH03-94]	HA722073	Mouse	WB,IHC-P,IF-Tissue,IP,IHC-Fr,FC,IF-Cell,mIHC
	CD19 Recombinant Rabbit mAb [JF099-9]	ET1702-74	Human, Mouse, Rat	WB,FC
	CD19 Recombinant Rabbit mAb	IRS018	Human	mIHC
	FITC Conjugated CD19 Recombinant Rabbit mAb [PSH13-44]	HA723762F	Human	FC
Relative expression	CD19 Recombinant Rabbit mAb [PSH03-94] - BSA and Azide free	HA750919	Mouse	WB,IHC-P,IF-Tissue,IP,IHC-Fr,FC,IF-Cell
Relative expression	APC Conjugated CD19 Recombinant Rabbit mAb [PSH13-44]	HA723761F	Human	FC
	PE Conjugated CD19 Recombinant Rabbit mAb [PSH13-44]	HA723763F	Human	FC
Relative expression	CD19 Recombinant Rabbit mAb [JF100-06]	ET1702-93	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC,mIHC
	FITC Conjugated CD19 Recombinant Mouse mAb [PSH14-97]	HA601518F	Human	FC
	CD19 Mouse mAb [A3-A8]	EM40308	Human, Mouse	WB,IHC-P,FC
	CD1a Recombinant Rabbit mAb [JE28-74]	ET7110-48	Human	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	CD1a Recombinant Rabbit mAb [JM21-33]	ET1705-17	Human	WB,IHC-P,IF-Tissue
	CD1a Mouse mAb [A3F9]	EM1901-89	Human	WB,IHC-P
Relative expression	CD1a Mouse mAb [A4B10]	HA600031	Human	IHC-P,FC
	CD1a Rabbit pAb	ER1902-73	Human	WB,FC
	CD1c Recombinant mAb	IRS111RB	Human	mIHC
	CD1c Recombinant Rabbit mAb [JE55-00]	HA722562	Human	WB,IHC-P
	CD1C Mouse mAb [3G1B3]	EM1710-74	Human	WB,IHC-P,FC
Relative expression	CD2 Recombinant Rabbit mAb [PSH05-85]	HA722484	Human	WB,IHC-P,IF-Cell,FC,IF-Tissue
Relative expression	CD2 Recombinant Rabbit mAb [PSH05-85] - BSA and Azide free	HA751033	Human	WB,IHC-P,IF-Cell,FC,IF-Tissue
	CD2 Mouse mAb [B4-D1]	EM1707-56	Human, Mouse	WB,IHC-P,FC
Relative expression	CD2 Recombinant Rabbit mAb [PSH05-84]	HA722483	Human	WB,IHC-P
	CD20 Recombinant Rabbit mAb	IRS003	Human	mIHC
Relative expression	CD20 Recombinant Rabbit mAb [PSH08-12]	HA722953	Human	WB,IF-Cell,IHC-P,FC,IP
	APC Conjugated CD20 Recombinant Rabbit mAb [PSH08-12]	HA723838F	Human	FC
	PE Conjugated CD20 Recombinant Rabbit mAb [PSH08-12]	HA723839F	Human	FC
	CD20 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS113RT	Human	mIHC
	CD20 Rabbit pAb	IRS114RB	Mouse	mIHC
	CD20 Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS112MS	Human	mIHC
Relative expression	CD20 Recombinant Rabbit mAb [PSH08-12] - BSA and Azide free	HA751207	Human	WB,IF-Cell,IHC-P,FC,IP
Relative expression	CD20 Recombinant Rabbit mAb [SY12-01]	ET1605-33	Human, Monkey	WB,IHC-P,IF-Tissue,IF-Cell,FC,IP
Relative expression	CD20 Recombinant Antibody [PD00-02] - Mouse IgG1 (Chimeric)	HA601544	Human, Mouse	WB,IHC-P,IHC-Fr
Relative expression	CD20 Recombinant Rabbit mAb [PD00-02]	HA721138	Human, Mouse	IHC-P,mIHC,WB,IF-Cell,IP,IHC-Fr
Relative expression	CD20 Recombinant Antibody [PD00-02] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610360	Human, Mouse	WB,IHC-P,IHC-Fr
Relative expression	CD200 Recombinant Rabbit mAb [PSH01-48]	HA721691	Human	WB,IHC-P,IF-Tissue,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	CD200 Recombinant Rabbit mAb [PSH01-49]	HA721692	Human	WB,IHC-P,FC,IF-Tissue
Cell treatment	CD204 Mouse mAb [A8G2]	HA601090	Human	WB,IHC-P
Cell treatment	CD204 Recombinant Mouse mAb [A8G2-R] - BSA and Azide free	HA610052	Human, Rat	WB,IHC-P,IF-Cell
Cell treatment	CD204 Recombinant Mouse mAb [A8G2-R]	HA601165	Human, Rat	WB,IHC-P,IF-Cell
Relative expression	CD21 Recombinant Rabbit mAb [SC0681]	ET1610-61	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC,IF-Tissue
	CD21 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS265RT	Human	mIHC
	CD21 Recombinant Rabbit mAb [PD00-23]	HA721163	Human	WB,IHC-P,IF-Cell,FC,mIHC,IF-Tissue
	CD21 Recombinant Rabbit mAb	IRS019	Human	mIHC
Relative expression	CD21 Mouse mAb [A4E5]	HA600001	Human	WB,IF-Cell,IHC-P
	CD21 Rabbit pAb	ER1706-56	Human	IF-Cell,IHC-P,FC
	CD22 Mouse mAb [A2G6]	EM1902-13	Human	WB,IHC-P,FC
	CD22 Recombinant Rabbit mAb [PO00-31]	HA721157	Human	WB,IHC-P,IF-Tissue
	CD223 Mouse mAb [4H6A3]	EM1710-97	Human	WB,FC
	CD23 Rabbit pAb	HA500386	Human	IHC-P,FC
	CD23 Recombinant Rabbit mAb	IRS021	Human	mIHC
	CD23 Recombinant Rabbit mAb [PD00-03]	HA721139	Human	IHC-P,mIHC,IF-Tissue
	CD24 Rabbit pAb	0804-4	Human	WB,IF-Cell,FC
Cell treatment, Relative expression	CD24 Recombinant Rabbit mAb [PSH10-85] - BSA and Azide free	HA751377	Human	IF-Cell,FC,WB
	CD24 Rabbit pAb	R1303-1	Human	WB,IHC-P
Cell treatment, Relative expression	CD24 Recombinant Rabbit mAb [PSH10-85]	HA723252	Human	IF-Cell,FC,WB
	CD25 / IL-2 Receptor alpha Recombinant Mouse mAb [PSH14-46]	HA601486	Human, Mouse	FC,IHC-P,IF-Tissue
	APC Conjugated CD25 / IL-2 Receptor alpha Recombinant Mouse mAb [PSH14-46]	HA601530F	Human	FC
Relative expression	CD27 Recombinant Rabbit mAb [PSH08-90] - BSA and Azide free	HA751255	Human	WB,IF-Cell,IHC-P,FC,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	CD27 Recombinant Rabbit mAb [JB40-98]	ET7107-73	Human, Mouse, Rat	WB,IHC-P
	CD27 Recombinant Rabbit mAb	IRS115RB	Human	mIHC
	APC Conjugated CD27 Recombinant Rabbit mAb [PSH08-90]	HA723828F	Human	FC
Relative expression	CD27 Recombinant Rabbit mAb [PSH08-90]	HA723042	Human	WB,IF-Cell,IHC-P,FC,IP
Relative expression, KD	CD276 Recombinant Rabbit mAb [PD00-36]	HA721245	Human	IHC-P,FC,WB
KD	CD276 Rabbit pAb	ER1803-36	Human	WB,IF-Cell
Relative expression, KD	CD276 Recombinant Rabbit mAb [PSH02-46] - BSA and Azide free	HA750811	Human	WB,IF-Cell,FC
Relative expression, KD	CD276 Recombinant Rabbit mAb [PSH02-46]	HA721823	Human	WB,IF-Cell,FC
	CD276 Mouse mAb [C10-A1]	EM1709-85	Human	WB,IHC-P,FC
	CD28 Rabbit pAb	ER2001-42	Human, Mouse, Rat	WB,IHC-P,FC
	CD3 Recombinant Rabbit mAb	IRS074	Mouse	mIHC
	CD3 Recombinant Rabbit mAb	IRS022	Human	mIHC
	CD3 Recombinant Rabbit mAb [JE80-02]	HA720082	Human	WB,IHC-P,IF-Tissue,IP,mIHC,IF-Cell
	CD3 Recombinant Mouse mAb [PD01-43]	HA601110	Human	WB,IHC-P,IF-Tissue
Relative expression	CD3 Recombinant Rabbit mAb [PSH19-70]	HA724094	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Relative expression	CD30 Recombinant Rabbit mAb [PSH04-10]	HA722091	Human	WB,IHC-P,FC,IF-Tissue
Relative expression	CD30 Recombinant Rabbit mAb [JE57-71]	HA721579	Human	WB,IHC-P
Relative expression	CD30 Recombinant Rabbit mAb [PSH04-11]	HA722092	Human	WB,IF-Cell,IHC-P
Relative expression	CD30 Recombinant Rabbit mAb [PSH04-12]	HA722093	Human	WB,IF-Cell,IHC-P,IF-Tissue
Cell treatment	CD32 Recombinant Rabbit mAb [JM10-70] - BSA and Azide free	HA750382	Human	WB,IHC-P,IF-Cell,FC,IP
	CD32 Rabbit pAb	ER62945	Human, Mouse, Rat	WB,IHC-P
Cell treatment	CD32 Recombinant Rabbit mAb [JM10-70]	ET1703-43	Human	WB,IHC-P,IF-Cell,FC,IP
	CD32A + CD32B + CD32C Recombinant Rabbit mAb [JE64-17]	HA721930	Human	WB,IHC-P
	CD32B Mouse mAb [A9C11]	HA601121	Human	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CD32-C Rabbit pAb	ER62946	Human	WB
Relative expression	CD33 Recombinant Rabbit mAb [PSH02-51]	HA721827	Human	WB,IF-Cell,FC
	CD33 Mouse mAb [B7-C6]	EM1709-20	Human	IHC-P,FC
Relative expression, Cell treatment	CD33 Recombinant Rabbit mAb [PSH02-50] - BSA and Azide free	HA750813	Human	WB,IF-Cell,FC
Relative expression	CD33 Recombinant Rabbit mAb [PSH02-51] - BSA and Azide free	HA750814	Human	WB,IF-Cell,FC
Relative expression, Cell treatment	CD33 Recombinant Rabbit mAb [PSH02-50]	HA721826	Human	WB,IF-Cell,FC
	CD33 Rabbit pAb	ER1905-89	Mouse	WB,IHC-P,IF-Tissue
	CD34 Recombinant Rabbit mAb	IRS050	Mouse	mIHC
	CD34 Recombinant Rabbit mAb	IRS024	Human	mIHC
	CD34 Mouse mAb [15H1]	EM1901-01	Human, Mouse, Rat	IHC-P,FC
	CD34 Recombinant Rabbit mAb [SI16-01] - BSA and Azide free	HA750092	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	CD34 Rabbit pAb	0108-2	Human, Mouse, Zebrafish	WB,IF-Cell,IHC-P
	CD34 Rabbit pAb	ER0802	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC
	CD34 Recombinant Mouse mAb [PDM0-12]	HA601109	Human	IHC-P,IF-Tissue
	CD34 Recombinant Rabbit mAb [SI16-01]	ET1606-11	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,mIHC,FC
	CD34 Rabbit pAb	ER1706-66	Human, Mouse	WB,IHC-P
Relative expression	CD38 Recombinant Rabbit mAb [PD01-49]	HA721268	Human	WB,IHC-P,IF-Cell,FC,mIHC,IF-Tissue
	CD38 Recombinant Antibody [PD01-49] - Mouse IgG1 (Chimeric)	HA601554	Human	WB,IHC-P,FC
Relative expression	CD38 Recombinant Antibody [PD01-49] - Rat IgG1 (Chimeric)	HA601417		
	CD38 Recombinant Rabbit mAb [JE29-35]	ET7110-53	Human	WB,IHC-P,FC,IF-Tissue
Relative expression	CD38 Recombinant Antibody [PD01-49] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610257		
	CD38 Recombinant Antibody [PD01-49] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610369	Human	WB,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CD38 Recombinant Rabbit mAb	IRS006	Human	mIHC
	CD38 Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS120MS	Human	mIHC
	CD38 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS119RT	Human	mIHC
	CD38 Recombinant Rabbit mAb [SN07-20]	ET1611-66	Human	WB,IHC-P,FC,IF-Tissue
	CD39 Recombinant Rabbit mAb [JA90-36]	ET1704-74	Human, Mouse, Rat	WB,IHC-P
	CD39 Recombinant Rabbit mAb [JA90-36] - BSA and Azide free	HA750421	Human, Mouse, Rat	WB,IHC-P
	CD3D Recombinant Rabbit mAb [JJ08-97]	ET1701-57	Human	WB,IF-Cell,IHC-P,FC,IP,IF-Tissue
	CD3E Rabbit pAb	HA500344	Human, Mouse, Rat	WB
	CD3 epsilon Recombinant Rabbit mAb [SY0239] - BSA and Azide free	HA750115	Human	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP
	CD3 epsilon Recombinant Rabbit mAb [JE17-74]	ET7110-63	Human	WB,IHC-P,FC
	CD3 epsilon Recombinant Rabbit mAb [SY0239]	ET1607-29	Human	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP
	CD3 epsilon Rabbit pAb	ER80501	Human, Mouse, Rat	WB
	CD3G Recombinant Rabbit mAb [JB38-29]	ET7107-55	Human, Mouse	WB,IHC-P,FC,IP,IF-Cell,IF-Tissue
Relative expression	CD3 zeta Recombinant Rabbit mAb [JE09-37]	HA722692	Human, Rat	WB,IHC-P,IF-Tissue
	CD3 zeta Rabbit pAb	0805-7	Human	WB,IF-Cell,FC
Relative expression	CD4 Recombinant Rabbit mAb [PSH08-22] - BSA and Azide free	HA751213	Mouse	WB,IF-Cell,IHC-P,IHC-Fr
	CD4 Rabbit pAb	ER1803-84	Human	WB,FC
	CD4 Recombinant Rabbit mAb	IRS005	Human	mIHC
	CD4 Recombinant Rabbit Multiclonal Antibody [PSH13-24] - BSA and Azide free	HA751460	Human, Mouse	WB,IF-Cell,IHC-P
	CD4 Recombinant Rabbit mAb [JE56-36]	ET7110-64	Human	WB,IHC-P
	CD4 Recombinant Rabbit mAb	IRS051	Mouse	mIHC
KD	CD4 Recombinant Rabbit mAb [ST0488]	ET1609-52	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC,mIHC
	CD4 Recombinant Antibody [PSH08-22] - Rat IgG1 (Chimeric)	HA601534	Mouse	WB,IF-Tissue,IHC-P,IHC-Fr
Relative expression	CD4 Recombinant Rabbit mAb [PSH08-22]	HA722966	Mouse	WB,IF-Cell,IHC-P,mIHC,IHC-Fr

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CD4 Recombinant Rabbit Multiclonal Antibody [PSH13-24]	HA723511	Human, Mouse	WB,IF-Cell,IHC-P
	CD4 Rabbit pAb	R1211-1	Human	WB,IF-Cell
Relative expression	CD4 Recombinant Antibody [PSH08-22] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610353	Mouse	WB,IF-Tissue,IHC-P,IF-Cell
	CD4 Recombinant Antibody [ST0488] - Mouse IgG1 (Chimeric)	HA601557	Human	WB,IF-Cell,IHC-P
KD	CD4 Rabbit pAb	0108-6	Human, Mouse, Rat	WB,IF-Cell,FC
	CD4 Recombinant Antibody [PSH08-22] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610354	Mouse	WB,IF-Tissue,IHC-P,IHC-Fr
Relative expression	CD4 Recombinant Antibody [PSH08-22] - Mouse IgG1 (Chimeric)	HA601533	Mouse	WB,IF-Tissue,IHC-P,IF-Cell
	CD40 Rabbit pAb	ER1803-54	Human, Rat, Mouse	WB,IF-Cell,IHC-P,FC
	CD40 Mouse mAb [A1E5]	EM1901-41	Human	WB,IHC-P
Relative expression	CD40 Mouse mAb [A1E6]	EM1901-40	Human	WB,IF-Cell,IHC-P
	CD40 Recombinant Rabbit mAb [PSH10-01] - BSA and Azide free	HA751331	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-P,IHC-Fr,IF-Cell,FC,IF-Tissue
Relative expression	CD40 Recombinant Rabbit mAb [PSH05-90]	HA722488	Human, Cynomolgus monkey, Pig	WB,IF-Cell,FC
	CD40 Recombinant Rabbit mAb [PSH10-01]	HA723170	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-P,IHC-Fr,IF-Cell,FC,IF-Tissue
Relative expression	CD40 Recombinant Rabbit mAb [PSH05-90] - BSA and Azide free	HA751037	Human, Cynomolgus monkey, Pig	WB,IF-Cell,FC
	CD40L Recombinant Rabbit mAb [JM11-34] - BSA and Azide free	HA750387	Human, Mouse, Rat	WB
	CD40L Recombinant Rabbit mAb [JM11-34]	ET1703-64	Human, Mouse, Rat	WB
Relative expression	CD45 Recombinant Rabbit Multiclonal Antibody [PSH14-19] - BSA and Azide free	HA751510	Human, Mouse	WB,IHC-P,IF-Cell
	CD45 Recombinant Rabbit mAb [JE03-05]	ET7111-03	Human	WB,IHC-P,FC,IF-Tissue,mIHC
	CD45 Recombinant Rabbit mAb	IRS025	Human	mIHC
	CD45 Rabbit pAb	ER1902-02	Human	WB,FC
	CD45 Rabbit pAb	ER1901-29	Human	WB,IHC-P,FC
	CD45 Rabbit pAb	ER1905-93	Human	FC
Relative expression	CD45 Recombinant Rabbit Multiclonal Antibody [PSH14-19]	HA723606	Human, Mouse	WB,IHC-P,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	PE Conjugated CD45RO Recombinant Mouse mAb [PSH13-94]	HA601531F	Human	FC
	CD45RO Recombinant Mouse mAb [PSH13-94]	HA601441	Human	IHC-P,FC
	CD45RO Recombinant Mouse mAb	IRS124MS	Human	mIHC
	CD5 Recombinant Rabbit mAb [JE25-41]	ET7110-56	Human	WB,IHC-P
	CD5 Recombinant Rabbit mAb [PD00-01]	HA721137	Human	IHC-P,IF-Tissue,WB
	CD62P Recombinant Rabbit mAb [JM10-95]	ET1703-49	Human	WB,IHC-P
	CD64 Mouse mAb [5F5F8]	EM1710-99	Human	WB,IHC-P,FC
	CD64 Mouse mAb [A7E4]	HA601003	Human	WB
	CD66b Rabbit pAb	HA500100	Human	WB,IHC-P,mIHC,IF-Tissue
	CD66b pAb	IRS127RB	Human	mIHC
Relative expression	CD68 Mouse mAb [A3C3]	EM1901-94	Human	WB,IHC-P,FC
Relative expression	CD68 Recombinant Mouse mAb [A3C2-R] - BSA and Azide free	HA610164	Human	WB,IF-Cell,IHC-P
Relative expression	CD68 Recombinant Mouse mAb [A3C3-R]	HA601290	Human	WB,IF-Cell,IHC-P,FC
	CD68 Mouse mAb [F3-D7]	EM1706-11	Human	WB,IHC-P,IF-Cell,FC
	CD68 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS264RT	Human	mIHC
Relative expression	CD68 Mouse mAb [A3C2]	EM1901-93	Human	WB,IHC-P,FC
	CD68 Recombinant Rabbit mAb [SN07-27]	ET1611-53	Human	WB
	CD68 Recombinant Antibody [PDM0-13] - Rat IgG1 (Chimeric)	HA601520	Human	IF-Tissue,IHC-P
Relative expression, Cell treatment	CD68 Recombinant Rabbit mAb [PSH06-65]	HA722737	Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,IHC-Fr
Relative expression, Cell treatment	CD68 Recombinant Rabbit mAb [PSH06-65] - BSA and Azide free	HA751094	Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,IHC-Fr
	CD68 Recombinant Antibody [PDM0-13] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610347	Human	IF-Tissue,IHC-P
Cell treatment	CD68 Recombinant Rabbit mAb [PSH05-47]	HA722285	Mouse, Rat	WB,IF-Cell,IHC-P,IHC-Fr,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CD68 Rabbit pAb	ER1901-32	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	CD68 Recombinant Mouse mAb	IRS004	Human	mIHC
	CD68 Recombinant Antibody - Rabbit IgG (Chimeric)	IRS266RB	Human	mIHC
Relative expression	CD68 Recombinant Mouse mAb [PDM0-13]	HA601115	Human	IHC-P,IF-Cell,FC,mIHC,IF-Tissue
	CD68 Recombinant Antibody [PDM0-13] - Rabbit IgG (Chimeric)	HA723708	Human	IF-Tissue,IHC-P
	CD68 Recombinant Antibody [PDM0-13] - Rabbit IgG (Chimeric) - BSA and Azide free	HA751544	Human	IF-Tissue,IHC-P
Relative expression	CD68 Recombinant Mouse mAb [A3C4-R]	HA601292	Human	WB,IF-Cell,IHC-P
Relative expression	CD68 Recombinant Mouse mAb [A3C2-R]	HA601291	Human	WB,IF-Cell,IHC-P
Relative expression	CD68 Recombinant Mouse mAb [A3C4-R] - BSA and Azide free	HA610165	Human	WB,IF-Cell,IHC-P
Relative expression	CD68 Recombinant Mouse mAb [PDM0-13] - BSA and Azide free	HA610034	Human	IHC-P,IF-Cell,FC,IF-Tissue
Cell treatment	CD68 Recombinant Rabbit mAb [PSH05-47] - BSA and Azide free	HA751006	Mouse, Rat	WB,IF-Cell,IHC-P,IHC-Fr,IF-Tissue
Relative expression	CD68 Mouse mAb [A3C4]	EM1901-95	Human	WB,IHC-P,FC,mIHC
	CD69 Mouse mAb [E6-B6]	EM1709-21	Human	WB,IHC-P,FC
	CD69 Rabbit pAb	HA500073	Human, Mouse, Rat	WB
Cell treatment	CD69 Mouse mAb [A8C4]	HA601013	Human	WB
Relative expression	CD7 Recombinant Rabbit mAb [JE57-72]	HA720023	Human	WB,IF-Cell,IHC-P
	CD74 Recombinant Rabbit mAb [JF40-10]	ET1702-51	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
Relative expression	CD74 Recombinant Mouse mAb [PDM0-22]	HA601117	Human	WB,IHC-P,IF-Cell,IF-Tissue
	CD79a Rabbit pAb	R1204-4	Human, Mouse	WB,IF-Cell
	CD79a Recombinant Rabbit mAb [JE32-14]	HA722896	Human	WB,IF-Cell,IHC-P,IF-Tissue
	CD79a Recombinant Rabbit mAb [JE43-50]	ET7109-33	Human	WB,IF-Tissue,IHC-P,FC
	CD79a Recombinant Mouse mAb [PDM0-23]	HA601118	Human	WB,IHC-P
	CD79a Mouse mAb [C1-E8]	EM1902-29	Human, Rat, Mouse	WB,FC
	CD79b Recombinant Rabbit mAb [JE58-34]	HA720043	Human	WB,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CD8 Rabbit pAb	0108-7	Mouse	WB,IF-Cell,FC
	CD8 Rabbit pAb	IRS083	Mouse	mIHC
Relative expression	CD80 Recombinant Rabbit mAb [PSH09-79] - BSA and Azide free	HA751321	Human	WB,IF-Cell,FC
	CD80 Rabbit pAb	ER1905-99	Human, Mouse, Rat	WB
Relative expression, Cell treatment	CD80 Recombinant Rabbit mAb [PSH18-13]	HA723976	Mouse	WB,IF-Cell,FC,IP
Relative expression	CD80 Recombinant Rabbit mAb [JF100-4]	ET1702-95	Human	WB,IHC-P
Relative expression	CD80 Recombinant Rabbit mAb [PSH09-79]	HA723144	Human	WB,IF-Cell,FC
	CD80 Recombinant Mouse mAb [8-E5-R] - BSA and Azide free	HA610122	Human	WB
Relative expression	CD80 Mouse mAb [8-E5]	M1007-10	Human, Mouse	WB,IHC-P,FC
	CD80 Recombinant Mouse mAb [8-E5-R]	HA601247	Human	WB
	CD80 Rabbit pAb	1007-8	Human	WB,IF-Cell
Relative expression	CD81 Recombinant Rabbit mAb [PSH10-26]	HA723193	Human, Rat	WB,IHC-P,FC
Cell treatment	CD82 Recombinant Rabbit mAb [JE32-17]	HA722533	Human	WB
	CD83 Rabbit pAb	ER62949	Human	WB,IHC-P
	CD83 Recombinant Rabbit mAb [JE33-61]	HA721556	Human	WB
	CD84 Recombinant Rabbit mAb [JE59-27]	HA720037	Human	WB,IF-Cell,IHC-P,FC
	CD86 Rabbit pAb	ER1906-01	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
	CD86 Recombinant Rabbit mAb [SJ20-00] - BSA and Azide free	HA750109	Human, Mouse, Rat	WB,IF-Cell,IP
	CD86 Recombinant Rabbit mAb [SJ20-00]	ET1606-50	Human, Mouse, Rat	WB,IF-Cell,IP
	CD8a Mouse mAb	IRS007	Human	mIHC
	CD8 alpha Recombinant Rabbit mAb [SI18-01] - BSA and Azide free	HA750101	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
Relative expression	CD8 alpha Recombinant Rabbit mAb [PSH19-81]	HA724105	Human, Rat	WB,IHC-P,IF-Tissue,FC
	CD8 alpha Recombinant Rabbit mAb [SI18-01]	ET1606-31	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
	CD90 / THY1 Recombinant Rabbit mAb [SU35-07]	ET1608-46	Human, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CD90 / THY1 Recombinant Rabbit mAb [SU35-07] - BSA and Azide free	HA750149	Human, Rat	WB,IHC-P
KD	CD90 / THY1 Recombinant Rabbit mAb [JF10-09] - BSA and Azide free	HA750364	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
KD	CD90 / THY1 Recombinant Rabbit mAb [JF10-09]	ET1702-92	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
	CD90 / THY1 Rabbit pAb	R1107-9	Human, Mouse, Rat	WB
	CD99 Recombinant Rabbit mAb [PD00-13]	HA721150	Human	IHC-P,IF-Tissue
	CD99 Rabbit pAb	ER1803-81	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	CD99 Recombinant Rabbit mAb [JF0991]	ET1702-35	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Relative expression	CFTR Recombinant Rabbit mAb [PSH09-99]	HA723168	Human	WB,IP
Relative expression	CFTR Recombinant Rabbit mAb [PSH09-99] - BSA and Azide free	HA751330	Human	WB,IP
	ChT1 Rabbit pAb	HA500314	Human, Mouse	WB,IF-Cell,IHC-P,FC
	ChT1 Rabbit pAb	HA500237	Human, Mouse, Rat	WB,IHC-P
	c-kit Rabbit pAb	0108-4	Human	WB,IF-Cell
Relative expression, KD	c-Kit Recombinant Rabbit mAb [ST04-99]	ET1609-60	Human	WB,IHC-P
	c-Kit Recombinant Rabbit mAb	IRS013	Human	mIHC
Relative expression	c-Kit Recombinant Rabbit mAb [PSH17-96]	HA723961	Mouse, Rat	WB,IHC-Fr,IHC-P
Relative expression, KD	c-Kit Recombinant Rabbit mAb [PD00-24]	HA721154	Human	IHC-P,WB,mIHC
Relative expression	c-Kit Recombinant Rabbit mAb [PSH17-79]	HA723948	Mouse	IHC-Fr,IF-Cell,FC
	COPS8 / CSN8 Rabbit pAb	ER2001-29	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Relative expression	CSF-1-R Recombinant Rabbit mAb [PSH09-77] - BSA and Azide free	HA751319	Human, Mouse, Rat	WB,IHC-P
Relative expression	CSF-1-R Recombinant Rabbit mAb [PSH09-77]	HA723142	Human, Mouse, Rat	WB,IHC-P
	CSF-1-R Recombinant Rabbit mAb [PD00-45]	HA721180	Human	IHC-P,FC
	CSF1R Mouse mAb [B9-B9]	EM1708-56	Human	WB,IF-Cell,IHC-P
	CTLA-4 Recombinant Rabbit mAb [PS01-34] - BSA and Azide free	HA750581	Human, Mouse	IHC-P
	CTLA-4 Recombinant Rabbit mAb [PS01-34]	HA721269	Human, Mouse	IHC-P,mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CTLA-4 Recombinant Rabbit mAb	IRS016	Human	mIHC
	CTLA4 Rabbit pAb	ER1906-81	Human, Mouse, Rat	WB,FC
	CX3CR1 Mouse mAb [A5C7]	HA600074	Human	IF-Cell
	CXCR4 Rabbit pAb	ER1802-28	Human, Rat	WB,IF-Cell,FC
	CXCR4 Recombinant Rabbit mAb [JE03-35]	HA722304	Human	WB,IF-Cell,IHC-P
	CXCR5 Recombinant Rabbit mAb [JB11-40]	ET7107-74	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell
	CXCR5 Recombinant Rabbit mAb [JB11-40] - BSA and Azide free	HA750463	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell
	DC-SIGN Recombinant Rabbit mAb [JE31-24]	HA721382	Human	WB
	DC-SIGN Rabbit pAb	ER1907-23	Human	WB,IHC-P,IF-Tissue,FC
KO	DPP4 / CD26 Recombinant Rabbit mAb [JM11-42] - BSA and Azide free	HA750396	Human, Mouse, Rat	WB
Cell treatment	DPP4 / CD26 Recombinant Rabbit mAb [PSH12-59] - BSA and Azide free	HA751446	Human	WB,IHC-P,IF-Cell,FC,IP
Cell treatment	DPP4 / CD26 Recombinant Rabbit mAb [PSH12-59]	HA723473	Human	WB,IHC-P,IF-Cell,FC,IP
KO	DPP4 / CD26 Recombinant Rabbit mAb [JM11-42]	ET1703-93	Human, Mouse, Rat	WB
Relative expression	F4/80 Recombinant Rabbit mAb [PSH01-87]	HA721745	Mouse	WB,IHC-P,IF-Tissue,IHC-Fr
	F4/80 Recombinant Antibody [PSH01-87] - Rat IgG1 (Chimeric)	HA601538	Mouse	WB,IHC-Fr,IHC-P
Relative expression	F4/80 Recombinant Antibody [PSH01-87] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610356	Mouse	WB,IHC-P
Relative expression	F4/80 Recombinant Antibody [PSH01-87] - Guinea pig IgG2 (Chimeric)	HA601567	Mouse	WB,IHC-Fr
	F4/80 Recombinant Rabbit mAb	IRS052	Mouse	mIHC
Relative expression	F4/80 Recombinant Rabbit mAb [PSH01-87] - BSA and Azide free	HA750761	Mouse	WB,IHC-P,IF-Tissue,IHC-Fr
Relative expression	F4/80 Recombinant Antibody [PSH01-87] - Mouse IgG1 (Chimeric)	HA601537	Mouse	WB,IHC-P
	F4/80 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS143RT	Mouse	mIHC
	F4/80 Recombinant Antibody [PSH0-85] - Rat IgG1 (Chimeric)	HA601485	Mouse, Rat	IHC-P
	F4/80 Recombinant Antibody [PSH0-85] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610319	Mouse, Rat	IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	F4/80 Recombinant Rabbit mAb [PSH0-85]	HA721520	Mouse, Rat	WB,IHC-P,IF-Tissue,IHC-Fr,mIHC
Relative expression	F4/80 Recombinant Rabbit mAb [PSH0-85] - BSA and Azide free	HA750665	Mouse, Rat	WB,IHC-P,IF-Tissue,IHC-Fr
	F4/80 Recombinant Antibody [PSH01-87] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610357	Mouse	WB,IHC-Fr,IHC-P
	FCGR1A/CD64 Rabbit pAb	ER1908-78	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue
	FIt3 Rabbit pAb	ER2001-25	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell
	FLT3 Mouse mAb [C3-B7]	EM1709-22	Human	WB,IHC-P
	FOXP3 Recombinant Rabbit mAb [JF0898] - BSA and Azide free	HA750335	Human, Mouse, Rat	WB,IF-Cell,FC,IHC-P
Relative expression	FOXP3 Recombinant Rabbit mAb [PSH07-47]	HA722835	Human, Mouse, Rat	WB,IHC-P,mIHC
Relative expression	FOXP3 Recombinant Rabbit mAb [PSH07-47] - BSA and Azide free	HA751168	Human, Mouse, Rat	WB,IHC-P
	FOXP3 Recombinant Rabbit mAb [JF0926]	ET1702-16	Human	WB,FC,IHC-P
	FOXP3 Recombinant Rabbit mAb [JF0898]	ET1702-12	Human, Mouse, Rat	WB,IF-Cell,FC,IHC-P
	FOXP3 Recombinant Rabbit mAb	IRS053	Mouse	mIHC
	FOXP3 Recombinant Rabbit mAb	IRS008	Human	mIHC
	FOXP3 Rabbit pAb	IRS147RB	Mouse	mIHC
Relative expression	GATA2 Recombinant Rabbit mAb [PSH04-53] - BSA and Azide free	HA750960	Human, Rat	WB,IF-Cell
Relative expression	GATA2 Recombinant Rabbit mAb [PSH04-53]	HA722152	Human, Rat	WB,IF-Cell
	GATA-2/3 Rabbit pAb	ER63507	Human, Mouse	WB,IHC-P,IF-Cell
Relative expression	GATA2 + GATA3 Recombinant Rabbit mAb [JE04-46]	HA722713	Human, Mouse, Rat	WB,IHC-P
	GATA3 Recombinant Mouse mAb [A3G3-R]	HA601189	Human	WB,IF-Cell,IHC-P,FC
	GATA3 Recombinant Mouse mAb [A3G3-R] - BSA and Azide free	HA610069	Human	WB,IF-Cell,IHC-P,FC
	GATA3 Mouse mAb [A3G6]	EM1902-23	Human	WB
	GATA3 Mouse mAb [A3G3]	EM1902-33	Human	WB,IF-Cell,IHC-P
	GATA3 Rabbit pAb	ER1902-69	Human, Rat	WB,IHC-P,FC
	GATA3 Rabbit pAb	ER1901-20	Human, Rat	WB,IF-Cell,IHC-P,FC
	GATA3 Recombinant Mouse mAb [PD00-41]	HA601132	Human	WB,IHC-P,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	GITR Recombinant Rabbit mAb [PSH04-06] - BSA and Azide free	HA750929	Human	WB,IF-Cell,FC
	GITR Recombinant Rabbit mAb [PSH04-06]	HA722087	Human	WB,IF-Cell,FC
Relative expression	Granzyme B Recombinant Rabbit mAb [PSH14-66]	HA723664	Human, Mouse	WB,IHC-P
Relative expression	Granzyme B Recombinant Rabbit mAb [PSH14-66] - BSA and Azide free	HA751527	Human, Mouse	WB,IHC-P
	Granzyme B Rabbit pAb	HA500252	Human	IHC-P,mIHC
	Granzyme B Mouse mAb [4G9]	RT1259	Human, Mouse, Rat	WB,IP,IF,IHC-P
	Granzyme B Rabbit pAb	IRS032	Human	mIHC
	Histamine H2 receptor Rabbit pAb	HA500154	Human	WB,IHC-P
KD	HNRNPA0 Recombinant Rabbit mAb [PSH02-77] - BSA and Azide free	HA750832	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC
KD	HNRNPA0 Recombinant Rabbit mAb [PSH02-77]	HA721855	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC
	FITC Conjugated Human CD10 Recombinant Rabbit mAb [SN75-07]	HA720209F	Human	FC
	FITC Conjugated Human CD11b Recombinant Mouse mAb [A2-9-R]	HA600130F	Human	FC
	Biotin Conjugated Human CD127 Recombinant Rabbit mAb [PSH18-28] - Detector	HA725346B	Human	ELISA(Det),ELISA
	Human CD127 Recombinant Rabbit mAb [PSH18-28] - BSA and Azide free (Detector)	HA725346	Human	ELISA(Det)
	Human CD127 Recombinant Rabbit mAb [PSH18-27] - BSA and Azide free (Capture)	HA725345	Human	ELISA(Cap)
	iFluor™ 647 Conjugated Human CD14 Recombinant Rabbit mAb [SC69-02]	HA720184F	Human	FC
	FITC Conjugated Human CD14 Recombinant Rabbit mAb [SC69-02]	HA720194F	Human	FC
	PE Conjugated Human CD16 Recombinant Mouse mAb [PSH04-95]	HA600129F	Human	FC
	Human CD16 Recombinant Mouse mAb [PSH04-95]	HA601340	Human	FC
	Human CD163 Recombinant Rabbit mAb [PSH13-79] - BSA and Azide free (Detector)	HA723567	Human	ELISA(Det)
	Human CD163 Recombinant Rabbit mAb [PSH13-80] - BSA and Azide free (Capture)	HA723569	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human CD163 Recombinant Rabbit mAb [PSH13-78] - BSA and Azide free (Capture)	HA723566	Human	ELISA(Cap)
	Biotin Conjugated Human CD163 Recombinant Rabbit mAb [PSH13-79] - Detector	HA723568B	Human	ELISA(Det),ELISA
	Human CD200 Recombinant Rabbit mAb [PSH05-59] - BSA and Azide free (Capture)	HA722379	Human	ELISA(Cap)
	Biotin Conjugated Human CD200 Recombinant Rabbit mAb [PSH05-60]	HA722591B	Human	ELISA(Det)
	Human CD200 Recombinant Rabbit mAb [PSH05-60] - BSA and Azide free (Detector)	HA722380	Human	ELISA(Det)
	FITC Conjugated Human CD200 Recombinant Rabbit mAb [PSH01-48]	HA720208F	Human	FC
	FITC Conjugated Human CD200R1 Recombinant Rabbit mAb [PSH05-02]	HA720212F	Human	FC
	Human CD200R1 Recombinant Rabbit mAb [PSH05-28] - BSA and Azide free (Detector)	HA722257	Human	ELISA(Det)
	Biotin Conjugated Human CD200R1 Recombinant Rabbit mAb [PSH05-28]	HA722598B	Human	ELISA(Det)
Cell treatment	Human CD200R1 Recombinant Rabbit mAb [PSH05-02]	HA722211	Human	IF-Cell,FC
	Human CD200R1 Recombinant Rabbit mAb [PSH05-27] - BSA and Azide free (Capture)	HA722256	Human	ELISA(Cap)
	Human CD276 Recombinant Rabbit mAb [PSH05-93] - BSA and Azide free (Detector)	HA722492	Human	ELISA(Det)
	Biotin Conjugated Human CD276 Recombinant Rabbit mAb [PSH05-93] - Detector	HA722588B	Human	ELISA(Det),ELISA
	Human CD276 Recombinant Rabbit mAb [PSH05-92] - BSA and Azide free (Capture)	HA722491	Human	ELISA(Cap)
	Human CD27/TNFRSF7 Recombinant Rabbit mAb [PSH08-31] - BSA and Azide free (Capture)	HA722973	Human	ELISA(Cap)
	Biotin Conjugated Human CD27/TNFRSF7 Recombinant Rabbit mAb [PSH08-33] - Detector	HA722976B	Human	ELISA(Det),ELISA
	Human CD27/TNFRSF7 Recombinant Rabbit mAb [PSH08-33] - BSA and Azide free (Detector)	HA722975	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human CD27/TNFRSF7 Recombinant Rabbit mAb [PSH08-32] - BSA and Azide free (Capture)	HA722974	Human	ELISA(Cap)
	Human CD28 Recombinant Rabbit mAb [PSH12-51] - BSA and Azide free (Detector)	HA723458	Human	ELISA(Det)
	Biotin Conjugated Human CD28 Recombinant Rabbit mAb [PSH12-51] - Detector	HA723459B	Human	ELISA(Det),ELISA
	Human CD28 Recombinant Rabbit mAb [PSH12-50] - BSA and Azide free (Capture)	HA723457	Human	ELISA(Cap)
	FITC Conjugated Human CD3 Recombinant Mouse mAb [PSH04-92]	HA600116F	Human	FC
	APC Conjugated Human CD3 Recombinant Mouse mAb [PSH04-92]	HA600118F	Human	FC
	Percp-Cy5.5 Conjugated Human CD3 Recombinant Mouse mAb [PSH04-92]	HA601548F	Human	FC
	PE Conjugated Human CD3 Recombinant Mouse mAb [PSH04-92]	HA600117F	Human	FC
	Human CD3 Recombinant Mouse mAb [PSH04-92]	HA601337	Human	FC
	Biotin Conjugated Human CD30 Recombinant Rabbit mAb [PSH04-82]	HA722601B	Human	ELISA(Det)
	Human CD30 Recombinant Rabbit mAb [PSH04-82] - BSA and Azide free (Detector)	HA722188	Human	ELISA(Det)
	Human CD30 Recombinant Rabbit mAb [PSH04-81] - BSA and Azide free (Capture)	HA722187	Human	ELISA(Cap)
	Human CD33 Recombinant Rabbit mAb [PSH05-26] - BSA and Azide free (Detector)	HA722255	Human	ELISA(Det)
	FITC Conjugated Human CD33 Recombinant Rabbit mAb [PSH04-96]	HA720197F	Human	FC
	Biotin Conjugated Human CD33 Recombinant Rabbit mAb [PSH05-26]	HA722599B	Human	ELISA(Det)
	Human CD33 Recombinant Rabbit mAb [PSH05-25] - BSA and Azide free (Capture)	HA722254	Human	ELISA(Cap)
	FITC Conjugated Human CD38 Recombinant Rabbit mAb [SN07-20]	HA720196F	Human	FC
	FITC Conjugated Human CD4 Recombinant Mouse mAb [PSH04-93]	HA600119F	Human	FC
	PE Conjugated Human CD4 Recombinant Mouse mAb [PSH04-93]	HA600120F	Human	FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	APC Conjugated Human CD4 Recombinant Mouse mAb [PSH04-93]	HA600121F	Human	FC
	Human CD4 Recombinant Mouse mAb [PSH04-93]	HA601338	Human	FC
	PE-Cy7 Conjugated Human CD4 Recombinant Mouse mAb [PSH04-93]	HA601542F	Human	FC
	FITC Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600112F	Human	FC
	PE Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600126F	Human	FC
	APC Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600127F	Human	FC
	FITC Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600125F	Human	FC
	Human CD45 Recombinant Mouse mAb [PSH03-64]	HA601285	Human	IF-Cell,FC
	iFluor™ 488 Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600109F	Human	FC
	Human CD45 Recombinant Mouse mAb [PSH03-64] - BSA and Azide free	HA610158	Human	IF-Cell,FC
	iFluor™ 555 Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA601556F	Human	FC
	Percp-Cy5.5 Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600128F	Human	FC
	Human CD45RA Recombinant Mouse mAb [PSH04-97]	HA601341	Human	FC,IF-Cell,IHC-P
	FITC Conjugated Human CD45RA Recombinant Mouse mAb [PSH04-97]	HA600131F	Human	FC
	PE Conjugated Human CD45RA Recombinant Mouse mAb [PSH04-97]	HA601525F	Human	FC
	FITC Conjugated Human CD5 Recombinant Rabbit mAb [PSH05-01]	HA720211F	Human	FC
Relative expression	Human CD5 Recombinant Rabbit mAb [PSH05-01]	HA722210	Human	IF-Cell,FC
	FITC Conjugated Human CD68 Recombinant Mouse mAb [A3C3-R]	HA600132F	Human	FC
	Biotin Conjugated Human CD80 Recombinant Rabbit mAb [PSH09-62] - Detector	HA723126B	Human	ELISA(Det),ELISA
	Human CD80 Recombinant Rabbit mAb [PSH09-63] - BSA and Azide free (Detector)	HA723125	Human	ELISA(Det)
	Human CD80 Recombinant Rabbit mAb [PSH09-61] - BSA and Azide free (Capture)	HA723123	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human CD80 Recombinant Rabbit mAb [PSH09-63] - Detector	HA723127B	Human	ELISA(Det),ELISA
	Human CD80 Recombinant Rabbit mAb [PSH09-62] - BSA and Azide free (Detector)	HA723124	Human	ELISA(Det)
	Human CD83 Recombinant Rabbit mAb [PSH20-08] - BSA and Azide free (Capture)	HA725305	Human	ELISA(Cap)
	Biotin Conjugated Human CD83 Recombinant Rabbit mAb [PSH20-09] - Detector	HA725306B	Human	ELISA(Det),ELISA
	Human CD83 Recombinant Rabbit mAb [PSH20-09] - BSA and Azide free (Detector)	HA725306	Human	ELISA(Det)
	APC-Cy7 Conjugated Human CD8 alpha Recombinant Mouse mAb [PSH04-94]	HA601541F	Human	FC
	FITC Conjugated Human CD8 alpha Recombinant Mouse mAb [PSH04-94]	HA600122F	Human	FC
	PE Conjugated Human CD8 alpha Recombinant Mouse mAb [PSH04-94]	HA600123F	Human	FC
	APC Conjugated Human CD8 alpha Recombinant Mouse mAb [PSH04-94]	HA600124F	Human	FC
	Human CD8 alpha Recombinant Mouse mAb [PSH04-94]	HA601339	Human	FC
	Biotin Conjugated Human CSF-1-R Recombinant Rabbit mAb [PSH12-67] - Detector	HA723476B	Human	ELISA(Det),ELISA
	Human CSF-1-R Recombinant Rabbit mAb [PSH12-66] - BSA and Azide free (Capture/Detector)	HA723474	Human	ELISA(Cap),ELISA(Det)
	Human CSF-1-R Recombinant Rabbit mAb [PSH12-67] - BSA and Azide free (Detector)	HA723475	Human	ELISA(Det)
	HRP Conjugated Human CSF-1-R Recombinant Rabbit mAb [PSH12-66] - Detector	HA723478H	Human	ELISA(Det),ELISA
	Human CSF-1-R Recombinant Rabbit mAb [PSH12-68] - BSA and Azide free (Capture)	HA723477	Human	ELISA(Cap)
	Biotin Conjugated Human CSF-1-R Recombinant Rabbit mAb [PSH12-66] - Detector	HA723478B	Human	ELISA(Det),ELISA
	Human CTLA-4 Recombinant Rabbit mAb [PSH08-13] - BSA and Azide free (Capture)	HA722954	Human	ELISA(Cap)
	Human CTLA-4 Recombinant Rabbit mAb [PSH08-14] - BSA and Azide free (Detector)	HA722955	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human CTLA-4 Recombinant Rabbit mAb [PSH08-14]	HA722956B	Human	ELISA(Det)
	Biotin Conjugated Human GITR Recombinant Rabbit mAb [PSH05-64]	HA722593B	Human	ELISA(Det)
	Human GITR Recombinant Rabbit mAb [PSH05-63] - BSA and Azide free (Capture)	HA722383	Human	ELISA(Cap)
	Human GITR Recombinant Rabbit mAb [PSH05-64] - BSA and Azide free (Detector)	HA722384	Human	ELISA(Det)
	Human Granzyme B Recombinant Rabbit mAb [PSH15-99] - BSA and Azide free (Detector)	HA723804	Human	ELISA(Det)
	Human Granzyme B Recombinant Rabbit mAb [PSH15-98] - BSA and Azide free (Capture)	HA723803	Human	ELISA(Cap)
	Biotin Conjugated Human Granzyme B Recombinant Rabbit mAb [PSH15-99] - Detector	HA723804B	Human	ELISA(Det),ELISA
	HRP Conjugated Human ICAM1 Recombinant Rabbit mAb [PSH10-16] - Detector	HA723474H	Human	ELISA(Det),ELISA
	Human ICAM1 Recombinant Rabbit mAb [PSH10-15] - BSA and Azide free (Capture)	HA723184	Human	ELISA(Cap)
	Biotin Conjugated Human ICAM1 Recombinant Rabbit mAb [PSH10-16] - Detector	HA723186B	Human	ELISA(Det),ELISA
	Human ICAM1 Recombinant Rabbit mAb [PSH10-16] - BSA and Azide free (Detector)	HA723185	Human	ELISA(Det)
	Human ICOSLG/B7-H2/CD275 Recombinant Rabbit mAb [PSH11-33] - BSA and Azide free (Detector)	HA723318	Human	ELISA(Det)
	Biotin Conjugated Human ICOSLG/B7-H2/CD275 Recombinant Rabbit mAb [PSH11-33] - Detector	HA723319B	Human	ELISA(Det),ELISA
	Human ICOSLG/B7-H2/CD275 Recombinant Rabbit mAb [PSH11-34] - BSA and Azide free (Capture)	HA723320	Human	ELISA(Cap)
	Human ICOSLG/B7-H2/CD275 Recombinant Rabbit mAb [PSH11-35] - BSA and Azide free (Detector)	HA723321	Human	ELISA(Det)
	Biotin Conjugated Human ICOSLG/B7-H2/CD275 Recombinant Rabbit mAb [PSH11-35] - Detector	HA723322B	Human	ELISA(Det),ELISA
	Human ICOSLG/B7-H2/CD275 Recombinant Rabbit mAb [PSH11-32] - BSA and Azide free (Capture)	HA723317	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human IDO1 Recombinant Rabbit mAb [PSH09-86] - BSA and Azide free (Capture)	HA723154	Human	ELISA(Cap)
	Human IDO1 Recombinant Rabbit mAb [PSH09-87] - BSA and Azide free (Detector)	HA723155	Human	ELISA(Det)
	Biotin Conjugated Human IDO1 Recombinant Rabbit mAb [PSH09-87] - Detector	HA723156B	Human	ELISA(Det),ELISA
	Human IgG4 Mouse mAb [2G1]	RT1315	Human	ELISA
	HRP Conjugated Mouse Anti-Human IgG heavy chain [54-76-R] monoclonal Antibody	HA601284	Human	WB
	Human IgG heavy chain Mouse mAb [54-76]	M0806-1	Human	WB,ELISA
	Human IgG heavy chain Recombinant Mouse mAb [54-76-R]	HA601283	Human	WB,ELISA,IHC-P
	human IgM Mouse mAb [B6-8]	M1510-6	Human	WB,ELISA
	human IgM Mouse mAb [9-C1]	M1510-9	Human	WB,IHC-P,ELISA
	human IgM Mouse mAb [2-B5]	M1510-7	Human	WB,ELISA
	human IgM Mouse mAb [B7-3]	M1510-5	Human	WB,IHC-P,ELISA
	Human IgM Mouse mAb [4-B2]	M1510-8	Human	WB,IHC-P,ELISA
Cell treatment	Human IgM Recombinant Rabbit mAb [JF84-09]	ET1702-81	Human	WB,IHC-P,IP,FC,IF-Cell
	Human IgM Rabbit pAb	R1511-23	Human	WB,IHC-P
	Human IL-10 Recombinant Rabbit mAb [PS00-98] - BSA and Azide free (Capture)	HA721236	Human	ELISA(Cap)
	Biotin Conjugated Human IL-10 Recombinant Rabbit mAb [PS00-81]	HA722652B	Human	ELISA(Det)
	Human IL-10 Recombinant Rabbit mAb [PS00-81] - BSA and Azide free (Detector)	HA721237	Human	ELISA(Det)
	Human IL-12B p40 Recombinant Rabbit mAb [PSH02-58] - BSA and Azide free (Detector)	HA721835	Human	ELISA(Det)
	Biotin Conjugated Human IL-12B p40 Recombinant Rabbit mAb [PSH03-04]	HA722629B	Human	ELISA(Det)
	Human IL-12B p40 Recombinant Rabbit mAb [PSH02-57] - BSA and Azide free (Capture)	HA721834	Human	ELISA(Cap)
	Biotin Conjugated Human IL-12B p40 Recombinant Rabbit mAb [PSH02-58]	HA722633B	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human IL-12B p40 Recombinant Rabbit mAb [PSH03-04] - BSA and Azide free (Detector)	HA721934	Human	ELISA(Det)
	Human IL-18 Recombinant Rabbit mAb [PSH02-08] - BSA and Azide free	HA750778	Human	WB,IHC-P,IF-Tissue
Relative expression	Human IL-18 Recombinant Rabbit mAb [PSH02-09]	HA721768	Human	WB
	Human IL-18 Recombinant Rabbit mAb [PSH01-71] - BSA and Azide free (Capture)	HA721722	Human	ELISA(Cap)
	Biotin Conjugated Human IL-18 Recombinant Rabbit mAb [PSH01-70]	HA722644B	Human	ELISA(Det)
	Human IL-18 Recombinant Rabbit mAb [PSH01-70] - BSA and Azide free (Detector)	HA721723	Human	ELISA(Det)
	Human IL-18 Recombinant Rabbit mAb [PSH02-08]	HA721767	Human	WB,IHC-P,IF-Tissue,IP
Relative expression	Human IL-18 Recombinant Rabbit mAb [PSH02-09] - BSA and Azide free	HA750779	Human	WB
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-87] - BSA and Azide free (Capture)	HA723883	Human	ELISA(Cap)
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS01-39]	HA721712	Human	WB
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-86] - BSA and Azide free (Detector)	HA723882	Human	ELISA(Det)
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS00-82] - BSA and Azide free	HA750745	Human	WB
	Human IL-1 beta Recombinant Rabbit mAb [PS00-82] - BSA and Azide free (Detector)	HA721279	Human	ELISA(Det)
	Biotin Conjugated Human IL-1 beta Recombinant Rabbit mAb [PS00-82]	HA722649B	Human	ELISA(Det)
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-85] - BSA and Azide free (Capture)	HA723881	Human	ELISA(Cap)
	Biotin Conjugated Human IL-1 beta Recombinant Rabbit mAb [PSH16-86] - Detector	HA723882B	Human	ELISA(Det),ELISA
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS01-39] - BSA and Azide free	HA750744	Human	WB
	Human IL-1 beta Recombinant Rabbit mAb [PS01-39] - BSA and Azide free (Capture)	HA721278	Human	ELISA(Cap)
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS00-82]	HA721713	Human	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human IL-2 Mouse mAb [A8B7] - BSA and Azide free (Detector)	HA601006	Human	WB,ELISA(Det)
	Human IL-2 Recombinant Rabbit mAb [PSH01-42] - BSA and Azide free (Capture)	HA721680	Human	ELISA(Cap)
	Biotin Conjugated Human IL-2 Mouse mAb [A8B7]	HA601361B	Human	ELISA(Det)
	Human IL-23A Recombinant Rabbit mAb [PSH03-03] - BSA and Azide free (Capture)	HA721933	Human	ELISA(Cap)
	Human IL-2 Receptor alpha Recombinant Rabbit mAb [PSH05-65] - BSA and Azide free (Capture)	HA722385	Human	ELISA(Cap)
	Human IL-2 Receptor alpha Recombinant Rabbit mAb [PS00-52] - BSA and Azide free (Capture)	HA721280	Human	ELISA(Cap)
	Biotin Conjugated Human IL-2 Receptor alpha Recombinant Rabbit mAb [PS01-03]	HA722650B	Human	ELISA(Det)
	Human IL-2 Receptor alpha Recombinant Rabbit mAb [PSH05-66] - BSA and Azide free (Detector)	HA722386	Human	ELISA(Det)
	Human IL-2 Receptor alpha Recombinant Rabbit mAb [PS01-03] - BSA and Azide free (Detector)	HA721281	Human	ELISA(Det)
	Biotin Conjugated Human IL-2 Receptor alpha Recombinant Rabbit mAb [PSH05-66]	HA722594B	Human	ELISA(Det)
	Biotin Conjugated Human IL-4 Recombinant Rabbit mAb [PS00-66] - Detector	HA721368	Human	ELISA(Det)
	Human IL-4 Recombinant Rabbit mAb [PS00-66] - BSA and Azide free (Detector)	HA721231	Human	ELISA(Det)
	Human IL-4 Recombinant Rabbit mAb [PS00-67] - BSA and Azide free (Capture)	HA721230	Human	ELISA(Cap)
	HRP Conjugated Human IL-4 Recombinant Rabbit mAb [PS00-66] - Detector	HA721295	Human	ELISA(Det)
	Human IL-5 Recombinant Rabbit mAb [PSH04-64] - BSA and Azide free (Capture)	HA722164	Human	ELISA(Cap)
	Human IL-5 Recombinant Rabbit mAb [PSH04-65] - BSA and Azide free (Detector)	HA722165	Human	ELISA(Det)
	Biotin Conjugated Human IL-5 Recombinant Rabbit mAb [PSH04-65] - Detector	HA722609B	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human IL-6 Recombinant Rabbit mAb [PS00-21] - BSA and Azide free (Capture)	HA721167	Human	ELISA(Cap)
	HRP Conjugated Human IL-6 Recombinant Rabbit mAb [PS00-15] - Detector	HA721293	Human	ELISA(Det)
	Human IL-6 Recombinant Rabbit mAb [PS00-15] - BSA and Azide free (Detector)	HA721168	Human	ELISA(Det)
	Biotin Conjugated Human IL-6 Recombinant Rabbit mAb [PS00-15] - Detector	HA722233B	Human	ELISA(Det)
	Human IL-6 Mouse mAb [A6D1]	HA601051	Human	WB,ELISA
	Human IL-6 Mouse mAb [A6D2]	HA601050	Human	WB,ELISA
	Human IL-7 Recombinant Rabbit mAb [PSH03-09] - BSA and Azide free (Capture)	HA721937	Human	ELISA(Cap)
	Human IL-7 Recombinant Rabbit mAb [PSH03-10] - BSA and Azide free (Detector)	HA721938	Human	ELISA(Det)
	Biotin Conjugated Human IL-7 Recombinant Rabbit mAb [PSH03-10]	HA722632B	Human	ELISA(Det)
	Human Interferon gamma Recombinant Rabbit mAb [PSH01-46] - BSA and Azide free (Capture)	HA721684	Human	ELISA(Cap)
	Human Interferon gamma Recombinant Rabbit mAb [PS01-75] - BSA and Azide free (Capture)	HA721574	Human	ELISA(Cap)
	Biotin Conjugated Human Interferon gamma Recombinant Rabbit mAb [PSH01-47]	HA722645B	Human	ELISA(Det)
	Human Interferon gamma Recombinant Rabbit mAb [PS01-76] - BSA and Azide free (Detector)	HA721575	Human	ELISA(Det)
	Human Interferon gamma Recombinant Rabbit mAb [PSH01-45] - BSA and Azide free (Capture)	HA721683	Human	ELISA(Cap)
	Biotin Conjugated Human Interferon gamma Recombinant Rabbit mAb [PS01-76] - Detector	HA721576	Human	ELISA(Det),ELISA
	Human Interferon gamma Recombinant Rabbit mAb [PSH01-47] - BSA and Azide free (Detector)	HA721685	Human	ELISA(Det)
	Biotin Conjugated Human KIM-1/HAVCR1 Recombinant Rabbit mAb [PSH18-81]	HA725250B	Human	ELISA(Det),ELISA
	Human KIM-1/HAVCR1 Recombinant Rabbit mAb [PSH18-81] - BSA and Azide free (Detector)	HA725250	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human KIM-1/HAVCR1 Recombinant Rabbit mAb [PSH18-80] - BSA and Azide free (Capture)	HA725249	Human	ELISA(Cap)
	Human LAG-3 Recombinant Rabbit mAb [PSH05-29] - BSA and Azide free (Capture)	HA722258	Human	ELISA(Cap)
	Biotin Conjugated Human LAG-3 Recombinant Rabbit mAb [PSH05-30]	HA722596B	Human	ELISA(Det)
	Human LAG-3 Recombinant Rabbit mAb [PSH05-30] - BSA and Azide free (Detector)	HA722259	Human	ELISA(Det)
	Human Lipocalin-2 Recombinant Rabbit mAb [PSH13-68] - BSA and Azide free (Detector)	HA725110	Human	ELISA(Det)
	HRP Conjugated Human Lipocalin-2 Recombinant Rabbit mAb [PSH13-68] - Detector	HA725111H	Human	ELISA(Det),ELISA
	Human Lipocalin-2 Recombinant Rabbit mAb [PSH13-67] - BSA and Azide free (Capture)	HA725109	Human	ELISA(Cap)
	Biotin Conjugated Human Lipocalin-2 Recombinant Rabbit mAb [PSH13-68] - Detector	HA725111B	Human	ELISA(Det),ELISA
	HRP Conjugated Human MERTK Recombinant Rabbit mAb [PSH12-70] - Detector	HA723481H	Human	ELISA(Det),ELISA
	Human MERTK Recombinant Rabbit mAb [PSH12-69] - BSA and Azide free (Capture)	HA723479	Human	ELISA(Cap)
	Biotin Conjugated Human MERTK Recombinant Rabbit mAb [PSH12-70] - Detector	HA723481B	Human	ELISA(Det),ELISA
	Human MERTK Recombinant Rabbit mAb [PSH12-70] - BSA and Azide free (Detector)	HA723480	Human	ELISA(Det)
	Biotin Conjugated Human MICA Recombinant Rabbit mAb [PSH11-40] - Detector	HA723329B	Human	ELISA(Det),ELISA
	Human MICA Recombinant Rabbit mAb [PSH11-39] - BSA and Azide free (Capture)	HA723327	Human	ELISA(Cap)
	Human MICA Recombinant Rabbit mAb [PSH11-40] - BSA and Azide free (Detector)	HA723328	Human	ELISA(Det)
	Human MICB Recombinant Rabbit mAb [PSH10-73] - BSA and Azide free (Capture)	HA723238	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human MICB Recombinant Rabbit mAb [PSH10-74] - BSA and Azide free (Detector)	HA723239	Human	ELISA(Det)
	Biotin Conjugated Human MICB Recombinant Rabbit mAb [PSH10-74]	HA723240B	Human	ELISA(Det),ELISA
	Human NCAM1 Recombinant Rabbit mAb [PSH08-08] - BSA and Azide free (Capture/Detector)	HA722947	Human	ELISA(Det),ELISA(Cap)
	Biotin Conjugated Human NCAM1 Recombinant Rabbit mAb [PSH08-09] - Detector	HA722950B	Human	ELISA(Det)
	Human NCAM1 Recombinant Rabbit mAb [PSH08-07] - BSA and Azide free (Capture)	HA722946	Human	ELISA(Cap)
	Human NCAM1 Recombinant Rabbit mAb [PSH08-09] - BSA and Azide free (Detector)	HA722949	Human	ELISA(Det)
	Biotin Conjugated Human NCAM1 Recombinant Rabbit mAb [PSH08-08] - Detector	HA722948B	Human	ELISA(Det)
	PE Conjugated Human NCAM1 / CD56 Recombinant Rabbit mAb [PSH06-53]	HA720215F	Human	FC
	APC Conjugated Human NCAM1 / CD56 Recombinant Rabbit mAb [PSH06-53]	HA723757F	Human	FC
	Human Neuropilin-1 Recombinant Rabbit mAb [PSH14-73] - BSA and Azide free (Capture)	HA725112	Human	ELISA(Cap)
	Human Neuropilin-1 Recombinant Rabbit mAb [PSH14-74] - BSA and Azide free (Detector)	HA725113	Human	ELISA(Det)
	HRP Conjugated Human Neuropilin-1 Recombinant Rabbit mAb [PSH14-74] - Detector	HA725114H	Human	ELISA(Det),ELISA
	Human OX40/TNFRSF4 Recombinant Rabbit mAb [PSH10-17] - BSA and Azide free (Capture)	HA723187	Human	ELISA(Cap)
	Human OX40/TNFRSF4 Recombinant Rabbit mAb [PSH10-18] - BSA and Azide free (Detector)	HA723188	Human	ELISA(Det)
	Biotin Conjugated Human OX40/TNFRSF4 Recombinant Rabbit mAb [PSH10-18] - Detector	HA723189B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human PD1 Recombinant Rabbit mAb [PSH13-39] - Detector	HA723525B	Human	ELISA(Det),ELISA
	Human PD1 Recombinant Rabbit mAb [SJ01-91] - BSA and Azide free (Capture)	HA723521	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human PD1 Recombinant Rabbit mAb [PSH13-38] - Detector	HA723523B	Human	ELISA(Det),ELISA
	Human PD1 Recombinant Rabbit mAb [PSH13-38] - BSA and Azide free (Detector)	HA723522	Human	ELISA(Det)
	Human PD1 Recombinant Rabbit mAb [PSH13-39] - BSA and Azide free (Detector)	HA723524	Human	ELISA(Det)
	Human RANKL Recombinant Rabbit mAb [PSH16-59] - BSA and Azide free (Detector)	HA723861	Human	ELISA(Det)
	Human RANKL Recombinant Rabbit mAb [PSH16-60] - BSA and Azide free (Capture)	HA723862	Human	ELISA(Cap)
	Biotin Conjugated Human RANKL Recombinant Rabbit mAb [PSH16-59] - Detector	HA723861B	Human	ELISA(Det),ELISA
	Human RANKL Recombinant Rabbit mAb [PSH16-58] - BSA and Azide free (Capture)	HA723860	Human	ELISA(Cap)
	Human TNF alpha Recombinant Rabbit mAb	IRS040	Human	mIHC
	Human TNF alpha Recombinant Rabbit mAb [PS00-40] - BSA and Azide free (Capture)	HA721265	Human	ELISA(Cap)
	Human TNF alpha Recombinant Rabbit mAb [PS00-40] - BSA and Azide free	HA750593	Human	IHC-P
	Human TNF alpha Recombinant Rabbit mAb [PS01-07] - BSA and Azide free (Detector)	HA721266	Human	ELISA(Det)
	Human TNF alpha Recombinant Rabbit mAb [PS00-40]	HA721286	Human	IHC-P,mIHC
	Biotin Conjugated Human TNF alpha Recombinant Rabbit mAb [PS01-07] - Detector	HA721559	Human	ELISA(Det)
	HRP Conjugated Human TNF alpha Recombinant Rabbit mAb [PS01-07] - Detector	HA721294	Human	ELISA(Det)
	Biotin Conjugated Human ULBP1 Recombinant Rabbit mAb [PSH11-97] - Detector	HA723410B	Human	ELISA(Det),ELISA
	Human ULBP1 Recombinant Rabbit mAb [PSH11-96] - BSA and Azide free (Detector)	HA723407	Human	ELISA(Det)
	Human ULBP1 Recombinant Rabbit mAb [PSH11-95] - BSA and Azide free (Capture)	HA723406	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human ULBP1 Recombinant Rabbit mAb [PSH11-96] - Detector	HA723408B	Human	ELISA(Det),ELISA
	Human ULBP1 Recombinant Rabbit mAb [PSH11-97] - BSA and Azide free (Detector)	HA723409	Human	ELISA(Det)
KD	ICAM1 Mouse mAb [F11-A4]	M1505-7	Human	WB,IF-Cell,FC
	ICAM1 Rabbit pAb	ER131207	Human	WB,IF-Cell,IHC-P,FC
Relative expression	ICAM1 Recombinant Rabbit mAb [ST0487] - BSA and Azide free	HA750176	Human	WB,IF-Tissue,IHC-P
	ICAM1 Mouse mAb [J0-F7]	M1511-6	Human	WB,IF-Cell,FC
Relative expression, KD	ICAM1 Recombinant Rabbit mAb [PSH09-31]	HA723089	Human	WB,IHC-P,IF-Tissue
Relative expression	ICAM1 Recombinant Rabbit mAb [ST0487]	ET1609-46	Human	WB,IF-Tissue,IHC-P
Cell treatment, Relative expression	ICOSLG / B7-H2 / CD275 Recombinant Rabbit mAb [JE52-85]	HA721299	Human	WB
Cell treatment	IDO Recombinant Rabbit mAb [PD00-62]	HA721331	Human	IHC-P,IF-Cell,WB,mIHC
Cell treatment, Relative expression	IDO Recombinant Rabbit mAb [PSH12-31] - BSA and Azide free	HA751445	Human	WB,IF-Cell,IHC-P,FC,IP
	IDO Recombinant Rabbit mAb	IRS045	Human	mIHC
	IDO Rabbit pAb	ER1706-94	Human, Mouse	WB,IF-Cell,IHC-P,FC
Cell treatment, Relative expression	IDO Recombinant Rabbit mAb [PSH12-31]	HA723456	Human	WB,IF-Cell,IHC-P,FC,IP
	IDO Rat mAb [1G1]	RT1318	Human, Mouse	WB,IP,IF,IHC-P
	IgA Mouse mAb [3-F5]	HA601137	Human	ELISA,WB
	IgA Recombinant Rabbit mAb [JM10-42]	ET1703-10	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP
	IgA Mouse mAb [5-C3]	M1510-30	Human	WB,ELISA
	IgA Mouse mAb [3-B6]	M1510-31	Human	WB,ELISA
	IgG4 Recombinant Rabbit mAb [ST05-12]	ET1609-56	Human	WB,IHC-P
	IL-10 Rabbit pAb	ER1911-19	Human, Mouse	WB,IHC-P,IF-Tissue
	IL-10 Rabbit pAb	ER65464	Human, Mouse, Rat	WB,IHC-P
Cell treatment	IL-10 Recombinant Rabbit mAb [JE59-71]	HA722032	Human, Mouse, Rat	WB,IF-Cell,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	IL12B Rabbit pAb	ER1803-31	Human, Mouse	WB
	IL-18 Recombinant Rabbit mAb [PS01-06] - BSA and Azide free	HA750677	Human	WB,FC
Relative expression	IL-18 Recombinant Rabbit mAb [PSH18-15]	HA723978	Human, Mouse, Rat	WB,IHC-P,IP
	IL-18 Recombinant Rabbit mAb [PS00-41] - BSA and Azide free	HA750676	Human	WB,IHC-P
	IL-18 Recombinant Rabbit mAb [PS00-41]	HA721536	Human	WB,IHC-P
	IL-18 Recombinant Rabbit mAb [PS01-06]	HA721537	Human	WB,FC
	IL-1 alpha Rabbit pAb	ER1802-36	Human, Mouse	IHC-P,FC,WB
Cell treatment	IL-1 alpha Recombinant Rabbit mAb [PSH15-65]	HA723772	Mouse	WB
	IL-1 alpha Recombinant Rabbit mAb [PSH15-66]	HA723773	Mouse	IF-Cell,FC
Cell treatment	IL-1 alpha Recombinant Rabbit mAb [PSH15-65] - BSA and Azide free	HA751576	Mouse	WB
	IL-1 alpha Recombinant Rabbit mAb [PSH15-66] - BSA and Azide free	HA751577	Mouse	IF-Cell,FC
	IL-1 beta Recombinant Rabbit mAb [PSH18-00]	HA723965	Mouse, Rat	WB
Cell treatment	IL-1 beta Recombinant Rabbit mAb [JJ087-3]	ET1701-39	Human	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	IL-1 beta Recombinant Mouse mAb [A7F7-R]	HA601036	Human, Mouse, Rat	WB,ELISA,IHC-P,IF-Tissue
Cell treatment	IL-1 beta Mouse mAb [A7A10]	HA601035	Human	WB,ELISA,IF-Cell
Cell treatment	IL-1 beta Mouse mAb [A7F9]	HA601002	Human, Mouse, Rat	WB,ELISA,IF-Cell
Cell treatment	IL-1 beta Recombinant Mouse mAb [A7F7-R] - BSA and Azide free	HA610004	Human, Mouse, Rat	WB,ELISA,IHC-P,IF-Tissue
	IL-2 Rabbit pAb	ER63748	Human, Rat	WB,IHC-P
	IL-2 Mouse mAb [A8B9]	HA601016	Human	WB,ELISA
	IL-2 Mouse mAb [A8B8]	HA601015	Human	WB,ELISA
	IL-23 Recombinant Rabbit mAb [JE38-99]	HA721037	Human	WB,IHC-P,IF-Tissue
	IL-2 receptor alpha Rabbit pAb	HA500076	Human, Rat	WB
	IL-4 Recombinant Rabbit mAb [JA86-10]	ET1704-28	Human	WB,IF-Cell,IF-Tissue,IHC-P
	IL-4 Mouse mAb [A7D7]	HA601034	Human	ELISA
	IL-4 Rabbit pAb	ER1706-55	Human	WB,IF-Cell,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	IL-5 Rabbit pAb	HA500081	Human	WB,IF-Cell,ELISA
	IL-6 Mouse mAb [1-6]	EM1701-45	Mouse, Rat	WB,IF-Cell
	IL-6 Mouse mAb	EM30301	Human	WB
	IL-6 Recombinant Rabbit mAb	IRS160RB	Human	mIHC
	IL-6 Rabbit pAb	R1203-2	Mouse	WB,ELISA
	IL-6 Mouse mAb [3-A4]	M1303-3	Mouse	WB,IF-Cell,FC,ELISA
	IL-6 Rabbit pAb	R1412-2	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	IL-7 Mouse mAb [A8B11]	HA601048	Human	WB
	IL-7 Rabbit pAb	HA500171	Human	WB,IF-Cell,IHC-P
	IL7 Mouse mAb [AH31-34]	EM170412	Human, Rat	ELISA,WB,ICC,IHC-P
	IL-7R alpha Rabbit pAb	ER1706-68	Human, Zebrafish	WB,IHC-P,FC
KD	Integrin alpha 2 Recombinant Rabbit mAb [SN0752] - BSA and Azide free	HA750263	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
KD	Integrin alpha 2 Recombinant Rabbit mAb [SN0752]	ET1611-57	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	Integrin alpha E Recombinant Rabbit mAb [SN62-06]	ET1611-27	Human	IHC-P
	Interferon gamma Recombinant Rabbit mAb [JM10-10]	ET1703-17	Human	WB,FC
	Interferon gamma Recombinant Rabbit mAb [JM10-10] - BSA and Azide free	HA750373	Human	WB,FC
	intestinal alkaline phosphatase Recombinant Rabbit mAb [PSH17-83]	HA723952	Human, Mouse, Rat	IHC-P
	IRF5 Recombinant Rabbit mAb [SN201-05] - BSA and Azide free	HA750277	Human	WB,IHC-P,FC
	IRF5 Recombinant Rabbit mAb [SD203-07]	ET1612-15	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,IP
	IRF5 Recombinant Rabbit mAb [SN201-05]	ET1611-94	Human	WB,IHC-P,FC
Cell treatment	LAG-3 Recombinant Rabbit mAb [PS01-33]	HA721358	Human	IHC-P,IF-Cell,WB,mIHC
Cell treatment	LAG-3 Recombinant Rabbit mAb [PSH0-43]	HA721345	Human, Mouse	WB,IHC-P
Cell treatment	LAG-3 Recombinant Rabbit mAb [PSH0-44]	HA721346	Human, Mouse	WB,IHC-P,IF-Cell
Cell treatment	LAG-3 Recombinant Rabbit mAb [PSH0-42]	HA721344	Mouse	WB,IF-Cell
Cell treatment	LAG-3 Recombinant Rabbit mAb [PSH0-43] - BSA and Azide free	HA750621	Human, Mouse	WB,IHC-P
	LAG-3 Recombinant Rabbit mAb	IRS020	Human	mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	LAG-3 Recombinant Rabbit mAb [PSH0-44] - BSA and Azide free	HA750622	Human, Mouse	WB,IHC-P,IF-Cell
Cell treatment	LAG-3 Recombinant Rabbit mAb [PSH0-42] - BSA and Azide free	HA750620	Mouse	WB,IF-Cell
Relative expression	LAIR1 Recombinant Mouse mAb [PSH04-50] - BSA and Azide free	HA610201	Human	WB,IHC-P
Relative expression	LAIR1 Recombinant Mouse mAb [PSH04-50]	HA601329	Human	WB,IHC-P
Relative expression	Lck Recombinant Rabbit mAb [PSH09-25]	HA723083	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC
	Lck Rabbit pAb	HA500280	Human, Rat	WB,FC
	Lck Rabbit pAb	ER63849	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Relative expression	Lck Recombinant Rabbit mAb [PSH09-25] - BSA and Azide free	HA751279	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC
	Lipocalin-2 Recombinant Rabbit mAb [JM10-67]	ET1703-39	Human	WB,IF-Tissue,IHC-P
Relative expression	Liver Arginase Recombinant Rabbit mAb [PSH09-78]	HA723143	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Cell treatment	Liver Arginase Recombinant Rabbit mAb [SY09-06]	ET1605-8	Human, Mouse, Rat	WB
	Liver Arginase Recombinant Rabbit mAb [PD00-05]	HA721147	Human	WB,IHC-P,IF-Tissue
Relative expression	Liver Arginase Recombinant Rabbit mAb [JB21-49]	ET7107-90	Human	WB,IHC-P,IP
	LSP1 Recombinant Rabbit mAb [SN06-25]	ET1611-9	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	LY6A Recombinant Rabbit mAb [JM53-11]	ET1703-67	Mouse	WB,IHC-P
	LY6A Recombinant Rabbit mAb	IRS166RB	Mouse	mIHC
Relative expression	Ly6c Recombinant Rabbit mAb [PSH04-51]	HA722134	Mouse	WB
	Ly6c Rabbit pAb	HA500087	Mouse	WB,IHC-P,FC
Relative expression	Ly6c Rabbit pAb	HA500088	Mouse	WB,IHC-P,FC
Relative expression	Ly6c Recombinant Rabbit mAb [PSH04-51] - BSA and Azide free	HA750958	Mouse	WB
	LY75 Recombinant Rabbit mAb [JB87-35]	ET7107-58	Human, Mouse	WB,IHC-P,IP
	Lysozyme Recombinant Mouse mAb [A5A8-R] - BSA and Azide free	HA610104	Human	WB,IHC-P
	Lysozyme Mouse mAb [A5A9]	EM1901-99	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Lysozyme Recombinant Mouse mAb [A5A8-R]	HA601229	Human	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Lysozyme Recombinant Rabbit mAb [ST50-02] - BSA and Azide free	HA750170	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Lysozyme Recombinant Rabbit mAb	IRS168RB	Mouse	mIHC
	Lysozyme Mouse mAb [A5A8]	EM1901-98	Human	WB,IHC-P
	Lysozyme Recombinant Rabbit mAb [ST50-02]	ET1609-35	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP,mIHC
	Lysozyme Rabbit pAb	ER1803-88	Human, Mouse, Rat	WB,IHC-P
	Mannose Receptor(CD206) Recombinant Antibody [PSH07-76] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610286	Mouse, Rat	IHC-Fr,WB
	Mannose Receptor(CD206) Recombinant Antibody [PSH07-76] - Guinea pig IgG2 (Chimeric)	HA601452	Mouse, Rat	IHC-Fr,WB
Cell treatment	Mannose Receptor(CD206) Recombinant Rabbit mAb [PSH07-76] - BSA and Azide free	HA751188	Mouse, Rat	WB,IHC-P,IF-Tissue,IHC-Fr,IF-Cell,IP
Cell treatment	Mannose Receptor(CD206) Recombinant Rabbit mAb [PSH07-76]	HA722892	Mouse, Rat	WB,IHC-P,IF-Tissue,IHC-Fr,IF-Cell,IP,mIHC
KD	Mannose Receptor(CD206) Recombinant Rabbit mAb [JF0953] - BSA and Azide free	HA750332	Human, Mouse	WB,IF-Cell,FC
	FITC Conjugated Mannose Receptor(CD206) Recombinant Rabbit mAb [JF0953]	HA720214F	Human	FC
	Mannose Receptor(CD206) Recombinant Antibody [PSH07-76] - Mouse IgG1 (Chimeric)	HA601453	Mouse, Rat	IHC-Fr,IF-Tissue,IHC-P,WB,IF-Cell
	Mannose Receptor(CD206) Recombinant Rabbit mAb	IRS049	Mouse	mIHC
KD	Mannose Receptor(CD206) Recombinant Rabbit mAb [JF0953]	ET1702-04	Human, Mouse	WB,IF-Cell,FC
	Mannose Receptor(CD206) Recombinant Antibody [PSH07-76] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610287	Mouse, Rat	IHC-Fr,IF-Tissue,IHC-P,WB,IF-Cell
KD	MERTK Recombinant Rabbit mAb [SR29-07] - BSA and Azide free	HA750036	Human	WB,IHC-P,IP,FC,IF-Tissue
KD	MERTK Recombinant Rabbit mAb [SR29-07]	ET1602-21	Human	WB,IHC-P,FC,IF-Tissue
	MHC Class II Recombinant Rabbit mAb [JA10-94] - BSA and Azide free	HA750404	Human	WB,IP,IF-Cell,IF-Tissue,IHC-P
	MHC Class II Recombinant Rabbit mAb [JA10-94]	ET1704-13	Human	WB,IP,IF-Cell,IF-Tissue,IHC-P,mIHC
Relative expression	MICA Recombinant Rabbit mAb [PSH13-60] - BSA and Azide free	HA751479	Human	WB,FC,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	MICA Recombinant Rabbit mAb [PSH13-60]	HA723549	Human	WB,FC,IP
	MICA Mouse mAb [A4C5]	HA600002	Human	WB
Relative expression	MICA + MICB Rabbit pAb	ER1803-95	Human	WB
Relative expression	MICA + MICB Recombinant Rabbit mAb [PSH10-88] - BSA and Azide free	HA751380	Human	WB,FC,IP
Relative expression	MICA + MICB Recombinant Rabbit mAb [PSH10-88]	HA723255	Human	WB,FC,IP
	Mouse anti-Human IgG Mouse mAb [A5D12]	HA600076	Human	WB,ELISA
	Mouse anti-Human IgG Recombinant Mouse mAb [A5D12-R]	HA601220	Human	WB,ELISA
	PE Conjugated Mouse CD16/32 Recombinant Rabbit mAb [PSH06-60]	HA720216F	Mouse	FC
	FITC Conjugated Anti-Mouse CD3 Recombinant Antibody [PSH04-98] - Rabbit IgG (Chimeric)	HA720199F	Mouse	FC
	PE Conjugated Anti-Mouse CD3 Recombinant Antibody [PSH04-98] - Rabbit IgG (Chimeric)	HA720200F	Mouse	FC
	APC Conjugated Anti-Mouse CD3 Recombinant Antibody [PSH04-98] - Rabbit IgG (Chimeric)	HA720201F	Mouse	FC
	Anti-Mouse CD3 Recombinant Antibody [PSH04-98] - Rabbit IgG (Chimeric)	HA723347	Mouse	FC
	PE Conjugated Anti-Mouse CD4 Recombinant Antibody [PSH04-99] - Rabbit IgG (Chimeric)	HA720203F	Mouse	FC
	FITC Conjugated Anti-Mouse CD4 Recombinant Antibody [PSH04-99] - Rabbit IgG (Chimeric)	HA720202F	Mouse	FC
	APC Conjugated Anti-Mouse CD4 Recombinant Antibody [PSH04-99] - Rabbit IgG (Chimeric)	HA720204F	Mouse	FC
	Mouse CD68 Recombinant Rabbit mAb [PSH12-26] - BSA and Azide free (Capture)	HA723449	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse CD68 Recombinant Rabbit mAb [PSH12-25] - Detector	HA723448B	Mouse	ELISA(Det),ELISA
	Biotin Conjugated Mouse CD68 Recombinant Rabbit mAb [PSH12-27] - Detector	HA723451B	Mouse	ELISA(Det),ELISA
	Mouse CD68 Recombinant Rabbit mAb [PSH12-24] - BSA and Azide free (Capture)	HA723446	Mouse	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse CD68 Recombinant Rabbit mAb [PSH12-27] - BSA and Azide free (Detector)	HA723450	Mouse	ELISA(Det)
	Mouse CD68 Recombinant Rabbit mAb [PSH12-25] - BSA and Azide free (Detector)	HA723447	Mouse	ELISA(Det)
	Mouse CD86 Recombinant Rabbit mAb [PSH13-65] - BSA and Azide free (Capture)	HA723555	Mouse	ELISA(Cap)
	Mouse CD86 Recombinant Rabbit mAb [PSH13-66] - BSA and Azide free (Detector)	HA723556	Mouse	ELISA(Det)
	Biotin Conjugated Mouse CD86 Recombinant Rabbit mAb [PSH13-66] - Detector	HA723557B	Mouse	ELISA(Det),ELISA
	Mouse CD86 Recombinant Rabbit mAb [PSH13-64] - BSA and Azide free (Detector)	HA723553	Mouse	ELISA(Det)
	Mouse CD86 Recombinant Rabbit mAb [PSH13-63] - BSA and Azide free (Capture)	HA723552	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse CD86 Recombinant Rabbit mAb [PSH13-64] - Detector	HA723554B	Mouse	ELISA(Det),ELISA
	Mouse CD8 alpha Recombinant Rat mAb [PSH16-78]	HA601543	Mouse	FC
	Mouse CD8 alpha Recombinant Rat mAb [PSH16-78] - BSA and Azide free	HA610359	Mouse	FC
	FITC Conjugated Anti-Mouse CD8 alpha Recombinant Antibody [PSH05-00] - Rabbit IgG (Chimeric)	HA720205F	Mouse	FC
	PE Conjugated Anti-Mouse CD8 alpha Recombinant Antibody [PSH05-00] - Rabbit IgG (Chimeric)	HA720206F	Mouse	FC
	Anti-Mouse CD8 alpha Recombinant Antibody [PSH05-00] - Rabbit IgG (Chimeric) - BSA and Azide free	HA751593	Mouse	FC,IHC-Fr
	APC Conjugated Anti-Mouse CD8 alpha Recombinant Antibody [PSH05-00] - Rabbit IgG (Chimeric)	HA720207F	Mouse	FC
	Anti-Mouse CD8 alpha Recombinant Antibody [PSH05-00] - Rabbit IgG (Chimeric)	HA723820	Mouse	FC,IHC-Fr
	Mouse IL-10 Recombinant Rabbit mAb [PS01-93] - BSA and Azide free (Detector)	HA721518	Mouse	ELISA(Det)
	Mouse IL-10 Recombinant Rabbit mAb [PS01-85] - BSA and Azide free (Capture)	HA721517	Mouse	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Mouse IL-10 Recombinant Rabbit mAb [PS01-93] - Detector	HA721519	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL12 p40 Recombinant Rabbit mAb [PSH04-88]	HA722603B	Mouse	ELISA(Det)
	Mouse IL12 p40 Recombinant Rabbit mAb [PSH04-88] - BSA and Azide free (Detector)	HA722194	Mouse	ELISA(Det)
	Mouse IL12 p40 Recombinant Rabbit mAb [PSH04-87] - BSA and Azide free (Capture)	HA722193	Mouse	ELISA(Cap)
	Mouse IL-1 alpha Recombinant Rabbit mAb [PSH13-41] - BSA and Azide free (Detector)	HA723527	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-1 alpha Recombinant Rabbit mAb [PSH13-41] - Detector	HA723528B	Mouse	ELISA(Det),ELISA
	Mouse IL-1 alpha Recombinant Rabbit mAb [PSH13-42] - BSA and Azide free (Capture)	HA723529	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-1 alpha Recombinant Rabbit mAb [PSH13-43] - Detector	HA723531B	Mouse	ELISA(Det),ELISA
	Mouse IL-1 alpha Recombinant Rabbit mAb [PSH13-43] - BSA and Azide free (Detector)	HA723530	Mouse	ELISA(Det)
	Mouse IL-1 alpha Recombinant Rabbit mAb [PSH13-40] - BSA and Azide free (Capture)	HA723526	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-1 beta Recombinant Rabbit mAb [PSH03-32] - Detector	HA722625B	Mouse	ELISA(Det),ELISA
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH03-32] - BSA and Azide free (Detector)	HA721977	Mouse	ELISA(Det)
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH09-59] - BSA and Azide free	HA723121	Mouse	ELISA
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH09-58] - BSA and Azide free	HA723120	Mouse	ELISA
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH04-66] - BSA and Azide free (Capture)	HA722166	Mouse	ELISA(Cap)
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH04-67] - BSA and Azide free (Detector)	HA722167	Mouse	ELISA(Det)
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH03-31] - BSA and Azide free (Capture)	HA721976	Mouse	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Mouse IL-1 beta Recombinant Rabbit mAb [PSH04-67]	HA722608B	Mouse	ELISA(Det)
	Mouse IL-6 Recombinant Rabbit mAb [PSH0-60] - BSA and Azide free (Detector)	HA721413	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PSH0-60] - Detector	HA722647B	Mouse	ELISA(Det)
	HRP Conjugated Mouse IL-6 Mouse mAb [3-A4]	M1303-2	Mouse	WB,ELISA
	Mouse IL-6 Recombinant Rabbit mAb [PSH0-58] - BSA and Azide free (Capture)	HA721411	Mouse	ELISA(Cap)
	Mouse IL-6 Recombinant Rabbit mAb [PS01-83] - BSA and Azide free (Detector)	HA721601	Mouse	ELISA(Det)
	Mouse IL-6 Recombinant Rabbit mAb [PS01-74] - BSA and Azide free (Capture/Detector)	HA721600	Mouse	ELISA(Cap),ELISA(Det)
	Mouse IL-6 Recombinant Rabbit mAb [PSH0-59] - BSA and Azide free (Detector)	HA721412	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PS01-83] - BSA and Azide free (Detector)	HA721602	Mouse	ELISA(Det)
	Mouse IL-6 Recombinant Rabbit mAb [PSH14-47] - BSA and Azide free (Capture)	HA723642	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PSH0-59] - Detector	HA722646B	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PS01-74] - Detector	HA723643B	Mouse	ELISA(Det)
Cell treatment	Biotin Conjugated Mouse RANKL Recombinant Rabbit mAb [PSH07-92]	HA722925B	Mouse	ELISA(Det)
Cell treatment	Mouse RANKL Recombinant Rabbit mAb [PSH07-90] - BSA and Azide free (Capture)	HA722922	Mouse	ELISA(Cap)
Cell treatment	Mouse RANKL Recombinant Rabbit mAb [PSH07-91] - BSA and Azide free (Capture)	HA722923	Mouse	ELISA(Cap)
Cell treatment	Mouse RANKL Recombinant Rabbit mAb [PSH07-92] - BSA and Azide free (Detector)	HA722924	Mouse	ELISA(Det)
	Biotin Conjugated Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-30]	HA722626B	Mouse, Rat	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-30] - BSA and Azide free (Detector)	HA721975	Mouse, Rat	ELISA(Det)
	Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-29] - BSA and Azide free (Capture)	HA721974	Mouse, Rat	ELISA(Cap)
	Mouse TNF alpha Recombinant Rabbit mAb [PSH03-77] - BSA and Azide free (Detector)	HA722041	Mouse	ELISA(Det)
	Biotin Conjugated Mouse TNF alpha Recombinant Rabbit mAb [PSH03-77]	HA722623B	Mouse	ELISA(Det)
	Mouse TNF alpha Recombinant Rabbit mAb [PSH03-76] - BSA and Azide free (Capture)	HA722040	Mouse	ELISA(Cap)
	MSR1 Rabbit pAb	ER1913-21	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
	MSR1 / CD204 Rabbit pAb	ER1902-71	Human	WB,FC
Relative expression	MUM1 Mouse mAb [A7D1]	HA601030	Human	IHC-P,WB
Relative expression	MUM1 Recombinant Rabbit mAb [JE52-60]	ET7110-44	Human	WB,IHC-P
	MUM1 Mouse mAb [A7C12]	HA601031	Human	IHC-P
	Munc13-4 Rabbit pAb	HA500052	Human, Rat	WB,IF-Cell,IHC-P
Relative expression	NCAM1 / CD56 Recombinant Mouse mAb [PDH0-07]	HA601126	Human	IHC-P,WB
	NCAM1 / CD56 Recombinant Rabbit mAb	IRS026	Human	mIHC
	NCAM1 / CD56 Recombinant Rabbit mAb [JF1021]	ET1702-43	Human, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP,mIHC
Relative expression	NCAM1 / CD56 Recombinant Rabbit mAb [PSH06-53]	HA722756	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
	NCAM1 / CD56 Recombinant Rabbit mAb [JF1021] - BSA and Azide free	HA750346	Human, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
Relative expression	NCAM1 / CD56 Recombinant Rabbit mAb [PSH06-82] - BSA and Azide free	HA751104	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Cell,FC,IF-Tissue
Relative expression	NCAM1 / CD56 Recombinant Rabbit mAb [PSH06-53] - BSA and Azide free	HA751105	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
	NCAM1 / CD56 Rabbit pAb	R1204-1	Human, Mouse, Rat	WB,IF-Cell
Relative expression	NCAM1 / CD56 Recombinant Rabbit mAb [PSH06-82]	HA722755	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Cell,FC,IF-Tissue
Relative expression	NCAM1 / CD56 Recombinant Rabbit mAb [PSH08-89]	HA723041	Human	WB,IHC-P,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression, KD	Neuropilin 1 Recombinant Rabbit mAb [ST05-30] - BSA and Azide free	HA750189	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
Relative expression, KD	Neuropilin 1 Recombinant Rabbit mAb [ST05-30]	ET1609-69	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Neuropilin 1 Mouse mAb [A3-D1]	M0910-5	Human, Mouse	WB,IF-Cell,IHC-P
	Neutrophil Elastase Recombinant Rabbit mAb [JF098-6] - BSA and Azide free	HA750360	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P
	Neutrophil Elastase Recombinant Rabbit mAb [JF098-6]	ET1702-78	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P
Relative expression	NFAT1 Recombinant Rabbit mAb [JA11-08] - BSA and Azide free	HA750405	Human	WB,IF-Tissue,IHC-P
Relative expression	NFAT1 Recombinant Rabbit mAb [JA11-08]	ET1704-14	Human	WB,IF-Tissue,IHC-P
	NFAT2 Recombinant Rabbit mAb [JA81-03]	ET1704-45	Human	WB,IHC-P,FC
	NK1.1 Recombinant Rabbit mAb [PSH08-47]	HA722993	Mouse	WB,IHC-P
	FITC Conjugated NK1.1 Recombinant Rabbit mAb [PSH08-46]	HA720219F	Mouse	FC
	NK1.1 Recombinant Rabbit mAb [PSH08-46]	HA722992	Mouse	FC
	NK1.1 Recombinant Rabbit mAb [PSH08-46] - BSA and Azide free	HA751221	Mouse	FC
	NK1.1 Recombinant Rabbit mAb [PSH08-47] - BSA and Azide free	HA751222	Mouse	WB,IHC-P
	NK-p44 Recombinant Rabbit mAb [JE61-04]	HA720070	Human	WB
	NLRC3 Rabbit pAb	ER1706-69	Human, Mouse	IF-Cell,IHC-P,FC
	NLRC3 Rabbit pAb	ER1706-45	Human, Mouse	WB,IF-Cell,IHC-P,FC
	NOX2 / gp91phox Recombinant Rabbit mAb [SN07-16] - BSA and Azide free	HA750255	Human, Mouse, Rat	WB
	NOX2 / gp91phox Recombinant Rabbit mAb [SN07-16]	ET1611-44	Human, Mouse, Rat	WB
Relative expression	NOX2 / gp91phox Recombinant Rabbit mAb [PSH10-59] - BSA and Azide free	HA751362	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IHC-Fr,IP
Relative expression	NOX2 / gp91phox Recombinant Rabbit mAb [PSH10-59]	HA723224	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IHC-Fr,IP
Relative expression	NOX2 / gp91phox Rabbit pAb	ER1511-34	Human, Mouse	IF-Cell,IHC-P,FC,WB
	NOXA2 Recombinant Rabbit mAb [JG86-31]	ET7108-72	Human, Mouse, Rat	WB,IHC-P,IP
	Orai1 Rabbit pAb	ER1803-11	Human, Mouse	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	PD1 Recombinant Rabbit mAb [SJ01-91]	ET1606-41	Human	WB,IHC-P
Cell treatment	PD1 Mouse mAb [1F2]	EM1707-60	Human	WB,IF-Cell,IHC-P,FC
	PD1 Mouse mAb	IRS190MS	Human	mIHC
	PD1 Recombinant Rabbit mAb [SJ01-91] - BSA and Azide free	HA750103	Human	WB,IHC-P
	PD1 Blocking Mouse mAb [3A11]	EM1701-17	Human	ELISA,FC
	Perforin Rabbit pAb	ER1803-77	Human	WB,IF-Cell,IHC-P,FC
Relative expression	Perforin-1 Recombinant Rabbit mAb [PSH15-92]	HA723797	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Relative expression	Perforin-1 Recombinant Rabbit mAb [PSH15-92] - BSA and Azide free	HA751583	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-GATA3 (S308) Recombinant Rabbit mAb [ST44-09]	ET1609-17	Human, Mouse, Rat	WB,IHC-P,IP,IF-Tissue
Cell treatment, Relative expression	Phospho-LYN (Y397)/LCK (Y394)/HCK (Y411)/BLK (Y389) Recombinant Rabbit mAb [PSH09-43]	HA723105	Human, Mouse, Rat	WB,IF-Cell
Cell treatment, Relative expression	Phospho-LYN (Y397)/LCK (Y394)/HCK (Y411)/BLK (Y389) Recombinant Rabbit mAb [PSH09-43] - BSA and Azide free	HA751290	Human, Mouse, Rat	WB,IF-Cell
	Phospho-SHIP (Y1020) Recombinant Rabbit mAb [JJ082-03]	ET1701-5	Human	WB,IF-Cell,IP
Cell treatment, Relative expression	Phospho-STAT1 (S727) Recombinant Rabbit mAb [SN67-04] - BSA and Azide free	HA750247	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment, Relative expression	Phospho-STAT1 (S727) Recombinant Rabbit mAb [SN67-04]	ET1611-20	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-STAT1 (Y701) Recombinant Rabbit mAb [PSH04-02]	HA722083	Human	WB,IF-Cell
Cell treatment	Phospho-STAT1 (Y701) Recombinant Rabbit mAb [PSH04-02] - BSA and Azide free	HA750925	Human	WB,IF-Cell
Cell treatment	Phospho-STAT6 (Y641) Recombinant Rabbit mAb [PSH06-98] - BSA and Azide free	HA751125	Human, Mouse	WB
Cell treatment	Phospho-STAT6 (Y641) Recombinant Rabbit mAb [PSH06-98]	HA722776	Human, Mouse	WB
	Pin1 Recombinant Rabbit mAb [PSH10-81]	HA723248	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
	Pin1 Recombinant Rabbit mAb [PSH10-81] - BSA and Azide free	HA751373	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	P-Selectin / CD62P Recombinant Rabbit mAb [PSH19-03]	HA724039	Human, Mouse, Rat	WB,IHC-P,FC,IHC-Fr,IP
	Biotin Conjugated P-Selectin / CD62P Recombinant Rabbit mAb [PSH16-71] - Detector	HA723873B	Human, Mouse, Rat	ELISA(Det),ELISA
	HRP Conjugated P-Selectin / CD62P Recombinant Rabbit mAb [PSH16-71] - Detector	HA723873H	Human, Mouse, Rat	ELISA(Det),ELISA
	P-Selectin / CD62P Recombinant Rabbit mAb [PSH16-71] - BSA and Azide free (Detector)	HA723873	Human, Mouse, Rat	ELISA(Det)
	P-Selectin / CD62P Recombinant Rabbit mAb [PSH16-70] - BSA and Azide free (Capture)	HA723872	Human, Mouse, Rat	ELISA(Cap)
	PSMB10 Rabbit pAb	HA500393	Human, Mouse	WB,IHC-P,IF-Cell
Cell treatment	PSMB10 Recombinant Rabbit mAb [PSH01-58] - BSA and Azide free	HA750734	Human, Mouse	WB,IHC-P
Cell treatment	PSMB10 Recombinant Rabbit mAb [PSH01-58]	HA721701	Human, Mouse	WB,IHC-P
	PTPRC Mouse mAb [C11-C10]	EM1708-69	Human	WB,IHC-P
	PU.1/Spi1 Recombinant Rabbit mAb [JB72-35]	ET7107-49	Human	WB,IHC-P,FC
	PU.1 / Spi1 Recombinant Rabbit mAb [PD00-14]	HA721162	Human	WB,IHC-P
	Rab 2A Rabbit pAb	ER64714	Human, Mouse, Rat	WB,IHC-P
	HRP Conjugated Rabbit anti-Human IgG polyclonal Antibody	HA1024	Human	WB,ELISA
	RANKL Rabbit pAb	HA500369	Human, Mouse, Rat	WB,IF-Cell
	Biotin Conjugated Rat IL-10 Recombinant Rabbit mAb [PSH02-22]	HA722637B	Rat	ELISA(Det)
	Rat IL-10 Recombinant Rabbit mAb [PSH02-21] - BSA and Azide free (Capture)	HA721786	Rat	ELISA(Cap)
	Rat IL-10 Recombinant Rabbit mAb [PSH02-22] - BSA and Azide free (Detector)	HA721787	Rat	ELISA(Det)
	Biotin Conjugated Rat IL-1 beta Recombinant Rabbit mAb [PSH02-60]	HA722634B	Rat	ELISA(Det)
	Rat IL-1 beta Recombinant Rabbit mAb [PSH02-60] - BSA and Azide free (Detector)	HA721837	Rat	ELISA(Det)
	Rat IL-1 beta Recombinant Rabbit mAb [PS01-87] - BSA and Azide free (Detector)	HA721610	Rat	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Rat IL-1 beta Recombinant Rabbit mAb [PSH02-59] - BSA and Azide free (Capture)	HA721836	Rat	ELISA(Cap)
	Biotin Conjugated Rat IL-1 beta Recombinant Rabbit mAb [PS01-87] - Detector	HA721611	Rat	ELISA(Det),ELISA
	Rat IL-1 beta Recombinant Rabbit mAb [PS01-86] - BSA and Azide free (Capture)	HA721609	Rat	ELISA(Cap)
	Rat KIM-1 / HAVCR1 Recombinant Rabbit mAb [PSH20-50] - BSA and Azide free (Detector)	HA725382	Rat	ELISA(Det)
	Biotin Conjugated Rat KIM-1 / HAVCR1 Recombinant Rabbit mAb [PSH20-50] - Detector	HA725382B	Rat	ELISA(Det),ELISA
	Rat KIM-1 / HAVCR1 Recombinant Rabbit mAb [PSH20-49] - BSA and Azide free (Capture)	HA725381	Rat	ELISA(Cap)
	Biotin Conjugated Rat TNF alpha Recombinant Rabbit mAb [PSH04-71]	HA722606B	Rat	ELISA(Det)
	Rat TNF alpha Recombinant Rabbit mAb [PSH04-70] - BSA and Azide free (Capture)	HA722170	Rat	ELISA(Cap)
	Rat TNF alpha Recombinant Rabbit mAb [PSH04-71] - BSA and Azide free (Detector)	HA722171	Rat	ELISA(Det)
	SCA-1 Rabbit pAb	0804-10	Mouse	WB,IF-Cell
	SET Recombinant Rabbit mAb [JB60-34]	ET7107-18	Human, Mouse	WB,IF-Cell,IHC-P,IP,IF-Tissue
	SHIP Recombinant Rabbit mAb [SY11-08]	ET1605-30	Human	WB,IHC-P,IP,FC
	STAT1 Recombinant Mouse mAb [G3-B11-R] - BSA and Azide free	HA610181	Human	WB,IHC-P
	STAT1 Rabbit pAb	R1408-2	Human, Mouse, Rat, Zebrafish	WB,IF-Cell
	STAT1 Recombinant Mouse mAb [G3-B11-R]	HA601309	Human	WB,IHC-P
	STAT1 Mouse mAb [G3-B11]	M1407-1	Human	WB,IF-Cell,IHC-P,IF-Tissue
KD	STAT1 Recombinant Rabbit mAb [SD20-75]	ET1612-22	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
KD	STAT1 alpha Recombinant Rabbit mAb [SJ01-89]	ET1606-39	Human, Mouse	WB,IP,FC,IHC-P,IF-Cell,IF-Tissue
KD	STAT1 alpha Recombinant Rabbit mAb [SJ01-89] - BSA and Azide free	HA750102	Human, Mouse	WB,IP,FC,IHC-P,IF-Cell,IF-Tissue
	STAT1 alpha/beta Rabbit pAb	ER31221	Human, Mouse, Rat	WB,IF-Cell,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	STAT6 Recombinant Rabbit mAb [SY02-72]	ET1607-61	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	STAT6 Recombinant Rabbit mAb [SY13-09]	ET1605-49	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC,IF-Tissue
KD	STAT6 Rabbit pAb	ER40121	Human	WB,IF-Cell,IHC-P,FC
KD	TAPA1/CD81 Recombinant Rabbit mAb [SN206-01]	ET1611-87	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Relative expression	T-bet / Tbx21 Recombinant Rabbit mAb [JE60-20]	HA721901	Human	WB,IHC-P,IF-Cell,FC
	T-bet / Tbx21 Rabbit pAb	ER1802-65	Human	WB,FC,IHC-P
	TBP-1 Rabbit pAb	HA500413	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
KO	TDG Recombinant Rabbit mAb [PSH10-11] - BSA and Azide free	HA751531	Human, Mouse, Rat	WB,IHC-P
	TDG Rabbit pAb	ER65129	Human, Mouse, Rat	WB,IHC-P,IF-Cell
KO	TDG Recombinant Rabbit mAb [PSH10-11]	HA723180	Human, Mouse, Rat	WB,IHC-P
	Telomerase reverse transcriptase Recombinant Rabbit mAb [JF0586]	ET1702-22	Human	WB,FC
	TERT Mouse mAb [G9-H2]	EM1707-13	Human	IHC-P,FC
	TFE3 Recombinant Rabbit mAb [JE60-60]	HA721896	Human, Mouse, Rat	WB,IHC-P,FC
	TFE3 Rabbit pAb	ER65143	Human, Mouse	WB,IHC-P
Relative expression	TGF beta 1 Recombinant Rabbit mAb [PD00-17]	HA721143	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Tissue
Relative expression	TGF beta 1 Recombinant Rabbit mAb [PD00-17] - BSA and Azide free	HA750485	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Tissue
	TGF beta 1 Mouse mAb [B2-F6]	EM1709-93	Human	WB,IHC-P,FC
	TGF beta 1 Rabbit pAb	HA500496	Human	WB
	Biotin Conjugated TGF beta 1 Recombinant Rabbit mAb [PSH15-22] - Detector	HA723727B	Human, Mouse, Rat	ELISA(Det),ELISA
	TGF beta 1 Recombinant Rabbit mAb	IRS205RB	Human	mIHC
	TGF beta 1 Recombinant Rabbit mAb [PSH15-21] - BSA and Azide free (Capture)	HA723725	Human, Mouse, Rat	ELISA(Cap)
	TGF beta 1 Recombinant Rabbit mAb [PSH18-59]	HA724012	Human, Mouse, Monkey	WB,IF-Cell,IHC-P,FC
	TGF beta 1 Recombinant Rabbit mAb [PSH15-22] - BSA and Azide free (Detector)	HA723726	Human, Mouse, Rat	ELISA(Det)
	Themis Recombinant Rabbit mAb [JE55-73]	ET7111-20	Human, Rat	WB,IF-Cell,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Thy-1 (OX7) Mouse mAb [1G2]	RT1615	Mouse, Rat	WB,IP,IF,IHC-P,FC
	TIA1 Recombinant Rabbit mAb [JM42-11]	ET1703-59	Human, Mouse	WB,IF-Tissue,IHC-P,IP
	TIA1 Recombinant Rabbit mAb [JM42-11] - BSA and Azide free	HA750385	Human, Mouse	WB,IF-Tissue,IHC-P,IP
	TIM 1 Recombinant Rabbit mAb [JE39-41]	HA721535	Human, Mouse, Rat	WB
	TLR4 Mouse mAb	RT1666	Human, Mouse, Rat	WB
	TLR4 Rabbit pAb	ER1706-43	Human, Mouse, Rat	IF-Cell,IHC-P,FC
Cell treatment	TNF alpha Recombinant Rabbit mAb [PS01-07]	HA722022	Human, Mouse	WB
	TNF alpha Rabbit pAb	ER1919-22	Human, Mouse	WB
Cell treatment	TNF alpha Recombinant Rabbit mAb [PS01-07] - BSA and Azide free	HA750897	Human, Mouse	WB
Cell treatment	TNF alpha Rabbit pAb	R1203-1	Mouse	WB,IF-Cell,ELISA
	TNF alpha Rabbit pAb	ER65189	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	TNFRSF18 Mouse mAb [4H2D6]	EM1711-16	Human	WB,IHC-P,FC
	TNFSF11 Mouse mAb [8A7B9]	EM1710-42	Human, Monkey	WB,IHC-P,IF-Cell,FC
	TNFSF14 Rabbit pAb	HA500384	Human, Mouse	WB,IHC-P
	TNR18 Rabbit pAb	ER62416	Human	WB
KD	TRK fused gene Recombinant Rabbit mAb [SD082-0]	ET1612-52	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	ULBP1 Recombinant Rabbit mAb [JE65-75]	HA720092	Human	WB
Relative expression	ZAP70 Recombinant Mouse mAb [A1B5-R]	HA601218	Human, Mouse, Rat	WB,IHC-P
Relative expression	ZAP70 Mouse mAb [A1B7]	EM1902-28	Human, Mouse, Rat	WB,IF-Cell,FC
Relative expression	ZAP70 Recombinant Rabbit mAb [PSH11-52]	HA723344	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Relative expression	ZAP70 Recombinant Rabbit mAb [JU08-39]	ET1706-42	Human	WB,IF-Cell,IHC-P,FC,IP
Relative expression	ZAP70 Mouse mAb [A1B5]	EM1901-27	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell
Relative expression	ZAP70 Recombinant Mouse mAb [A1B5-R] - BSA and Azide free	HA610096	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层