

脓毒症研究-产品清单

高级验证	产品名	货号	反应种属	验证应用
	6CKine Rabbit pAb	ER65800	Human	IHC-P
Cell treatment, KD	Active+Pro Caspase-9 Rabbit pAb	R1308-12	Human	WB,IHC-P
Cell treatment, KD	Active+Pro Caspase-9 Recombinant Rabbit mAb [SC65-05]	ET1610-95	Human, Mouse	WB,IP
Cell treatment, KD	Active+Pro Caspase-9 Recombinant Rabbit mAb [JJ08-05]	ET1701-22	Human	WB,IHC-P,FC
Cell treatment, KD	Active+Pro Caspase-9 Recombinant Rabbit mAb [SC65-05] - BSA and Azide free	HA750238	Human, Mouse	WB,IP
KD	ALIX Recombinant Rabbit mAb [JM85-31]	ET1705-74	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell
KD	ALIX Recombinant Rabbit mAb [JM85-31] - BSA and Azide free	HA750439	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell
	Angiotensinogen Recombinant Rabbit mAb [SD201-02]	ET1612-48	Human	WB,IF-Cell,IF-Tissue,IP
KD	Annexin V Recombinant Rabbit mAb [JF50-11]	ET1702-62	Human, Mouse, Rat	WB,FC
	Apolipoprotein A II Recombinant Rabbit mAb [JB70-35]	ET7107-47	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
KO	ATM Recombinant Rabbit mAb [SI70-01] - BSA and Azide free	HA750096	Human, Rat	WB,IF-Cell,IF-Tissue,IHC-P
KO	ATM Recombinant Rabbit mAb [SI70-01]	ET1606-20	Human, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	Bad Recombinant Rabbit mAb [SA32-01]	ET1601-7	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	Bad Rabbit pAb	R1511-32	Human	WB,IF-Cell,IHC-P,FC
	Bad Recombinant Rabbit mAb [SA32-01] - BSA and Azide free	HA750025	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	Bad Mouse mAb [2H4]	RT1067	Human, Mouse, Rat	WB,IP,IF,IHC-P
	BAFF Mouse mAb [A2-C3]	EM1707-99	Human	WB,IF-Cell,FC
	Bcl10 Recombinant Rabbit mAb [PSH13-49]	HA723538	Human, Mouse, Rat	WB,IF-Cell,FC,IP
	Bcl10 Recombinant Rabbit mAb [SN74-04]	ET1611-79	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Bcl10 Recombinant Rabbit mAb [PSH13-49] - BSA and Azide free	HA751468	Human, Mouse, Rat	WB,IF-Cell,FC,IP
	Bcl-2 Mouse mAb [9F3]	EM1701-83	Human	WB,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	Bcl-2 Recombinant Rabbit mAb [JF104-8]	ET1702-53	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	Bcl-2 Recombinant Rabbit mAb [SZ10-03] - BSA and Azide free	HA750057	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Bcl-2 Rabbit pAb	ER0602	Human	WB,IF-Cell,IHC-P,FC
	Bcl-2 Recombinant Rabbit mAb [PSH09-49]	HA723111	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
	Bcl-2 Recombinant Mouse mAb [9F3-R]	HA601256	Human	WB,IHC-P,FC
	Bcl-2 Recombinant Rabbit mAb [PSH09-49] - BSA and Azide free	HA751296	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
KD	Bcl-2 Recombinant Rabbit mAb [JF104-8] - BSA and Azide free	HA750350	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	FITC Conjugated Bcl-2 Recombinant Rabbit mAb [JF104-8]	HA720210F	Human, Mouse	FC
	Bcl-2 Recombinant Mouse mAb [9F3-R] - BSA and Azide free	HA610131	Human	WB,IHC-P,FC
	Bcl-2 Mouse mAb [9F1]	EM1701-82	Human	WB,IHC-P
KD	Bcl-2 Recombinant Rabbit mAb [SZ10-03]	ET1603-11	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Bcl-2 Mouse mAb [F7-1-A1]	M1206-4	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	Bcl-2 Recombinant Rabbit mAb [PD01-41]	HA721235	Human	WB,IHC-P,IF-Cell,FC
	Bcl-2 Rabbit pAb	ER1706-47	Human	WB,IF-Cell,IHC-P,FC
	Bcl-2 Rabbit pAb	ER1802-97	Human, Mouse	WB,IF-Cell,FC
	BCL2A1 Recombinant Rabbit mAb [SC05-42]	ET1610-20	Human, Mouse	WB,IHC-P
	BCL2A1 Recombinant Rabbit mAb [SC05-42] - BSA and Azide free	HA750205	Human, Mouse	WB,IHC-P
	BCL2L1 Rabbit pAb	R1308-16	Human, Mouse, Rat	WB,IHC-P
	Bcl-x Recombinant Rabbit mAb	RE6009	Human, Mouse, Rat	WB,ICC,IHC-P,FC,IP
	Bcl-XL Recombinant Rabbit mAb [SZ3-03]	ET1603-28	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	Bcl-XL Recombinant Rabbit mAb [SZ3-03] - BSA and Azide free	HA750063	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	BIRC2 Rabbit pAb	ER62503	Human, Mouse	WB
	BLNK Recombinant Rabbit mAb [JE63-50]	HA721028	Human, Mouse, Rat	WB
	BLNK Mouse mAb [G1-F9]	EM1709-94	Human, Mouse	WB,IF-Cell,IHC-P,FC
	BMP-9 Rabbit pAb	ER65787	Human, Mouse	IHC-P
KO, Relative expression	BTK Recombinant Rabbit mAb [SN69-07]	ET1611-4	Human, Mouse	WB,IHC-P,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KO, Relative expression	BTK Recombinant Rabbit mAb [SN69-07] - BSA and Azide free	HA750252	Human, Mouse	WB,IHC-P,IP
	CARD 10 Rabbit pAb	ER62894	Human, Mouse	WB,IHC-P,IF-Cell
KD	CARD11 Recombinant Rabbit mAb [JE37-28]	HA722309	Human, Mouse, Rat	WB,IHC-P
	CARD14 Rabbit pAb	ER1905-39	Human, Mouse	WB,IF-Cell
	Caspase-1 Rabbit pAb	R1510-22	Human	IHC-P,FC
	Caspase-1 Rabbit pAb	R1510-23	Human	IF-Cell,IHC-P,FC
	Caspase-1 Recombinant Rabbit mAb [SU40-07]	ET1608-69	Human, Mouse, Rat	WB,FC,IF-Cell
	Caspase-1 Recombinant Rabbit mAb [SU40-07] - BSA and Azide free	HA750161	Human, Mouse, Rat	WB,FC,IF-Cell
Relative expression	Caspase-1 + p10 + p12 Recombinant Rabbit mAb [JE56-35]	HA722222	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Caspase-1 p20 Rabbit pAb	ER1905-47	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
Cell treatment, KD	Caspase-3 p12 Recombinant Rabbit mAb [JE07-51]	HA722235	Human, Mouse, Rat	WB,IHC-P
	Caspase-4 Mouse mAb [A9G2]	HA601431	Human, Cynomolgus monkey	WB
	Caspase-4 Rabbit pAb	ER62908	Human	WB,IHC-P
	Caspase-5 Recombinant Rabbit mAb [SD203-2]	ET1612-29	Human	WB,IF-Tissue,IHC-P,IP
Cell treatment	Caspase-7 Recombinant Rabbit mAb [PSH13-69] - BSA and Azide free	HA751482	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Caspase-7 Recombinant Rabbit mAb [PSH13-69]	HA723558	Human, Mouse, Rat	WB,IHC-P
	Caspase-9 Recombinant Rabbit mAb [SZ29-01]	ET1603-27	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Caspase-9 Recombinant Rabbit mAb [SZ29-01] - BSA and Azide free	HA750062	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	CCL13 Rabbit pAb	ER61222	Human	WB
	CCL2 Mouse mAb [E11-D8]	EM1710-22	Human, Mouse, Monkey	WB,IF-Cell,IHC-P,FC
Cell treatment	CCL2 / MCP-1 Recombinant Rabbit mAb [PSH18-16]	HA723979	Mouse, Rat	WB,IF-Cell,FC
	CCL2 / MCP-1 Recombinant Rabbit mAb [PSH12-44]	HA723691	Mouse	WB,IF-Cell,FC
	CCL2 / MCP-1 Rabbit pAb	HA500042	Human	WB,IHC-P,IF-Cell
	CCL2 / MCP-1 Rabbit pAb	HA500267	Mouse	WB,IF-Cell,FC
	CCL2 / MCP-1 Recombinant Rabbit mAb [PSH12-44] - BSA and Azide free	HA751533	Mouse	WB,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	CCL4 / MIP-1 beta Recombinant Rabbit mAb [JE33-21]	HA721904	Human, Mouse	WB
	PE Conjugated CD123 / IL3RA Recombinant Mouse mAb [PSH14-51]	HA601528F	Human	FC
	CD123 / IL3RA Recombinant Mouse mAb [PSH14-51]	HA601490	Human	FC
KD	CD14 Recombinant Rabbit mAb [SC69-02] - BSA and Azide free	HA750232	Human, Mouse	WB,IF-Cell,IHC-P,FC
	CD14 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS106RT	Human	mIHC
KD	CD14 Recombinant Rabbit mAb [SC69-02]	ET1610-85	Human, Mouse	WB,IF-Cell,IHC-P,mIHC,FC
	CD14 Recombinant Rabbit mAb	IRS015	Human	mIHC
KD	CD14 Recombinant Antibody [SC69-02] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610260	0	
KD	CD14 Recombinant Rabbit mAb [PD01-63]	HA721321	Human, Mouse	WB,IHC-P,IF-Tissue
KD	CD14 Recombinant Antibody [SC69-02] - Rat IgG1 (Chimeric)	HA601420		
	CD14 Rabbit pAb	R1107-2	Human, Mouse	WB,IF-Cell,FC
	CD354 Mouse mAb [6B12D1]	EM1711-44	Human, Rat, Monkey	WB,FC
KD	CD4 Rabbit pAb	0108-6	Human, Mouse, Rat	WB,IF-Cell,FC
	CD4 Rabbit pAb	ER1803-84	Human	WB,FC
Relative expression	CD4 Recombinant Rabbit mAb [PSH08-22]	HA722966	Mouse	WB,IF-Cell,IHC-P,mIHC,IHC-Fr
Relative expression	CD4 Recombinant Rabbit mAb [PSH08-22] - BSA and Azide free	HA751213	Mouse	WB,IF-Cell,IHC-P,IHC-Fr
	CD4 Rabbit pAb	R1211-1	Human	WB,IF-Cell
	CD4 Recombinant Rabbit mAb	IRS005	Human	mIHC
	CD4 Recombinant Rabbit mAb	IRS051	Mouse	mIHC
	CD4 Recombinant Antibody [PSH08-22] - Rat IgG1 (Chimeric)	HA601534	Mouse	WB,IF-Tissue,IHC-P,IHC-Fr
	CD4 Recombinant Antibody [PSH08-22] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610354	Mouse	WB,IF-Tissue,IHC-P,IHC-Fr
Relative expression	CD4 Recombinant Antibody [PSH08-22] - Mouse IgG1 (Chimeric)	HA601533	Mouse	WB,IF-Tissue,IHC-P,IF-Cell
	CD4 Recombinant Rabbit mAb [JE56-36]	ET7110-64	Human	WB,IHC-P
	CD4 Recombinant Antibody [ST0488] - Mouse IgG1 (Chimeric)	HA601557	Human	WB,IF-Cell,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	CD4 Recombinant Antibody [PSH08-22] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610353	Mouse	WB,IF-Tissue,IHC-P,IF-Cell
	CD4 Recombinant Rabbit Multiclonal Antibody [PSH13-24] - BSA and Azide free	HA751460	Human, Mouse	WB,IF-Cell,IHC-P
	CD4 Recombinant Rabbit Multiclonal Antibody [PSH13-24]	HA723511	Human, Mouse	WB,IF-Cell,IHC-P
KD	CD4 Recombinant Rabbit mAb [ST0488]	ET1609-52	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC,mIHC
	CD40 Recombinant Rabbit mAb [PSH10-01]	HA723170	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-P,IHC-Fr,IF-Cell,FC,IF-Tissue
	CD40 Mouse mAb [A1E5]	EM1901-41	Human	WB,IHC-P
Relative expression	CD40 Mouse mAb [A1E6]	EM1901-40	Human	WB,IF-Cell,IHC-P
	CD40 Rabbit pAb	ER1803-54	Human, Rat, Mouse	WB,IF-Cell,IHC-P,FC
	CD40 Recombinant Rabbit mAb [PSH10-01] - BSA and Azide free	HA751331	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-P,IHC-Fr,IF-Cell,FC,IF-Tissue
Relative expression	CD40 Recombinant Rabbit mAb [PSH05-90]	HA722488	Human, Cynomolgus monkey, Pig	WB,IF-Cell,FC
Relative expression	CD40 Recombinant Rabbit mAb [PSH05-90] - BSA and Azide free	HA751037	Human, Cynomolgus monkey, Pig	WB,IF-Cell,FC
	CD40L Recombinant Rabbit mAb [JM11-34]	ET1703-64	Human, Mouse, Rat	WB
	CD40L Recombinant Rabbit mAb [JM11-34] - BSA and Azide free	HA750387	Human, Mouse, Rat	WB
Relative expression	CD74 Recombinant Mouse mAb [PDM0-22]	HA601117	Human	WB,IHC-P,IF-Cell,IF-Tissue
	CD74 Recombinant Rabbit mAb [JF40-10]	ET1702-51	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
Relative expression, KD	CD9 Recombinant Rabbit mAb [PSH0-95] - BSA and Azide free	HA750675	Human	WB,IHC-P
Relative expression, KD	CD9 Recombinant Rabbit mAb [PSH0-95]	HA721533	Human	WB,IHC-P
	CD9 Rabbit pAb	ER80402	Human, Mouse	WB
	CD9 Recombinant Rabbit mAb [SA35-08]	ET1601-9	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC
	CD9 Recombinant Rabbit mAb [SA35-08] - BSA and Azide free	HA750027	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC
	CDKN1A/P21 Rabbit pAb	ER1906-07	Human, Rat	IHC-P,IF-Cell,IF-Tissue,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	CDKN1A/P21 Rabbit pAb	HA500156	Human, Mouse	WB,IHC-P
	CFLAR Mouse mAb [D12-A3]	EM1708-94	Human, Mouse	WB,IHC-P,FC
KD	cIAP1 Recombinant Rabbit mAb [PSH06-71]	HA722743	Human, Mouse, Rat	WB
KD	cIAP1 Recombinant Rabbit mAb [PSH06-71] - BSA and Azide free	HA751100	Human, Mouse, Rat	WB
	cIAP1 Recombinant Rabbit mAb [PSH06-72]	HA722744	Human, Mouse, Rat	WB,IF-Cell,FC
	cIAP1 Recombinant Rabbit mAb [PSH06-71] - BSA and Azide free (Detector)	HA723282	Human, Mouse, Rat	ELISA(Det)
	cIAP1 Recombinant Rabbit mAb [PSH06-72] - BSA and Azide free	HA751101	Human, Mouse, Rat	WB,IF-Cell,FC
KD	cIAP1 Recombinant Rabbit mAb [PSH02-10]	HA721769	Human	WB
	cIAP1 Recombinant Rabbit mAb [PSH06-72] - BSA and Azide free (Capture)	HA723280	Human, Mouse, Rat	ELISA(Cap)
	cIAP1 Recombinant Rabbit mAb [PSH02-10] - BSA and Azide free (Capture)	HA723281	Human	ELISA(Cap)
KD	cIAP1 Recombinant Rabbit mAb [PSH02-10] - BSA and Azide free	HA750780	Human	WB
	Biotin Conjugated cIAP1 Recombinant Rabbit mAb [PSH06-71] - Detector	HA723283B	Human, Mouse, Rat	ELISA(Det),ELISA
	cIAP1/2 Rabbit pAb	HA500535	Human, Mouse	WB,IF-Cell,IHC-P,FC
	cIAP2 Recombinant Rabbit mAb [PSH06-92]	HA722767	Human	WB
	cIAP2 Recombinant Rabbit mAb [PSH06-92] - BSA and Azide free	HA751116	Human	WB
	CK II alpha Rabbit pAb	ER1906-34	Mouse, Rat	WB
	Cleaved-Bad (D71) Rabbit pAb	ER60015	Mouse, Rat	WB,IHC-P
Cell treatment	Cleaved Caspase-3 Recombinant Rabbit mAb [SR01-02] - BSA and Azide free	HA750053	Human, Pig	WB,IF-Cell,IHC-P
	Cleaved Caspase-3 Recombinant Rabbit mAb	IRS136RB	Human	mIHC
Cell treatment	Cleaved Caspase-3 Recombinant Rabbit mAb [SR01-02]	ET1602-47	Human, Pig	WB,IF-Cell,IHC-P
	Cleaved Caspase-3 Rabbit pAb	IRS227RB	Mouse	mIHC
	Cleaved-Caspase-4/5 p20 (D270/D311) Rabbit pAb	ER60012	Human	WB
	Cleaved-Caspase-4 p20 (Q81) Rabbit pAb	ER60011	Human, Mouse	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Cleaved Caspase-7 (Ser199) Rabbit pAb	ER60002	Human, Mouse	WB,IHC-P,IF-Cell
	Cleaved-Caspase-9 (D353) Rabbit pAb	ER60008	Mouse, Rat, Human	IF-Cell,WB,IHC-P
Cell treatment, KD	Cleaved PARP Recombinant Rabbit mAb [SU0314] - BSA and Azide free	HA750136	Human	WB,IF-Cell,IP
Cell treatment, KD	Cleaved PARP Recombinant Rabbit mAb [SU0314]	ET1608-10	Human	WB,IF-Cell,IP
Cell treatment	Cleaved PARP1 Recombinant Rabbit mAb [JE00-15]	HA722218	Human, Mouse, Rat, Monkey	WB
Cell treatment, KD	Cleaved+pro Caspase-3 Recombinant Rabbit mAb [JE75-05]	HA722367	Human, Mouse	WB
Cell treatment, KO	Cleaved+pro Caspase-3 Recombinant Rabbit mAb [SU38-04] - BSA and Azide free	HA750159	Human, Mouse, Rat	WB,IHC-P,IP
Cell treatment, KO	Cleaved+pro Caspase-3 Recombinant Rabbit mAb [SU38-04]	ET1608-64	Human, Mouse, Rat	WB,IHC-P,IP
	COX2 Mouse mAb	RT1159	Human, Mouse, Rat	WB,IP,IF,IHC-P,ELISA
	COX2 Mouse mAb [A3F8]	EM1902-11	Human, Rat	WB,IHC-P,FC
	COX2 Recombinant Rabbit mAb [6J8]	RE6064	Human, Mouse	WB,ICC
KD	COX2 Mouse mAb [A3F7]	EM1902-12	Human, Rat	WB,IHC-P,FC
KD	COX2/Cyclooxygenase 2 Recombinant Rabbit mAb [SC56-06]	ET1610-23	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
KD	COX2/Cyclooxygenase 2 Recombinant Rabbit mAb [SC56-06] - BSA and Azide free	HA750207	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	CSNK2A1 Rabbit pAb	HA500528	Human, Mouse	WB,IHC-P
	CSNK2A1 Recombinant Rabbit mAb [PSH04-54]	HA722153	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	CSNK2A1 Recombinant Rabbit mAb [PSH04-54] - BSA and Azide free	HA750961	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	CSNK2A2 Rabbit pAb	HA500514	Human, Mouse, Rat	WB
	CSNK2B Recombinant Rabbit mAb [JU37-62]	ET7106-67	Human, Mouse, Rat	WB,IHC-P
	CXCL12/SDF1 Recombinant Rabbit mAb [JJ09-20]	ET1701-77	Human	WB,IF-Cell,IF-Tissue,IHC-P
	CXCL12/SDF1 Rabbit pAb	ER1902-35	Human, Mouse, Rat	WB,IHC-P
	CXCL12/SDF1 Rabbit pAb	ER1802-6	Human	WB,IHC-P,FC
Cell treatment	CXCL1 / GRO alpha Recombinant Rabbit mAb [PSH07-53]	HA722841	Human	WB
Cell treatment	CXCL1 / GRO alpha Recombinant Rabbit mAb [PSH07-53] - BSA and Azide free	HA751174	Human	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CYLD Rabbit pAb	ER1906-97	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
KO	CYLD Recombinant Rabbit mAb [PSH11-21]	HA723303	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell
KO	CYLD Recombinant Rabbit mAb [PSH11-22]	HA723304	Human, Mouse, Monkey	WB,IF-Cell
Relative expression	Cystatin F Recombinant Rabbit mAb [PSH09-00] - BSA and Azide free	HA751261	Human, Mouse	WB,IHC-Fr
Relative expression	Cystatin F Recombinant Rabbit mAb [PSH09-00]	HA723053	Human, Mouse	WB,IHC-Fr
	DDX58 Recombinant Rabbit mAb [JE60-05]	HA720090	Human	WB
	DDX58 Rabbit pAb	HA500374	Human, Mouse	WB,IHC-P,IF-Cell,FC
	EDA Rabbit pAb	ER1907-93	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue
	ELKS Recombinant Rabbit mAb [PSH17-33]	HA723916	Human, Mouse, Rat	WB
	ELKS Recombinant Rabbit mAb [PSH17-33] - BSA and Azide free	HA751621	Human, Mouse, Rat	WB
	ELKS Rabbit pAb	HA500376	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	EpCAM Recombinant Rabbit mAb	IRS254RB	Human	mIHC
	EpCAM Recombinant Rabbit mAb [SU03-32]	ET1608-26	Human, Mouse, Rat	WB,IHC-P
	EpCAM Recombinant Rabbit mAb [PS01-69] - BSA and Azide free	HA750703	Human	IHC-P,IF-Cell,FC
	iFluor™ 555 Conjugated EpCAM Recombinant Rabbit mAb [PS01-69]	HA720218F	Human	FC
	iFluor™ 488 Conjugated EpCAM Recombinant Rabbit mAb [PS01-69]	HA720217F	Human	IF-Cell
	EpCAM Rabbit pAb	ER1802-80	Human, Mouse, Rat	WB,IHC-P
	EpCAM Recombinant Rabbit mAb [PS01-69]	HA721644	Human	IHC-P,IF-Cell,FC
	EpCAM Mouse mAb [3-E1-E2]	EM1111	Human, Mouse	WB,IF-Cell
	iFluor™ 647 Conjugated EpCAM Recombinant Rabbit mAb [PS01-69]	HA720185F	Human	FC,IF-Cell
	EpCAM Mouse mAb [7E11]	EM1712-43	Human	WB,IHC-P,ELISA
KD	ERK1 Recombinant Rabbit mAb [SP00-74]	ET1604-32	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	ERK1 Recombinant Rabbit mAb [SP05-09]	ET1604-16	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	ERK1 Mouse mAb [7-D6-E5]	M1505-9	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Erk1/2 Mouse mAb [2G1]	RT1484	Human, Mouse, Rat	WB,IP,IF

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	ERK1/2 Recombinant Rabbit mAb [SA43-03]	ET1601-29	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IP,FC,IHC-P,IHC-Fr
KD	ERK1/2 Recombinant Rabbit mAb [SA43-03] - BSA and Azide free	HA750019	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IP,FC,IHC-P,IHC-Fr
	ERK1/2 Rabbit pAb	ER131011	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	ERK2 Recombinant Mouse mAb [C3B12-R] - BSA and Azide free	HA610079	Human, Mouse, Rat	WB
	ERK2 Mouse mAb [C3B12]	EM1901-55	Human, Mouse, Rat	WB,IHC-P,FC
	ERK2 Mouse mAb [5-D2]	M1407-3	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	ERK2 Recombinant Mouse mAb [5-D2-R]	HA601202	Human, Mouse, Rat	WB
KD	ERK2 Mouse mAb [B10-8]	EM1901-53	Human, Mouse, Rat	WB,IHC-P,FC
KD	ERK2 Mouse mAb [A10D2]	EM1901-54	Human, Mouse, Rat	WB,IHC-P,FC
KD	ERK2 Rabbit pAb	ER131218	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	ERK2 Recombinant Rabbit mAb [SZ25-01]	ET1603-23	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
KD	ERK2 Recombinant Mouse mAb [5-D2-R] - BSA and Azide free	HA610080	Human, Mouse, Rat	WB
KD	ERK2 Recombinant Mouse mAb [C3B12-R]	HA601201	Human, Mouse, Rat	WB
	ERK 3 Rabbit pAb	ER63397	Human, Mouse, Rat	WB,IHC-P,IF-Cell
KD	ERK5 Recombinant Rabbit mAb [SD2084] - BSA and Azide free	HA750299	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IP,FC
KD	ERK5 Recombinant Rabbit mAb [SD2084]	ET1612-7	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IP,FC
	F13A1 Mouse mAb [16D1]	EM1901-39	Human	WB,IF-Cell,IHC-P,FC
	F13A1 Mouse mAb [16D3]	EM1901-38	Human	WB,IHC-P,FC
	Factor XIIIa Rabbit pAb	ER1902-34	Human	WB,IHC-P,FC
	Factor XIIIa Recombinant Rabbit mAb [PD00-08]	HA721148	Human	WB,IHC-P,IF-Tissue
	FLIP Recombinant Rabbit mAb [PSH01-95]	HA721753	Human	WB
	FLIP Recombinant Rabbit mAb [PSH01-95] - BSA and Azide free	HA750769	Human	WB
	Flk-1/Flt-4 Rabbit pAb	ER63443	Human, Mouse, Rat	IHC-P,IF-Cell
	Flotillin 1 Recombinant Rabbit mAb [JB19-45]	ET7107-82	Human, Mouse, Rat	WB,IHC-P
	Flt-4 Rabbit pAb	ER65750	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	GADD 45 gamma Rabbit pAb	ER63492	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Cell treatment	Gasdermin D (N terminal) Recombinant Rabbit mAb [PD00-18] - BSA and Azide free	HA750486	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Gasdermin D (N terminal) Rabbit pAb	ER1901-37	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell
Cell treatment	Gasdermin D (N terminal) Recombinant Rabbit mAb [PSH10-87]	HA723254	Human	WB,IP
Cell treatment	Gasdermin D (N terminal) Recombinant Rabbit mAb [PD00-18]	HA721144	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Cell treatment	Gasdermin D (N terminal) Recombinant Rabbit mAb [PSH10-87] - BSA and Azide free	HA751379	Human	WB,IP
	GM130 Recombinant Rabbit mAb [JE42-53]	HA721282	Human, Mouse	WB,IHC-P,IF-Cell
	GM130 (cis-Golgi Marker) Rabbit pAb	R1608-7	Human, Mouse, Rat	WB,IHC-P
	GRO alpha Rabbit pAb	ER63645	Human, Mouse, Rat	IHC-P
	GRO gamma Rabbit pAb	ER65798	Human	IHC-P
	Group VI iPLA2 Rabbit pAb	ER63644	Human, Mouse, Rat	WB
	GSDMD Mouse mAb [A8D8]	HA601047	Human	WB
	GSDMD Mouse mAb [A8D7]	HA601046	Human, Mouse	WB
	HLA-DR Recombinant Mouse mAb	IRS154MS	Human	mIHC
	HLA-DR Recombinant Mouse mAb [10-D8-R]	HA601252	Human	WB,IF-Cell,IHC-P,IF-Tissue
	HLA-DR Rabbit pAb	R1510-34	Human, Mouse	WB,IF-Cell,IHC-P,FC
	HLA-DR Recombinant Rabbit mAb [SC06-78]	ET1610-66	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	HLA-DR Recombinant Rabbit mAb [SC06-78] - BSA and Azide free	HA750230	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	HLA-DR Mouse mAb [10-D8]	M1701-5	Human	WB,IHC-P,IF-Cell
	HLA-DR Rabbit pAb	R1510-35	Human	WB,IF-Cell,IHC-P,FC
	HLA-DR Recombinant Mouse mAb [10-D8-R] - BSA and Azide free	HA610127	Human	WB,IF-Cell,IHC-P,IF-Tissue
KO	HMGB1 Recombinant Rabbit mAb [SA39-03]	ET1601-2	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IHC-Fr
KO	HMGB1 Rabbit pAb	R1309-12	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KO	HMGB1 Recombinant Rabbit mAb [SA39-03] - BSA and Azide free	HA750011	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IHC-Fr
KD	HMGB1 Rabbit pAb	ER0913	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Hsp70 Recombinant Rabbit mAb [PD01-65]	HA721497	Human	WB,IHC-P,IF-Tissue
	Hsp70 Rabbit pAb	ER50802	Human, Mouse, Rat, Cow	WB,IF-Cell,IHC-P,FC
KD	Hsp70 Recombinant Rabbit mAb [SA0379]	ET1601-11	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	HSP70 Mouse mAb [3G10]	EM60062	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	HSP70 (Acetyl K246) Rabbit pAb	ER60056	Human	WB
	Human Caspase-9 Recombinant Rabbit mAb [PSH17-19] - BSA and Azide free (Capture)	HA723906	Human	ELISA(Cap)
	Human Caspase-9 Recombinant Rabbit mAb [PSH17-20] - BSA and Azide free (Detector)	HA723907	Human	ELISA(Det)
	Biotin Conjugated Human Caspase-9 Recombinant Rabbit mAb [PSH17-20] - Detector	HA723907B	Human	ELISA(Det),ELISA
	Human CCL21 Recombinant Rabbit mAb [PSH13-37] - BSA and Azide free (Detector)	HA723519	Human	ELISA(Det)
	Human CCL21 Recombinant Rabbit mAb [PSH13-36] - BSA and Azide free (Capture)	HA723518	Human	ELISA(Cap)
	Biotin Conjugated Human CCL21 Recombinant Rabbit mAb [PSH13-37] - Detector	HA723520B	Human	ELISA(Det),ELISA
	Human CCL2 / MCP-1 Recombinant Rabbit mAb [PSH18-21] - BSA and Azide free (Capture)	HA725195	Human	ELISA(Cap)
	Human CCL2 / MCP-1 Recombinant Rabbit mAb [PSH18-22] - BSA and Azide free (Detector)	HA725196	Human	ELISA(Det)
	Biotin Conjugated Human CCL2 / MCP-1 Recombinant Rabbit mAb [PSH18-22] - Detector	HA725196B	Human	ELISA(Det)
	Biotin Conjugated Human CCL4 / MIP-1 beta Recombinant Rabbit mAb [PSH19-49] - Detector	HA725356B	Human	ELISA(Det),ELISA
	Human CCL4 / MIP-1 beta Recombinant Rabbit mAb [PSH19-48] - BSA and Azide free (Capture)	HA725355	Human	ELISA(Cap)
	Human CCL4 / MIP-1 beta Recombinant Rabbit mAb [PSH19-49] - BSA and Azide free (Detector)	HA725356	Human	ELISA(Det)
	Biotin Conjugated Human CCL5 Recombinant Rabbit mAb [PSH13-17] - Detector	HA725057B	Human	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human CCL5 Recombinant Rabbit mAb [PSH13-16] - BSA and Azide free (Capture)	HA725055	Human	ELISA(Cap)
	HRP Conjugated Human CCL5 Recombinant Rabbit mAb [PSH13-17] - Detector	HA725057H	Human	ELISA(Det),ELISA
	Human CCL5 Recombinant Rabbit mAb [PSH13-17] - BSA and Azide free (Detector)	HA725056	Human	ELISA(Det)
	iFluor™ 647 Conjugated Human CD14 Recombinant Rabbit mAb [SC69-02]	HA720184F	Human	FC
	FITC Conjugated Human CD14 Recombinant Rabbit mAb [SC69-02]	HA720194F	Human	FC
	PE-Cy7 Conjugated Human CD4 Recombinant Mouse mAb [PSH04-93]	HA601542F	Human	FC
	APC Conjugated Human CD4 Recombinant Mouse mAb [PSH04-93]	HA600121F	Human	FC
	Human CD4 Recombinant Mouse mAb [PSH04-93]	HA601338	Human	FC
	PE Conjugated Human CD4 Recombinant Mouse mAb [PSH04-93]	HA600120F	Human	FC
	FITC Conjugated Human CD4 Recombinant Mouse mAb [PSH04-93]	HA600119F	Human	FC
	Biotin Conjugated Human cIAP2 Recombinant Rabbit mAb [PSH09-39] - Detector	HA723101B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human cIAP2 Recombinant Rabbit mAb [PSH09-40] - Detector	HA723102B	Human	ELISA(Det),ELISA
	Human cIAP2 Recombinant Rabbit mAb [PSH09-38] - BSA and Azide free (Capture)	HA723098	Human	ELISA(Cap)
	Human cIAP2 Recombinant Rabbit mAb [PSH09-40] - BSA and Azide free (Detector)	HA723100	Human	ELISA(Det)
	Human cIAP2 Recombinant Rabbit mAb [PSH09-39] - BSA and Azide free (Capture/Detector)	HA723099	Human	ELISA(Det),ELISA(Cap)
	Human CXCL10 / IP-10 Recombinant Rabbit mAb [PSH12-37] - BSA and Azide free (Detector)	HA725008	Human	ELISA(Det)
	Biotin Conjugated Human CXCL10 / IP-10 Recombinant Rabbit mAb [PSH12-37] - Detector	HA725009B	Human	ELISA(Det),ELISA
	Human CXCL10 / IP-10 Recombinant Rabbit mAb [PSH12-36] - BSA and Azide free (Capture)	HA725007	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human CXCL1/GRO alpha Recombinant Rabbit mAb [PSH05-94] - BSA and Azide free (Capture)	HA722493	Human	ELISA(Cap)
	Biotin Conjugated Human CXCL1/GRO alpha Recombinant Rabbit mAb [PSH05-95]	HA722587B	Human	ELISA(Det)
	Human CXCL1/GRO alpha Recombinant Rabbit mAb [PSH05-95] - BSA and Azide free (Detector)	HA722494	Human	ELISA(Det)
	Human ICAM1 Recombinant Rabbit mAb [PSH10-15] - BSA and Azide free (Capture)	HA723184	Human	ELISA(Cap)
	Human ICAM1 Recombinant Rabbit mAb [PSH10-16] - BSA and Azide free (Detector)	HA723185	Human	ELISA(Det)
	HRP Conjugated Human ICAM1 Recombinant Rabbit mAb [PSH10-16] - Detector	HA723474H	Human	ELISA(Det),ELISA
	Biotin Conjugated Human ICAM1 Recombinant Rabbit mAb [PSH10-16] - Detector	HA723186B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human IFN beta Recombinant Rabbit mAb [PSH17-30] - Detector	HA725054B	Human	ELISA(Det),ELISA
	Human IFN beta Recombinant Rabbit mAb [PSH17-29] - BSA and Azide free (Capture)	HA725052	Human	ELISA(Cap)
	Human IFN beta Recombinant Rabbit mAb [PSH17-30] - BSA and Azide free (Detector)	HA725053	Human	ELISA(Det)
	Human IL-10 Recombinant Rabbit mAb [PS00-98] - BSA and Azide free (Capture)	HA721236	Human	ELISA(Cap)
	Biotin Conjugated Human IL-10 Recombinant Rabbit mAb [PS00-81]	HA722652B	Human	ELISA(Det)
	Human IL-10 Recombinant Rabbit mAb [PS00-81] - BSA and Azide free (Detector)	HA721237	Human	ELISA(Det)
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-87] - BSA and Azide free (Capture)	HA723883	Human	ELISA(Cap)
	Human IL-1 beta Recombinant Rabbit mAb [PS01-39] - BSA and Azide free (Capture)	HA721278	Human	ELISA(Cap)
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS01-39] - BSA and Azide free	HA750744	Human	WB
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS00-82]	HA721713	Human	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human IL-1 beta Recombinant Rabbit mAb [PSH16-86] - BSA and Azide free (Detector)	HA723882	Human	ELISA(Det)
	Biotin Conjugated Human IL-1 beta Recombinant Rabbit mAb [PSH16-86] - Detector	HA723882B	Human	ELISA(Det),ELISA
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS01-39]	HA721712	Human	WB
	Biotin Conjugated Human IL-1 beta Recombinant Rabbit mAb [PS00-82]	HA722649B	Human	ELISA(Det)
	Human IL-1 beta Recombinant Rabbit mAb [PS00-82] - BSA and Azide free (Detector)	HA721279	Human	ELISA(Det)
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-85] - BSA and Azide free (Capture)	HA723881	Human	ELISA(Cap)
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS00-82] - BSA and Azide free	HA750745	Human	WB
	Biotin Conjugated Human IL-2 Mouse mAb [A8B7]	HA601361B	Human	ELISA(Det)
	Human IL-2 Mouse mAb [A8B7] - BSA and Azide free (Detector)	HA601006	Human	WB,ELISA(Det)
	Human IL-2 Recombinant Rabbit mAb [PSH01-42] - BSA and Azide free (Capture)	HA721680	Human	ELISA(Cap)
	Biotin Conjugated Human IL3RA/CD123 Recombinant Rabbit mAb [PSH11-58] - Detector	HA723353B	Human	ELISA(Det),ELISA
	Human IL3RA/CD123 Recombinant Rabbit mAb [PSH11-57] - BSA and Azide free (Capture)	HA723351	Human	ELISA(Cap)
	Human IL3RA/CD123 Recombinant Rabbit mAb [PSH11-58] - BSA and Azide free (Detector)	HA723352	Human	ELISA(Det)
	Human IL-4 Recombinant Rabbit mAb [PS00-66] - BSA and Azide free (Detector)	HA721231	Human	ELISA(Det)
	Biotin Conjugated Human IL-4 Recombinant Rabbit mAb [PS00-66] - Detector	HA721368	Human	ELISA(Det)
	HRP Conjugated Human IL-4 Recombinant Rabbit mAb [PS00-66] - Detector	HA721295	Human	ELISA(Det)
	Human IL-4 Recombinant Rabbit mAb [PS00-67] - BSA and Azide free (Capture)	HA721230	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human IL-5 Recombinant Rabbit mAb [PSH04-65] - Detector	HA722609B	Human	ELISA(Det)
	Human IL-5 Recombinant Rabbit mAb [PSH04-65] - BSA and Azide free (Detector)	HA722165	Human	ELISA(Det)
	Human IL-5 Recombinant Rabbit mAb [PSH04-64] - BSA and Azide free (Capture)	HA722164	Human	ELISA(Cap)
	Human IL-6 Mouse mAb [A6D1]	HA601051	Human	WB,ELISA
	HRP Conjugated Human IL-6 Recombinant Rabbit mAb [PS00-15] - Detector	HA721293	Human	ELISA(Det)
	Human IL-6 Recombinant Rabbit mAb [PS00-21] - BSA and Azide free (Capture)	HA721167	Human	ELISA(Cap)
	Human IL-6 Mouse mAb [A6D2]	HA601050	Human	WB,ELISA
	Human IL-6 Recombinant Rabbit mAb [PS00-15] - BSA and Azide free (Detector)	HA721168	Human	ELISA(Det)
	Biotin Conjugated Human IL-6 Recombinant Rabbit mAb [PS00-15] - Detector	HA722233B	Human	ELISA(Det)
	Human IL-8 Recombinant Rabbit mAb [PS01-22] - BSA and Azide free (Detector)	HA722099	Human	ELISA(Det)
	Biotin Conjugated Human IL-8 Recombinant Rabbit mAb [PS01-22]	HA722611B	Human	ELISA(Det)
	Human IL-8 Recombinant Rabbit mAb [PS00-83] - BSA and Azide free (Capture)	HA722098	Human	ELISA(Cap)
	Biotin Conjugated Human MCP4/CCL13 Recombinant Rabbit mAb [PSH15-70] - Detector	HA725200B	Human	ELISA(Det),ELISA
	Human MCP4/CCL13 Recombinant Rabbit mAb [PSH15-70] - BSA and Azide free (Detector)	HA725200	Human	ELISA(Det)
	Human MCP4/CCL13 Recombinant Rabbit mAb [PSH15-69] - BSA and Azide free (Capture)	HA725199	Human	ELISA(Cap)
	Human MMP-3 Recombinant Rabbit mAb [PSH07-62] - BSA and Azide free (Detector)	HA722878	Human	ELISA(Det)
	Human MMP-3 Recombinant Rabbit mAb [PSH07-61] - BSA and Azide free (Capture)	HA722877	Human	ELISA(Cap)
	Biotin Conjugated Human MMP-3 Recombinant Rabbit mAb [PSH07-62]	HA722879B	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human NLRP3 Recombinant Rabbit mAb [PSH11-08] - Detector	HA723277B	Human	ELISA(Det),ELISA
	Human NLRP3 Recombinant Rabbit mAb [PSH11-07] - BSA and Azide free (Capture)	HA723275	Human	ELISA(Cap)
	Human NLRP3 Recombinant Rabbit mAb [PSH11-08] - BSA and Azide free (Detector)	HA723276	Human	ELISA(Det)
	Biotin Conjugated Human NLRP3 Recombinant Rabbit mAb [PSH11-09] - Detector	HA723279B	Human	ELISA(Det),ELISA
	Human NLRP3 Recombinant Rabbit mAb [PSH11-09] - BSA and Azide free (Detector)	HA723278	Human	ELISA(Det)
	Human PAI1 Recombinant Rabbit mAb [PSH12-19] - BSA and Azide free (Detector)	HA723436	Human	ELISA(Det)
	Biotin Conjugated Human PAI1 Recombinant Rabbit mAb [PSH12-18] - Detector	HA723439B	Human	ELISA(Det),ELISA
	Human PAI1 Recombinant Rabbit mAb [PSH12-18] - BSA and Azide free (Capture/Detector)	HA723435	Human	ELISA(Cap),ELISA(Det)
	Biotin Conjugated Human PAI1 Recombinant Rabbit mAb [PSH12-19] - Detector	HA723437B	Human	ELISA(Det),ELISA
	Human PAI1 Recombinant Rabbit mAb [PSH12-20] - BSA and Azide free (Capture)	HA723438	Human	ELISA(Cap)
	HRP Conjugated Human PAI1 Recombinant Rabbit mAb [PSH12-18] - Detector	HA723439H	Human	ELISA(Det),ELISA
	Human PD-L1 Recombinant Rabbit mAb [PSH15-36] - BSA and Azide free (Capture)	HA723740	Human	ELISA(Cap)
	Human PD-L1 Recombinant Rabbit mAb [PSH15-37] - BSA and Azide free (Detector)	HA723741	Human	ELISA(Det)
	Biotin Conjugated Human PD-L1 Recombinant Rabbit mAb [PSH15-37] - Detector	HA723742B	Human	ELISA(Det),ELISA
	Human RANKL Recombinant Rabbit mAb [PSH16-59] - BSA and Azide free (Detector)	HA723861	Human	ELISA(Det)
	Human RANKL Recombinant Rabbit mAb [PSH16-60] - BSA and Azide free (Capture)	HA723862	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human RANKL Recombinant Rabbit mAb [PSH16-58] - BSA and Azide free (Capture)	HA723860	Human	ELISA(Cap)
	Biotin Conjugated Human RANKL Recombinant Rabbit mAb [PSH16-59] - Detector	HA723861B	Human	ELISA(Det),ELISA
	Human sTNF RI Recombinant Rabbit mAb [PSH09-89] - BSA and Azide free (Capture)	HA723158	Human	ELISA(Cap)
	HRP Conjugated Human sTNF RI Recombinant Rabbit mAb [PSH09-90] - Detector	HA723374H	Human	ELISA(Det)
	Biotin Conjugated Human sTNF RI Recombinant Rabbit mAb [PSH09-90] - BSA and Azide free (Detector)	HA723160B	Human	ELISA(Det),ELISA
	Human sTNF RI Recombinant Rabbit mAb [PSH09-90] - BSA and Azide free (Detector)	HA723159	Human	ELISA(Det)
	Human TNF alpha Recombinant Rabbit mAb [PS00-40]	HA721286	Human	IHC-P,mIHC
	Human TNF alpha Recombinant Rabbit mAb [PS01-07] - BSA and Azide free (Detector)	HA721266	Human	ELISA(Det)
	Human TNF alpha Recombinant Rabbit mAb [PS00-40] - BSA and Azide free	HA750593	Human	IHC-P
	Human TNF alpha Recombinant Rabbit mAb	IRS040	Human	mIHC
	Human TNF alpha Recombinant Rabbit mAb [PS00-40] - BSA and Azide free (Capture)	HA721265	Human	ELISA(Cap)
	HRP Conjugated Human TNF alpha Recombinant Rabbit mAb [PS01-07] - Detector	HA721294	Human	ELISA(Det)
	Biotin Conjugated Human TNF alpha Recombinant Rabbit mAb [PS01-07] - Detector	HA721559	Human	ELISA(Det)
	Human TNF-beta Recombinant Rabbit mAb [PS01-05] - BSA and Azide free (Capture)	HA721242	Human	ELISA(Cap)
	Human TNF-beta Recombinant Rabbit mAb [PS01-19] - BSA and Azide free (Detector)	HA721243	Human	ELISA(Det)
	Biotin Conjugated Human TNF-beta Recombinant Rabbit mAb [PS01-19]	HA722651B	Human	ELISA(Det)
	Human TRAF1 Recombinant Rabbit mAb [PSH14-27] - BSA and Azide free (Capture)	HA723618	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human TRAF1 Recombinant Rabbit mAb [PSH14-28] - BSA and Azide free (Detector)	HA723619	Human	ELISA(Det)
	Biotin Conjugated Human TRAF1 Recombinant Rabbit mAb [PSH14-28] - Detector	HA723620B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human TREM1 Recombinant Rabbit mAb [PSH17-86] - Detector	HA725234B	Human	ELISA(Det),ELISA
	Human TREM1 Recombinant Rabbit mAb [PSH17-86] - BSA and Azide free (Detector)	HA725234	Human	ELISA(Det)
	Human TREM1 Recombinant Rabbit mAb [PSH17-85] - BSA and Azide free (Capture)	HA725233	Human	ELISA(Cap)
	Biotin Conjugated Human VCAM-1 Recombinant Rabbit mAb [PSH12-47] - Detector	HA725066B	Human	ELISA(Det),ELISA
	HRP Conjugated Human VCAM-1 Recombinant Rabbit mAb [PSH12-47] - Detector	HA725066H	Human	ELISA(Det),ELISA
	Human VCAM-1 Recombinant Rabbit mAb [PSH12-47] - BSA and Azide free (Detector)	HA725065	Human	ELISA(Det)
	Human VCAM-1 Recombinant Rabbit mAb [PSH12-46] - BSA and Azide free (Capture)	HA725064	Human	ELISA(Cap)
	Human VEGFA Recombinant Rabbit mAb [PSH12-91] - BSA and Azide free (Capture)	HA723493	Human	ELISA(Cap)
	Human VEGFA Recombinant Rabbit mAb [PSH12-92] - BSA and Azide free (Detector)	HA723494	Human	ELISA(Det)
	Human VEGFA Recombinant Rabbit mAb [PSH12-93] - BSA and Azide free (Detector)	HA723496	Human	ELISA(Det)
	Biotin Conjugated Human VEGFA Recombinant Rabbit mAb [PSH12-93] - Detector	HA723497B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human VEGFA Recombinant Rabbit mAb [PSH12-92] - Detector	HA723495B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human VEGF Receptor 1 Recombinant Rabbit mAb [PSH19-23] - BSA and Azide free (Detector)	HA725366B	Human	ELISA(Det),ELISA
	Human VEGF Receptor 1 Recombinant Rabbit mAb [PSH19-22] - BSA and Azide free (Capture)	HA725365	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human VEGF Receptor 1 Recombinant Rabbit mAb [PSH19-23] - BSA and Azide free (Detector)	HA725366	Human	ELISA(Det)
	Human VEGF Receptor 2 Recombinant Rabbit mAb [PSH12-84] - BSA and Azide free (Detector)	HA725122	Human	ELISA(Det)
	Human VEGF Receptor 2 Recombinant Rabbit mAb [PSH12-83] - BSA and Azide free (Capture)	HA725121	Human	ELISA(Cap)
	Biotin Conjugated Human VEGF Receptor 2 Recombinant Rabbit mAb [PSH12-84] - Detector	HA725123B	Human	ELISA(Det),ELISA
Relative expression	ICAM1 Recombinant Rabbit mAb [ST0487]	ET1609-46	Human	WB,IF-Tissue,IHC-P
Relative expression	ICAM1 Recombinant Rabbit mAb [ST0487] - BSA and Azide free	HA750176	Human	WB,IF-Tissue,IHC-P
KD	ICAM1 Mouse mAb [F11-A4]	M1505-7	Human	WB,IF-Cell,FC
	ICAM1 Mouse mAb [J0-F7]	M1511-6	Human	WB,IF-Cell,FC
	ICAM1 Rabbit pAb	ER131207	Human	WB,IF-Cell,IHC-P,FC
Relative expression, KD	ICAM1 Recombinant Rabbit mAb [PSH09-31]	HA723089	Human	WB,IHC-P,IF-Tissue
	IFN beta Rabbit pAb	ER65794	Human	IHC-P
	IKB alpha Recombinant Rabbit mAb [SZ00-07]	ET1603-6	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	IKB alpha Recombinant Rabbit mAb [SZ00-07] - BSA and Azide free	HA750068	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	IKK alpha Rabbit pAb	ER30911	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	IKK alpha Recombinant Rabbit mAb [SN06-99] - BSA and Azide free	HA750246	Human, Mouse, Rat	WB,IHC-P,IP,FC
	IKK alpha Recombinant Rabbit mAb [SN06-99]	ET1611-15	Human, Mouse, Rat	WB,IHC-P,IP,FC
KD	IKK alpha + IKK beta Recombinant Rabbit mAb [SN63-02] - BSA and Azide free	HA750248	Human, Mouse, Rat	WB,IF-Cell,IP
KD	IKK alpha + IKK beta Recombinant Rabbit mAb [SN63-02]	ET1611-23	Human, Mouse, Rat	WB,IF-Cell,IP
	IKK beta Rabbit pAb	R1706-13	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	IKK beta Recombinant Rabbit mAb [JE32-81]	HA721715	Human	WB,FC
KD	IKK gamma Recombinant Rabbit mAb [JA11-83]	ET1704-40	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,IP,FC
	IL-10 Rabbit pAb	ER1911-19	Human, Mouse	WB,IHC-P,IF-Tissue
	IL-10 Rabbit pAb	ER65464	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	IL-10 Recombinant Rabbit mAb [JE59-71]	HA722032	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	IL-13 Rabbit pAb	ER1911-23	Human	IHC-P,IF-Tissue
	IL-1 alpha Rabbit pAb	ER1802-36	Human, Mouse	IHC-P,FC,WB
Cell treatment	IL-1 alpha Recombinant Rabbit mAb [PSH15-65]	HA723772	Mouse	WB
	IL-1 alpha Recombinant Rabbit mAb [PSH15-66] - BSA and Azide free	HA751577	Mouse	IF-Cell,FC
Cell treatment	IL-1 alpha Recombinant Rabbit mAb [PSH15-65] - BSA and Azide free	HA751576	Mouse	WB
	IL-1 alpha Recombinant Rabbit mAb [PSH15-66]	HA723773	Mouse	IF-Cell,FC
	IL-1 beta Recombinant Rabbit mAb [PSH18-00]	HA723965	Mouse, Rat	WB
Cell treatment	IL-1 beta Recombinant Rabbit mAb [JJ087-3]	ET1701-39	Human	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	IL-1 beta Mouse mAb [A7A10]	HA601035	Human	WB,ELISA,IF-Cell
Cell treatment	IL-1 beta Mouse mAb [A7F9]	HA601002	Human, Mouse, Rat	WB,ELISA,IF-Cell
Cell treatment	IL-1 beta Recombinant Mouse mAb [A7F7-R]	HA601036	Human, Mouse, Rat	WB,ELISA,IHC-P,IF-Tissue
Cell treatment	IL-1 beta Recombinant Mouse mAb [A7F7-R] - BSA and Azide free	HA610004	Human, Mouse, Rat	WB,ELISA,IHC-P,IF-Tissue
	IL1R1 Mouse mAb [D12-D2]	EM1708-33	Human	WB,FC
	IL-2 Rabbit pAb	ER63748	Human, Rat	WB,IHC-P
	IL-2 Mouse mAb [A8B8]	HA601015	Human	WB,ELISA
	IL-2 Mouse mAb [A8B9]	HA601016	Human	WB,ELISA
	IL24 Rabbit pAb	ER1911-33	Human	WB
	IL3RA Mouse mAb [C5-E11]	EM1709-19	Human	WB,IF-Cell,IHC-P,FC
Relative expression	IL3RA / CD123 Recombinant Rabbit mAb [PSH09-18]	HA723076	Human	WB,IF-Cell,IHC-P,FC,IF-Tissue
	IL3RA / CD123 Rabbit pAb	HA500426	Human	WB
	IL3RA / CD123 Recombinant mAb	IRS104RB	Human	mIHC
	IL3RA / CD123 Rabbit pAb	HA500291	Human	WB,IF-Cell,IHC-P,FC
	IL-4 Mouse mAb [A7D7]	HA601034	Human	ELISA
	IL-4 Recombinant Rabbit mAb [JA86-10]	ET1704-28	Human	WB,IF-Cell,IF-Tissue,IHC-P
	IL-4 Rabbit pAb	ER1706-55	Human	WB,IF-Cell,IHC-P
	IL-5 Rabbit pAb	ER1911-38	Human, Rat	WB,IHC-P,IF-Tissue
	IL-5 Rabbit pAb	HA500081	Human	WB,IF-Cell,ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	IL-6 Mouse mAb [1-6]	EM1701-45	Mouse, Rat	WB,IF-Cell
	IL-6 Mouse mAb [3-A4]	M1303-3	Mouse	WB,IF-Cell,FC,ELISA
	IL-6 Mouse mAb	EM30301	Human	WB
	IL-6 Recombinant Rabbit mAb	IRS160RB	Human	mIHC
	IL-6 Rabbit pAb	R1412-2	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	IL-6 Rabbit pAb	R1203-2	Mouse	WB,ELISA
	IL-8 Rabbit pAb	ER1706-67	Human	ELISA,IF-Cell,IHC-P
Cell treatment	IL-8 Recombinant Rabbit mAb [PSH10-96]	HA723263	Human	WB,FC
Cell treatment	IL-8 Recombinant Rabbit mAb [PSH10-97] - BSA and Azide free	HA751387	Human	WB
Cell treatment	IL-8 Recombinant Rabbit mAb [PSH10-97]	HA723264	Human	WB
Cell treatment	IL-8 Recombinant Rabbit mAb [PSH10-96] - BSA and Azide free	HA751386	Human	WB,FC
Cell treatment	IL-8 Rabbit pAb	ER1901-61	Human	WB,IHC-P,IF-Cell,FC
Cell treatment	IL-8 Rabbit pAb	R1511-15	Human	WB,IHC-P,IF-Cell,FC
Cell treatment	IP10 Recombinant Rabbit mAb [JA10-82]	ET1704-27	Human	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	IP10 Recombinant Rabbit mAb [JA10-82] - BSA and Azide free	HA750411	Human	WB,IF-Cell,IF-Tissue,IHC-P
KD	IRAK-1 Recombinant Rabbit mAb [JE60-03]	HA720079	Human, Mouse	WB,IHC-P,IP
	IRAK-1 Rabbit pAb	HA500377	Human, Mouse	WB,IHC-P,IF-Cell
	IRAK-1 Rabbit pAb	ER63778	Human, Mouse, Rat	WB,IF-Cell
	IRAK4 Recombinant Rabbit mAb [JE41-77]	HA721347	Human, Mouse, Monkey	WB,IHC-P
	IRAK4 Mouse mAb [B9-H2]	EM1709-86	Human, Mouse, Monkey	WB,IHC-P,FC
KD	IRF3 Recombinant Rabbit mAb [SD2062]	ET1612-14	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
KD	IRF3 Recombinant Rabbit mAb [SD2062] - BSA and Azide free	HA750283	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
KD	JAK2 Mouse mAb [6-D3]	M1501-8	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	JAK2 Recombinant Rabbit mAb [SY0245]	ET1607-35	Human, Mouse, Rat	WB,IF-Cell
	JAK2 Rabbit pAb	RT1343	Human, Mouse, Rat	WB,IP,IF,IHC-P
	JAK2 Rabbit pAb	ER1706-58	Human, Mouse, Rat	IF-Cell,IHC-P
KD	JAK2 Recombinant Rabbit mAb [SY0245] - BSA and Azide free	HA750118	Human, Mouse, Rat	WB,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	JNK1 Rabbit pAb	R1309-1	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC
	JNK1 Recombinant Rabbit mAb [JE63-04]	HA721899	Human, Mouse, Rat	WB,IHC-P
	JNK1+JNK2+JNK3 Recombinant Rabbit mAb [SA43-06]	ET1601-28	Human, Mouse, Rat, Cow, Monkey, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP
	JNK1+JNK2+JNK3 Recombinant Rabbit mAb [SA43-06] - BSA and Azide free	HA750018	Human, Mouse, Rat, Cow, Monkey, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP
	JNK1+JNK3 Recombinant Rabbit mAb [JJ08-86]	ET1701-47	Human, Mouse, Rat	WB,IP,IHC-P
KD	JNK2 Recombinant Rabbit mAb [SC05-85]	ET1610-11	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
	JNK2 Rabbit pAb	HA500110	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	JNK3 Recombinant Rabbit mAb [SD082-09]	ET1612-68	Human, Mouse, Rat	WB,IF-Tissue,IHC-P
	KDR / VEGFR-2 Mouse mAb [1G1]	RT1347	Human, Mouse, Rat	WB,IP,IF,IHC-P
Relative expression	LAT Recombinant Rabbit mAb [PSH09-76]	HA723141	Human	WB,IF-Cell,FC
Relative expression	LAT Recombinant Rabbit mAb [PSH09-76] - BSA and Azide free	HA751318	Human	WB,IF-Cell,FC
	LAT Rabbit pAb	ER63843	Human, Mouse, Rat	WB
	LBP Rabbit pAb	ER63848	Human, Mouse, Rat	WB,IHC-P
Relative expression	Lck Recombinant Rabbit mAb [PSH09-25]	HA723083	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC
	Lck Rabbit pAb	HA500280	Human, Rat	WB,FC
	Lck Rabbit pAb	ER63849	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Relative expression	Lck Recombinant Rabbit mAb [PSH09-25] - BSA and Azide free	HA751279	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC
	LT-alpha Mouse mAb [6-C10-10]	M0807-15	Human, Rat	WB
	LT-alpha Rabbit pAb	0807-5	Human, Rat	WB
	LT-alpha Mouse mAb [C2-F2]	M0807-12	Human	WB
Relative expression	LT-alpha (C-term) Rabbit pAb	0807-6	Human, Mouse, Rat	WB,IF-Cell
	LT beta Rabbit pAb	ER63881	Human, Mouse	WB,IF-Cell
	LTBR Rabbit pAb	ER1912-51	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	LTBR Mouse mAb [1D1C5C8]	EM1711-57	Human, Mouse	WB,FC
	LY96 Rabbit pAb	ER61807	Human, Mouse	WB
KD	Lyn Recombinant Rabbit mAb [JE44-75]	HA721991	Human, Mouse, Rat	WB,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Lyn Rabbit pAb	HA500142	Human, Mouse, Rat	WB
	Lyn Rabbit pAb	ER1912-57	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
KD	MALT1 Recombinant Rabbit mAb [PSH02-31] - BSA and Azide free	HA750563	Human, Mouse, Rat	WB,IHC-P
	MALT1 Rabbit pAb	ER63902	Human, Mouse	WB,IHC-P,IF-Cell
	MALT1 Recombinant Rabbit mAb [PSH02-32] - BSA and Azide free	HA750722	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	MALT1 Recombinant Rabbit mAb [PSH02-32]	HA721689	Human, Mouse, Rat	WB,IHC-P,IF-Cell
KD	MALT1 Recombinant Rabbit mAb [PSH02-31]	HA721244	Human, Mouse, Rat	WB,IHC-P
	MAP3K14 Mouse mAb [2G2B4]	EM1712-23	Human	WB,IHC-P,FC
	MAP3K7 Mouse mAb [1G5]	RT1377	Human, Mouse, Rat	WB,IP,IF
KD	MAPK6 / ERK3 Recombinant Rabbit mAb [PSH03-39] - BSA and Azide free	HA750873	Human, Mouse, Rat	WB,IF-Cell,FC
KD	MAPK6 / ERK3 Recombinant Rabbit mAb [PSH03-39]	HA721984	Human, Mouse, Rat	WB,IF-Cell,FC
	MD2 Recombinant Rabbit mAb [PSH07-72] - BSA and Azide free	HA751184	Human, Mouse, Rat	WB
	MD2 Recombinant Rabbit mAb [PSH07-72]	HA722888	Human, Mouse, Rat	WB
	MIP-1b Rabbit pAb	ER65660	Human	WB
	MIP-3 beta Rabbit pAb	ER63961	Human	IHC-P
	MMP-3 Recombinant Rabbit mAb [JM46-22]	ET1705-98	Human, Mouse, Rat	WB,IHC-P
	MMP-3 Rabbit pAb	ER1706-77	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue
	MNDA Rabbit pAb	ER63983	Human	WB,IHC-P
	Mouse CCL2 / MCP-1 Recombinant Rabbit mAb [PSH12-44] - BSA and Azide free (Capture)	HA725040	Mouse	ELISA(Cap)
	Mouse CCL2 / MCP-1 Recombinant Rabbit mAb [PSH12-45] - BSA and Azide free (Detector)	HA725041	Mouse	ELISA(Det)
	Biotin Conjugated Mouse CCL2 / MCP-1 Recombinant Rabbit mAb [PSH12-45] - Detector	HA725042B	Mouse	ELISA(Det),ELISA
	Mouse CCL5 Recombinant Rat mAb	HA601395	Human, Mouse	ELISA
	Mouse CCL5/RANTES Recombinant Rabbit mAb [PSH14-15] - BSA and Azide free (Capture)	HA725046	Mouse	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse CCL5/RANTES Recombinant Rabbit mAb [PSH14-16] - BSA and Azide free (Detector)	HA725047	Mouse	ELISA(Det)
	Biotin Conjugated Mouse CCL5/RANTES Recombinant Rabbit mAb [PSH14-16] - Detector	HA725048B	Mouse	ELISA(Det),ELISA
	APC Conjugated Anti-Mouse CD4 Recombinant Antibody [PSH04-99] - Rabbit IgG (Chimeric)	HA720204F	Mouse	FC
	FITC Conjugated Anti-Mouse CD4 Recombinant Antibody [PSH04-99] - Rabbit IgG (Chimeric)	HA720202F	Mouse	FC
	PE Conjugated Anti-Mouse CD4 Recombinant Antibody [PSH04-99] - Rabbit IgG (Chimeric)	HA720203F	Mouse	FC
	Mouse IL-10 Recombinant Rabbit mAb [PS01-93] - BSA and Azide free (Detector)	HA721518	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-10 Recombinant Rabbit mAb [PS01-93] - Detector	HA721519	Mouse	ELISA(Det)
	Mouse IL-10 Recombinant Rabbit mAb [PS01-85] - BSA and Azide free (Capture)	HA721517	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-13 Recombinant Rabbit mAb [PSH02-05]	HA722641B	Mouse	ELISA(Det)
	Mouse IL-13 Recombinant Rabbit mAb [PSH02-05] - BSA and Azide free (Detector)	HA721785	Mouse	ELISA(Det)
	Mouse IL-13 Recombinant Rabbit mAb [PSH02-04] - BSA and Azide free (Capture)	HA721784	Mouse	ELISA(Cap)
	Mouse IL-1 alpha Recombinant Rabbit mAb [PSH13-42] - BSA and Azide free (Capture)	HA723529	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-1 alpha Recombinant Rabbit mAb [PSH13-43] - Detector	HA723531B	Mouse	ELISA(Det),ELISA
	Mouse IL-1 alpha Recombinant Rabbit mAb [PSH13-43] - BSA and Azide free (Detector)	HA723530	Mouse	ELISA(Det)
	Mouse IL-1 alpha Recombinant Rabbit mAb [PSH13-40] - BSA and Azide free (Capture)	HA723526	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-1 alpha Recombinant Rabbit mAb [PSH13-41] - Detector	HA723528B	Mouse	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse IL-1 alpha Recombinant Rabbit mAb [PSH13-41] - BSA and Azide free (Detector)	HA723527	Mouse	ELISA(Det)
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH03-32] - BSA and Azide free (Detector)	HA721977	Mouse	ELISA(Det)
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH03-31] - BSA and Azide free (Capture)	HA721976	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-1 beta Recombinant Rabbit mAb [PSH03-32] - Detector	HA722625B	Mouse	ELISA(Det),ELISA
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH09-58] - BSA and Azide free	HA723120	Mouse	ELISA
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH04-66] - BSA and Azide free (Capture)	HA722166	Mouse	ELISA(Cap)
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH04-67] - BSA and Azide free (Detector)	HA722167	Mouse	ELISA(Det)
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH09-59] - BSA and Azide free	HA723121	Mouse	ELISA
	Biotin Conjugated Mouse IL-1 beta Recombinant Rabbit mAb [PSH04-67]	HA722608B	Mouse	ELISA(Det)
	Mouse IL-5 Recombinant Rabbit mAb [PSH12-78] - BSA and Azide free (Detector)	HA725035	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-5 Recombinant Rabbit mAb [PSH12-78] - Detector	HA725036B	Mouse	ELISA(Det),ELISA
	Mouse IL-5 Recombinant Rabbit mAb [PSH12-77] - BSA and Azide free (Capture)	HA725034	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PSH0-60] - Detector	HA722647B	Mouse	ELISA(Det)
	Mouse IL-6 Recombinant Rabbit mAb [PS01-74] - BSA and Azide free (Capture/Detector)	HA721600	Mouse	ELISA(Cap),ELISA(Det)
	Mouse IL-6 Recombinant Rabbit mAb [PSH0-59] - BSA and Azide free (Detector)	HA721412	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PS01-83] - BSA and Azide free (Detector)	HA721602	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PSH0-59] - Detector	HA722646B	Mouse	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse IL-6 Recombinant Rabbit mAb [PS01-83] - BSA and Azide free (Detector)	HA721601	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PS01-74] - Detector	HA723643B	Mouse	ELISA(Det)
	Mouse IL-6 Recombinant Rabbit mAb [PSH14-47] - BSA and Azide free (Capture)	HA723642	Mouse	ELISA(Cap)
	Mouse IL-6 Recombinant Rabbit mAb [PSH0-60] - BSA and Azide free (Detector)	HA721413	Mouse	ELISA(Det)
	HRP Conjugated Mouse IL-6 Mouse mAb [3-A4]	M1303-2	Mouse	WB,ELISA
	Mouse IL-6 Recombinant Rabbit mAb [PSH0-58] - BSA and Azide free (Capture)	HA721411	Mouse	ELISA(Cap)
	Mouse PD-L1 Recombinant Rabbit mAb [PSH04-83] - BSA and Azide free (Capture)	HA722189	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse PD-L1 Recombinant Rabbit mAb [PSH04-84]	HA722602B	Mouse	ELISA(Det)
	Mouse PD-L1 Recombinant Rabbit mAb [PSH04-84] - BSA and Azide free (Detector)	HA722190	Mouse	ELISA(Det)
Cell treatment	Mouse PD-L1 Recombinant Rabbit mAb [PSH05-04]	HA722209	Mouse	FC
Cell treatment	FITC Conjugated Mouse PD-L1 Recombinant Rabbit mAb [PSH05-04]	HA720213F	Mouse	FC
Cell treatment	Mouse PD-L1 Recombinant Rabbit mAb [PSH05-03]	HA722208	Mouse	FC
Cell treatment	Mouse RANKL Recombinant Rabbit mAb [PSH07-90] - BSA and Azide free (Capture)	HA722922	Mouse	ELISA(Cap)
Cell treatment	Mouse RANKL Recombinant Rabbit mAb [PSH07-92] - BSA and Azide free (Detector)	HA722924	Mouse	ELISA(Det)
Cell treatment	Biotin Conjugated Mouse RANKL Recombinant Rabbit mAb [PSH07-92]	HA722925B	Mouse	ELISA(Det)
Cell treatment	Mouse RANKL Recombinant Rabbit mAb [PSH07-91] - BSA and Azide free (Capture)	HA722923	Mouse	ELISA(Cap)
	Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-30] - BSA and Azide free (Detector)	HA721975	Mouse, Rat	ELISA(Det)
	Biotin Conjugated Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-30]	HA722626B	Mouse, Rat	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-29] - BSA and Azide free (Capture)	HA721974	Mouse, Rat	ELISA(Cap)
	Mouse TNF alpha Recombinant Rabbit mAb [PSH03-76] - BSA and Azide free (Capture)	HA722040	Mouse	ELISA(Cap)
	Mouse TNF alpha Recombinant Rabbit mAb [PSH03-77] - BSA and Azide free (Detector)	HA722041	Mouse	ELISA(Det)
	Biotin Conjugated Mouse TNF alpha Recombinant Rabbit mAb [PSH03-77]	HA722623B	Mouse	ELISA(Det)
	Mouse VCAM1 Recombinant Rabbit mAb [PSH16-32] - BSA and Azide free (Detector)	HA723830	Mouse	ELISA(Det)
	HRP Conjugated Mouse VCAM1 Recombinant Rabbit mAb [PSH16-32] - Detector	HA723830H	Mouse	ELISA(Det),ELISA
	Mouse VCAM1 Recombinant Rabbit mAb [PSH16-33] - BSA and Azide free (Capture)	HA723831	Mouse	ELISA(Cap)
	Mouse VCAM1 Recombinant Rabbit mAb [PSH16-31] - BSA and Azide free (Capture)	HA723829	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse VCAM1 Recombinant Rabbit mAb [PSH16-32] - Detector	HA723830B	Mouse	ELISA(Det),ELISA
	Mouse VEGF Recombinant Rabbit mAb [PSH02-72] - BSA and Azide free (Capture)	HA721850	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse VEGF Recombinant Rabbit mAb [PSH02-73]	HA722636B	Mouse	ELISA(Det)
	Mouse VEGF Recombinant Rabbit mAb [PSH02-73] - BSA and Azide free (Detector)	HA721851	Mouse	ELISA(Det)
	MyD88 Recombinant Rabbit mAb [PSH05-58]	HA722378	Human, Monkey	WB,IF-Cell,IHC-P,FC,IP
	MyD88 Recombinant Rabbit mAb [PSH09-74]	HA723139	Human, Mouse, Rat, Monkey	WB
	MyD88 Recombinant Rabbit mAb [PSH09-74] - BSA and Azide free	HA751316	Human, Mouse, Rat, Monkey	WB
	MyD88 Recombinant Rabbit mAb [SC65-04] - BSA and Azide free	HA750231	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
	MyD88 Recombinant Rabbit mAb [PSH05-58] - BSA and Azide free	HA751021	Human, Monkey	WB,IF-Cell,IHC-P,FC,IP
	MyD88 Recombinant Rabbit mAb [SC65-04]	ET1610-81	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	NAK / TBK1 Recombinant Rabbit mAb [PSH11-69]	HA723365	Human, Mouse, Rat	WB,IF-Cell,FC,IP
	NAK / TBK1 Recombinant Rabbit mAb [PSH11-69] - BSA and Azide free	HA751425	Human, Mouse, Rat	WB,IF-Cell,FC,IP
KD	NAK / TBK1 Mouse mAb [A8D9]	HA601045	Human, Mouse, Rat	WB
	NAK / TBK1 Mouse mAb [2G1]	RT1413	Human, Mouse, Rat	WB,IP
KD	NAK / TBK1 Mouse mAb [A8D10]	HA601044	Human, Mouse, Rat	WB
	NFkB Inducing Kinase NIK Rabbit pAb	ER1913-80	Human, Mouse	WB,IHC-P,IF-Tissue,FC
KD	NF-kB p100 / NFKB2 Mouse mAb [16B2]	EM1901-78	Human	WB,IHC-P,FC
KD	NF-kB p100 / NFKB2 Mouse mAb [16B3]	EM1901-79	Human, Mouse, Rat	WB,IHC-P,FC
KD	NF-kB p100 / NFKB2 Recombinant Rabbit mAb [JM82-03]	ET1705-45	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
KD	NF-kB p105 / p50 Recombinant Rabbit mAb [SZ20-01]	ET1603-18	Human, Mouse, Rat	WB,IHC-P,ICC,FC,IF-Tissue
	NF-kB p105 / p50 Rabbit pAb	R1309-9	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	NF-kB p105 / p50 Mouse mAb [A2-7]	M1505-8	Human, Mouse	WB,IF-Cell,IHC-P
KD	NF-kB p105 / p50 Recombinant Rabbit mAb [SZ20-01] - BSA and Azide free	HA750060	Human, Mouse, Rat	WB,IHC-P,ICC,FC,IF-Tissue
KO, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [PSH0-27] - BSA and Azide free	HA750601	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
KO, KD, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [SZ10-04]	ET1603-12	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KO, KD, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [SZ10-04] - BSA and Azide free	HA750058	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KO, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [PSH0-27]	HA721307	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
	NF-kB p65 Rabbit pAb	ER0815	Human, Mouse, Rat, Zebrafish	WB,IHC-P,FC
	NF-kB p65 Rabbit pAb	HA500402	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
	NLRP3 Rabbit pAb	ER1706-72	Human, Mouse	WB,IF-Cell,IHC-P,FC
Cell treatment	NLRP3 Recombinant Rabbit mAb [SC06-23] - BSA and Azide free	HA750236	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Cell treatment	NLRP3 Recombinant Rabbit mAb [SC06-23]	ET1610-93	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Relative expression	NLRP3 Recombinant Rabbit mAb [PSH19-43]	HA724075	Human, Mouse	WB
Cell treatment	Nrf2 Rabbit pAb	ER1706-41	Human, Mouse, Rat	IHC-P,IF-Cell,WB
Cell treatment	Nrf2 Recombinant Rabbit mAb [PSH11-20]	HA723302	Human, Mouse, Rat	WB,IF-Cell,ChIP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Nrf2 Rabbit pAb	R1312-8	Human, Mouse	WB,IHC-P,IF-Tissue
Cell treatment	Nrf2 Recombinant Rabbit mAb [PSH11-20] - BSA and Azide free	HA751401	Human, Mouse, Rat	WB,IF-Cell,ChIP
Cell treatment	Nrf2 Recombinant Rabbit mAb [JE30-14]	HA721432	Human, Mouse	WB,IF-Cell
	p21 Recombinant Rabbit mAb [PSH03-93]	HA722065	Mouse	WB,IF-Cell,IHC-P,FC
	p21 Mouse mAb [A8C11]	HA601018	Human	IHC-P
	p21 Recombinant Rabbit mAb [PSH03-93] - BSA and Azide free	HA750918	Mouse	WB,IF-Cell,IHC-P,FC
	p21 Recombinant Mouse mAb [A8C11-R]	HA601221	Human	IHC-P
	p21 Recombinant Rabbit mAb	IRS260RB	Mouse	mIHC
	p21 Recombinant Rabbit mAb [JE00-63]	HA722685	Human	WB,IF-Cell,IHC-P,IP
	p21 Recombinant Mouse mAb [A8C11-R] - BSA and Azide free	HA610098	Human	IHC-P
	p21 Mouse mAb [A8C12]	HA601017	Human	IHC-P
	PAF1 Recombinant Rabbit mAb [JE38-80]	HA722401	Human, Mouse, Rat	WB
	PAI1 Rabbit pAb	HA500124	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	PAI1 Mouse mAb [A9G8]	HA601152	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
	PARP1 Rabbit pAb	ER1802-67	Human, Mouse	WB,IF-Cell,FC
Cell treatment, KD	PARP1 Recombinant Rabbit mAb [SU03-68]	ET1608-56	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
Cell treatment, KD	PARP1 Recombinant Rabbit mAb [SU03-68] - BSA and Azide free	HA750154	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	PD-L1 Recombinant Rabbit mAb [JJ08-95]	ET1701-41	Human	WB,IHC-P
Cell treatment, Relative expression, KD	PD-L1 Recombinant Rabbit mAb [PD01-02]	HA721176	Human	IHC-P,mIHC,WB
	PD-L1 Recombinant Rabbit mAb [PD00-97]	HA722300	Human	IHC-P
	PD-L1 Recombinant Rabbit mAb	IRS059	Mouse	mIHC
	PD-L1 Recombinant Rabbit mAb	IRS009	Human	mIHC
Cell treatment	PD-L1 Recombinant Rabbit mAb [PSH04-80]	HA722184	Mouse	WB,IHC-P,mIHC
Cell treatment	PD-L1 Recombinant Rabbit mAb [PSH04-80] - BSA and Azide free	HA750974	Mouse	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	PD-L1 Recombinant Rabbit mAb [JJ08-95] - BSA and Azide free	HA750309	Human	WB,IHC-P
	PD-L1 Recombinant Rabbit Multiclonal Antibody [PSH13-23]	HA723510	Human, Mouse	WB
	PD-L1 Recombinant Rabbit mAb [PD00-97] - BSA and Azide free	HA751015	Human	IHC-P
	PD-L1 Recombinant Rabbit Multiclonal Antibody [PSH13-23] - BSA and Azide free	HA751459	Human, Mouse	WB
Cell treatment	Phospho-ATM (S1981) Recombinant Rabbit mAb [JM93-25]	ET1705-50	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Cell treatment	Phospho-ATM (S1981) Recombinant Rabbit mAb [JM93-25] - BSA and Azide free	HA750432	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Cell treatment	Phospho-Bcl-2 (S70) Recombinant Rabbit mAb [PSH11-02] - BSA and Azide free	HA751392	Human	WB,IF-Cell,FC
Cell treatment	Phospho-Bcl-2 (S70) Recombinant Rabbit mAb [PSH11-02]	HA723269	Human	WB,IF-Cell,FC
Cell treatment	Phospho-BTK (Y223) Recombinant Rabbit mAb [ST51-05]	ET1609-55	Human, Mouse	WB,IP
Cell treatment	Phospho-Erk1 (T202)+Erk2 (T185) Recombinant Rabbit mAb [SZ2-4] - BSA and Azide free	HA750061	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
Cell treatment	Phospho-Erk1 (T202)+Erk2 (T185) Recombinant Rabbit mAb [SZ2-4]	ET1603-22	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
Cell treatment	Phospho-Erk1 (T202 + Y204) + Erk2 (T185 + Y187) Recombinant Rabbit mAb [SC58-01] - BSA and Azide free	HA750199	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP
Cell treatment	Phospho-Erk1 (T202 + Y204) + Erk2 (T185 + Y187) Recombinant Rabbit mAb [SC58-01]	ET1610-13	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP
Cell treatment	Phospho-IKB alpha (S32) Recombinant Rabbit mAb [ST53-05] - BSA and Azide free	HA750196	Human, Mouse	WB,IF-Cell,IP
Cell treatment	Phospho-IKB alpha (S32) Recombinant Rabbit mAb [ST53-05]	ET1609-78	Human, Mouse	WB,IF-Cell,IP
Cell treatment	Phospho-IKB alpha (S32/S36) Recombinant Rabbit mAb [PSH06-76] - BSA and Azide free	HA751119	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Cell treatment	Phospho-IKB alpha (S32/S36) Recombinant Rabbit mAb [PSH06-76]	HA722770	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Cell treatment	Phospho-IKB alpha (S36) Recombinant Rabbit mAb [JE56-65]	HA721802	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-IRF3 (S386) Recombinant Rabbit mAb [SU03-28] - BSA and Azide free	HA750138	Human	WB,IF-Cell,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-IRF3 (S386) Recombinant Rabbit mAb [SU03-28]	ET1608-22	Human	WB,IF-Cell,IF-Tissue
Cell treatment	Phospho-IRF3 (S396) Recombinant Rabbit mAb [PSH06-75] - BSA and Azide free	HA751121	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-IRF3 (S396) Recombinant Rabbit mAb [PSH06-75]	HA722772	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-JAK2 (Y1007 + Y1008) Recombinant Rabbit mAb [SY24-03]	ET1607-34	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
Cell treatment	Phospho-JAK2 (Y1007 + Y1008) Recombinant Rabbit mAb [SY24-03] - BSA and Azide free	HA750117	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
Cell treatment	Phospho-JNK1/2/3 (T183 + T183 + T221) Recombinant Rabbit mAb [ST500] - BSA and Azide free	HA750173	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
Cell treatment	Phospho-JNK1/2/3 (T183 + T183 + T221) Recombinant Rabbit mAb [ST500]	ET1609-42	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
Cell treatment, Relative expression	Phospho-LYN (Y397)/LCK (Y394)/HCK (Y411)/BLK (Y389) Recombinant Rabbit mAb [PSH09-43]	HA723105	Human, Mouse, Rat	WB,IF-Cell
Cell treatment, Relative expression	Phospho-LYN (Y397)/LCK (Y394)/HCK (Y411)/BLK (Y389) Recombinant Rabbit mAb [PSH09-43] - BSA and Azide free	HA751290	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-NF-kB p65 (S529) Recombinant Rabbit mAb [SP07-00]	ET1604-27	Human	WB,IHC-P,FC,IP
Cell treatment	Phospho-NF-kB p65 (S536) Recombinant Rabbit mAb [PSH10-58]	HA723223	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-NF-kB p65 (S536) Recombinant Rabbit mAb [PSH10-58] - BSA and Azide free	HA751361	Human, Mouse, Rat	WB,IF-Cell
Cell treatment, Relative expression	Phospho-Nrf2 (S40) Recombinant Rabbit mAb [SU0334]	ET1608-28	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
Cell treatment, Relative expression	Phospho-Nrf2 (S40) Recombinant Rabbit mAb [SU0334] - BSA and Azide free	HA750142	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
Cell treatment	Phospho-PKC alpha/beta II (T638/T641) Recombinant Rabbit mAb [PSH16-27] - BSA and Azide free	HA751595	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,Dot Blot
Cell treatment	Phospho-PKC alpha/beta II (T638/T641) Recombinant Rabbit mAb [PSH16-27]	HA723825	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,Dot Blot

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-PKC alpha/beta II (T638/T641) Recombinant Rabbit mAb [PSH15-91]	HA723796	Human, Mouse, Rat	WB,IF-Cell,Dot Blot
Cell treatment	Phospho-PKC alpha/beta II (T638/T641) Recombinant Rabbit mAb [PSH15-91] - BSA and Azide free	HA751582	Human, Mouse, Rat	WB,IF-Cell,Dot Blot
Cell treatment	Phospho-PKC alpha (T497) Recombinant Rabbit mAb [JE44-81]	HA722278	Human, Mouse, Rat	WB,IF-Cell,FC
Cell treatment	Phospho-PKC alpha (T638) Recombinant Rabbit mAb [JF0964]	ET1702-17	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,Dot Blot
Cell treatment	Phospho-PKC delta (Y311) Recombinant Rabbit mAb [JE44-20]	HA722458	Human, Mouse, Rat	WB
Cell treatment	Phospho-PKC (S729) Recombinant Rabbit mAb [JE77-54]	HA722846	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-PKC zeta (T560) Recombinant Rabbit mAb [SN0753]	ET1611-78	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-PLC-gamma-1 (S1248) Recombinant Rabbit mAb [PSH13-62]	HA723551	Human	WB,IP
Cell treatment	Phospho-PLC-gamma-1 (S1248) Recombinant Rabbit mAb [PSH13-62] - BSA and Azide free	HA751481	Human	WB,IP
Cell treatment	Phospho-PLC-gamma-1 (Y783) Recombinant Rabbit mAb [PSH06-29] - BSA and Azide free	HA751066	Human, Mouse, Rat	WB
Cell treatment	Phospho-PLC-gamma-1 (Y783) Recombinant Rabbit mAb [PSH06-29]	HA722661	Human, Mouse, Rat	WB
Cell treatment, KD	Phospho-STAT3 (S727) Recombinant Rabbit mAb [SY24-09] - BSA and Azide free	HA750121	Human, Rat, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP
Cell treatment, KD	Phospho-STAT3 (S727) Recombinant Rabbit mAb [SY24-09]	ET1607-39	Human, Rat, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP
Cell treatment	Phospho-STAT3 (Y705) Recombinant Rabbit mAb [SZ43-01]	ET1603-40	Human, Mouse, Rat	WB,IHC-P,IP
Cell treatment	Phospho-STAT3 (Y705) Recombinant Rabbit mAb [SZ43-01] - BSA and Azide free	HA750067	Human, Mouse, Rat	WB,IHC-P,IP
Cell treatment	Phospho-STING (S366) Recombinant Rabbit mAb [PSH09-72]	HA723137	Human	WB,Dot Blot,IP
Cell treatment	Phospho-STING (S366) Recombinant Rabbit mAb [PSH09-72] - BSA and Azide free	HA751314	Human	WB,Dot Blot,IP
Cell treatment	Phospho-Syk (Y352) Recombinant Rabbit mAb [PSH11-66]	HA723362	Human	WB,IF-Cell,FC
Cell treatment	Phospho-Syk (Y352) Recombinant Rabbit mAb [PSH11-66] - BSA and Azide free	HA751422	Human	WB,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-Syk (Y525 + Y526) Recombinant Rabbit mAb [PSH03-70]	HA722034	Human	WB,IF-Cell
Cell treatment	Phospho-Syk (Y525 + Y526) Recombinant Rabbit mAb [PSH03-70] - BSA and Azide free	HA750904	Human	WB,IF-Cell
Cell treatment	Phospho-TAK1 (S412) Recombinant Rabbit mAb [PSH14-64] - BSA and Azide free	HA751525	Human, Mouse, Rat	WB,FC
Cell treatment	Phospho-TAK1 (S412) Recombinant Rabbit mAb [PSH14-64]	HA723662	Human, Mouse, Rat	WB,FC
Cell treatment	Phospho-TAK1 (S439) Recombinant Rabbit mAb [JE44-93]	HA722242	Human, Mouse, Rat	WB
	PI3 Kinase p110 beta Recombinant Rabbit mAb [JE65-38]	HA722474	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	PI3 Kinase p110 α Recombinant Rabbit mAb [SJ0186]	ET1606-36	Human, Mouse	WB,IF-Cell,IF- Tissue,IP,IHC-P
	PIAS4 Rabbit pAb	HA500411	Human, Mouse, Rat	WB,IHC-P
	PKC Rabbit pAb	ER64620	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
KD	PKC alpha Recombinant Rabbit mAb [SU31-08]	ET1608-15	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC-P,FC,IP
	PKC beta Rabbit pAb	ER1802-14	Human, Mouse	WB,IF-Cell,IHC-P
	PKC beta 1 Recombinant Rabbit mAb [JM43-26]	ET1705-30	Human, Mouse, Rat	WB,IHC-P,FC
	PKC beta 2 Recombinant Rabbit mAb [ST48-06]	ET1609-44	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC-P,IP,FC
	PKC delta Recombinant Rabbit mAb [JJ093-02]	ET1701-85	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC-P
	PKC delta Recombinant Rabbit mAb [JJ093-02] - BSA and Azide free	HA750325	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC-P
	PKC epsilon Recombinant Rabbit mAb [JA38-21]	ET1704-85	Human, Mouse, Rat	WB,IHC-P
	PKC eta Recombinant Rabbit mAb [JE32-41]	HA722389	Human, Mouse, Rat, Monkey	WB
	PKC gamma Recombinant Rabbit mAb [JG69-32]	ET7108-42	Human, Mouse, Rat	WB,IHC-P
	PKC theta Rabbit pAb	ER64621	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Relative expression	PKC zeta Recombinant Rabbit mAb [PSH06-74] - BSA and Azide free	HA751122	Human	WB
Relative expression	PKC zeta Recombinant Rabbit mAb [PSH06-74]	HA722773	Human	WB
	PLCG1 Mouse mAb [D11-C2]	EM1707-02	Human, Rat, Monkey	WB,IF-Cell,IHC-P,FC
	PLCG 2 Recombinant Rabbit mAb [JE39-12]	HA721477	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	PLCG2 Rabbit pAb	R1512-4	Human	WB,IF-Cell,IHC-P,FC
	PLC-gamma-1 Recombinant Rabbit mAb [SN06-10]	ET1611-25	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
Cell treatment, Relative expression, KO	Pro Caspase-3 Recombinant Rabbit mAb [SR03-01]	ET1602-39	Human, Monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Pro Caspase-3 Rabbit pAb	ER30804	Human, Mouse	WB,IF-Cell,IHC-P,FC,IHC-Fr
Cell treatment, Relative expression, KO	Pro Caspase-3 Recombinant Rabbit mAb [SR03-01] - BSA and Azide free	HA750048	Human, Monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	Pro Caspase-3 Recombinant Rabbit mAb [SZ02-08]	ET1603-26	Human, Monkey, Mouse	WB,IF-Tissue,IHC-P
	pro Caspase-7 Recombinant Rabbit mAb [SD20-57]	ET1612-28	Human	WB,IF-Cell,IF-Tissue,IHC-P
	RANK Rabbit pAb	ER1915-70	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	RANKL Rabbit pAb	HA500369	Human, Mouse, Rat	WB,IF-Cell
	RANTES Recombinant Rabbit mAb [JM03-45]	ET1705-70	Human, Mouse	WB,IHC-P,IF-Cell,IF-Tissue
	Rat IL-10 Recombinant Rabbit mAb [PSH02-22] - BSA and Azide free (Detector)	HA721787	Rat	ELISA(Det)
	Biotin Conjugated Rat IL-10 Recombinant Rabbit mAb [PSH02-22]	HA722637B	Rat	ELISA(Det)
	Rat IL-10 Recombinant Rabbit mAb [PSH02-21] - BSA and Azide free (Capture)	HA721786	Rat	ELISA(Cap)
	Biotin Conjugated Rat IL-1 beta Recombinant Rabbit mAb [PSH02-60]	HA722634B	Rat	ELISA(Det)
	Biotin Conjugated Rat IL-1 beta Recombinant Rabbit mAb [PS01-87] - Detector	HA721611	Rat	ELISA(Det),ELISA
	Rat IL-1 beta Recombinant Rabbit mAb [PSH02-59] - BSA and Azide free (Capture)	HA721836	Rat	ELISA(Cap)
	Rat IL-1 beta Recombinant Rabbit mAb [PS01-86] - BSA and Azide free (Capture)	HA721609	Rat	ELISA(Cap)
	Rat IL-1 beta Recombinant Rabbit mAb [PS01-87] - BSA and Azide free (Detector)	HA721610	Rat	ELISA(Det)
	Rat IL-1 beta Recombinant Rabbit mAb [PSH02-60] - BSA and Azide free (Detector)	HA721837	Rat	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	RelB Recombinant Rabbit mAb [SD07-39]	ET1612-18	Human	WB,IF-Tissue,IHC-P,IP
	RIP Recombinant Rabbit mAb [PSH19-57]	HA724083	Human, Mouse, Rat, Monkey	WB
	RIP Recombinant Rabbit mAb [JJ092-04]	ET1701-79	Human	WB,FC,IP
	ROCK1 Recombinant Rabbit mAb [ST05-19]	ET1609-53	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,IP
	ROCK1 Recombinant Rabbit mAb [ST05-19] - BSA and Azide free	HA750181	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,IP
KD	ROCK1 Rabbit pAb	ER40122	Human, Mouse, Rat	WB,IHC-P
	SERPINE1 Mouse mAb [H9-D5]	EM1709-36	Human, Mouse, Rat, Rabbit	IHC-P,FC
	STAT3 Rabbit pAb	RT1591	Human, Mouse, Rat	WB,IP,IF,IHC-P
KD	STAT3 Recombinant Rabbit mAb [SY34-01] - BSA and Azide free	HA750087	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	STAT3 Recombinant Rabbit mAb [SY24-08]	ET1607-38	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC
KD	STAT3 Recombinant Rabbit mAb [SY24-08] - BSA and Azide free	HA750120	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC
KD	STAT3 Recombinant Rabbit mAb [SY34-01]	ET1605-45	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	STING Recombinant Rabbit mAb [PSH07-44] - BSA and Azide free	HA751165	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue,IP
	STING Recombinant Rabbit mAb [JM03-47] - BSA and Azide free	HA750437	Human, Rat	WB,IHC-P,IF-Cell,FC,IF-Tissue,IP
	STING Recombinant Rabbit mAb [PSH07-44]	HA722832	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue,IP
	STING Recombinant Rabbit mAb [JM03-47]	ET1705-68	Human, Rat	WB,IHC-P,IF-Cell,FC,IF-Tissue,IP
KD	Syk Recombinant Rabbit mAb [PSH0-76] - BSA and Azide free	HA750655	Human	WB,IF-Cell
KD	Syk Recombinant Rabbit mAb [PSH0-76]	HA721494	Human	WB,IF-Cell
	Syk Rabbit pAb	HA500080	Human, Rat	WB,IF-Cell,IHC-P
KD	TAB1 Rabbit pAb	HA500400	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	TAB1 Mouse mAb [1G2]	RT1603	Human	WB,IP,IF
	TAB3 Rabbit pAb	HA500421	Human, Mouse	WB,IHC-P
	TAK1 Rabbit pAb	ER1802-90	Human, Mouse, Rat	IF-Cell,IHC-P,FC
KD	TAK1 Recombinant Rabbit mAb [JM73-19] - BSA and Azide free	HA750424	Human, Mouse, Rat, Monkey	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	TAK1 Recombinant Rabbit mAb [JM73-19]	ET1705-14	Human, Mouse, Rat, Monkey	WB
	TCAM2 Rabbit pAb	ER61885	Human, Mouse	WB
Relative expression	TGF beta 1 Recombinant Rabbit mAb [PD00-17] - BSA and Azide free	HA750485	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Tissue
	TGF beta 1 Recombinant Rabbit mAb [PSH18-59]	HA724012	Human, Mouse, Monkey	WB,IF-Cell,IHC-P,FC
	TGF beta 1 Mouse mAb [B2-F6]	EM1709-93	Human	WB,IHC-P,FC
	TGF beta 1 Rabbit pAb	HA500496	Human	WB
	Biotin Conjugated TGF beta 1 Recombinant Rabbit mAb [PSH15-22] - Detector	HA723727B	Human, Mouse, Rat	ELISA(Det),ELISA
	TGF beta 1 Recombinant Rabbit mAb [PSH15-21] - BSA and Azide free (Capture)	HA723725	Human, Mouse, Rat	ELISA(Cap)
Relative expression	TGF beta 1 Recombinant Rabbit mAb [PD00-17]	HA721143	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Tissue
	TGF beta 1 Recombinant Rabbit mAb	IRS205RB	Human	mIHC
	TGF beta 1 Recombinant Rabbit mAb [PSH15-22] - BSA and Azide free (Detector)	HA723726	Human, Mouse, Rat	ELISA(Det)
	TGF beta 3 Recombinant Rabbit mAb [PSH09-26] - BSA and Azide free	HA751280	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	TGF beta 3 Rabbit pAb	ER1917-64	Human, Mouse, Rat, Pig	WB,IHC-P,IF-Tissue,FC
	TGF beta 3 Recombinant Rabbit mAb [PSH09-26]	HA723084	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	TIRAP Recombinant Rabbit mAb [JE46-35]	HA722001	Human, Mouse	WB
	TLR4 Rabbit pAb	ER1706-43	Human, Mouse, Rat	IF-Cell,IHC-P,FC
	TLR4 Mouse mAb	RT1666	Human, Mouse, Rat	WB
	TNFAIP3 Recombinant Rabbit mAb [SN07-31]	ET1611-40	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
	TNFAIP3 Recombinant Rabbit mAb [SN07-31] - BSA and Azide free	HA750253	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
	TNFAIP3 Rabbit pAb	ER1706-32	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	TNF alpha Rabbit pAb	ER1919-22	Human, Mouse	WB
	TNF alpha Rabbit pAb	ER65189	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Cell treatment	TNF alpha Recombinant Rabbit mAb [PS01-07]	HA722022	Human, Mouse	WB
Cell treatment	TNF alpha Recombinant Rabbit mAb [PS01-07] - BSA and Azide free	HA750897	Human, Mouse	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	TNF alpha Rabbit pAb	R1203-1	Mouse	WB,IF-Cell,ELISA
	TNFR1 Mouse mAb [2G2]	RT1624	Human, Mouse, Rat	WB,IP,IF,IHC-P
	TNF Receptor I Rabbit pAb	HA500140	Human	WB,IHC-P,FC
	TNFSF11 Mouse mAb [8A7B9]	EM1710-42	Human, Monkey	WB,IHC-P,IF-Cell,FC
	TNFSF14 Rabbit pAb	HA500384	Human, Mouse	WB,IHC-P
	TNR27 Rabbit pAb	ER61889	Human, Mouse	WB
	TR13C Rabbit pAb	ER61707	Human	WB
	TRADD Recombinant Rabbit mAb [JE32-60]	HA722000	Human, Mouse, Rat	WB,IF-Cell
	TRADD Rabbit pAb	ER1901-79	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	TRADD Rabbit pAb	ER1917-98	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Relative expression	TRAF1 Recombinant Rabbit mAb [PSH13-88] - BSA and Azide free	HA751490	Human, Rat	WB,IP
	TRAF1 Rabbit pAb	ER65717	Human, Mouse	WB
Relative expression	TRAF1 Recombinant Rabbit mAb [PSH13-88]	HA723581	Human, Rat	WB,IP
	TRAF2 Recombinant Rabbit mAb [SD205-06] - BSA and Azide free	HA750297	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,IP
	TRAF2 Recombinant Rabbit mAb [SD205-06]	ET1612-5	Human, Mouse, Rat	WB,IHC-P,IP
	TRAF3 Recombinant Rabbit mAb [PSH01-30] - BSA and Azide free	HA750708	Human, Mouse, Rat	WB
	TRAF3 Rabbit pAb	ER65205	Human, Mouse	WB
	TRAF3 Recombinant Rabbit mAb [PSH01-30]	HA721654	Human, Mouse, Rat	WB
	TRAF6 Recombinant Rabbit mAb [JE01-42]	HA722857	Human, Mouse, Rat	WB,IF-Cell,FC
	TRAF6 Rabbit pAb	R1311-2	Human, Mouse, Rat, Chicken, Zebrafish	IF-Cell,WB
	TREM1 Rabbit pAb	ER1918-03	Human, Rat	WB,IHC-P,IF-Tissue,FC
	TRIF Recombinant Rabbit mAb [PSH05-39]	HA722268	Human	WB
	TRIF Rabbit pAb	ER65721	Human	WB
	TRIF Recombinant Rabbit mAb [PSH05-39] - BSA and Azide free	HA750998	Human	WB
	TRIM25 Mouse mAb [B12-B5]	EM1707-34	Human	WB,IF-Cell,IHC-P,FC
KD	TRIM25 Recombinant Rabbit mAb [JE57-43]	HA721506	Human	WB,IF-Cell
KD	TRIM25 Rabbit pAb	HA500104	Human, Rat	WB,IHC-P
	UBC9 Rabbit pAb	RT1639	Human, Mouse, Rat	WB,IP,IF,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	UBC9 Recombinant Rabbit mAb [SC0534]	ET1610-21	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	UBE2I Rabbit pAb	HA500401	Human, Mouse, Rat	WB
	Urokinase Rabbit pAb	HA500429	Human, Rat	WB,IF-Cell
	Urokinase Recombinant Rabbit mAb [JM106-09]	ET1703-26	Human, Mouse, Rat	WB,IHC-P
Relative expression, KD	VCAM1 Recombinant Rabbit mAb [SA05-04]	ET1601-18	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Relative expression, KD	VCAM1 Recombinant Rabbit mAb [SA05-04] - BSA and Azide free	HA750009	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	VCAM1 Rabbit pAb	R1512-13	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	VEGF Rabbit pAb	ER30607	Human	WB,IF-Cell,IHC-P,FC
	VEGF Recombinant Rabbit mAb [SP07-01] - BSA and Azide free	HA750074	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	VEGF Recombinant Rabbit mAb [SP07-01]	ET1604-28	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	VEGFA Recombinant Rabbit mAb [PSH19-17]	HA724051	Human, Mouse, Rat, Monkey	WB
KD	VEGF Receptor 1 Recombinant Rabbit mAb [SY09-09] - BSA and Azide free	HA750081	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	VEGF Receptor 1 Rabbit pAb	ER50803	Human, Mouse, Rat	WB,IHC-P
KD	VEGF Receptor 1 Recombinant Rabbit mAb [SY09-09]	ET1605-11	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	VEGF Receptor 1 Rabbit pAb	RT1651	Human, Mouse, Rat	WB,IP,IF,IHC-P
Relative expression	VEGF Receptor 2 Recombinant Rabbit mAb [SU03-42]	ET1608-33	Human, Rat	WB,IF-Cell,IHC-P,IP
	XIAP Rabbit pAb	ER65320	Human, Mouse, Rat	WB,IHC-P
KD	XIAP Recombinant Rabbit mAb [PSH04-30] - BSA and Azide free	HA750939	Human, Mouse, Rat, Monkey	WB,IHC-P
KD	XIAP Recombinant Rabbit mAb [PSH04-30]	HA722113	Human, Mouse, Rat, Monkey	WB,IHC-P
	XIAP/BIRC4 Rabbit pAb	ER1918-75	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC
Relative expression	ZAP70 Recombinant Rabbit mAb [JU08-39]	ET1706-42	Human	WB,IF-Cell,IHC-P,FC,IP
Relative expression	ZAP70 Recombinant Mouse mAb [A1B5-R]	HA601218	Human, Mouse, Rat	WB,IHC-P
Relative expression	ZAP70 Mouse mAb [A1B7]	EM1902-28	Human, Mouse, Rat	WB,IF-Cell,FC
Relative expression	ZAP70 Mouse mAb [A1B5]	EM1901-27	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	ZAP70 Recombinant Mouse mAb [A1B5-R] - BSA and Azide free	HA610096	Human, Mouse, Rat	WB,IHC-P
Relative expression	ZAP70 Recombinant Rabbit mAb [PSH11-52]	HA723344	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC