

迁移体研究-产品清单

高级验证	产品名	货号	反应种属	验证应用
	AER61 Recombinant Rabbit mAb [JE64-87]	HA721973	Human, Mouse, Rat	WB
	alpha Tubulin Rabbit pAb	ER130905	Human, Mouse, Rat, Zebrafish	WB,IHC-P,FC,IF-Cell,IF-Tissue
	alpha Tubulin Mouse mAb [A8-6]	M1501-1	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC
	alpha Tubulin Recombinant Rabbit mAb [PSH02-94]	HA721913	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
	alpha Tubulin Recombinant Rabbit mAb [PSH0-14]	HA721272	Human, Mouse, Rat	WB,IHC-P
	alpha Tubulin Recombinant Rabbit mAb [PSH02-95]	HA721914	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
	alpha Tubulin Recombinant Rabbit mAb [PSH0-14] - BSA and Azide free	HA750584	Human, Mouse, Rat	WB,IHC-P
	alpha Tubulin Recombinant Rabbit mAb [PSH02-95] - BSA and Azide free	HA750850	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
	alpha Tubulin Recombinant Rabbit mAb [PSH02-94] - BSA and Azide free	HA750849	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
	alpha Tubulin Rabbit pAb	HA500366	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	alpha Tubulin Rabbit pAb	R1308-14	Human, Rat, Mouse	WB,IF-Cell
KD	alpha Tubulin 4A Recombinant Rabbit mAb [JM73-24]	ET1705-31	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	alpha Tubulin (Acetyl K112) Rabbit pAb	ER60054	Human, Mouse, Rat	WB
Cell treatment	alpha Tubulin (acetyl K40) Recombinant Rabbit mAb [PSH0-13] - BSA and Azide free	HA750583	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	alpha Tubulin (acetyl K40) Recombinant Rabbit mAb [PSH19-53]	HA724079	Human, Mouse, Rat	WB,IP,IF-Cell
Cell treatment	alpha Tubulin (acetyl K40) Recombinant Rabbit mAb [PSH0-13]	HA721271	Human, Mouse, Rat	WB,IF-Cell
	Arp2 Recombinant Rabbit mAb [JE40-62]	HA721348	Human, Rat, Mouse	WB,IHC-P,IF-Cell
	Arp3 Recombinant Rabbit mAb [JB33-44]	ET7107-97	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell,IF-Tissue
	ATG4B Recombinant Rabbit mAb [JE59-97]	HA721334	Human, Mouse, Rat, Monkey	WB
KD	ATG4B Recombinant Rabbit mAb [PSH02-53] - BSA and Azide free	HA750816	Human, Mouse, Rat, Monkey	WB
	ATG4B Rabbit pAb	ER62759	Human, Mouse, Rat, Monkey	WB,IHC-P

KD	ATG4B Recombinant Rabbit mAb [PSH02-53]	HA721829	Human, Mouse, Rat, Monkey	WB
	beta II Tubulin Recombinant Rabbit mAb [ST52-04]	ET1609-48	Mouse, Rat, Human	WB,IF-Tissue,IHC-P
	HRP Conjugated beta Tubulin Recombinant Rabbit mAb [JF41-50]	ET1702-68	Human, Mouse, Rat	WB
	beta Tubulin Recombinant Rabbit mAb [SR25-04]	ET1602-4	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	beta Tubulin Recombinant Mouse mAb [1-B11-R] - BSA and Azide free	HA610049	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	iFluor™ 594 Conjugated beta Tubulin Recombinant Mouse mAb [A1-A4-R]	HA600110F	Human, Mouse	IF-Cell,IF-Tissue,FC
	beta Tubulin Recombinant Rabbit mAb [SR25-04] - BSA and Azide free	HA750049	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	HRP Conjugated beta Tubulin Mouse mAb [1-B11]	EM1701-59	Human, Mouse, Rat, Zebrafish	WB,IHC-P
	beta Tubulin Recombinant Mouse mAb [A1-A4-R]	HA601187	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IHC-P
	beta Tubulin Recombinant Mouse mAb [A1-A4-R] - BSA and Azide free	HA610067	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IHC-P
	beta Tubulin Mouse mAb [A1-A4]	M1305-2	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	beta Tubulin Mouse mAb [1-B11]	EM0103	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC,IF-Tissue
	beta Tubulin Recombinant Mouse mAb [1-B11-R]	HA601147	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	CD231 Rabbit pAb	ER65779	Human, Mouse, Rat	IHC-P
	CD34 Recombinant Rabbit mAb	IRS024	Human	mIHC
	CD34 Rabbit pAb	ER1706-66	Human, Mouse	WB,IHC-P
	CD34 Rabbit pAb	ER0802	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC
	CD34 Recombinant Rabbit mAb	IRS050	Mouse	mIHC
	CD34 Recombinant Rabbit mAb [SI16-01]	ET1606-11	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,mIHC,FC
	CD34 Recombinant Mouse mAb [PDM0-12]	HA601109	Human	IHC-P,IF-Tissue
	CD34 Mouse mAb [15H1]	EM1901-01	Human, Mouse, Rat	IHC-P,FC
	CD34 Rabbit pAb	0108-2	Human, Mouse, Zebrafish	WB,IF-Cell,IHC-P
	CD34 Recombinant Rabbit mAb [SI16-01] - BSA and Azide free	HA750092	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	CD5L Rabbit pAb	ER61441	Human, Mouse	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	CD63 Recombinant Rabbit mAb [PSH06-64]	HA722731	Human	WB,IF-Cell,FC,IP,IHC-P
KD	CD63 Recombinant Rabbit mAb [SY21-02] - BSA and Azide free	HA750112	Human	WB,IHC-P
KD	CD63 Recombinant Rabbit mAb [SY21-02]	ET1607-2	Human	WB,IHC-P
Relative expression	CD63 Recombinant Rabbit mAb [PSH06-64] - BSA and Azide free	HA751088	Human	WB,IF-Cell,FC,IP,IHC-P
Relative expression	CD81 Recombinant Rabbit mAb [PSH10-26]	HA723193	Human, Rat	WB,IHC-P,FC
Cell treatment	CD82 Recombinant Rabbit mAb [JE32-17]	HA722533	Human	WB
	CXCL12/SDF1 Rabbit pAb	ER1902-35	Human, Mouse, Rat	WB,IHC-P
	CXCL12/SDF1 Recombinant Rabbit mAb [JJ09-20]	ET1701-77	Human	WB,IF-Cell,IF-Tissue,IHC-P
	CXCL12/SDF1 Rabbit pAb	ER1802-6	Human	WB,IHC-P,FC
	DNM1L Mouse mAb [C7-A4]	EM1708-45	Human, Mouse, Rat	WB,IF-Cell,FC
	DRP1 Rabbit pAb	HA500487	Human, Mouse, Rat	WB,IF-Cell,FC,IHC-P
	Human CD63 Recombinant Rabbit mAb [PSH08-43] - BSA and Azide free (Capture)	HA722988	Human	ELISA(Cap)
	Human CD63 Recombinant Rabbit mAb [PSH08-44] - BSA and Azide free (Detector)	HA722989	Human	ELISA(Det)
	Biotin Conjugated Human CD63 Recombinant Rabbit mAb [PSH08-44] - Detector	HA722990B	Human	ELISA(Det),ELISA
	Human VCAM-1 Recombinant Rabbit mAb [PSH12-46] - BSA and Azide free (Capture)	HA725064	Human	ELISA(Cap)
	Human VCAM-1 Recombinant Rabbit mAb [PSH12-47] - BSA and Azide free (Detector)	HA725065	Human	ELISA(Det)
	Biotin Conjugated Human VCAM-1 Recombinant Rabbit mAb [PSH12-47] - Detector	HA725066B	Human	ELISA(Det),ELISA
	HRP Conjugated Human VCAM-1 Recombinant Rabbit mAb [PSH12-47] - Detector	HA725066H	Human	ELISA(Det),ELISA
	Human VEGFA Recombinant Rabbit mAb [PSH12-91] - BSA and Azide free (Capture)	HA723493	Human	ELISA(Cap)
	Human VEGFA Recombinant Rabbit mAb [PSH12-93] - BSA and Azide free (Detector)	HA723496	Human	ELISA(Det)
	Biotin Conjugated Human VEGFA Recombinant Rabbit mAb [PSH12-93] - Detector	HA723497B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human VEGFA Recombinant Rabbit mAb [PSH12-92] - Detector	HA723495B	Human	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human VEGFA Recombinant Rabbit mAb [PSH12-92] - BSA and Azide free (Detector)	HA723494	Human	ELISA(Det)
	ILK Rabbit pAb	RT1330	Human, Mouse, Rat	WB,IP,IF
KD	Integrin alpha 2 Recombinant Rabbit mAb [SN0752]	ET1611-57	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
KD	Integrin alpha 2 Recombinant Rabbit mAb [SN0752] - BSA and Azide free	HA750263	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	Integrin alpha 3 Recombinant Mouse mAb [A7H2-R]	HA601273	Human	WB,IF-Cell
	Integrin alpha 3 Recombinant Mouse mAb [A7H2-R] - BSA and Azide free	HA610148	Human	WB,IF-Cell
	Integrin alpha 3 Mouse mAb [A7H2]	HA600100	Human	WB,IF-Cell
	Integrin alpha 3 Rabbit pAb	HA500111	Human, Rat	WB,IHC-P
Cell treatment	Integrin alpha 4 Recombinant Rabbit mAb [JA09-36]	ET1704-31	Human, Mouse, Rat	WB,IHC-P,FC
	Integrin alpha 4 Rabbit pAb	ER65508	Human, Mouse, Rat	WB,IHC-P
KD	Integrin alpha 5 Recombinant Rabbit mAb [JJ08-94]	ET1701-58	Human, Mouse, Monkey, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Integrin alpha 6 Recombinant Rabbit mAb [SR45-00] - BSA and Azide free	HA750028	Human, Mouse, Rat	WB,IHC-P
	Integrin alpha 6 Recombinant Rabbit mAb [SR45-00]	ET1602-1	Human, Mouse, Rat	WB,IHC-P
	Integrin alpha 9 Recombinant Rabbit mAb [JE41-58]	HA722406	Human, Mouse, Rat	WB
	Integrin alpha D Rabbit pAb	ER65586	Human	WB,IHC-P
	Integrin alpha E Recombinant Rabbit mAb [SN62-06]	ET1611-27	Human	IHC-P
KD	Integrin alpha V Recombinant Rabbit mAb [SC56-07]	ET1610-15	Human, Mouse, Rat	WB,IF-Cell,FC,IHC-P
KD	Integrin alpha V Recombinant Rabbit mAb [SC56-07] - BSA and Azide free	HA750201	Human, Mouse, Rat	WB,IF-Cell,FC,IHC-P
	Integrin Alpha V Rabbit pAb	ER1911-52	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,IP
KD	Integrin beta 1 Rabbit pAb	ER31001	Human, Mouse	WB,IHC-P,FC
	Integrin beta 1 Mouse mAb [B10-A5]	EM1901-77	Human, Mouse	WB,IHC-P,FC
KD	Integrin beta 1 Recombinant Rabbit mAb [SA40-08] - BSA and Azide free	HA750008	Human, Mouse, Rat, Zebrafish	WB,IHC-P
KD	Integrin beta 1 Recombinant Rabbit mAb [SA40-08]	ET1601-17	Human, Mouse, Rat, Zebrafish	WB,IHC-P
KD	Integrin beta 1 Recombinant Rabbit mAb [SR30-03]	ET1602-10	Human, Mouse, Rat	WB,IHC-P,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	Integrin beta 3 Recombinant Rabbit mAb [PSH09-27]	HA723085	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Integrin beta 3 Rabbit pAb	ER1601-4	Human, Mouse, Rat	WB,IHC-P
KD	Integrin beta 3 Recombinant Rabbit mAb [SJ19-09]	ET1606-49	Human, Mouse	WB,IHC-P,IF-Tissue
	Integrin beta 3 Recombinant Rabbit mAb [JE22-64]	ET7110-61	Human, Mouse	WB,IHC-P
	Integrin beta 4 Recombinant Rabbit mAb [JM11-06]	ET1703-52	Human, Mouse, Rat	WB,IHC-P
	Integrin beta 5 Recombinant Rabbit mAb [PSH11-53]	HA723345	Human	IF-Cell,FC
	Integrin beta 5 Recombinant Rabbit mAb [PSH10-13]	HA723182	Human, Rat, Monkey	WB,IHC-P
	Integrin beta 5 Recombinant Rabbit mAb [PSH10-13] - BSA and Azide free	HA751342	Human, Rat, Monkey	WB,IHC-P
	Integrin beta 5 Recombinant Rabbit mAb [PSH11-53] - BSA and Azide free	HA751415	Human	IF-Cell,FC
	Integrin beta 5 Rabbit pAb	ER63768	Human, Mouse, Rat, Monkey	WB
	Integrin beta 7 Rabbit pAb	ER65852	Human, Mouse	WB
	Integrin linked ILK Rabbit pAb	ER1803-78	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	Integrin linked ILK Recombinant Rabbit mAb [SC68-04]	ET1610-76	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC
	KIF5B Recombinant Rabbit mAb [JE64-16]	HA721045	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Cell treatment	LAMP1 Recombinant Rabbit mAb [PSH07-39]	HA722827	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IHC-Fr,IP
Cell treatment	LAMP1 Recombinant Rabbit mAb [PSH05-55]	HA722302	Mouse	WB,IF-Cell,IHC-P,FC,IP,IF-Tissue,IHC-Fr
KD	LAMP1 Recombinant Rabbit mAb [JJ0940]	ET1701-94	Human	WB,IHC-P,IF-Tissue
KD	LAMP1 Recombinant Rabbit mAb [JJ0940] - BSA and Azide free	HA750328	Human	WB,IHC-P,IF-Tissue
Cell treatment	LAMP1 Recombinant Rabbit mAb [PSH07-39] - BSA and Azide free	HA751160	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IHC-Fr,IP
Cell treatment	LAMP1 Recombinant Rabbit mAb [PSH05-55] - BSA and Azide free	HA751017	Mouse	WB,IF-Cell,IHC-P,FC,IP,IF-Tissue,IHC-Fr
Cell treatment	LC3B Recombinant Antibody [JJ090-6] - Mouse IgG1 (Chimeric)	HA601532	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	LC3B Recombinant Antibody [JJ090-6] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610352	Human, Mouse, Rat	WB,IF-Cell,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment, KD	LC3B Recombinant Rabbit mAb [JJ090-6] - BSA and Azide free	HA750319	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,IP
Cell treatment, KD	LC3B Recombinant Rabbit mAb [JJ090-6]	ET1701-65	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,IP,mIHC
	L-Lactyl-alpha Tubulin (K40) Rabbit pAb	HA500554	Human, Mouse, Rat, Zebrafish	WB,IHC-P
	MLP3C Rabbit pAb	ER62516	Human	WB
	Biotin Conjugated Mouse VCAM1 Recombinant Rabbit mAb [PSH16-32] - Detector	HA723830B	Mouse	ELISA(Det),ELISA
	Mouse VCAM1 Recombinant Rabbit mAb [PSH16-32] - BSA and Azide free (Detector)	HA723830	Mouse	ELISA(Det)
	HRP Conjugated Mouse VCAM1 Recombinant Rabbit mAb [PSH16-32] - Detector	HA723830H	Mouse	ELISA(Det),ELISA
	Mouse VCAM1 Recombinant Rabbit mAb [PSH16-33] - BSA and Azide free (Capture)	HA723831	Mouse	ELISA(Cap)
	Mouse VCAM1 Recombinant Rabbit mAb [PSH16-31] - BSA and Azide free (Capture)	HA723829	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse VEGF Recombinant Rabbit mAb [PSH02-73]	HA722636B	Mouse	ELISA(Det)
	Mouse VEGF Recombinant Rabbit mAb [PSH02-73] - BSA and Azide free (Detector)	HA721851	Mouse	ELISA(Det)
	Mouse VEGF Recombinant Rabbit mAb [PSH02-72] - BSA and Azide free (Capture)	HA721850	Mouse	ELISA(Cap)
	MYO19 Recombinant Rabbit mAb [JE64-30]	HA721051	Human, Mouse, Rat	WB,IF-Cell
	NET-5 Rabbit pAb	ER64155	Human, Mouse	WB
Cell treatment	Phospho-DRP1 (S616) Recombinant Rabbit mAb [PSH06-77]	HA722760	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell,FC
Cell treatment	Phospho-DRP1 (S616) Recombinant Rabbit mAb [PSH06-77] - BSA and Azide free	HA751109	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell,FC
	PIP5K1A Rabbit pAb	HA500201	Human, Mouse, Rat	IHC-P,IF-Cell,FC
KD	Syntenin Recombinant Rabbit mAb [JE40-72]	ET7109-14	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
KD	TAPA1/CD81 Recombinant Rabbit mAb [SN206-01]	ET1611-87	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Tetraspanin-8 Rabbit pAb	ER65141	Human	WB,IHC-P,IF-Cell
Relative expression, KD	TOMM20 Recombinant Rabbit mAb [ST04-72]	ET1609-25	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell,FC,IP,IHC-Fr
	TOMM20 Recombinant Antibody [ST04-72] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610288	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression, KD	TOMM20 Recombinant Rabbit mAb [ST04-72] - BSA and Azide free	HA750166	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell,FC,IP,IHC-Fr
	TOMM20 Recombinant Antibody [ST04-72] - Mouse IgG1 (Chimeric)	HA601454	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	TSG101 Recombinant Rabbit mAb [JJ0900]	ET1701-59	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	TSG101 Recombinant Rabbit mAb [JJ0900] - BSA and Azide free	HA750316	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	TSPAN1 Rabbit pAb	ER1918-23	Human	WB
	Tubulin alpha 1/3/4 Rabbit pAb	ER65239	Human, Mouse, Rat	WB
	Tubulin alpha 3C/D/E Rabbit pAb	ER65240	Human, Mouse, Rat	WB
	VCAM1 Rabbit pAb	R1512-13	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Relative expression, KD	VCAM1 Recombinant Rabbit mAb [SA05-04]	ET1601-18	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Relative expression, KD	VCAM1 Recombinant Rabbit mAb [SA05-04] - BSA and Azide free	HA750009	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	VEGF Rabbit pAb	ER30607	Human	WB,IF-Cell,IHC-P,FC
	VEGF Recombinant Rabbit mAb [SP07-01]	ET1604-28	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	VEGF Recombinant Rabbit mAb [SP07-01] - BSA and Azide free	HA750074	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	VEGFA Recombinant Rabbit mAb [PSH19-17]	HA724051	Human, Mouse, Rat, Monkey	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层