

乳腺癌研究-产品清单

高级验证	产品名	货号	反应种属	验证应用
	Actin alpha 3 Rabbit pAb	ER62586	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	AKT1 Mouse mAb [D9-9-C9]	EM40507	Human, Mouse, Rat	WB,IHC-P
KD	AKT1 Recombinant Rabbit mAb [ST05-09] - BSA and Azide free	HA750177	Human, Mouse, Rat	WB,IHC-P,IP,FC,IHC-Fr,IF-Cell
	AKT1 Recombinant Rabbit mAb [PSH0-50]	HA721364	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC
	AKT1 Recombinant Rabbit mAb [PSH0-50] - BSA and Azide free	HA750627	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC
KD	AKT1 Recombinant Rabbit mAb [ST05-09]	ET1609-47	Human, Mouse, Rat	WB,IHC-P,IP,FC,IHC-Fr,IF-Cell
KD	AKT1/2/3 Recombinant Rabbit mAb [JE75-09]	HA721870	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IP
KD	AKT1/2/3 Recombinant Rabbit mAb [JE75-09] - BSA and Azide free	HA750840	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IP
	AKT1/2/3 Rabbit pAb	ER1901-33	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	AKT1/2/3 Recombinant Rabbit mAb [ST48-09] - BSA and Azide free	HA750179	Human, Mouse, Rat, Monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
	AKT1/2/3 Recombinant Rabbit mAb [ST48-09]	ET1609-51	Human, Mouse, Rat, Monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
	Akt1/3 Rabbit pAb	ER62638	Human, Mouse, Rat	WB,IHC-P
	Akt2 Mouse mAb [2A1]	RT1027	Human, Mouse, Rat	WB,IP,IF,IHC-P
	AKT2 Rabbit pAb	HA500091	Human, Mouse, Rat	WB,IHC-P
	AP-1/Jun D Rabbit pAb	ER62677	Human, Mouse, Rat	WB,IHC-P
	APC Rabbit pAb	ER1802-39	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue,IHC-Fr
	APC Recombinant Rabbit mAb [SC61-03]	ET1610-80	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IP,IHC-P,IHC-Fr
KD	A-RAF Recombinant Rabbit mAb [JJ08-68]	ET1701-18	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Axin2 Recombinant Rabbit mAb [JM11-30]	ET1703-96	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Bak Rabbit pAb	ER1904-10	Human, Mouse	WB,IHC-P,IF-Tissue,FC
KO	Bax Recombinant Rabbit mAb [SZ3-07] - BSA and Azide free	HA750065	Human, Mouse, Rat, Goat	WB,IF-Cell,FC
KO	Bax Recombinant Rabbit mAb [SZ3-07]	ET1603-34	Human, Mouse, Rat, Goat	WB,IF-Cell,FC
KD	BAX Rabbit pAb	ER0907	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	BAX Recombinant Rabbit mAb [PD01-44]	HA721324	Human	IHC-P,IF-Tissue
	BAX Recombinant Mouse mAb [PSH0-20]	HA601123	Human	WB,IHC-P
	BAX Recombinant Mouse mAb [PSH0-20] - BSA and Azide free	HA610038	Human	WB,IHC-P
	BAX Mouse mAb [3-D2]	EM1203	Human, Mouse	WB,IF-Cell,IHC-P,FC
	Beta Catenin Rabbit pAb	ER0805	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell,IF-Tissue
	Beta Catenin Recombinant Mouse mAb [A6-F8-R] - BSA and Azide free	HA610062	Human, Mouse, Rat	WB,IF-Cell,FC
	Beta Catenin Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS277MS	Human	mIHC
	Beta Catenin Recombinant Rabbit mAb	IRS043	Human	mIHC
KD	Beta Catenin Recombinant Rabbit mAb [SA30-04]	ET1601-5	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IP,mIHC,IF-Cell,IHC-Fr,FC
KD	Beta Catenin Recombinant Rabbit mAb [SA30-04] - BSA and Azide free	HA750023	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IP,IF-Cell,IHC-Fr,FC
	Beta Catenin Recombinant Mouse mAb [A6-F8-R]	HA601179	Human, Mouse, Rat	WB,IF-Cell,FC
	Beta Catenin Recombinant Rabbit mAb	IRS061	Mouse	mIHC
	Beta Catenin Mouse mAb [10-C0-B7]	M1405-6	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Beta Catenin Rabbit pAb	0407-16	Human, Mouse, Rat, Zebrafish	WB,IF-Cell
	iFluor™ 488 Conjugated Beta Catenin Recombinant Rabbit mAb [SA30-04]	HA720182F	Human, Mouse	IF,ICC
	Beta Catenin Mouse mAb [A6-F8]	EM0306	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue,ChIP
	Beta Catenin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS278RT	Human	mIHC
	bFGF Rabbit pAb	ER1904-22	Human, Mouse	WB,FC
KD	BRAF Recombinant Rabbit mAb [SU34-04]	ET1608-36	Human, Mouse, Rat	WB,IHC-P,FC,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	BRAF Mouse mAb [F1-F4-E10]	EM1701-32	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
	BRCA1 Rabbit pAb	HA500015	Human, Mouse, Rat	WB,IHC-P
	BRCA1 Mouse mAb [A7E3]	HA601009	Human, Mouse, Rat	IHC-P
	BRCA1 Recombinant Mouse mAb [A7E3-R] - BSA and Azide free	HA610177	Human, Mouse, Rat	IHC-P,IF-Cell
	BRCA1 Mouse mAb [D2-C6]	EM1710-25	Human	IHC-P
	BRCA1 Recombinant Mouse mAb [A7E3-R]	HA601304	Human, Mouse, Rat	IHC-P,IF-Cell
	BRCA2 Recombinant Rabbit mAb [PSH07-55] - BSA and Azide free	HA751177	Human	WB,FC
	BRCA2 Recombinant Rabbit mAb [PSH07-55]	HA722850	Human	WB,FC
	Calponin Recombinant Rabbit mAb [SI67-01] - BSA and Azide free	HA750094	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	Calponin Recombinant Rabbit mAb [SI67-01]	ET1606-17	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	Calponin-1 Rabbit pAb	R1510-7	Human, Mouse, Rat	IF-Cell,IHC-P,FC
	Casein Kinase 1 alpha Recombinant Rabbit mAb [PSH04-46] - BSA and Azide free	HA750955	Human, Mouse, Rat	WB,IF-Cell,FC,IP
	Casein Kinase 1 alpha Recombinant Rabbit mAb [PSH04-46]	HA722131	Human, Mouse, Rat	WB,IF-Cell,FC,IP
	Casein Kinase I alpha Rabbit pAb	ER62901	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Casein Kinase I alpha 1/1L Rabbit pAb	ER62900	Human, Mouse, Rat	IHC-P,IF-Cell,WB
	Caveolin-1 Recombinant Rabbit mAb	IRS100RB	Mouse	mIHC
	Caveolin-1 Mouse mAb [A3E4]	EM1902-01	Human, Mouse	WB,IHC-P,FC
	Caveolin-1 Rabbit pAb	0407-6	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Caveolin-1 Mouse mAb [19-D6]	EM40728	Human, Mouse, Rat	WB,IHC-P
	Caveolin-1 Mouse mAb [A3E6]	HA601005	Human	WB,IHC-P
	Caveolin-1 Recombinant Rabbit mAb [SZ02-01]	ET1603-1	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	Caveolin-1 Mouse mAb [A3E7]	EM1901-87	Human	WB,IF-Cell,IHC-P,FC
	Caveolin-1 Mouse mAb [A3E8]	EM1901-88	Human	WB,IF-Cell,IHC-P,FC
	Cdk4 Recombinant Mouse mAb [5-C3-G3-R]	HA601227	Human, Mouse, Rat	WB
	Cdk4 Recombinant Rabbit mAb [SD205-1]	ET1612-23	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Cdk4 Recombinant Mouse mAb [5-C3-G3-R] - BSA and Azide free	HA610102	Human, Mouse, Rat	WB
	Cdk4 Recombinant Rabbit mAb [SD20-42]	ET1612-1	Human	WB,IF-Tissue,IHC-P
	Cdk4 Recombinant Rabbit mAb [SD205-1] - BSA and Azide free	HA750286	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Cdk4 Mouse mAb [5-C3-G3]	M1505-5	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	CDKN1A/P21 Rabbit pAb	ER1906-07	Human, Rat	IHC-P,IF-Cell,IF-Tissue,FC
Cell treatment	CDKN1A/P21 Rabbit pAb	HA500156	Human, Mouse	WB,IHC-P
Cell treatment	c-Fos Recombinant Rabbit mAb [PSH05-77]	HA722666	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Cell,IHC-Fr,mIHC,IF-Tissue
Cell treatment	c-Fos Recombinant Antibody [PSH05-77] - Rat IgG1 (Chimeric)	HA601373	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,IF-Tissue,IF-Cell,WB
	c-Fos Mouse mAb [A10-G2]	EM1710-29	Human	WB,IHC-P,FC
Cell treatment	c-Fos Recombinant Rabbit mAb [JJ0938]	ET1701-95	Human, Rat	WB,IF-Cell,IHC-P
Cell treatment	c-Fos Recombinant Antibody [PSH07-51] - Mouse IgG1 (Chimeric)	HA601440	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,IF-Cell,WB
	c-Fos Recombinant Rabbit mAb [PSH07-51] - BSA and Azide free	HA751172	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-Fr,IF-Cell,IHC-P,ChIP,IF-Tissue,IP
Cell treatment	c-Fos Recombinant Antibody [PSH05-77] - Guinea pig IgG2 (Chimeric)	HA601374	Mouse, Rat, Human, Cynomolgus monkey, Pig	IHC-Fr,WB
Cell treatment	c-Fos Recombinant Antibody [PSH07-51] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610275	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,IF-Cell,WB
Cell treatment	c-Fos Recombinant Rabbit mAb [JJ0938] - BSA and Azide free	HA750329	Human, Rat	WB,IF-Cell,IHC-P
Cell treatment	c-Fos Recombinant Antibody [PSH05-77] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610220	Mouse, Rat, Human, Cynomolgus monkey, Pig	IHC-Fr,WB
	c-Fos Recombinant Rabbit mAb	IRS131RB	Mouse	mIHC
Cell treatment	c-Fos Recombinant Antibody [PSH05-77] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610219	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,IF-Tissue,IF-Cell,WB
Cell treatment	c-Fos Recombinant Rabbit mAb [PSH05-77] - BSA and Azide free	HA751071	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Cell,IHC-Fr,IF-Tissue
	c-Fos Recombinant Rabbit mAb [PSH07-51]	HA722839	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-Fr,IF-Cell,IHC-P,ChIP,IF-Tissue,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	c-Jun Recombinant Rabbit mAb [SY0298] - BSA and Azide free	HA750143	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
	c-Jun Recombinant Rabbit mAb [SY0298]	ET1608-3	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
	c-kit Rabbit pAb	0108-4	Human	WB,IF-Cell
Relative expression, KD	c-Kit Recombinant Rabbit mAb [ST04-99]	ET1609-60	Human	WB,IHC-P
	c-Kit Recombinant Rabbit mAb	IRS013	Human	mIHC
Relative expression, KD	c-Kit Recombinant Rabbit mAb [PD00-24]	HA721154	Human	IHC-P,WB,mIHC
	Cleaved-Notch 2 (A1734) Rabbit pAb	ER60031	Human, Mouse, Rat	WB
	Cleaved-Notch 2 (D1733) Rabbit pAb	ER60032	Human, Mouse, Rat	IF-Cell,WB,IHC-P
	Cleaved-Notch 2 (V1697) Rabbit pAb	ER60033	Human, Mouse, Rat	WB
	c-Myc Recombinant Rabbit mAb [JE01-20]	HA722895	Human, Mouse, Rat, Monkey	WB,IF-Cell
	HRP Conjugated c-Myc Rabbit pAb	0912-4	Human	WB
	c-Myc Recombinant Rabbit mAb [PS00-11]	HA721182	Human, Mouse, Rat	WB,IHC-P,IP,ChIP
	c-Myc Rabbit pAb	0912-2	Human	WB
	c-Myc Mouse mAb [6G6]	RT1149	Human, Mouse, Rat, Monkey	WB,IP,IF,IHC-P
	c-Myc Recombinant Rabbit mAb [PS00-11] - BSA and Azide free	HA750519	Human, Mouse, Rat	WB,IHC-P,IP,ChIP
	COX2 Recombinant Rabbit mAb [6J8]	RE6064	Human, Mouse	WB,ICC
	COX2 Mouse mAb [A3F8]	EM1902-11	Human, Rat	WB,IHC-P,FC
KD	COX2 Mouse mAb [A3F7]	EM1902-12	Human, Rat	WB,IHC-P,FC
	COX2 Mouse mAb	RT1159	Human, Mouse, Rat	WB,IP,IF,IHC-P,ELISA
KD	COX2/Cyclooxygenase 2 Recombinant Rabbit mAb [SC56-06]	ET1610-23	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
KD	COX2/Cyclooxygenase 2 Recombinant Rabbit mAb [SC56-06] - BSA and Azide free	HA750207	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	CSNK1A1 + CSNK1A1L Recombinant Rabbit mAb [PSH04-58]	HA722157	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
	CSNK1A1 + CSNK1A1L Recombinant Rabbit mAb [PSH04-58] - BSA and Azide free	HA750965	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
Relative expression, KD	Cyclin D1 Recombinant Rabbit mAb [SA38-08]	ET1601-31	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Cyclin D1 Rabbit pAb	ER0722	Human, Mouse, Rat	WB,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression, KD	Cyclin D1 Recombinant Rabbit mAb [PD01-64]	HA721322	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	Cyclin D1 Rabbit pAb	0407-27	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,FC
	iFluor™ 488 Conjugated Cytokeratin 10 Recombinant Rabbit mAb [ST05-43]	HA720134F	Human, Rat	IF-Tissue
KD	Cytokeratin 10 Recombinant Rabbit mAb [ST05-43]	ET1609-75	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	iFluor™ 647 Conjugated Cytokeratin 10 Recombinant Rabbit mAb [ST05-43]	HA720146F	Human, Rat, Mouse	IF-Tissue
	iFluor™ 594 Conjugated Cytokeratin 10 Recombinant Rabbit mAb [ST05-43]	HA720119F	Human, Rat	IF-Tissue
KD	Cytokeratin 10 Recombinant Rabbit mAb [ST05-43] - BSA and Azide free	HA750193	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	Cytokeratin 14 Mouse mAb [A2C12]	EM1901-31	Human	WB,IHC-P,FC
Relative expression	Cytokeratin 14 Recombinant Mouse mAb [A2C11-R]	HA601259	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Cytokeratin 14 Rabbit pAb	ER1901-08	Human, Rat	WB,IHC-P
	iFluor™ 594 Conjugated Cytokeratin 14 Recombinant Rabbit mAb [SC65-06]	HA720117F	Rat, Human, Mouse	IF-Tissue
Relative expression	Cytokeratin 14 Recombinant Mouse mAb [A2C10-R]	HA601323	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Relative expression	Cytokeratin 14 Mouse mAb [A2C10]	EM1901-33	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue
	iFluor™ 647 Conjugated Cytokeratin 14 Recombinant Rabbit mAb [SC65-06]	HA720147F	Human, Mouse, Rat	IF-Cell,IF-Tissue
Relative expression	Cytokeratin 14 Recombinant Rabbit mAb [SC65-06]	ET1610-42	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	iFluor™ 488 Conjugated Cytokeratin 14 Recombinant Rabbit mAb [SC65-06]	HA720135F	Mouse, Rat, Human	IF-Tissue,IF-Cell
Relative expression	Cytokeratin 14 Recombinant Rabbit mAb [SC65-06] - BSA and Azide free	HA750216	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Cytokeratin 14 Recombinant Rabbit mAb	IRS162RB	Human	mIHC
Relative expression	Cytokeratin 14 Mouse mAb [A2C11]	EM1901-32	Human	WB,IHC-P,IF-Cell
	iFluor™ 594 Conjugated Cytokeratin 15 Recombinant Rabbit mAb [ST04-85]	HA720122F	Human, Mouse, Rat	IF-Tissue
	Cytokeratin 15 Recombinant Rabbit mAb [ST04-85]	ET1609-54	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell
	Cytokeratin 15 Recombinant Rabbit mAb [ST04-85] - BSA and Azide free	HA750182	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell
	iFluor™ 647 Conjugated Cytokeratin 15 Recombinant Rabbit mAb [ST04-85]	HA720148F	Human, Mouse, Rat	IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	iFluor™ 488 Conjugated Cytokeratin 15 Recombinant Rabbit mAb [ST04-85]	HA720136F	Human, Mouse, Rat	IF-Tissue,IF-Cell,FC
	iFluor™ 488 Conjugated Cytokeratin 16 Recombinant Rabbit mAb [SC52-09]	HA720137F	Human	IF-Tissue,IF-Cell,FC
	Cytokeratin 16 Recombinant Rabbit mAb [SC52-09] - BSA and Azide free	HA750203	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
	iFluor™ 594 Conjugated Cytokeratin 16 Recombinant Rabbit mAb [SC52-09]	HA720114F	Human	IF-Cell,IF-Tissue
	iFluor™ 647 Conjugated Cytokeratin 16 Recombinant Rabbit mAb [SC52-09]	HA720149F	Human	IF-Tissue
	Cytokeratin 16 Recombinant Rabbit mAb [SC52-09]	ET1610-17	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Cytokeratin 19 Recombinant Antibody [SA30-06] - Rat IgG1 (Chimeric)	HA601555	Human, Rat	WB,IF-Cell,IHC-P,FC
	Cytokeratin 19 Recombinant Antibody [SA30-06] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610370	Human, Rat	WB,IF-Cell,IHC-P,FC
KD	Cytokeratin 19 Mouse mAb [A3D1]	EM1901-75	Human	WB,IHC-P,IF-Cell,FC
	iFluor™ 647 Conjugated Cytokeratin 19 Recombinant Rabbit mAb [SA30-06]	HA720152F	Human, Mouse	IF-Cell,IF-Tissue,FC
	iFluor™ 488 Conjugated Cytokeratin 19 Recombinant Rabbit mAb [SA30-06]	HA720140F	Human, Mouse	IF-Cell,IF-Tissue,FC
KO	Cytokeratin 19 Recombinant Rabbit mAb [SA30-06]	ET1601-6	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,mIHC
KD	Cytokeratin 19 Recombinant Mouse mAb [A3D1-R] - BSA and Azide free	HA610133	Human	WB,IF-Cell,IHC-P,FC
	Cytokeratin 19 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS134RT	Human	mIHC
KD	Cytokeratin 19 Mouse mAb [A3D5]	EM1901-76	Human	WB,IHC-P,IF-Cell,FC
	iFluor™ 594 Conjugated Cytokeratin 19 Recombinant Rabbit mAb [SA30-06]	HA720115F	Human, Mouse	IF-Cell,IF-Tissue,FC
	Cytokeratin 19 Rabbit pAb	ER1803-79	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Cytokeratin 19 Recombinant Rabbit mAb	IRS011	Human	mIHC
KD	Cytokeratin 19 Recombinant Mouse mAb [A3D1-R]	HA601258	Human	WB,IF-Cell,IHC-P,FC
	Cytokeratin 20 Rabbit pAb	ER1901-76	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	Cytokeratin 20 Recombinant Rabbit mAb [SA35-03] - BSA and Azide free	HA750026	Human, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	Cytokeratin 20 Recombinant mAb	IRS135RB	Human	mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Cytokeratin 20 Rabbit pAb	ER1901-97	Human, Mouse	WB,IHC-P,FC
	Cytokeratin 20 Mouse mAb [A3D12]	EM1901-96	Human	WB,FC,IHC-P
	iFluor™ 488 Conjugated Cytokeratin 20 Recombinant Rabbit mAb [SA35-03]	HA720189F	Human	IF-Cell,IF-Tissue
	Cytokeratin 20 Recombinant Rabbit mAb [JE20-39]	ET7110-54	Human, Mouse	WB,IHC-P
	Cytokeratin 20 Recombinant Mouse mAb [A3H7-R]	HA601077	Human	IHC-P,IF-Tissue
	Cytokeratin 20 Mouse mAb [A3H7]	EM1901-97	Human	IHC-P,IF-Tissue
	Cytokeratin 20 Recombinant Rabbit mAb [SA35-03]	ET1601-8	Human, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	iFluor™ 594 Conjugated Cytokeratin 5 Recombinant Rabbit mAb [SC62-04]	HA720121F	Human	IF-Tissue
	iFluor™ 488 Conjugated Cytokeratin 5 Recombinant Rabbit mAb [SC62-04]	HA720129F	Human	IF-Tissue,IF-Cell
	Cytokeratin 5 Recombinant Rabbit mAb	IRS249RB	Human	mIHC
	iFluor™ 647 Conjugated Cytokeratin 5 Recombinant Rabbit mAb [SC62-04]	HA720126F	Human	IF-Cell,IF-Tissue
Relative expression, KD	Cytokeratin 5 Recombinant Rabbit mAb [SC62-04]	ET1610-43	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell,FC
	Cytokeratin 5 Mouse mAb [6D-4D]	EM1706-36	Human, Mouse, Monkey	WB,IHC-P,IF-Cell,FC
	Cytokeratin 5 Recombinant Rabbit mAb [SY88-19]	ET1605-17	Human, Mouse	WB,IHC-P
	Cytokeratin 5+6 Recombinant Rabbit mAb [PSH0-68]	HA721456	Human	WB,IF-Cell,FC
	Cytokeratin 5+6 Recombinant Mouse mAb [A2B2-R] - BSA and Azide free	HA610108	Human, Rat	WB,IHC-P
	Cytokeratin 5+6 Rabbit pAb	ER1901-03	Human, Mouse, Rat	WB,IHC-P
	Cytokeratin 5+6 Mouse mAb [A2B2]	EM1901-08	Human, Rat	WB,IHC-P,FC
	Cytokeratin 5+6 Rabbit pAb	HA500508	Human	WB,IHC-P
	Cytokeratin 5+6 Recombinant Mouse mAb [A2B2-R]	HA601233	Human, Rat	WB,IHC-P
	Cytokeratin 5+6 Recombinant Rabbit mAb [PSH0-67]	HA721455	Human, Mouse, Rat	IHC-P,IF-Tissue,WB
	Cytokeratin 5+6 Rabbit pAb	ER1901-86	Human, Mouse, Rat	WB,IHC-P,FC
	Cytokeratin 5+6 Recombinant Rabbit mAb [PSH0-68] - BSA and Azide free	HA750645	Human	WB,IF-Cell,FC
	Cytokeratin 5+6 Mouse mAb [A2A12]	EM1901-07	Human, Rat	WB,IHC-P,FC
	Cytokeratin 6 Recombinant Rabbit mAb [SN71-07]	ET1611-70	Human, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	iFluor™ 647 Conjugated Cytokeratin 6 Recombinant Rabbit mAb [SN71-07]	HA720143F	Human, Rat	IF-Cell,IF-Tissue
	Cytokeratin 6 Recombinant Rabbit mAb [SN71-07] - BSA and Azide free	HA750270	Human, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	iFluor™ 488 Conjugated Cytokeratin 6 Recombinant Rabbit mAb [SN71-07]	HA720130F	Human, Rat	IF-Cell
	Cytokeratin 7 Rabbit pAb	ER1803-89	Human, Mouse, Rat	WB,IHC-P
KD, Relative expression	Cytokeratin 7 Recombinant Rabbit mAb [ST50-05]	ET1609-62	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,mIHC
KD	Cytokeratin 7 Rabbit pAb	R1309-4	Human	WB,IF-Cell,IHC-P,FC
KO	Cytokeratin 7 Mouse mAb [3-G3-D8]	EM0702	Human	WB,IHC-P,IF-Cell,IF-Tissue
KD, Relative expression	Cytokeratin 7 Recombinant Rabbit mAb [ST50-05] - BSA and Azide free	HA750185	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Cytokeratin 7 Recombinant Mouse mAb [3-G3-D8-R]	HA601065	Human	WB,IHC-P,IF-Cell,FC,IF-Tissue
	iFluor™ 594 Conjugated Cytokeratin 7 Recombinant Rabbit mAb [ST50-05]	HA720120F	Human, Mouse	IF-Cell,IF-Tissue,FC
	iFluor™ 647 Conjugated Cytokeratin 7 Recombinant Rabbit mAb [ST50-05]	HA720144F	Human, Mouse	IF-Cell,IF-Tissue,FC
	Cytokeratin 7 Recombinant Rabbit mAb	IRS029	Human	mIHC
	iFluor™ 488 Conjugated Cytokeratin 7 Recombinant Rabbit mAb [ST50-05]	HA720131F	Human, Mouse	IF-Tissue,IF-Cell,FC
	Cytokeratin (Pan) Mouse mAb [7H8C4]	EM1712-42	Human	WB,IHC,IF-Cell,FC
	delta 1 Catenin/CAS Rabbit pAb	ER1803-71	Human, Rat	WB,IF-Cell
	delta 1 Catenin/CAS Rabbit pAb	ER1901-07	Human, Mouse	WB,IF-Cell,IHC-P
	delta 1 Catenin/CAS Recombinant Rabbit mAb [PD00-43]	HA721193	Human	WB,IHC-P,IF-Tissue
	delta 1 Catenin/CAS Rabbit pAb	ER1901-80	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
KD	delta 1 Catenin/CAS Rabbit pAb	ER1901-92	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Dishevelled 2 Recombinant Rabbit mAb [PSH04-07]	HA722088	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Dishevelled 2 Recombinant Rabbit mAb [PSH04-07] - BSA and Azide free	HA750930	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	DLL1 Rabbit pAb	ER61747	Human, Rat, Mouse	WB
	DLL3 Rabbit pAb	ER1907-51	Human, Mouse	WB
	DLL4 Rabbit pAb	ER1706-29	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	DLL4 Rabbit pAb	ER1902-20	Human, Mouse	WB,FC
	DLL4 Rabbit pAb	ER1706-63	Human, Rat	WB,IF-Cell,IHC-P,FC
	DNA Ligase I Recombinant Rabbit mAb [JB43-39]	ET7107-14	Human, Rat	WB,IHC-P
	Dvl-2 Rabbit pAb	ER63295	Human, Mouse	WB,IHC-P
	E2F1 Recombinant Rabbit mAb [JJ092-02]	ET1701-73	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	E2F1 Recombinant Rabbit mAb [JJ092-02] - BSA and Azide free	HA750322	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	E2F2 Recombinant Rabbit mAb [SN201-04]	ET1611-88	Human, Mouse, Rat	WB,FC
	iFluor™ 488 Conjugated E-Cadherin Recombinant Rabbit mAb [SY0287]	HA720159F	Human	IF-Cell,FC
	E-Cadherin Recombinant Mouse mAb	IRS229MS	Mouse	mIHC
	E-Cadherin Rabbit pAb	0407-25	Human, Mouse	WB,IF-Cell,FC
	E-Cadherin Recombinant Mouse mAb [A0-G11-2-R] - BSA and Azide free	HA610047	Human, Mouse, Rat	WB,IHC
KD	E-Cadherin Mouse mAb [A0-G11-2]	EM0502	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Relative expression, KD	E-Cadherin Recombinant Rabbit mAb [SY0287] - BSA and Azide free	HA750133	Human, Cynomolgus monkey	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	E-Cadherin Recombinant Rabbit mAb	IRS030	Human	mIHC
	E-Cadherin Recombinant Mouse mAb [A0-G11-2-R]	HA601143	Human, Mouse, Rat	WB,IHC,mIHC
	E-Cadherin Recombinant Mouse mAb	IRS031	Human	mIHC
Relative expression	E-Cadherin Recombinant Rabbit mAb [PSH13-75] - BSA and Azide free	HA751488	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Cell,FC,IP
	iFluor™ 594 Conjugated E-Cadherin Recombinant Rabbit mAb [SY0287]	HA720164F	Human	IF-Tissue
Relative expression, KD	E-Cadherin Recombinant Rabbit mAb [SY0287]	ET1607-75	Human, Cynomolgus monkey	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP,mIHC
Relative expression	E-Cadherin Recombinant Rabbit mAb [PSH13-75]	HA723564	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Cell,FC,IP
	E-Cadherin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS284RT	Mouse	mIHC
Relative expression	E-Cadherin Rabbit pAb	ER63312	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IHC-Fr,FC
	E-Cadherin Mouse mAb [3-F9]	M1405-3	Human, Mouse	IHC-P
	EGF Recombinant Rabbit mAb [JA14-15]	ET1704-33	Human, Mouse, Rat	WB,IHC-P,IP
	EGFR Mouse mAb [B6-E5-D9]	EM1701-22	Human, Mouse, Monkey, Rat	WB,IHC-P,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression, KD	EGFR Recombinant Rabbit mAb [PD00-75]	HA721220	Human, Monkey	WB,IHC-P,IF-Cell,FC
Relative expression, KD	EGFR Mouse mAb [A8E12]	HA601073	Human, Monkey	WB,IHC-P
	EGFR Recombinant Rabbit mAb	IRS046	Human	mIHC
Relative expression, KO	EGFR Mouse mAb [7-F2-F8]	EM1901-67	Human, Monkey	WB,IHC-P,FC
Relative expression, KO	EGFR Recombinant Rabbit mAb [SP00-86]	ET1604-44	Human, Mouse, Rat, Monkey	WB,IF-Cell,IP,FC
	EGFR Rabbit pAb	R1608-3	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC
Relative expression, KD	EGFR Recombinant Rabbit mAb [SZ40-19]	ET1603-37	Human, Mouse, Monkey, Rat	WB,IF-Cell,IP,FC,mIHC
Relative expression, KO	EGFR Recombinant Rabbit mAb [SP00-86] - BSA and Azide free	HA750078	Human, Mouse, Rat, Monkey	WB,IF-Cell,IP,FC
	EGFR Rabbit pAb	ER1512-6	Human, Mouse, Monkey	WB,IF-Cell,IHC-P,FC
Relative expression, KD	EGFR Recombinant Rabbit mAb [PD00-75] - BSA and Azide free	HA750550	Human, Monkey	WB,IHC-P,IF-Cell,FC
	EGFR Rabbit pAb	0407-21	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P
	EGFR Rabbit pAb	R1511-18	Human, Monkey	WB,IF-Cell,IHC-P
KD	ERK1 Recombinant Rabbit mAb [SP05-09]	ET1604-16	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	ERK1 Recombinant Rabbit mAb [SP00-74]	ET1604-32	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	ERK1 Mouse mAb [7-D6-E5]	M1505-9	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Erk1/2 Mouse mAb [2G1]	RT1484	Human, Mouse, Rat	WB,IP,IF
KD	ERK1/2 Recombinant Rabbit mAb [SA43-03] - BSA and Azide free	HA750019	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IP,FC,IHC-P,IHC-Fr
KD	ERK1/2 Recombinant Rabbit mAb [SA43-03]	ET1601-29	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IP,FC,IHC-P,IHC-Fr
	ERK1/2 Rabbit pAb	ER131011	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	ERK2 Mouse mAb [C3B12]	EM1901-55	Human, Mouse, Rat	WB,IHC-P,FC
KD	ERK2 Recombinant Mouse mAb [C3B12-R]	HA601201	Human, Mouse, Rat	WB
KD	ERK2 Recombinant Mouse mAb [5-D2-R] - BSA and Azide free	HA610080	Human, Mouse, Rat	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	ERK2 Mouse mAb [A10D2]	EM1901-54	Human, Mouse, Rat	WB,IHC-P,FC
KD	ERK2 Rabbit pAb	ER131218	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	ERK2 Mouse mAb [5-D2]	M1407-3	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	ERK2 Recombinant Mouse mAb [5-D2-R]	HA601202	Human, Mouse, Rat	WB
KD	ERK2 Recombinant Mouse mAb [C3B12-R] - BSA and Azide free	HA610079	Human, Mouse, Rat	WB
KD	ERK2 Mouse mAb [B10-8]	EM1901-53	Human, Mouse, Rat	WB,IHC-P,FC
KD	ERK2 Recombinant Rabbit mAb [SZ25-01]	ET1603-23	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
	ESR1 Rabbit pAb	0803-1	Human	WB,IF-Cell,IHC-P
	ESR1 Mouse mAb [E8-A6]	EM1708-51	Human, Mouse, Rat, Rabbit	WB,IF-Cell,IHC-P,FC
Relative expression	Estrogen Receptor alpha Recombinant Rabbit mAb [JE26-14]	ET7110-60	Human, Mouse	WB,IHC-P
	Estrogen Receptor alpha Recombinant mAb	IRS142RB	Human	mIHC
	Estrogen Receptor alpha Rabbit pAb	ER1803-83	Human, Mouse, Rat	WB,FC
	Estrogen Receptor alpha Recombinant Rabbit mAb [PD00-04]	HA721140	Human	IHC-P,IF-Tissue
Relative expression	Estrogen Receptor beta Recombinant Rabbit mAb [JJ09-02]	ET1701-69	Human	WB,IHC-P
	FGF-1 Rabbit pAb	ER63425	Human, Mouse, Rat	IHC-P
	FGF1 Recombinant Rabbit mAb [PSH18-61]	HA724014	Human, Mouse, Rat	WB,FC
	FGF1 Recombinant Rabbit mAb [JE33-93]	HA722695	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	FGF10 Rabbit pAb	ER60274	Human, Mouse, Rat	WB
	FGF-16 Rabbit pAb	ER65578	Human, Mouse, Rat	WB,IHC-P
	FGF-17 Rabbit pAb	ER65622	Human, Mouse, Rat	WB
	FGF-18 Rabbit pAb	ER60001	Human, Mouse, Rat	WB,IHC-P
Cell treatment	FGF18 Recombinant Rabbit mAb [PSH01-25]	HA721642	Mouse	WB
Cell treatment	FGF18 Recombinant Rabbit mAb [PSH01-24]	HA721641	Human, Mouse	WB
Cell treatment	FGF18 Recombinant Rabbit mAb [PSH01-25] - BSA and Azide free	HA750701	Mouse	WB
Cell treatment	FGF18 Recombinant Rabbit mAb [PSH01-24] - BSA and Azide free	HA750700	Human, Mouse	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	FGF-19 Rabbit pAb	ER65959	Human, Pig	WB
	FGF2 Recombinant Rabbit mAb [JM28-10]	ET1703-18	Human, Mouse	WB,IP,FC,IHC-P
	FGF-20 Rabbit pAb	ER65477	Human, Mouse, Rat	WB,IHC-P
	FGF21 Recombinant Rabbit mAb [JA10-97]	ET1704-04	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	FGF21 Recombinant Rabbit mAb [JA10-97] - BSA and Azide free	HA750399	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	FGF-22 Rabbit pAb	ER63427	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	FGF23 Rabbit pAb	ER1908-82	Human, Rat	WB
Relative expression	FGF23 Recombinant Rabbit mAb [PSH19-00]	HA724036	Human, Mouse, Rat	WB
	FGF3 Rabbit pAb	ER61973	Human, Mouse	WB
	FGF5 Rabbit pAb	ER1908-83	Human, Mouse, Rat	IHC-P,IF-Tissue
	FGF-6 Rabbit pAb	ER65475	Human, Mouse, Rat	WB
	FGF-8 Rabbit pAb	ER65576	Human, Mouse, Rat	WB
	FGF-9 Rabbit pAb	ER65478	Human, Mouse, Rat	WB,IHC-P
Relative expression, KD	FGFR1 Recombinant Rabbit mAb [SD08-25]	ET1612-62	Human	WB,IF-Cell
	FGFR2 Recombinant Rabbit mAb [JM10-60]	ET1703-78	Human, Mouse, Rat	WB,IP
	FGFR2/CD332 Recombinant Mouse mAb [A5-E10-G7-R]	EM50103	Human	WB
	FGFR2/CD332 Rabbit pAb	R1311-7	Human	WB,IHC-P
	FGFR2/CD332 Rabbit pAb	R1310-10	Human	WB,IF-Cell,IHC-P,FC
	FGFR2/CD332 Recombinant Mouse mAb [A5-E10-G7-R] - BSA and Azide free	HA610003	Human	WB
	FOXA1 Recombinant Mouse mAb [A2E8-R]	HA601298	Human	WB,IF-Cell,IHC-P
	FOXA1 Recombinant Rabbit mAb [JF10-02] - BSA and Azide free	HA750363	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	FOXA1 Recombinant Rabbit mAb [JF10-02]	ET1702-89	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	FOXA1 Mouse mAb [A2E8]	EM1902-34	Human	WB,IHC-P,FC,IF-Cell
	FOXA1 Recombinant Mouse mAb [A2E8-R] - BSA and Azide free	HA610171	Human	WB,IF-Cell,IHC-P
	FOXA1 Mouse mAb [A2E9]	HA600022	Human	WB,IHC-P,IF-Cell
	Frizzled-2 Rabbit pAb	ER63470	Human, Mouse, Rat	WB,IF-Cell
	Frizzled-4 Rabbit pAb	ER65833	Human, Mouse, Rat	IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Frizzled 8 Recombinant Rabbit mAb [JB40-16]	ET7107-68	Human, Mouse	WB,IHC-P
	FZD2 Rabbit pAb	HA500406	Human, Mouse	WB,IHC-P
Relative expression	GATA2 Recombinant Rabbit mAb [PSH04-53]	HA722152	Human, Rat	WB,IF-Cell
Relative expression	GATA2 Recombinant Rabbit mAb [PSH04-53] - BSA and Azide free	HA750960	Human, Rat	WB,IF-Cell
	GATA-2/3 Rabbit pAb	ER63507	Human, Mouse	WB,IHC-P,IF-Cell
Relative expression	GATA2 + GATA3 Recombinant Rabbit mAb [JE04-46]	HA722713	Human, Mouse, Rat	WB,IHC-P
	GATA3 Rabbit pAb	ER1902-69	Human, Rat	WB,IHC-P,FC
	GATA3 Recombinant Mouse mAb [A3G3-R] - BSA and Azide free	HA610069	Human	WB,IF-Cell,IHC-P,FC
	GATA3 Mouse mAb [A3G3]	EM1902-33	Human	WB,IF-Cell,IHC-P
	GATA3 Rabbit pAb	ER1901-20	Human, Rat	WB,IF-Cell,IHC-P,FC
	GATA3 Mouse mAb [A3G6]	EM1902-23	Human	WB
	GATA3 Recombinant Mouse mAb [A3G3-R]	HA601189	Human	WB,IF-Cell,IHC-P,FC
	GATA3 Recombinant Mouse mAb [PD00-41]	HA601132	Human	WB,IHC-P,IF-Cell
	GRB2 Mouse mAb [11-B11]	EM1701-41	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	GRB2 Recombinant Rabbit mAb [JE49-73]	HA721923	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC
	GRB2 Rabbit pAb	ER31203	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	GSK3 alpha Rabbit pAb	R1312-1	Human	WB,IF-Cell,IHC-P,FC
	GSK3 alpha Recombinant Rabbit mAb [JE33-27]	HA722217	Human, Mouse, Rat, Monkey	WB,IHC-P
	GSK3 alpha/beta Rabbit pAb	ER63648	Human, Mouse, Rat	WB,IHC-P
	GSK3 beta Rabbit pAb	R1510-4	Human, Mouse, Rat	IF-Cell,IHC-P,FC
KD	GSK3 beta Recombinant Rabbit mAb [SY28-03]	ET1607-71	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IHC-Fr
KD	GSK3 beta Recombinant Rabbit mAb [SY28-03] - BSA and Azide free	HA750131	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IHC-Fr
	GTPase HRAS Recombinant Rabbit mAb [JE42-15]	HA721355	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	HER2 / ErbB2 Recombinant Rabbit mAb [PD00-99]	HA721210	Human	WB,IHC-P,IF-Tissue
	HER2 / ErbB2 Recombinant Rabbit mAb [PD00-53]	HA721178	Human	IHC-P,WB,IF-Cell,FC,mIHC
	HER2 / ErbB2 Rabbit pAb	R1511-33	Human	WB,IF-Cell,IHC-P,FC
	HER2 / ErbB2 Rabbit pAb	ER0106	Human	WB,IHC-P
	HER2 / ErbB2 Recombinant Rabbit mAb	IRS034	Human	mIHC
	Hes1 Recombinant Rabbit mAb [SC06-21] - BSA and Azide free	HA750240	Human, Mouse, Rat	WB,IHC-P,IHC-Fr,IF-Tissue,IF-Cell,FC
	Hes1 Recombinant Rabbit mAb [SC06-21]	ET1610-97	Human, Mouse, Rat	WB,IHC-P,IHC-Fr,IF-Tissue,IF-Cell,FC
	HES5 Rabbit pAb	ER1910-51	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	HES5 Recombinant Rabbit mAb [JE35-96]	HA722699	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P
	Human E-Cadherin Recombinant Rabbit mAb [SY0287] - BSA and Azide free (Detector)	HA723056	Human	ELISA(Det)
	Human E-Cadherin Recombinant Rabbit mAb [PSH09-02] - BSA and Azide free (Capture)	HA723055	Human	ELISA(Cap)
	Biotin Conjugated Human E-Cadherin Recombinant Rabbit mAb [SY0287] - BSA and Azide free (Detector)	HA723057B	Human	ELISA(Det),ELISA
	Human EGF Recombinant Rabbit mAb [PSH02-20] - BSA and Azide free (Detector)	HA721783	Human	ELISA(Det)
	Human EGF Recombinant Rabbit mAb [PSH02-19] - BSA and Azide free (Capture)	HA721782	Human	ELISA(Cap)
	Biotin Conjugated Human EGF Recombinant Rabbit mAb [PSH02-20] - Detector	HA722638B	Human	ELISA(Det)
	Human FGF-19 Recombinant Rabbit mAb [PSH13-30] - BSA and Azide free (Detector)	HA725101	Human	ELISA(Det)
	Human FGF-19 Recombinant Rabbit mAb [PSH13-29] - BSA and Azide free (Capture)	HA725100	Human	ELISA(Cap)
	Biotin Conjugated Human FGF-19 Recombinant Rabbit mAb [PSH13-30] - Detector	HA725102B	Human	ELISA(Det),ELISA
	Human FGF-2 Recombinant Rabbit mAb [PSH02-09] - BSA and Azide free (Detector)	HA721781	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human FGF-2 Recombinant Rabbit mAb [PSH02-09]	HA722639B	Human	ELISA(Det)
	Human FGF-2 Recombinant Rabbit mAb [PSH02-08] - BSA and Azide free (Capture)	HA721780	Human	ELISA(Cap)
	Human/Mouse/Rat FGF1 Recombinant Rabbit mAb [PSH18-34] - BSA and Azide free (Capture)	HA723987	Human, Mouse, Rat	ELISA(Cap)
	Human/Mouse/Rat FGF1 Recombinant Rabbit mAb [PSH18-35] - BSA and Azide free (Detector)	HA723988	Human, Mouse, Rat	ELISA(Det)
	Biotin Conjugated Human/Mouse/Rat FGF1 Recombinant Rabbit mAb [PSH18-35] - Detector	HA723988B	Human, Mouse, Rat	ELISA(Det),ELISA
	Biotin Conjugated Human RANKL Recombinant Rabbit mAb [PSH16-59] - Detector	HA723861B	Human	ELISA(Det),ELISA
	Human RANKL Recombinant Rabbit mAb [PSH16-58] - BSA and Azide free (Capture)	HA723860	Human	ELISA(Cap)
	Human RANKL Recombinant Rabbit mAb [PSH16-60] - BSA and Azide free (Capture)	HA723862	Human	ELISA(Cap)
	Human RANKL Recombinant Rabbit mAb [PSH16-59] - BSA and Azide free (Detector)	HA723861	Human	ELISA(Det)
	Biotin Conjugated Human VE Cadherin Recombinant Rabbit mAb [PSH02-62] - Detector	HA722635B	Human	ELISA(Det)
	Human VE Cadherin Recombinant Rabbit mAb [PSH02-62] - BSA and Azide free (Detector)	HA721839	Human	ELISA(Det)
	Human VE Cadherin Recombinant Rabbit mAb [PSH02-61] - BSA and Azide free (Capture)	HA721838	Human	ELISA(Cap)
Cell treatment	IGF1 Recombinant Rabbit mAb [JE30-01]	HA722527	Human	WB
	IGF1 (E peptide) Rabbit pAb	HA500107	Human	WB,IHC-P,FC
	IGF1R Rabbit pAb	R1310-6	Human	WB,IHC-P
	IGF1R Mouse mAb [D5-G2]	HA600033	Human	WB,FC
Relative expression	IGF1 Receptor beta Recombinant Rabbit mAb [PSH11-72]	HA723368	Human, Mouse, Rat	WB,IF-Cell
Relative expression	IGF1 Receptor beta Recombinant Rabbit mAb [PSH11-72] - BSA and Azide free	HA751427	Human, Mouse, Rat	WB,IF-Cell
	IGF I Rabbit pAb	ER1911-10	Human, Mouse, Rat, Dog	WB,IHC-P,IF-Tissue
	IGF-IR Rabbit pAb	ER63734	Human, Mouse, Rat	IF-Cell,WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD, Relative expression	Insulin Receptor Beta Recombinant Rabbit mAb [SN20-13]	ET1611-90	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC
KD, Relative expression	Insulin Receptor Beta Recombinant Rabbit mAb [SN20-13] - BSA and Azide free	HA750276	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC
	Jagged 1 Rabbit pAb	R1706-10	Human, Mouse, Rat	IF-Cell,IHC-P,FC
	Jagged1 Recombinant Rabbit mAb [JF96-06]	ET1702-63	Human, Mouse, Rat	WB,IHC-P
	JunD Recombinant Rabbit mAb [PSH08-85]	HA723037	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP,ChIP
	JunD Recombinant Rabbit mAb [PSH08-85] - BSA and Azide free	HA751250	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP,ChIP
	JunD Recombinant Rabbit mAb [SD0830]	ET1612-92	Human	WB,IF-Cell,IF-Tissue,IP,IHC-P,FC
	JunD Recombinant Rabbit mAb [JE59-51]	HA722913	Human, Mouse, Rat	WB,IF-Cell,IHC-P,ChIP,IP
	KGF Recombinant Rabbit mAb [JE40-14]	HA722404	Human, Mouse, Rat, Monkey	WB
	Ki67 Recombinant Rabbit mAb	IRS036	Human	mIHC
	Ki67 Rabbit pAb	ER1706-46	Human	WB,IF-Cell,IHC-P,FC
	Ki67 Recombinant Rabbit mAb [SR00-02]	HA721115	Human, Mouse, Rat, Cynomolgus monkey	IHC-P,IF-Tissue,mIHC,IF-Cell,FC
	iFluor™ 488 Conjugated Ki67 Recombinant Rabbit mAb [SR00-02]	HA720162F	Human	IF-Cell,IF-Tissue,FC
	Ki67 Recombinant Rabbit mAb [PSH0-02] - BSA and Azide free	HA750570	Human	IHC-P,FC
	Ki67 Rabbit pAb	ER1902-75	Human	WB,IF-Cell,FC
	Ki67 Recombinant Mouse mAb [PD00-10]	HA601053	Human, Mouse	IHC-P,IF-Tissue
	iFluor™ 647 Conjugated Ki67 Recombinant Rabbit mAb [SR00-02]	HA720163F	Human, Mouse, Rat	IF-Cell,IF-Tissue,FC,IHC-Fr
	Ki67 Recombinant Rabbit mAb	IRS057	Mouse	mIHC
	Ki67 Recombinant Rabbit mAb [ST50-01]	ET1609-34	Human	WB,IF-Cell,IHC-P,FC
	Ki67 Rabbit pAb	RT1349	Human	WB,IP,IF,IHC-P
	Ki67 Recombinant Rabbit mAb [PSH0-02]	HA721254	Human	IHC-P,FC
	K-Ras Rabbit pAb	ER65706	Human, Mouse, Rat	WB,IHC-P
	KRAS Recombinant Mouse mAb [A8E6] - BSA and Azide free	HA610008	Human, Mouse, Rat	WB,IHC-P
	KRAS Mouse mAb [A8E5]	HA601059	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	KRAS Recombinant Mouse mAb [A8E6]	HA601060	Human, Mouse, Rat	WB,IHC-P
	KRAS Mouse mAb [A8E4]	HA601058	Human, Mouse, Rat	WB,IF-Cell,FC
Relative expression, KD	LEF1 Recombinant Rabbit mAb [JE52-28]	HA721992	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	LEF1 Mouse mAb [B11-A4]	EM1708-21	Human	WB,IHC-P
	LEF1 Rabbit pAb	HA500273	Human, Mouse, Rat	WB,IHC-P,IF-Cell
KD	MEK1 Recombinant Rabbit mAb [SZ22-01]	ET1603-20	Human, Mouse, Rat, Cow, Dog	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	MEK1/2 Recombinant Rabbit mAb [SR13-07]	ET1602-3	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	MEK1/2 Recombinant Rabbit mAb [SR13-07] - BSA and Azide free	HA750042	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	MEK2 Recombinant Rabbit mAb [SD207-08]	ET1612-6	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Mouse FGF-21 Recombinant Rabbit mAb [PSH15-27] - BSA and Azide free (Detector)	HA723732	Mouse	ELISA(Det)
	Mouse FGF-21 Recombinant Rabbit mAb [PSH15-26] - BSA and Azide free (Capture)	HA723731	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse FGF-21 Recombinant Rabbit mAb [PSH15-27] - Detector	HA723733B	Mouse	ELISA(Det),ELISA
Cell treatment	Biotin Conjugated Mouse RANKL Recombinant Rabbit mAb [PSH07-92]	HA722925B	Mouse	ELISA(Det)
Cell treatment	Mouse RANKL Recombinant Rabbit mAb [PSH07-92] - BSA and Azide free (Detector)	HA722924	Mouse	ELISA(Det)
Cell treatment	Mouse RANKL Recombinant Rabbit mAb [PSH07-90] - BSA and Azide free (Capture)	HA722922	Mouse	ELISA(Cap)
Cell treatment	Mouse RANKL Recombinant Rabbit mAb [PSH07-91] - BSA and Azide free (Capture)	HA722923	Mouse	ELISA(Cap)
	HRP Conjugated Mouse/Rat E-Cadherin Recombinant Rabbit mAb [PSH17-25] - Detector	HA723912H	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat E-Cadherin Recombinant Rabbit mAb [PSH17-24] - BSA and Azide free (Capture)	HA723911	Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Mouse/Rat E-Cadherin Recombinant Rabbit mAb [PSH17-25] - Detector	HA723912B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat E-Cadherin Recombinant Rabbit mAb [PSH17-25] - BSA and Azide free (Detector)	HA723912	Mouse, Rat	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	mTOR Recombinant Rabbit mAb [SU30-00] - BSA and Azide free	HA750153	Human, Mouse, Rat	WB,IHC-P,IP,IF-Cell
	mTOR Rabbit pAb	R1510-21	Human, Mouse, Rat	IF-Cell,IHC-P
	mTOR Rabbit pAb	HA500126	Human, Rat, Mouse	WB,IF-Cell,IHC-P
	mTOR Recombinant Rabbit mAb [SU30-00]	ET1608-5	Human, Mouse, Rat	WB,IHC-P,IP,IF-Cell
	HRP Conjugated Myc tag Rabbit pAb	R1208-2	Species independent	WB
	Myc tag Rabbit pAb	R1208-1	Species independent	WB,IP,IF-Cell,FC
	N Cadherin Recombinant Antibody [SY02-46] - Mouse IgG1 (Chimeric)	HA601464	Human, Mouse, Rat	IHC-Fr,IHC-P
Relative expression, KD	N Cadherin Recombinant Rabbit mAb [SY02-46]	ET1607-37	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue
	N Cadherin Recombinant Antibody [SY02-46] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610298	Human, Mouse, Rat	IHC-Fr,IHC-P
Relative expression, KD	N Cadherin Recombinant Rabbit mAb [SY02-46] - BSA and Azide free	HA750119	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue
	N-Cadherin Rabbit pAb	ER0503	Human, Mouse, Rat	WB,IHC,IF
	N-Cadherin Mouse mAb [A10-4]	M1304-1	Human, Mouse	WB,IHC-P
	NCoA-3 Rabbit pAb	ER65603	Human, Mouse, Rat	WB,IHC-P
KD	NF-kB p100 / NFKB2 Mouse mAb [16B3]	EM1901-79	Human, Mouse, Rat	WB,IHC-P,FC
KD	NF-kB p100 / NFKB2 Mouse mAb [16B2]	EM1901-78	Human	WB,IHC-P,FC
KD	NF-kB p100 / NFKB2 Recombinant Rabbit mAb [JM82-03]	ET1705-45	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
	N/H/K-Ras Rabbit pAb	ER64095	Human, Mouse, Rat	WB,IHC-P
KD	Notch 1 Recombinant Rabbit mAb [SJ205]	ET1606-55	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	Notch 1 Recombinant Rabbit mAb [SJ205] - BSA and Azide free	HA750111	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Notch2 Recombinant Rabbit mAb [PSH0-55]	HA721402	Human, Mouse, Rat	WB
	Notch2 Recombinant Rabbit mAb [PSH0-55] - BSA and Azide free	HA750632	Human, Mouse, Rat	WB
	Notch2 Recombinant Rabbit mAb [PSH0-54] - BSA and Azide free	HA750631	Human, Mouse, Rat	WB
	Notch2 Recombinant Rabbit mAb [PSH0-54]	HA721401	Human, Mouse, Rat	WB
	NOTCH3 Rabbit pAb	HA500519	Human	WB,FC
	NOTCH4 Recombinant Rabbit mAb [JM101-05]	ET1703-20	Mouse	WB,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	N-Shc Rabbit pAb	ER64229	Human, Mouse, Rat	IHC-P
	p21 Recombinant Rabbit mAb [JE00-63]	HA722685	Human	WB,IF-Cell,IHC-P,IP
	p21 Recombinant Rabbit mAb	IRS260RB	Mouse	mIHC
	p21 Mouse mAb [A8C12]	HA601017	Human	IHC-P
	p21 Recombinant Mouse mAb [A8C11-R] - BSA and Azide free	HA610098	Human	IHC-P
	p21 Mouse mAb [A8C11]	HA601018	Human	IHC-P
	p21 Recombinant Mouse mAb [A8C11-R]	HA601221	Human	IHC-P
	p21 Recombinant Rabbit mAb [PSH03-93]	HA722065	Mouse	WB,IF-Cell,IHC-P,FC
	p21 Recombinant Rabbit mAb [PSH03-93] - BSA and Azide free	HA750918	Mouse	WB,IF-Cell,IHC-P,FC
Cell treatment, KD, Relative expression	p53 Recombinant Mouse mAb [6B5-R]	HA601315	Human, Mouse, Monkey	WB,IF-Cell
Cell treatment	p53 Recombinant Rabbit mAb [PSH03-95] - BSA and Azide free	HA750920	Mouse	WB,IF-Cell,IP
Cell treatment, Relative expression, KD	p53 Recombinant Mouse mAb [6B5-R] - BSA and Azide free	HA610187	Human, Mouse, Monkey	WB,IF-Cell
KD	p53 Recombinant Rabbit mAb [JB42-26]	ET7107-33	Human	WB,IHC-P,FC,IP,IF-Cell,IF-Tissue
	p53 Recombinant Rabbit mAb [PD01-32] - BSA and Azide free	HA750545	Human	IHC-P
KD	p53 Rabbit pAb	R1308-7	Human	WB,IF-Cell
KO	p53 Mouse mAb [7-5]	EM20603	Human	WB,IF-Cell,IHC-P,FC
Relative expression, KD, Cell treatment	p53 Recombinant Rabbit mAb [SA39-07] - BSA and Azide free	HA750005	Human	WB,IF-Cell,ChIP
KD	p53 Recombinant Mouse mAb [6B2-R] - BSA and Azide free	HA610194	Human, Monkey	WB,IF-Cell
KD	p53 Recombinant Mouse mAb [6B2-R]	HA601322	Human, Monkey	WB,IF-Cell
Relative expression	p53 Mouse mAb [6B2]	EM1701-91	Human, Mouse	WB,IF-Cell
Relative expression, KD	p53 Recombinant Rabbit mAb [SY010-6]	ET1605-16	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	p53 Mouse mAb [6B5]	EM1701-92	Human	WB
	p53 Recombinant Rabbit mAb [PD01-32]	HA721213	Human	IHC-P
Cell treatment	p53 Recombinant Rabbit mAb [PSH03-95]	HA722074	Mouse	WB,IF-Cell,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	p53 Mouse mAb [C8-A11]	M1312-2	Human	WB,IF-Cell,IHC-P
Relative expression, KD, Cell treatment	p53 Recombinant Rabbit mAb [SA39-07]	ET1601-13	Human	WB,IF-Cell,ChIP
Cell treatment	p53 (acetyl K370) Recombinant Rabbit mAb [SR40-09] - BSA and Azide free	HA750047	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,FC
Cell treatment	p53 (acetyl K370) Recombinant Rabbit mAb [SR40-09]	ET1602-38	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,FC
Cell treatment	p53 (acetyl K382) Recombinant Rabbit mAb [PSH06-87] - BSA and Azide free	HA751111	Human, Mouse	WB,IF-Cell
Cell treatment	p53 (acetyl K382) Recombinant Rabbit mAb [PSH06-88] - BSA and Azide free	HA751112	Human	WB,IF-Cell
Cell treatment	p53 (acetyl K382) Recombinant Rabbit mAb [PSH06-88]	HA722763	Human	WB,IF-Cell
Cell treatment	p53 (acetyl K382) Recombinant Rabbit mAb [PSH06-87]	HA722762	Human, Mouse	WB,IF-Cell
Cell treatment	p53 (acetyl K382) Recombinant Rabbit mAb [SD0801]	ET1612-71	Human, Mouse	WB,IF-Cell,IHC-P
	p70 S6 Kinase Rabbit pAb	ER1916-13	Human, Rat	WB,IHC-P,IF-Tissue
	p70 S6 Kinase Rabbit pAb	HA500095	Human, Mouse, Rat	WB,IHC-P
	p70 S6 Kinase Recombinant Rabbit mAb [PSH13-50] - BSA and Azide free	HA751469	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IP
	p70 S6 Kinase Recombinant Rabbit mAb [PSH13-50]	HA723539	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IP
	p70 S6 Kinase Recombinant Rabbit mAb [JE03-20]	HA722520	Human, Mouse, Rat	WB,IHC-P,IP
	p70 S6 Kinase beta Rabbit pAb	ER31205	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	P70 S6 Kinase beta Recombinant Rabbit mAb [JE35-69]	HA721354	Human, Mouse	WB,IF-Cell,FC,IHC-P
	Pan Cadherin Recombinant Rabbit mAb	RE6045	Human, Mouse, Rat	WB,ICC,IHC-P,FC,IP
	Pan-Cadherin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS288RT	Mouse	mIHC
	Pan-Cadherin Recombinant Rabbit mAb	IRS185RB	Mouse	mIHC
	Pan-Cadherin Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS285MS	Human	mIHC
	Pan-Cadherin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS287RT	Human	mIHC
	Pan-Cadherin Recombinant Rabbit mAb [ST54-01] - BSA and Azide free	HA750190	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,IHC-Fr

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Pan-Cadherin Recombinant Rabbit mAb [ST54-01]	ET1609-70	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,IHC-Fr
	Pan-Cadherin Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS286MS	Mouse	mIHC
	pan Cytokeratin Recombinant Antibody - Rabbit IgG (Chimeric)	IRS186RB	Human	mIHC
	iFluor™ 488 Conjugated pan Cytokeratin Recombinant Mouse mAb [PDH0-10]	HA600111F	Human	IF-Tissue
	pan Cytokeratin Recombinant Antibody [PDH09-10] - Rat IgG1 (Chimeric)	HA601519	Human	IHC-P
	pan Cytokeratin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS187RT	Human	mIHC
	pan Cytokeratin Recombinant Mouse mAb [PDH09-10]	HA601138	Human	IHC-P,mIHC,IF-Cell,IF-Tissue,WB,FC
Cell treatment	Phospho-AKT1 (S124) Recombinant Rabbit mAb [JJ08-46]	ET1701-36	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
Cell treatment	Phospho-AKT1 (S129) Recombinant Rabbit mAb [JE59-50]	HA721677	Human, Mouse, Rat	WB
Cell treatment	Phospho-AKT1 (T450) Recombinant Rabbit mAb [SD08-12]	ET1612-73	Human, Mouse, Rat	WB,IHC-P,IP,IF-Cell,FC
Cell treatment	Phospho-AKT2 (S474) Recombinant Rabbit mAb [PSH09-46] - BSA and Azide free	HA751293	Human, Mouse, Rat	WB,FC
Cell treatment	Phospho-AKT2 (S474) Rabbit pAb	HA500116	Human, Mouse	WB,IHC-P
Cell treatment	Phospho-AKT2 (S474) Recombinant Rabbit mAb [PSH09-46]	HA723108	Human, Mouse, Rat	WB,FC
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [SY28-05]	ET1607-73	Human, Mouse, Rat, Dog	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [PSH04-44] - BSA and Azide free	HA750953	Human, Mouse, Rat	WB,IHC-P,FC
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [SY28-05] - BSA and Azide free	HA750132	Human, Mouse, Rat, Dog	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [PSH04-44]	HA722129	Human, Mouse, Rat	WB,IHC-P,FC
Cell treatment	Phospho-AKT (T308) Recombinant Rabbit mAb [PSH08-10] - BSA and Azide free	HA751205	Human	WB,IF-Cell,FC
Cell treatment	Phospho-AKT (T308) Recombinant Rabbit mAb [PSH08-10]	HA722951	Human	WB,IF-Cell,FC
Cell treatment	Phospho-Beta Catenin (S33 + S37) Recombinant Rabbit mAb [JE59-59]	HA721580	Human, Rat, Mouse	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-Beta Catenin (S552) Recombinant Rabbit mAb [PSH08-72] - BSA and Azide free	HA751237	Human, Mouse, Rat	WB,IF-Cell,IHC- P,FC
Cell treatment	Phospho-Beta Catenin (S552) Recombinant Rabbit mAb [PSH08-72]	HA723024	Human, Mouse, Rat	WB,IF-Cell,IHC- P,FC
Cell treatment	Phospho-Beta Catenin (T41/S45) Recombinant Rabbit mAb [JE54-02]	HA722316	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-BRAF (T401) Recombinant Rabbit mAb [JJ08-72]	ET1701-20	Human, Mouse, Rat	WB,IF-Cell,IP,IHC-P
Cell treatment	Phospho-c-Jun (S63) Recombinant Rabbit mAb [SY0297]	ET1608-4	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-c-Jun (S63) Recombinant Rabbit mAb [SY0297] - BSA and Azide free	HA750146	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-c-Jun (S73)+JunD (S100) Recombinant Rabbit mAb [PSH07-12] - BSA and Azide free	HA751133	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-c-Jun (S73)+JunD (S100) Recombinant Rabbit mAb [JE65-46]	HA722475	Human, Mouse, Rat	WB,IF-Cell,FC
Cell treatment	Phospho-c-Jun (S73)+JunD (S100) Recombinant Rabbit mAb [PSH07-12]	HA722800	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-c-Jun (T91) Recombinant Rabbit mAb [JJ080-9]	ET1701-32	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC-P,IP
Cell treatment	Phospho-c-Myc (S62) Recombinant Rabbit mAb [ST49-08] - BSA and Azide free	HA750186	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC-P,IP
Cell treatment	Phospho-c-Myc (S62) Recombinant Rabbit mAb [ST49-08]	ET1609-64	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC-P,IP
Cell treatment	Phospho-c-Myc (T58) Recombinant Rabbit mAb [SN60-01]	ET1611-24	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,FC,IHC-P
Cell treatment	Phospho-c-Myc (T58 + S62) Recombinant Rabbit mAb [SZ02-06]	ET1603-24	Human, Mouse	WB,IP,FC
	Phospho-EGFR (S695) Recombinant Rabbit mAb [SN07-36]	ET1611-50	Human	WB,IF-Cell,IF-Tissue
Cell treatment	Phospho-EGFR (T669) Recombinant Rabbit mAb [PSH10-03] - BSA and Azide free	HA751333	Human, Mouse, Rat	WB,IP
Cell treatment	Phospho-EGFR (T669) Recombinant Rabbit mAb [PSH10-03]	HA723172	Human, Mouse, Rat	WB,IP
Cell treatment	Phospho-EGFR (Y1045) Recombinant Rabbit mAb [PSH03-28] - BSA and Azide free	HA750868	Human	WB,IF-Cell,FC
Cell treatment	Phospho-EGFR (Y1045) Recombinant Rabbit mAb [PSH03-28]	HA721970	Human	WB,IF-Cell,FC
Cell treatment	Phospho-EGFR (Y1068) Recombinant Rabbit mAb [SD2055]	ET1612-30	Human	WB,IF-Cell,IF- Tissue,IP,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-EGFR (Y1068) Recombinant Rabbit mAb [SD2055] - BSA and Azide free	HA750290	Human	WB,IF-Cell,IF-Tissue,IP,FC
Cell treatment	Phospho-EGFR (Y1092) Recombinant Rabbit mAb [SJ0194]	ET1606-44	Human, Mouse	WB,IF-Tissue,IHC-P
	Phospho-EGFR (Y1148) Rabbit pAb	HA500580	Human	WB,IF-Cell
Cell treatment	Phospho-EGFR (Y1173) Recombinant Rabbit mAb [SC57-04]	ET1610-4	Human, Mouse	WB,IF-Cell,IHC-P,IP
Cell treatment	Phospho-EGFR (Y845) Recombinant Rabbit mAb [JE46-77]	HA721488	Human	WB
Cell treatment	Phospho-Erk1 (T202)+Erk2 (T185) Recombinant Rabbit mAb [SZ2-4] - BSA and Azide free	HA750061	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
Cell treatment	Phospho-Erk1 (T202)+Erk2 (T185) Recombinant Rabbit mAb [SZ2-4]	ET1603-22	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
Cell treatment	Phospho-Erk1 (T202 + Y204) + Erk2 (T185 + Y187) Recombinant Rabbit mAb [SC58-01]	ET1610-13	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP
Cell treatment	Phospho-Erk1 (T202 + Y204) + Erk2 (T185 + Y187) Recombinant Rabbit mAb [SC58-01] - BSA and Azide free	HA750199	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP
Cell treatment	Phospho-Estrogen Receptor alpha (S118) Recombinant Rabbit mAb [SC05-87]	ET1610-32	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	Phospho-GATA3 (S308) Recombinant Rabbit mAb [ST44-09]	ET1609-17	Human, Mouse, Rat	WB,IHC-P,IP,IF-Tissue
Cell treatment	Phospho-GSK3 (alpha + beta) (Y216 + Y279) Recombinant Rabbit mAb [SY26-05] - BSA and Azide free	HA750127	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,IF-Tissue,FC
Cell treatment	Phospho-GSK3 (alpha + beta) (Y216 + Y279) Recombinant Rabbit mAb [SY26-05]	ET1607-54	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,IF-Tissue,FC
Cell treatment	Phospho-GSK3 alpha (S21) Recombinant Rabbit mAb [PSH09-11] - BSA and Azide free	HA751265	Human, Mouse	WB,IHC-P,Dot Blot
Cell treatment	Phospho-GSK3 alpha (S21) Recombinant Rabbit mAb [PSH09-11]	HA723069	Human, Mouse	WB,IHC-P,Dot Blot
Cell treatment	Phospho-GSK3 beta (S9) Recombinant Rabbit mAb [SY02-71]	ET1607-60	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
Cell treatment	Phospho-GSK3 beta (S9) Recombinant Rabbit mAb [SY02-71] - BSA and Azide free	HA750128	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Phospho-HER2 / ErbB2 (T686) Mouse mAb [4G3]	RT1275	Human, Mouse, Rat	WB,IP
Cell treatment	Phospho-HER2 / ErbB2 (Y1139) Recombinant Rabbit mAb [JE60-90]	HA722201	Human	WB,IF-Cell,FC
Cell treatment	Phospho-HER2 / ErbB2 (Y1221 + Y1222) Recombinant Rabbit mAb [JE44-12]	HA721433	Human	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-HER2 / ErbB2 (Y1248) Recombinant Rabbit mAb [PSH04-03]	HA722084	Human	WB,IF-Cell
Cell treatment	Phospho-HER2 / ErbB2 (Y1248) Recombinant Rabbit mAb [PSH04-03] - BSA and Azide free	HA750926	Human	WB,IF-Cell
Cell treatment	Phospho-JunD (S255) Recombinant Rabbit mAb [JJ08-21]	ET1701-35	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC-P
Cell treatment	Phospho-MEK1/2 (S218 + S222) Recombinant Rabbit mAb [ST0490]	ET1609-50	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC-P,IP
	Phospho-MEK1 (S298) Recombinant Rabbit mAb [SD206-7]	ET1612-40	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC-P
Cell treatment	Phospho-MEK1 (T292) Recombinant Rabbit mAb [SD2088]	ET1612-42	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-mTOR (S2448) Mouse mAb [A5D5]	HA600094	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-mTOR (S2481) Recombinant Rabbit mAb [JE59-62]	HA721632	Human	WB
Cell treatment	Phospho-p53 (S15) Recombinant Rabbit mAb [PSH01-98] - BSA and Azide free	HA750772	Human	WB,IF-Cell
Cell treatment	Phospho-p53 (S15) Recombinant Rabbit mAb [PSH01-98]	HA721756	Human	WB,IF-Cell
Cell treatment	Phospho-p53 (S20) Recombinant Rabbit mAb [PSH06-86]	HA722761	Human	WB
Cell treatment	Phospho-p53 (S20) Recombinant Rabbit mAb [PSH06-86] - BSA and Azide free	HA751110	Human	WB
Cell treatment	Phospho-p53 (S33) Recombinant Rabbit mAb [JE41-97]	HA722703	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-p53 (S376) Recombinant Rabbit mAb [ST0440]	ET1609-14	Human, Mouse	WB
Cell treatment	Phospho-p53 (S392) Recombinant Rabbit mAb [SI17-04] - BSA and Azide free	HA750097	Human, Mouse, Rat	WB,IHC-P,IP
Cell treatment	Phospho-p53 (S392) Recombinant Rabbit mAb [SI17-04]	ET1606-24	Human, Mouse, Rat	WB,IHC-P,IP
Cell treatment	Phospho-p53 (T55) Recombinant Rabbit mAb [ST0439]	ET1609-13	Human	WB,IHC-P
	Phospho-p70 S6 Kinase (S434) Mouse mAb [2G1]	RT1456	Human, Mouse, Rat	WB,IP,IF,IHC-P
Cell treatment	Phospho-p70 S6 Kinase (T389) Recombinant Rabbit mAb [PSH02-23] - BSA and Azide free	HA750791	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-p70 S6 Kinase (T389) Recombinant Rabbit mAb [PSH02-23]	HA721803	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-PI3K p85 (Y467) + PI3K p55 (Y199) Recombinant Rabbit mAb [PSH01-39]	HA721673	Human, Mouse	WB,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-PI3K p85 (Y467) + PI3K p55 (Y199) Recombinant Rabbit mAb [PSH01-39] - BSA and Azide free	HA750717	Human, Mouse	WB,IF-Cell
Cell treatment	Phospho-PI3K p85 (Y467) + PI3K p55 (Y199) Recombinant Rabbit mAb [PSH01-38] - BSA and Azide free	HA750716	Human, Mouse	WB,IF-Cell
Cell treatment	Phospho-PI3K p85 (Y467) + PI3K p55 (Y199) Recombinant Rabbit mAb [PSH01-38]	HA721672	Human, Mouse	WB,IF-Cell
	Phospho-PTEN (S380) Recombinant Rabbit mAb [JJ08-81]	ET1701-46	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-Raf1 (S259) Recombinant Rabbit mAb [SD85-07] - BSA and Azide free	HA750304	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-Raf1 (S259) Recombinant Rabbit mAb [SD85-07]	ET1612-87	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	Phospho-Raf1 (S43) Recombinant Rabbit mAb [SD087-05]	ET1612-75	Human, Mouse, Rat	WB
	Phospho-Raf1 (S621) Recombinant Rabbit mAb [JJ085-05]	ET1701-3	Human, Rat, Mouse	WB,IHC-P
Cell treatment	Phospho-Rb (S608) Recombinant Rabbit mAb [PSH11-04] - BSA and Azide free	HA751394	Human, Monkey	WB
Cell treatment	Phospho-Rb (S608) Recombinant Rabbit mAb [PSH11-04]	HA723271	Human, Monkey	WB
Cell treatment	Phospho-Rb (S780) Recombinant Rabbit mAb [JE44-61]	HA721513	Human, Mouse, Rat	WB
Cell treatment	Phospho-Rb (S795) Recombinant Rabbit mAb [PSH06-89]	HA722764	Human	WB,IF-Cell
Cell treatment	Phospho-Rb (S795) Recombinant Rabbit mAb [PSH06-89] - BSA and Azide free	HA751113	Human	WB,IF-Cell
Cell treatment	Phospho-Rb (S807) Recombinant Rabbit mAb [SR3-82]	ET1602-36	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	Phospho-Rb (S807) Recombinant Rabbit mAb [SR3-82] - BSA and Azide free	HA750045	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	Phospho-Rb (S807/S811) Recombinant Rabbit mAb [PSH11-01] - BSA and Azide free	HA751391	Human, Mouse	WB,IHC-P,IF-Cell,FC,IP
Cell treatment	Phospho-Rb (S807/S811) Recombinant Rabbit mAb [PSH11-01]	HA723268	Human, Mouse	WB,IHC-P,IF-Cell,FC,IP
	PI3 Kinase p110 beta Recombinant Rabbit mAb [JE65-38]	HA722474	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	PI 3 Kinase p110 delta Recombinant Rabbit mAb [JM83-24]	ET1705-46	Human, Mouse, Rat	WB,IHC-P
KD	PI3-kinase p110 subunit beta Recombinant Rabbit mAb [SC0580]	ET1610-36	Human	WB,IP,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	PI3 Kinase p110 α Recombinant Rabbit mAb [SJ0186]	ET1606-36	Human, Mouse	WB,IF-Cell,IF-Tissue,IP,IHC-P
	PI 3 Kinase p55 gamma Recombinant Rabbit mAb [SR44-06]	ET1602-27	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
KD	PI 3 Kinase p85 alpha Mouse mAb [L3-D6]	EM1701-62	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue
	PI 3 Kinase p85 alpha Recombinant Rabbit mAb [SU04-07]	ET1608-70	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	PI 3 Kinase p85 alpha Recombinant Mouse mAb [A3-D0-R] - BSA and Azide free	HA610084	Human	WB,IF-Cell
	PI 3 Kinase p85 alpha Rabbit pAb	R1312-6	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	PI 3 Kinase p85 alpha Recombinant Rabbit mAb [SU04-07] - BSA and Azide free	HA750162	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	PI 3 Kinase p85 alpha Recombinant Mouse mAb [A3-D0-R]	HA601206	Human	WB,IF-Cell
KD	PI 3 Kinase p85 alpha Mouse mAb [A3-D0]	M1510-2	Human	WB,IF-Cell
	PI3 Kinase p85 alpha Recombinant Rabbit mAb [JE04-11]	HA722522	Human, Mouse, Rat	WB
	PI 3 Kinase p85 alpha/gamma Rabbit pAb	ER64588	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Tissue
	PI 3 Kinase p85 beta Recombinant Rabbit mAb [ST04-77]	ET1609-30	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	Progesterone receptor Rabbit pAb	ER1901-96	Human	WB,IF-Cell,FC
Relative expression	Progesterone Receptor Recombinant Rabbit mAb [JF0549]	ET1702-24	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Progesterone Receptor Rabbit pAb	ER1802-19	Human, Mouse	IF-Cell,IHC-P
	Progesterone Receptor Rabbit pAb	ER1802-37	Human, Mouse	WB,IF-Cell,FC
	PTEN Rabbit pAb	ER1706-51	Human, Mouse	WB,IF-Cell,IHC-P,FC
Relative expression, KD	PTEN Recombinant Rabbit mAb [SJ19-03]	ET1606-43	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
Relative expression, KD	PTEN Recombinant Rabbit mAb [JE51-96]	HA721962	Human, Mouse, Rat	WB,IHC-P,FC
	PTEN Rabbit pAb	ER1902-27	Human, Mouse, Rat	WB,IHC-P,FC
	PTEN Mouse mAb [4D8]	RT1519	Human, Mouse, Rat	WB,IP,IF
KD	PTEN Rabbit pAb	ER1902-33	Human, Mouse, Rat	WB,IHC-P,FC
Relative expression, KD	PTEN Recombinant Rabbit mAb [SJ19-03] - BSA and Azide free	HA750105	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	Raf1 Recombinant Rabbit mAb [JJ086-06]	ET1701-21	Human, Mouse, Rat	WB
	RANKL Rabbit pAb	HA500369	Human, Mouse, Rat	WB,IF-Cell
	Ras Recombinant Rabbit mAb [SA39-05]	ET1601-16	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,FC,IP
	Ras Rabbit pAb	ER40115	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Ras Recombinant Rabbit mAb [JE44-13]	HA721883	Human, Mouse, Rat	WB,IF-Cell,FC
	Ras Recombinant Rabbit mAb [JF10-11]	ET1702-94	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IP,FC
	Rb Recombinant Rabbit mAb [SY63-03]	ET1607-9	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Rb Recombinant Rabbit mAb [SY63-03] - BSA and Azide free	HA750135	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP
	R-cadherin Rabbit pAb	ER64773	Human, Mouse	IHC-P,IF-Cell,WB
	Sck Rabbit pAb	ER64915	Human, Rat	WB
	SHC Recombinant Rabbit mAb [JE57-45]	HA722564	Human, Mouse, Rat	WB,IHC-P
	SHC1 Mouse mAb [7C-7F]	EM1706-40	Human, Mouse	WB,IF-Cell,FC
	SHC4 Rabbit pAb	ER1916-54	Mouse	WB
	smooth muscle Myosin heavy chain 11 Recombinant Rabbit mAb [JA03-35] - BSA and Azide free	HA750417	Human, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC
	smooth muscle Myosin heavy chain 11 Recombinant Rabbit mAb [JE60-21]	HA722863	Human, Mouse, Rat	WB,IHC-P
	smooth muscle Myosin heavy chain 11 Recombinant Rabbit mAb [JA03-35]	ET1704-61	Human, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC
	smooth muscle Myosin heavy chain 11 Recombinant Rabbit mAb	IRS176RB	Human	mIHC
	SOS1 Recombinant Rabbit mAb [JE30-06]	HA720077	Human, Mouse, Rat	WB,IHC-P
	Sos 2 Rabbit pAb	ER64998	Human, Mouse	WB,IHC-P
	SP1 Recombinant Rabbit mAb [JF0950]	ET1702-02	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,ChIP,IP
	SP1 Recombinant Rabbit mAb [JF0950] - BSA and Azide free	HA750330	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,ChIP,IP
	SRC1 Recombinant Rabbit mAb [JB66-31]	ET7107-44	Human	WB,IHC-P,IP,IF-Cell
	TCF-4 / TCF7L2 Recombinant Rabbit mAb [SC06-90]	ET1610-73	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	TNFSF11 Mouse mAb [8A7B9]	EM1710-42	Human, Monkey	WB,IHC-P,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	VE Cadherin Recombinant Rabbit mAb [PSH03-15]	HA721943	Human	WB,IF-Cell,FC
KD	VE Cadherin Recombinant Rabbit mAb [PSH03-15] - BSA and Azide free	HA750856	Human	WB,IF-Cell,FC
	Wnt-1 Rabbit pAb	ER65317	Human, Mouse	WB,IHC-P,IF-Cell
	Wnt10a Rabbit pAb	HA500356	Human, Mouse	WB
	Wnt16 Rabbit pAb	HA500145	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Wnt2b Recombinant Rabbit mAb [JB48-34]	ET7107-20	Human, Mouse, Rat	WB,IF-Cell,FC
	WNT4 Mouse mAb [4G1]	RT1662	Human, Mouse, Rat	WB,IP,IF,IHC-P
	WNT8A Recombinant Rabbit mAb [PSH0-39]	HA721341	Human, Mouse, Rat	WB,IF-Cell,FC
	WNT8A Recombinant Rabbit mAb [PSH0-39] - BSA and Azide free	HA750617	Human, Mouse, Rat	WB,IF-Cell,FC
	WNT9B Rabbit pAb	HA500412	Human, Mouse	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层