

上皮间质转化 (EMT) 研究-产品清单

高级验证	产品名	货号	反应种属	验证应用
	Beta Catenin Recombinant Mouse mAb [A6-F8-R]	HA601179	Human, Mouse, Rat	WB,IF-Cell,FC
KD	Beta Catenin Recombinant Rabbit mAb [SA30-04] - BSA and Azide free	HA750023	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IP,IF-Cell,IHC-Fr,FC
	Beta Catenin Rabbit pAb	ER0805	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell,IF-Tissue
	Beta Catenin Recombinant Rabbit mAb	IRS043	Human	mIHC
	iFluor™ 488 Conjugated Beta Catenin Recombinant Rabbit mAb [SA30-04]	HA720182F	Human, Mouse	IF,ICC
	Beta Catenin Rabbit pAb	0407-16	Human, Mouse, Rat, Zebrafish	WB,IF-Cell
	Beta Catenin Recombinant Rabbit mAb	IRS061	Mouse	mIHC
	Beta Catenin Mouse mAb [A6-F8]	EM0306	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue,ChIP
	Beta Catenin Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS277MS	Human	mIHC
KD	Beta Catenin Recombinant Rabbit mAb [SA30-04]	ET1601-5	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IP,mIHC,IF-Cell,IHC-Fr,FC
	Beta Catenin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS278RT	Human	mIHC
	Beta Catenin Recombinant Mouse mAb [A6-F8-R] - BSA and Azide free	HA610062	Human, Mouse, Rat	WB,IF-Cell,FC
	Beta Catenin Mouse mAb [10-C0-B7]	M1405-6	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Relative expression	E-Cadherin Recombinant Rabbit mAb [PSH13-75]	HA723564	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Cell,FC,IP
KD	E-Cadherin Mouse mAb [A0-G11-2]	EM0502	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	E-Cadherin Rabbit pAb	0407-25	Human, Mouse	WB,IF-Cell,FC
	E-Cadherin Recombinant Mouse mAb [A0-G11-2-R] - BSA and Azide free	HA610047	Human, Mouse, Rat	WB,IHC
	E-Cadherin Recombinant Rabbit mAb	IRS030	Human	mIHC
Relative expression, KD	E-Cadherin Recombinant Rabbit mAb [SY0287]	ET1607-75	Human, Cynomolgus monkey	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP,mIHC
	iFluor™ 488 Conjugated E-Cadherin Recombinant Rabbit mAb [SY0287]	HA720159F	Human	IF-Cell,FC

Relative expression, KD	E-Cadherin Recombinant Rabbit mAb [SY0287] - BSA and Azide free	HA750133	Human, Cynomolgus monkey	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	E-Cadherin Mouse mAb [3-F9]	M1405-3	Human, Mouse	IHC-P
Relative expression	E-Cadherin Recombinant Rabbit mAb [PSH13-75] - BSA and Azide free	HA751488	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Cell,FC,IP
	E-Cadherin Recombinant Mouse mAb	IRS229MS	Mouse	mIHC
	E-Cadherin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS284RT	Mouse	mIHC
	iFluor™ 594 Conjugated E-Cadherin Recombinant Rabbit mAb [SY0287]	HA720164F	Human	IF-Tissue
Relative expression	E-Cadherin Rabbit pAb	ER63312	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IHC-Fr,FC
	E-Cadherin Recombinant Mouse mAb [A0-G11-2-R]	HA601143	Human, Mouse, Rat	WB,IHC,mIHC
	E-Cadherin Recombinant Mouse mAb	IRS031	Human	mIHC
KD	Fibronectin Recombinant Rabbit mAb [JF0582] - BSA and Azide free	HA750340	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P
	Fibronectin Mouse mAb [3G4]	RT1224	Human, Mouse, Rat	WB,IP,IF,IHC-P
	Fibronectin Recombinant Rabbit mAb [PSH14-22] - BSA and Azide free	HA751515	Human, Mouse	WB,IF-Cell,FC
KD	Fibronectin Recombinant Rabbit mAb [JF0582]	ET1702-25	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P
	Fibronectin Recombinant Rabbit mAb [PSH14-22]	HA723611	Human, Mouse	WB,IF-Cell,FC
	Fibronectin Recombinant Rabbit mAb [PSH14-21]	HA723610	Human, Mouse	WB,IHC-P,IF-Cell,FC
	Fibronectin Recombinant mAb	IRS144RB	Mouse	mIHC
	Fibronectin Recombinant Rabbit mAb [PSH14-21] - BSA and Azide free	HA751514	Human, Mouse	WB,IHC-P,IF-Cell,FC
	FOXC1 Recombinant Rabbit mAb [PSH03-89] - BSA and Azide free	HA750914	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	FOXC1 Recombinant Rabbit mAb [PSH03-89]	HA722060	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	FoxC1/2 Rabbit pAb	ER63452	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	GSC Recombinant Rabbit mAb [JE51-35]	HA722460	Human, Rat	WB
	GSK3 alpha Recombinant Rabbit mAb [JE33-27]	HA722217	Human, Mouse, Rat, Monkey	WB,IHC-P
KD	GSK3 alpha Rabbit pAb	R1312-1	Human	WB,IF-Cell,IHC-P,FC
	GSK3 alpha/beta Rabbit pAb	ER63648	Human, Mouse, Rat	WB,IHC-P
	GSK3 beta Rabbit pAb	R1510-4	Human, Mouse, Rat	IF-Cell,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	GSK3 beta Recombinant Rabbit mAb [SY28-03] - BSA and Azide free	HA750131	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IHC-Fr
KD	GSK3 beta Recombinant Rabbit mAb [SY28-03]	ET1607-71	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IHC-Fr
	Human E-Cadherin Recombinant Rabbit mAb [PSH09-02] - BSA and Azide free (Capture)	HA723055	Human	ELISA(Cap)
	Biotin Conjugated Human E-Cadherin Recombinant Rabbit mAb [SY0287] - BSA and Azide free (Detector)	HA723057B	Human	ELISA(Det),ELISA
	Human E-Cadherin Recombinant Rabbit mAb [SY0287] - BSA and Azide free (Detector)	HA723056	Human	ELISA(Det)
	Biotin Conjugated Human MMP-8 Recombinant Rabbit mAb [PSH03-34]	HA722624B	Human	ELISA(Det)
	Human MMP-8 Recombinant Rabbit mAb [PSH03-33] - BSA and Azide free (Capture)	HA721978	Human	ELISA(Cap)
	Human MMP-8 Recombinant Rabbit mAb [PSH03-34] - BSA and Azide free (Detector)	HA721979	Human	ELISA(Det)
	Human MMP-9 Recombinant Rabbit mAb [PSH04-23] - BSA and Azide free (Capture)	HA722106	Human	ELISA(Cap)
	Human MMP-9 Recombinant Rabbit mAb [PSH04-62] - BSA and Azide free (Capture)	HA722162	Human	ELISA(Cap)
	Biotin Conjugated Human MMP-9 Recombinant Rabbit mAb [PSH04-24]	HA722615B	Human	ELISA(Det)
	Human MMP-9 Recombinant Rabbit mAb [PSH04-24] - BSA and Azide free (Detector)	HA722107	Human	ELISA(Det)
	Human MMP-9 Recombinant Rabbit mAb [PSH04-63] - BSA and Azide free (Detector)	HA722163	Human	ELISA(Det)
	Human VE Cadherin Recombinant Rabbit mAb [PSH02-62] - BSA and Azide free (Detector)	HA721839	Human	ELISA(Det)
	Human VE Cadherin Recombinant Rabbit mAb [PSH02-61] - BSA and Azide free (Capture)	HA721838	Human	ELISA(Cap)
	Biotin Conjugated Human VE Cadherin Recombinant Rabbit mAb [PSH02-62] - Detector	HA722635B	Human	ELISA(Det)
	ILK Rabbit pAb	RT1330	Human, Mouse, Rat	WB,IP,IF
KD	Integrin alpha V Recombinant Rabbit mAb [SC56-07]	ET1610-15	Human, Mouse, Rat	WB,IF-Cell,FC,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	Integrin alpha V Recombinant Rabbit mAb [SC56-07] - BSA and Azide free	HA750201	Human, Mouse, Rat	WB,IF-Cell,FC,IHC-P
	Integrin Alpha V Rabbit pAb	ER1911-52	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,IP
	Integrin beta 3 Recombinant Rabbit mAb [JE22-64]	ET7110-61	Human, Mouse	WB,IHC-P
	Integrin beta 3 Rabbit pAb	ER1601-4	Human, Mouse, Rat	WB,IHC-P
KD	Integrin beta 3 Recombinant Rabbit mAb [SJ19-09]	ET1606-49	Human, Mouse	WB,IHC-P,IF-Tissue
Relative expression	Integrin beta 3 Recombinant Rabbit mAb [PSH09-27]	HA723085	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Integrin linked ILK Rabbit pAb	ER1803-78	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	Integrin linked ILK Recombinant Rabbit mAb [SC68-04]	ET1610-76	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC
	MMP-2 Mouse mAb [B11-E6]	EM1706-99	Human, Mouse	WB,FC
Relative expression	MMP-2 Recombinant Rabbit mAb [PSH11-26]	HA723308	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
Relative expression	MMP-2 Recombinant Rabbit mAb [PSH11-26] - BSA and Azide free	HA751405	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
	MMP-2 Recombinant Rabbit mAb [SI15-04]	ET1606-4	Human, Mouse, Rat	WB,IHC-P,IP,IF-Tissue
	MMP-8 Rabbit pAb	ER63980	Human, Mouse, Rat	WB,IHC-P
	MMP-8 Recombinant Rabbit mAb [JE59-38]	HA721686	Human, Mouse, Rat	IHC-P
Cell treatment	MMP-9 Recombinant Rabbit mAb [PSH06-33] - BSA and Azide free	HA751070	Human	WB,IHC-P
	MMP-9 Recombinant Mouse mAb [10A1]	EM1801-22	Human	ELISA
Cell treatment	MMP-9 Recombinant Rabbit mAb [PSH06-33]	HA722665	Human	WB,IHC-P
	MMP-9 Mouse mAb [2G3]	RT1401	Human	WB,IP,IF,IHC-P
	MMP-9 Recombinant Rabbit mAb	IRS172RB	Human	mIHC
Cell treatment, KD	MMP-9 Recombinant Rabbit mAb [JA80-73] - BSA and Azide free	HA750420	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,FC
Cell treatment, KD	MMP-9 Recombinant Rabbit mAb [JA80-73]	ET1704-69	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,FC
	Mouse MMP-9 Recombinant Rabbit mAb [PSH15-81] - BSA and Azide free (Detector)	HA723787	Mouse	ELISA(Det)
	Mouse MMP-9 Recombinant Rabbit mAb [PSH15-80] - BSA and Azide free (Capture/Detector)	HA723785	Mouse	ELISA(Det),ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	HRP Conjugated Mouse MMP-9 Recombinant Rabbit mAb [PSH15-81] - Detector	HA723788H	Mouse	ELISA(Det),ELISA
	Biotin Conjugated Mouse MMP-9 Recombinant Rabbit mAb [PSH15-80] - Detector	HA723786B	Mouse	ELISA(Det),ELISA
	Biotin Conjugated Mouse MMP-9 Recombinant Rabbit mAb [PSH15-81] - Detector	HA723788B	Mouse	ELISA(Det),ELISA
	Mouse MMP-9 Recombinant Rabbit mAb [PSH15-79] - BSA and Azide free (Capture)	HA723784	Mouse	ELISA(Cap)
	Mouse/Rat E-Cadherin Recombinant Rabbit mAb [PSH17-25] - BSA and Azide free (Detector)	HA723912	Mouse, Rat	ELISA(Det)
	HRP Conjugated Mouse/Rat E-Cadherin Recombinant Rabbit mAb [PSH17-25] - Detector	HA723912H	Mouse, Rat	ELISA(Det),ELISA
	Biotin Conjugated Mouse/Rat E-Cadherin Recombinant Rabbit mAb [PSH17-25] - Detector	HA723912B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat E-Cadherin Recombinant Rabbit mAb [PSH17-24] - BSA and Azide free (Capture)	HA723911	Mouse, Rat	ELISA(Cap)
	N Cadherin Recombinant Antibody [SY02-46] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610298	Human, Mouse, Rat	IHC-Fr,IHC-P
Relative expression, KD	N Cadherin Recombinant Rabbit mAb [SY02-46]	ET1607-37	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue
Relative expression, KD	N Cadherin Recombinant Rabbit mAb [SY02-46] - BSA and Azide free	HA750119	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue
	N Cadherin Recombinant Antibody [SY02-46] - Mouse IgG1 (Chimeric)	HA601464	Human, Mouse, Rat	IHC-Fr,IHC-P
	N-Cadherin Mouse mAb [A10-4]	M1304-1	Human, Mouse	WB,IHC-P
	N-Cadherin Rabbit pAb	ER0503	Human, Mouse, Rat	WB,IHC,IF
	NF-kB p65 Rabbit pAb	ER0815	Human, Mouse, Rat, Zebrafish	WB,IHC-P,FC
KO, KD, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [SZ10-04] - BSA and Azide free	HA750058	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KO, KD, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [SZ10-04]	ET1603-12	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	NF-kB p65 Rabbit pAb	HA500402	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KO, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [PSH0-27]	HA721307	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
KO, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [PSH0-27] - BSA and Azide free	HA750601	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
	Occludin Recombinant Rabbit mAb	IRS183RB	Mouse	mIHC
KD	Occludin Recombinant Rabbit mAb [JJ091-08]	ET1701-76	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IP,FC,IF-Cell,IHC-P,IHC-Fr,IF-Tissue
	Occludin Rabbit pAb	R1510-33	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Pan Cadherin Recombinant Rabbit mAb	RE6045	Human, Mouse, Rat	WB,ICC,IHC-P,FC,IP
	Pan-Cadherin Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS285MS	Human	mIHC
	Pan-Cadherin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS287RT	Human	mIHC
	Pan-Cadherin Recombinant Rabbit mAb [ST54-01] - BSA and Azide free	HA750190	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,IHC-Fr
	Pan-Cadherin Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS286MS	Mouse	mIHC
	Pan-Cadherin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS288RT	Mouse	mIHC
	Pan-Cadherin Recombinant Rabbit mAb [ST54-01]	ET1609-70	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,IHC-Fr
	Pan-Cadherin Recombinant Rabbit mAb	IRS185RB	Mouse	mIHC
Cell treatment	Phospho-Beta Catenin (S33 + S37) Recombinant Rabbit mAb [JE59-59]	HA721580	Human, Rat, Mouse	WB
Cell treatment	Phospho-Beta Catenin (S552) Recombinant Rabbit mAb [PSH08-72]	HA723024	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Cell treatment	Phospho-Beta Catenin (S552) Recombinant Rabbit mAb [PSH08-72] - BSA and Azide free	HA751237	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Cell treatment	Phospho-Beta Catenin (T41/S45) Recombinant Rabbit mAb [JE54-02]	HA722316	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-GSK3 (alpha + beta) (Y216 + Y279) Recombinant Rabbit mAb [SY26-05] - BSA and Azide free	HA750127	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,IF-Tissue,FC
Cell treatment	Phospho-GSK3 (alpha + beta) (Y216 + Y279) Recombinant Rabbit mAb [SY26-05]	ET1607-54	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,IF-Tissue,FC
Cell treatment	Phospho-GSK3 alpha (S21) Recombinant Rabbit mAb [PSH09-11] - BSA and Azide free	HA751265	Human, Mouse	WB,IHC-P,Dot Blot
Cell treatment	Phospho-GSK3 alpha (S21) Recombinant Rabbit mAb [PSH09-11]	HA723069	Human, Mouse	WB,IHC-P,Dot Blot
Cell treatment	Phospho-GSK3 beta (S9) Recombinant Rabbit mAb [SY02-71]	ET1607-60	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-GSK3 beta (S9) Recombinant Rabbit mAb [SY02-71] - BSA and Azide free	HA750128	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
Cell treatment	Phospho-NF-kB p65 (S529) Recombinant Rabbit mAb [SP07-00]	ET1604-27	Human	WB,IHC-P,FC,IP
Cell treatment	Phospho-NF-kB p65 (S536) Recombinant Rabbit mAb [PSH10-58] - BSA and Azide free	HA751361	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-NF-kB p65 (S536) Recombinant Rabbit mAb [PSH10-58]	HA723223	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-Smad2 (S250) Recombinant Rabbit mAb [SD207-1]	ET1612-32	Human, Mouse, Rat	WB,FC,IP,IHC-P
Cell treatment	Phospho-Smad2 (S250) Recombinant Rabbit mAb [SD207-1] - BSA and Azide free	HA750291	Human, Mouse, Rat	WB,FC,IP,IHC-P
Cell treatment	Phospho-Smad2 (S255) Recombinant Rabbit mAb [JF0882]	ET1702-34	Human, Mouse, Rat	WB,IHC-P,IP
Cell treatment	Phospho-SMAD2 (S465+467) Recombinant Rabbit mAb [PSH05-76]	HA722443	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-SMAD2 (S465+467) Recombinant Rabbit mAb [PSH05-76] - BSA and Azide free	HA751026	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-Smad3 (S423 + S425) Recombinant Rabbit mAb [ST0493]	ET1609-41	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	Phospho-Smad3 (S423 + S425) Recombinant Rabbit mAb [ST0493] - BSA and Azide free	HA750172	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	Phospho-Vimentin (S39) Recombinant Rabbit mAb [JE43-26]	HA721442	Human, Mouse, Rat	WB,IF-Cell
	R-cadherin Rabbit pAb	ER64773	Human, Mouse	IHC-P,IF-Cell,WB
	RhoA Recombinant Rabbit mAb [SN0612] - BSA and Azide free	HA750242	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	RhoA Recombinant Rabbit mAb [SN0612]	ET1611-10	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Rho A+B+C Recombinant Rabbit mAb [SR38-00]	ET1602-23	Human, Mouse, Rat	WB,IF-Cell,FC
	Rho B Rabbit pAb	ER64810	Human, Mouse, Rat	WB
	SIP1 Rabbit pAb	ER64964	Human, Mouse, Rat	WB,IHC-P
Relative expression	SLUG Recombinant Rabbit mAb [PSH07-40]	HA722828	Human	WB,IF-Cell
Relative expression	SLUG Recombinant Rabbit mAb [PSH07-40] - BSA and Azide free	HA751161	Human	WB,IF-Cell
	Smad2 Mouse mAb [C9-B0]	EM1701-44	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Smad2 Rabbit pAb	R1511-29	Human, Mouse, Rat	IF-Cell,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KO	Smad2 Recombinant Rabbit mAb [SP06-05]	ET1604-22	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC
KO	Smad2 Recombinant Rabbit mAb [SP06-05] - BSA and Azide free	HA750073	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC
	Smad2/3 Rabbit pAb	ER64981	Human, Mouse, Rat, Monkey	WB,IHC-P
	Smad2/3 Mouse mAb [4D3]	RT1566	Human, Mouse, Rat	WB,IP,IF,IHC-P
	Smad2/3 (Acetyl K19) Rabbit pAb	ER60061	Human, Mouse, Rat	IHC-P,WB
	Smad3 Rabbit pAb	R1510-25	Human	WB,IF-Cell,IHC-P,FC
	Smad3 Mouse mAb [8E1]	RT1567	Human, Mouse, Rat	WB,IP,IF,IHC-P
KO	Smad3 Recombinant Rabbit mAb [SY25-01] - BSA and Azide free	HA750123	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KO	Smad3 Recombinant Rabbit mAb [SY25-01]	ET1607-41	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Smad3 Rabbit pAb	R1510-24	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC
	SNAI2 Mouse mAb [D5-B6]	EM1706-65	Human	WB,IF-Cell,IHC-P,FC
	SNAIL Rabbit pAb	ER1706-22	Human, Mouse, Rat	WB,IF-Cell
KD	SOX10 Recombinant Rabbit mAb [SD204-04]	ET1612-79	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Relative expression	SOX10 Recombinant Rabbit mAb [PDH0-04]	HA721241	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC
Relative expression	SOX10 Recombinant Rabbit mAb [PDH0-03]	HA721240	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Cell,IF-Tissue,FC,IHC-Fr,IP
	SOX10 Recombinant Antibody [PDH0-03] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610253	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,IF-Tissue
	SOX10 Recombinant Rabbit mAb [SN20-41]	ET1611-89	Human, Mouse, Rat	WB
	SOX10 Rabbit pAb	HA500503	Human, Mouse	WB,IHC-P
	SOX10 Mouse mAb [1-11]	M1311-1	Mouse	WB,IHC-P
	SOX10 Recombinant Antibody [PDH0-03] - Rat IgG1 (Chimeric)	HA601413	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,IF-Tissue
Relative expression	SOX10 Recombinant Rabbit mAb [PDH0-05]	HA721239	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP
	Twist Mouse mAb [1G3]	RT1635	Human, Mouse, Rat	WB,IP,IF,FC
	Twist1 Recombinant Rabbit mAb [PSH11-98] - BSA and Azide free	HA751431	Human, Mouse, Rat	WB
	Twist1 Recombinant Rabbit mAb [PSH11-98]	HA723411	Human, Mouse, Rat	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	TWIST1 Mouse mAb [E6-F4]	EM1710-24	Human, Mouse	WB,IF-Cell,IHC-P,FC
KD	VE Cadherin Recombinant Rabbit mAb [PSH03-15] - BSA and Azide free	HA750856	Human	WB,IF-Cell,FC
KD	VE Cadherin Recombinant Rabbit mAb [PSH03-15]	HA721943	Human	WB,IF-Cell,FC
KO	Vimentin Mouse mAb [D4-B11]	EM0401	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue,IHC-Fr,IP
KO, Relative expression	Vimentin Recombinant Rabbit mAb [SC60-05]	ET1610-39	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue,IF-Cell,FC,IP
KO, Relative expression	Vimentin Recombinant Rabbit mAb [SC60-05] - BSA and Azide free	HA750215	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue,IF-Cell,FC,IP
	iFluor™ 488 Conjugated Vimentin Recombinant Rabbit mAb [SC60-05]	HA720166F	Human, Mouse, Rat	IF-Cell,IF-Tissue,FC
Relative expression	Vimentin Recombinant Mouse mAb [A6-C1-R]	HA601251	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Vimentin Recombinant Antibody - Rat IgG1 (Chimeric)	HA601564	Human, Mouse, Rat	IHC,ICC
	Vimentin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS214RT	Human	mIHC
	Vimentin Recombinant Rabbit mAb	IRS060	Mouse	mIHC
KO	Vimentin Rabbit pAb	R1308-6	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Relative expression, KO	Vimentin Mouse mAb [A6-C1]	M1412-1	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	iFluor™ 647 Conjugated Vimentin Recombinant Rabbit mAb [SC60-05]	HA720165F	Human, Mouse, Rat	IF-Cell,IF-Tissue,FC
	Vimentin Recombinant Rabbit mAb [PDH0-01]	HA721174	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IF-Tissue,mIHC
Relative expression	Vimentin Recombinant Mouse mAb [A6-C1-R] - BSA and Azide free	HA610126	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Vimentin Recombinant Rabbit mAb	IRS044	Human	mIHC
	Vimentin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS263RT	Mouse	mIHC
	ZEB Recombinant Rabbit mAb	IRS245RB	Mouse	mIHC
	ZEB Recombinant Rabbit mAb [PSH0-61]	HA721438	Human, Mouse, Rat	WB,IHC-P,FC
	ZEB Recombinant Rabbit mAb	IRS216RB	Human	mIHC
	ZEB Recombinant Rabbit mAb [PSH0-61] - BSA and Azide free	HA750638	Human, Mouse, Rat	WB,IHC-P,FC
	ZEB1 Mouse mAb [A2-D8]	EM1710-27	Human	WB,IHC-P,FC
	ZEB1 Recombinant Rabbit mAb [JE57-64]	HA721318	Human	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

