

纳米抗体



产品优势



结构简单，组织穿透力强



亲和力高，特异性强




批次间稳定性高



生物相容性高


检测数据精选

DYKDDDDK tag, AlpHcAbs® Mouse antibody(HRP) (HA710083)




1. HeLa+N-flag
2. L-929+C-flag

GST, AlpHcAbs® Rabbit antibody(HRP) (HA710105)




1. 293T+OE-GST
2. 293T+OE-GST-Histone H3.1

HA tag Nanoselector Agarose (HA710007)




IP of C-HA protein with HA tag nanoselector
Lane 1: Input
Lane 2: Flow-through
Lane 3: PH=3 Elution
Lane 4: PH=2 Elution
Lane 5: Eluted beads

Rabbit IgG(H+L), AlpSdAbs® VHH(HRP) (HA710051)




1. 293T+OE-Claudin18.2-Myc unbioled
2. 293T+OE-Myc-Histone H3.1

Mouse IgG kappa, AlpHcAbs® Goat antibody(HRP) (HA710127)



A titer ELISA of mouse IgG-Anti-Mouse IgG kappa, AlpHcAbs® Goat antibody(HRP)



ELISA analysis of specificity for IgG-Anti-Mouse IgG kappa, AlpHcAbs® Goat antibody(HRP)

杭州华安生物技术有限公司

Hangzhou HuaAn Biotechnology Co.,Ltd.

地址：杭州市钱塘新区白杨街道6号大街452号高科技企业

孵化园区2号楼16层D区

邮编：310012

邮箱：sales@huabio.cn

电话：0571-88062880

传真：0571-88060400



聚焦抗体，专业服务！www.huabio.cn

售后技术支持电话：0571-86052305



优质产品推荐

扫码查看更多
“纳米抗体”产品信息



Target name	Product name	Catalog#	Tested applications
Nanoselector			
C-His	C-His tag Nanoselector Agarose	HA710174	IP,CHIP,MS,Purification
GFP	GFP Nanoselector Agarose	HA710003	IP,Co-IP,Mass Spectrometry,ChIP
	GFP Nanoselector Magnetic beads	HA710169	IP,CHIP,MS,Purification
GST	GST Nanoselector Agarose	HA710006	IP,Co-IP,Mass Spectrometry,ChIP
HA	HA tag Nanoselector Agarose	HA710007	IP,Co-IP,Mass Spectrometry,ChIP
Halo	Halo Nanoselector Agarose	HA710173	IP,CHIP,MS,Purification
MBP	MBP Nanoselector Agarose	HA710171	IP,CHIP,MS,Purification
mNeogreen	mNeogreen Nanoselector Agarose	HA710172	IP,CHIP,MS,Purification
Mouse IgG	Mouse IgG Nanoselector Agarose	HA710175	IP,CHIP,MS,Purification
mWasabi	mWasabi Nanoselector Agarose	HA710164	IP,CHIP,MS,Purification
Rabbit IgG	Rabbit IgG Nanoselector Agarose	HA710167	IP,CHIP,MS,Purification
RFP	RFP Nanoselector Agarose	HA710004	IP,Co-IP,Mass Spectrometry,ChIP
	RFP Nanoselector Magnetic beads	HA710168	IP,CHIP,MS,Purification
SNAP	SNAP tag Nanoselector Agarose	HA710165	IP,CHIP,MS,Purification
TagFP	TagFP Nanoselector Agarose	HA710170	IP,CHIP,MS,Purification
turboGFP	turboGFP Nanoselector Agarose	HA710005	IP,Co-IP,Mass Spectrometry,ChIP
V5	V5 tag Nanoselector Agarose	HA710166	IP,CHIP,MS,Purification
	V5 tag Nanoselector Magnetic beads	HA710002	IP,ChIP
Primary Antibodies			
BCMA	BCMA, AlpHcAbs® Human	HA710324	ELISA,FC
CD154	CD154, AlpHcAbs® Human	HA710332	ELISA,FC
CD16a	CD16a, AlpHcAbs® Human	HA710339	ELISA,FC
CD19	CD19, AlpHcAbs® Human	HA710342	ELISA,FC
CD33	CD33, AlpHcAbs® Human	HA710344	ELISA,FC
CD38	CD38, AlpHcAbs® Human	HA710331	ELISA,FC
CD7	CD7, AlpHcAbs® Human	HA710330	ELISA,FC
CD70	CD70, AlpHcAbs® Human	HA710335	ELISA,FC
Claudin-18	Claudin-18, AlpHcAbs® Human	HA710341	ELISA,FC
Claudin-18.2	Claudin-18.2, AlpHcAbs® Human	HA710340	ELISA,FC
CTLA4	CTLA4, AlpHcAbs® Human	HA710327	ELISA,FC
DLL3	DLL3, AlpHcAbs® Human	HA710326	ELISA,FC
DYKDDDDK	DYKDDDDK tag, AlpHcAbs® Mouse antibody(HRP)	HA710083	WB,ELISA
	DYKDDDDK tag, AlpHcAbs® Rabbit	HA710226	WB,ICC/IF,ELISA,IP,FC
FGFR2(IIIb)	FGFR2(IIIb), AlpHcAbs® Human	HA710328	ELISA,FC
FOLR1	FOLR1, AlpHcAbs® Human	HA710333	ELISA,FC