



细胞器标志物研究


线粒体标志物




IF-Cell staining with NIH/3T3 cells
TOMM20 rabbit mAb | ET1609-25
Beta tubulin | M1305-2



IF-Cell staining with HeLa cells
OPA1 rabbit mAb | HA722673
Beta tubulin | M1305-2




IF-Cell staining with HeLa cells
VDAC1 rabbit mAb | ET1601-20
Beta tubulin | M1305-2




IF-Cell staining with HeLa cells
Mitofillin rabbit mAb | HA721636
Beta tubulin | M1305-2


其他细胞器标志物




IF-Cell staining with MCF7 cells
GM130 rabbit mAb | HA721282
Beta tubulin | M1305-2



IF-Cell staining with HeLa cells
Calreticulin rabbit mAb | ET1608-60
Beta tubulin | M1305-2




IF-Cell staining with C2C12 cells
LC3B rabbit mAb | ET1701-65
Beta tubulin | M1305-2




IF-Cell staining with HeLa cells
Lamin A/C rabbit mAb | ET7110-12
Beta tubulin | M1305-2

客户文献引用案例



(Nucleic Acids Res)

IF: 16.1004 PMID:34570240
产品名称: Histone H3 Mouse mAb | M1309-1
科研单位: 浙江大学



(Sci Immunol)

IF: 30.63 PMID:35239449
产品名称: TOMM20 rabbit mAb | ET1609-25
科研单位: 浙江大学

优质产品推荐

扫码了解更多
“细胞器标志物” 抗体产品信息



Organelle type	Product name	Catalog#	Reactivity	Tested applications	Antibody type
Mitochondrion	ACAA2	HA720006	H,Ms,R	WB,IHC-P	Rabbit mAb
	ACAT-1	HA721093	H,Ms,R	WB,IHC-P,IF-Cell	Rabbit mAb
	ALDH2	HA601277	H,R	WB,IF-Cell,IHC-P	Mouse mAb
	Mitofilin	HA721636	H,Ms,R	WB,IHC-P,IF-Cell,IF-Tissue	Rabbit mAb
	NDUFV2	HA721869	H,Ms,R	WB,IHC-P,IF-Cell	Rabbit mAb
	TOMM20	ET1609-25	H,Ms,R	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC	Rabbit mAb
	VDAC1	ET1601-20	H,Ms,R	WB,IF-Cell,IHC-P,FC,IF-Tissue	Rabbit mAb
	DRP1	HA500487	H,Ms,R	WB,IF-Cell,FC	Rabbit pAb
	OPA1	HA722673	H,Ms,R	WB, IHC-P, IF-Cell, FC, IP	Rabbit mAb
	OPA1	ET1705-9	H,Ms,R	WB, IF-Tissue, IHC-P, FC	Rabbit mAb
	ATP5A1	ET1703-53	H,Ms,R	WB,IF-Cell,IF-Tissue,IHC-P	Rabbit mAb
	CPS1	ET7107-69	H,Ms,R	WB,IF-Cell,IHC-P,IP	Rabbit mAb
	HSP60	ET1609-45	H,Ms,R	WB,IHC-P,IF-Cell,IP	Rabbit mAb
	MTCO1	HA722838	H,Ms,R	WB,IF-Cell,IHC-P,FC,IF-Tissue	Rabbit mAb
	Prohibitin	ET1602-31	H,Ms,R,Z	WB,IF-Cell,IHC-P,IP,FC	Rabbit mAb
	Cytochrome C	ET1610-60	H,Ms,R	WB,IF-Tissue,IHC-P,IP,IF-Cell	Rabbit mAb
Endoplasmic Reticulum	Calreticulin	ET1608-60	H,Ms,R	WB,IF-Cell,IHC-P,FC,IP	Rabbit mAb
	FADS1	ET7111-19	H,Ms,R	WB,IF-Cell,IHC-P,FC	Rabbit mAb
	Glucosidase 2 subunit beta	HA720042	H	WB,IHC-P,IF-Cell,IF-Tissue	Rabbit mAb
	GRP94	M1701-12	H,Ms,R	WB,IHC-P,IF-Cell	Mouse mAb
	HSPA14	ET1612-93	H,Ms,R	WB,IF-Cell,IF-Tissue,IHC-P,IP	Rabbit mAb
Autophagosome & Lysosome	ATG5	ET1611-38	H,Ms,R,Mk	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC	Rabbit mAb
	Beclin 1	HA721216	H,Ms,R	WB,IHC-P,IF-Cell	Rabbit mAb
	LAMP1	HA722827	H,Ms,R	WB,IHC-P,IF-Tissue,IHC-Fr	Rabbit mAb
	LAMP1	ET1701-94	H	WB,IHC-P,IF-Tissue	Rabbit mAb
	LAMP2	EM1901-62	H	WB,IF-Cell,IHC-P,FC	Mouse mAb
	MAP1LC3A	ET1609-26	H,Ms,R	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC	Rabbit mAb
	LC3B	ET1701-65	H,Ms,R	WB,IF-Cell,IHC-P,IF-Tissue,IP,mIHC	Rabbit mAb
Golgi Apparatus	ALIX	ET1705-74	H,Ms,R	WB,IHC-P,FC,IF-Cell	Rabbit mAb
	Fukutin	HA720031	H,Ms,R	WB,IHC-P	Rabbit mAb
	Furin	ET7107-37	H,Ms,R	WB,IHC-P	Rabbit mAb
	Giantin	HA722052	H,Ms,R,Mk	WB,IF-Cell,IHC-P,FC	Rabbit mAb
	GM130	HA721282	H,Ms,R	WB,IHC-P,IF-Cell	Rabbit mAb
Nucleus	Histone H3	ET1701-64	H,Ms,R	WB,IF-Cell,IF-Tissue,IHC-P,ChIP,IP	Rabbit mAb
	Histone H3	M1309-1	H,Ms,R,Z	WB,IF-Cell,IHC-P,IF-Tissue	Mouse mAb
	Lamin A/C	ET7110-12	H,Ms,R	WB,IF-Cell,IHC-P,FC	Rabbit mAb
	Lamin B1	ET1606-27	H,Ms,R	WB,IF-Tissue,IHC-P,CUT&Tag-seq	Rabbit mAb

*H: Human; Ms: Mouse; R: Rat; Z: Zebrafish; Mk: Monkey

杭州华安生物技术有限公司

Hangzhou HuaAn Biotechnology Co.,Ltd.

地址：杭州市钱塘新区白杨街道6号大街452号高科技企业
孵化园区2号楼16层D区
邮编：310012 邮箱：sales@huabio.cn
电话：0571-88062880 传真：0571-88060400



聚焦抗体，专业服务！ www.huabio.cn
售后技术支持电话：0571-86052305

