

Apoptosis相关抗体

在Cell Biology—Apoptosis领域，家族蛋白主要分为凋亡促进蛋白和凋亡抑制蛋白两大类，它们通过调控caspase级联反应和线粒体途径精确控制细胞凋亡过程。

优质抗体推荐

Product name	Catalog#	Reactivity	Tested applications	Antibody type
Granzyme K	HA723826	H	WB,IHC-P	Rabbit mAb
AMID	HA723767	H,Ms,R	WB,IF-Cell,IHC-P,FC,IP	Rabbit mAb
Acid sphingomyelinase	HA723768	H,Ms,R,Mk	WB,IHC-P,IP	Rabbit mAb
WWOX	HA723750	H,Ms,R,Mk	WB,IP	Rabbit mAb
WWOX	HA723749	H,Ms,R,Mk	WB,IF-Cell,IP	Rabbit mAb
JMJD6	HA723703	H,Ms,R,Mk	WB,IF-Cell,IHC-P,FC,IP,CHIP	Rabbit mAb
MFF	HA723630	H,Ms,R,Mk	WB,IHC-P,IF-Cell,FC,IP	Rabbit mAb
DR3 / LARD	HA723583	H,Ms,R	WB,IF-Cell	Rabbit mAb

*H: Human; Ms: Mouse; R: Rat; Mk: Monkey

产品特点

高灵敏性

高稀释比，信号强，背景低。

01



高特异性

高低表达和诱导材料验证，确保抗体的特异性。

02



种属/应用全面

满足不同种属、不同应用的实验需求。

03



杭州华安生物技术有限公司

Hangzhou HuaAn Biotechnology Co.,Ltd.


地址：杭州市钱塘新区白杨街道6号大街452号高科技企业
孵化园区2号楼16层D区
邮编：310012 邮箱：sales@huabio.cn
电话：0571-88062880 传真：0571-88060400




聚焦抗体，专业服务！ www.huabio.cn
售后技术支持电话：0571-86052305




Acid sphingomyelinase 重组兔单抗 | HA723768




WB analysis on different lysates
(1/5,000 dilution)
The expression is up-regulated
in response to TPA then LPS and
BFA treatment.



IP analysis on HepG2 lysates
(1/5,000 dilution)




IHC-P staining with human liver
tissue
(1/200 dilution)




IHC-P staining with mouse kidney
tissue
(1/100 dilution)

WVVOX重组兔单抗 | HA723749




WB analysis on different lysates
(1/10,000 dilution)




IF-cell analysis of NIH/3T3 cells
(1/100 dilution)

Granzyme K重组兔单抗 | HA723826




WB analysis on different lysates
(1/10,000 dilution)
Relative expression (RE)




IHC-P staining with human NK/T
lymphoma tissue
(1/1,000 dilution)


AMID重组兔单抗 | HA723767




WB analysis on different lysates
(1/5,000 dilution)



IF-Tissue analysis of rat liver tissue
(1/200 dilution)



IHC-P staining with mouse liver
tissue
(1/200 dilution)



FC analysis of PC-12 cells
(1/1,000 dilution)