







扫码了解更多
“HIF-1 α信号通路图”
相关研究产品信息

Product name	Catalog#	Reactivity	Application
AKT1/2/3	HA721870	H,Ms,R,Mk	WB,IF-Cell,FC,IP
Phospho-AKT (S473)	ET1607-73	H,Ms,R,D	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
Aldolase A/B/C	HA721733	H,Ms,R,Cmk,P	WB,IHC-P,IHC-Fr,IF-Tissue,IF-Cell
AMPK alpha 1	ET1608-40	H,Ms,R	WB,IF-Cell,IF-Tissue,FC,IP,IHC-Fr,IHC-P
Phospho-AMPK alpha 1 (S496)	ET1612-72	H,Ms	WB,IP,IHC-Fr
AMPK alpha 2	HA600078	H,Ms,R	WB,IF-Cell,IHC-P
Phospho-AMPK alpha 2 (S345)	ET1701-37	H,Ms,R	WB,IF-Tissue,IHC-P,FC
Phospho-AMPK alpha 2 (S491)	HA722567	H	WB
Angiotensin 1	ET1611-28	H,Ms,R	WB
Angiotensin 2	ET1705-6	H,Ms,R	WB,IHC-P,FC
ANP	ET1705-75	H	WB,IP
Bcl-2	ET1702-53	H,Ms,R	WB,IF-Cell,IF-Tissue,IHC-P,FC
Bcl-2	HA721235	H	WB,IHC-P,IF-Cell,FC
BNIP3	ET1704-08	H,Ms,R	WB,IF-Cell,IF-Tissue,IHC-P
CaMKII	ET1608-47	H,Ms,R,Cmk,P	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr,FC
Phospho-CaMKII (T286)	HA721794	H,Ms,R	WB
CCL2	HA500267	Ms,H	WB,FC
CD73	HA601011	H,Ms,R	WB,IHC-P,IF-Tissue,FC,IF-Cell
CDKN1A/P21	HA500156	H,Ms	WB,IHC-P
CREBBP	HA722297	H,Ms	IHC-P
EGFR	ET1603-37	H,Ms,Mk,R	WB,IF-Cell,IP,FC,mIHC
Phospho-EGFR (Y1068)	ET1612-30	H	WB,IF-Cell,IF-Tissue,IP,FC
eIF4E	HA721192	H,Ms,R	WB,IHC-P,IF-Cell,FC
Phospho-eIF4E (S209)	ET1608-66	H,Ms,R	WB,IF-Cell,IHC-P,IP,FC
eIF4EBP1	ET1701-83	H,Ms,R	WB,IF-Cell,IF-Tissue,IHC-P,IP
Phospho-eIF4EBP1 (S65)	HA723106	H	WB,IF-Cell,IP
Phospho-eIF4EBP1 (T70)	HA723222	H,Ms	WB,IHC-P,IP
ERK1/2	ET1601-29	H,Ms,R,Z	WB,IF-Cell,IF-Tissue,IP,FC,IHC-P,IHC-Fr
Phospho-Erk1 (T202)+Erk2 (T185)	ET1603-22	H,Ms,R	WB,IF-Cell,IHC-P,FC,IF-Tissue
Glucose Transporter GLUT1	ET1601-10	H,Ms,R	WB,IF-Cell,IF-Tissue,IHC-P,FC
Heme Oxygenase 1 (HO-1)	ER1802-73	H,Ms	WB,IF-Cell,IHC-P,FC
Heme Oxygenase 1 (HO-1)	ET1604-45	H	WB,IHC-P,IP,FC,IF-Tissue
HER2 / ErbB2	HA721178	H	IHC-P,WB,IF-Cell,FC,mIHC
Hexokinase 1	ET1609-28	H,Ms,R	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
HIF-1 alpha	HA721997	H,Ms,R	WB,IF-Cell
HIF-1 beta	ET1609-67	H,Ms,R	WB,IHC-P
HKDC1	ER1803-32	H,Ms,R	WB,IF-Cell,IHC-P,FC
IFNGR1	HA721952	H,Ms,R	WB,IHC-P
IL-6	R1412-2	H,Ms,R	WB,IF-Cell,IHC-P,FC
IL-6	EM1701-45	Ms,R	WB,IF-Cell
IL-6R	HA500438	H	WB,IF-Cell,IHC-P
iNOS	ER1706-89	H,Ms,R,C	IF-Cell,FC,WB,IHC-Fr
Insulin	EM1701-33	H	IHC-P,ELISA
Insulin	EM80714	H,Ms,R	IHC-P,ELISA,IF-Tissue,IHC-Fr
Interferon gamma	ET1703-17	H	WB,FC
Lactate Dehydrogenase	ET1608-57	H,Ms,R,Z	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
LDHA	ER00702	H,Ms,R	WB,IHC-P,FC,IF-Cell
MEK1/2	ET1602-3	H,Ms,R,Z	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Phospho-MEK1/2 (S218 + S222)	ET1609-50	H,Ms,R	WB,IF-Cell,IF-Tissue,IHC-P,IP



Product name	Catalog#	Reactivity	Application
mTOR	ET1608-5	H,Ms,R	WB,IHC-P,IP,IF-Cell
Phospho-mTOR (S2448)	HA600094	H,Ms,R	WB,IHC-P
Phospho-mTOR (S2481)	HA721632	H	WB
NF-kB p105 / p50	ET1603-18	H,Ms,R	WB,IHC-P,ICC,FC,IF-Tissue
NF-kB p65	ET1603-12	H,Ms,R,Z	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Phospho-NF-kB p65 (S529)	ET1604-27	H	WB,IHC-P,FC,IP
Phospho-NF-kB p65 (S536)	HA723223	H,Ms	WB,IF-Cell
NOX2 / gp91phox	HA723224	H,Ms,R	WB,IF-Cell,IHC-P,IHC-Fr,IP
NSE	ET1610-96	H,Ms,R,Z	WB,IF-Cell,IHC-P
p27 KIP 1	ET1608-61	H,Ms,R	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
p70 S6 Kinase	HA722520	H,Ms,R	WB,IHC-P,IP
Phospho-p70 S6 Kinase (T389)	HA721803	H,Ms,R	WB,IF-Cell,IHC-P
PDHA1	ET1702-75	H,Ms,R	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Phospho-PDHA1 (S293)	HA721443	H,Ms,R	WB,IHC-P,IP
PDHB	ET1705-57	H,Ms,R	WB,IHC-P,FC,IP
PDK1	ET1704-66	H,Ms,R	WB,IHC-P,IP,FC,IF-Cell,IF-Tissue
PFKFB3	ET1705-66	H,Ms,R	WB,IF-Cell,IHC-P,FC,IP,IF-Tissue
PFKM	ET1706-97	H,Ms,R	WB,IF-Cell,IF-Tissue,IHC-P,FC
PFKP	HA500472	H,Ms,R	WB,IHC-P
PGK1	ET1609-63	H,Ms,R	WB
PHD1 / prolyl hydroxylase	ET1604-5	H,Ms,R	WB,IF-Cell,IF-Tissue,IHC-P,FC
Phospho-PKC alpha (T638)	ET1702-17	H,Ms,R	WB,IF-Cell,IF-Tissue,IHC-P,IP
Phospho-RPS6 (S235 + S236)	HA721275	H,Ms,R	WB,IHC-P,IF-Cell
Phospho-STAT3 (S727)	ET1607-39	H,R,Ms	WB,IF-Cell,IF-Tissue,IHC-P,IP
Phospho-STAT3 (Y705)	ET1603-40	H,Ms,R	WB,IHC-P,IP
PI 3 Kinase p55 gamma	ET1602-27	H,Ms,R	WB,IF-Cell,IF-Tissue,IHC-P
PI 3 Kinase p85 alpha	ET1608-70	H,Ms,R	WB,IF-Cell,IF-Tissue,IHC-P,FC
Phospho-PI3K p85 (Y467) + PI3K p55 (Y199)	HA721673	H,Ms	WB,IF-Cell
PI3 Kinase p110 beta	HA722474	H,Ms,R	WB,IHC-P,IF-Tissue
PI3 Kinase p110α	ET1606-36	H,Ms	WB,IF-Cell,IF-Tissue,IP,IHC-P
PI 3 Kinase p110 delta	ET1705-46	H,Ms,R	WB,IHC-P
PI3-kinase p110 subunit beta	ET1610-36	H	WB,IP,IHC-P,FC
PKC alpha	ET1608-15	H,Ms,R	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
PLCG 2	HA721477	H,Ms,R	WB,IHC-P,IF-Cell,FC
RBX1	HA722779	H,Ms,R	WB,IHC-P,IF-Cell,IP,IF-Tissue
RPS6	HA600084	H,Ms,R	WB,IF-Cell,IHC-P,FC
Smad4	ET1604-12	H,Ms,R	WB,IF-Cell,IF-Tissue,IHC-P,FC
SOD2 (acetyl K68)	HA722463	H,Ms,R	WB,IHC-P
STAT3	ET1605-45	H,Ms,R	WB,IF-Cell,IF-Tissue,IHC-P,FC
STAT3	ET1607-38	H,Ms,R,Z	WB,IF-Cell,IHC-P,FC
TGF beta 1	HA721143	H,Ms,R	WB,IHC-P,IF-Tissue
TIE2	HA722732	H	WB,FC
TLR4	ER1706-43	H,Ms,R	IF-Cell,IHC-P,FC
Transferrin Receptor (CD71)	ET1702-06	H,Ms,R	WB,IP,FC
VEGF	ET1604-28	H,Ms,R	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
VEGF Receptor 1	ET1605-11	H,Ms,R	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
VEGF Receptor 2	ET1608-33	H,R	WB,IF-Cell,IHC-P,IP

*H: Human; Ms: Mouse; R: Rat; Z: Zebrafish; MK: Monkey; D: Dog; C:Chicken;
Cmk:Cynomolgus monkey; P:Pig

