




Inflammasome Signaling Pathway



扫码了解更多
“Inflammasome信号通路图”
相关研究产品信息

一抗产品				
Product name	Catalog#	Reactivity	Application	Antibody type
AIM2	HA500270	H	WB, IF-Cell, FC	Rabbit pAb
ATG4B	HA721829	H, Ms, R, Mk	WB	Rabbit mAb
ATG5	ET1611-38	H, Ms, R, Mk	WB, IF-Cell, IF-Tissue, IHC-P, IP, FC	Rabbit mAb
Bcl-2	HA721235	H	WB, IHC-P, IF-Cell, FC	Rabbit mAb
BRCC36	HA722226	H, Ms, R	WB, IF-Cell, IHC-P, FC	Rabbit mAb
Caspase-1	ET1608-69	H, Ms, R	WB, FC, IF-Cell	Rabbit mAb
Caspase-1 + p10 + p12	HA722222	H, Ms, R	WB, IF-Cell, IHC-P	Rabbit mAb
Caspase-12	HA500144	H, Ms, R	WB, IHC-P	Rabbit pAb
Caspase-12	HA722837	H, Ms, R	WB, IHC-P	Rabbit mAb
Caspase-5	ET1612-29	H	WB, IF-Cell, IF-Tissue, IHC-P, IP	Rabbit mAb
Caspase-8	ET1603-16	H	WB, IF-Tissue, IHC-P	Rabbit mAb
Caspase-8	ET1612-70	H	WB, IF-Tissue, IHC-P	Rabbit mAb
Caspase-8	HA722181	H	WB, IF-Cell, IHC-P, FC, IP	Rabbit mAb
Caspase-8 (p18)	HA722482	Ms	WB, IF-Cell	Rabbit mAb
GSDMD	HA601046	H, Ms	WB	Mouse mAb
Gasdermin D (N terminal)	ER1901-37	H, Ms, R	WB, IHC-P, FC, IF-Cell	Rabbit pAb
Gasdermin D (N terminal)	HA721144	H, Ms, R	WB, IHC-P, IF-Cell	Rabbit mAb
DFNA5 / GSDME	EM1901-48	H	WB, IF-Cell, IHC-P, FC	Mouse mAb
DFNA5 / GSDME	HA723113	H, Ms	WB, IP	Rabbit mAb
Human IL-1 beta	HA721712	H	WB	Rabbit mAb
IL-1 beta	HA601002	H, Ms, R	WB, ELISA, IF-Cell	Mouse mAb
IL-1 beta	HA601036	H, Ms, R	WB, ELISA, IHC-P	Mouse mAb
IL-1 beta	ET1701-39	H	WB, IF-Cell, IF-Tissue, IHC-P, FC	Rabbit mAb
IL-6	R1412-2	H, Ms, R	WB, IF-Cell, IHC-P, FC	Rabbit pAb
IL-6	EM1701-45	Ms, R	WB, IF-Cell	Mouse mAb
IL-18	HA721536	H	WB, IHC-P	Rabbit mAb
Human IL-18	HA721767	H	WB, IHC-P	Rabbit mAb
Human TNF alpha	IRS040	H	mIHC	Rabbit mAb
IKB alpha	ET1603-6	H, Ms, R	WB, IF-Cell, IF-Tissue, IHC-P, IP, FC	Rabbit mAb
Phospho-IKB alpha (S36)	HA721802	H, Ms, R	WB	Rabbit mAb
Phospho-IKB alpha (S32)	ET1609-78	H, Ms	WB, IF-Cell, IP	Rabbit mAb
Pyrin	HA722509	H	WB, FC	Rabbit mAb
RIP	ET1701-79	H	WB, FC, IP	Rabbit mAb
RIP3	ER1901-27	H, Ms, R	WB, IF-Cell, IHC-P	Rabbit pAb
Phospho-RIP3 (S232)	HA721428	Ms, R	WB, IF-Cell	Rabbit mAb
RIP3	HA722182	H	WB, IHC-P	Rabbit mAb
IKB beta	ET1611-1	H	WB, IF-Cell, IHC-P, IP	Rabbit mAb
IKK alpha	ET1611-15	H, Ms, R	WB, IHC-P, IP, FC	Rabbit mAb
IKK alpha + IKK beta	ET1611-23	H, Ms, R	WB, IF-Cell, IP	Rabbit mAb
IKK beta	HA721715	H	WB, FC	Rabbit mAb
IKK gamma	ET1704-40	H, Ms, R	WB, IHC-P, IF-Cell, IF-Tissue, IP, FC	Rabbit mAb
MyD88	ET1610-81	H	WB, IF-Cell, IF-Tissue, IHC-P, FC	Rabbit mAb
NF-kB p100 / NFKB2	ET1705-45	H, Ms, R	WB, IF-Cell, IHC-P, IP	Rabbit mAb
NF-kB p105 / p50	ET1603-18	H, Ms, R	WB, IHC-P, ICC, FC, IF-Tissue	Rabbit mAb
NF-kB p65	ET1603-12	H, Ms, R, Z	WB, IF-Cell, IF-Tissue, IHC-P, IP, FC	Rabbit mAb
Phospho-NF-kB p65 (S529)	ET1604-27	H	WB, IHC-P, FC, IP	Rabbit mAb
NOX2 / gp91phox	ET1611-44	H, Ms, R	WB	Rabbit mAb
NOX2 / gp91phox	HA723224	H, Ms, R	WB, IF-Cell, IHC-P, IHC-Fr, IP	Rabbit mAb
TNF alpha	HA722022	H, Ms	WB	Rabbit mAb
Human TNF alpha	IRS040	H	mIHC	Rabbit mAb
TNF alpha	R1203-1	Ms	WB, ELISA	Rabbit pAb
TNFAIP3	ET1611-40	H	WB, IF-Cell, IF-Tissue, IHC-P, FC	Rabbit mAb



一抗产品				
Product name	Catalog#	Reactivity	Application	Antibody type
TXNIP	ET1705-72	H, Ms, R	WB, IF-Cell, IF-Tissue, IHC-P, FC	Rabbit mAb
XIAP	HA722113	H, Ms, R, Mk	WB, IHC-P	Rabbit mAb
配对抗体				
Product name	Catalog#	Reactivity	Application	Antibody type
Human IL-1 beta	HA721278	H	ELISA (Cap)	Rabbit mAb
Human IL-1 beta	HA721279	H	ELISA (Det)	Rabbit mAb
Human IL-8	HA722098	H	ELISA (Cap)	Rabbit mAb
Human IL-8	HA722099	H	ELISA (Det)	Rabbit mAb
Human IL-18	HA721722	H	ELISA (Cap)	Rabbit mAb
Human IL-18	HA721723	H	ELISA (Det)	Rabbit mAb
Human NLRP3	HA723275	H	ELISA (Cap)	Rabbit mAb
Human NLRP3	HA723276	H	ELISA (Det)	Rabbit mAb
Human NLRP3	HA723278	H	ELISA (Det)	Rabbit mAb
Human TNF alpha	HA721294	H	ELISA (Det)	Rabbit mAb
Human TNF alpha	HA721559	H	ELISA (Det)	Rabbit mAb
Mouse IL-1 beta	HA722166	Ms	ELISA (Cap)	Rabbit mAb
Mouse IL-1 beta	HA722167	Ms	ELISA (Det)	Rabbit mAb
Mouse TNF alpha	HA722040	Ms	ELISA (Cap)	Rabbit mAb
Mouse TNF alpha	HA722041	Ms	ELISA (Det)	Rabbit mAb
Rat IL-1 beta	HA721836	R	ELISA (Cap)	Rabbit mAb
Rat IL-1 beta	HA721837	R	ELISA (Det)	Rabbit mAb
Rat TNF alpha	HA722170	R	ELISA (Cap)	Rabbit mAb
Rat TNF alpha	HA722171	R	ELISA (Det)	Rabbit mAb
ELISA Kits				
Product name	Catalog#	Sensitivity	Type	
Human CCL2 ELISA kit	EH0010	0.04 ng/ml (50 µl); 0.03 ng/ml (10 µl)	完全即用型	
Mouse CCL2/MCP-1 ELISA kit	EM0025	1.78 pg/ml (10 µl)	完全即用型	
Mouse CCL5 / RANTES ELISA kit	EM0032	1.86 pg/ml (10 µl)	完全即用型	
Human CCL5/RANTES ELISA kit	EH0027	2.78 pg/ml (50 µl); 4.12 pg/ml (10 µl)	完全即用型	
Mouse CXCL1/KC ELISA kit	EM0017	2.11 pg/ml (10 µl)	完全即用型	
Human CXCL1/GRO alpha ELISA kit	EH0054	14.89 pg/ml (50 µl); 32.67 pg/ml (10 µl)	完全即用型	
Mouse CXCL2/MIP-2 ELISA kit	EM0018	10.2 pg/ml (10 µl)	完全即用型	
Human IFN-β ELISA kit	EH0017	0.35 pg/ml (50 µl)	完全即用型	
Human IFN-α ELISA kit	EH0085	0.03 pg/ml (50 µl); 0.24 pg/ml (10 µl)	即用型	
Mouse IFN-β ELISA kit	EM0004	7.69 pg/ml (10 µl)	完全即用型	
Human IL-18 ELISA kit	EH0053	0.46 pg/ml (50 µl); 2.13 pg/ml (10 µl)	完全即用型	
Mouse IL-18 ELISA kit	EM0034	1.46 pg/ml (10 µl)	完全即用型	
Human IL-1β ELISA kit	EH0006	0.77 pg/ml (10 µl)	完全即用型	
Rat IL-1β ELISA kit	ER0001	17.11 pg/ml (50 µl)	完全即用型	
Mouse IL-1β ELISA kit	EM0029	0.91 pg/ml (10 µl)	完全即用型	
Human IL-33 ELISA kit	EH0047	0.76 pg/ml (50 µl); 1.25 pg/ml (10 µl)	完全即用型	
Human IL-6 ELISA kit	EH0001	0.26 pg/ml (50 µl); 0.43 pg/ml (10 µl)	完全即用型	
Rat IL-6 ELISA kit	EH0002	0.21 pg/ml (50 µl); 0.26 pg/ml (10 µl)	完全即用型	
Mouse IL-6 ELISA kit	EM0004	0.29 pg/ml (10 µl)	完全即用型	
Human IL-8/CXCL8 ELISA kit	EH0014	1.45 pg/ml (50 µl); 3.63 pg/ml (10 µl)	完全即用型	
Human TNF-α ELISA kit	EH0002	0.21 pg/ml (50 µl); 0.26 pg/ml (10 µl)	完全即用型	
Mouse TNF-α ELISA kit	EM0010	0.097 pg/ml (10 µl)	完全即用型	
Rat TNF-α ELISA kit	ER0007	0.16 pg/ml (50 µl); 0.08 pg/ml (10 µl)	完全即用型	

*H: Human; Ms: Mouse; R: Rat; Z: Zebrafish; MK: Monkey

