








扫码了解更多
“p53信号通路”
相关研究产品信息

Product name	Catalog#	Reactivity	Application
14-3-3 sigma	ET1612-41	H,M,S,R	WB,IHC-P,IF-Tissue
Active+Pro Caspase-9	ET1610-95	H,M,S	WB,IP
Caspase-9	ET1603-27	H,M,S	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Cleaved-Caspase-9 (D353)	ER60008	Ms,R,H	IF-Cell,WB,IHC-P
Apaf-1	ET1607-12	H,M,S	WB,IF-Cell,IF-Tissue,IHC-P
ATM	ET1606-20	H,R	WB,IF-Cell,IF-Tissue,IHC-P
Phospho-ATM (S1981)	ET1705-50	H,M,S	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
ATR	HA500176	H,R	WB,IHC-P
Phospho-ATR (T1989)	HA500405	H	WB
Phospho-ATR (T1989)	HA721190	H	WB,IF-Cell
Bax	ET1603-34	H,M,S,R,G	WB,IF-Cell,FC
Bcl-2	ET1603-11	H	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Bcl-XL	ET1603-28	H,M,S,R	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
Caspase-8	ET1603-16	H	WB,IF-Tissue,IHC-P
CD82	HA722533	H	WB
CD95	ER1902-56	H	WB,FC
CDK1	ET1605-54	H,M,S,R	WB,IF-Cell,IHC-P,FC
CDK1	ET1607-51	H,M,S,R	WB,IF-Cell,IF-Tissue,IHC-P,IP
Phospho-CDK1 (T161)	HA721987	H,M,S,R	WB,IF-Cell,IHC-P,FC
Cdk2	ET1602-6	H,M,S,R	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Phospho-Cdk2 (Y15)	ET1611-52	H,M,S,R	WB,IHC-P,IP,IF-Tissue
Cdk4	ET1612-23	H,M,S,R	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Cdk6	ET1612-3	H	IF-Cell,IF-Tissue,IHC-P,FC,WB
CDKN1A/p21	HA500156	H,M,S	WB,IHC-P
CDKN2A/p16INK4a	ER1803-53	H	WB,IF-Cell,IHC-P,FC
Chk1	HA722855	H,M,S,R	WB,IF-Cell,IP
Phospho-Chk1 (S280)	HA721471	H,M,S,R	WB
Phospho-Chk1 (S296)	ET1611-76	H,M,S,R	WB,IHC-P,IF-Cell
Phospho-Chk1 (S317)	HA722050	H,M,S,R	WB,IHC-P
Phospho-Chk1 (S345)	HA721189	H,M,S	WB,IF-Cell
Chk2	ET1610-52	H	WB,IF-Cell,IF-Tissue,IHC-P,IP
Phospho-Chk2 (T68)	HA721633	H,M,S,R	WB
Cleaved Caspase-3	ET1602-47	H,P	WB,IF-Cell,IF-Tissue,IHC-P
Cleaved+pro Caspase-3	ET1608-64	H,M,S,R	WB,IHC-P,IP
Pro Caspase-3	ET1602-39	H	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Cyclin B1	ET1608-27	H	WB,IF-Cell,IF-Tissue,IHC-P,IP
Cyclin B2	ET1612-21	H,R	WB,IF-Cell,IP
Cyclin D1	ET1601-31	H,M,S,R	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Cyclin D3	ET1612-4	H,M,S,R	WB
Cyclin E1	ET1612-16	H,M,S	WB,IF-Cell,IF-Tissue,IHC-P
Phospho-Cyclin E1 (T77)	ET1612-31	H,M,S,R	WB,IF-Cell,IF-Tissue,IHC-P
Cyclin E2	M0407-15	H,M,S,R	WB,IF-Cell,IHC-P
Cyclin E2	HA601146	H	WB,IHC-P,IF-Cell
Cytochrome C	ET1610-60	H,M,S,R	WB,IF-Tissue,IHC-P,IP,IF-Cell



Product name	Catalog#	Reactivity	Application
DDB2	ET7109-29	H	WB,IHC-P,FC
DR5	ET1704-55	H	WB,IP,IHC-P,FC
Fas(CD95)	ET1701-92	H	WB,IF-Tissue,IHC-P
GADD45B	HA723114	H,M,S,R	WB
GTSE-1	ER1910-29	H,M,S	WB,IHC-P,IF-Tissue
IGF1	HA722527	H	WB
IGFBP3	ER1911-12	H,R	IHC-P,IF-Tissue
Maspin	ER1912-78	H,M,S,R	WB,IHC-P,IF-Tissue
MDM2	RT1382	H,M,S,R	WB,IP,IF,IHC-P
MDM2	HA601310	H,Mk	WB,IF-Cell,IHC-P
Phospho-MDM2 (S166)	HA721791	H,M,S,R	WB
Noxa	ET1704-35	H,M,S,R	WB,IHC-P
p16	HA601131	H	IHC-P,WB,IF-Cell
p16INK4a	ET1608-62	H	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
p21	HA722065	Ms	WB,IF-Cell,IHC-P,FC
p21	HA722685	H	WB,IF-Cell,IHC-P,IP
p53	ET1601-13	H	WB,IF-Cell,ChIP
p53	ET1605-16	H	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
p53	HA722074	Ms	WB,IF-Cell,IP
p53 (acetyl K370)	ET1602-38	H,M,S,R	WB,IF-Cell,IF-Tissue,IP,FC
p53 (acetyl K382)	ET1612-71	H,M,S	WB,IF-Cell,IHC-P
Phospho-p53 (S15)	HA721756	H	WB,IF-Cell
Phospho-p53 (S20)	HA722761	H	WB
Phospho-p53 (S33)	HA722703	H,M,S,R	WB,IF-Cell,IHC-P
Phospho-p53 (S376)	ET1609-14	H,M,S	WB
Phospho-p53 (S392)	ET1606-24	H,M,S,R	WB,IHC-P,IP
Phospho-p53 (T55)	ET1609-13	H	WB,IHC-P
p73	ET1609-80	H,M,S	WB,IF-Cell,IF-Tissue,IHC-P
PAI1	HA601152	H,M,S,R	WB,IHC-P,IF-Tissue,FC
Pirh2	ET7111-16	H,M,S	WB
PTEN	ET1606-43	H,M,S,R	WB,IF-Cell,IF-Tissue,IHC-P
Phospho-PTEN (S380)	ET1701-46	H,M,S,R	WB,IF-Cell
PUMA	ET1602-24	H,M,S,R	WB,IF-Cell,IF-Tissue,IHC-P,FC
RRM2	ET1705-62	H,M,S,R	WB,IF-Cell,IF-Tissue,IHC-P
SERPINE1	EM1709-36	H,M,S	WB,IHC-P,FC
SESN1	HA722307	H	WB,IF-Cell,IHC-P,FC
SESN2 / Sestrin-2	HA722489	H,M,S,R	WB
SIAH1	HA721363	H,M,S,R	WB,IHC-P,IF-Cell,FC
Thrombospondin 1	HA721916	H,M,S,R	WB,IHC-P,IF-Cell,IF-Tissue,FC
Tuberin	ET1610-10	H,M,S,R	WB,IHC-P,FC
Phospho-Tuberin (S939)	HA721762	H,M,S	WB
Phospho-Tuberin (T1462)	HA721792	H,M,S,R	WB

*H: Human; Ms: Mouse; R: Rat; MK: Monkey; D: Dog; G: Goat; P: Pig

