






PI3K-AKT Signaling Pathway



扫码了解更多
“PI3K-AKT信号通路图”
相关研究产品信息

Product name	Catalog#	Reactivity	Application
Caspase-9	ET1603-27	H,Ms	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Active+Pro Caspase-9	R1308-12	H	WB,IHC-P
Active+Pro Caspase-9	ET1701-22	H	WB,IHC-P,FC
AKT1	ET1609-47	H,Ms,R	WB,IHC-P,IP,FC,IHC-Fr,IF-Cell
Phospho-AKT1 (S124)	ET1701-36	H,Ms,R	WB,IF-Cell,IF-Tissue,IHC-P,IP
Phospho-AKT1 (S129)	HA721677	H,Ms,R	WB
Phospho-AKT1 (T450)	ET1612-73	H,Ms,R	WB,IHC-P,IP,IF-Cell,FC
AKT1/2/3	ET1609-51	H,Ms,R,Mk	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
AKT1/2/3	HA721870	H,Ms,R,Mk	WB,IF-Cell,FC,IP
Phospho-AKT (S473)	ET1607-73	H,Ms,R,D	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
Phospho-AKT (T308)	HA722951	H	WB,IF-Cell,FC
AKT2	HA500091	H,Ms,R	WB,IHC-P
Phospho-AKT2 (S474)	HA500116	H,Ms	WB,IHC-P
Phospho-AKT2 (S474)	HA723108	H,Ms,R	WB,FC
Bad	ET1601-7	H,Ms,R	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
Cdk2	ET1602-6	H,Ms,R	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Cdk4	ET1612-1	H	WB,IF-Tissue,IHC-P
Cdk4	ET1612-23	H,Ms,R	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Phospho-CDK1 (T161) + CDK2/CDK3 (T160)	HA722582	H	WB, Dot Blot
Phospho-Cdk1/2/3 (T14)	HA722272	H,Ms,R	WB,IHC-P
Phospho-Cdk2 (Y15)	ET1611-52	H,Ms,R	WB,IHC-P,IP,IF-Tissue
Cdk6	ET1612-3	H	IF-Cell,IF-Tissue,IHC-P,FC,WB
Cyclin D1	ET1601-31	H,Ms,R	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Cyclin D2	HA500325	H,Ms,R	WB
Cyclin D3	ET1612-4	H,Ms,R	WB
Cyclin E1	ET1612-16	H,Ms	WB,IF-Cell,IF-Tissue,IHC-P
Phospho-Cyclin E1 (T77)	ET1612-31	H,Ms,R	WB,IF-Cell,IF-Tissue,IHC-P
Cyclin E2	ET1612-17	H,Ms	IF-Cell,IF-Tissue,IHC-P,FC
Cyclin E2	M0407-15	H,Ms,R	WB,IF-Cell,IHC-P
EGFR	ET1604-44	H,Ms,R,Mk	WB,IF-Cell,IP,FC
EGFR	IRS046	H	mIHC
eIF4EBP1	ET1701-83	H,Ms,R	WB,IF-Cell,IF-Tissue,IHC-P,IP
Phospho-eIF4EBP1 (S65)	HA723106	H	WB,IF-Cell,IP
FOXO3A	ER1706-79	H,Ms	WB,IF-Cell,IHC-P,FC
FOXO3A	ET1604-11	H,Ms,R	WB,IF-Cell,IF-Tissue
Phospho-FOXO1A (T24)+ FOXO3A (T32)+ FOXO4 (T28)+ FOXO6 (T26)	HA722499	H,Ms	WB,IF-Cell,FC
Phospho-FOXO3a (S253)	ET1609-49	H,Ms,R	WB,IHC-P,IF-Cell,FC
GDNF	ET1704-46	H,Ms,R	WB,IF-Cell
Growth Hormone	ET1703-47	H	WB,IHC-P,IP
Growth hormone receptor	ET1706-49	H,Ms,R	WB,IF-Cell,IHC-P,FC,IF-Tissue
GSK3 beta	ET1607-71	H,Ms,R	WB,IF-Cell,IF-Tissue,IHC-P,FC,IHC-Fr
Phospho-GSK3 (alpha + beta) (Y216 + Y279)	ET1607-54	H,Ms,R	WB,IF-Cell,IHC-P,IP,IF-Tissue,FC
Phospho-GSK3 beta (S9)	ET1607-60	H,Ms,R	WB,IF-Cell,IF-Tissue,IHC-P,FC
IGF1 (E peptide)	HA500107	H	WB,IHC-P,FC
IGF1R	HA600033	H	WB,FC
IGF2	ET1701-88	H	WB,IF-Cell,IF-Tissue,IHC-P



Product name	Catalog#	Reactivity	Application
IKK alpha	ET1611-15	H,Ms,R	WB,IHC-P,IP,FC
IKK alpha + IKK beta	ET1611-23	H,Ms,R	WB,IF-Cell,IP
IKK beta	HA721715	H	WB,FC
IKK gamma	ET1704-40	H,Ms,R	WB,IHC-P,IF-Cell,IF-Tissue,IP,FC
MDM2	HA601310	H,Mk	WB,IF-Cell,IHC-P
Phospho-MDM2 (S166)	HA721791	H,Ms,R	WB
mTOR	ET1608-5	H,Ms,R	WB,IHC-P,IP,IF-Cell
Phospho-mTOR (S2448)	HA600094	H,Ms,R	WB,IHC-P
Phospho-mTOR (S2481)	HA721632	H	WB
NF-kB p105 / p50	ET1603-18	H,Ms,R	WB,IHC-P,ICC,FC,IF-Tissue
NF-kB p65	ET1603-12	H,Ms,R,Z	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Phospho-NF-kB p65 (S529)	ET1604-27	H	WB,IHC-P,FC,IP
NGF	ET1606-29	H,Ms,R,Z	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
p21	HA722685	H	WB,IF-Cell,IHC-P,IP
p21	HA722065	Ms	WB,IF-Cell,IHC-P,FC
p27 KIP 1	ET1608-61	H,Ms,R	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
p53	EM20603	H	WB,IF-Cell,IHC-P,FC
p53	ET1601-13	H	WB,IF-Cell,ChIP
p53	ET1605-16	H	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
p53 (acetyl K370)	ET1602-38	H,Ms,R	WB,IF-Cell,IF-Tissue,IP,FC
p53 (acetyl K382)	ET1612-71	H,Ms	WB,IF-Cell,IHC-P
Phospho-p53 (S15)	HA721756	H	WB,IF-Cell
Phospho-p53 (S20)	HA722761	H	WB
Phospho-p53 (S33)	HA722703	H,Ms,R	WB,IF-Cell,IHC-P
Phospho-p53 (S376)	ET1609-14	H,Ms	WB
Phospho-p53 (S392)	ET1606-24	H,Ms,R	WB,IHC-P,IP
Phospho-p53 (T55)	ET1609-13	H	WB,IHC-P
p70 S6 Kinase	HA722520	H,Ms,R	WB,IHC-P,IP
P70 S6 Kinase beta	HA721354	H,Ms	WB,IF-Cell,FC,IHC-P
Phospho-p70 S6 Kinase (T389)	HA721803	H,Ms,R	WB,IF-Cell,IHC-P
p75 NGF Receptor	ET1601-22	H,Ms,R	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
PDGF Receptor beta	ET1605-20	H,Ms,R	WB,IHC-P,IP,FC,IF-Tissue,IF-Cell,mIHC,IHC-Fr
PDGF Receptor beta	IRS058	Ms	mIHC
PDGF Receptor beta	IRS038	H	mIHC
PDGFR alpha	ET1702-49	H,Ms,R,Cmk,P	WB,IHC-P,IHC-Fr
Phospho-RPS6 (S235 + S236)	HA721275	H,Ms,R	WB,IHC-P,IF-Cell
Phospho-RPS6 (S240 + S244)	HA721589	H,Ms,R	WB,IHC-P,IF-Cell
PI 3 Kinase p85 alpha	EM1701-62	H,Ms,R	WB,IHC-P,IF-Cell,IF-Tissue
PI 3 Kinase p85 alpha	ET1608-70	H,Ms,R	WB,IF-Cell,IF-Tissue,IHC-P,FC
PI3 Kinase p110 beta	HA722474	H,Ms,R	WB,IF-Cell,IF-Tissue
PI3 Kinase p110α	ET1606-36	H,Ms	WB,IF-Cell,IF-Tissue,IP,IHC-P
PI3-kinase p110 subunit beta	ET1610-36	H	WB,IP,IHC-P,FC
Phospho-PI3K p85 (Y467) + PI3K p55 (Y199)	HA721672	H,Ms	WB,IF-Cell
PTEN	ET1606-43	H,Ms,R	WB,IF-Cell,IF-Tissue,IHC-P
Phospho-PTEN (S380)	ET1701-46	H,Ms,R	WB,IF-Cell
VEGF	ET1604-28	H,Ms,R	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
VEGF Receptor 1	ET1605-11	H,Ms,R	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
VEGF Receptor 2	ET1608-33	H,R	WB,IF-Cell,IHC-P,IP

*H: Human; Ms: Mouse; R: Rat; Z: Zebrafish; Mk: Monkey; D: Dog; Cmk: Cynomolgus monkey; P: Pig

