



Hangzhou HuaAn Biotechnology Co.,Ltd.

ELISA 试剂盒 配对抗体

- Quality ELISA 试剂盒
- DTSet ELISA 试剂盒
- 配对抗体

让研究变得简单

杭州华安生物技术有限公司

☎ 0571-88062880

🖨 0571-88060400

✉ sales@huabio.cn

📍 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园区2号楼16层

 华安生物
H U A B I O

HUABIO

完全即用型： Quality ELISA 试剂盒

我们基于双抗夹心法开发了一系列性能稳定，品质卓越的完全即用型Quality ELISA 试剂盒产品。为满足客户对产品质量的深层次需求，在生产、研发和质控环节，我们严格筛选兔单克隆抗体对，并经过多种生物样本验证，确保Quality ELISA试剂盒的最优性能。

目录

CONTENTS

1. Quality ELISA 试剂盒	01-08
2. DTSet ELISA 试剂盒	09-11
3. 配对抗体	12-28

01



产品优势

开盒即用

实验流程简化高效，一抗孵育仅需90 min，一天可以做多轮实验

高灵敏度

准确检测pg级别蛋白，适用于低表达蛋白检测

高精密度

CV值低于8 %，保证检测结果的一致性

高稳定性

2-8℃保存，有效期高达18个月

CE认证

进口品质，畅销国外，深受国内、外客户认可

Quality ELISA 试剂盒

应用领域广，适用于热门领域相关靶点的定性、定量检测

肿瘤/免疫治疗-热销靶点

BCMA	Chitinase 3-like 1	GDF-15	IFN-β	IL-8	PD-1	TNF-α
CCL3	Clusterin	GM-CSF	IL-18	Lipocalin-2	PDGF-AB	TNF-β
CD147	CXCL1	HGF	IL-2	MCP-1	PDGF-BB	VCAM-1
CD25	EGF	HGF	IL-22	NSE	PD-L1	VEGF
CD28	G-CSF	IDO	IL-6	Osteoprotegerin(OPG)	TGF-β	VEGFR2
CD47						

神经科研-热销靶点

BDNF	C5a	Enolase 2/NSE	Lipocalin-2	Neuropilin-1	PDGF-BB	VCAM-1
b-NGF	CRP					

心血管研究-热销靶点

ACE2	C-peptide	GH	IL-1β	Insulin	PCSK9	TNF-α
Angiogenin(ANG)	CRP	HGF	IL-2	Leptin	PDGF-BB	tPA
ANGPTL3	FGF-21	IFN-γ	IL-6	MCP-1	ST2	VEGF
Clusterin	GDF-15	IL-18	IL-8	MMP-9	TGF-β	

免疫/炎症/自身免疫性疾病研究-热销靶点

C5a	GM-CSF	IL-10	IL-18	IL-22	IL-5	MMP-9
CRP	HIV-1 Gag p24	IL-12p70	IL-1β	IL-23	IL-6	TNF-α
FGF-19	IFN-β	IL-13	IL-2	IL-4	MCP-1	TNF-β
FGF-21	IFN-γ	IL-17	IL-21			

血液系统疾病研究-热销靶点

BCMA	Granzyme B	IDO	IL-15	IL-21	ST2	TPO
Fibronectin	HGF	IFN-γ	IL-2	IL-7	TGF-β	

肥胖/糖尿病/代谢性疾病研究-热销靶点

ACE2	CRP	GDF-15	HGF	IL-1β	Leptin	TNF-α
ANGPTL3	FGF-19	GH	IGF-1	IL-6	Lipocalin-2	VEGF
C-peptide	FGF-21	Glucagon	IL-10	Insulin	MCP-1	

高灵敏度，适用于pg/mL级别蛋白的检测

靶点名称	华安生物 & 攸克生物		竞品 S
	Cat Number	Sensitivity	Sensitivity
Human IL-6	EH0001	0.26 pg/ml	0.7 pg/ml
Human TNF-alpha	EH0002	0.21 pg/ml	6.23 pg/ml
Human IFN-gamma	EH0003	0.36 pg/ml	5.69 pg/ml
Human IL-17A	EH0004	0.26 pg/ml	15 pg/ml
Human IL-4	EH0007	0.03 pg/ml	10 pg/ml

在相同实验操作条件下，对人源IL-6，TNF-alpha等ELISA 试剂盒和竞品S相应靶标进行最低检测限评价实验。结果表明，我们的ELISA 试剂盒灵敏度显著优于竞品，高达竞品的300倍。

精密度优秀，CV值低于8 %，保证检测值的一致性

Human IL-6	华安生物 & 攸克生物			竞品 S		
Sample Number	S1	S2	S3	S1	S2	S3
	22	22	22	20	20	20
Average (pg/ml)	24	115.4	263.2	15.8	95.6	179
Standard Deviation	1.1	4.5	9.3	0.7	3	3.1
Coefficient of variation(%)	5.6	3.9	3.6	4.4	3.1	1.7

*精密度评价标准：及格，CV≤15%；良好，CV≤10%；优秀，CV≤8%；卓越，CV≤5%

在相同实验操作条件下，使用人源IL-6 ELISA 试剂盒和竞品S检测3个批次样本中 IL-6 的含量，并进行精密度统计分析。结果表明，我们的ELISA 试剂盒CV值低于8%，严格控制批内/批间差，与竞品S精密度相当。


准确度优秀，线性范围广，确保实验结果的准确可靠

靶点名称	华安生物 & 攸克生物		竞品 S
	Cat Number	Recovery (%)	Recovery (%)
Human IL-6	EH0001	84-115	86-99
Human TNF-alpha	EH0002	95-105	95-108
Human IFN-gamma	EH0003	89-113	90-109
Human IL-17A	EH0004	80-112	90-105
Human IL-4	EH0007	88-105	82-111

*准确度评价标准：及格，回收率在50%~150%；优秀，回收率在80%~120%

在相同实验操作条件下，对人源IL-6，TNF-alpha等ELISA 试剂盒和竞品S相应靶标进行准确度评价分析。结果表明，我们的ELISA 试剂盒回收率严格控制在80%-120%之间，较为优秀，与竞品准确度相当。

客户引用案例精选




使用人源IFN-β（货号EH0017），CXCL10（货号EH0026），IL-6（货号EH0001）及TNF-α（货号EH0002）ELISA 试剂盒检测转染TFEB siRNAs和TFE3 siRNAs的THP-1细胞在S. Typhimurium感染16h后细胞上清中相应蛋白的含量。

数据来源：Immunity. 2025 Feb 11;58(2):309-325.e6

IF: 25.5

科研单位：湖南第一师范学院



Quality ELISA 试剂盒的实验流程图

热销Quality ELISA 试剂盒清单精选

白介素类			
Product name	Catalog#	Size	Sensitivity
Human IL-6 ELISA Kit	EH0001	96T	0.26 pg/ml (50 µl) ; 0.43 pg/ml (10 µl)
Human IL-10 ELISA Kit	EH0009	96T	0.89 pg/ml (10 µl)
Human IL-12/IL-23 p40 ELISA Kit	EH0072	96T	13.59 pg/ml (50 µl); 37.7 pg/ml (10 µl)
Human IL-12B p40 ELISA Kit	EH0098	96T	0.29 pg/ml (50 µl); 1.84 pg/ml (10 µl)
Human IL-12p70 ELISA Kit	EH0021	96T	3.42 pg/ml (50 µl) ; 11.63 pg/ml (10 µl)
Human IL-13 ELISA Kit	EH0049	96T	101.81 pg/ml (50 µl); 223.66 pg/ml (10 µl)
Human IL-15 ELISA Kit	EH0031	96T	0.14 pg/ml (50 µl) ; 0.79 pg/ml (10 µl)
Human IL-17A ELISA Kit	EH0004	96T	0.26 pg/ml (50 µl) ; 0.65 pg/ml (10 µl)
Human IL-17E ELISA Kit	EH0066	96T	0.69 pg/ml (50 µl); 1.45 pg/ml (10 µl)
Human IL-18 ELISA Kit	EH0053	96T	0.46pg/ml (50 µl); 2.13 pg/ml (10 µl)
Human IL-1β ELISA Kit	EH0006	96T	0.77 pg/ml (10 µl)
Human IL-2 ELISA Kit	EH0013	96T	1.67 pg/ml (50 µl) ; 3.08 pg/ml (10 µl)
Human IL-21 ELISA Kit	EH0018	96T	21.7 pg/ml (50 µl) ; 67.9 pg/ml (10 µl)
Human IL-22 ELISA Kit	EH0086	96T	0.03 pg/ml (50 µl); 0.21 pg/ml (10 µl)
Human IL-23 p19 ELISA Kit	EH0109	96T	0.32 pg/ml (50 µl); 1.3 pg/ml (10 µl)
Human IL-27 ELISA Kit	EH0094	96T	5.51 pg/ml (50 µl); 17.28 pg/ml (10 µl)
Human IL-3 ELISA Kit	EH0050	96T	3.74 pg/ml (50 µl); 20.64 pg/ml (10 µl)
Human IL-31 ELISA Kit	EH0097	96T	0.4 pg/ml (50 µl); 2.68 pg/ml (10 µl)
Human IL-33 ELISA Kit	EH0047	96T	0.76 pg/ml (50 µl); 1.25 pg/ml (10 µl)
Human IL-4 ELISA Kit	EH0007	96T	0.03 pg/ml (50 µl) ; 0.28 pg/ml (10 µl)
Human IL-5 ELISA Kit	EH0008	96T	0.15 pg/ml (50 µl) ; 0.73 pg/ml (10 µl)
Human IL-7 ELISA Kit	EH0045	96T	0.62 pg/ml (50 µl)
Human IL-8/CXCL8 ELISA Kit	EH0014	96T	1.45 pg/ml (50 µl) ; 3.63 pg/ml (10 µl)
Mouse IL-10 ELISA Kit	EM0005	96T	5.75 pg/ml (10 µl)
Mouse IL-12/IL-23 p40 ELISA Kit	EM0012	96T	1.50 pg/ml (10 µl)
Mouse IL-12p70 ELISA Kit	EM0006	96T	5.50 pg/ml (10 µl)
Mouse IL-13 ELISA Kit	EM0007	96T	0.37 pg/ml (10 µl)
Mouse IL-17A ELISA Kit	EM0008	96T	0.25 pg/ml (10 µl)
Mouse IL-17E ELISA Kit	EM0028	96T	12.33 pg/ml (50 µl)
Mouse IL-18 ELISA Kit	EM0034	96T	1.46 pg/ml (10 µl)
Mouse IL-1β ELISA Kit	EM0029	96T	0.91 pg/ml (10 µl)

白介素类			
Product name	Catalog#	Size	Sensitivity
Mouse IL-2 ELISA Kit	EM0002	96T	0.16 pg/ml (10 µl)
Mouse IL-21 ELISA Kit	EM0020	96T	10.79 pg/ml (10 µl)
Mouse IL-22 ELISA Kit	EM0021	96T	11.2 pg/ml (10 µl)
Mouse IL-33 ELISA Kit	EM0033	96T	0.33 pg/ml (10 µl)
Mouse IL-4 ELISA Kit	EM0003	96T	0.35 pg/ml (10 µl)
Mouse IL-5 ELISA Kit	EM0023	96T	0.24 pg/ml (10 µl)
Mouse IL-6 ELISA Kit	EM0004	96T	0.29 pg/ml (10 µl)
Mouse IL-9 ELISA Kit	EM0022	96T	0.10 pg/ml (10 µl)
Rat IL-10 ELISA Kit	ER0005	96T	5.76 pg/ml (50 µl)
Rat IL-1β ELISA Kit	ER0001	96T	17.11 pg/ml (50 µl)
Rat IL-2 ELISA Kit	ER0002	96T	8.13 pg/ml (50 µl)
Rat IL-4 ELISA Kit	ER0003	96T	17.29pg/ml (50 µl)
Rat IL-6 ELISA Kit	ER0004	96T	7.05pg/ml (50 µl)
集落刺激因子			
Product name	Catalog#	Size	Sensitivity
Mouse GM-CSF/CSF2 ELISA Kit	EM0011	96T	0.19 pg/ml (10 µl)
Human GM-CSF ELISA Kit	EH0051	96T	3.57 pg/ml (50 µl); 5.45 pg/ml (10 µl)
肿瘤坏死因子			
Product name	Catalog#	Size	Sensitivity
Human HVEM/TNFRSF14 ELISA Kit	EH0107	96T	1.68 pg/ml (50 µl); 1.86 pg/ml (10 µl)
Human OX40/TNFRSF4 ELISA Kit	EH0120	96T	1.15 pg/ml (50 µl); 8.63 pg/ml (10 µl)
Human TNF-α ELISA Kit	EH0002	96T	0.21 pg/ml (50 µl) ; 0.26 pg/ml (10 µl)
Mouse TNF-α ELISA Kit	EM0010	96T	0.097 pg/ml (10 µl)
Rat TNF-α ELISA Kit	ER0007	96T	0.16 pg/ml (50 µl); 0.08pg/ml (10 µl)
趋化因子及其受体			
Product name	Catalog#	Size	Sensitivity
Human CCL2/MCP-1 ELISA Kit	EH0010	96T	0.04 ng/ml (50 µl); 0.03 ng/ml (10 µl)
Human CCL5/RANTES ELISA Kit	EH0027	96T	2.78 pg/ml (50 µl) ; 4.12 pg/ml (10 µl)
Human CCL20/ MIP3 ELISA Kit	EH0028	96T	8.58 pg/ml (50 µl) ; 12.19 pg/ml (10 µl)
Human CCL25 ELISA Kit	EH0117	96T	0.06 pg/ml (50 µl); 0.32 pg/ml (10 µl)
Human CXCL1/GRO alpha ELISA Kit	EH0054	96T	14.89 pg/ml (50 µl); 32.67 pg/ml (10 µl)
Human CXCL10/IP-10 ELISA Kit	EH0026	96T	0.33 pg/ml (50 µl) ; 1.14 pg/ml (10 µl)

趋化因子及其受体			
Product name	Catalog#	Size	Sensitivity
Human CXCL9/MIG ELISA Kit	EH0015	96T	2.24 pg/ml (50 µl) ; 5.32 pg/ml (10 µl)
Mouse CCL2/MCP-1 ELISA Kit	EM0025	96T	1.78 pg/ml (10 µl)
Mouse CCL3/MIP-1 alpha ELISA Kit	EM0031	96T	0.51 pg/ml (10 µl)
Mouse CCL5/RANTES ELISA Kit	EM0032	96T	1.86 pg/ml (10 µl)
Mouse CXCL1/KC ELISA Kit	EM0017	96T	2.11 pg/ml (10 µl)
Mouse CXCL10/IP-10/CRG-2 ELISA Kit	EM0014	96T	213.35 pg/ml (10 µl)
Mouse CXCL2 ELISA Kit	EM0018	96T	10.2 pg/ml (10 µl)
Mouse CXCL9/MIG ELISA Kit	EM0015	96T	0.29 pg/ml (10 µl)
干扰素			
Product name	Catalog#	Size	Sensitivity
Human IFN-α ELISA Kit	EH0085	96T	0.03 pg/ml (50 µl) ; 0.24 pg/ml (10 µl)
Human IFN-β ELISA Kit	EH0017	96T	0.35 pg/ml (50 µl)
Human IFN-γ ELISA Kit	EH0039	96T	0.36 pg/ml (50 µl); 0.85 pg/ml (10 µl)
Human IFN-gamma R1 ELISA Kit	EH0104	96T	2.28 pg/ml (50 µl); 11.4 pg/ml (10 µl)
Mouse IFN-beta ELISA Kit	EM0040	96T	7.69 pg/ml (10 µl)
Mouse IFN-γ ELISA Kit	EM0009	96T	3.24 pg/ml (10 µl)
Rat IFN-γ ELISA Kit	ER0006	96T	7.39 pg/ml (50 µl)
生长因子 & 转化因子			
Product name	Catalog#	Size	Sensitivity
Human beta-NGF ELISA Kit	EH0065	96T	1.05 pg/ml (50 µl); 2.53 pg/ml (10 µl)
Human EGF ELISA Kit	EH0023	96T	3.27 pg/ml (50 µl); 1.91 pg/ml (10 µl)
Human FGF basic/FGF2/bFGF ELISA Kit	EH0024	96T	0.12 pg/ml (50 µl) ; 0.59 pg/ml (10 µl)
Human HGF ELISA Kit	EH0019	96T	1.51 pg/ml (50 µl) ; 9.93 pg/ml (10 µl)
Human Leptin ELISA Kit	EH0056	96T	4.07 pg/ml (50 µl) ; 6.01 pg/ml (10 µl)
Human Prolactin ELISA Kit	EH0089	96T	1.76 pg/ml (50 µl); 3.43 pg/ml (10 µl)
Human TGF alpha ELISA Kit	EH0059	96T	9.45 pg/ml (50 µl); 18.83 pg/ml (10 µl)
Human VEGF ELISA Kit	EH0075	96T	19.30 pg/ml (50 µl)
Human/Mouse BDNF ELISA Kit	EH0069	96T	3.47 pg/ml (50 µl); 11.20 pg/ml (10 µl)
Human/Mouse/Rat/Porcine/Canine TGF-β1 ELISA Kit	EH0012	96T	1.91 pg/ml (50 µl); 3.89 pg/ml (10 µl)
Mouse Leptin ELISA Kit	EM0027	96T	33.32 pg/ml (50 µl)
Mouse VEGF ELISA Kit	EM0013	96T	0.01 pg/ml (50 µl); 1.02 pg/ml (10 µl)
Mouse/Rat IGF1/IGF-I/IGF-1 ELISA Kit	EM0024	96T	0.57 pg/ml (10 µl)

生长因子 & 转化因子			
Product name	Catalog#	Size	Sensitivity
Rat VEGF ELISA Kit	ER0008	96T	2.60 pg/ml (50 µl); 3.77 pg/ml (10 µl)
免疫&炎症			
Product name	Catalog#	Size	Sensitivity
Human B7-H6/NCR3LG1 ELISA Kit	EH0116	96T	3.98 pg/ml (50 µl); 4.61 pg/ml (10 µl)
Human CD200R1 ELISA Kit	EH0087	96T	0.1 pg/ml (50 µl); 3.09 pg/ml (10 µl)
Human CD27/TNFRSF7 ELISA Kit	EH0114	96T	0.10 pg/ml (50 µl); 18.37 pg/ml (10 µl)
Human CD30 ELISA Kit	EH0105	96T	0.8 pg/ml (50 µl); 3.19 pg/ml (10 µl)
Human CD33 ELISA Kit	EH0106	96T	0.83 pg/ml (50 µl); 0.31 pg/ml (10 µl)
Human CD48 ELISA Kit	EH0119	96T	0.25 pg/ml (50 µl); 1.05 pg/ml (10 µl)
Human CD80 ELISA Kit	EH0121	96T	0.08 pg/ml (50 µl); 0.52 pg/ml (10 µl)
Human CD85d ELISA Kit	EH0046	96T	10.40 pg/ml (50 µl): 14.80 pg/ml (10 µl)
Human CDCP1 ELISA Kit	EH0083	96T	1.41pg/ml (50 µl)
Human CRP ELISA Kit	EH0076	96T	4.26 pg/ml (50 µl); 21.23 pg/ml (10 µl)
Human Granzyme B ELISA Kit	EH0052	96T	1.12 pg/ml (50 µl); 5.96 pg/ml (10 µl)
Human GITR/TNFRSF18 ELISA Kit	EH0095	96T	0.07 pg/ml (50 µl); 1.20 pg/ml (10 µl)
Human ICAM-1/CD54 ELISA Kit	EH0064	96T	0.49 pg/ml(50 µl); 6.36 pg/ml (10 µl)
Human IDO1 ELISA Kit	EH0118	96T	3.04pg/ml (50 µl); 3.61 pg/ml (10 µl)
Human LIF ELISA Kit	EH0096	96T	0.10 pg/ml (50 µl); 0.36 pg/ml (10 µl)
Human MBL2 ELISA Kit	EH0090	96T	1.92 pg/ml (50 µl)
Human MICA ELISA Kit	EH0044	96T	8.3 pg/ml (50 µl)
Human NCAM 1/CD56 ELISA Kit	EH0102	96T	3.69 pg/ml (50 µl); 37.50 pg/ml (10 µl)
Human PD-L1/B7-H1 ELISA Kit	EH0048	96T	8.50 pg/ml (50 µl); 26.02 pg/ml (10 µl)
Human TSLP ELISA Kit	EH0061	96T	1.56 pg/ml (50 µl); 4.58 pg/ml (10 µl)
Human VCAM-1/CD106 ELISA Kit	EH0067	96T	7.07 pg/ml (50 µl)
Mouse CRP/C-Reactive Protein ELISA Kit	EM0035	96T	0.01 ng/ml (50 µl)
Mouse Granzyme B ELISA Kit	EM0016	96T	9.91 pg/ml (10 µl)
Mouse LIF ELISA Kit	EM0038	96T	2.77 pg/ml (10 µl)
Mouse Myeloperoxidase ELISA Kit	EM0030	96T	0.10 ng/ml (10 µl)
Mouse PD-1 ELISA Kit	EM0037	96T	7.65 pg/ml (10 µl)
Mouse PD-L1 ELISA Kit	EM0036	96T	0.16 pg/ml (10 µl)
Mouse TREM2 ELISA Kit	EM0039	96T	0.60 pg/ml (10 µl)
Rat CRP/C-Reactive Protein ELISA Kit	ER0009	96T	27.24 pg/ml (50 µl); 63.91 pg/ml (10 µl)

心血管研究			
Product name	Catalog#	Size	Sensitivity
Human Adiponectin ELISA Kit	EH0099	96T	0.15 pg/ml (50 µl); 3.56 pg/ml (10 µl)
Human MMP-8 ELISA Kit	EH0103	96T	1.59 pg/ml (50 µl); 9.04 pg/ml (10 µl)
Human MMP-9 ELISA Kit	EH0020	96T	2.7 pg/ml (50 µl) ; 7.9 pg/ml (10 µl)
Human Osteopontin (OPN) ELISA Kit	EH0057	96T	0.04 ng/ml (50 µl); 0.12ng/ml (10 µl)
Human Osteoprotegerin/TNFRSF11B ELISA Kit	EH0070	96T	0.37 pg/ml (50 µl); 2.53 pg/ml (10 µl)
Human PDGF-AA ELISA Kit	EH0058	96T	0.10 ng/ml (50 µl); 0.17 ng/ml (10 µl)
Human PDGF-BB ELISA Kit	EH0068	96T	2.79 pg/ml (50 µl); 7.56 pg/ml (10 µl)
Human uPAR ELISA Kit	EH0077	96T	8.79 pg/ml (50 µl)
细胞生物学			
Product name	Catalog#	Size	Sensitivity
Human Asprosin ELISA Kit	EH0032	96T	0.78 ng/ml (50 µl)
Human E-Cadherin ELISA Kit	EH0115	96T	25.58 pg/ml (50 µl)
Human Total HSP60 ELISA Kit	EH0060	96T	3.22 pg/ml (50 µl); 4.45 pg/ml (10 µl)
Human VE Cadherin ELISA Kit	EH0091	96T	0.01 ng/ml (50 µl); 0.12 ng/ml (10 µl)
Human/Rat/Mouse Irisin ELISA Kit	EH0081	96T	2.60 pg/ml (50 µl); 9.39 pg/ml (10 µl)
微生物学			
Product name	Catalog#	Size	Sensitivity
HIV-1 p24 ELISA Kit	EH0078	96T	0.64 pg/ml (50 µl); 2.62 pg/ml (10 µl)
HIV-1 p24 ELISA Kit	EHY0002	96T	0.12 ng/ml (10 µl)



扫码了解更多
“Quality ELISA试剂盒”
产品信息



自主灵活： DTSet ELISA 试剂盒

我们提供的DTSet ELISA 试剂盒包含夹心法免疫分析的必需组分，具有显著的自主灵活性和高性价比的特征。同时，能够兼顾应用场景的通用性及实验结果的稳定可靠。

试剂盒组分		Volume
捕获抗体	CA (capture antibody 100×)	1 vial
标准品	S (Standard)	1 vial
检测抗体	DA (Detect Antibody 100×)	1 vial
链霉亲和素-HRP	SH (Streptavidin-HRP 100×)	1 vial



热销DTSet ELISA 试剂盒清单精选

Product name	Catalog#	Size	Catalog#	Size	Sensitivity
Human IL-5 DTSet ELISA Kit	DSEH000805	96T*5	DSEH000815	96T*15	1.69 pg/ml(50 µl)
Human IL-2 DTSet ELISA Kit	DSEH001305	96T*5	DSEH001315	96T*15	0.25 pg/ml(50 µl)
Human IL-8 DTSet ELISA Kit	DSEH001405	96T*5	DSEH001415	96T*15	0.11 pg/ml(50 µl)
Human IFN-beta DTSet ELISA Kit	DSEH001705	96T*5	DSEH001715	96T*15	6.26 pg/ml(50 µl)
Human HGF DTSet ELISA Kit	DSEH001905	96T*5	DSEH001915	96T*15	0.45 pg/ml(50 µl)
Human MMP-9 DTSet ELISA Kit	DSEH002005	96T*5	DSEH002015	96T*15	1.27 pg/ml(50 µl)
Human EGF DTSet ELISA Kit	DSEH002305	96T*5	DSEH002315	96T*15	0.12 pg/ml(50 µl)
Human FGF-2 DTSet ELISA Kit	DSEH002405	96T*5	DSEH002415	96T*15	1.66 pg/ml(50 µl)
Human CXCL10/IP-10 DTSet ELISA Kit	DSEH002605	96T*5	DSEH002615	96T*15	0.008 pg/ml(50 µl)
Human IL-15 DTSet ELISA Kit	DSEH003105	96T*5	DSEH003115	96T*15	0.06 pg/ml(50 µl)
Human IFN-γ DTSet ELISA Kit	DSEH003905	96T*5	DSEH003915	96T*15	0.14 pg/ml(50 µl)
Human IL-7 DTSet ELISA Kit	DSEH004505	96T*5	DSEH004515	96T*15	0.86 pg/ml(50 µl)
Human CD85d DTSet ELISA Kit	DSEH004605	96T*5	DSEH004615	96T*15	77.03 pg/ml(50 µl)
Human IL-33 DTSet ELISA Kit	DSEH004705	96T*5	DSEH004715	96T*15	2.21 pg/ml(50 ul)
Human IL-13 DTSet ELISA Kit	DSEH004905	96T*5	DSEH004915	96T*15	10.26pg/ml(50 ul)
Human IL-3 DTSet ELISA Kit	DSEH005005	96T*5	DSEH005015	96T*15	0.02 pg/ml(50 µl)
Human GM-CSF DTSet ELISA Kit	DSEH005105	96T*5	DSEH005115	96T*15	0.06 pg/ml(50 µl)
Human IL-18 DTSet ELISA Kit	DSEH005305	96T*5	DSEH005315	96T*15	1.26 pg/ml(50 µl)
Human CXCL1/GRO alpha DTSet ELISA Kit	DSEH005405	96T*5	DSEH005415	96T*15	16.04 pg/ml(50 µl)
Human ICAM-1/CD54 DTSet ELISA Kit	DSEH006405	96T*5	DSEH006415	96T*15	2.48pg/ml(50 µl)
Human PDGF-BB DTSet ELISA Kit	DSEH006805	96T*5	DSEH006815	96T*15	0.10 pg/ml(50 µl)
Human CRP DTSet ELISA Kit	DSEH007605	96T*5	DSEH007615	96T*15	8.06 pg/ml(50 µl)
Human IFN-α DTSet ELISA Kit	DSEH008505	96T*5	DSEH008515	96T*15	0.06 pg/ml(50 µl)
Human IL-22 DTSet ELISA Kit	DSEH008605	96T*5	DSEH008615	96T*15	0.41pg/ml(50 µl)
Human CD200R1 DTSet ELISA Kit	DSEH008705	96T*5	DSEH008715	96T*15	0.24 pg/ml(50 µl)
Human Prolactin DTSet ELISA Kit	DSEH008905	96T*5	DSEH008915	96T*15	3.25 pg/ml(50 µl)
Human MBL2 DTSet ELISA Kit	DSEH009005	96T*5	DSEH009015	96T*15	7.53 pg/ml(50 µl)
Human VE Cadherin DTSet ELISA Kit	DSEH009105	96T*5	DSEH009115	96T*15	48.52pg/ml(50 µl)
Human IL-27 DTSet ELISA Kit	DSEH009405	96T*5	DSEH009415	96T*15	0.40 pg/ml(50 µl)
Human GITR/TNFRSF18 DTSet ELISA Kit	DSEH009505	96T*5	DSEH009515	96T*15	1.49 pg/ml(50 µl)
Human LIF DTSet ELISA Kit	DSEH009605	96T*5	DSEH009615	96T*15	0.15 pg/ml(50 µl)

Product Name	Catalog#	Size	Catalog#	Size	Sensitivity
Human IL-31 DTSet ELISA Kit	DSEH009705	96T*5	DSEH009715	96T*15	0.07 pg/ml(50 µl)
Human IL-12B p40 DTSet ELISA Kit	DSEH009805	96T*5	DSEH009815	96T*15	0.03 pg/ml(50 µl)
Human Adiponectin DTSet ELISA Kit	DSEH009905	96T*5	DESH009915	96T*15	0.23pg/ml(50 µl)
Human CD276 DTSet ELISA Kit	DSEH010105	96T*5	DSEH010115	96T*15	0.12pg/ml(50 µl)
Human MMP-8 DTSet ELISA Kit	DSEH010305	96T*5	DSEH010315	96T*15	3.33 pg/ml(50 µl)
Human IFN-γ R1 DTSet ELISA Kit	DSEH010405	96T*5	DSEH010415	96T*15	0.12pg/ml(50 µl)
Human CD30 DTSet ELISA Kit	DSEH010505	96T*5	DSEH010515	96T*15	0.74 pg/ml(50 µl)
Human CD33 DTSet ELISA Kit	DSEH010605	96T*5	DSEH010615	96T*15	0.17 pg/ml(50 µl)
Human HVEM/TNFRSF14 DTSet ELISA Kit	DSEH010705	96T*5	DSEH010715	96T*15	2.08pg/ml(50 µl)
Human Complemenet C5 DTSet ELISA Kit	DSEH010805	96T*5	DSEH010815	96T*15	0.80pg/ml(50 µl)
Human IL-23 p19 DTSet ELISA Kit	DSEH010905	96T*5	DSEH010915	96T*15	0.28pg/ml(50 µl)
Human beta Amyloid 1-40 DTSet ELISA Kit	DSEH011105	96T*5	DSEH011115	96T*15	0.55pg/ml(50 µl)
Human beta Amyloid 1-42 DTSet ELISA Kit	DSEH011205	96T*5	DSEH011215	96T*15	15.61 pg/ml(50 µl)
Mammal VEGF-C DTSet ELISA Kit	DSEH011305	96T*5	DSEH011315	96T*15	4.23 pg/ml(50 µl)
Human CD27/TNFRSF7 DTSet ELISA Kit	DSEH011405	96T*5	DSEH011415	96T*15	0.00pg/ml(50 µl)
Human E-Cadherin DTSet ELISA Kit	DSEH011505	96T*5	DSEH011515	96T*15	14.75pg/ml(50 µl)
Human B7-H6/NCR3LG1 ELISA Kit	DSEH011605	96T*5	DSEH011615	96T*15	13.41 pg/ml(50 µl)
Human CCL25 ELISA Kit	DSEH011705	96T*5	DSEH011715	96T*15	4.37 pg/ml(50 µl)
Human IDO1 DTSet ELISA Kit	DSEH011805	96T*5	DSEH011815	96T*15	3.78 pg/ml(50 µl)
Human CD48 DTSet ELISA Kit	DSEH011905	96T*5	DSEH011915	96T*15	0.67pg/ml(50 µl)
Human OX40/TNFRSF4 DTSet ELISA Kit	DSEH012005	96T*5	DSEH012015	96T*15	1.08 pg/ml(50 µl)
Human Human B7-1/CD80 DTSet ELISA Kit	DSEH012105	96T*5	DSEH012115	96T*15	0.22pg/ml(50 µl)
Human TREM2 DTSet ELISA Kit	DSEH012205	96T*5	DSEH012215	96T*15	0.46 pg/ml(50 µl)
Human cIAP2 DTSet ELISA Kit	DSEH012305	96T*5	DSEH012315	96T*15	0.11 pg/ml(50 µl)
Human CD137 DTSet ELISA Kit	DSEH012405	96T*5	DSEH012415	96T*15	1.55 pg/ml(10 µl)
Human MMP-7 DTSet ELISA Kit	DSEH012505	96T*5	DSEH012515	96T*15	10.53 pg/ml(50 µl)
Human M-CSF DTSet ELISA Kit	DSEH012605	96T*5	DSEH012615	96T*15	0.81 pg/ml(50 µl)
Human MICB DTSet ELISA Kit	DSEH012705	96T*5	DSEH012715	96T*15	6.98 pg/ml(50 µl)
Human TNF R1 DTSet ELISA Kit	DSEH012805	96T*5	DSEH012815	96T*15	0.00pg/ml(50 µl)
Human Granzyme A/GZMA DTSet ELISA Kit	DSEH012905	96T*5	DSEH012915	96T*15	0.08 pg/ml(50 µl)
Human CCL23/MPIF-1 ELISA Kit	DSEH013005	96T*5	DSEH013015	96T*15	0.06 pg/ml(50 µl)
Mouse IL-1 beta DTSet ELISA Kit	DSEM000105	96T*5	DSEM000115	96T*15	0.54pg/ml(10 µl)

Product Name	Catalog#	Size	Catalog#	Size	Sensitivity
Mouse IL-4 DTSet ELISA Kit	DSEM000305	96T*5	DSEM000315	96T*15	0.81 pg/ml(10 µl)
Mouse IL-6 DTSet ELISA Kit	DSEM000405	96T*5	DSEM000415	96T*15	0.88pg/ml(10 µl)
Mouse IL-13 DTSet ELISA Kit	DSEM000705	96T*5	DSEM000715	96T*15	0.19 pg/ml(10 µl)
Mouse IL-17A DTSet ELISA Kit	DSEM000805	96T*5	DSEM000815	96T*15	0.35 pg/ml(10 µl)
Mouse TNF alpha DTSet ELISA Kit	DSEM001005	96T*5	DSEM001015	96T*15	1.10 pg/ml(10 µl)
Mouse IL-12/IL-23 p40 DTSet ELISA Kit	DSEM001205	96T*5	DSEM001215	96T*15	0.02pg/ml(10 µl)
Mouse VEGF DTSet ELISA Kit	DSEM001305	96T*5	DSEM001315	96T*15	3.76pg/ml(10 µl)
Mouse IL-9 DTSet ELISA Kit	DSEM002205	96T*5	DSEM002215	96T*15	0.02 pg/ml(10 µl)
Mouse IL-33 DTSet ELISA Kit	DSEM003305	96T*5	DSEM003315	96T*15	1.84pg/ml(10 µl)
Mouse PD-L1 DTSet ELISA Kit	DSEM003605	96T*5	DSEM003615	96T*15	0.03 pg/ml(10 µl)
Mouse PD1 DTSet ELISA Kit	DSEM003705	96T*5	DSEM003715	96T*15	7.33 pg/ml(10 µl)
Mouse LIF DTSet ELISA Kit	DSEM003805	96T*5	DSEM003815	96T*15	1.56 pg/ml(10 µl)
Mouse TREM2 DTSet ELISA Kit	DSEM003905	96T*5	DSEM003915	96T*15	0.10 pg/ml(10 µl)
Rat IL-1beta DTSet ELISA Kit	DSER000105	96T*5	DSER000115	96T*15	1.92 pg/ml(50 µl)
Rat IL-2 DTSet ELISA Kit	DSER000205	96T*5	DSER000215	96T*15	3.43 pg/ml(50 µl)
Rat IL-4 DTSet ELISA Kit	DSER000305	96T*5	DSER000315	96T*15	0.40 pg/ml(50 µl)
Rat IL-6 DTSet ELISA Kit	DSER000405	96T*5	DSER000415	96T*15	15.08 pg/ml(50 µl)
Rat IL-10 DTSet ELISA Kit	DSER000505	96T*5	DSER000515	96T*15	3.44 pg/ml(50 µl)
Rat IFN gamma DTSet ELISA Kit	DSER000605	96T*5	DSER000615	96T*15	6.05 pg/ml(50 µl)



扫码了解更多
“DTSet ELISA试剂盒”
产品信息



配对抗体

抗体对筛选“万里挑一”，品质保障，成就高质量的ELISA试剂盒

我们严格遵循夹心法ELISA筛选和验证捕获抗体和检测抗体。为保证实际应用场景中检测结果的特异性和一致性，我们使用生物样本进行筛选高品质和性能稳定的抗体对。同时，我们提供的抗体具有高稳定性和批次间一致性，确保实验结果的稳定可靠和可重复。



产品优势

生物样本验证

使用血清、血浆及细胞样本进行筛选和验证，保证实验结果的特异性和一致性

高灵敏度

对目标抗原的亲和力高，准确检测pg级别蛋白，适用于低表达天然蛋白的检测

高稳定性

在标准曲线范围内，严格控制CV值低于 10%，保证实验结果的准确可靠

高批次间一致性

采用HiMab重组兔单抗研发生产，确保稳定持续供应

自行搭配灵活兼容

提供PBS-only形式，实现捕获和检测抗体自行搭配，随时偶联


配对抗体

生物样本筛选，保证实验结果的特异性和一致性

人/小鼠/大鼠 NF-L配对抗体


Capture antibody: HA723330; Detect antibody: HA723331

NF-L是一种是分子量最小（68kDa）、含量最为丰富的神经丝亚单位，负责维持神经系统功能。NF-L通常作为多种神经系统疾病的生物标志物，在神经科学研究和临床诊断中具有重要价值。




使用Human/Mouse/Rat NF-L配对抗体检测不同样本中NF-L的含量。结果表明：小鼠脑组织提取物上清中NF-L含量的平均值是643,742 pg/mL；大鼠脑组织提取物上清中NF-L含量的平均值是127,718 pg/m；U-87 MG细胞提取物上清中NF-L含量低于检测限。

检测灵敏度高，有效捕获pg/mL级别蛋白





在相同实验操作条件下，使用Human IL-3配对抗体（Capture antibody: HA601153; Detect antibody: HA601155）和竞品S检测人源IL-3重组蛋白，并进行标准曲线统计分析。结果表明，较于竞品，我们的IL-3配对抗体最低检测限为0.0140 pg/ml，检测灵敏度明显优于竞品。

高批次间一致性，保证实验数据的可重复性




使用双抗夹心法对Human TREM2 配对抗体（Capture antibody: HA723149; Detect antibody: HA723151）的2次批次进行ELISA检测分析。结果表明，Human TREM2 配对抗体具有稳定的批次间一致性。

高稳定性，CV值低于10%，确保实验结果的稳定可靠



使用夹心ELISA方法对Human IL-1RA 配对抗体（**Capture antibody: HA722788**;
Detect antibody: HA722789）分别对冻融稳定性（左图）和37°C热稳定性（右图）
进行评价分析。
结果表明，Human IL-1RA 配对抗体CV值严格控制在低于10 %，具有高稳定性。



夹心法ELISA筛选配对抗体的实验流程

热销配对抗体清单精选

白介素及其受体			
Product name	Tested applications	Catalog#	Antibody type
Human IL-1 beta	ELISA(Cap)	HA721278	Rabbit mAb
Human IL-1 beta	ELISA(Det)	HA721279	Rabbit mAb
Mouse IL-1 beta	ELISA(Cap)	HA721976	Rabbit mAb
Mouse IL-1 beta	ELISA(Det)	HA721977	Rabbit mAb
Mouse IL-1 beta	ELISA(Cap)	HA722166	Rabbit mAb
Mouse IL-1 beta	ELISA(Det)	HA722167	Rabbit mAb
Rat IL-1 beta	ELISA(Cap)	HA721609	Rabbit mAb
Rat IL-1 beta	ELISA(Det)	HA721610	Rabbit mAb
Rat IL-1 beta	ELISA(Det)	HA721611	Rabbit mAb
Rat IL-1 beta	ELISA(Cap)	HA721836	Rabbit mAb
Rat IL-1 beta	ELISA(Det)	HA721837	Rabbit mAb
Human IL-1RA	ELISA(Cap)	HA722788	Rabbit mAb
Human IL-1RA	ELISA(Det)	HA722789	Rabbit mAb
Human IL-2 Receptor alpha	ELISA(Cap)	HA721280	Rabbit mAb
Human IL-2 Receptor alpha	ELISA(Det)	HA721281	Rabbit mAb
Human IL-2 Receptor alpha	ELISA(Cap)	HA722385	Rabbit mAb
Human IL-2 Receptor alpha	ELISA(Det)	HA722386	Rabbit mAb
Human IL-2	WB,ELISA(Det)	HA601006	Mouse mAb
Human IL-2	ELISA(Cap)	HA721680	Rabbit mAb
Mouse IL-2	ELISA(Cap)	HA722110	Rabbit mAb
Mouse IL-2	ELISA(Det)	HA722111	Rabbit mAb
Rat IL-2	ELISA(Cap)	HA721612	Rabbit mAb
Rat IL-2	ELISA(Det)	HA721613	Rabbit mAb
Rat IL-2	ELISA(Det)	HA721614	Rabbit mAb
Rat IL-2	ELISA(Cap)	HA722423	Rabbit mAb
Rat IL-2	ELISA(Det)	HA722424	Rabbit mAb
Human IL-3	ELISA(Cap)	HA601153	Mouse mAb
Human IL-3	ELISA(Cap)	HA601154	Mouse mAb
Human IL-3	ELISA(Det)	HA601155	Mouse mAb
Human IL3RA/CD123	ELISA(Cap)	HA723351	Rabbit mAb
Human IL3RA/CD123	ELISA(Det)	HA723352	Rabbit mAb

白介素及其受体			
Product name	Tested applications	Catalog#	Antibody type
Human IL-4	ELISA(Cap)	HA721230	Rabbit mAb
Human IL-4	ELISA(Det)	HA721231	Rabbit mAb
Human IL-4	ELISA(Det)	HA721295	Rabbit mAb
Human IL-4	ELISA(Det)	HA721368	Rabbit mAb
Mouse IL-4	ELISA(Cap)	HA722168	Rabbit mAb
Mouse IL-4	ELISA(Det)	HA722169	Rabbit mAb
Mouse IL-4	ELISA(Cap)	HA722262	Rabbit mAb
Mouse IL-4	ELISA(Det)	HA722263	Rabbit mAb
Rat IL-4	ELISA(Cap)	HA722104	Rabbit mAb
Rat IL-4	ELISA(Det)	HA722105	Rabbit mAb
Human IL-4R	ELISA(Cap)	HA723336	Rabbit mAb
Human IL-4R	ELISA(Det)	HA723337	Rabbit mAb
Human IL-5	ELISA(Cap)	HA722164	Rabbit mAb
Human IL-5	ELISA(Det)	HA722165	Rabbit mAb
Mouse IL-5	ELISA(Cap)	HA725034	Rabbit mAb
Mouse IL-5	ELISA(Det)	HA725035	Rabbit mAb
Human IL-6	ELISA(Cap)	HA721167	Rabbit mAb
Human IL-6	ELISA(Det)	HA721168	Rabbit mAb
Human IL-6	ELISA(Det)	HA721293	Rabbit mAb
Mouse IL-6	ELISA(Cap)	HA721411	Rabbit mAb
Mouse IL-6	ELISA(Det)	HA721412	Rabbit mAb
Mouse IL-6	ELISA(Det)	HA721413	Rabbit mAb
Mouse IL-6	ELISA(Cap)	HA721600	Rabbit mAb
Mouse IL-6	ELISA(Det)	HA721601	Rabbit mAb
Mouse IL-6	ELISA(Det)	HA721602	Rabbit mAb
Rat IL-6	ELISA(Cap)	HA721615	Rabbit mAb
Rat IL-6	ELISA(Det)	HA721616	Rabbit mAb
Rat IL-6	ELISA(Det)	HA721617	Rabbit mAb
Human IL-7	ELISA(Cap)	HA721937	Rabbit mAb
Human IL-7	ELISA(Det)	HA721938	Rabbit mAb
Human IL-8	ELISA(Cap)	HA722098	Rabbit mAb
Human IL-8	ELISA(Det)	HA722099	Rabbit mAb

白介素及其受体			
Product name	Tested applications	Catalog#	Antibody type
Mouse IL-9	ELISA(Cap)	HA721603	Rabbit mAb
Mouse IL-9	ELISA(Det)	HA721604	Rabbit mAb
Mouse IL-9	ELISA(Det)	HA721605	Rabbit mAb
Human IL-10	ELISA(Cap)	HA721236	Rabbit mAb
Human IL-10	ELISA(Det)	HA721237	Rabbit mAb
Mouse IL-10	ELISA(Cap)	HA721517	Rabbit mAb
Mouse IL-10	ELISA(Det)	HA721518	Rabbit mAb
Mouse IL-10	ELISA(Det)	HA721519	Rabbit mAb
Rat IL-10	ELISA(Cap)	HA721786	Rabbit mAb
Rat IL-10	ELISA(Det)	HA721787	Rabbit mAb
Mouse IL12 p40	ELISA(Cap)	HA722193	Rabbit mAb
Mouse IL12 p40	ELISA(Det)	HA722194	Rabbit mAb
Human IL-12B p40	ELISA(Cap)	HA721834	Rabbit mAb
Human IL-12B p40	ELISA(Det)	HA721835	Rabbit mAb
Human IL-13R alpha 2	ELISA(Cap)	HA722880	Rabbit mAb
Human IL-13R alpha 2	ELISA(Det)	HA722881	Rabbit mAb
Human IL-13	ELISA(Cap)	HA725013	Rabbit mAb
Human IL-13	ELISA(Det)	HA725014	Rabbit mAb
Mouse IL-13	ELISA(Cap)	HA721784	Rabbit mAb
Mouse IL-13	ELISA(Det)	HA721785	Rabbit mAb
Human IL-15	ELISA(Cap)	HA722102	Rabbit mAb
Human IL-15	ELISA(Det)	HA722103	Rabbit mAb
Human IL-17A	ELISA(Cap)	HA721232	Rabbit mAb
Human IL-17A	ELISA(Det)	HA721233	Rabbit mAb
Human IL-17A	ELISA(Det)	HA721296	Rabbit mAb
Human IL-17A	ELISA(Det)	HA721631	Rabbit mAb
Mouse IL-17A	ELISA(Cap)	HA721606	Rabbit mAb
Mouse IL-17A	ELISA(Det)	HA721607	Rabbit mAb
Mouse IL-17A	ELISA(Det)	HA721608	Rabbit mAb
Mouse IL-17A	ELISA(Cap)	HA722264	Rabbit mAb
Mouse IL-17A	ELISA(Det)	HA722265	Rabbit mAb
Human IL-18	ELISA(Cap)	HA721722	Rabbit mAb

白介素及其受体			
Product name	Tested applications	Catalog#	Antibody type
Human IL-18	ELISA(Det)	HA721723	Rabbit mAb
Human IL-18BP	ELISA(Cap)	HA723235	Rabbit mAb
Human IL-18BP	ELISA(Det)	HA723236	Rabbit mAb
Human IL-21	ELISA(Cap)	HA601195	Mouse mAb
Human IL-21	ELISA(Det)	HA601196	Mouse mAb
Mouse IL-21	ELISA(Cap)	HA722174	Rabbit mAb
Mouse IL-21	ELISA(Det)	HA722175	Rabbit mAb
Human IL-22	ELISA(Cap)	HA721939	Rabbit mAb
Human IL-22	ELISA(Det)	HA721940	Rabbit mAb
Human IL-23A	ELISA(Cap)	HA721933	Rabbit mAb
Human IL-12B p40	ELISA(Det)	HA721934	Rabbit mAb
Human IL-27	ELISA(Cap)	HA722046	Rabbit mAb
Human IL-27	ELISA(Det)	HA722047	Rabbit mAb
Human IL-27	ELISA(Cap)	HA722048	Rabbit mAb
Human IL-27	ELISA(Det)	HA722049	Rabbit mAb
Human IL-27 p28	ELISA(Cap)	HA722044	Rabbit mAb
Human IL-27 p28	ELISA(Det)	HA722045	Rabbit mAb
Human IL-31	ELISA(Cap)	HA721931	Rabbit mAb
Human IL-31	ELISA(Det)	HA721932	Rabbit mAb
Human IL-33	ELISA(Cap)	HA721935	Rabbit mAb
Human IL-33	ELISA(Det)	HA721936	Rabbit mAb
Mouse IL-33	ELISA(Cap)	HA722985	Rabbit mAb
Mouse IL-33	ELISA(Det)	HA722986	Rabbit mAb
肿瘤坏死因子			
Product name	Tested applications	Catalog#	Antibody type
Human TNF alpha	ELISA(Cap)	HA721265	Rabbit mAb
Human TNF alpha	ELISA(Det)	HA721266	Rabbit mAb
Human TNF alpha	ELISA(Det)	HA721294	Rabbit mAb
Human TNF alpha	ELISA(Det)	HA721559	Rabbit mAb
Mouse TNF alpha	ELISA(Cap)	HA722040	Rabbit mAb
Mouse TNF alpha	ELISA(Det)	HA722041	Rabbit mAb
Mouse/Rat TNF alpha	ELISA(Cap)	HA721974	Rabbit mAb

肿瘤坏死因子			
Product name	Tested applications	Catalog#	Antibody type
Mouse/Rat TNF alpha	ELISA(Det)	HA721975	Rabbit mAb
Rat TNF alpha	ELISA(Cap)	HA722170	Rabbit mAb
Rat TNF alpha	ELISA(Det)	HA722171	Rabbit mAb
Human TNF-beta	ELISA(Cap)	HA721242	Rabbit mAb
Human TNF-beta	ELISA(Det)	HA721243	Rabbit mAb
Human sTNF RI	ELISA(Cap)	HA723158	Rabbit mAb
Human sTNF RI	ELISA(Det)	HA723159	Rabbit mAb
Human TNF Receptor II	ELISA(Cap)	HA723402	Rabbit mAb
Human TNF Receptor II	ELISA(Det)	HA723403	Rabbit mAb
Human TNF Receptor II	ELISA(Cap)	HA723405	Rabbit mAb
干扰素			
Product name	Tested applications	Catalog#	Antibody type
Human IFN-alpha	ELISA(Cap)	HA721290	Rabbit mAb
Human IFN-alpha	ELISA(Det)	HA721291	Rabbit mAb
Human IFN-alpha	ELISA(Det)	HA721369	Rabbit mAb
Mouse IFN beta	ELISA(Cap)	HA722172	Rabbit mAb
Mouse IFN beta	ELISA(Det)	HA722173	Rabbit mAb
Human Interferon gamma	ELISA(Cap)	HA721574	Rabbit mAb
Human Interferon gamma	ELISA(Det)	HA721575	Rabbit mAb
Human Interferon gamma	ELISA(Det),ELISA	HA721576	Rabbit mAb
Human Interferon gamma	ELISA(Cap)	HA721683	Rabbit mAb
Human Interferon gamma	ELISA(Cap)	HA721684	Rabbit mAb
Human Interferon gamma	ELISA(Det)	HA721685	Rabbit mAb
Rat IFN gamma	ELISA(Cap)	HA721618	Rabbit mAb
Rat IFN gamma	ELISA(Det)	HA721619	Rabbit mAb
Rat IFN gamma	ELISA(Det)	HA721620	Rabbit mAb
Rat IFN gamma	ELISA(Cap)	HA722421	Rabbit mAb
Rat IFN gamma	ELISA(Det)	HA722422	Rabbit mAb
Human IFNGR1	ELISA(Cap)	HA722679	Rabbit mAb
Human IFNGR1	ELISA(Det)	HA722680	Rabbit mAb

集落刺激因子			
Product name	Tested applications	Catalog#	Antibody type
Human GM-CSF	ELISA(Cap)	HA721288	Rabbit mAb
Human GM-CSF	ELISA(Det)	HA721289	Rabbit mAb
Mouse GM-CSF	ELISA(Cap)	HA725019	Rabbit mAb
Mouse GM-CSF	ELISA(Det)	HA725020	Rabbit mAb
Human M-CSF	ELISA(Cap)	HA723323	Rabbit mAb
Human M-CSF	ELISA(Cap)	HA723324	Rabbit mAb
Human M-CSF	ELISA(Det)	HA723325	Rabbit mAb
生长因子			
Product name	Tested applications	Catalog#	Antibody type
Human Beta Nerve Growth Factor	ELISA(Cap)	HA722902	Rabbit mAb
Human Beta Nerve Growth Factor	ELISA(Cap)	HA722903	Rabbit mAb
Human Beta Nerve Growth Factor	ELISA(Det)	HA722904	Rabbit mAb
Human EGF	ELISA(Cap)	HA721782	Rabbit mAb
Human EGF	ELISA(Det)	HA721783	Rabbit mAb
Human FGF-2	ELISA(Cap)	HA721780	Rabbit mAb
Human FGF-2	ELISA(Det)	HA721781	Rabbit mAb
Human HGF	ELISA(Cap)	HA721778	Rabbit mAb
Human HGF	ELISA(Det)	HA721779	Rabbit mAb
Human HGF	ELISA(Cap)	HA722042	Rabbit mAb
Human HGF	ELISA(Det)	HA722043	Rabbit mAb
Human Insulin	ELISA(Cap)	HA722096	Rabbit mAb
Human Insulin	ELISA(Det)	HA722097	Rabbit mAb
Human IGFBP3	ELISA(Cap)	HA723432	Rabbit mAb
Human IGFBP3	ELISA(Det)	HA723433	Rabbit mAb
Mouse/Rat IGF1	ELISA(Cap)	HA725037	Rabbit mAb
Mouse/Rat IGF1	ELISA(Det)	HA725038	Rabbit mAb
Human Prolactin	ELISA(Cap)	HA721941	Rabbit mAb
Human Prolactin	ELISA(Det)	HA721942	Rabbit mAb
Human/Mouse VEGFC	ELISA(Cap)	HA723241	Rabbit mAb
Human/Mouse VEGFC	ELISA(Det)	HA723242	Rabbit mAb
Mouse VEGF	ELISA(Cap)	HA721850	Rabbit mAb
Mouse VEGF	ELISA(Det)	HA721851	Rabbit mAb

生长因子			
Product name	Tested applications	Catalog#	Antibody type
Human Leptin	ELISA(Cap)	HA723421	Rabbit mAb
Human Leptin	ELISA(Det)	HA723422	Rabbit mAb
Mouse/Rat Leptin	ELISA(Cap)	HA723046	Rabbit mAb
Mouse/Rat Leptin	ELISA(Det)	HA723047	Rabbit mAb
Human GH	ELISA(Cap)	HA725106	Rabbit mAb
Human GH	ELISA(Det)	HA725107	Rabbit mAb
Human SCF	ELISA(Cap)	HA722871	Rabbit mAb
Human SCF	ELISA(Det)	HA722872	Rabbit mAb
趋化因子及其受体			
Product name	Tested applications	Catalog#	Antibody type
Human CCL3 / MIP-1 alpha	ELISA(Cap)	HA723377	Rabbit mAb
Human CCL3 / MIP-1 alpha	ELISA(Det)	HA723378	Rabbit mAb
Human CCL3 / MIP-1 alpha	ELISA(Cap)	HA723380	Rabbit mAb
Mouse CCL2 / MCP-1	ELISA(Cap)	HA725040	Rabbit mAb
Mouse CCL2 / MCP-1	ELISA(Det)	HA725041	Rabbit mAb
Mouse CXCL2	ELISA(Cap)	HA725070	Rabbit mAb
Mouse CXCL2	ELISA(Det)	HA725071	Rabbit mAb
Mouse CCL3	ELISA(Cap)	HA722982	Rabbit mAb
Mouse CCL3	ELISA(Det)	HA722983	Rabbit mAb
Human CCL20/MIP-3 alpha	ELISA(Cap)	HA725010	Rabbit mAb
Human CCL20/MIP-3 alpha	ELISA(Det)	HA725011	Rabbit mAb
Human CCL23	ELISA(Cap)	HA723061	Rabbit mAb
Human CCL23	ELISA(Cap)	HA723062	Rabbit mAb
Human CCL23	ELISA(Det)	HA723063	Rabbit mAb
Human CCL25	ELISA(Cap)	HA723015	Rabbit mAb
Human CCL25	ELISA(Det)	HA723016	Rabbit mAb
Human CCL25	ELISA(Det)	HA723017	Rabbit mAb
Human CXCL1/GRO alpha	ELISA(Cap)	HA722493	Rabbit mAb
Human CXCL1/GRO alpha	ELISA(Det)	HA722494	Rabbit mAb
Human CXCL6/GCP-2	ELISA(Cap)	HA723314	Rabbit mAb
Human CXCL6/GCP-2	ELISA(Det)	HA723315	Rabbit mAb
Human CXCL9	ELISA(Cap)	HA725004	Rabbit mAb

趋化因子及其受体			
Product name	Tested applications	Catalog#	Antibody type
Human CXCL9	ELISA(Det)	HA725005	Rabbit mAb
Human CXCL10 / IP-10	ELISA(Cap)	HA725007	Rabbit mAb
Human CXCL10 / IP-10	ELISA(Det)	HA725008	Rabbit mAb
Mouse CXCL9	ELISA(Cap)	HA725022	Rabbit mAb
Mouse CXCL9	ELISA(Det)	HA725023	Rabbit mAb
Human CXCL13	ELISA(Cap)	HA722791	Rabbit mAb
Human CXCL13	ELISA(Det)	HA722792	Rabbit mAb
Human CXCL13	ELISA(Det)	HA722794	Rabbit mAb
神经科学			
Product name	Tested applications	Catalog#	Antibody type
Human Alpha-Synuclein	ELISA(Cap)	HA723467	Rabbit mAb
Human Alpha-Synuclein	ELISA(Det)	HA723468	Rabbit mAb
beta Amyloid 1-40	ELISA(Cap)	HA722077	Rabbit mAb
beta Amyloid 1-40	ELISA(Det)	HA722078	Rabbit mAb
beta Amyloid 1-42	ELISA(Cap)	HA722079	Rabbit mAb
beta Amyloid 1-42	ELISA(Det)	HA722080	Rabbit mAb
beta Amyloid 1-42	ELISA(Det)	HA722293	Rabbit mAb
Human CNTF	ELISA(Cap)	HA723381	Rabbit mAb
Human CNTF	ELISA(Det)	HA723382	Rabbit mAb
Human CNTF	ELISA(Cap)	HA723384	Rabbit mAb
Human CNTF	ELISA(Det)	HA723385	Rabbit mAb
GFAP	ELISA(Cap)	HA722234	Rabbit mAb
GFAP	ELISA(Det)	HA601342	Mouse mAb
Human CSPG4	ELISA(Cap)	HA723244	Rabbit mAb
Human CSPG4	ELISA(Cap)	HA723245	Rabbit mAb
Human CSPG4	ELISA(Det)	HA723246	Rabbit mAb
Human Mesothelin	ELISA(Cap)	HA723287	Rabbit mAb
Human Mesothelin	ELISA(Det)	HA723288	Rabbit mAb
Human NCAM1	ELISA(Cap)	HA722946	Rabbit mAb
Human NCAM1	ELISA(Det),ELISA(Cap)	HA722947	Rabbit mAb
Human NCAM1	ELISA(Det)	HA722949	Rabbit mAb
Human Nectin-2	ELISA(Cap)	HA723460	Rabbit mAb

神经科学			
Product name	Tested applications	Catalog#	Antibody type
Human Nectin-2	ELISA(Det)	HA723461	Rabbit mAb
Human Nectin-2	ELISA(Cap)	HA723463	Rabbit mAb
Human Nectin-4	ELISA(Cap)	HA723470	Rabbit mAb
Human Nectin-4	ELISA(Det)	HA723471	Rabbit mAb
Human Oxytocin-neurophysin 1	ELISA(Cap)	HA723452	Rabbit mAb
Human Oxytocin-neurophysin 1	ELISA(Det)	HA723453	Rabbit mAb
Human/Mouse/Rat NF-L	ELISA(Cap)	HA723330	Rabbit mAb
Human/Mouse/Rat NF-L	ELISA(Det)	HA723331	Rabbit mAb
Human/Mouse/Rat NF-L	ELISA(Cap)	HA723333	Rabbit mAb
Human/Mouse/Rat NF-L	ELISA(Det)	HA723334	Rabbit mAb
免疫学			
Product name	Tested applications	Catalog#	Antibody type
Human B7-H6	ELISA(Cap)	HA723190	Rabbit mAb
Human B7-H6	ELISA(Det)	HA723191	Rabbit mAb
Human CD27/TNFRSF7	ELISA(Cap)	HA722973	Rabbit mAb
Human CD27/TNFRSF7	ELISA(Cap)	HA722974	Rabbit mAb
Human CD27/TNFRSF7	ELISA(Det)	HA722975	Rabbit mAb
Human CD28	ELISA(Cap)	HA723457	Rabbit mAb
Human CD28	ELISA(Det)	HA723458	Rabbit mAb
Human CD30	ELISA(Cap)	HA722187	Rabbit mAb
Human CD30	ELISA(Det)	HA722188	Rabbit mAb
Human CD31	ELISA(Cap)	HA722906	Rabbit mAb
Human CD31	ELISA(Det)	HA722907	Rabbit mAb
Human CD33	ELISA(Cap)	HA722254	Rabbit mAb
Human CD33	ELISA(Det)	HA722255	Rabbit mAb
Human CD47	ELISA(Cap)	HA723272	Rabbit mAb
Human CD47	ELISA(Det)	HA723273	Rabbit mAb
Human CD48	ELISA(Cap)	HA722977	Rabbit mAb
Human CD48	ELISA(Cap),ELISA(Det)	HA722978	Rabbit mAb
Human CD48	ELISA(Det)	HA722979	Rabbit mAb
Human CD63	ELISA(Cap)	HA722988	Rabbit mAb
Human CD63	ELISA(Det)	HA722989	Rabbit mAb

免疫学			
Product name	Tested applications	Catalog#	Antibody type
Mouse CD68	ELISA(Cap)	HA723446	Rabbit mAb
Mouse CD68	ELISA(Det)	HA723447	Rabbit mAb
Mouse CD68	ELISA(Cap)	HA723449	Rabbit mAb
Mouse CD68	ELISA(Det)	HA723450	Rabbit mAb
Human CD73	ELISA(Cap)	HA722957	Rabbit mAb
Human CD73	ELISA(Det)	HA722958	Rabbit mAb
Human CD80	ELISA(Cap)	HA723123	Rabbit mAb
Human CD80	ELISA(Det)	HA723124	Rabbit mAb
Human CD80	ELISA(Det)	HA723125	Rabbit mAb
Human CD85d	ELISA(Cap)	HA722381	Rabbit mAb
Human CD85d	ELISA(Det)	HA722382	Rabbit mAb
Human CD137	ELISA(Cap)	HA723009	Rabbit mAb
Human CD137	ELISA(Det)	HA723010	Rabbit mAb
Human CD137	ELISA(Cap)	HA723012	Rabbit mAb
Human CD137	ELISA(Det)	HA723013	Rabbit mAb
Human CD200	ELISA(Cap)	HA722379	Rabbit mAb
Human CD200	ELISA(Det)	HA722380	Rabbit mAb
Human CD200R1	ELISA(Cap)	HA722256	Rabbit mAb
Human CD200R1	ELISA(Det)	HA722257	Rabbit mAb
Human CD276	ELISA(Cap)	HA722491	Rabbit mAb
Human CD276	ELISA(Det)	HA722492	Rabbit mAb
Human Complement C5	ELISA(Cap)	HA722682	Rabbit mAb
Human Complement C5	ELISA(Det)	HA722683	Rabbit mAb
Human CRP	ELISA(Cap)	HA722108	Rabbit mAb
Human CRP	ELISA(Det)	HA722109	Rabbit mAb
Mouse CRP	ELISA(Cap)	HA725073	Rabbit mAb
Mouse CRP	ELISA(Det)	HA725074	Rabbit mAb
Rat CRP	ELISA(Cap)	HA725049	Rabbit mAb
Rat CRP	ELISA(Det)	HA725050	Rabbit mAb
Human CTGF	ELISA(Cap)	HA723296	Rabbit mAb
Human CTGF	ELISA(Det)	HA723297	Rabbit mAb
Human CTLA-4	ELISA(Cap)	HA722954	Rabbit mAb

免疫学			
Product name	Tested applications	Catalog#	Antibody type
Human CTLA-4	ELISA(Det)	HA722955	Rabbit mAb
Human FOLR1	ELISA(Cap)	HA723354	Rabbit mAb
Human FOLR1	ELISA(Det)	HA723355	Rabbit mAb
Human GITR	ELISA(Cap)	HA722383	Rabbit mAb
Human GITR	ELISA(Det)	HA722384	Rabbit mAb
Human Granzyme A	ELISA(Cap)	HA723092	Rabbit mAb
Human Granzyme A	ELISA(Det)	HA723093	Rabbit mAb
Human Granzyme A	ELISA(Cap)	HA723095	Rabbit mAb
Human Granzyme A	ELISA(Det)	HA723096	Rabbit mAb
Human ICAM1	ELISA(Cap)	HA723184	Rabbit mAb
Human ICAM1	ELISA(Det)	HA723185	Rabbit mAb
Human ICOSLG/B7-H2/CD275	ELISA(Cap)	HA723317	Rabbit mAb
Human ICOSLG/B7-H2/CD275	ELISA(Det)	HA723318	Rabbit mAb
Human ICOSLG/B7-H2/CD275	ELISA(Cap)	HA723320	Rabbit mAb
Human ICOSLG/B7-H2/CD275	ELISA(Det)	HA723321	Rabbit mAb
Human LAG-3	ELISA(Cap)	HA722258	Rabbit mAb
Human LAG-3	ELISA(Det)	HA722259	Rabbit mAb
Human LIF	ELISA(Cap)	HA722495	Rabbit mAb
Human LIF	ELISA(Det)	HA722496	Rabbit mAb
Mouse LIF	ELISA(Cap)	HA721758	Rabbit mAb
Mouse LIF	ELISA(Det)	HA721759	Rabbit mAb
Human MBL2	ELISA(Cap)	HA722191	Rabbit mAb
Human MBL2	ELISA(Det)	HA722192	Rabbit mAb
Human MICA	ELISA(Cap)	HA723327	Rabbit mAb
Human MICA	ELISA(Det)	HA723328	Rabbit mAb
Human MICB	ELISA(Cap)	HA723238	Rabbit mAb
Human MICB	ELISA(Det)	HA723239	Rabbit mAb
Human NLRP3	ELISA(Cap)	HA723275	Rabbit mAb
Human NLRP3	ELISA(Det)	HA723276	Rabbit mAb
Human NLRP3	ELISA(Det)	HA723278	Rabbit mAb
Human OX40/TNFRSF4	ELISA(Cap)	HA723187	Rabbit mAb
Human OX40/TNFRSF4	ELISA(Det)	HA723188	Rabbit mAb

免疫学			
Product name	Tested applications	Catalog#	Antibody type
Human PD-L2	ELISA(Cap)	HA723145	Rabbit mAb
Human PD-L2	ELISA(Cap)	HA723146	Rabbit mAb
Human PD-L2	ELISA(Det)	HA723147	Rabbit mAb
Human SAA1	ELISA(Cap)	HA601333	Mouse mAb
Human SAA1	ELISA(Det)	HA601334	Mouse mAb
Human Siglec 7	ELISA(Cap)	HA723392	Rabbit mAb
Human Siglec 7	ELISA(Cap),ELISA(Det)	HA723393	Rabbit mAb
Human Siglec 7	ELISA(Det)	HA723395	Rabbit mAb
Human Siglec 9	ELISA(Cap)	HA723397	Rabbit mAb
Human Siglec 9	ELISA(Det)	HA723398	Rabbit mAb
Human Siglec 9	ELISA(Det)	HA723400	Rabbit mAb
Human TREM2	ELISA(Cap),ELISA(Det)	HA723149	Rabbit mAb
Human TREM2	ELISA(Cap)	HA723150	Rabbit mAb
Human TREM2	ELISA(Det)	HA723151	Rabbit mAb
Human TREM2	ELISA(Det)	HA723372	Rabbit mAb
Mouse Trem2	ELISA(Cap)	HA721760	Rabbit mAb
Mouse Trem2	ELISA(Det)	HA721761	Rabbit mAb
Mouse PD1	ELISA(Cap)	HA722100	Rabbit mAb
Mouse PD1	ELISA(Det)	HA722101	Rabbit mAb
Mouse PD-L1	ELISA(Cap)	HA722189	Rabbit mAb
Mouse PD-L1	ELISA(Det)	HA722190	Rabbit mAb
Mouse RANKL	ELISA(Cap)	HA722922	Rabbit mAb
Mouse RANKL	ELISA(Cap)	HA722923	Rabbit mAb
Mouse RANKL	ELISA(Det)	HA722924	Rabbit mAb
TNFRSF14/HVEM	ELISA(Cap)	HA722676	Rabbit mAb
TNFRSF14/HVEM	ELISA(Det)	HA722677	Rabbit mAb
Human ULBP1	ELISA(Cap)	HA723406	Rabbit mAb
Human ULBP1	ELISA(Det)	HA723407	Rabbit mAb
Human ULBP1	ELISA(Det)	HA723409	Rabbit mAb
心血管研究			
Product name	Tested applications	Catalog#	Antibody type
Human Adiponectin	ELISA(Cap)	HA722260	Rabbit mAb

心血管研究			
Product name	Tested applications	Catalog#	Antibody type
Human Adiponectin	ELISA(Det)	HA722261	Rabbit mAb
Human Angiogenin	ELISA(Cap)	HA725076	Rabbit mAb
Human Angiogenin	ELISA(Det)	HA725077	Rabbit mAb
Human CKMB	ELISA(Cap)	HA723348	Rabbit mAb
Human CKMB	ELISA(Det)	HA723349	Rabbit mAb
Human E-Cadherin	ELISA(Cap)	HA723055	Rabbit mAb
Human E-Cadherin	ELISA(Det)	HA723056	Rabbit mAb
Mouse Factor XII	ELISA(Cap)	HA723440	Rabbit mAb
Mouse Factor XII	ELISA(Det)	HA723441	Rabbit mAb
Mouse/Rat Factor XII	ELISA(Cap)	HA723443	Rabbit mAb
Mouse/Rat Factor XII	ELISA(Det)	HA723444	Rabbit mAb
Human IDO1	ELISA(Cap)	HA723154	Rabbit mAb
Human IDO1	ELISA(Det)	HA723155	Rabbit mAb
Human MMP-1	ELISA(Cap)	HA722874	Rabbit mAb
Human MMP-1	ELISA(Det)	HA722875	Rabbit mAb
Human MMP-3	ELISA(Cap)	HA722877	Rabbit mAb
Human MMP-3	ELISA(Det)	HA722878	Rabbit mAb
Human MMP-7	ELISA(Cap)	HA722926	Rabbit mAb
Human MMP-7	ELISA(Det)	HA722927	Rabbit mAb
Human MMP-8	ELISA(Cap)	HA721978	Rabbit mAb
Human MMP-8	ELISA(Det)	HA721979	Rabbit mAb
Human MMP-9	ELISA(Cap)	HA722106	Rabbit mAb
Human MMP-9	ELISA(Det)	HA722107	Rabbit mAb
Human MMP-9	ELISA(Cap)	HA722162	Rabbit mAb
Human MMP-9	ELISA(Det)	HA722163	Rabbit mAb
Human MMP13	ELISA(Cap)	HA723293	Rabbit mAb
Human MMP13	ELISA(Det)	HA723294	Rabbit mAb
Human Osteopontin	ELISA(Cap)	HA723004	Rabbit mAb
Human Osteopontin	ELISA(Det)	HA723005	Rabbit mAb
Human Osteopontin	ELISA(Det)	HA723006	Rabbit mAb
Human PDGF-BB	ELISA(Cap)	HA722497	Rabbit mAb
Human PDGF-BB	ELISA(Det)	HA722498	Rabbit mAb

心血管研究			
Product name	Tested applications	Catalog#	Antibody type
Human PAI1	ELISA(Cap),ELISA(Det)	HA723435	Rabbit mAb
Human PAI1	ELISA(Det)	HA723436	Rabbit mAb
Human PAI1	ELISA(Cap)	HA723438	Rabbit mAb
Human Pentraxin 3 / PTX3	ELISA(Cap)	HA723387	Rabbit mAb
Human Pentraxin 3 / PTX3	ELISA(Det)	HA723388	Rabbit mAb
Human Pentraxin 3 / PTX3	ELISA(Det)	HA723390	Rabbit mAb
Human TRAIL	ELISA(Cap)	HA722943	Rabbit mAb
Human TRAIL	ELISA(Det)	HA722944	Rabbit mAb
Human VE Cadherin	ELISA(Cap)	HA721838	Rabbit mAb
Human VE Cadherin	ELISA(Det)	HA721839	Rabbit mAb
Human VCAM-1	ELISA(Cap)	HA725064	Rabbit mAb
Human VCAM-1	ELISA(Det)	HA725065	Rabbit mAb
Human Von Willebrand Factor	ELISA(Cap)	HA723058	Rabbit mAb
Human Von Willebrand Factor	ELISA(Det)	HA723059	Rabbit mAb
细胞生物学			
Product name	Tested applications	Catalog#	Antibody type
clAP1	ELISA(Cap)	HA723280	Rabbit mAb
clAP1	ELISA(Cap)	HA723281	Rabbit mAb
clAP1	ELISA(Det)	HA723282	Rabbit mAb
Human clAP2	ELISA(Cap)	HA723098	Rabbit mAb
Human clAP2	ELISA(Det),ELISA(Cap)	HA723099	Rabbit mAb
Human clAP2	ELISA(Det)	HA723100	Rabbit mAb
Human LAP3	ELISA(Cap)	HA723284	Rabbit mAb
Human LAP3	ELISA(Det)	HA723285	Rabbit mAb
Mouse Pro-Collagen I alpha 1	ELISA(Cap)	HA722934	Rabbit mAb
Mouse Pro-Collagen I alpha 1	ELISA(Det)	HA722935	Rabbit mAb
NuMA	ELISA(Cap)	HA601223	Mouse mAb
NuMA	ELISA(Det)	HA601224	Mouse mAb
微生物学			
Product name	Tested applications	Catalog#	Antibody type
SARS-CoV-2 Nucleocapsid protein	WB,ELISA(Det)	HA600005	Mouse mAb
SARS-CoV-2 Nucleocapsid protein	WB,ELISA(Cap)	HA600006	Mouse mAb

微生物学			
Product name	Tested applications	Catalog#	Antibody type
SARS-CoV-2 Nucleocapsid protein	WB,ELISA(Cap)	HA600007	Mouse mAb
SARS-CoV-2 Nucleocapsid protein	WB,ELISA(Cap)	HA600008	Mouse mAb
HIV1 p24	ELISA(Det),WB	HA600009	Mouse mAb
HIV1 p24	WB,ELISA	HA600010	Mouse mAb
HIV1 p24	ELISA(Cap)	HA601104	Mouse mAb
HIV1 p24	ELISA(Det)	HA601105	Mouse mAb
HIV1 p24	ELISA(Cap)	HA721527	Mouse mAb
HIV1 p24	ELISA(Det)	HA721528	Mouse mAb
HIV1 p24	ELISA(Det)	HA721529	Rabbit mAb
HIV1 p24	ELISA(Cap)	HA601159	Rabbit mAb
HIV1 p24	ELISA(Det),WB	HA601160	Rabbit mAb
Human Poliovirus Receptor/PVR	ELISA(Cap)	HA723464	Mouse mAb
Human Poliovirus Receptor/PVR	ELISA(Det)	HA723465	Mouse mAb
			Rabbit mAb
			Rabbit mAb



扫码了解更多
“配对抗体研究”
产品信息