

流式抗体


(Flow Cytometry)

高 特异性 质量稳定可靠


检测数据案例

对人外周血T免疫细胞进行特异性分型与鉴定

华安生物




进口品牌



在人外周血中,对CD3⁺, CD4⁺, CD8⁺ T免疫细胞的分型结果与进口品牌一致 (CD45: HA600128F; CD3: HA600116F; CD4: HA600121F; CD8: HA600123F)。

稳定性实验案例

产品特性:质量稳定可靠



37°C FITC anti-human CD3 antibody (HA600116F)

CD3流式抗体 (HA600116F) 在37°C稳定性实验21天后,对CD3⁺ T细胞依然呈现特异性染色,且阳性信号较强。

热门流式抗体推荐

扫码查看全线
流式抗体产品清单



靶点名称	产品名称	货号	检测应用	抗体类型
Mouse IgG1 isotype control	FITC mouse IgG1 isotype control antibody	HA600113F	FC	Mouse mAb
	PE mouse IgG1 isotype control antibody	HA600114F	FC	Mouse mAb
	APC mouse IgG1 isotype control antibody	HA600115F	FC	Mouse mAb
Rabbit isotype control	FITC rabbit isotype control antibody	HA720191F	FC	Rabbit mAb
	PE rabbit isotype control antibody	HA720192F	FC	Rabbit mAb
	APC rabbit isotype control antibody	HA720193F	FC	Rabbit mAb
CD45	FITC anti-human CD45 antibody	HA600125F	FC	Mouse mAb
	PE anti-human CD45 antibody	HA600126F	FC	Mouse mAb
	APC anti-human CD45 antibody	HA600127F	FC	Mouse mAb
	percp-Cy5.5 anti-human CD45 antibody	HA600128F	FC	Mouse mAb
	iFluor 488 anti-human CD45 antibody	HA600109F	FC	Mouse mAb
CD3	FITC anti-human CD3 antibody	HA600116F	FC	Mouse mAb
	PE anti-human CD3 antibody	HA600117F	FC	Mouse mAb
	APC anti-human CD3 antibody	HA600118F	FC	Mouse mAb
CD4	FITC anti-human CD4 antibody	HA600119F	FC	Mouse mAb
	PE anti-human CD4 antibody	HA600120F	FC	Mouse mAb
	APC anti-human CD4 antibody	HA600121F	FC	Mouse mAb
	FITC anti-mouse CD4 antibody	HA720202F	FC	Ch mAb
	PE anti-mouse CD4 antibody	HA720203F	FC	Ch mAb
CD8a	APC anti-mouse CD4 antibody	HA720204F	FC	Ch mAb
	FITC anti-human CD8a antibody	HA600122F	FC	Mouse mAb
	PE anti-human CD8a antibody	HA600123F	FC	Mouse mAb
	APC anti-human CD8a antibody	HA600124F	FC	Mouse mAb
	FITC anti-mouse CD8 antibody	HA720205F	FC	Ch mAb
	PE anti-mouse CD8 antibody	HA720206F	FC	Ch mAb
CD11b	APC anti-mouse CD8 antibody	HA720207F	FC	Ch mAb
	FITC anti-human CD11b antibody	HA600130F	FC	Mouse mAb
CD14	FITC anti-human CD14 antibody	HA720194F	FC	Rabbit mAb
CD14	iFluor™ 647 anti-human CD14 antibody	HA720184F	FC	Rabbit mAb
CD15	FITC anti-human CD15 antibody	HA720195F	FC	Rabbit mAb
CD16	PE anti-human CD16 antibody	HA600129F	FC	Mouse mAb
CD33	FITC anti-human CD33 antibody	HA720197F	FC	Rabbit mAb
CD5	FITC anti-human CD5 antibody	HA720211F	FC	Rabbit mAb
CD58	FITC anti-human CD58 antibody	HA720198F	FC	Rabbit mAb
CD45RA	FITC anti-human CD45RA antibody	HA600131F	FC	Rabbit mAb
CD200	FITC anti-human CD200 antibody	HA720208F	FC	Rabbit mAb
BCL-2	FITC anti-human BCL-2 antibody	HA720210F	FC	Rabbit mAb
CD274 / PD-L1	FITC anti-mouse CD274 (PD-L1) antibody	HA720213F	FC	Rabbit mAb

杭州华安生物技术有限公司

Hangzhou HuaAn Biotechnology Co.,Ltd

地址：杭州市钱塘新区白杨街道6号大街452号高科技企业
孵化园区2号楼16层D区
邮编：310012 邮箱：sales@huabio.cn
电话：0571-88062880 传真：0571-88060400



聚焦抗体，专业服务！ www.huabio.cn
国内免费销售热线：4008-123-175

